

HLA-C08 2020 7K8

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42		
HLA-C allele^{1,2,3}																																												
*08:01:01:01 -08:01:06, 08:01:08-08:01:20, 08:01:22-08:01:29, 08:46, 08:50, 08:58, 08:60, 08:79, 08:81, 08:83-08:85, 08:91, 08:96-08:98, 08:105-08:106, 08:109, 08:119, 08:122, 08:124, 08:131, 08:133, 08:136, 08:145, 08:147-08:148, 08:153, 08:155, 08:157, 08:162, 08:164, 08:173N, 08:175-08:176, 08:186-08:187, 08:190, 08:192-08:194, 08:196-08:197, 08:199, 08:203	1	2										12									21												34									41		
*08:01:07	1	2										12									21							28					34								41			
*08:01:21, 08:163	1	2										12									21																				41			
*08:02:01:01 -08:02:05, *08:02:08 -08:02:24, 08:02:26-08:02:29, 08:45, 08:48, 08:67-08:68, 08:71, 08:76, 08:100, 08:108, 08:111, 08:114, 08:116, 08:120, 08:123, 08:125, 08:132, 08:134, 08:140, 08:146, 08:150-08:151, 08:156, 08:158-08:159, 08:161N, 08:166-08:167, 08:169-08:172, 08:179-08:181N, 08:184, 08:191, 08:195, 08:198, 08:200-08:201, 08:206			3									12	13																				34											
*08:02:06			3									12	13																	31			34											
*08:02:07			3									12	13															28					34											
*08:02:25, 08:183			3									12	13																															
*08:03:01 -08:03:04, 08:117	1	2		4								12									21												34											
*08:04:01:01 -08:04:03, 08:93, 08:139	1		3									12																					34											
*08:05			3		5							12	13				17																34											
*08:06	1	2		4		6						12									21												34											
*08:07			3				7					12	13																				34											
*08:08:01-08:08:02	1	2										12		14							21												34								41			
*08:09	1	2								10	11	12									21				25								34											
*08:10	1	2							9												21																				41			
*08:11 , 08:59	1	2									11	12									21				25								34								41			
*08:12:01:01-08:12:01:02			3								11	12	13																				34											
*08:13	1		3											15																			34											
*08:14	1	2		4		6						12									21				26								34											
*08:15:01-08:15:02 , 08:51					5							12																					33	34										
*08:16:01-08:16:02	1	2												15							21													34							41			
*08:17			3							10		12	13																					34										
*08:18			3										13	14																				34										
*08:19:01			3									12	13															27					31		34									
*08:19:02			3									12	13														27							34										
*08:20 , 08:99 ⁴	1	2										12									21											30			34						41			
*08:21	1	2			5							12					17				21													34								41		
*08:22:01:01-08:22:01:02 , 08:154, 08:178	1	2										12				16					21													34								41		
*08:23			3			6						12	13																					34										
*08:24, 08:88N ⁴	1	2							9			12									21													34								41		
*08:25			3										13	15			17																	34										
*08:26N	1	2										12						18			21													34								41		
*08:27			3									12							19		21													34										
*08:28			3									12	13				17													30				34										
*08:29			3									12							19															34										

five Control

[illegible]

HLA-C08 2020 7K8

*08:101	1	2	4							12					21			27				34						
*08:102	1	2			6					12					21							34					41	
*08:103			3		6					12	13							26				34						
*08:104	1		3			7				12												34						
*08:107			3							12	13											34				40		
*08:110			3							12	13								30			34						
*08:112			3			8				12	13											34						
*08:113	1		3						11	12							25					34						
*08:115			3			7				12	13				22							34						
*08:118		2								12	13				21							34					41	
*08:121N	1	2								12					21							34	36				41	
*08:126			3							12	13									32		34						
*08:127N	1	2								12					21							34			37		41	
*08:128	1	2							11	12					21			29				34					41	
*08:129N	1	2								12					21							34			38		41	
*08:130N	1	2								12					21							34				39	41	
*08:137	1	2			5					12			17		21				30			34					41	
*08:138	1		3			7				12					21	22						34						
*08:141Q	1	2								12					21							34			37	38	41	
*08:143	1	2								12					21			27				34					41	
*08:144:01	1	2		4						12					21			25		31		34						
*08:144:02	1	2		4						12					21			25				34						
*08:152			3						11	12	13							25				34						
*08:160, 16:139 ⁵										12					21							34						
*08:165	1		4							12					21			25				34						
*08:168	1		3							12			17						30			34						
*08:174	1	2		4						12												34						
*08:182	1		3			7				12					22							34						
*08:185			3							12	13		18									34						
*08:188	1					7				12												34						
*08:189	1	2								12					21							w					41	
*08:202			3							12	13		16									34						
*08:205	1	2								12					21							34						

[illegible]

[illegible]

[illegible]

[illegible]

*05:215		3					11	13						25	28							Negative Co
*05:224		3				9		13														
*06:67			5									21										
*06:129		3					11				19	21		25								
*06:262, 16:90						9						21										
*07:01:48, 07:02:35				7												30		34				
*07:01:74, 07:02:09				6									24	26				34				
*07:04:01:01-07:04:04, 07:04:06-07:04:10, 07:04:12-07:04:19, 07:11-07:12, 07:63, 07:139, 07:142, 07:181, 07:199:01-07:199:02, 07:272, 07:323-07:324, 07:329N, 07:336, 07:338, 07:354, 07:358, 07:361, 07:365, 07:394-07:395, 07:420, 07:426, 07:428, 07:459, 07:466-07:467, 07:480, 07:501, 07:523, 07:534-07:535, 07:552, 07:562-07:563, 07:569, 07:585-07:586, 07:600:01N-07:600:02N, 07:622, 07:626, 07:651, 07:655-07:656, 07:664, 07:672N, 07:674, 07:693, 07:698, 07:742, 07:751N, 07:780, 07:797N, 07:831												22					34					
*07:04:05, 07:328, 07:355, 07:357, 07:378, 07:403, 07:406												22										
*07:41		3						13								30		34				
*07:68, 07:260:01-07:260:02, 07:302							11					22		25				34				
*07:101			5				12					22						34				
*07:148, 07:161, 07:583			5				12											34				
*07:364							11					22						34				
*07:413, 07:422								14										34				
*07:447		3								19								34				
*07:487												22		25				34				
*12:16:01, 12:147, 16:80				6								21		24	26			34				
*12:109												21				30		34				
*12:127, 12:203			5				12					21						34				
*12:131, 16:42, 16:56												21	23					34				
*12:149								14				21						34				
*12:169	1											21						34	35			
*12:195:02, 12:217, 16:21												21		24				34				
*12:290, 12:304								16				21						34				
*14:87	1											21							35			
*15:23:01-15:23:02, 15:138, 15:158	1											21		24								
*15:29, 15:87	1							16				21										
*15:63	1		6									21		24	26							
*15:113	1		6									21			26							
*15:130	1	3										21					32					
*15:156N	1											21							36			
*15:200	1											21	23									
*16:33			6									21						34				
*16:100						9						21						34				
*17:05	1		5									21										
*17:47	1		6									21										
A*02:425, A*02:519, A*29:10:01-29:10:02, A*68:69														25								
A*32:72																	31					
B*14:32, B*15:337		3					11					22		25								

HLA-C08 2020 7K8

[illegible]