

INTERPRETATION TABLE**HLA-C*08 Add-on SSP typing**

	Well	
	1	2
Length of spec.	115	545
PCR product		
Length of int.	800	1070
pos. control ¹		
5'-primer(s) ²	527	972
	5' -TAC 3'	5' -CTA 3'
3'-primer(s) ³	601	1034
	5' -CTT 3'	5' -AgT 3'
Well No.	1	2
HLA-C allele		
*08:01:01-08:01:09, 08:03:01-08:03:02, 08:06, 08:08-08:11, 08:14, 08:16:01-08:16:02, 08:20-08:21, 08:24, 08:26N, 08:36N, 08:38, 08:40-08:42, 08:44, 08:46, 08:50, 08:58-08:61, 05:79	1	
*08:22, 08:56	1	2
*15:29		2
HLA-C allele ⁴		
Well No.	1	2