

| Well No. | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 | 25 | 26 | 27 | 28 | 29 |
|--|---|---|---|---|---|---|---|---|---|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| HLA-C allele ^{1,2,3} | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| *14:02:01, 14:02:04-14:02:16, 14:30, 14:36-14:37, 14:39-14:40, 14:42, 14:46, 14:50-14:52, 14:56, 14:59, 14:62, 14:65 | 1 | | | | 5 | | | | 9 | | | | | | | | | | | | | | | | | | | | |
| *14:02:02 | 1 | | | | 5 | | | | | | | | | | | | | | | | | | | | | | | | |
| *14:02:03 | | 2 | | 4 | 5 | | | | 9 | | | | | | | | | | | | | | | | | | | | |
| *14:03, 14:38, 14:61 | | 2 | | 4 | 5 | | | 8 | 9 | | | | | | | | | | | | | | | | | | | | |
| *14:04 | 1 | | 3 | | 5 | | 7 | | 9 | | | | | | | | | | | | | | | | | | | | |
| *14:05 | 1 | | | | | 6 | | | 9 | | | | | | | | | | | | | | | | | | | | |
| *14:06 | 1 | | | | | | | | 9 | 10 | | | | | | | | | | | | | | | | | | | |
| *14:07N | 1 | | | | 5 | | | | 9 | | 11 | | | | | | | | | | | | | | | | | | |
| *14:08 | | 2 | | 4 | | | | | 9 | | | 12 | | | | | | | | | | | | | | | | | |
| *14:09 | 1 | | | | | 6 | | | | | | | | | | | | | | | | | | | | | | | |
| *14:10 | | 2 | | 4 | 5 | | | 8 | 9 | | | | 13 | | | | | | | | | | | | | | | | |
| *14:11 | 1 | | | | 5 | | | | 9 | | | | | 14 | | | | | | | | | | | | | | | |
| *14:12 | 1 | | 3 | | 5 | | | | 9 | | | | | | 15 | | | | | | | | | | | | | | |
| *14:13 | 1 | | | | 5 | | | | 9 | | | | | | | 16 | | | | | | | | | | | | | |
| *14:14 | 1 | | | | 5 | | | | 9 | | | | 13 | | | | | | | | | | | | | | | | |
| *14:15 | 1 | | | | | | | | 9 | 10 | | | | 14 | | | | | | | | | | | | | | | |
| *14:16 | 1 | | | | | | | | 9 | | | | | | | | 17 | | | | | | | | | | | | |
| *14:17 | 1 | | | | 5 | | | | | | | | | | | | | 18 | | | | | | | | | | | |
| *14:18, 14:29 ⁴ | 1 | | | | 5 | | | | 9 | | | | | | | | | | 19 | | | | | | | | | | |
| *14:19 | 1 | | | | 5 | | | | 9 | | | | | | | | | | | | 20 | | | | | | | | |
| *14:20 | 1 | | | | 5 | | | | 9 | | | | | | | | | | | | | 21 | | | | | | | |
| *14:21N | 1 | | | | 5 | | | | 9 | 10 | | | | | | | | | | | | | | | | | | | |
| *14:22 | 1 | | | 4 | 5 | | | 8 | 9 | | | | | | | | 17 | | | | | | | | | | | | |
| *14:23 | 1 | | | | 5 | | | | 9 | | | 12 | | | | | | | | | | | | | | | | | |
| *14:24:01-14:24:02, 14:31 ⁴ | 1 | | | | 5 | | | | 9 | | | | | | | | | | | | | 22 | | | | | | | |
| *14:25 | 1 | | | | 5 | | | | 9 | | | | | | | | | | | | | | 23 | | | | | | |
| *14:26 | 1 | | | | 5 | | | | 9 | | | | | | | | | 18 | | | | | | | | | | | |
| *14:27 | 1 | | | | 5 | | | | 9 | | | | | | | | 17 | | | | | | | | | | | | |
| *14:28 | 1 | | | | | | | | | | | | | | | | | | | | 20 | | | | | | | | |
| *14:32 | 1 | | | | 5 | | | | 9 | | | | | | | | | | | | | 21 | | 23 | | | | | |
| *14:33 | 1 | | | | 5 | | | | 9 | | | | | | | | | | | | | | | | 24 | | | | |
| *14:34 | | | | | 5 | | | | 9 | | | | | | | | | | | | | | | | 24 | | | | |
| *14:35N | | 2 | | 4 | | | | 8 | 9 | | | | | | | | | | | | 20 | | | | | | | | |
| *14:41 | | 2 | | 4 | 5 | | | 8 | 9 | | | | | | | | | | | | | | | 24 | | | | | |
| *14:43 | 1 | | | | 5 | | | | 9 | | | | | | | | | | | | | | | | | | | 27 | |

Negative Control

