

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
HLA-C allele^{1,2,3}																																
*14:02:01:01 -14:02:01:03, 14:02:04-14:02:09, 14:02:11- 14:02:20, 14:30, 14:36-14:37, 14:39-14:40, 14:42, 14:50-14:52, 14:56, 14:59, 14:62, 14:65, 14:67- 14:69, 14:71-14:72, 14:75-14:76, 14:80, 14:83-14:84	1				5				9																							31
*14:02:02, 14:82	1				5																											31
*14:02:03		2		4	5				9																							31
*14:02:10, 14:85	1				5				9																							
*14:03 , 14:38, 14:61		2		4	5			8	9																							
*14:04	1		3		5		7		9																							31
*14:05	1					6			9																							31
*14:06	1								9	10																						31
*14:07N	1				5				9		11																					31
*14:08		2		4					9			12																				31
*14:09	1					6			9																							31
*14:10		2		4	5			8	9				13																			
*14:11	1				5				9				14																			31
*14:12	1		3		5				9					15																		31
*14:13	1				5				9						16																	31
*14:14 , 14:46:01-14:46:02 ⁴	1				5				9				13																			31
*14:15	1								9	10			14																			31
*14:16	1								9								17															31
*14:17	1				5													18														31
*14:18	1				5				9										19													31
*14:19	1				5				9											20												31
*14:20	1				5				9												21											31
*14:21N	1				5				9	10																						31
*14:22	1			4	5			8	9								17															
*14:23	1				5				9			12																				31
*14:24:01- 14:24:02 , 14:31 ⁴	1				5				9													22										31
*14:25	1				5				9														23									31
*14:26	1				5				9								18															31
*14:27	1				5				9								17															31
*14:28:01	1																			20												31

Negative Control

<p>*01:17, 01:23, 01:128, 06:11, 06:82, 07:01:01:01-07:01:40, 07:01:42-07:02:07, 07:02:10-07:02:28, 07:02:30-07:02:31, 07:02:33-07:02:40, 07:02:42-07:06, 07:10-07:17:02, 07:18-07:19, 07:22-07:33N, 07:35-07:38:02, 07:41-07:63, 07:65-07:68, 07:70-07:71, 07:74-07:75, 07:77-07:91, 07:93-07:95, 07:97-07:107, 07:109-07:126, 07:128-07:151, 07:153-07:155, 07:157-07:171, 07:173-07:176, 07:178-07:209, 07:211-07:222, 07:224-07:237, 07:239-07:241, 07:243-07:262, 07:264N-07:266, 07:268-07:294, 07:296-07:314:03, 07:316-07:321, 07:324-07:327, 07:329N-07:334, 07:336-07:356, 07:358-07:360, 07:362-07:389, 07:391-07:401, 07:403-07:404, 07:407-07:414, 07:417, 07:419-07:492, 07:494Q-07:558, 07:560-07:561, 07:564-07:577, 07:579-07:580, B*07:13, B*07:15, B*07:160, B*08:123, B*08:145, B*42:18</p>		7						Negative Control
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--	---	--	--	--	--	--	------------------

