

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
HLA-C allele ^{1,2,3}																																		
*15:02:01:01-15:02:20, 15:49-15:51, 15:53, 15:57, 15:60, 15:62, 15:67-15:68, 15:71, 15:73, 15:75, 15:78-15:80, 15:83, 15:86, 15:88-15:89, 15:91, 15:93-15:94, 15:98-15:101, 15:106	1	2									11				15																			
*15:03	1	2	3								11				15																			
*15:04, 16:70 ⁴		2		4											15																			
*15:05:01-15:05:11, 15:54, 15:59, 15:70, 15:90, 15:104	1	2			5										15																			
*15:06:01-15:06:03, 15:55	1	2				6									15																			
*15:07	1	2					7				11	12																				31		
*15:08, 15:74	1	2						8			11				15																			
*15:09, 15:61	1	2		4											15																			
*15:10:01-15:10:03	1									10	11				15																			
*15:11	1								9					14																				
*15:12	1	2									11		13		15																			
*15:13	1	2									11			14	15																			
*15:15	1	2									11				15	16																		
*15:16	1		3								11			14																				
*15:17	1										11			14																				
*15:18	1	2							9		11				15																			
*15:19	1	2		4				8							15																			
*15:21	1	2					7				11	12			15																			
*15:22	1	2			5	6									15																			
*15:23	1	2			5				9						15																			
*15:24	1	2											13		15																			
*15:25			w				7					12																				31		
*15:26	1	2				6					11				15	16																		
*15:27	1			4											16																			
*15:28	1	2	3								11	12			15																			
*15:29	1	2			5										15							23												
*15:30		2		4											15												27							
*15:31	1	2									11				15		18																	
*15:32Q, 15:41 ⁵	1	2									11				15							24												
*15:33, 15:84Q ⁵	1	2									11				15													28						
*15:34, 15:39 ⁵	1	2									11				15	17																		
*15:35, 15:47 ⁵	1	2									11				15						21													
*15:36	1	2			5										15	17																		
*15:37, 15:58, 15:72	1	2				6					11				15																			
*15:38	1	2									11				15							22												
*15:40, 15:102		2				6									15																			
*15:42	1	2									11				15			19																
*15:43	1										11	12		14																				
*15:44-15:45 ⁵	1	2									11				15					20														
*15:46	1	2			5										15			19																
*15:48	1	2									11											23												
*15:52	1	2									11				15																			
*15:56	1	2									11				15												25							
*15:63	1								9		11				15												26							
*15:64, 15:66, 15:76		2													15																			
*15:65		2		4		6									15																			
*15:69	1	2			5	6									15	16																		
*15:77		2		4											15	16																		
*15:81, 15:96Q ⁵	1	2									11				15																29			
*15:82	1	2									11				15																30			
*15:85	1	2									11																						31	
*15:87	1	2									11				15								23											
*15:92N	1	2			5						11				15						21													
*15:95N	1	2									11	12									20													
*15:97	1	2									11				15						20		22											
*15:103	1	2									11			14	15																			
*15:105Q		2													15									24						27				
*01:02:34, 08:14, 08:58, 08:80, 08:103, 08:105, 14:54, B*35:08:02, B*35:08:05, B*67:02							7																										31	
Well No.	1	2	3	4	5	6	7																											

