

INTERPRETATION TABLE**HLA-C*15 Add-on SSP typing**

	Well	
	1	2
Length of spec.	100	545
PCR product		
Length of int.	800	1070
pos. control ¹		
5'-primer(s) ²	420	972
	5' -TTC 3'	5' -CTA 3'
3'-primer(s) ³	477	1034
	5' -gCg 3'	5' -AgT 3'
Well No.	1	2
HLA-C allele ⁴		
*15:05:01-15:05:06, 15:22-15:23, 15:36, 15:46, 15:54	1	
*15:29	1	2
*08:22, 08:56		2
HLA-C allele		
Well No.	1	2