

INTERPRETATION TABLE		
HLA-C*15 Add-on SSP typing		
	Well	
	1	2
Length of spec.	100	545
PCR product		
Length of int.	100	545
pos. control¹		
5'-primer(s)²	420	972
	5' -TTC 3'	5' -CTA 3'
3'-primer(s)³	477	1034
	5' -gCg 3'	5' -AgT 3'
Well No.	1	2
HLA-C allele⁴		
*04:108, 15:05:01-15:05:08, 15:22-15:23, 15:36, 15:46, 15:54, 15:59, B*35:03:11, B*35:205	1	
*15:29	1	2
*08:22, 08:56		2
HLA-C allele		
Well No.	1	2