

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-C allele^{1,2,3}																								
*16:01:01:01-16:01:01:05, 16:01:03-16:01:15, 16:01:17-16:01:24, 16:34, 16:36, 16:41, 16:44, 16:51, 16:54, 16:59, 16:62, 16:65, 16:72-16:73, 16:75-16:76, 16:79, 16:81, 16:83, 16:86-16:87, 16:92-16:93, 16:95-16:96, 16:98, 16:100, 16:105, 16:111-16:112, 16:114, 16:116, 16:118-16:119, 16:122	1			4	5																			
*16:01:02, 16:01:16, 16:02:15, 16:70, 16:85, 16:94, 16:106	1				5																			
*16:02:01-16:02:14, 16:02:16, 16:47, 16:57, 16:60, 16:63, 16:69, 16:74, 16:84, 16:99, 16:101-16:102, 16:107-16:108, 16:120-16:121 ⁶	1	2			5																			
*16:04:01:01-16:04:01:02, 16:04:04-16:04:05, 16:33, 16:66, 16:78, 16:82, 16:109, 12:176 ⁴	1		3	4																				
*16:04:03	1		3	4											15									
*16:06				4	5	6										16								
*16:07:01				4	5		7									16								
*16:07:02, 16:117				4	5											16								
*16:08	1			4	5			8																
*16:09	1	2			5		7																	
*16:10	1			4		6																		
*16:11	1			4	5	6																		
*16:12	1	2			5				9															
*16:13	1			4	5					10														
*16:14	1			4	5								13											
*16:15:01-16:15:02, 16:20 ⁵	1			4	5						11													
*16:16Q-16:17 ⁵	1			4	5							12												
*16:18	1			4	5									14										
*16:19	1	2			5					10														
*16:21, 16:80	1			4	5										15									
*16:22				4	5							12												
*16:23	1			4	5									14						20				
*16:24, 16:58 ⁵	1			4	5																	22		
*16:25	1	2			5						11													
*16:26	1			4	5												17	18						
*16:27, 16:32 ⁵	1			4	5															20				
*16:28, 16:31, 16:50 ⁵	1			4	5															19				
*16:29	1		3	4																19				
*16:30N, 16:56 ⁵	1			4	5																		23	
*16:35				4					9															
*16:37, 16:71	1			w	5																			
*16:38	1			4	5			8										18						
*16:39:01-16:39:02	1			4	5	6								15										
*16:40, 16:110				4																	21			
*16:42	1		3	4																			23	
*16:43				4	5																			
*16:45	1			4			7																	
*16:46	1	2			5												17	18						
*16:48		2							9															
*16:49	1			4	5																	21		
*16:52	1			4	5				9															
*16:53				4				8														21		
*16:55	1		3	4													17	18						
*16:61	1		3	4																				
*16:64	1			4	5						11						17	18						
*16:67	1			4								12							19					
*16:68				4				8																
*16:77N	1	2			5								13											
*16:88		2						8																
*16:89N	1	2			5			8																
*16:90	1	2			5											16								
*16:91	1	2	3																					
*16:97	1			4	5												17							
*16:103	1	2																						
*16:104	1	2			5									14							20			
*16:113				4	5																	21		
*16:115	1	2			5	6																		
*16:123N	1	2			5			8					13											
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Negative Control

*17:21		2																19								
A*24:96, A*24:146							7																			
A*24:106						6																				
A*24:174					5																					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24		