

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-C allele^{1,2,3}																								
*16:01:01:01-16:01:01:37, 16:01:03-16:01:15, 16:01:17-16:01:25, 16:01:27-16:01:42, 16:34, 16:36, 16:41, 16:44, 16:49:02, 16:51, 16:54, 16:59, 16:62, 16:65, 16:72-16:73, 16:75-16:76, 16:79, 16:81, 16:83, 16:86-16:87, 16:92-16:93, 16:95-16:96, 16:98, 16:100, 16:105, 16:111-16:112, 16:114, 16:116, 16:118-16:119, 16:122, 16:125, 16:127-16:131, 16:134-16:135, 16:137-16:139, 16:141-16:142, 16:146-16:148, 16:151-16:152, 16:154, 16:157-16:162, 16:164-16:165, 16:168, 16:171-16:175, 16:177-16:178, 16:183, 16:185-16:187, 16:190, 16:192-16:193, 16:196-16:197	1			4	5																			
*16:01:02, 16:01:16, 16:01:26, 16:02:18, 16:85, 16:94, 16:106, 16:121, 16:126, 16:191	1				5																			
*16:02:01:01-16:02:17, 16:02:19-16:02:20, 16:47, 16:57, 16:60, 16:63, 16:69-16:70, 16:84, 16:99, 16:101-16:102, 16:107-16:108, 16:120, 16:132N-16:133, 16:136, 16:143-16:145, 16:155, 16:163, 16:166, 16:176, 16:179, 16:181, 16:194-16:195N	1	2			5																			
*16:04:01:01-16:04:01:03, 16:04:04-16:04:05, 16:66, 16:78, 16:82, 16:109, 16:124, 16:149-16:150, 16:182, 12:176, 12:333	1		3	4																				
*16:04:03	1		3	4											15									
*16:06				4	5	6										16								
*16:07:01				4	5		7									16								
*16:07:02, 16:117				4	5											16								
*16:08	1			4	5			8																
*16:09	1	2			5		7																	
*16:10	1			4		6																		
*16:11	1			4	5	6																		
*16:12	1	2			5				9															
*16:13	1			4	5					10														
*16:14	1			4	5								13											
*16:15:01-16:15:02, 16:20 ⁴	1			4	5						11													
*16:16Q-16:17 ⁴	1			4	5							12												
*16:18	1			4	5									14										
*16:19	1	2			5					10														
*16:21, 16:80	1			4	5										15									
*16:22				4	5							12												
*16:23	1			4	5									14						20				
*16:24, 16:58 ⁴	1			4	5																	22		
*16:25, 16:189	1	2			5						11													
*16:26	1			4	5												17	18						
*16:27, 16:32 ⁴	1			4	5															20				
*16:28, 16:31, 16:50 ⁴	1			4	5																19			
*16:29	1		3	4																	19			
*16:30N, 16:56 ⁴	1			4	5																		23	
*16:33	1		3	4					9															
*16:35				4					9															
*16:37, 16:71, 16:170	1			w	5																			
*16:38	1			4	5			8													18			
*16:39:01-16:39:02	1			4	5	6									15									
*16:40, 16:110				4																		21		
*16:42	1		3	4																			23	
*16:43				4	5																			
*16:45	1			4			7																	
*16:46	1	2			5												17	18						
*16:48		2							9															
*16:49:01	1			4	5																	21		
*16:52	1			4	5				9															
*16:53				4				8														21		
*16:55	1		3	4													17	18						
*16:61	1		3	4						10														
*16:64	1			4	5						11						17	18						
*16:67	1			4								12									19			
*16:68				4				8																
*16:74	1	2			5									14										
*16:77N	1	2			5								13											
*16:88		2						8																
*16:89N	1	2			5			8																
*16:90	1	2			5										16									
*16:91	1	2	3																					
*16:97	1			4	5												17							
*16:103, 16:140	1	2																						
*16:104	1	2			5									14							20			
*16:113				4	5																	21		
*16:115	1	2			5	6																		

Negative Control

*16:123N	1	2			5			8					13														
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24			
*16:153, 16:184		2			5																						
*16:156	1	2			5														19								
*16:167	1	2			5			8										18									
*16:169	1			4							11																
*16:180	1			4	5											16											
*16:188	1			4																							
*01:02:34, 03:03:64, 03:04:25, 08:01:11, 08:02:36, 08:12:01:01-08:12:01:02, 08:34, 08:44, 08:61, 08:82, 08:126, 08:171, 08:180N, 08:209, 08:219, 12:02:12, 12:02:43, 12:03:47, 12:19-12:20, 12:24, 12:44, 12:134, 12:153, 12:175, 12:183, 12:200, 12:251, 12:258, 12:270N, 12:309, 12:335, 12:344, 14:02:03, 14:03:01:01-14:03:01:03, 14:03:03-14:03:05, 14:08, 14:10, 14:22, 14:35N, 14:38, 14:41, 14:53-14:54, 14:61, 14:70, 14:79, 14:86, 14:112-14:113, 14:118, 14:126, 14:141N, 15:07:01:01-15:07:01:02, 15:144, B*07:02:77, B*08:01:36, B*15:01:50, B*15:436, B*18:03:02, B*35:08:02, B*35:08:05, B*55:01:25, B*67:02:01:01-67:02:01:02				4																							
*01:04, 01:09:01-01:09:02, 06:217, 12:03:04, 12:03:09, 12:03:34:01-12:03:34:02, 12:195:01, 12:195:03, 12:358, 14:16, B*15:510, B*35:475, B*46:84				3																							
*01:05, 01:36:01-01:36:02, 01:55, 01:79:01-01:79:02, 01:120, 02:02:01-02:02:03, 02:02:06-02:02:08, 02:02:10-02:02:30, 02:02:32-02:02:48, 02:02:50-02:04, 02:06:01-02:11, 02:13, 02:15-02:16:02:02, 02:18-02:20, 02:22-02:26:03, 02:28-02:33, 02:35-02:36:02, 02:38:01N-02:40:02, 02:42-02:48, 02:50-02:56, 02:58-02:61:02, 02:63-02:73, 02:76-02:78, 02:80, 02:82-02:86, 02:88-02:106, 02:108-02:114, 02:116-02:120, 02:122, 02:124-02:125, 02:127-02:130, 02:132-02:158, 02:160, 02:162, 02:165N-02:209, 02:211Q-02:212, 03:05:01:01-03:05:01:03, 03:13:01:01-03:13:02, 03:25, 03:27, 03:35:01-03:35:02:02, 03:135, 03:167, 03:178, 03:198, 03:267, 03:292, 03:296:01-03:296:02, 03:335, 03:386, 03:407, 03:474, 03:482, 03:494, 03:555, 05:129, 07:01:23, 07:02:11, 07:02:142, 07:03, 07:07, 07:16, 07:20:01, 07:36, 07:51, 07:56:01-07:56:02, 07:59-07:60, 07:64, 07:73:02, 07:79, 07:88, 07:96:01, 07:109:01-07:109:02, 07:127:01-07:127:02, 07:172:02, 07:227N, 07:263, 07:314:03, 07:327, 07:390, 07:402, 07:438, 07:441:02, 07:561, 07:578, 07:580, 07:598, 07:668, 07:756, 07:914, 07:919, 07:922, 07:970N, 07:993, 08:01:30, 08:02:11:01-08:02:11:02, 12:02:11, 12:02:29, 12:16:02, 12:27, 12:123, 14:09, 14:28:02, 15:22, 15:65, 15:72, 15:190, 15:246, 17:16:02															16												
*01:10, 01:108, 14:25, 14:119, 14:123, 17:21, B*07:239, B*14:46, B*14:52, B*40:243																			19								
*01:14, 01:59, 01:118, 03:07:01:01-03:07:02, 03:15, 03:45, 03:130, 03:140:01:01-03:140:01:02, 03:163, 03:243, 03:268, 03:297, 03:461, 05:01:21, 05:07N, 05:10, 05:16, 05:31, 05:85, 05:96, 05:104, 05:107, 05:148, 05:152, 05:175N, 05:214, 05:241, 06:204, 12:09, 14:04, 14:12, 14:49, 14:64, 14:77, 14:108, 14:140, 15:205, 17:02		2																									
*01:21, 02:27:01-02:27:02, 02:87, 02:131, 04:11, 04:55, 04:172, 07:723, 08:01:01:01-08:01:07, 08:01:09-08:01:10, 08:01:12-08:01:29, 08:01:31-08:02:10, 08:02:12-08:02:35, 08:03:01:01-08:09, 08:11, 08:13-08:27, 08:29-08:33:05, 08:35-08:36N, 08:38-08:43, 08:45-08:60, 08:62:01-08:63, 08:65-08:68, 08:70Q-08:81, 08:83-08:94, 08:97-08:125, 08:127N-08:134, 08:136, 08:138-08:167, 08:169-08:170, 08:172-08:179, 08:181N-08:208N, 08:210-08:218, 08:220, 08:222-08:243, 12:02:01-12:02:10, 12:02:13-12:02:22, 12:02:24-12:02:27, 12:02:30-12:02:42, 12:02:44-12:02:46, 12:08, 12:10:01-12:10:03, 12:14:01-12:14:02, 12:17-12:18:02, 12:22, 12:30, 12:36, 12:49, 12:56, 12:64, 12:67-12:69, 12:73-12:74, 12:80N, 12:83-12:86, 12:96, 12:103-12:106, 12:112, 12:114, 12:117-12:118, 12:124, 12:126-12:128, 12:130, 12:136, 12:137, 12:142, 12:145, 12:148N, 12:151, 12:155Q, 12:161-12:162, 12:164, 12:168-12:169, 12:177, 12:179, 12:181, 12:193, 12:196, 12:198, 12:204, 12:207-12:208, 12:212, 12:214, 12:219N, 12:221-12:222, 12:224, 12:226, 12:228, 12:231, 12:233-12:234, 12:236N, 12:239-12:241, 12:243, 12:247, 12:250, 12:252, 12:255, 12:261, 12:263, 12:268, 12:275, 12:279-12:281, 12:285, 12:294, 12:296, 12:298, 12:303-12:304, 12:307-12:308, 12:315, 12:327N, 12:329N, 12:336, 12:338, 12:351Q, 12:359, 12:363				4											16												

Negative Control

*04:80	2									17	18						
*04:101	2					10				16	17						
*04:214, 04:449			w							16	17						
*04:273	2									16	17					23	
*04:293	2									16	17					22	
*04:360	2	3									17						
*05:43	2									15	16						
*05:56	2										16						23
*05:81, 05:117	2					10					16						
*05:112	2				7						16				20		
*05:169N	2				8						16						
*06:02:01:01-06:02:01:92, 06:02:03-06:02:09, 06:02:11, 06:02:13-06:02:15, 06:02:17-06:02:42, 06:02:44-06:02:64, 06:02:66-06:02:71, 06:02:73, 06:02:75-06:03:02, 06:07, 06:09:01-06:10, 06:12-06:13, 06:15-06:16N, 06:18-06:21, 06:23-06:25, 06:27-06:29, 06:33-06:34:02, 06:36-06:39, 06:41-06:51, 06:54-06:71, 06:73-06:76:01, 06:77-06:78, 06:80, 06:82-06:86, 06:88-06:89, 06:91-06:100, 06:102:01-06:117, 06:119-06:121, 06:126, 06:128N-06:131, 06:134N-06:135, 06:137-06:142, 06:145-06:146, 06:148-06:152N, 06:154-06:160, 06:162-06:168, 06:170, 06:172-06:196, 06:199-06:202, 06:205-06:209:02, 06:211:01:01N-06:216, 06:218-06:235, 06:237-06:243, 06:245-06:247, 06:249-06:251, 06:253-06:263N, 06:265-06:270, 06:272-06:279, 06:281N-06:298, 06:300-06:302, 06:305-06:316N, 06:318-06:326Q, 06:328, 06:331-06:351	2	3										18					
*06:02:10, 06:02:74, 06:26, 06:30, 06:53:01-06:53:02, 06:125, 06:132:01-06:132:02, 06:169, 06:198, 06:210, 06:303, 12:04:01-12:05:02, 12:33, 12:54, 12:60, 12:353	2	3															
*06:02:12, 06:02:65, 06:11, 06:17, 06:52, 06:122, 06:124, 06:147, 06:236, 06:244, 06:248, 06:252, 06:280, 06:329:01:01-06:329:01:02			3									18					
*06:02:16, 06:02:43, 06:04:01-06:04:02:02, 06:06:01:01-06:06:01:02, 06:14, 06:35, 06:40, 06:72, 06:79N, 06:81, 06:123, 06:143, 06:153:01-06:153:02, 06:197, 06:203, 06:271, 06:317	2											18					
*06:02:72	2	3								15	17	18					
*06:05	2									15	17	18					
*06:08, 06:22, 06:90	2	3										18	19				
*06:31	1	2	3									18					
*06:32, 06:171:01:01N-06:171:01:02N	2	3							13			18					
*06:76:02	2	3									17	18					
*06:87, 06:330	2	3						10				18					
*06:118	2					9						18					
*06:133, 06:161, 06:299	2	3				9						18					
*06:327	2	3									16	18					

Negative Control

