

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-C allele^{1,2,3}																								
*16:01:01:01-16:01:01:38, 16:01:03-16:01:15, 16:01:17-16:01:25, 16:01:27-16:01:43, 16:34, 16:36, 16:41, 16:44, 16:49:02, 16:51, 16:54, 16:59, 16:62, 16:65, 16:72-16:73, 16:75-16:76, 16:79, 16:81, 16:83, 16:86-16:87, 16:92-16:93, 16:95-16:96, 16:98, 16:100, 16:105, 16:111-16:112, 16:114, 16:116, 16:118-16:119, 16:122, 16:125, 16:127-16:131, 16:134-16:135, 16:137-16:139, 16:141-16:142, 16:146-16:148, 16:151-16:152, 16:154, 16:157-16:162, 16:164-16:165, 16:168, 16:171-16:175, 16:177-16:178, 16:183, 16:185-16:187, 16:190, 16:192-16:193, 16:196-16:198, 16:200	1			4	5																			
*16:01:02, 16:01:16, 16:01:26, 16:02:18, 16:85, 16:94, 16:106, 16:121, 16:126, 16:191	1				5																			
*16:02:01:01-16:02:17, 16:02:19-16:02:20, 16:47, 16:57, 16:60, 16:63, 16:69-16:70, 16:84, 16:99, 16:101-16:102, 16:107-16:108, 16:120, 16:132N-16:133, 16:136, 16:143-16:145, 16:155, 16:163, 16:166, 16:176, 16:179, 16:181, 16:194-16:195N, 16:202	1	2			5																			
*16:04:01:01-16:04:01:03, 16:04:04-16:04:05, 16:66, 16:78, 16:82, 16:109, 16:124, 16:149-16:150, 16:182, 16:203, 12:176, 12:333 ⁵	1		3	4																				
*16:04:03	1		3	4											15									
*16:06				4	5	6										16								
*16:07:01				4	5		7									16								
*16:07:02, 16:117				4	5											16								
*16:08	1			4	5			8																
*16:09	1	2			5		7																	
*16:10	1			4		6																		
*16:11	1			4	5	6																		
*16:12	1	2			5				9															
*16:13	1			4	5					10														
*16:14	1			4	5								13											
*16:15:02, 16:20, 16:201 ⁴	1			4	5						11													
*16:16Q-16:17 ⁴	1			4	5							12												
*16:18	1			4	5									14										
*16:19	1	2			5					10														
*16:21, 16:80	1			4	5										15									
*16:22				4	5							12												
*16:23	1			4	5									14						20				
*16:24, 16:58 ⁴	1			4	5																	22		
*16:25, 16:189	1	2			5						11													
*16:26	1			4	5												17	18						
*16:27, 16:32 ⁴	1			4	5															20				
*16:28, 16:31, 16:50 ⁴	1			4	5														19					
*16:29	1		3	4															19					
*16:30N, 16:56 ⁴	1			4	5																		23	
*16:33	1		3	4					9															
*16:35				4					9															
*16:37, 16:71, 16:170	1			w	5																			
*16:38	1			4	5			8										18						
*16:39:01-16:39:02	1			4	5	6									15									
*16:40, 16:110				4																	21			
*16:42	1		3	4																			23	
*16:43, 16:204				4	5																			
*16:45	1			4			7																	
*16:46	1	2			5												17	18						
*16:48		2							9															
*16:49:01	1			4	5																21			
*16:52	1			4	5				9															
*16:53				4				8													21			
*16:55	1		3	4													17	18						
*16:61	1		3	4						10														
*16:64	1			4	5						11						17	18						
*16:67	1			4								12							19					
*16:68				4				8																
*16:74	1	2			5										14									
*16:77N	1	2			5								13											
*16:88		2						8																
*16:89N	1	2			5			8																
*16:90	1	2			5											16								
*16:91	1	2	3																					
*16:97	1			4	5												17							
*16:103, 16:140	1	2																						
*16:104	1	2			5									14						20				
*16:113				4	5																21			

Negative Control

[illegible]

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*04:69		2							9							16	17							
*04:80		2															17	18						
*04:101, 04:487		2								10						16	17							
*04:214, 04:449				w												16	17							
*04:273		2														16	17						23	
*04:293, 04:488		2														16	17					22		
*04:360, 04:477		2	3														17							
*05:43		2													15	16								
*05:56		2														16							23	
*05:81, 05:117		2								10						16								
*05:112		2					7									16			20					
*05:169N		2						8								16								
*06:02:01:01-06:02:01:93, 06:02:03-06:02:09, 06:02:11, 06:02:13-06:02:15, 06:02:17-06:02:42, 06:02:44-06:02:64, 06:02:66-06:02:71, 06:02:73, 06:02:75-06:03:02, 06:07, 06:09:01-06:10, 06:12-06:13, 06:15-06:16N, 06:18-06:21, 06:23-06:25, 06:27-06:29, 06:33-06:34:02, 06:36-06:39, 06:41-06:51, 06:54-06:71, 06:73-06:76:01, 06:77-06:78, 06:80, 06:82-06:86, 06:88-06:89, 06:91-06:100, 06:102:01-06:117, 06:119-06:121, 06:126, 06:128N-06:131, 06:134N-06:135, 06:137-06:142, 06:145-06:146, 06:148-06:152N, 06:154-06:160, 06:162-06:168, 06:170, 06:172-06:196, 06:199-06:202, 06:205-06:209:02, 06:211:01:01N-06:216, 06:218-06:235, 06:237-06:243, 06:245-06:247, 06:249-06:251, 06:253-06:263N, 06:265-06:270, 06:272-06:279, 06:281N-06:298, 06:300-06:302, 06:305-06:316N, 06:318-06:326Q, 06:328, 06:331-06:357N, 06:360-06:362N			2	3													18							
*06:02:10, 06:02:74, 06:26, 06:30, 06:53:01-06:53:02, 06:125, 06:132:01-06:132:02, 06:169, 06:198, 06:210, 06:303, 12:04:01-12:05:02, 12:33, 12:54, 12:60, 12:353		2	3																					
*06:02:12, 06:02:65, 06:11, 06:17, 06:52, 06:122, 06:124, 06:147, 06:236, 06:244, 06:248, 06:252, 06:280, 06:329:01:01-06:329:01:02			3														18							
*06:02:16, 06:02:43, 06:04:01-06:04:02:02, 06:06:01:01-06:06:01:02, 06:14, 06:35, 06:40, 06:72, 06:79N, 06:81, 06:123, 06:143, 06:153:01-06:153:02, 06:197, 06:203, 06:271, 06:317, 06:358		2															18							
*06:02:72		2	3												15	17	18							
*06:05		2													15	17	18							
*06:08, 06:22, 06:90		2	3														18	19						
*06:31	1	2	3														18							
*06:32, 06:171:01:01N-06:171:01:02N		2	3										13				18							
*06:76:02		2	3														17	18						
*06:87, 06:330		2	3							10							18							
*06:118		2							9								18							
*06:133, 06:161, 06:299		2	3						9								18							
*06:327		2	3													16	18							
*06:359N		2	3														18				21			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*07:01:01:01-07:01:02:14, 07:01:04-07:01:10, 07:01:12-07:01:22, 07:01:24-07:01:27, 07:01:29-07:01:73, 07:01:75-07:02:08, 07:02:10, 07:02:12-07:02:97, 07:02:99-07:02:126, 07:02:128-07:02:137, 07:02:139-07:02:141, 07:02:143-07:02:145, 07:05-07:06:01:12, 07:06:03-07:06:06, 07:08-07:09, 07:13:01:01-07:15, 07:17:01:01-07:19, 07:20:02-07:23, 07:26:01:01-07:30, 07:32N-07:33N, 07:35, 07:37-01:02-07:42, 07:44, 07:46-07:50, 07:52, 07:54-07:55N, 07:57																								

Negative Control

[illegible]