

Well No.	1	2	3	4	5	6	7	8
<b>HLA-C allele<sup>1,2,3</sup></b>								
*18:01	1	2		4				
*18:02	1		3	4				
*18:03	1	?	?	4	5			
*18:04	1	?	?	4		6		
*18:05	1	?	?				7	
*18:06	1	?	?	4	5	6		
*18:07N	1	?	?	4			7	
*18:08-18:09	1	?	?	4				
*02:22, 05:08, 05:52, 05:89, 05:106, 06:09, 06:144, 08:27, 08:29, 08:31, 12:31, 12:144, B*15:137					5			
*03:08, 03:29, 03:31, 03:246, 15:77							7	
*04:01:01:01-04:01:69, 04:03:01-04:20, 04:23- 04:93N, 04:95N-04:111, 04:113-04:122, 04:124- 04:168, 04:170N-04:217N			3					
*04:94:01-04:94:02			3		5			
*04:112, 04:123N, 04:169			3				7	
*06:02:08, 06:34:02				4				
*06:17	1							
*07:01:01:01-07:06, 07:08, 07:10-07:19, 07:21- 07:33N, 07:35-07:48, 07:50-07:63, 07:65-07:72, 07:74-07:75, 07:77-07:91, 07:93-07:95, 07:97- 07:171, 07:173-07:209, 07:211-07:237, 07:239- 07:246, 07:248-07:294, 07:296-07:314:02, 07:316- 07:327, 07:329N-07:347N, 07:349-07:389, 07:391- 07:402, 07:404-07:405, 07:407-07:459		2						
*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:210, 07:238, 07:247, 07:315, 07:328, 07:403, 07:406	1	2						
*07:20, 07:64, 07:73, 07:92, 07:96:01-07:96:02, 07:172:01-07:172:02, 07:390		2				6		
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>

Negative Control