

<b>INTERPRETATION TABLE</b>		
<b>HLA-Cw*0409N SSP typing</b>		
	<b>Well</b>	
	<b>1</b>	<b>2</b>
<b>Length of spec.</b>	<b>220</b>	<b>255</b>
<b>PCR product</b>		<b>360</b>
<b>Length of int.</b>	<b>800</b>	<b>1070</b>
<b>pos. control<sup>1</sup></b>		
<b>5'-primer(s)<sup>2</sup></b>	<b>1018</b>	<b>28</b>
	5'-gTg <sup>3'</sup>	5'-TCA <sup>3'</sup>
		<b>368</b>
		5'-gTT <sup>3'</sup>
<b>3'-primer(s)<sup>3</sup></b>	<b>1092</b>	<b>112</b>
	5'-TTA <sup>3'</sup>	5'-CCA <sup>3'</sup>
		<b>118</b>
		5'-gCT <sup>3'</sup>
		<b>218</b>
		5'-gCT <sup>3'</sup>
<b>Well No.</b>	<b>1</b>	<b>2</b>
<b>HLA-Cw allele</b>		
<b>*0409N</b>	<b>1</b>	<b>2</b>
<b>*04010101-040105, 0403-0408, 0410-0420, 0423-0432</b>		<b>2</b>
<b>HLA-Cw allele</b>		
<b>Well No.</b>	<b>1</b>	<b>2</b>