

INTERPRETATION TABLE		
HLA-Cw*0409N SSP typing		
	Well	
	1	2
Length of spec.	205	255
PCR product		360
Length of int.	800	1070
pos. control¹		
5'-primer(s)²	1018	28
	5' -gTg 3'	5' -TCA 3'
3'-primer(s)³	1092	112
	5' -TTA 3'	5' -CCA 3'
		118
		5' -gCT 3'
		218
		5' -gCT 3'
Well No.	1	2
HLA-Cw allele		
*0409N	1	2
*04010101-040106, 0403-0408, 0410-0420, 0423-0438, 0440		2
HLA-Cw allele		
Well No.	1	2