

INTERPRETATION TABLE

HLA-Cw*17 SSP subtyping

Amplification patterns of the HLA-Cw*17 alleles

	Tube					
	1	2	3	4	5	
Length of spec.	350	70	300	90	155	Length of spec.
PCR product						PCR product
Length of int.	800	1070	1070	1070	1070	Length of int.
pos. control ¹						pos. control ¹
5'-primer ²	20	28	70	20	126	5'-primer ²
	5'-CCA ^{3'}	5'-TCA ^{3'}	5'-ggA ^{3'}	5'-CCA ^{3'}	5'-ggA ^{3'}	
3'-primer ³	201	59	201	70	239	3'-primer ³
	5'-CTC ^{3'}	5'-CgA ^{3'}	5'-CTC ^{3'}	5'-ggC ^{3'}	5'-gCg ^{3'}	
Tube No.	1	2	3	4	5	Tube No.
HLA-Cw allele						HLA-Cw allele
*1701	+			+		*1701
*1702		+				*1702
*1703	+		+			*1703
*1704	?	?	?	?	+	*1704
HLA-Cw allele						HLA-Cw allele
Tube No.	1	2	3	4	5	Tube No.