

INTERPRETATION TABLE

HLA-Cw*17 SSP subtyping

Amplification patterns of the HLA-Cw*17 alleles

| | Well | | | | | |
|---------------------------|----------------------|----------------------|----------------------|----------------------|----------------------|---------------------------|
| | 1 | 2 | 3 | 4 | 5 | |
| Length of spec. | 350 | 70 | 300 | 90 | 155 | Length of spec. |
| PCR product | | | | | | PCR product |
| Length of int. | 800 | 1070 | 1070 | 1070 | 1070 | Length of int. |
| pos. control ¹ | | | | | | pos. control ¹ |
| 5'-primer ² | 20 | 28 | 70 | 20 | 126 | 5'-primer ² |
| | 5'-CCA ^{3'} | 5'-TCA ^{3'} | 5'-ggA ^{3'} | 5'-CCA ^{3'} | 5'-ggA ^{3'} | |
| 3'-primer ³ | 201 | 59 | 201 | 70 | 239 | 3'-primer ³ |
| | 5'-CTC ^{3'} | 5'-CgA ^{3'} | 5'-CTC ^{3'} | 5'-ggC ^{3'} | 5'-gCg ^{3'} | |
| Well No. | 1 | 2 | 3 | 4 | 5 | Well No. |
| HLA-Cw allele | | | | | | HLA-Cw allele |
| *1701 | 1 | | | 4 | | *1701 |
| *1702 | | 2 | | | | *1702 |
| *1703 | 1 | | 3 | | | *1703 |
| *1704 | ? | ? | ? | ? | 5 | *1704 |
| HLA-Cw allele | | | | | | HLA-Cw allele |
| Well No. | 1 | 2 | 3 | 4 | 5 | Well No. |