

INTERPRETATION TABLE					
HLA-Cw*18 SSP subtyping					
Amplification patterns of the Cw*18 alleles					
	Well				
	1	2	3	4	5
Length of spec.	425	560	225	535	165
PCR product					
Length of int.	800	1070	1070	1070	1070
pos. control¹					
5'-primer²	47	956	895	213	412
	5'-Agg 3'	5'-ggT 3'	5'-ggA 3'	5'-CCC 3'	5'-ATA 3'
3'-primer³	302	1034	956	459	538
	5'-ggT 3'	5'-AgC 3'	5'-CAg 3'	5'-AgA 3'	5'-CCA 3'
Well No.	1	2	3	4	5
HLA-Cw allele					
*1801	1	2		4	
*1802	1		3	4	
*1803	1	?	?	4	5
*0222, 0508, 0609, 0827, 0829, B*9537					5
*04010101-040112, 0403-0420, 0423-0448			3		
*0617	1				
*070101-0706, 0708, 0710-0733N, 0735-0748, 0750-0775, 0777-0791		2			
*0707, 0709, 0749, 0776	1	2			
HLA-Cw allele					
Well No.	1	2	3	4	5