

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26
KIR allele ¹																										
2DL1*001-026N	1																									
2DL2*0010101-013 ²		2																								
2DL3*0010101-032 ³			3																							
2DL4*00101-00602, 010, 01201-01202, 014-016, 018, 021-026				4																				24		
2DL4*007-00902, 011, 013, 017, 019, 020, 027				4																						
2DL5A*0010101-00105, 0050101-0050104, 01201-01202, 014-015					5	6																				
2DL5B*0020101-0020105, 0020107, 004, 00601, 00603-0070101, 0080101-00802, 009-01302					5		7																			25
2DL5B*0020106, 00202, 0070102, 00803, 01303, 016-018					5		7																			?
2DL5B*003, 00602					5		7																			
2DS1*001								8														21				
2DS1*0020101-008 ⁴								8															22			
2DS2*0010101-008									9																	
2DS3*00101-007										10																
2DS4*0010101-00104, 01101-01102, 014-016											11															
2DS4*0030101-0030104, 0040101-0040102, 0060101-0060102, 007-010, 012, 013												12														
2DS5*001-012													13													
3DL1*0010101-002, 0050101-009, 01501-018, 020010-0200102, 022-035, 038, 040-044, 051-054, 057, 059-077, 079														14												
3DL1*00401-00403, 019, 021, 036, 037, 039, 056, 072														14									23			
3DL2*0010101-063															15											
3DL3*00101-057																16										
3DS1*010-014, 045-049N, 050, 055, 058, 078, 082-085																	17									
2DP1*00101-014																		18								
3DP1*001-002, 004, 007, 0090101-00902 ⁴																			19				22			25
3DP1*0030101-0030402, 005, 006, 008, 010																			19	20						
3DP1*011-012																				19			22			?
3DP1*013-014																				19	20					?
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26

Negative Control