

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

Olerup SSP[®] HLA-A-B-DR SSP Combi Tray

Product number:	101.701-24/06 – including <i>Taq</i> pol. 101.701-24u/06u – without <i>Taq</i> pol.
Lot number:	9F4
Expiry date:	2020-07-01
Number of tests:	24 tests – Product No. 101.701-24/24u 6 tests – Product No. 101.701-06/06u
Number of wells per test:	95 +1
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
- Adhesive PCR seals	RT
- Product Insert	RT

This Product Description is only valid for Lot No. 9F4.

Complete product documentation consists of generic Instructions for Use (IFU), lot specific Product Insert, Worksheet and Certificate.

CHANGES COMPARED TO THE PREVIOUS OLERUP SSP[®] HLA-A-B-DR SSP COMBI TRAY LOT (6E8)

The **HLA-A low resolution** specificity and interpretation tables have been updated for the HLA-A alleles described since the previous Olerup SSP[®] HLA-A-B-DR SSP Combi Tray lot was made (**Lot No. 6E8**). The kit design is based on IMGT/HLA database 3.29.0.

The primers of the wells detailed below have been exchanged, modified or added compared to the previous lot.

Well	5'-primer	3'-primer	rationale
6	-	Modified	3'-primer modified for increased yield of the HLA-A*23 alleles.
8	Removed	-	5'-primer removed for the HLA-A*66:23 allele for improved resolution of the HLA-A*66 alleles.
9	Added	-	5'-primer added for the A*25:39 and A*32:93 alleles.
10	-	Added	3'-primer added for the A*26:137 allele.
12	-	Added	3'-primer added for the A*66:26Q allele.
13	Removed	-	5'-primer removed.
17	-	Added	3'-primer added for the A*03:273 allele.
18	-	Exchanged	3'-primer exchanged for the A*31:08 allele.
21	Exchanged	Modified, removed	5'-primer exchanged, 3'-primer modified, and 3'-primer removed for improved HLA-specific amplification.
22	Added	-	5'-primer added for increased yield of the HLA-A*36 alleles.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

The **HLA-B low resolution** specificity and interpretation tables have been updated for the HLA-B alleles described since the previous *Olerup* SSP[®] HLA-A-B-DR SSP Combi Tray lot was made (**Lot No. 6E8**). The kit design is based on IMGT/HLA database 3.29.0.

The primers of the wells detailed below have been exchanged, modified or added compared to the previous lot.

Well	5'-primer	3'-primer	rationale
41	-	Added	3'-primer added for improved allelic resolution.
42	-	Added	3'-primer added for the B*35:329 allele.
46	-	Added	3'-primer added for the B*38:61 allele.
52	-	Added	Modified 3'-primer added for increased yield.
56	-	Added	Modified 3'-primer added for increased yield.

The **HLA-DR low resolution** specificity and interpretation tables have been updated for the HLA-DRB alleles described since the previous *Olerup* SSP[®] HLA-A-B-DR SSP Combi Tray lot was made (**Lot No. 6E8**). The kit design is based on IMGT/HLA database 3.29.0.

The primers of the wells detailed below have been exchanged, modified or added compared to the previous lot.

Well	5'-primer	3'-primer	rationale
76	Removed	Removed, added	5'-primer removed, 3'-primers removed, 3'-primer added for improved HLA-specific amplification.
77	-	Added	3'-primer added for the DRB1*03:15:02 allele.
87	-	Added	3'-primer added for the DRB1*13:220 allele.
91	Added	-	5'-primer added for the DRB1*14:185 allele.
95	-	Added	5'-primer added for increased yield of the DRB1*01:12 allele.

Change in revision R01 compared to R00:

1. Primer mix 51 does not amplify the B*08:31, 41:06, 41:15 and 42:14 alleles, and primer mix 72 does not amplify the B*08:31 allele. This has been corrected in the Specificity and Interpretation Tables.

Change in revision R02 compared to R01:

1. The reactivities of primer mix 17 have been corrected in the Specificity and Interpretation Tables.

Change in revision R03 compared to R02:

1. Primer mix 4 amplifies the A*01:15N allele. This has been corrected in the Specificity and Interpretation Tables.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

Well **96** contains Negative Control primer pairs, that will amplify more than 95% of the *Olerup* SSP[®] HLA Class I, DRB, DQB1, DPB1 and DQA1 amplicons as well as all the amplicons generated by the control primer pairs matching the human growth hormone gene.

HLA-specific PCR product sizes range from 75 to 200 base pairs.

The PCR product generated by the positive control primer pair is 430 base pairs.

Length of PCR product	105	200	105	80	75	80	85
5'-primer¹	164	340	440	45	45	43	36
	5'-CAC ^{3'}	5'-Agg ^{3'}	5'-TTA ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-TAC ^{3'}
							36
							5'-TAT ^{3'}
3'-primer²	231	2nd I	507	59	58	57	47
	5'-TgC ^{3'}	5'-AAA ^{3'}	5'-TTg ^{3'}	5'-CTC ^{3'}	5'-ggC ^{3'}	5'-CTC ^{3'}	5'-ACA ^{3'}
							48
							5'-gCA ^{3'}
							48
							5'-gCC ^{3'}
							52
							5'-TgT ^{3'}
A*	+	+	+				
B*	+	+	+				
C*	+	+	+				
DRB1				+	+		
DRB3				+	+		
DRB5				+			
DQB1					+		
DPB1						+	
DQA1							+

¹The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

²The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

PRODUCT DESCRIPTION

HLA-A-B-DR SSP Combi Tray

CONTENT

The primer set contains 5'- and 3'-primers for grouping the HLA-A*01:01 to A*80:03 alleles into the corresponding serological groups A1 to A80.

The primer set contains 5'- and 3'-primers for grouping the B*07:02 to B*83:01 alleles into the corresponding serological groups B7 to B81 as well as primer pairs for recognizing the Bw4 and Bw6 sequence motifs.

The primer set contains 5'- and 3'-primers for grouping the DRB1*01:01 to DRB1*10:22 alleles into the corresponding serological groups DR1 to DR18 as well as primer pairs for recognizing the DRB3, DRB4 and DRB5 groups of alleles.

PLATE LAYOUT

Each test consists of 96 PCR reactions in a 96 well PCR plate.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48
49	50	51	52	53	54	55	56
57	58	59	60	61	62	63	64
65	66	67	68	69	70	71	72
73	74	75	76	77	78	79	80
81	82	83	84	85	86	87	88
89	90	91	92	93	94	95	NC

The 96 well PCR plate is marked with 'A-B-DR' in silver/gray ink.

Well No. 1 is marked with the Lot No. '9F4'.

Wells 1 to 24 – HLA-A low resolution primers.

Wells 25 to 72 – HLA-B low resolution primers.

Wells 73 to 95 – HLA-DR low resolution primers.

Well 96 – Negative Control.

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded. The PCR plates are covered with a PCR-compatible foil.

INTERPRETATION

Only HLA-A alleles will be amplified by the 24 wells of the HLA-A low resolution primer set, **wells 1 to 24**, except that a few HLA-B and HLA-C alleles will be amplified by some primer mixes. For further details see Specificity Table.

Only HLA-B alleles will be amplified by the 48 wells of the HLA-B low resolution, primer set, **wells 25 to 72**, except that a few HLA-A and HLA-C alleles will be amplified by some primer mixes. For further details see Specificity Table.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

Only HLA-DRB alleles will be amplified by the 23 wells of the DR low resolution primer set, **wells 73 to 95**. Thus, the interpretation of DR low resolution typings is not influenced by other HLA class II genes. For further details see Specificity Table.

UNIQUELY IDENTIFIED ALLELES

All the HLA-A alleles, i.e. **A*01:01 to A*80:03**, recognized by the HLA Nomenclature Committee in August 2017^{1,2} will be amplified by the primers in the HLA-A low resolution primer set, **wells 1 to 24**. The HLA-A alleles will be grouped into their corresponding serological specificities^{3,4}.

All the HLA-B alleles, i.e. **B*07:02 to B*83:01**, recognized by the HLA Nomenclature Committee in August 2017^{1,2} will be amplified by the primers in the HLA-B low resolution primer set, **wells 25 to 72**. The HLA-B alleles will be grouped into their corresponding serological specificities³.

All the HLA-DRB1, -DRB3, -DRB4⁴ and -DRB5 alleles, i.e. **DRB1*01:01:01 to DRB1*10:22, DRB3*01:01:02:01 to DRB3*03:13, DRB4*01:01:01:01 to DRB4*03:01N, DRB5*01:01:01 to DRB5*02:17**, recognized by the HLA Nomenclature Committee in August 2017^{1,2} will be amplified by the primers in the DR low resolution primer set, **wells 73 to 95**. The HLA-DRB alleles will be grouped into their corresponding serological specificities⁵.

¹HLA-A, HLA-B and HLA-DRB alleles listed on the IMGT/HLA web page 2017-August-10, release 3.29.0, www.ebi.ac.uk/imgt/hla.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

³The HLA-A, HLA-B and the HLA-DRB1, -DRB3, -DRB4 and -DRB5 alleles will be grouped into their corresponding serological specificities, except that following alleles give rise to identical amplification patterns:

Alleles	Alleles
A*11:116, 11:140, 11:199:01, 11:222, A*66:23	B*55:04, 55:08, 55:13, 55:27, 55:46, 55:49, 55:81, B*56:01:05, 56:15, 56:19N, 56:22
A*23:14:01-23:14:02, A*24:24, 24:71, 24:315	B*55:23, 55:32, B*56:18, 56:31-56:32, 56:50
A*30:01:01-30:02:11, 30:02:13-30:04:02, 30:06-30:07, 30:09-30:20, 30:22-30:54, 30:56-30:88, 30:90-30:119, B*07:260	B*57:67:01, B*58:36
A*31:08, A*33:53	DRB1*03:126, DRB1*14:179
A*31:109, A*33:125	DRB1*08:20, DRB1*13:18, 13:47, 13:55, 13:158, 13:164, 13:232
B*08:26:01-08:26:02, 08:50, 08:62, 08:85, 08:94, 08:146, B*42:07, 42:24	DRB1*08:31, 08:41, 08:75, DRB1*11:67, 11:193:01-11:193:02, 11:209
B*13:35, 13:59, 13:71, B*44:135, 44:158, 44:184	DRB1*12:57, DRB1*13:67, 13:103, 13:160, 13:195
B*14:08:01-14:08:02, 14:55, B*39:01:19, 39:25N, 39:30, 39:32-39:34, 39:47, 39:50, 39:74, 39:82, 39:102, 39:107, 39:112, 39:128	DRB1*13:13, 13:119, 13:154, 13:156, DRB1*14:84, 14:116, 14:144
B*18:29, 18:72:01-18:72:03, 18:92, 18:102, B*35:32:01-35:32:02, 35:37, 35:53N, 35:64:01-35:64:02, 35:68:01-35:68:02, 35:99, 35:118-35:119, 35:174	DRB1*13:44, 13:86, 13:206, DRB1*14:183
B*53:30, B*57:45, 57:51, 57:69	

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

⁴The A*01:26, 01:136 and 01:192 and the A*11:94, 11:112 and 11:211 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mixes 1 and 11.

The A*23:66, A*24:14:01:01-24:15, 24:51-24:53, 24:57, 24:64, 24:94, 24:114, 24:138, 24:188, 24:222N, 24:228, 24:291, 24:296, 24:304, 24:316 and 24:324 and the C*04:01:03 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 5.

The A*68:08:02, 68:105 and 68:157 alleles may be distinguished by the different sizes of the specific PCR products generated by primer mix 21.

⁵The DRB4*02:01N and DRB4*03:01N null alleles will not be amplified by the DR low resolution primer set.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information
SPECIFICITY TABLE

HLA-A low resolution primer set

Specificities and sizes of the PCR products of the 24 primer mixes used for HLA-A low resolution SSP typing

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	HLA-A serology ³	Amplified HLA-A alleles ⁴
1 ⁵	120 bp, 145 bp, 225 bp	800 bp	A1, A36	*01:01:01:01-01:04N, 01:06-01:33, 01:35-01:240N, 03:18, 03:135, 03:260, 11:94, 11:112, 11:211, 11:226, 26:120, 36:01-36:05
2	175 bp, 215 bp, 255 bp, 365 bp, 415 bp	800 bp	A2, A19, A28, A203, A210	*02:01:01:01-02:01:15, 02:01:17-02:01:19, 02:01:21-02:01:81, 02:01:83-02:22:02, 02:24:01-02:35:01, 02:35:03-02:47, 02:49-02:77, 02:78 ^w , 02:79:01-02:97:02, 02:99, 02:101:01-02:128, 02:130-02:570, 02:572-02:643N, 02:645-02:694
3 ⁶	205 bp, 235 bp	1070 bp	A1, A3, A11, A32, A34, A36	*01:12, 01:19, 01:21, 01:126, 01:200, 02:338, 03:01:01:01-03:17:02, 03:19-03:74, 03:76-03:94, 03:96-03:134, 03:136-03:176, 03:178N-03:186, 03:188-03:193, 03:195-03:214, 03:216-03:224, 03:226-03:248, 03:250-03:282, 11:03, 11:20, 11:25:01-11:25:02, 11:60, 11:130, 11:158, 11:175, 11:183, 11:209, 24:92, 32:04, 34:02:01-34:04, 34:07-34:10N, 34:13, 34:15, 36:02, 74:23
4	190 bp	800 bp	A1, A3, A11, A30, A36, A68	*01:01:01:01-01:01:22, 01:01:24-01:01:47, 01:01:49-01:01:64, 01:01:67-01:01:78, 01:01:80-01:04N, 01:06-01:07, 01:09:01-01:11N, 01:13, 01:15N-01:18N, 01:20-01:29, 01:31N-01:33, 01:35-01:78, 01:80-01:98, 01:100-01:144, 01:146, 01:148, 01:150-01:158, 01:160N-01:166, 01:168-01:177, 01:179N-01:199, 01:201-01:204, 01:206-01:207, 01:209-01:213, 01:215-01:227, 01:229-01:235, 01:237-01:240N, 02:78, 02:169, 03:12, 03:18, 03:88, 03:135, 11:01:01:01-11:27, 11:29-11:52Q, 11:54-11:264, 25:43, 26:19, 26:72, 29:67, 30:08, 32:64, 36:04, 66:23, 68:13:01, 68:66:01-68:66:02, 68:134, 74:19, C*03:365, C*12:131
5	160 bp, 335 bp, 505 bp	1070 bp	A3, A9, A23, A24, A2403, A29, A31, A32	*03:15, 03:19, 03:30, 03:152, 03:273, 11:139, 23:01:01:01-23:68, 23:70-23:83, 24:02:01:01-24:11N, 24:13:01-24:15, 24:17-24:64, 24:66-24:210, 24:212-24:384, 29:07, 29:49, 31:08, 31:29, 32:05, 32:79, 33:19, 33:53, C*04:01:03
6	135 bp, 200 bp	800 bp	A9, A23, A24, A29, A80	*11:166, 23:01:01:01-23:56, 23:58-23:65, 23:67-23:68, 23:70-23:83, 24:24, 24:71, 24:315, 29:07, 29:49, 31:29, 31:85, 32:72, 80:01:01:01-80:03, B*18:27
7	175 bp, 205 bp	1070 bp	A2, A3, A9, A23, A24, A2403, A26	*02:17:01 ^w -02:17:04 ^w , 11:139, 23:14:01-23:14:02, 24:02:01:01-24:11N, 24:13:01-24:13:02, 24:17-24:50, 24:54-24:56, 24:58-24:63, 24:66-24:91, 24:93, 24:95-24:113, 24:115-24:137, 24:139-24:187, 24:189-24:210, 24:212-24:221, 24:223-24:227, 24:229-24:290, 24:292-24:295, 24:297-24:303N, 24:305-24:315, 24:317-24:323N, 24:325-24:384, 26:16, 33:19, 33:119, 68:45, 68:117

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

Lot No.	Fragment Size	Probe Size	Alleles	SNPs
8	160 bp, 200 bp	800 bp	A2, A3, A10, A11, A25, A26, A28, A32, A34, A66, A68, A69	*01:51, 02:55, 02:527, 02:582, 02:644, 03:24, 03:50, 11:10, 11:183, 11:191, 25:01:01:01-25:16, 25:18-25:33, 25:35-25:45, 26:01:01:01-26:06, 26:08- 26:15, 26:17-26:18, 26:20-26:43:02, 26:45-26:63, 26:65-26:71N, 26:73-26:88, 26:90-26:91, 26:93- 26:145N, 29:28, 32:15, 33:51, 34:01:01-34:17, 66:01:01:01-66:22, 66:24-66:28N, 68:01:01:01- 68:166, 69:01:01:01-69:03
9^{5,7}	75 bp	800 bp	A2, A3, A25, A32	*02:81, 02:124, 25:01:01:01-25:45, 32:01:01:01- 32:02, 32:04, 32:06-32:37, 32:39-32:59, 32:61-32:70, 32:72-32:77, 32:81-32:103, B*07:81, B*08:52, B*18:67, B*38:41, B*38:67, B*51:185, B*53:05:01- 53:05:02, B*53:16, B*53:33, B*57:60
10^{5,9}	80 bp, 240 bp	800 bp	A10, A26, A43	*01:43, 01:51, 02:644, 11:17, 11:40, 11:223, 26:01:01:01-26:02:02, 26:04, 26:07:01-26:20, 26:22- 26:29, 26:31-26:43:02, 26:45-26:77, 26:79-26:91, 26:93-26:110, 26:112-26:145N, 33:13, 33:48, 43:01, 68:84
11⁵	80 bp, 175 bp, 500 bp	1070 bp	A1, A9, A10, A11, A26, A31, A34, A66	*01:13, 01:26, 01:28, 01:136, 01:176, 01:192, 01:229, 02:662, 03:63, 03:88, 11:01:01:01-11:27, 11:29-11:52Q, 11:54-11:264, 24:19, 24:44, 26:03:01, 26:06, 26:21, 26:78, 26:92, 26:111, 34:01:01-34:08, 34:10N-34:17, 66:01:01:01-66:01:03, 66:04-66:11, 66:13-66:14, 66:17-66:20, 66:22-66:24, 66:27N, 69:02, 80:02
12⁵	125 bp, 160 bp, 195 bp	800 bp	A1, A3, A10, A25, A26, A31, A34, A43, A66, A74	*01:01:56, 02:309, 02:454, 03:01:19, 11:11, 25:05- 25:06, 26:09, 26:54, 26:91, 31:03-31:04, 31:123, 34:01:01-34:17, 43:01, 66:02-66:03, 66:16, 66:19, 66:21, 66:25-66:26Q, 66:28N, 68:130:01, 74:01:03
13	180 bp, 225 bp	1070 bp	A1, A2, A3, A10, A25, A26, A34, A43, A66	*01:13, 01:176, 01:194, 02:34-02:35:03, 02:56:01- 02:56:02, 02:62, 02:103, 02:135, 02:580, 03:01:01:01-03:01:22, 03:01:24-03:07:02, 03:09- 03:11N, 03:13-03:31, 03:33-03:35, 03:37-03:40, 03:42-03:56, 03:58, 03:60-03:71, 03:73-03:87, 03:90- 03:106, 03:109-03:110, 03:112-03:141, 03:143- 03:151, 03:153-03:171, 03:174-03:175, 03:177, 03:179-03:193, 03:195-03:197N, 03:201-03:202, 03:204, 03:206-03:210, 03:212-03:218, 03:220- 03:251, 03:253-03:259, 03:261-03:266N, 03:268- 03:272, 03:274-03:282, 11:116, 11:140, 11:199:01- 11:199:02, 11:222, 25:01:01:01-25:05, 25:07-25:35, 25:37-25:42N, 25:44-25:45, 26:01:01:01-26:01:42, 26:02:01 ^w -26:02:02 ^w , 26:03:01, 26:05-26:08, 26:10- 26:28, 26:29 ^w , 26:30-26:43:02, 26:45-26:48, 26:49 ^w , 26:50-26:74, 26:76-26:77, 26:79-26:90, 26:92- 26:145N, 30:55, 31:24-31:25, 32:26:01-32:26:02, 33:61, 34:08, 43:01, 66:01:01:01-66:01:03, 66:04- 66:09, 66:10 ^w , 66:11-66:15, 66:17-66:20, 66:22- 66:24, 66:27N, 68:71, 74:13
14^{5,9}	80 bp, 115 bp, 200 bp, 240 bp, 460 bp	1070 bp	A3, A26, A29, A31, A34, A36, A68	*02:237, 02:309, 02:454, 03:01:38, 03:95, 03:123:02, 03:171, 03:231:01, 11:130, 26:19, 26:22, 29:01:01:01-29:103, 31:03-31:04, 31:123, 32:42, 34:04, 36:02, 66:09, 68:19, 68:130:02

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

15 ^{5,10,11}	90 bp, 135 bp, 205 bp	1070 bp	A1, A30	*01:13, 01:28, 01:176, 01:194, 01:229, 03:43, 03:82, 03:186, 03:260, 11:113, 11:162, 30:01:01-30:02:11, 30:02:13-30:04:02, 30:06-30:20, 30:22-30:119, 31:35, B*07:260
16	240 bp, 380 bp, 410 bp	1070 bp	A24, A31, A32	*02:237, 03:95, 29:14, 31:01:02:01-31:127, 32:05, 32:79, 33:53, 33:125, 74:13
17 ¹¹	140 bp, 180 bp, 235 bp, 260 bp	1070 bp	A1, A3, A25, A32, A74	*01:07, 01:95, 03:32, 03:43, 03:82, 03:152, 03:186, 03:219, 03:273, 23:64, 24:104, 24:243, 25:03, 25:30, 29:13, 29:39, 29:98, 30:89, 31:21, 31:35, 32:01:01:01-32:103, 74:07
18	200 bp, 390 bp	1070 bp	A24, A31, A32, A33, A68, A74	*02:243:01-02:243:02, 24:82, 29:48, 31:08, 31:109, 32:05, 32:15, 32:79, 32:93, 33:01:01:01-33:01:10, 33:03:01-33:37, 33:39-33:130, 68:29, 74:04, 74:21
19	340 bp, 375 bp	800 bp	A2, A19, A68, A74	*01:121, 02:65, 02:407, 02:449, 03:246, 03:282, 32:62, 68:25, 74:01:01-74:28
20 ^{9,11}	210 bp, 240 bp	800 bp	A2, A210, A25, A68	*02:34-02:35:03, 02:46, 02:48, 02:56:01-02:56:02, 02:62, 02:70, 02:78, 02:103, 02:129, 02:571, 02:580, 02:651, 11:199:02, 23:01:13, 24:340, 25:05, 26:54, 26:122, 34:02:04, 68:01:01:01-68:166
21 ^{5,12}	65 bp, 240 bp, 375 bp, 545 bp	800 bp	A2, A23, A24, A26, A28, A32, A66, A68, A69	*02:55, 02:149, 02:243:01-02:243:02, 02:309, 02:644, 23:57, 24:82, 25:30, 26:22, 29:80, 31:41, 32:06, 33:22, 66:06, 66:09, 68:08:02, 68:29, 68:105, 68:157, 69:01:01:01-69:03
22 ^{5,8}	85 bp, 240 bp, 400 bp	800 bp	A2, A36	*02:34-02:35:03, 02:46, 02:48, 02:56:01-02:56:02, 02:62, 02:70, 02:78, 02:103, 02:129, 02:571, 02:576, 02:580, 02:682, 03:187, 11:155, 11:226, 31:62, 36:01-36:05, 68:11N, 68:41, B*40:359, B*57:65, C*04:31, C*06:137, C*07:569
23 ^{5,9}	85 bp, 160 bp, 240 bp, 495 bp	800 bp	A2, A3, A24, A26, A28, A36, A68, A80	*02:55, 02:237, 02:671, 03:41, 03:63, 03:75, 03:88, 03:95, 03:177, 11:130, 24:18, 24:204, 24:213, 26:03:01, 26:05-26:06, 26:21, 26:30, 26:78, 26:111, 33:24, 36:02, 68:05, 68:15, 68:20, 68:109, 68:136, 80:01:01:01-80:01:01:02, 80:03, C*06:187
24 ^{10,11}	360 bp, 445 bp	1070 bp	A3, A10, A26, A31, A66, A74	*02:135, 02:309, 02:454, 03:01:19, 25:13, 26:07:01- 26:07:02, 26:30, 26:65, 26:92, 31:04, 31:123, 34:09, 66:02-66:03, 66:12, 66:16, 66:21, 66:25-66:26Q, 66:28N, 74:01:03

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-A low resolution SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length. Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings. PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings. Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: 9F4

Lot-specific information

1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³The serological reactivity of all HLA-A alleles is not known. In this table we use the expert-assigned serological grouping in Tissue Antigens (2009) 73:95-170 and the serological grouping of the sequence-defined allele.

⁴For several HLA Class I alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁵HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁶Primer mix 3 may faintly amplify the A*30:04:01-30:04:02, 30:06, 30:09, 30:17, 30:29, 30:46, 30:77, 30:90, 30:99, 30:103, 30:105 and 30:117 alleles.

⁷Primer mix 9 may weakly amplify the A*34 alleles.

⁸Primer mix 22 might faintly amplify most A*11 alleles.

⁹Primer mixes 10, 14, 20 and 23 may have a tendency of giving rise to primer oligomer formation.

¹⁰Primer mixes 15 and 24 may give rise to a lower yield of HLA-specific PCR product than the other HLA-A low primer mixes.

¹¹Primer mixes 15, 17, 20 and 24 may have tendencies of unspecific amplifications.

¹²Primer mix 21 may generate a false positive band of about 560 base pairs. This band should be disregarded when interpreting HLA-A low resolution typings.

‘w’, might be weakly amplified.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information
SPECIFICITY TABLE

HLA-B low resolution primer set

Specificities and sizes of the PCR products of the 48 primer mixes used for HLA-B low resolution SSP typing

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	HLA-B serol. ³	Amplified HLA-B alleles ⁴
25^{5,9}	110 bp	800 bp	7, 703, 35, 40, 41, 42, 48, 61	*07:02:01:01-07:18:02, 07:20-07:32, 07:34-07:39, 07:41-07:47, 07:49N-07:52, 07:54-07:59, 07:61-07:99, 07:101-07:121, 07:123-07:138, 07:140-07:163, 07:166-07:208, 07:210-07:222, 07:224-07:225, 07:227-07:234, 07:236-07:296, 08:79, 08:156, 15:138, 15:230, 15:241, 15:297, 35:66, 37:07, 40:15-40:16, 40:23, 40:32, 40:98, 40:128, 40:136, 40:158, 40:183, 40:198, 40:231, 40:279, 40:355, 41:08, 42:05:01-42:05:02, 44:150, 48:05, 48:08, 48:15, C*02:23, C*03:260, C*04:77
26	215 bp	1070 bp	8, 44	*08:01:01:01-08:05, 08:07-08:25, 08:27-08:49, 08:51-08:61, 08:63-08:64, 08:66-08:84, 08:86N-08:93, 08:95-08:122, 08:124-08:131, 08:133-08:145, 08:147-08:184, 13:93, 15:142, 15:180, 15:429, 35:218, 35:256, 37:62, 40:328, 44:49, 44:156, 49:25, 51:68, 51:176, C*06:147
27	140 bp, 235 bp	1070 bp	7, 8, 13, 15, 35, 4005, 44, 49, 61, 62, 77	*07:20, 07:24, 07:60, 07:100, 07:131, 07:220, 07:223, 08:21, 08:25, 08:175, 13:01:01:01-13:04, 13:06-13:08, 13:10-13:23, 13:25-13:38, 13:40-13:57, 13:59-13:66, 13:68-13:96, 13:99-13:103N, 15:07:01:01-15:07:03, 15:36 ^w , 15:45, 15:55, 15:68, 15:89 ^w , 15:126, 15:197:02, 15:207, 15:324, 15:329, 15:331, 15:360, 15:383, 15:405, 35:05:01:01-35:05:03, 35:16-35:17:02, 35:22, 35:30, 35:51, 35:58, 35:72, 35:89, 35:97, 35:113-35:114, 35:125, 35:164, 35:199, 35:232, 35:242, 35:318, 35:346, 40:05:01:01-40:05:01:02, 40:71, 40:174, 44:08 ^w , 44:54, 44:57 ^w , 44:60, 44:106, 44:110, 44:135, 44:158, 44:184, 44:211, 44:213, 46:12, 51:64, 51:148, 53:14, 55:51, 58:18, C*03:05, C*03:25, C*03:27, C*03:143, C*03:167, C*03:296, C*03:335, C*03:346
28^{6,7}	140 bp, 265 bp	800 bp	12, 13, 14, 15, 17, 18, 21, 27, 35, 39, 40, 41, 44, 45, 47, 48, 49, 50, 57, 58, 60, 61, 62, 64, 65	*07:68:01-07:68:03, 07:100, 07:218, 07:220, 07:237, 07:262, 08:71, 08:118, 13:01:01:01-13:01:09, 13:01:11-13:04, 13:06-13:13:02, 13:15-13:23, 13:25-13:57, 13:59-13:82, 13:84-13:92, 13:94-13:96, 13:99-13:101, 13:102 ^w , 13:103N, 14:01:01:01-14:02:08, 14:02:10-14:05, 14:07N, 14:09, 14:11-14:36, 14:38-14:54, 14:57, 15:07:01:01-15:07:03, 15:45-15:46, 15:53, 15:55, 15:68, 15:106, 15:126, 15:197:02, 15:207, 15:212, 15:324, 15:329, 15:331, 15:336, 15:360, 15:383, 15:399, 15:405, 18:14, 18:44:01-18:44:02, 18:107, 27:32, 27:138, 35:05:01:01-35:05:03, 35:22, 35:46-35:47, 35:51, 35:58, 35:63, 35:66, 35:72, 35:89, 35:97, 35:154, 35:199, 35:207, 35:217, 35:232, 35:317, 38:19, 39:03:01:01-39:03:01:02, 39:14, 39:24:01-39:24:02, 39:29, 39:37, 39:76, 39:79, 39:120, 40:01:01-40:01:39, 40:01:42-40:01:49, 40:07, 40:10:01:01-40:10:02, 40:14:01-40:16, 40:22N-40:23, 40:25, 40:30-40:36, 40:38, 40:42-40:43, 40:45, 40:47-40:49, 40:51-40:55, 40:58-40:63, 40:65-40:67, 40:69, 40:71-40:73, 40:76-40:77, 40:79-40:81, 40:84, 40:87:01-40:88, 40:92, 40:100-40:102, 40:106, 40:108, 40:110, 40:112-40:114:02, 40:116-40:118N, 40:121,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				40:123-40:126, 40:128-40:130:02, 40:132, 40:134-40:136, 40:139-40:141, 40:146-40:147, 40:150-40:156, 40:158, 40:160:01-40:160:02, 40:163-40:164, 40:166, 40:168, 40:171-40:172, 40:174-40:175, 40:178-40:179, 40:182-40:183, 40:186:01-40:188, 40:191-40:199, 40:204, 40:207-40:208, 40:212-40:213, 40:215-40:219, 40:221-40:223, 40:227-40:228, 40:231, 40:233-40:242, 40:245, 40:247, 40:249-40:250, 40:252-40:253, 40:257-40:264, 40:272-40:273, 40:277-40:282, 40:285-40:286N, 40:288, 40:298-40:301, 40:308, 40:310, 40:312, 40:315-40:316, 40:319, 40:321, 40:323-40:327, 40:329, 40:332-40:333, 40:336, 40:338N-40:339, 40:342, 40:344, 40:346, 40:352-40:353, 40:355, 40:357-40:358, 40:360, 41:01:01-41:40, 41:42-41:52, 42:10, 42:17, 44:02:01:01-44:02:31, 44:02:33-44:02:39, 44:02:41-44:03:06, 44:03:08-44:03:15, 44:03:17-44:03:18, 44:03:21-44:03:29, 44:03:32-44:05:04, 44:09-44:39, 44:41:01-44:43:02, 44:45-44:56N, 44:58N-44:80, 44:82, 44:84:01-44:102, 44:104-44:110, 44:112-44:129, 44:132-44:133, 44:135-44:144, 44:146-44:155, 44:157-44:193, 44:195N-44:199, 44:201-44:209, 44:211, 44:213-44:223, 44:226-44:230, 44:231 ^w , 44:232-44:245, 44:247-44:271, 45:01:01-45:20, 46:12, 47:01:01:01-47:10, 48:12, 48:14, 48:29, 49:01:01:01-49:02, 49:04:01-49:14, 49:16-49:24, 49:27-49:30, 49:32-49:41, 49:43-49:48, 50:01:01:01-50:02, 50:04-50:06, 50:08-50:20, 50:31-50:33, 50:35-50:42, 50:44-50:53, 51:64, 53:14, 53:30, 54:38, 57:11, 57:45, 57:51, 57:69, 57:71, 57:88, 58:02:01-58:02:02, 58:06-58:07, 58:18, 58:25, 58:38, 58:43, 58:60, C*01:30, C*01:32:01-01:32:02, C*06:20, C*12:87, C*14:76, C*14:82, C*16:85, C*16:98, C*16:102, C*16:110
29⁷	185 bp, 235 bp	800 bp	7, 17, 27, 2708, 37, 38, 39, 3901, 3902, 64, 67	*07:02:32, 07:197, 14:01:01:01-14:01:05, 14:07N-14:08:02, 14:10, 14:12, 14:14, 14:19, 14:26, 14:32, 14:40, 14:46-14:47, 14:49, 14:53-14:55, 27:01-27:05:15, 27:05:17-27:05:19, 27:05:21-27:06, 27:07:04, 27:08-27:10, 27:12-27:13, 27:16-27:18, 27:20, 27:23, 27:26-27:27, 27:29, 27:31, 27:35-27:37, 27:39-27:42, 27:44-27:46, 27:48-27:61, 27:64N-27:69, 27:72-27:75, 27:77-27:80, 27:82-27:124, 27:126, 27:128-27:129, 27:131-27:137, 27:140-27:149, 27:151-27:152, 27:154-27:163, 37:02, 38:01:01-38:01:07, 38:01:09-38:02:02, 38:02:04-38:03, 38:07-38:24, 38:26-38:32, 38:34N-38:68Q, 39:01:01:01-39:01:01:07, 39:01:03:01-39:01:08, 39:01:10-39:02:01, 39:03:01:01-39:03:01:02, 39:05:01-39:06:02:04, 39:06:04-39:09:03, 39:11, 39:14-39:15, 39:18, 39:19:02, 39:22, 39:24:01-39:39:01, 39:40:01N-39:44, 39:46-39:48, 39:50-39:62, 39:64-39:71, 39:73-39:75, 39:77, 39:79-39:82, 39:84-39:88, 39:90-39:91, 39:93-39:95N, 39:97N-39:98, 39:100, 39:102-39:104, 39:106-39:109, 39:111-39:125, 39:127-39:128, 47:04-47:05, 48:21, 48:26, 67:01:01, 67:03-67:04, 67:06-67:07
30	190 bp	800 bp	14, 35, 38, 39, 65	*07:28, 08:87, 14:02:01:01-14:02:02, 14:02:04-14:02:05, 14:02:07-14:06:02, 14:09, 14:11, 14:13, 14:15-14:18, 14:20, 14:22-14:23, 14:25, 14:27, 14:29-14:31, 14:33-14:39, 14:41N, 14:43-14:45, 14:48, 14:50-14:52, 14:56-14:57, 15:77, 15:189, 15:233, 35:03:17, 35:26, 38:05, 38:33, 39:04, 44:16, 44:37:01-44:37:02, 44:64:01-44:64:02, 44:91, 44:132, 44:150, 44:190, 44:268, 52:26, 57:04:01-57:04:02, A*23:31, A*24:106, C*01:115, C*02:97, C*07:231, C*14:79, C*16:10
31	180 bp, 290 bp	1070 bp	15, 22, 46, 62, 63, 71,	*15:01:01:01-15:01:03, 15:01:06-15:02:11, 15:03:03, 15:04:01:01-15:08:02, 15:10:02, 15:11:01-15:11:03,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

			72, 75, 76, 77	15:11:05-15:15, 15:17:01:01-15:17:06, 15:19-15:21, 15:24:01-15:28, 15:30-15:36, 15:38:01-15:40:02, 15:42-15:46, 15:48, 15:50, 15:55-15:58, 15:60, 15:63, 15:65-15:66, 15:70-15:71, 15:73, 15:75-15:79N, 15:81-15:89, 15:92, 15:94N, 15:96-15:97, 15:101-15:102, 15:104-15:107, 15:109-15:113, 15:116-15:118, 15:120-15:122, 15:125-15:126, 15:128-15:129, 15:135-15:150, 15:152, 15:154-15:155, 15:157, 15:159-15:160, 15:162-15:172, 15:174-15:175, 15:177-15:179:02, 15:181N-15:185, 15:187-15:196, 15:199, 15:201-15:209N, 15:211, 15:213-15:219, 15:223-15:225, 15:227-15:228, 15:230-15:234, 15:236-15:237, 15:239-15:241, 15:244-15:251, 15:256-15:262N, 15:264-15:265, 15:267-15:273, 15:276-15:280, 15:283-15:289, 15:291, 15:295-15:305, 15:308-15:310, 15:315-15:322, 15:324-15:328, 15:330-15:334, 15:339-15:346, 15:348-15:350, 15:352-15:353, 15:355-15:361, 15:363:01-15:368, 15:370-15:375N, 15:377Q-15:378, 15:381, 15:383-15:387, 15:390-15:396, 15:398-15:407, 15:409-15:413, 15:415-15:417, 15:419-15:425, 15:428-15:429, 46:01:01-46:70, 46:72, C*01:73, C*05:01:28, C*07:02:30, C*08:16:02
32⁸	165 bp, 225 bp, 285 bp, 330 bp	1070 bp	5, 8, 12, 21, 22, 37, 41, 42, 44, 45, 48, 51, 56, 57, 60, 62, 70, 71, 72, 82	*08:01:01:01-08:01:22, 08:01:24-08:05, 08:08N-08:13, 08:15-08:27, 08:29-08:36, 08:38-08:48, 08:50-08:73, 08:75-08:76, 08:78, 08:80-08:86N, 08:88, 08:90-08:106, 08:109-08:110, 08:112-08:155, 08:157-08:163, 08:165-08:184, 13:46, 15:03:01:01-15:03:06, 15:09:01-15:10:04, 15:18:01:01-15:18:06, 15:23, 15:29, 15:37, 15:46-15:47:02, 15:49, 15:51-15:54, 15:61-15:62, 15:64:01-15:64:02, 15:69, 15:72, 15:74, 15:80, 15:83, 15:90-15:91, 15:93, 15:98-15:99, 15:103, 15:106, 15:108, 15:114-15:115, 15:119, 15:123-15:124, 15:127, 15:131-15:134, 15:143, 15:151, 15:153, 15:156, 15:158, 15:161, 15:173, 15:176, 15:186, 15:197:01-15:198, 15:200, 15:210, 15:212, 15:220-15:221, 15:226N, 15:229, 15:235, 15:238, 15:242:01-15:243, 15:252-15:253, 15:255, 15:259, 15:263, 15:266, 15:274-15:275:02, 15:281-15:282, 15:290, 15:292-15:294N, 15:306-15:307, 15:311-15:314, 15:323, 15:329, 15:335, 15:337-15:338, 15:351, 15:354, 15:369, 15:379-15:380N, 15:382, 15:388-15:389, 15:397, 15:414, 15:426-15:427, 35:87, 37:12, 38:30, 40:12, 40:136, 40:231, 40:279, 41:01:01-41:03:02, 41:05-41:09, 41:11-41:17, 41:19-41:48, 41:50, 42:01:01-42:02:01:02, 42:04-42:08, 42:10-42:15, 42:17-42:19, 42:21-42:24, 44:14-44:15, 44:18, 44:20, 44:62, 44:100, 44:166, 44:213, 45:01:01-45:20, 51:08:01-51:08:02, 51:20, 51:36, 51:44N, 51:97, 51:141, 51:153, 52:19, 55:20, 55:56, 56:13, 57:09, 82:01-82:03, C*03:342, C*07:46, C*07:356, C*07:546, C*12:184, C*15:99, C*16:85
33^{7,8,10}	165 bp, 190 bp, 390 bp	1070 bp	5, 17, 21, 35, 51, 5102, 5103, 52, 56, 62, 63, 70, 78	*13:62, 15:01:02, 15:09:01-15:09:02, 15:16:01:01-15:17:06, 15:67, 15:95, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:222, 15:230, 15:243, 15:254, 15:268, 15:273, 15:277-15:278, 15:356, 15:361-15:362, 15:396, 15:403, 15:408, 15:411, 15:423-15:424, 27:153, 35:01:10, 35:04:02, 35:251, 40:26, 40:28, 44:03:40, 44:62, 49:18:01, 50:14, 51:01:01:01-51:01:03, 51:01:05-51:02:03, 51:02:05-51:09:03, 51:11N-51:12, 51:13:02-51:14, 51:16-51:23, 51:24:03-51:24:05, 51:26-51:33, 51:37-51:41N, 51:43-51:44N, 51:46, 51:48-51:61:02, 51:63:01, 51:64-51:80, 51:82-51:91, 51:92:02, 51:94-51:98N, 51:100-51:130, 51:132-51:142, 51:145-51:147, 51:149N-51:156, 51:158:01-51:196, 51:198, 51:200-51:219,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				51:221-51:224, 51:226, 52:01:01:01-52:01:04, 52:01:06-52:01:14, 52:01:16-52:13, 52:15-52:24, 52:25:02-52:57, 52:59-52:66, 52:68-52:71, 53:28, 56:05:01-56:06, 56:21, 58:08:01-58:08:02, 78:01:01:01-78:09
34	180 bp	800 bp	13, 22, 55, 60	*07:78, 13:01:01:01-13:02:21, 13:07N-13:09, 13:11, 13:14-13:20, 13:22:01-13:23, 13:25, 13:27-13:34, 13:36-13:45, 13:47, 13:49N-13:50, 13:52, 13:54-13:58, 13:60-13:61, 13:63N-13:70, 13:73-13:78, 13:80-13:85, 13:87-13:103N, 40:48, 45:10, 49:07, 50:31, 54:26, 55:09, 55:22, 55:24, 59:09
35⁵	100 bp, 195 bp	1070 bp	8, 12, 27, 38, 39, 3902, 44, 45, 48, 60, 70, 71, 72, 82	*07:27, 07:50, 07:219, 08:04:01-08:04:02, 08:17, 08:54:01-08:54:02, 08:110, 13:97, 15:03:01:01-15:03:06, 15:47:01-15:47:02, 15:49, 15:54, 15:61-15:62, 15:64:01-15:64:02, 15:68-15:69, 15:91, 15:98, 15:103, 15:123, 15:127, 15:131-15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220, 15:235, 15:242:01-15:242:02, 15:253, 15:259, 15:266, 15:274, 15:281-15:282, 15:369, 15:379, 15:389, 15:397, 18:37, 18:125, 27:18, 27:29, 37:28, 38:03, 39:02:01-39:02:03, 39:08, 39:13:01-39:13:02, 39:23, 39:39:01-39:39:02, 39:49, 39:88, 39:101, 39:105, 39:110, 40:12, 40:46, 40:93, 42:11, 44:10, 44:15, 44:18, 44:40, 44:44, 44:117, 44:130, 44:140, 44:156, 44:210, 44:251, 45:01:01-45:03, 45:05-45:07, 45:11-45:20, 48:01:01:01-48:05, 48:07-48:30, 48:32-48:40, 49:20, 50:02, 52:16, 55:18, 56:44, 58:44, 82:01-82:03, C*05:102, C*08:93, C*15:25
36	280 bp	1070 bp	5, 7, 8, 13, 15, 16, 17, 18, 22, 35, 40, 42, 44, 46, 48, 49, 51, 53, 55, 56, 57, 58, 60, 61, 62, 63, 70, 71, 72, 75, 76, 77, 78	*07:09, 07:11, 07:17, 07:162, 07:235, 07:237, 07:255, 07:287, 08:28, 08:35, 08:37, 08:89, 08:107, 13:04, 13:10, 13:26:01-13:26:02, 13:72, 14:56, 15:01:01:01-15:01:04, 15:01:06-15:08:02, 15:11:01-15:16:03, 15:18:01:01-15:21, 15:23-15:29, 15:31-15:36, 15:38:01-15:40:02, 15:43-15:44, 15:46-15:47:02, 15:49-15:57, 15:60-15:62, 15:64:01-15:72, 15:74-15:76, 15:78:01-15:82, 15:84-15:85, 15:87-15:89, 15:91-15:98, 15:101-15:129, 15:131-15:132, 15:134-15:136, 15:138-15:149N, 15:151-15:161, 15:163-15:167, 15:169-15:176, 15:178-15:187, 15:189-15:195, 15:197:01-15:207, 15:209N-15:215, 15:217-15:223, 15:225-15:228, 15:231-15:232, 15:234-15:242:02, 15:244-15:247, 15:249-15:251, 15:254-15:268, 15:270-15:272N, 15:274-15:286, 15:288-15:302N, 15:304N-15:311, 15:314-15:322, 15:325-15:328, 15:330-15:333, 15:335-15:337, 15:340-15:341, 15:343-15:355, 15:357-15:360, 15:362-15:395, 15:397-15:399, 15:401-15:402, 15:404-15:410, 15:412-15:415, 15:417-15:419, 15:421-15:422, 15:425-15:426, 15:428-15:429, 18:01:01:01-18:15, 18:17N-18:25, 18:27-18:40, 18:42-18:60, 18:62-18:121, 18:123-18:139, 27:41, 27:107, 35:01:01:01-35:01:47, 35:05:01:01-35:05:03, 35:07-35:08:09, 35:10-35:11:03, 35:14:01-35:17:02, 35:19-35:21, 35:23-35:30, 35:32:01-35:32:02, 35:35, 35:37, 35:40N-35:43:01, 35:45-35:54, 35:57-35:58, 35:61:01-35:64:02, 35:66-35:69, 35:71-35:72, 35:76-35:80, 35:82, 35:86, 35:89-35:94, 35:97, 35:99-35:105, 35:107-35:108:02, 35:110-35:126, 35:130N-35:135, 35:137-35:148, 35:158-35:159, 35:161, 35:165N-35:166, 35:168, 35:170-35:171, 35:173N-35:178, 35:180, 35:185-35:192, 35:194, 35:196-35:197, 35:200, 35:202-35:203, 35:206-35:210, 35:212, 35:215-35:219, 35:222, 35:224-35:229, 35:232, 35:238-35:241, 35:244-35:245, 35:247-35:250, 35:252-35:255, 35:259-35:260, 35:262-35:265, 35:268-35:269, 35:271-35:272, 35:275, 35:277, 35:280, 35:283-35:284, 35:286-35:289, 35:293-35:295, 35:297, 35:299-35:308, 35:310,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				35:313-35:315, 35:317-35:319, 35:322, 35:324-35:325, 35:327-35:328, 35:330-35:334, 35:336-35:338, 35:340-35:343, 35:346-35:348, 37:35, 37:56, 38:52, 39:07, 39:43, 39:91, 40:03:01:01-40:03:01:02, 40:20, 40:38, 40:52, 40:59-40:60, 40:105, 40:158-40:159, 40:222, 40:267, 40:284, 40:352, 41:24, 41:30, 41:51, 42:09, 44:17, 44:43:01-44:43:02, 44:144, 44:230, 44:256, 45:09, 46:01:01-46:10, 46:12-46:17, 46:19-46:60, 46:62-46:69, 46:71-46:72, 48:02:01-48:02:03, 48:14, 48:23, 48:25, 49:04:01-49:05, 50:33, 51:37, 51:45, 51:63:01-51:63:02, 51:97, 51:216, 51:220, 52:39, 53:01:01-53:03, 53:05:01-53:06, 53:08:01-53:18, 53:20-53:30, 53:32, 53:34-53:35, 53:37-53:41, 53:43-53:48N, 54:06, 54:09, 54:14, 55:14, 55:23, 55:32, 55:68, 56:03, 56:09, 56:18, 56:31-56:32, 56:50, 57:01:01-57:01:24, 57:06, 57:08, 57:10-57:11, 57:13-57:16, 57:18-57:27, 57:29-57:31, 57:33-57:38, 57:40-57:41, 57:43-57:45, 57:47-57:56, 57:58-57:62, 57:64-57:65, 57:67:01-57:69, 57:71-57:75, 57:77-57:79N, 57:81, 57:83, 57:85-57:90, 58:01:01:01-58:02:02, 58:04-58:07, 58:09-58:19, 58:21-58:26, 58:29, 58:31N-58:33, 58:35-58:64, 58:66-58:72N, 58:74-58:91
37	195 bp	1070 bp	27, 35, 37, 44, 58, 76	*07:27, 07:236, 07:240, 08:126, 08:131, 08:139, 15:12, 15:14, 15:19, 15:91, 15:131, 15:161, 15:270, 15:298, 15:304N, 18:54, 18:56, 18:108, 18:132, 35:45, 35:71, 37:01:01-37:09, 37:12-37:13, 37:15-37:21, 37:23-37:28, 37:30N-37:33N, 37:35-37:36, 37:38-37:66, 37:68, 38:17, 44:17, 44:43:01-44:43:02, 44:144, 44:230, 44:256, 45:09, 46:17, 53:22, 58:07, A*24:174, C*02:82, C*06:165, C*14:61
38^{5,8}	105 bp, 395 bp, 435 bp	1070 bp	5, 7, 15, 41, 42, 62, 63, 70, 71, 75, 77	*07:04, 07:25, 07:146, 15:09:01-15:10:04, 15:13:01-15:13:02, 15:16:01:01-15:18:06, 15:21, 15:23-15:24:02, 15:37, 15:44, 15:51-15:52, 15:66-15:67, 15:72, 15:76, 15:80, 15:87, 15:90, 15:93, 15:95, 15:99, 15:101, 15:108, 15:114-15:115, 15:119, 15:124, 15:133-15:134, 15:153, 15:157, 15:161-15:162, 15:168, 15:176-15:177, 15:186, 15:189, 15:196-15:198, 15:200, 15:208, 15:216, 15:221-15:222, 15:226N, 15:229-15:230, 15:238, 15:243, 15:252, 15:254-15:255, 15:263, 15:268, 15:273, 15:275:01-15:275:02, 15:290, 15:292-15:294N, 15:306-15:307, 15:311-15:314, 15:323, 15:329, 15:335, 15:337-15:338, 15:351, 15:354, 15:356, 15:361-15:362, 15:376, 15:380N, 15:382, 15:388, 15:394, 15:396, 15:403-15:404, 15:408, 15:411, 15:414, 15:418, 15:420, 15:423-15:424, 15:426-15:427, 40:136, 40:231, 40:279, 41:08, 42:05:01-42:05:02, C*07:335, C*14:62
39⁵	115 bp, 150 bp	1070 bp	18	*08:176, 18:01:01:01-18:15, 18:17N-18:28, 18:30-18:71, 18:73-18:91, 18:93-18:101, 18:103-18:106, 18:108-18:139, 35:289, 48:28
40⁵	80 bp	1070 bp	7, 27, 2708, 44, 60	*07:73, 08:88, 08:102, 27:01-27:05:22, 27:05:24-27:21, 27:23-27:51, 27:53-27:66N, 27:68-27:74, 27:76-27:139, 27:141-27:164, 38:22, 40:46, 40:93, 44:40, 44:44, 44:130, 44:210, 49:26
41⁵	100 bp, 150 bp	800 bp	12, 13, 17, 18, 22, 27, 35, 37, 39, 44, 45, 47, 48, 51, 53, 56, 57, 58, 60, 61, 62,	*07:137, 07:235, 08:49, 08:60, 08:76, 08:129, 08:181, 13:01:01:01-13:01:12, 13:06-13:07N, 13:12:01-13:13:02, 13:17, 13:20-13:23, 13:25-13:26:02, 13:28-13:29, 13:36, 13:39, 13:43, 13:50-13:52, 13:57, 13:60-13:61, 13:63N, 13:73, 13:76N-13:80, 13:83, 13:86-13:87, 13:90, 13:92, 13:95, 14:10, 15:02:01-15:02:11, 15:13:01-15:13:02, 15:20-15:21, 15:25:01-15:25:03, 15:36, 15:44, 15:62, 15:77, 15:80, 15:85, 15:88-15:89, 15:106, 15:112, 15:121, 15:139, 15:144, 15:154, 15:165,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

			70, 75, 77, 78	15:170, 15:194, 15:204, 15:213-15:214, 15:223, 15:240, 15:250, 15:265, 15:271, 15:283, 15:289, 15:291, 15:297, 15:301-15:302N, 15:308, 15:319, 15:325, 15:328, 15:330, 15:341, 15:345, 15:357-15:358, 15:366, 15:374, 15:378, 15:393-15:394, 15:399, 15:402, 15:404, 15:407, 15:418, 15:420, 15:425, 18:22, 18:69, 18:105, 27:19, 27:30, 27:127, 27:150, 35:01:01:01-35:04:03, 35:06-35:21, 35:23-35:30, 35:33-35:36, 35:38-35:42:02, 35:45-35:50, 35:52, 35:54-35:57, 35:59:01-35:59:02, 35:61:01-35:63, 35:65Q, 35:69-35:71, 35:74, 35:76-35:78, 35:80-35:85, 35:88, 35:90-35:96, 35:98, 35:100-35:101:02, 35:103-35:113, 35:115-35:116, 35:120-35:126, 35:128-35:134N, 35:136-35:150:02, 35:152-35:173N, 35:175-35:184, 35:186-35:198, 35:200-35:204, 35:206-35:231, 35:233, 35:236-35:250, 35:252-35:259, 35:261-35:264, 35:266-35:272, 35:274-35:275, 35:277-35:291, 35:293-35:330, 35:332-35:349, 37:01:01-37:01:07, 37:01:09, 37:01:11, 37:03N-37:06:02, 37:08, 37:10-37:43, 37:45-37:50, 37:52-37:68, 38:20, 39:42, 40:04, 40:28, 40:30, 40:34, 40:64, 40:68, 40:99, 40:129, 40:137, 40:160:01-40:160:02, 40:343, 41:48, 42:13, 44:02:01:01-44:02:33, 44:02:35-44:14, 44:16-44:17, 44:19N, 44:21-44:30, 44:32-44:40, 44:42-44:46, 44:48-44:52N, 44:55-44:64:02, 44:66-44:98, 44:101-44:105, 44:107-44:134, 44:136-44:137, 44:139-44:157, 44:159-44:165, 44:167-44:181, 44:183, 44:185-44:191, 44:194-44:196, 44:198N-44:212, 44:214-44:229, 44:231-44:235, 44:237N-44:245, 44:247-44:271, 46:33, 46:67, 48:02:01-48:02:03, 48:25, 50:17, 51:04, 51:42, 51:46, 51:56:01-51:56:03, 51:139, 51:186, 53:01:01-53:13, 53:15-53:43, 53:45-53:46, 53:48N, 55:14, 56:09, 56:11-56:12, 57:01:01-57:01:04, 57:01:06-57:10, 57:12, 57:14:01-57:20, 57:22-57:30, 57:32-57:46, 57:48-57:81, 57:84-57:90, 58:01:01:01-58:01:02, 58:01:04-58:01:15, 58:01:17-58:01:22, 58:04-58:05, 58:09-58:15, 58:17N, 58:19, 58:21-58:24, 58:28:01-58:29, 58:31N-58:37, 58:39N-58:42, 58:44-58:45:02, 58:47-58:59:02, 58:62-58:91, 83:01, C*03:102, C*03:263:01-03:263:02, C*03:317, C*15:24, C*15:126
42^{7,8}	140 bp	1070 bp	8, 18, 35, 37, 39, 56, 78	*07:65, 07:134, 08:32, 08:133, 08:180, 15:202, 15:239, 15:308, 15:376, 18:01:01:01-18:08, 18:10-18:15, 18:17N-18:36, 18:38-18:47, 18:49-18:53, 18:55-18:65, 18:68-18:106, 18:108-18:124, 18:126-18:135, 18:137-18:139, 35:01:01:01-35:18, 35:20:01-35:24:02, 35:27-35:45, 35:48, 35:50-35:62, 35:64:01-35:72, 35:74-35:153, 35:155-35:197, 35:199-35:206, 35:208-35:209, 35:211-35:216N, 35:218-35:253, 35:255-35:270, 35:272-35:273, 35:275-35:310, 35:312-35:316, 35:318-35:349, 37:05, 37:11, 37:14, 37:29, 37:37, 37:67, 39:19:01-39:19:02, 40:185, 40:209, 40:308, 40:350, 46:40:01-46:40:02, 56:06, 56:45, 58:64, 78:01:01:01-78:09
43^{5,9}	55 bp, 245 bp, 390 bp, 415 bp	1070 bp	5, 15, 44, 48, 60, 62, 63, 70, 71, 72, 75, 76, 77, 82	*15:01:01:01-15:01:04, 15:01:06-15:21, 15:23-15:40:02, 15:42-15:58, 15:60-15:67, 15:69-15:99, 15:102-15:129, 15:131-15:179:02, 15:181N-15:235, 15:237-15:254, 15:256-15:429, 40:12, 46:50, 82:01-82:03, C*03:179, C*06:147, C*14:61, C*16:85
44⁹	210 bp	800 bp	27, 37	*07:27, 07:236, 07:273, 08:49, 08:126, 18:54, 37:01:01-37:68, 38:17, C*07:302
45	170 bp	800 bp	16, 35, 38, 39, 3901, 3902, 67, 72	*08:55, 15:69, 15:186, 35:35, 35:276, 38:01:01-38:02:05, 38:02:07-38:09, 38:11-38:21, 38:23-38:29, 38:31-38:35, 38:38-38:44, 38:46-38:49, 38:51, 38:53-38:62, 38:64, 38:66-38:68Q, 39:01:01:01-39:01:01:07, 39:01:03:01-39:01:18,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				39:01:20-39:20, 39:22-39:24:02, 39:26-39:29, 39:31, 39:35, 39:37-39:42, 39:44-39:46, 39:48-39:49, 39:51-39:73, 39:75-39:81, 39:83-39:101, 39:103-39:106, 39:108-39:111, 39:113-39:127, 40:184, 50:49, 51:101, 52:71, 58:20, 58:84, 67:01:01-67:07, C*15:51
46⁵	110 bp	1070 bp	38	*27:07:04, 38:01:01-38:01:07, 38:01:09-38:02:02, 38:02:04-38:03, 38:05, 38:08-38:24, 38:26-38:34N, 38:36-38:68Q, 49:26, 52:01:10, A*23:01:18, A*25:01:08
47	395 bp	1070 bp	14, 16, 38, 39, 3901, 64, 65	*07:69, 07:85, 07:180, 08:65, 14:01:01:01-14:20, 14:22-14:41N, 14:43-14:48, 14:50-14:56, 38:01:01-38:02:07, 38:04-38:05, 38:08-38:25, 38:27-38:53, 38:55Q-38:68Q, 39:01:01:01-39:01:01:07, 39:01:03:01-39:01:24, 39:03:01:01-39:07, 39:09:01:01-39:09:03, 39:11-39:12, 39:14-39:15, 39:18-39:19:02, 39:22, 39:24:01-39:38Q, 39:40:01N-39:48, 39:50-39:57, 39:59-39:62, 39:64-39:82, 39:84-39:87N, 39:90-39:95N, 39:97N-39:98, 39:100, 39:102-39:104, 39:106-39:109, 39:111-39:125, 39:127-39:128, 40:322 [?] , 49:26
48	160 bp, 425 bp	1070 bp	7, 703, 8, 14, 15, 16, 21, 39, 3901, 3902, 40, 41, 42, 45, 48, 50, 60, 61, 64, 65, 67, 70, 73, 81	*07:02:01:01-07:26, 07:28-07:35, 07:37, 07:39-07:64, 07:66-07:80, 07:82-07:148, 07:150-07:178, 07:181N-07:191, 07:193-07:194, 07:196-07:218, 07:220-07:235, 07:237-07:256, 07:258-07:270, 07:272N, 07:274-07:296, 08:01:01:01-08:01:36, 08:04:01-08:05, 08:07-08:35, 08:37-08:51, 08:53:01-08:77, 08:79-08:116, 08:118-08:125, 08:127-08:130, 08:132-08:135, 08:137-08:138, 08:140-08:146, 08:148N-08:184, 14:01:01:01-14:01:04, 14:01:06-14:02:04, 14:02:06-14:02:07, 14:02:09-14:29, 14:31-14:57, 27:14, 27:81, 27:130, 39:01:01:01-39:01:01:07, 39:01:03:01-39:06:04, 39:07-39:20, 39:22-39:24:01, 39:25N-39:103, 39:105-39:128, 40:01:01-40:01:49, 40:06:01:01-40:06:04:02, 40:06:06-40:07, 40:09 [?] , 40:11:01 [?] -40:11:02 [?] , 40:14:01 [?] -40:15 [?] , 40:16, 40:18 [?] , 40:21 [?] -40:22N [?] , 40:23, 40:24 [?] -40:31 [?] , 40:32, 40:33 [?] -40:38 [?] , 40:42-40:44, 40:45 [?] -40:46 [?] , 40:48 [?] -40:49 [?] , 40:51, 40:52 [?] , 40:53, 40:54 [?] -40:58 [?] , 40:59-40:60, 40:61 [?] -40:63 [?] , 40:65 [?] , 40:66, 40:67 [?] -40:69 [?] , 40:70:01-40:70:02, 40:71 [?] , 40:72:01-40:73, 40:74 [?] , 40:75, 40:77, 40:78 [?] , 40:79, 40:80 [?] -40:82 [?] , 40:83, 40:84 [?] -40:85 [?] , 40:86, 40:87:01 [?] -40:88 [?] , 40:90 [?] -40:92 [?] , 40:93, 40:95 [?] , 40:96, 40:98 [?] , 40:100 [?] -40:102 [?] , 40:103, 40:104 [?] -40:108 [?] , 40:109-40:110, 40:111 [?] -40:112 [?] , 40:113-40:114:01, 40:114:02 [?] -40:116 [?] , 40:118N [?] -40:121 [?] , 40:123 [?] , 40:124:01, 40:124:02 [?] -40:126 [?] , 40:127, 40:128 [?] -40:130:01 [?] , 40:130:02-40:131, 40:132 [?] -40:134 [?] , 40:136 [?] -40:137 [?] , 40:138-40:140, 40:141 [?] , 40:145 [?] , 40:146-40:148, 40:149 [?] , 40:150, 40:151 [?] , 40:152-40:153, 40:154 [?] , 40:155:01N-40:155:02N, 40:156 [?] , 40:158 [?] , 40:159, 40:160:01 [?] -40:160:02 [?] , 40:161-40:162, 40:163 [?] -40:164 [?] , 40:165, 40:166 [?] , 40:167, 40:168 [?] -40:170 [?] , 40:171-40:172, 40:173 [?] -40:176 [?] , 40:177, 40:178 [?] -40:185 [?] , 40:186:01, 40:186:02 [?] -40:187 [?] , 40:189 [?] , 40:190, 40:191 [?] -40:200 [?] , 40:202 [?] -40:213 [?] , 40:215 [?] -40:228 [?] , 40:230, 40:231 [?] -40:236 [?] , 40:237, 40:238 [?] -40:243 [?] , 40:244, 40:245 [?] -40:248 [?] , 40:249-40:253, 40:256N, 40:257 [?] , 40:258-40:265N, 40:268-40:270, 40:272-40:273, 40:274 [?] , 40:275, 40:276 [?] , 40:277-40:278, 40:279 [?] -40:280 [?] , 40:282 [?] -40:283 [?] , 40:285-40:286N, 40:287 [?] -40:291N [?] , 40:293 [?] , 40:294-40:295, 40:296 [?] -40:297 [?] , 40:298, 40:299 [?] -40:300 [?] , 40:301, 40:306-40:307, 40:308 [?] , 40:310 [?] , 40:311, 40:312 [?] , 40:313, 40:315 [?] , 40:316-40:318, 40:319 [?] , 40:321, 40:322 [?] , 40:323, 40:324 [?] -40:329 [?] , 40:330, 40:331 [?] , 40:332, 40:333 [?] , 40:336 [?] -40:339 [?] , 40:340-40:342,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				40:343 [?] , 40:346, 40:348-40:350, 40:351 [?] -40:352 [?] , 40:354, 40:355 [?] , 40:357 [?] -40:360 [?] , 41:01:01-41:45N, 41:47-41:52, 42:01:01-42:02:01:02, 42:04-42:21, 42:23-42:24, 45:01:01-45:20, 48:01:01:01-48:17, 48:19-48:41, 50:01:01:01-50:02, 50:04-50:20, 50:31-50:53, 51:01:09, 51:01:21, 51:01:25, 51:10, 51:24:02-51:24:04, 51:197, 51:217, 52:54, 54:11, 56:20:02, 67:01:01-67:07, 73:01-73:02, 81:02, C*05:152, C*07:463
49¹¹	180 bp	1070 bp	15, 40, 44, 48, 60, 62, 71	*15:116, 15:124, 15:352, 15:412, 27:125, 35:230, 35:331, 38:63, 40:01:01-40:01:06, 40:01:08-40:01:17, 40:01:19-40:01:49, 40:07, 40:10:01:01-40:10:02, 40:12, 40:21-40:23, 40:25, 40:30, 40:33-40:34, 40:36, 40:38, 40:42-40:43, 40:46-40:49, 40:51-40:52, 40:54-40:55, 40:59-40:63, 40:65-40:67, 40:69, 40:73-40:74, 40:76, 40:79, 40:81, 40:84, 40:87:01-40:88, 40:92, 40:100-40:102, 40:106, 40:108, 40:112-40:114:02, 40:116-40:118N, 40:123-40:126, 40:128, 40:130:01-40:130:02, 40:132, 40:134-40:135, 40:137-40:141, 40:146-40:147, 40:149-40:156, 40:160:01-40:160:02, 40:163, 40:166, 40:168, 40:170-40:172, 40:175, 40:178-40:179, 40:182-40:183, 40:185-40:188, 40:191, 40:193-40:196, 40:198-40:199, 40:204, 40:207-40:208, 40:210, 40:212-40:213, 40:215-40:218, 40:221-40:223, 40:227-40:228, 40:231, 40:234-40:236, 40:238-40:242, 40:245, 40:247, 40:249-40:253, 40:257-40:265N, 40:272-40:273, 40:277-40:282, 40:285-40:286N, 40:288, 40:299-40:301, 40:308, 40:310, 40:312, 40:315, 40:319, 40:321, 40:323-40:329, 40:332-40:333, 40:336, 40:338N-40:339, 40:344, 40:346, 40:352-40:353, 40:355, 40:357-40:358, 40:360, 44:31, 44:169, 44:182, 46:06, 46:43, 48:03:01-48:03:02, 48:17, 48:23, 48:40, 57:59, C*03:129, C*03:234, C*05:136, C*17:16
50⁹	290 bp, 465 bp	1070 bp	13, 17, 21, 35, 37, 40, 4005, 41, 46, 47, 50, 57, 60, 61, 62, 63	*13:02:06, 13:26:01, 13:72, 15:16:03, 15:20, 15:85, 15:194, 15:393, 15:407, 18:48, 35:46, 35:207, 37:06:01-37:06:02, 40:02:01:01-40:06:15, 40:08, 40:09 [?] , 40:10:02 [?] -40:11:02 [?] , 40:13, 40:14:01 [?] -40:15 [?] , 40:18 [?] -40:19 [?] , 40:20, 40:22N [?] , 40:24 [?] -40:31 [?] , 40:33 [?] -40:38 [?] , 40:39, 40:44 [?] -40:45 [?] , 40:47 [?] -40:49 [?] , 40:50, 40:52 [?] -40:58 [?] , 40:61 [?] -40:63 [?] , 40:64, 40:65 [?] , 40:67 [?] -40:69 [?] , 40:70:01-40:70:02, 40:71 [?] , 40:72:02 [?] , 40:74 [?] -40:76 [?] , 40:78 [?] , 40:80 [?] -40:82 [?] , 40:84 [?] -40:88 [?] , 40:89, 40:90 [?] -40:92 [?] , 40:94, 40:95 [?] -40:96 [?] , 40:97, 40:98 [?] , 40:99, 40:100 [?] -40:112 [?] , 40:114:02 [?] -40:116 [?] , 40:118N [?] -40:121 [?] , 40:122, 40:123 [?] , 40:124:02 [?] -40:130:01 [?] , 40:131 [?] -40:136 [?] , 40:141 [?] , 40:142N-40:144N, 40:145 [?] , 40:148 [?] -40:149 [?] , 40:151 [?] , 40:154 [?] , 40:156 [?] -40:169 [?] , 40:173 [?] -40:176 [?] , 40:178 [?] -40:184 [?] , 40:186:02 [?] -40:208 [?] , 40:210 [?] -40:213 [?] , 40:214, 40:215 [?] -40:228 [?] , 40:229, 40:230 [?] -40:236 [?] , 40:238 [?] -40:248 [?] , 40:254-40:255, 40:256N [?] -40:257 [?] , 40:266-40:271, 40:273 [?] -40:276 [?] , 40:279 [?] -40:283 [?] , 40:284, 40:287 [?] -40:293 [?] , 40:295 [?] -40:300 [?] , 40:302-40:305, 40:306 [?] -40:308 [?] , 40:309, 40:310 [?] -40:313 [?] , 40:314, 40:315 [?] , 40:317 [?] -40:319 [?] , 40:320, 40:322 [?] , 40:324 [?] -40:331 [?] , 40:333 [?] , 40:334, 40:335 [?] -40:345N [?] , 40:347 [?] , 40:356 [?] , 40:357-40:360, 46:01:07, 46:33, 47:01:01:01-47:01:02, 47:02-47:03, 47:05-47:10, 57:05, 57:67:01, 58:36, C*01:73
51⁵	105 bp	1070 bp	7, 8, 22, 41, 42	*07:04, 07:25, 07:146, 08:01:01:01-08:01:24, 08:01:26-08:05, 08:07-08:12:03, 08:14-08:19N, 08:21-08:24, 08:26:01-08:30N, 08:32-08:39, 08:41-08:54:02, 08:56:01-08:69, 08:71-08:78, 08:80-08:102, 08:104-08:112, 08:114-08:115, 08:117-08:133, 08:135-08:142, 08:144-08:155, 08:157-08:168, 08:170, 08:172-

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				08:174, 08:176-08:184, 35:87, 40:136, 40:231, 40:279, 41:01:01-41:05, 41:07-41:08, 41:10-41:14, 41:16-41:28, 41:30-41:31, 41:33-41:52, 42:01:01-42:02:01:02, 42:04-42:07, 42:09-42:13, 42:15-42:18, 42:20-42:24, 44:166, 55:20
52	325 bp	800 bp	7, 22, 27, 2708, 35, 42, 45, 46, 54, 55, 56, 67, 73, 78, 81, 82	*07:02:01:01-07:02:52, 07:04-07:07, 07:09-07:15, 07:17-07:26, 07:28-07:31, 07:33:01-07:36, 07:39-07:46, 07:47 ^w , 07:48-07:49N, 07:51-07:68:03, 07:70-07:84, 07:86-07:179, 07:181N-07:206, 07:208-07:218, 07:220-07:227, 07:229-07:231N, 07:233-07:296, 08:123, 08:132, 14:21, 14:57, 15:76, 15:101, 15:255, 27:01-27:11, 27:13-27:15, 27:17, 27:19-27:21, 27:24-27:28, 27:30-27:38, 27:40-27:58, 27:60-27:76, 27:78-27:91, 27:93-27:118, 27:120-27:128, 27:130-27:136, 27:138-27:152, 27:154-27:164, 35:76, 35:186, 35:229, 37:60, 38:26, 40:166, 42:01:01-42:02:01:02, 42:04-42:06, 42:08-42:10, 42:12-42:23, 44:90, 44:97, 45:06, 46:01:01-46:01:22, 46:02 ^w , 46:03-46:40:02, 46:42-46:48, 46:50-46:72, 54:01:01-54:29, 54:31-54:38, 55:01:01-55:05, 55:07-55:17, 55:19-55:77, 55:79-55:85, 56:01:01:01-56:16, 56:18-56:22, 56:24-56:32, 56:34-56:41, 56:43, 56:45-56:53, 67:01:01-67:07, 73:01-73:02, 81:01-81:08, 82:01-82:03, 83:01
53^{5,7}	80 bp, 115 bp, 160 bp, 195 bp, 225 bp, 260 bp	1070 bp	12, 35, 37, 44, 45, 47, 48, 57	*08:49, 08:60, 08:76, 08:87, 08:129, 08:181, 27:127, 35:38, 35:45, 35:71, 35:115, 35:169, 35:260, 35:309, 37:01:01-37:01:07, 37:01:09, 37:01:11, 37:03N-37:06:02, 37:08, 37:10-37:38, 37:41-37:43, 37:45-37:50, 37:52-37:65, 37:67-37:68, 40:132, 40:178, 40:231, 40:279, 41:48, 42:13, 44:02:01:01-44:14, 44:16-44:17, 44:19N-44:36, 44:38-44:63, 44:65-44:134, 44:136-44:149N, 44:151-44:157, 44:159-44:165, 44:167-44:183, 44:185-44:212, 44:214-44:229, 44:231-44:235, 44:237N-44:245, 44:247-44:267N, 44:269-44:271, 48:11, 51:42, 53:22, 57:07, 57:09, 57:24, 83:01, C*15:25
54	215 bp	800 bp	5, 8, 12, 17, 21, 22, 41, 42, 44, 45, 49, 51, 5102, 5103, 52, 53, 56, 57, 62	*08:09, 08:84, 13:46, 15:83, 41:01:01-41:01:05, 41:05-41:07, 41:09, 41:12, 41:14, 41:16-41:17, 41:20-41:22, 41:25-41:26, 41:28-41:30, 41:32-41:35, 41:37, 41:44, 41:50, 42:04, 44:06, 44:15, 44:18, 44:20, 44:100, 45:01:01-45:20, 49:03, 51:01:01:01-51:01:01:21, 51:01:03-51:01:08, 51:01:10-51:01:39, 51:01:41-51:02:01:02, 51:02:03-51:04, 51:06:01-51:07:01, 51:08:01-51:09:01, 51:09:03-51:14, 51:16-51:24:05, 51:26-51:46, 51:48-51:53, 51:55-51:77, 51:79-51:88, 51:90-51:92:02, 51:94-51:102, 51:104-51:115, 51:117-51:118N, 51:120-51:131, 51:133-51:143, 51:145-51:147, 51:149N-51:175, 51:177-51:178N, 51:181-51:217, 51:219-51:225, 52:01:02:01-52:01:02:03, 52:01:04, 52:01:09, 52:02:01, 52:03, 52:06:01-52:06:02, 52:09, 52:10:03, 52:19, 52:21, 52:31:02, 53:06, 55:20, 55:56, 56:13, 57:01:20, 58:08:01, A*23:01:14
55	130 bp, 280 bp	800 bp	7, 17, 18, 22, 35, 37, 40, 46, 51, 53, 60, 62, 78	*07:12, 07:14, 07:18:01-07:18:02, 07:137, 07:162-07:163, 07:199, 07:225, 07:227, 07:248, 14:10, 15:57 ^w , 15:289, 15:297, 15:318, 18:22, 18:69, 18:105, 35:21, 35:24:01-35:24:02, 35:81, 35:96, 35:109, 35:157, 35:188, 35:190, 35:233, 35:287, 37:04:01-37:04:02, 40:28, 40:30, 40:34, 40:59, 40:137, 40:160:01-40:160:02, 44:150, 44:190, 44:203, 46:01:01-46:39, 46:41N-46:46, 46:48-46:71, 48:17, 51:04, 51:46, 51:56:01-51:56:03, 51:139, 53:02, 53:06, 53:28, 57:14:01-57:14:02, 58:09, 58:76, C*01:73, C*03:260, C*15:39, C*15:130
56^{5,12}	90 bp, 410 bp	1070 bp	7, 15, 22, 27, 2708, 35, 47, 54,	*07:65, 27:01 ^w , 27:02:01:01-27:11, 27:13-27:15, 27:17, 27:19-27:21, 27:24-27:28, 27:30-27:38, 27:40-27:58, 27:60-27:76, 27:78-27:91, 27:93-27:118, 27:120-27:128, 27:130-27:136, 27:138-27:152, 27:154-27:164, 35:76, 35:186, 35:229, 37:60,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				55, 56, 78, 81, 82	44:90, 44:97, 47:01:01:01-47:03, 47:06-47:10, 54:01:01-54:29, 54:31-54:38, 55:01:01-55:05, 55:07-55:17, 55:19-55:77, 55:79-55:85, 56:01:01:01-56:16, 56:18-56:22, 56:24-56:32, 56:34-56:41, 56:43, 56:45-56:53, 81:01, 82:01-82:03, 83:01, C*14:79
57⁵	90 bp, 175 bp	1070 bp	15, 27, 35, 48, 60, 62, 75, 77, 81	*07:174, 07:202, 07:222, 07:279, 13:36, 15:02:01-15:02:11, 15:08:01-15:08:02, 15:11:01-15:11:07, 15:13:01-15:13:02, 15:15, 15:21, 15:31, 15:44, 15:55, 15:76 ^w , 15:88-15:89, 15:112, 15:121, 15:139, 15:144, 15:148, 15:170, 15:189, 15:191, 15:194, 15:209N, 15:213-15:215, 15:223, 15:265, 15:283-15:284, 15:291, 15:297, 15:301-15:303, 15:305, 15:319, 15:325, 15:328, 15:330, 15:341, 15:345, 15:350, 15:357-15:358, 15:367, 15:374, 15:378, 15:384, 15:394, 15:399, 15:404, 15:418, 15:420, 15:425, 27:24, 35:46, 35:311, 40:31, 40:45, 40:80, 48:01:01:01-48:01:06, 48:04:01-48:04:02, 48:06-48:07, 48:09, 48:11, 48:15-48:16, 48:18-48:20, 48:22, 48:24, 48:27-48:39, 48:41, 56:39 ^w , 81:01-81:08	
58	145 bp, 430 bp	1070 bp	44, 49, 59, 61	*40:13, 40:19, 40:109, 40:117, 40:292, 40:340, 41:46, 44:18, 44:25, 44:50:01, 44:95, 49:01:01:01-49:01:10, 49:04:01-49:14, 49:16-49:25, 49:27-49:48, 51:112, 53:30, 53:39, 54:12, 56:21, 57:45, 57:51, 57:69, 59:01:01:01-59:07, 59:09-59:10N	
59⁵	120 bp, 210 bp	800 bp	5, 13, 15, 17, 22, 35, 45, 49, 50, 51, 5102, 5103, 52, 55, 56, 62, 63, 78	*07:78, 07:84, 13:16, 13:31 ^w , 13:48, 13:62, 14:37, 15:04:01:01 ^w -15:04:04 ^w , 15:16:01:01-15:16:03, 15:42, 15:67, 15:95, 15:137 ^w , 15:254, 15:293, 15:303 ^w , 15:310, 15:362, 15:408, 27:153, 35:265 ^w , 40:95, 40:148, 40:161, 46:32 ^w , 49:01:01:01-49:01:01:02, 49:01:03-49:10, 49:12-49:14, 49:16-49:44, 49:46, 49:48, 50:01:01:01-50:02, 50:04-50:11, 50:13-50:16, 50:18-50:20, 50:31-50:46, 50:47 ^w , 50:48, 50:50, 50:52-50:53, 51:01:01:01-51:03, 51:05, 51:07:01-51:09:03, 51:11N-51:14, 51:16-51:24:05, 51:26-51:41N, 51:43-51:44N, 51:48-51:55, 51:57-51:58, 51:60-51:61:02, 51:63:01-51:63:02, 51:65-51:80, 51:82-51:134, 51:136-51:138, 51:140-51:147, 51:149N-51:156, 51:158:01-51:193, 51:195-51:212, 51:214-51:219, 51:221-51:226, 52:01:01:01-52:71, 54:20, 54:33 ^w , 55:01:01-55:01:06, 55:01:08-55:01:16, 55:03, 55:05, 55:09, 55:11, 55:15, 55:17, 55:21 ^w , 55:24-55:25, 55:28-55:29, 55:31, 55:33, 55:36, 55:38, 55:40, 55:44-55:45, 55:52-55:55N, 55:58-55:60, 55:64, 55:66, 55:68, 55:73-55:76, 55:78-55:79, 55:84-55:85, 56:05:01-56:06, 56:21, 56:25, 56:37, 58:08:01-58:08:02, 78:01:01:01-78:09, C*05:152, C*12:183^w	
60	430 bp	1070 bp	5, 17, 27, 44, 51, 5102, 5103, 52, 53, 57, 58, 61	*18:67, 18:136, 27:02:01:01-27:02:03, 27:30, 27:53, 27:57, 27:62, 27:65N, 27:75, 27:77, 27:83, 27:95, 27:102, 27:119, 27:126, 27:134, 27:157, 27:163, 37:34, 40:13, 40:19, 40:109, 40:292, 40:340 [?] , 44:06, 44:25, 44:50:01, 44:95, 51:01:01:01-51:01:39, 51:01:41-51:24:05, 51:26-51:46, 51:48-51:53, 51:55-51:77, 51:79-51:226, 52:01:01:01-52:01:28, 52:02:01-52:19, 52:21-52:71, 53:01:01-53:02, 53:04-53:08:02, 53:10, 53:14-53:30, 53:32-53:35, 53:37, 53:39-53:48N, 57:01:01-57:11, 57:13-57:90, 58:01:01:01-58:02:02, 58:04-58:16:02, 58:18-58:29, 58:31N-58:63, 58:65-58:91	
61	145 bp	1070 bp	12, 21, 35, 40, 4005, 41, 45, 47, 50, 60, 61	*07:133, 08:155, 15:46, 15:53, 15:106, 15:143, 15:212, 18:48, 18:107, 35:19, 35:47, 35:63, 35:154, 35:217, 35:274, 35:317, 39:69, 40:01:01-40:11:02, 40:14:01-40:16, 40:18, 40:20, 40:22N-40:40, 40:42-40:45, 40:48-40:75, 40:77-40:92, 40:94-40:95, 40:97-40:108, 40:111-40:116, 40:118N-40:136, 40:138-40:156, 40:158-40:184, 40:186:01-40:187, 40:189-40:200, 40:202-40:207, 40:210-40:264, 40:266-40:291N, 40:293-	

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				40:307, 40:309-40:339, 40:341-40:349, 40:351-40:360, 41:01:01-41:04, 41:05 [?] , 41:06-41:25, 41:27-41:45N, 41:47-41:52, 44:09, 44:46, 44:75, 44:90, 44:129, 44:131, 44:254, 45:01:01-45:15, 45:17-45:20, 47:02, 47:03 ^w , 48:41, 50:01:01:01-50:02, 50:04-50:05, 50:07-50:20, 50:31-50:53
62⁸	300 bp	1070 bp	7, 12, 13, 15, 17, 18, 21, 27, 2708, 35, 37, 38, 39, 3902, 40, 4005, 41, 44, 45, 47, 48, 49, 50, 52, 60, 61, 62, 70, 72, 76	*07:54, 07:123, 08:17, 08:38, 08:54:01-08:54:02, 08:101, 13:01:01:01-13:04, 13:06, 13:08-13:23, 13:25-13:35, 13:37-13:103N, 15:01:01:01-15:01:04, 15:01:06-15:01:18, 15:01:20-15:01:40, 15:03:01:01-15:03:04, 15:03:06-15:07:03, 15:12, 15:14, 15:19-15:20, 15:24:01-15:27:04, 15:30, 15:32:01-15:36, 15:38:01-15:40:02, 15:42-15:43, 15:45-15:50, 15:53-15:54, 15:56-15:58, 15:60-15:63, 15:65-15:66, 15:68-15:71, 15:73-15:75, 15:77-15:79N, 15:81-15:87, 15:91-15:92, 15:94N, 15:96-15:98, 15:101-15:107, 15:109-15:111N, 15:113, 15:116-15:118, 15:122-15:123, 15:125-15:129, 15:131-15:132, 15:135-15:138, 15:140-15:142, 15:145-15:147, 15:150-15:152, 15:154-15:160, 15:163-15:167, 15:169, 15:171-15:175, 15:178-15:179:02, 15:181N-15:185, 15:187-15:188, 15:190N, 15:192-15:193, 15:195, 15:199, 15:201-15:207, 15:210-15:212, 15:217-15:220, 15:224-15:225, 15:227-15:228, 15:231-15:237, 15:239-15:245Q, 15:247-15:251, 15:253, 15:255-15:262N, 15:264, 15:266-15:267, 15:269-15:272N, 15:274, 15:276-15:282, 15:286-15:289, 15:295-15:296, 15:298-15:300, 15:304N, 15:309-15:310, 15:315-15:318, 15:320-15:322, 15:324, 15:326-15:327, 15:331-15:334, 15:336, 15:339-15:340, 15:342-15:344, 15:346-15:349:02, 15:352-15:353, 15:355, 15:359-15:360, 15:363:01-15:366, 15:368-15:373, 15:375N-15:377Q, 15:379, 15:381, 15:383, 15:385-15:387, 15:389-15:391, 15:393, 15:395, 15:397-15:398, 15:400N-15:402, 15:405-15:407, 15:409-15:410, 15:412-15:413, 15:415-15:417, 15:419, 15:421-15:422, 15:429, 18:12:01-18:12:02, 27:01-27:21, 27:24-27:91, 27:93-27:108, 27:110-27:156, 27:158-27:164, 35:10, 35:13, 35:16, 35:28, 35:63, 35:69, 35:80, 35:207, 35:226, 35:270, 35:317-35:318, 37:01:01-37:07, 37:09-37:66, 37:68, 38:03-38:04, 38:25, 39:02:01-39:02:03, 39:08, 39:13:01-39:13:02, 39:22-39:23, 39:49, 39:101, 39:105, 39:110, 40:01:01-40:07, 40:09-40:16, 40:18-40:24, 40:26-40:40, 40:42-40:67, 40:70:01-40:105, 40:107-40:165, 40:167-40:231, 40:233-40:312, 40:314-40:356, 40:358-40:360, 41:01:01-41:52, 42:11, 44:02:01:01-44:02:21, 44:02:23-44:03:35, 44:03:37-44:05:04, 44:07-44:11, 44:13-44:170, 44:172-44:231, 44:233-44:245, 44:247-44:262, 44:264-44:271, 45:01:01-45:20, 47:01:01:01-47:10, 48:01:01:01-48:05, 48:07-48:15, 48:17-48:38, 48:40-48:41, 49:01:01:01-49:14, 49:16-49:25, 49:27-49:48, 50:01:01:01-50:02, 50:04-50:14, 50:16-50:20, 50:31-50:53, 52:01:01:01-52:71, 53:17:01-53:17:02, 53:28, 53:38, 54:35, 55:18, 55:34, 55:40, 78:05-78:06
63^{7,13}	160 bp, 370 bp	1070 bp	22, 44, 45, 47, 54	*08:87, 27:127, 44:02:01:01-44:02:47, 44:06, 44:08-44:09, 44:11-44:12, 44:16, 44:19N-44:24, 44:27:01-44:27:02, 44:33-44:34:02, 44:41:01-44:41:02, 44:44, 44:48-44:49, 44:52N-44:53:02, 44:55, 44:59:01-44:59:02, 44:63, 44:66-44:68, 44:71-44:74, 44:80, 44:83-44:84:02, 44:86-44:87, 44:89-44:91, 44:93, 44:95, 44:97, 44:99-44:102, 44:104, 44:106, 44:112-44:113, 44:116, 44:118-44:119, 44:121, 44:126:01-44:127, 44:131-44:132, 44:137-44:138Q, 44:142, 44:145, 44:148-44:149N, 44:151-44:152, 44:162, 44:168-44:173, 44:176-44:177, 44:179, 44:185, 44:187, 44:190-44:191, 44:195N-44:196, 44:200-

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				44:201, 44:206, 44:208, 44:211-44:212, 44:214, 44:216, 44:218-44:221, 44:225, 44:229, 44:235, 44:238, 44:240-44:244, 44:249, 44:253-44:255, 44:257, 44:260-44:262, 44:265, 44:267N, 44:269-44:270, 54:01:01-54:38, 83:01
64	180 bp, 210 bp	1070 bp	13, 15, 18, 35, 39, 40, 54, 55, 56, 59, 78	*07:65, 07:206, 08:26:03, 13:06, 13:53, 14:37, 15:32:02, 15:42, 15:86, 15:224, 18:04:01, 18:125, 35:08:05, 35:42:02, 35:60, 35:185, 39:17, 39:63, 40:01:11, 40:58, 45:08, 46:18, 46:40:01, 50:45, 50:49, 51:157, 54:01:01-54:02, 54:05N, 54:07-54:08N, 54:10, 54:12-54:13, 54:16-54:25, 54:27-54:32, 54:34-54:38, 55:01:01-55:03, 55:05, 55:07, 55:10-55:12, 55:15-55:16, 55:18-55:19, 55:21, 55:25-55:26, 55:29-55:31, 55:33-55:45, 55:47-55:48, 55:50, 55:52-55:57, 55:59-55:67, 55:69-55:80, 55:82-55:85, 56:05:01, 56:06, 56:10, 56:23, 59:01:01:01-59:01:01:02, 59:05-59:08, 59:10N, 78:01:01:01-78:01:02, 78:02:02-78:03, 78:07-78:09, C*01:90, C*07:536, C*15:02:04, C*15:91, C*15:128
65	180 bp, 205 bp	1070 bp	7, 12, 15, 21, 22, 44, 45, 49, 50, 56, 61, 62, 82	*07:29, 07:186, 08:123, 08:132, 13:03, 13:48, 15:73, 15:303, 40:71, 44:10, 44:15, 44:18, 44:140, 45:01:01-45:01:02, 45:04-45:07, 45:11-45:20, 46:11, 49:01:01:01-49:01:07, 49:01:09-49:03, 49:06, 49:08-49:14, 49:16-49:17, 49:19N-49:48, 50:01:01:01-50:02, 50:04-50:08, 50:10-50:13, 50:15-50:19, 50:32, 50:34-50:38, 50:40-50:44, 50:46-50:48, 50:50-50:53, 51:15, 51:62, 51:106, 51:223, 51:225, 52:25:01-52:25:02, 54:03, 54:33, 56:01:01:01-56:01:04, 56:01:06-56:02:02, 56:04, 56:07-56:08, 56:13-56:14, 56:16-56:17, 56:20:01-56:20:02, 56:24-56:30, 56:33-56:49, 56:51-56:53, 58:73, 59:04, 82:01-82:03, C*03:19, C*03:102
66⁵	90 bp, 240 bp	800 bp	57	*08:181, 55:14, 57:01:01-57:01:22, 57:01:24-57:04:01, 57:05-57:44, 57:46-57:50, 57:52-57:68, 57:70-57:88, 57:90, 58:36, C*06:72
67^{5,12}	90 bp	1070 bp	5, 13, 17, 39, 55, 56, 58, 60, 61, 62, 71	*07:182N, 13:11, 15:18:03, 15:73, 15:224, 39:17, 40:48, 40:71, 51:62, 56:02:01-56:02:02, 56:04, 56:10, 56:51, 58:01:01:01-58:02:02, 58:04-58:29, 58:31N-58:35, 58:37-58:43, 58:45:01-58:91, C*02:06:01^w, C*02:47^w
68⁵	95 bp	1070 bp	8, 18, 22, 35, 37, 38, 39, 44, 51, 5102, 5103, 53, 78	*07:65 ^w , 07:134 ^w , 08:32, 08:133, 08:180, 15:308, 18:01:01:01-18:11, 18:13-18:15, 18:17N-18:36, 18:38-18:47, 18:49-18:65, 18:67-18:106, 18:108-18:124, 18:126-18:139, 35:01:01:01-35:08:04, 35:08:06-35:09:03, 35:11:01-35:12:04, 35:14:01-35:15:02, 35:17:01-35:18, 35:20:01-35:24:02, 35:27, 35:29:01-35:45, 35:48, 35:50-35:62, 35:64:01-35:68:02, 35:70-35:72, 35:74-35:75, 35:76 ^w , 35:77-35:79, 35:81-35:153, 35:155-35:185, 35:187-35:190, 35:192-35:197, 35:199-35:206, 35:208-35:209, 35:211-35:216N, 35:218-35:225, 35:227-35:253, 35:255-35:269, 35:272-35:273, 35:275-35:310, 35:312-35:316, 35:319-35:349, 37:08, 37:67, 38:06-38:07, 39:19:01-39:19:02, 44:06, 51:01:01:01-51:24:05, 51:26-51:46, 51:48-51:103, 51:105-51:111, 51:113-51:117, 51:119-51:146, 51:148-51:226, 53:01:01-53:16, 53:18-53:27, 53:29, 53:31-53:37, 53:40-53:48N, 56:06 ^w , 56:45 ^w , 78:01:01:01-78:04, 78:07-78:09
69^{5,9}	115 bp, 150 bp	1070 bp	15, 18,22, 27, 35, 39, 46, 55, 62, 72, 73, 75, 76	*07:100, 13:31, 13:41, 15:01:01:01-15:01:04, 15:01:06-15:01:16, 15:01:18-15:01:30, 15:01:32-15:01:40, 15:04:01:01-15:04:04, 15:07:01:01-15:08:02, 15:11:01-15:12, 15:14-15:15, 15:19, 15:24:01-15:24:02, 15:26N-15:28, 15:30, 15:32:01-15:32:02, 15:34-15:35, 15:38:01-15:38:02, 15:43, 15:45-15:46, 15:50, 15:53-15:54, 15:56-15:58, 15:60, 15:63, 15:66, 15:68,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				15:70-15:71, 15:73, 15:75-15:77, 15:79N, 15:81-15:82, 15:85, 15:87, 15:92, 15:94N, 15:96-15:97, 15:101-15:102, 15:104-15:105, 15:109-15:111N, 15:113, 15:117-15:118, 15:120, 15:122, 15:125-15:126, 15:128-15:129, 15:135, 15:137, 15:140, 15:142-15:149N, 15:152, 15:154, 15:157, 15:159-15:160, 15:163-15:167, 15:169, 15:171-15:172, 15:174-15:175, 15:178, 15:180-15:184, 15:187, 15:189-15:193, 15:201-15:203, 15:205-15:207, 15:209N, 15:211-15:212, 15:215, 15:217, 15:225, 15:227-15:228, 15:231-15:234, 15:236, 15:239, 15:241, 15:244-15:247, 15:251, 15:256-15:262N, 15:264, 15:267, 15:269-15:270, 15:272N, 15:276-15:280, 15:284-15:286, 15:296, 15:298-15:299, 15:303-15:305, 15:309, 15:315-15:318, 15:320-15:321Q, 15:324, 15:326-15:327, 15:331, 15:333-15:334, 15:336, 15:340, 15:342-15:344, 15:346-15:347, 15:349:01-15:350, 15:353, 15:355, 15:359, 15:363:01-15:363:02, 15:365, 15:367-15:368, 15:370-15:373, 15:375N, 15:377Q, 15:381, 15:385-15:387, 15:390-15:392, 15:398, 15:400N-15:401, 15:405-15:406, 15:409-15:410, 15:412-15:413, 15:415-15:417, 15:419, 15:422, 15:428-15:429, 18:19, 27:25, 27:75, 27:115, 27:161, 35:14:01-35:14:02, 35:43:01-35:44, 35:62, 35:67, 35:79, 35:86, 35:102, 35:117, 35:135, 35:163, 35:185, 35:213, 35:265, 35:343, 37:60, 37:66, 39:18, 39:36, 44:146, 46:01:01-46:02, 46:04-46:05, 46:07N, 46:09-46:10, 46:12, 46:14-46:17, 46:20, 46:22-46:24, 46:26-46:32, 46:34-46:42, 46:44, 46:46-46:53, 46:55N-46:57, 46:60-46:61, 46:63-46:72, 50:47, 51:61:01-51:61:02, 51:220, 52:21, 54:06, 54:33, 55:21, 56:03, 56:43, 58:78, 73:01-73:02, A*26:68, A*26:100, A*68:56, C*02:56, C*06:20, C*12:50
70 ^{8,14}	360 bp	1070 bp	B ^w 4	*07:27, 07:36, 07:38, 07:81, 07:149, 07:180, 07:219, 07:236, 07:273, 08:02-08:03, 08:52, 08:78, 08:117, 08:126, 13:01:01:01-13:04, 13:06-13:08, 13:10-13:23, 13:25-13:38, 13:40-13:66, 13:68-13:103N, 15:13:01-15:13:02, 15:16:01:01-15:17:06, 15:23-15:24:02, 15:36, 15:43, 15:67, 15:87, 15:89, 15:95, 15:115, 15:157, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:222, 15:230, 15:254, 15:256, 15:268, 15:273, 15:339, 15:356, 15:361-15:362, 15:396, 15:403, 15:408, 15:411, 15:418, 15:423-15:424, 18:09, 18:54, 18:67, 18:136, 27:01-27:07:06, 27:09-27:11, 27:13-27:17, 27:19-27:21, 27:23-27:25, 27:27-27:32, 27:34-27:39, 27:41, 27:43, 27:45-27:88, 27:90:01-27:152, 27:154-27:164, 35:36, 35:329, 37:01:01-37:04:02, 37:06:01-37:10, 37:12-37:36, 37:38-37:66, 37:68, 38:01:01-38:68Q, 40:13, 40:19, 40:47, 40:96, 40:109-40:110, 40:117, 40:157, 40:188, 40:201, 40:207, 40:246, 40:292, 40:340, 41:46, 44:02:01:01-44:08, 44:10-44:45, 44:47-44:74, 44:76-44:89, 44:91-44:128:02, 44:130, 44:132-44:245, 44:247-44:253, 44:255-44:271, 47:01:01:01-47:01:03, 47:03-47:10, 48:18, 49:01:01:01-49:14, 49:16-49:20, 49:22-49:48, 51:01:01:01-51:24:05, 51:26-51:46, 51:48-51:226, 52:01:01:01-52:71, 53:01:01-53:48N, 54:12, 56:07, 56:21, 56:46, 57:01:01-57:11, 57:13-57:90, 58:01:01:01-58:02:02, 58:04-58:29, 58:31N-58:63, 58:65-58:91, 59:01:01:01-59:10N
71	350 bp	1070 bp	B ^w 6	*07:02:01:01-07:26, 07:28-07:35, 07:37, 07:39-07:80, 07:82-07:148, 07:150-07:179, 07:181N-07:191, 07:193-07:218, 07:220-07:235, 07:237-07:272N, 07:274-07:296, 08:01:01:01-08:01:36, 08:04:01-08:05, 08:07-08:51, 08:53:01-08:77, 08:79-08:116, 08:118-08:125, 08:127-08:130, 08:132-08:135, 08:137-08:138, 08:140-08:152, 08:154-08:156, 08:158-08:184, 13:09,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				13:39, 13:67, 14:01:01:01-14:57, 15:01:01:01-15:01:04, 15:01:06-15:12, 15:14-15:15, 15:18:01:01-15:21, 15:25:01-15:35, 15:37-15:40:02, 15:42, 15:44-15:58, 15:60-15:66, 15:68-15:86, 15:88, 15:90-15:94N, 15:96-15:99, 15:101-15:114, 15:116-15:129, 15:131-15:156, 15:158-15:161, 15:163-15:167, 15:169-15:176, 15:178-15:195, 15:197:01-15:207, 15:209N-15:215, 15:217-15:221, 15:223-15:229, 15:231-15:253, 15:255, 15:257-15:267, 15:269-15:272N, 15:274-15:276, 15:278-15:338, 15:340-15:355, 15:357-15:360, 15:363:01-15:376, 15:378-15:395, 15:397-15:402, 15:404-15:407, 15:409-15:410, 15:412-15:417, 15:419-15:422, 15:425-15:429, 18:01:01:01-18:08, 18:10-18:15, 18:17N-18:53, 18:55-18:66, 18:68-18:135, 18:137-18:139, 27:08, 27:12, 27:18, 27:26, 27:33, 27:40, 27:42, 27:44, 27:89, 27:153, 35:01:01:01-35:01:38, 35:01:40-35:72, 35:74-35:221, 35:223-35:247, 35:249-35:292, 35:294-35:328, 35:330-35:349, 37:05, 37:11, 37:14 ^w , 37:37, 37:67, 39:01:01:01-39:01:01:07, 39:01:03:01-39:01:14, 39:01:16-39:20, 39:22-39:103, 39:105-39:128, 40:01:01-40:12, 40:14:01-40:16, 40:18, 40:20-40:40, 40:42-40:46, 40:48-40:75, 40:77-40:95, 40:97-40:108, 40:111-40:116, 40:118N-40:156, 40:158-40:187, 40:189-40:200, 40:202-40:291N, 40:293-40:339, 40:341-40:360, 41:01:01-41:25, 41:27-41:45N, 41:47-41:52, 42:01:01-42:02:01:02, 42:04-42:21, 42:23-42:24, 44:09, 44:46, 44:75, 44:90, 44:129, 44:131, 44:254, 45:01:01-45:15, 45:17-45:20, 46:01:01-46:72, 47:02, 47:03 ^w , 48:01:01:01-48:17, 48:19-48:41, 50:01:01:01-50:02, 50:04-50:20, 50:31-50:53, 54:01:01-54:11, 54:13-54:38, 55:01:01-55:05, 55:07-55:85, 56:01:01:01-56:06, 56:08-56:20:02, 56:22-56:45, 56:47-56:53, 57:12, 58:64, 67:01:01-67:07, 73:01-73:02, 78:01:01:01-78:09, 81:01-81:08, 82:01-82:03, 83:01
72	285 bp	800 bp	7, 703, 8, 35, 40, 41, 42, 48, 53, 60, 61, 81	*07:02:01:01-07:02:38, 07:02:40-07:06:02, 07:08:01-07:18:02, 07:20-07:32, 07:34-07:39, 07:41-07:52, 07:54-07:59, 07:61-07:76, 07:79-07:83, 07:85-07:99, 07:101-07:121, 07:123-07:138, 07:140-07:164, 07:166-07:192, 07:194-07:195, 07:199-07:213, 07:215-07:234, 07:236-07:238, 07:240-07:254, 07:256-07:296, 08:01:01:01-08:01:20, 08:01:22-08:01:24, 08:01:26-08:05, 08:07-08:08N, 08:10-08:11, 08:13-08:15, 08:17-08:30N, 08:32-08:83, 08:85-08:88, 08:90-08:102, 08:104-08:106, 08:108-08:112, 08:114-08:115, 08:117-08:133, 08:135-08:168, 08:170-08:174, 08:176-08:183, 15:297, 27:164, 35:66, 35:87, 37:07, 40:15-40:16, 40:30-40:32, 40:34, 40:45, 40:59, 40:80, 40:98, 40:137, 40:160:01-40:160:02, 41:02:01-41:02:07, 41:04, 41:10-41:11, 41:13, 41:18-41:19, 41:23-41:24, 41:27, 41:31, 41:36, 41:38-41:43, 41:45N-41:49, 41:51-41:52, 42:01:01-42:02:01:02, 42:05:01-42:07, 42:09-42:13, 42:15-42:24, 44:150, 44:166, 48:01:01:01-48:01:06, 48:05-48:12, 48:14-48:20, 48:22, 48:27-48:39, 48:41, 53:15, 81:01-81:08, C*03:260

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-B low resolution SSP typings. When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length. Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

primers. PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings. PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings. Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 25 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³The serological reactivity of all HLA-B alleles is not known. In this table we use the expert-assigned serological grouping in Tissue Antigens (2009) 73:95-170 and the serological grouping of the sequence-defined allele.

⁴For several HLA Class I alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁵HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁶Primer mix 28 may give rise to a lower yield of HLA-specific PCR product than the other HLA-B low resolution primer mixes in B*40, B*41, B*45, B*49 and B*50 alleles.

⁷Primer mixes 28, 29, 33, 42, 53 and 63 may give rise to a lower yield of HLA-specific PCR product than the other HLA-B low resolution primer mixes.

⁸Primer mixes 32, 33, 38, 42, 62 and 70 have tendencies of unspecific amplifications, most pronounced in primer mixes 32 and 33.

⁹Primer mixes 25, 43, 44, 50, and 69 have a tendency giving rise to primer oligomer formation.

¹⁰The B*57 and B*58 alleles might be faintly amplified by primer mix 33.

¹¹The C*17:01 to C*17:04 alleles might be faintly amplified by primer mix 49.

¹²Primer mixes 56 and 67 may generate a false positive band of about 800 base pairs. This band should be disregarded when interpreting HLA-B low resolution typings.

¹³Primer mix 63 may give rise to a lower yield of B*54 alleles than the other B low primer mixes.

¹⁴The Bw4-associated HLA-A specificities A9, A23, A24, A2403, A25 and A32 are not amplified by the primer pair in primer mix 70.

‘w’, might be weakly amplified.

‘?’, nucleotide sequence information not available for the primer matching sequence.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information
SPECIFICITY TABLE

Specificities and sizes of the PCR products of the 23+1 primer mixes used for DR low resolution SSP typing

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	DR serology ³	Amplified HLA-DRB alleles ^{4,5}
73⁷	210 bp, 235 bp, 260 bp	515 bp	1	*01:01:01-01:02:12, 01:04-01:38, 01:40N-01:83
74	200 bp	430 bp	103	*01:03:01-01:03:02, 01:39N, 01:42, 01:61
75^{7,9}	200 bp, 230 bp	430 bp	2, 15	*15:01:01-01-15:147
76	210 bp	430 bp	16	*16:01:01-16:05:02, 16:07-16:47
77^{5,6}	125 bp, 225 bp	515 bp	3, 11, 13, 14, 17, 18	*03:01:01-01-03:125, 03:127-03:140, 11:07:01-11:07:02, 11:27:02-11:27:03, 11:84:02-11:84:03, 11:103:01- 11:103:02, 11:105, 11:107, 11:125, 11:136, 11:173, 13:33:01, 13:61:02, 13:94:01, 13:96:01, 14:38:02, 14:104, 14:111, 14:171, 15:25, 15:37:01, 15:100, 15:104:01
78^{5,6}	75 bp, 205 bp	430 bp	3, 6, 11, 13, 14, 17	*03:01:01-01-03:01:23, 03:01:25-03:01:26, 03:04:01-03:06, 03:08-03:16, 03:18-03:20, 03:22-03:23, 03:25-03:26, 03:28, 03:30-03:31, 03:33-03:34, 03:36-03:37, 03:43- 03:48, 03:50-03:52, 03:54-03:68N, 03:70-03:73, 03:75- 03:86, 03:89, 03:91-03:93, 03:95-03:96, 03:98-03:100:02, 03:104, 03:106-03:110, 03:112-03:114, 03:116-03:118, 03:121-03:130, 03:132-03:140, 08:40, 11:02:01-11:03:02, 11:11:01, 11:11:03, 11:14:01-11:14:02, 11:16, 11:20-11:21, 11:36, 11:40-11:41, 11:48, 11:59, 11:63:01-11:63:02, 11:65:01-11:65:02, 11:68, 11:70, 11:73, 11:76, 11:79- 11:80, 11:83, 11:85-11:87, 11:93, 11:118, 11:122, 11:124, 11:127, 11:131-11:132, 11:135, 11:138-11:139, 11:142, 11:151, 11:153, 11:161, 11:168, 11:171, 11:176, 11:182, 11:184, 11:213, 13:01:01-01-13:04, 13:08, 13:10, 13:15- 13:17, 13:19-13:20, 13:22-13:24, 13:27-13:29, 13:31- 13:41, 13:43, 13:45, 13:48, 13:51-13:54, 13:57, 13:59, 13:61:01-13:61:02, 13:63-13:66:02, 13:68-13:76, 13:78- 13:81, 13:83-13:85, 13:87-13:99, 13:101-13:102, 13:104- 13:107, 13:109, 13:111-13:117, 13:120-13:131, 13:133, 13:135, 13:137N-13:145, 13:147-13:149, 13:151-13:153, 13:155, 13:159, 13:162, 13:165-13:168, 13:170-13:180, 13:182, 13:184-13:188, 13:190-13:191, 13:193-13:194, 13:196, 13:198, 13:200N-13:202, 13:204-13:205, 13:207- 13:209, 13:212-13:219, 13:221, 13:223, 13:225-13:227, 13:230, 13:233-13:234, 13:236-13:240, 14:16, 14:19, 14:21, 14:45, 14:82, 14:95, 14:109, 14:120, 14:122, 14:132, 14:178-14:179
79⁶	85 bp, 210 bp	430 bp	3, 6, 11, 13, 14, 18	*03:02:01-03:03, 03:27, 03:29, 03:38, 03:53, 03:74, 03:88, 03:90, 03:102-03:103, 03:115, 03:119, 03:126, 03:131, 11:13:01 ^w -11:13:02 ^w , 11:26, 11:34, 13:15, 13:19, 13:26:01- 13:26:02, 13:44, 13:53, 13:57, 13:85-13:86, 13:104, 13:193, 13:198, 13:206, 13:235, 14:02:01-14:03:02, 14:06:01-14:06:03, 14:09, 14:12:01-14:13, 14:17-14:21, 14:24, 14:27:01-14:27:02, 14:29-14:30, 14:32:01 ^w - 14:32:02 ^w , 14:33, 14:40-14:41, 14:47-14:49, 14:51, 14:63, 14:65 ^w , 14:67, 14:77-14:78, 14:80-14:81, 14:83, 14:85,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				14:89, 14:94, 14:98, 14:102, 14:106, 14:108-14:109, 14:115, 14:119, 14:121, 14:135, 14:146, 14:154, 14:156, 14:159-14:160, 14:165, 14:170, 14:174, 14:176-14:177, 14:179, 14:181, 14:183
80^{6,7,9}	100 bp, 175 bp	430 bp	3, 4	*04:01:01:01-04:05:11, 04:05:13-04:233
81⁹	170 bp, 210 bp, 235 bp	430 bp	7, 13, 14	*07:01:01:01-07:01:21, 07:03-07:80, 08:67, 10:11, 12:22, 13:17, 13:116, 13:175, 13:214, 14:50
82	170 bp, 215 bp, 250 bp	515 bp	8, 11, 12, 14	*08:01:01-08:01:02, 08:01:04-08:19, 08:21-08:56, 08:58-08:66, 08:68-08:79, 08:81-08:82, 11:67, 11:193:01-11:193:02, 11:209, 12:04, 12:16:01-12:16:03, 12:22, 12:39, 12:49, 14:11, 14:15, 14:68:01-14:68:02, 14:93, 14:148
83⁶	85 bp, 100 bp, 135 bp, 165 bp, 190 bp	430 bp	3, 9, 11	*03:08, 03:65, 03:140, 09:01:02:01-09:31, 11:07:01-11:07:02, 11:103:01-11:103:02, 11:105, 11:107, 11:125, 11:173, DRB3*01:42
84	175 bp	430 bp	9, 10,11,13	*03:76, 03:125-03:126, 09:05, 10:01:01:01-10:22, 11:59, 11:80, 11:83, 11:87, 11:135, 11:142, 11:182, 13:27, 13:41, 13:71, 13:129, 13:176, 14:82, 14:95, 14:132, 14:179
85⁶	100 bp, 175 bp	430 bp	3, 8, 11, 13, 14	*03:08, 03:65, 03:140, 08:31, 08:41, 08:63, 08:75, 11:01:01:01-11:11:01, 11:11:03-11:70, 11:72-11:217N
86⁶	90 bp, 110 bp	430 bp	12, 14	*08:32, 08:53, 12:01:01:01-12:56, 12:58-12:67, 14:52
87⁷	220 bp	430 bp	6, 8, 11, 13, 14, 1403	*03:76, 03:125, 08:20-08:21, 11:01:01:01-11:01:27, 11:01:28 [?] , 11:01:29-11:04:15, 11:06:01-11:06:03, 11:08:01-11:11:01, 11:11:03-11:12:03, 11:14:01-11:16, 11:18-11:21, 11:23:01-11:25, 11:27:01-11:33, 11:35-11:51, 11:54:01-11:54:02, 11:56-11:66, 11:68, 11:70, 11:72-11:81, 11:83-11:88, 11:90-11:97, 11:99-11:102:02, 11:106, 11:108-11:124, 11:126-11:135, 11:137-11:142, 11:144-11:158, 11:160-11:172, 11:174-11:192, 11:194-11:201, 11:203-11:208, 11:210-11:217N, 12:57, 13:01:01:01-13:02:01:03, 13:02:03-13:08, 13:10-13:16, 13:18-13:43, 13:45-13:85, 13:87-13:115, 13:117-13:128, 13:130-13:145, 13:147-13:174, 13:176-13:182, 13:184-13:191, 13:194-13:205, 13:207-13:209, 13:212-13:213, 13:215-13:221, 13:223-13:240, 14:03:01-14:03:02, 14:12:01-14:12:02, 14:16, 14:19, 14:21-14:22, 14:25, 14:27:01-14:27:02, 14:40, 14:53, 14:63, 14:67, 14:69, 14:74, 14:77-14:78, 14:84-14:85, 14:98, 14:102, 14:105, 14:109, 14:115-14:116, 14:128, 14:135, 14:144, 14:156, 14:173-14:174, 14:181, DRB3*02:27
88⁷	195 bp, 225 bp	430 bp	6, 8, 11, 12, 13, 14	*03:125, 08:01:01-08:01:02, 08:01:04-08:01:05, 08:02:01:01-08:02:04, 08:04:01-08:09, 08:11, 08:16-08:17, 08:20-08:22, 08:24:01-08:24:02, 08:26, 08:28, 08:31, 08:39, 08:41-08:44, 08:50, 08:52, 08:54-08:55, 08:57, 08:59, 08:64, 08:67, 08:70, 08:72-08:73, 08:75, 08:77-08:80, 08:82, 11:01:01:01-11:01:17, 11:01:20-11:06:03, 11:09-11:11:01, 11:11:03-11:12:03, 11:14:01-11:16, 11:20-11:21, 11:23:01-11:25, 11:27:01-11:30, 11:32-11:33, 11:35-11:41, 11:43-11:44, 11:46:01-11:51, 11:54:01-11:56, 11:58:01-11:63:02, 11:65:01-11:70, 11:72, 11:74:01-11:78, 11:80-11:88, 11:90-11:97, 11:99-11:102:02, 11:106,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				11:108-11:118, 11:120-11:124, 11:126-11:129, 11:133-11:135, 11:137-11:142, 11:144-11:152, 11:154-11:158, 11:160-11:172, 11:174-11:183, 11:185-11:190, 11:192-11:201, 11:203-11:217N, 12:02:01:01-12:02:07, 12:13, 12:15-12:16:03, 12:18-12:21, 12:23, 12:26-12:27, 12:31N-12:33, 12:37, 12:42-12:45, 12:49-12:52, 12:55-12:57, 12:60N-12:61, 12:64-12:65, 12:67, 13:01:01:01-13:02:01:03, 13:02:03-13:02:13, 13:04-13:05:03, 13:07:01-13:09, 13:11:01-13:11:02, 13:14:01-13:24, 13:26:01-13:29, 13:31-13:32, 13:34-13:36, 13:38-13:43, 13:45-13:55, 13:57, 13:59, 13:61:01-13:65, 13:67-13:76, 13:78-13:80, 13:83-13:84, 13:87, 13:91-13:93, 13:96:01-13:100, 13:102-13:109, 13:111-13:114, 13:116-13:117, 13:121, 13:123-13:132, 13:135-13:136, 13:138-13:150, 13:153, 13:155, 13:158-13:160, 13:162, 13:164-13:166, 13:168-13:169, 13:171:01-13:171:02, 13:173, 13:175, 13:177, 13:179, 13:182-13:187, 13:189-13:192, 13:195, 13:197-13:205, 13:207-13:215, 13:217-13:218, 13:221-13:223, 13:225-13:226, 13:229-13:230, 13:232-13:234, 13:236-13:240, 14:15-14:16, 14:22, 14:24-14:25, 14:27:01-14:27:02, 14:37, 14:53, 14:73, 14:105, 14:128
89⁵	175 bp, 240 bp	430 bp	3, 4, 6, 8, 11, 13, 14, 1403,1404, 17, 18	*03:01:01:01-03:01:05, 03:01:07-03:01:08, 03:01:10-03:07, 03:09, 03:11:01-03:41:02, 03:43-03:45, 03:47-03:63, 03:66-03:68N, 03:70-03:86, 03:88-03:91, 03:93-03:110, 03:112-03:129, 03:131-03:139, 08:20, 11:13:01-11:13:02, 11:149, 11:202, 12:20, 13:01:01:01-13:16, 13:18-13:42, 13:44, 13:46-13:66:02, 13:68-13:102, 13:104-13:115, 13:117-13:121, 13:123-13:158, 13:161-13:164, 13:166-13:170, 13:171:02-13:174, 13:176-13:178, 13:180-13:190, 13:192-13:194, 13:196, 13:198-13:213, 13:215-13:216, 13:218-13:227, 13:229-13:240, 14:01:01, 14:01:02 [?] -14:01:04 [?] , 14:02:01-14:04:01, 14:04:02 [?] -14:04:03 [?] , 14:05:01:01-14:07:01, 14:07:02 [?] , 14:08-14:14, 14:15 [?] -14:16 [?] , 14:17-14:21, 14:22 [?] , 14:23:01, 14:23:02 [?] , 14:23:03-14:24, 14:25 [?] -14:26 [?] , 14:27:01-14:27:02, 14:28 [?] , 14:29-14:30, 14:31 [?] -14:32:03 [?] , 14:33, 14:34 [?] -14:35 [?] , 14:36-14:37, 14:38:01 [?] -14:39 [?] , 14:40-14:45, 14:47-14:48, 14:49 [?] -14:50 [?] , 14:51, 14:52 [?] , 14:54:01:01, 14:54:01:02 [?] , 14:54:02-14:54:06, 14:55 [?] , 14:56-14:57, 14:58 [?] , 14:59, 14:60 [?] -14:62 [?] , 14:63-14:65, 14:67, 14:68:01 [?] -14:76 [?] , 14:77-14:78, 14:79 [?] , 14:80-14:85, 14:86 [?] -14:88 [?] , 14:89, 14:90 [?] , 14:91, 14:92N [?] -14:93 [?] , 14:94-14:96, 14:97 [?] , 14:98, 14:99 [?] , 14:100, 14:101 [?] , 14:102-14:103, 14:104 [?] -14:105 [?] , 14:106, 14:107 [?] , 14:108-14:109, 14:110 [?] -14:114 [?] , 14:115-14:116, 14:117 [?] -14:120 [?] , 14:121, 14:122 [?] , 14:123, 14:124 [?] -14:126:02 [?] , 14:127:01-14:127:02, 14:128 [?] -14:133 [?] , 14:134-14:137N, 14:138 [?] -14:140 [?] , 14:141, 14:142 [?] -14:143 [?] , 14:144, 14:145 [?] -14:151 [?] , 14:152N, 14:153 [?] , 14:154-14:156, 14:157 [?] -14:158 [?] , 14:159-14:161, 14:162 [?] -14:164 [?] , 14:165, 14:166N [?] , 14:167, 14:168 [?] , 14:169-14:171, 14:172 [?] -14:173 [?] , 14:174, 14:175 [?] , 14:176-14:177, 14:178 [?] , 14:179, 14:180 [?] , 14:181, 14:182 [?] , 14:183, 14:184 [?] -14:186 [?]
90^{6,7}	100 bp, 150 bp, 195 bp, 240 bp	430 bp	4, 6, 8,13, 14	*04:62, 04:69, 04:73, 04:105:01-04:105:02, 04:122, 04:146, 08:08, 11:69, 11:82, 13:45, 13:197, 14:01:01-14:01:02, 14:01:04, 14:04:01-14:04:03, 14:07:01-14:07:02, 14:10, 14:16, 14:22, 14:25-14:26, 14:28, 14:31-14:32:03,

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

				14:35, 14:37-14:39, 14:49-14:50, 14:53-14:54:01:02, 14:54:03-14:55, 14:57-14:58, 14:60-14:62, 14:68:01-14:71, 14:73-14:76, 14:79, 14:82, 14:86-14:88, 14:90, 14:93, 14:99, 14:101, 14:104-14:105, 14:107, 14:110-14:114, 14:117-14:120, 14:122, 14:124-14:125, 14:128-14:129, 14:131, 14:137N-14:140, 14:142-14:143, 14:145-14:147, 14:149-14:153, 14:157-14:158, 14:163-14:164, 14:166N, 14:168 ^w , 14:171-14:173, 14:178, 14:180, 14:182, 14:184, 14:186, DRB4*01:03:01:02N
91^{6,8}	110 bp, 145 bp, 170 bp	430 bp	3, 9, 11, 13, 14	*03:10, 03:126, 04:211, 09:01:02:01-09:01:05, 09:01:07-09:02:02, 09:04-09:22, 09:24-09:28, 09:30-09:31, 11:13:01-11:13:02, 11:17, 11:52, 11:202, 12:66, 13:43, 13:159, 13:171:01, 13:179, 13:191, 13:193, 13:217, 14:01:01-14:02:03, 14:04:01-14:11, 14:13-14:14, 14:16-14:18, 14:19 ^w , 14:20, 14:21 ^w , 14:22-14:23:04, 14:26, 14:28-14:36, 14:38:01-14:39, 14:41, 14:43-14:52, 14:54:01:01-14:54:03, 14:54:04 ^w , 14:54:05-14:57, 14:59-14:62, 14:64-14:65, 14:68:01-14:68:02, 14:70-14:76, 14:79-14:83, 14:86-14:88, 14:90-14:97, 14:99-14:101, 14:103-14:108, 14:109 ^w , 14:110-14:114, 14:117-14:127:02, 14:129-14:134, 14:137N-14:140, 14:142-14:143, 14:145-14:155, 14:157-14:160, 14:162-14:173, 14:175-14:180, 14:182, 14:184-14:185, 15:27, 15:34, 15:54, 15:66:01-15:66:02, 15:122
92^{6,7,9}	110 bp, 150 bp, 180 bp, 220 bp	430 bp	2, 3, 6, 8, 11, 13, 14 15, 16	*03:10, 08:09, 08:20-08:21, 08:32, 08:35, 08:36:02, 08:53, 08:68 ^w , 08:70 ^w , 11:13:01-11:13:02, 11:17, 11:23:01-11:23:02, 11:25, 11:31, 11:45, 11:52, 11:55, 11:64, 11:89, 11:96, 11:119, 11:148, 11:159, 11:202, 13:13, 13:18, 13:43, 13:45, 13:47, 13:55, 13:119, 13:144, 13:146, 13:154, 13:156, 13:158-13:159, 13:164, 13:171:01, 13:179, 13:191, 13:197, 13:217, 13:232, 14:01:01-14:01:04, 14:03:01-14:05:04, 14:07:01-14:08, 14:10-14:12:02, 14:14-14:16, 14:18, 14:22-14:23:04, 14:25-14:28, 14:31-14:32:03, 14:34-14:36, 14:38:01-14:40, 14:42-14:45, 14:49-14:50, 14:53-14:65, 14:67-14:79, 14:81-14:82, 14:84-14:93, 14:95-14:97, 14:99-14:105, 14:107, 14:110-14:120, 14:122-14:140, 14:142-14:158, 14:160-14:164, 14:166N-14:169, 14:171-14:175, 14:178, 14:180-14:182, 14:184-14:186, 15:21 ^w , 16:04:01 ^w -16:04:02 ^w , 16:18 ^w , DRB4*01:31
93^{5,6,7}	80 bp, 160 bp, 240 bp	430 bp	52	*14:141, DRB3*01:01:02:01-01:45, DRB3*02:01-02:65, DRB3*03:01:01-03:13
94⁷	130 bp, 215 bp	430 bp	53	*10:19, 15:141, DRB4*01:01:01:01-01:56N, DRB4*02:01N, DRB4*03:01N
95	175 bp	430 bp	51	DRB5*01:01:01-01:35, DRB5*02:02:01-02:17
96¹⁰	-	-	-	Negative control

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of DR low resolution SSP subtypings. When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length. Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers, e.g. the primers in wells 75 and 90 to 92. PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings. PCR

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: 9F4

Lot-specific information

fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings. Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 73 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases. In the presence of a specific amplification the intensity of the control band often decreases.

³The serological split of all DRB1 alleles is not known. In this table we use the expert-assigned serological grouping in Tissue Antigens (2009) 73:95-170 and the serological grouping of the sequence-defined allele.

⁴For several DRB1 alleles 1st and/or 3rd exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁵Due to sharing of sequence motifs in codon 38 and 47, DRB3*01:14 will also be amplified in primer mixes 77, 78 and 89, and DRB3*01:23 and 02:32 in mix 77, in addition to primer mix 93.

⁶HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁷Primer mix 73, 80, 87, 88, 90, and 92 to 94 may have a tendency giving rise to primer oligomer formation.

⁸Primer mix 91 has a tendency of primer oligomer formation and also has an intense primer cloud due to the high number of primers present in the primer mix.

⁹Primer mixes 75, 80, 81 and 92 may have tendencies of unspecific amplifications.

¹⁰Primer mix 96 contains a negative control, which will amplify more than 95% of HLA amplicons as well as the amplicons generated by the control primer pairs matching the human growth hormone gene. HLA-specific PCR product sizes range from 75 to 200 base pairs and the PCR product generated by the HGH positive control primer pair is 430 base pairs.

'w', might be weakly amplified.

'?', nucleotide sequence information not available for the primer matching sequence.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

HLA-A LOW PRIMER SPECIFICATION

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec.	120	175	205	190	160	135	175	165	80	80	80	125
PCR product	145	215	235		335	200	205	200		240	175	160
	225	255			505						500	190
		365										
		415										
Length of int.	800	800	1070	800	1070	800	1070	800	800	800	1070	800
pos. control ¹												
5'-primer(s) ²	98	48	357	98	144	176	98	98	257	98	301	103
	5'-CTT 3'	5'-gCT 3'	5'-ATg 3'	5'-CTA 3'	5'-gCC 3'	5'-gCA 3'	5'-CTC 3'	5'-CTA 3'	5'-Cgg 3'	5'-CTA 3'	5'-Cgg 3'	5'-CCT 3'
	103	78	363	413	317	368	368	102	259	261	302	415
	5'-CCT 3'	5'-TCT 3'	5'-ATA 3'	5'-CCg 3'	5'-gCT 3'	5'-gTT 3'	5'-gTT 3'	5'-ACA 3'	5'-AgA 3'	5'-AAC 3'	5'-ggA 3'	5'-ggT 3'
	123	106	363					413	261		385	423
	5'-AgT 3'	5'-CCA 3'	5'-ATA 3'					5'-CCg 3'	5'-AAC 3'		5'-ggC 3'	5'-gCT 3'
	363								266			
	5'-ATA 3'								5'-ACg 3'			
3'-primer(s) ³	203	240	527	256	265	270	259	259	302	299	341	257
	5'-TCT 3'	5'-ggA 3'	5'-CCA 3'	5'-CTg 3'	5'-CCC 3'	5'-ACA 3'	5'-gTT 3'	5'-gTT 3'	5'-ggC 3'	5'-TCg 3'	5'-CgT 3'	5'-gCA 3'
	545	292	527	559	368	521	502	259	303	308	521	506
	5'-AgA 3'	5'-gTg 3'	5'-CCT 3'	5'-CCg 3'	5'-CAA 3'	5'-ggg 3'	5'-CTT 3'	5'-gTT 3'	5'-AgA 3'	5'-TCT 3'	5'-ggg 3'	5'-TgT 3'
			527			534	539	538				538
			5'-CCT 3'			5'-CgT 3'	5'-TCT 3'	5'-CCA 3'				5'-CCg 3'
			555									559
			5'-CCA 3'									5'-CTC 3'
			555									559
			5'-gCA 3'									5'-CgT 3'
												559
												5'-CgC 3'
												563
												5'-CgA 3'
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

Well No.	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec.	180	80	90	240	140	200	340	210	65	85	80	360
PCR product	225	115	135	380	180	390	375	240	240	240	160	445
		200	205	410	240				375	400	240	
		240			260				545		495	
		460										
Length of int. pos. control ¹	1070	1070	1070	1070	1070	1070	800	800	800	800	800	1070
5'-primer(s) ²	98	98	203	41	98	41	302	78	28	47	176	257
	5'-CTT ^{3'}	5'-CAC ^{3'}	5'-gAA ^{3'}	5'-CTT ^{3'}	5'-CTT ^{3'}	5'-CTT ^{3'}	5'-ggA ^{3'}	5'-TCT ^{3'}	5'-TCg ^{3'}	5'-g.T ^{3'}	5'-gCA ^{3'}	5'-Cgg ^{3'}
	414	219	362	355	180	98	302	106	261	78	261	341
	5'-CAg ^{3'}	5'-gCA ^{3'}	5'-ggT ^{3'}	5'-CCg ^{3'}	5'-TTT ^{3'}	5'-CAC ^{3'}	5'-ggA ^{3'}	5'-CCA ^{3'}	5'-AAC ^{3'}	5'-TCT ^{3'}	5'-AAC ^{3'}	5'-ggC ^{3'}
	423	238	363		203		341	2 nd I	502	521	341	
	5'-gCT ^{3'}	5'-AgA ^{3'}	5'-ATA ^{3'}		5'-gAA ^{3'}		5'-ggA ^{3'}	5'-CCT ^{3'}	5'-CCC ^{3'}	5'-ggT ^{3'}	5'-ggA ^{3'}	
		355	363		418					527	355	
		5'-CCg ^{3'}	5'-ATA ^{3'}		5'-AgC ^{3'}					5'-TgC ^{3'}	5'-CCC ^{3'}	
		489	363								362	
		5'-gCA ^{3'}	5'-ATA ^{3'}								5'-gAg ^{3'}	
			369								362	
			5'-TAC ^{3'}								5'-gAg ^{3'}	
3'-primer(s) ³	282	180	299	238	290	256	397	265	97	265	292	418
	5'-gAC ^{3'}	5'-TCA ^{3'}	5'-CCA ^{3'}	5'-CCT ^{3'}	5'-CAA ^{3'}	5'-CCC ^{3'}	5'-gAg ^{3'}	5'-CCC ^{3'}	5'-ggT ^{3'}	5'-CCC ^{3'}	5'-gTg ^{3'}	5'-gTC ^{3'}
	282	257	411	238	299	259		282	355	282	292	
	5'-gAC ^{3'}	5'-gCA ^{3'}	5'-TCA ^{3'}	5'-CCT ^{3'}	5'-TCT ^{3'}	5'-gTT ^{3'}		5'-gAC ^{3'}	5'-gAC ^{3'}	5'-gAC ^{3'}	5'-gTT ^{3'}	
	559	418	526	243	317	265		282	524	282	299	
	5'-CCC ^{3'}	5'-gTC ^{3'}	5'-CCA ^{3'}	5'-TCA ^{3'}	5'-ggA ^{3'}	5'-CCC ^{3'}		5'-gAC ^{3'}	5'-CAT ^{3'}	5'-gAC ^{3'}	5'-TCT ^{3'}	
	560	555		265	555	265		502		570	555	
	5'-ACC ^{3'}	5'-CCA ^{3'}		5'-CCC ^{3'}	5'-CCA ^{3'}	5'-CCC ^{3'}		5'-CTT ^{3'}		5'-CAC ^{3'}	5'-CCA ^{3'}	
				282				506				
				5'-gAC ^{3'}				5'-TgT ^{3'}				
				555								
				5'-CCA ^{3'}								
Well No.	13	14	15	16	17	18	19	20	21	22	23	24

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

HLA-B LOW PRIMER SPECIFICATION

Well No.	25	26	27	28	29	30	31	32	33	34	35	36
Length of spec. PCR product	110	215	140	140	185	190	180	165	165	180	100	280
			235	265	235		290	225	190		195	
								285	390			
								330				
Length of int. pos. control ¹	800	1070	1070	800	800	800	1070	1070	1070	800	1070	1070
5'-primer(s) ²	527	97	209	103	103	103	45	45	45	420	206	419
	5'-TgA 3'	5'-TCg 3'	5'-ggC 3'	5'-CCg 3'	5'-CCT 3'	5'-CCg 3'	5'-ggA 3'	5'-ggA 3'	5'-ggA 3'	5'-TTA 3'	5'-AgA 3'	5'-gTC 3'
	539		363	103	103	418		357	435		416	
	5'-gCg 3'		5'-AgC 3'	5'-CCg 3'	5'-CCT 3'	5'-Agg 3'		5'-Tgg 3'	5'-AAA 3'		5'-CCg 3'	
	540		363	361	363			412			420	
	5'-gAC 3'		5'-AgC 3'	5'-AgT 3'	5'-AAT 3'			5'-ATA 3'			5'-TTA 3'	
				363							420	
				5'-AgC 3'							5'-TTC 3'	
3'-primer(s) ³	605	272	309	193	246	246	1 st I	206	266	559	256	3 rd I
	5'-gCT 3'	5'-Tgg 3'	5'-gTg 3'	5'-CgT 3'	5'-TAT 3'	5'-TAT 3'	5'-gAC 3'	5'-CCT 3'	5'-TCC 3'	5'-CTC 3'	5'-CCC 3'	5'-TAT 3'
		272	312	193	246	572	165	538	559		272	
		5'-TgA 3'	5'-gCC 3'	5'-CgT 3'	5'-TAT 3'	5'-gCC 3'	5'-Tgg 3'	5'-gTC 3'	5'-CAG 3'		5'-Tgg 3'	
			559	463	559			603	583		572	
			5'-CAG 3'	5'-gCT 3'	5'-CTC 3'			5'-gTg 3'	5'-gTg 3'		5'-gCg 3'	
				583								
				5'-gTg 3'								
Well No.	25	26	27	28	29	30	31	32	33	34	35	36

Well No.	37	38	39	40	41	42	43	44	45	46	47	48
Length of spec. PCR product	195	105	115	80	100	140	55	210	170	110	395	160
		395	150		150		245					425
		435					390					
							415					
Length of int. pos. control ¹	1070	1070	1070	1070	800	1070	1070	800	800	1070	1070	1070
5'-primer(s) ²	142	45	161	167	349	206	45	142	409	246	44	44
	5'-TCT 3'	5'-ggA 3'	5'-Cgg 3'	5'-gCT 3'	5'-CTT 3'	5'-gAC 3'	5'-ggA 3'	5'-TCT 3'	5'-ggC 3'	5'-gAA 3'	5'-ggC 3'	5'-ggC 3'
	419	540			355		368	368	420			357
	5'-gTC 3'	5'-gAC 3'			5'-TCA 3'		5'-gTT 3'	5'-gTC 3'	5'-TTA 3'			5'-Tgg 3'
							557					
							5'-ggA 3'					
3'-primer(s) ³	301	272	234	204	412	301	259	311	544	317	272	302
	5'-gTC 3'	5'-TgC 3'	5'-TCT 3'	5'-TCT 3'	5'-gTT 3'	5'-gCT 3'	5'-gTT 3'	5'-ggg 3'	5'-ggT 3'	5'-ggA 3'	5'-TgC 3'	5'-ggC 3'
	301	272	272		463	302	259	311		320		302
	5'-gTC 3'	5'-TgT 3'	5'-TgA 3'		5'-gCT 3'	5'-ggC 3'	5'-gTT 3'	5'-ggg 3'		5'-AgT 3'		5'-ggT 3'
	570	309	272			311	266	538				477
	5'-CCg 3'	5'-ATC 3'	5'-Tgg 3'			5'-ggT 3'	5'-TCC 3'	5'-gTC 3'				5'-gCg 3'
	572	605				312	272					
	5'-gCg 3'	5'-gCT 3'				5'-AgT 3'	5'-Tgg 3'					
							292					
							5'-gTA 3'					
							572					
							5'-gCg 3'					
Well No.	37	38	39	40	41	42	43	44	45	46	47	48

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

Well No.	49	50	51	52	53	54	55	56	57	58	59	60
Length of spec. PCR product	180	290	105	325	80	215	130	90	90	145	120	430
		465			115		280	410	175	430	210	
					160							
					195							
					225							
					260							
Length of int. pos. control ¹	1070	1070	1070	800	1070	800	800	1070	1070	1070	800	1070
5'-primer(s) ²	355	72	540	1 st	355	141	209	41	209	48	357	49
	5'-TCA 3'	5'-gCT 3'	5'-gAC 3'	5'-CAg 3'	5'-TCA 3'	5'-ATT 3'	5'-ggC 3'	5'-CTg 3'	5'-ggC 3'	5'-gCC 3'	5'-Tgg 3'	5'-CAg 3'
	363	209			418	357	209	368	363	206		
	5'-Agg 3'	5'-ggC 3'			5'-Agg 3'	5'-Tgg 3'	5'-ggg 3'	5'-gTT 3'	5'-AgC 3'	5'-gAA 3'		
					499		355					
					5'-TCT 3'		5'-TCA 3'					
3'-primer(s) ³	499	193	603	282	538	309	299	282	259	309	435	309
	5'-ggA 3'	5'-CgT 3'	5'-gTg 3'	5'-gCC 3'	5'-gTC 3'	5'-ATC 3'	5'-TCA 3'	5'-gCC 3'	5'-gTT 3'	5'-ATC 3'	5'-TCT 3'	5'-ATC 3'
		193	605		572	538	583	418	499		527	
		5'-TgT 3'	5'-gCT 3'		5'-gCA 3'	5'-gTC 3'	5'-gTg 3'	5'-gTC 3'	5'-ggA 3'		5'-CCT 3'	
		387			572		605					
		5'-TCg 3'			5'-gCg 3'		5'-gCT 3'					
					573							
					5'-AgT 3'							
Well No.	49	50	51	52	53	54	55	56	57	58	59	60

Well No.	61	62	63	64	65	66	67	68	69	70	71	72
Length of spec. PCR product	145	300	160	180	180	90	90	95	115	360	350	285
			370	210	205	240			150			
Length of int. pos. control ¹	1070	1070	1070	1070	1070	800	1070	1070	1070	1070	1070	800
5'-primer(s) ²	206	1 st	15	141	97	209	207	206	165	1 st	1 st	355
	5'-gAA 3'	5'-CAg 3'	5'-gCA 3'	5'-ATT 3'	5'-TCg 3'	5'-ggC 3'	5'-ACC 3'	5'-gAC 3'	5'-ACC 3'	5'-CAg 3'	5'-CAg 3'	5'-TCA 3'
			418	420	420	209	208		463			363
			5'-Agg 3'	5'-TTA 3'	5'-TTA 3'	5'-ggA 3'	5'-CgT 3'		5'-TgA 3'			5'-AgC 3'
						362	209					
						5'-ggT 3'	5'-ggA 3'					
							369					
							5'-TAC 3'					
3'-primer(s) ³	311	259	207	311	259	256	256	259	272	317	311	603
	5'-ggT 3'	5'-CTC 3'	5'-TCC 3'	5'-ggT 3'	5'-gTT 3'	5'-CCC 3'	5'-CTC 3'	5'-gTT 3'	5'-TgC 3'	5'-ggA 3'	5'-ggT 3'	5'-gTg 3'
		259	226	559	559	259	256		538	320		605
		5'-CTC 3'	5'-CAC 3'	5'-CgT 3'	5'-CAg 3'	5'-CTT 3'	5'-CCC 3'		5'-CCA 3'	5'-AgT 3'		5'-gCT 3'
		262	538			263	420					
		5'-TgC 3'	5'-gTC 3'			5'-gTT 3'	5'-gCT 3'					
						559						
						5'-CgT 3'						
Well No.	61	62	63	64	65	66	67	68	69	70	71	72

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 25 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

DR LOW PRIMER SPECIFICATION

Well No.	73	74	75	76	77	78	79	80	81	82	83	84
Length of spec.	210	200	200	210	125	75	85	100	170	170	85	175
PCR product	235		230		225	205	210	175	210	215	100	
	260								235	250	135	
											160	
											190	
Length of int.	515	430	430	430	515	430	430	430	430	515	430	430
pos. control ¹												
5'-primer(s) ²	13(124) 5'-A.T 3'	14(129) 5'-gAA 3'	12(121) 5'-CTg 3'	13(126) 5'-Agg 3'	13(125) 5'-gTC 3'	13(125) 5'-gTC 3'	13(125) 5'-gTC 3'	13(125) 5'-ACA 3'	11(118) 5'-AgA 3'	16(133) 5'-gTT 3'	26(165) 5'-TAT 3'	26(164) 5'-gTA 3'
	14(129) 5'-gAA 3'		13(126) 5'-AgA 3'	13(126) 5'-AAg 3'	47(227) 5'-gTT 3'	16(133) 5'-gTT 3'		13(125) 5'-ACC 3'	14(127) 5'-gTA 3'	16(133) 5'-gCT 3'	37(196) 5'-AgA 3'	31(178) 5'-gCg 3'
			13(126) 5'-AAg 3'					13(125) 5'-ATA 3'	14(127) 5'-ATA 3'	16(133) 5'-gTT 3'	58(261) 5'-gAg 3'	
			13(126) 5'-Agg 3'					13(125) 5'-gTC 3'	14(127) 5'-ATA 3'			
									16(133) 5'-gTT 3'			
									31(178) 5'-gCg 3'			
3'-primer(s) ³	67(286) 5'-gAg 3'	67(286) 5'-gAT 3'	67(286) 5'-gAT 3'	67(286) 5'-gAA 3'	73(305) 5'-ggC 3'	26(164) 5'-ggT 3'	28(171) 5'-CTC 3'	33(184) 5'-gTg 3'	71(298) 5'-CTC 3'	58(260) 5'-CCT 3'	57(257) 5'-CgA 3'	72(303) 5'-gCA 3'
	67(286) 5'-gAg 3'		67(286) 5'-gAT 3'	67(286) 5'-gAg 3'	73(305) 5'-ggC 3'	67(286) 5'-gAT 3'	70(295) 5'-CTg 3'	58(260) 5'-Cgg 3'	73(305) 5'-ggC 3'	74(307) 5'-CAg 3'	73(305) 5'-ggC 3'	74(307) 5'-CgC 3'
	67(286) 5'-gAT 3'		67(286) 5'-gAT 3'	70(296) 5'-TgT 3'	73(305) 5'-ggC 3'	71(299) 5'-gCT 3'			77(317) 5'-AAT 3'	86(344) 5'-CAC 3'	78(319) 5'-CAC 3'	74(308) 5'-CCT 3'
	70(297) 5'-CTg 3'		70(295) 5'-Tg 3'	72(302) 5'-Cgg 3'	73(305) 5'-ggC 3'				78(319) 5'-CgC 3'			
	71(299) 5'-gCg 3'		70(295) 5'-CTg 3'		73(306) 5'-Tgg 3'				78(319) 5'-CAA 3'			
	77(317) 5'-AgT 3'		71(299) 5'-gCT 3'		74(308) 5'-CCC 3'				78(319) 5'-CAC 3'			
	86(344) 5'-CCA 3'		71(299) 5'-gCg 3'		75(310) 5'-CAA 3'							
			72(301) 5'-Cg 3'		77(317) 5'-AgT 3'							
			73(305) 5'-ggC 3'									
			77(317) 5'-AgT 3'									
Well No.	73	74	75	76	77	78	79	80	81	82	83	84

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

Well No.	85	86	87	88	89	90	91	92	93	94	95
Length of spec. PCR product	100	90	220	195	175	100	110	110	80	130	175
	175	110		225	240	150	145	150	160	215	
						195	170	180	240		
						240		220			
Length of int. pos. control ¹	430	430	430	430	430	430	430	430	430	430	430
5'-primer(s) ²	13(125) 5'-gTC 3'	13(124) 5'-Cgg 3'	10(116) 5'-gCT 3'	10(116) 5'-gCT 3'	13(125) 5'-gTC 3'	1 st I 5'-CAA 3'	26(164) 5'-gTA 3'	13(125) 5'-gTC 3'	10(116) 5'-gCT 3'	28(170) 5'-gAT 3'	13(125) 5'-gTA 3'
	16(133) 5'-gTC 3'	16(133) 5'-gTT 3'	12(122) 5'-TAT 3'	12(122) 5'-TAT 3'	13(125) 5'-gTg 3'	37(197) 5'-gTT 3'	34(189) 5'-CAg 3'	34(189) 5'-CAg 3'	10(116) 5'-gCT 3'	28(170) 5'-gAg 3'	
	38(200) 5'-CgT 3'		13(125) 5'-gTC 3'	12(123) 5'-ACA 3'	114(429) 5'-CTg 3'	37(197) 5'-gTA 3'	34(189) 5'-CAg 3'	37(196) 5'-AgC 3'	38(199) 5'-TCC 3'	105(401) 5'-AAA 3'	
				13(125) 5'-gTC 3'							
				13(125) 5'-gTg 3'							
				16(133) 5'-gTC 3'							
				16(133) 5'-gTT 3'							
3'-primer(s) ³	58(260) 5'-CCT 3'	28(171) 5'-CTC 3'	70(295) 5'-gTC 3'	67(286) 5'-gAA 3'	58(260) 5'-Cgg 3'	42(213) 5'-TCA 3'	57(257) 5'-CAg 3'	57(257) 5'-CAg 3'	51(239) 5'-CCC 3'	87(346) 5'-CTC 3'	57(256) 5'-gTC 3'
	58(260) 5'-CCT 3'	30(175) 5'-gTg 3'	71(299) 5'-gCT 3'	71(298) 5'-CTC 3'	58(260) 5'-CAg 3'	57(257) 5'-CAg 3'	70(295) 5'-CTg 3'	60(265) 5'-gTg 3'	77(317) 5'-AAT 3'	87(346) 5'-CTT 3'	57(258) 5'-gCg 3'
	58(260) 5'-CCT 3'	38(199) 5'-CAg 3'	71(299) 5'-ACT 3'	71(298) 5'-CgC 3'	181(630) 5'-CTT 3'	71(298) 5'-CgC 3'	70(296) 5'-TCC 3'	70(296) 5'-TCC 3'		135(490) 5'-gCT 3'	58(260) 5'-CCT 3'
	58(261) 5'-TCT 3'		74(307) 5'-CgT 3'					74(307) 5'-CAg 3'			
			78(319) 5'-CAC 3'								
Well No.	85	86	87	88	89	90	91	92	93	94	95

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 73 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

CELL LINE VALIDATION SHEET																				
HLA-A low resolution primer set ²																				
				Well																
				1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
				Prod. No.:	201787801	201785202	201787803	201785204	201674905	201785206	201674907	201785208	201785209	201787810	201785211	201785212	201785213	201674914	201674915	201785216
	IHWC cell line ¹	A*	A*																	
1	9001 SA	*24:02		-	-	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-
2	9280 LK707	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011 E4181324	*01:01		+	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275 GU373	*30:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-
5	9009 KAS011	*01:01		+	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM	*02:01	*26:03	-	+	-	-	-	-	-	-	+	-	-	+	-	+	-	-	-
7	9020 QBL	*26:01		-	-	-	-	-	-	-	-	+	-	+	-	-	+	-	-	-
8	9025 DEU	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
9	9026 YAR	*26:01		-	-	-	-	-	-	-	-	+	-	+	-	-	+	-	-	-
10	9107 LKT3	*24:02		-	-	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT	*29:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-
12	9052 DBB	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
15	9075 DKB	*24:02		-	-	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007	*29:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
17	9282 CTM3953540	*03:01	*80:01	-	-	+	-	-	+	-	-	-	-	-	-	+	-	-	-	-
18	9257 32367	*33:03	*74:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA	*02:17		-	+	-	-	-	-	W	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE	*02:01		-	+	-	-	-	-	W	-	-	-	-	-	-	-	-	-	-
23	9124 IHL	*02:01	*34:01	-	+	-	-	-	-	-	+	-	-	+	+	-	-	-	-	-
24	9035 JBUSH	*32:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
25	9049 IBW9	*33:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49	*02:05		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007	*24:10	*29:01	-	-	-	-	+	-	+	-	-	-	-	-	-	+	-	-	-
28	9320 BEL5GB	*02:01	*29:02	-	+	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
29	9050 MOU	*29:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
30	9021 RSH	*30:01	*68:02	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	+	-
31	9019 DUCAF	*30:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
32	9297 HAG	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
34	9104 DHIF	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
35	9302 SSTO	*32:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
36	9024 KT17	*02:06	*11:01	-	+	-	+	-	-	-	-	-	-	+	-	-	-	-	-	-
37	9065 HHKB	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	+	-	-	-	-
38	9099 LZL	*02:17		-	+	-	-	-	-	W	-	-	-	-	-	-	-	-	-	-
39	9315 CML	*01:01	*03:01	+	-	+	+	-	-	-	-	-	-	-	-	+	-	-	-	-
40	9134 WHONP199	*02:07	*30:01	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
41	9055 H0301	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	+	-	-	-	-
42	9066 TAB089	*02:07		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526	*02:06	*02:07	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM	*66:01		-	-	-	-	-	-	-	+	-	-	+	-	+	-	-	-	-
45	9239 SHJO	*23:01	*24:02	-	-	-	-	+	+	+	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	+	-	-	-	-
47	9045 TUBO	*02:16	*03:01	-	+	+	-	-	-	-	-	-	-	-	-	+	-	-	-	-
48	9303 TER-ND	*02:01	*11:01	-	+	-	+	-	-	-	-	-	-	+	-	-	-	-	-	-

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

CELL LINE VALIDATION SHEET													
HLA-A low resolution primer set ²													
					Well								
					17	18	19	20	21	22	23	24	
					Prod. No.:	201787817	201785218	201674919	201674920	201787821	201787822	201779223	201787824
	IHCW cell line ¹		A*	A*									
1	9001	SA	*24:02		-	-	-	-	-	-	-	-	-
2	9280	LK707	*02:01		-	-	-	-	-	-	-	-	-
3	9011	E4181324	*01:01		-	-	-	-	-	-	-	-	-
4	9275	GU373	*30:01		-	-	-	-	-	-	-	-	-
5	9009	KAS011	*01:01		-	-	-	-	-	-	-	-	-
6	9353	SM	*02:01	*26:03	-	-	-	-	-	-	-	+	-
7	9020	QBL	*26:01		-	-	-	-	-	-	-	-	-
8	9025	DEU	*31:01		-	-	-	-	-	-	-	-	-
9	9026	YAR	*26:01		-	-	-	-	-	-	-	-	-
10	9107	LKT3	*24:02		-	-	-	-	-	-	-	-	-
11	9051	PITOUT	*29:02		-	-	-	-	-	-	-	-	-
12	9052	DBB	*02:01		-	-	-	-	-	-	-	-	-
13	9004	JESTHOM	*02:01		-	-	-	-	-	-	-	-	-
14	9071	OLGA	*31:01		-	-	-	-	-	-	-	-	-
15	9075	DKB	*24:02		-	-	-	-	-	-	-	-	-
16	9037	SWEIG007	*29:02		-	-	-	-	-	-	-	-	-
17	9282	CTM3953540	*03:01	*80:01	-	-	-	-	-	-	-	+	-
18	9257	32367	*33:03	*74:01	-	+	+	-	-	-	-	-	-
19	9038	BM16	*02:01		-	-	-	-	-	-	-	-	-
20	9059	SLE005	*02:01		-	-	-	-	-	-	-	-	-
21	9064	AMALA	*02:17		-	-	-	-	-	-	-	-	-
22	9056	KOSE	*02:01		-	-	-	-	-	-	-	-	-
23	9124	IHL	*02:01	*34:01	-	-	-	-	-	-	-	-	-
24	9035	JBUSH	*32:01		+	-	-	-	-	-	-	-	-
25	9049	IBW9	*33:01		-	+	-	-	-	-	-	-	-
26	9285	WT49	*02:05		-	-	-	-	-	-	-	-	-
27	9191	CH1007	*24:10	*29:01	-	-	-	-	-	-	-	-	-
28	9320	BEL5GB	*02:01	*29:02	-	-	-	-	-	-	-	-	-
29	9050	MOU	*29:02		-	-	-	-	-	-	-	-	-
30	9021	RSH	*30:01	*68:02	-	-	-	+	-	-	-	-	-
31	9019	DUCAF	*30:02		-	-	-	-	-	-	-	-	-
32	9297	HAG	*02:01		-	-	-	-	-	-	-	-	-
33	9098	MT14B	*31:01		-	-	-	-	-	-	-	-	-
34	9104	DHIF	*31:01		-	-	-	-	-	-	-	-	-
35	9302	SSTO	*32:01		+	-	-	-	-	-	-	-	-
36	9024	KT17	*02:06	*11:01	-	-	-	-	-	-	-	-	-
37	9065	HHKB	*03:01		-	-	-	-	-	-	-	-	-
38	9099	LZL	*02:17		-	-	-	-	-	-	-	-	-
39	9315	CML	*01:01	*03:01	-	-	-	-	-	-	-	-	-
40	9134	WHONP199	*02:07	*30:01	-	-	-	-	-	-	-	-	-
41	9055	H0301	*03:01		-	-	-	-	-	-	-	-	-
42	9066	TAB089	*02:07		-	-	-	-	-	-	-	-	-
43	9076	T7526	*02:06	*02:07	-	-	-	-	-	-	-	-	-
44	9057	TEM	*66:01		-	-	-	-	-	-	-	-	-
45	9239	SHJO	*23:01	*24:02	-	-	-	-	-	-	-	-	-
46	9013	SCHU	*03:01		-	-	-	-	-	-	-	-	-
47	9045	TUBO	*02:16	*03:01	-	-	-	-	-	-	-	-	-
48	9303	TER-ND	*02:01	*11:01	-	-	-	-	-	-	-	-	-

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

Additional 5'- and 3'-primers in primer solutions 4, 5, 7, 8, 12, 14 to 17 and 20 to 23 were tested by separately adding one or more 3'-primers, respectively one or more 5'-primers.

Additional 5'-primers in primer solution 1, 9, 13, 19 and 24 were tested by separately adding one or two 3'-primers. Additional 3'-primers in primer solutions 3, 10 and 18 were tested by separately adding one or two 5'-primers. In primer solutions 2, 3, 9, 11,12, 15 and 22 one or two 5'-primers were not possible to test, and in primer solutions 3, 6, 8, 9, 12, 13, 16 and 18 one or two 3'-primers were not possible to test.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

CELL LINE VALIDATION SHEET																			
HLA-B low resolution SSP kit²																			
			Prod. No.:	Well															
				25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40
				201673601	201673602	201673603	201673604	201673605	201673606	201673607	201673608	201785309	201776710	201776711	201673612	201673613	201673614	201673615	201673616
	IHWC cell line ¹	B*																	
1	9001 SA	*07:02		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707	*52:01	*73:01	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
3	9011 E4181324	*52:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
4	9275 GU373	*15:10	*53:01	-	-	-	-	-	-	-	+	-	-	-	+	-	+	-	-
5	9009 KAS011	*37:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
6	9353 SM	*39:01	*51:01	-	-	-	-	+	-	-	-	+	-	-	-	-	-	-	-
7	9020 QBL	*18:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-
8	9025 DEU	*35:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
9	9026 YAR	*38:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3	*54:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT	*44:03		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB	*57:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
13	9004 JESTHOM	*27:05		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	+
14	9071 OLGA	*15:01	*15:20	-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	-
15	9075 DKB	*40:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007	*40:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282 CTM3953540	*08:01	*55:01	-	+	-	-	-	-	-	+	-	-	-	-	-	-	-	-
18	9257 32367	*14:01	*56:01	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16	*18:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-
20	9059 SLE005	*40:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA	*15:01		-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	-
22	9056 KOSE	*35:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL	*40:02	*56:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035 JBUSH	*38:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9	*14:02		-	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-
26	9285 WT49	*58:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
27	9191 CH1007	*07:05	*51:01	+	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
28	9320 BEL5GB	*44:02	*44:03	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU	*44:03		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH	*42:01		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
31	9019 DUCAF	*18:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-
32	9297 HAG	*41:02		-	-	-	+	-	-	-	+	-	-	-	-	-	-	-	-
33	9098 MT14B	*40:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF	*38:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO	*44:02		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17	*15:01	*35:01	-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	-
37	9065 HHKB	*07:02		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL	*15:01		-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	-
39	9315 CML	*08:01	*27:05	-	+	-	-	+	-	-	+	-	-	-	-	-	-	-	+
40	9134 WHONP199	*13:02	*46:01	-	-	+	+	-	-	+	-	-	+	-	+	-	-	-	-
41	9055 H0301	*14:02		-	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089	*46:01		-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	-
43	9076 T7526	*46:01		-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	-
44	9057 TEM	*38:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO	*42:01	*50:01	-	-	-	+	-	-	-	+	-	-	-	-	-	-	-	-
46	9013 SCHU	*07:02		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO	*51:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
48	9303 TER-ND	*35:01	*44:03	-	-	-	+	-	-	-	-	-	-	-	+	-	-	-	-

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

CELL LINE VALIDATION SHEET																				
HLA-B low resolution SSP kit ²																				
				Well																
				41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	
				Prod. No.:	201787917	201785318	201673619	201673620	201776721	201785322	201673623	201673624	201673625	201673626	201673627	201787928	201673629	201673630	201673631	201787932
	IHC cell line ¹	B*																		
1	9001 SA	*07:02		-	-	-	-	-	-	-	-	+	-	-	-	+	-	-	-	
2	9280 LK707	*52:01 *73:01		-	-	-	-	-	-	-	-	+	-	-	-	+	-	-	-	
3	9011 E4181324	*52:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
4	9275 GU373	*15:10 *53:01		+	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
5	9009 KAS011	*37:01		+	-	-	+	-	-	-	-	-	-	-	-	+	-	-	-	
6	9353 SM	*39:01 *51:01		-	-	-	-	+	-	+	+	-	-	-	-	-	+	-	-	
7	9020 QBL	*18:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
8	9025 DEU	*35:01		+	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
9	9026 YAR	*38:01		-	-	-	-	+	+	+	-	-	-	-	-	-	-	-	-	
10	9107 LKT3	*54:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	+	
11	9051 PITOUT	*44:03		+	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
12	9052 DBB	*57:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
13	9004 JESTHOM	*27:05		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	+	
14	9071 OLGA	*15:01 *15:20		+	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
15	9075 DKB	*40:01		-	-	-	-	-	-	-	+	+	-	-	-	-	-	-	-	
16	9037 SWEIG007	*40:02		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
17	9282 CTM3953540	*08:01 *55:01		-	-	-	-	-	-	-	+	-	-	+	+	-	-	-	+	
18	9257 32367	*14:01 *56:01		-	-	-	-	-	-	+	+	-	-	+	-	-	-	-	+	
19	9038 BM16	*18:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
20	9059 SLE005	*40:01		-	-	-	-	-	-	-	+	+	-	-	-	-	-	-	-	
21	9064 AMALA	*15:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
22	9056 KOSE	*35:03		+	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
23	9124 IHL	*40:02 *56:02		-	-	-	-	-	-	-	-	-	+	-	+	-	-	-	+	
24	9035 JBUSH	*38:01		-	-	-	-	+	+	+	-	-	-	-	-	-	-	-	-	
25	9049 IBW9	*14:02		-	-	-	-	-	-	+	+	-	-	-	-	-	-	-	-	
26	9285 WT49	*58:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
27	9191 CH1007	*07:05 *51:01		-	-	-	-	-	-	-	+	-	-	+	-	+	-	-	-	
28	9320 BEL5GB	*44:02 *44:03		+	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
29	9050 MOU	*44:03		+	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
30	9021 RSH	*42:01		-	-	-	-	-	-	-	+	-	-	+	+	-	-	-	-	
31	9019 DUCAF	*18:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
32	9297 HAG	*41:02		-	-	-	-	-	-	-	+	-	-	+	-	-	-	-	-	
33	9098 MT14B	*40:01		-	-	-	-	-	-	-	+	+	-	-	-	-	-	-	-	
34	9104 DHIF	*38:01		-	-	-	-	+	+	+	-	-	-	-	-	-	-	-	-	
35	9302 SSTO	*44:02		+	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
36	9024 KT17	*15:01 *35:01		+	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
37	9065 HHKB	*07:02		-	-	-	-	-	-	-	+	-	-	-	+	-	-	-	-	
38	9099 LZL	*15:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
39	9315 CML	*08:01 *27:05		-	-	-	-	-	-	-	+	-	-	+	+	-	-	-	+	
40	9134 WHONP199	*13:02 *46:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
41	9055 H0301	*14:02		-	-	-	-	-	-	+	+	-	-	-	-	-	-	-	-	
42	9066 TAB089	*46:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
43	9076 T7526	*46:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
44	9057 TEM	*38:01		-	-	-	-	+	+	+	-	-	-	-	-	-	-	-	-	
45	9239 SHJO	*42:01 *50:01		-	-	-	-	-	-	-	+	-	-	+	+	-	-	-	-	
46	9013 SCHU	*07:02		-	-	-	-	-	-	-	+	-	-	-	+	-	-	-	-	
47	9045 TUBO	*51:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	
48	9303 TER-ND	*35:01 *44:03		+	+	-	-	-	-	-	-	-	-	-	-	+	-	-	-	

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

CELL LINE VALIDATION SHEET																				
HLA-B low resolution SSP kit ²																				
				Well																
				57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	
				Prod. No.:	201673633	201673634	201673635	201673636	201673637	201673638	201673639	201673640	201776741	201673642	201673643	201673644	201673645	201785346	201673647	201673648
	IHC cell line ¹	B*																		
1	9001 SA	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	+	
2	9280 LK707	*52:01	*73:01	-	-	+	+	-	+	-	-	-	-	-	-	+	+	+	-	
3	9011 E4181324	*52:01		-	-	+	+	-	+	-	-	-	-	-	-	-	+	-	-	
4	9275 GU373	*15:10	*53:01	-	-	-	+	-	-	-	-	-	-	-	+	-	+	+	-	
5	9009 KAS011	*37:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
6	9353 SM	*39:01	*51:01	-	-	+	+	-	-	-	-	-	-	-	+	-	+	+	-	
7	9020 QBL	*18:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
8	9025 DEU	*35:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
9	9026 YAR	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	
10	9107 LKT3	*54:01		-	-	-	-	-	-	-	+	+	-	-	-	-	-	+	-	
11	9051 PITOUT	*44:03		-	-	-	-	-	+	-	-	-	-	-	-	-	+	-	-	
12	9052 DBB	*57:01		-	-	-	+	-	-	-	-	-	-	+	-	-	-	-	-	
13	9004 JESTHOM	*27:05		-	-	-	-	-	+	-	-	-	-	-	-	-	+	-	-	
14	9071 OLGA	*15:01	*15:20	-	-	-	-	-	+	-	-	-	-	-	-	+	-	+	-	
15	9075 DKB	*40:01		-	-	-	-	-	+	+	-	-	-	-	-	-	-	+	-	
16	9037 SWEIG007	*40:02		-	-	-	-	-	+	+	-	-	-	-	-	-	-	+	-	
17	9282 CTM3953540	*08:01	*55:01	-	-	+	-	-	-	-	-	+	-	-	-	-	-	+	+	
18	9257 32367	*14:01	*56:01	-	-	-	-	-	-	-	-	-	+	-	-	-	-	+	-	
19	9038 BM16	*18:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
20	9059 SLE005	*40:01		-	-	-	-	-	+	+	-	-	-	-	-	-	-	+	-	
21	9064 AMALA	*15:01		-	-	-	-	-	-	+	-	-	-	-	-	+	-	+	-	
22	9056 KOSE	*35:03		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
23	9124 IHL	*40:02	*56:02	-	-	-	-	-	+	+	-	-	+	-	+	-	-	+	-	
24	9035 JBUSH	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	
25	9049 IBW9	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
26	9285 WT49	*58:01		-	-	-	+	-	-	-	-	-	-	+	-	-	+	-	-	
27	9191 CH1007	*07:05	*51:01	-	-	+	+	-	-	-	-	-	-	-	+	-	+	+	+	
28	9320 BEL5GB	*44:02	*44:03	-	-	-	-	-	-	+	+	-	-	-	-	-	+	-	-	
29	9050 MOU	*44:03		-	-	-	-	-	-	+	-	-	-	-	-	-	+	-	-	
30	9021 RSH	*42:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	+	
31	9019 DUCAF	*18:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
32	9297 HAG	*41:02		-	-	-	-	+	+	-	-	-	-	-	-	-	-	+	+	
33	9098 MT14B	*40:01		-	-	-	-	+	+	-	-	-	-	-	-	-	-	+	-	
34	9104 DHIF	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	
35	9302 SSTO	*44:02		-	-	-	-	-	+	+	-	-	-	-	-	-	-	+	-	
36	9024 KT17	*15:01	*35:01	-	-	-	-	-	-	+	-	-	-	-	+	+	-	+	-	
37	9065 HHKB	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	+	
38	9099 LZL	*15:01		-	-	-	-	-	+	-	-	-	-	-	-	+	-	+	-	
39	9315 CML	*08:01	*27:05	-	-	-	-	-	-	+	-	-	-	-	-	-	+	+	+	
40	9134 WHONP199	*13:02	*46:01	-	-	-	-	-	-	+	-	-	-	-	-	+	+	+	-	
41	9055 H0301	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
42	9066 TAB089	*46:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	
43	9076 T7526	*46:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	
44	9057 TEM	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
45	9239 SHJO	*42:01	*50:01	-	-	+	-	+	+	-	-	+	-	-	-	-	-	+	+	
46	9013 SCHU	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	+	
47	9045 TUBO	*51:01		-	-	+	+	-	-	-	-	-	-	-	+	-	+	-	-	
48	9303 TER-ND	*35:01	*44:03	-	-	-	-	-	-	+	-	-	-	-	+	-	+	+	-	

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

Additional 5'-primers and 3'-primers in primer solutions 27, 30, 35, 38, 43, 53, 55, 64 and 65 were tested by separately adding one or two additional 3'-primers, respectively one or two additional 5'-primers. Additional 3'-primers in primer solutions 26, 32, 37, 39 and 42 were tested by separately adding one, two or three additional 5'-primers. Additional 5'-primers in primer solutions 25, 45, 49 and 72 were tested by separately adding one or two additional 3'-primers.

In primer mixes 35, 41, 55, 66 and 67 one or two 5'-primers were not possible to test, and in primer mixes 27, 31, 46, 53, 62, 63, 66 and 70 one or two 3'-primer were not possible to test.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

CELL LINE VALIDATION SHEET																				
DR low resolution primer set ²																				
				Well																
				73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	
				Prod. No.:	201673401	201673402	201673403	201785404	201788005	201673406	201673407	201673408	201673409	201675110	201675111	201675112	201673413	201673414	201782715	201673416
IHCW cell line ¹		DRB1																		
1	9001 SA	*01:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
2	9280 LK707	*15:02	*04:05	-	-	+	-	-	-	-	-	+	-	-	-	-	-	-	-	
3	9011 E4181324	*15:02		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
4	9275 GU373	*03:01		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-	
5	9009 KAS011	*16:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
6	9353 SM	*04:07	*08:03	-	-	-	-	-	-	-	-	+	-	+	-	-	-	-	-	
7	9020 QBL	*03:01		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-	
8	9025 DEU	*04:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	
9	9026 YAR	*04:02		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	
10	9107 LKT3	*04:05		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	
11	9051 PITOUT	*07:01		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
12	9052 DBB	*07:01		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
13	9004 JESTHOM	*01:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
14	9071 OLGA	*08:02		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	+	
15	9075 DKB	*09:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	
16	9037 SWEIG007	*11:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	+	+	
17	9282 CTM3953540	*03:01	*13:01	-	-	-	-	+	+	-	-	-	-	-	-	-	-	+	+	
18	9257 32367	*09:01	*11:01	-	-	-	-	-	-	-	-	-	-	+	-	+	-	+	+	
19	9038 BM16	*12:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	
20	9059 SLE005	*13:02		-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	+	
21	9064 AMALA	*14:02		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
22	9056 KOSE	*13:02	*14:54	-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	+	
23	9124 IHL	*08:03	*14:14	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
24	9035 JBUSH	*11:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	+	+	
25	9049 IBW9	*07:01		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
26	9285 WT49	*03:01		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-	
27	9191 CH1007	*04:05	*10:01	-	-	-	-	-	-	-	-	+	-	+	-	+	-	-	-	
28	9320 BEL5GB	*04:16	*07:01	-	-	-	-	-	-	-	-	+	+	-	-	-	-	-	-	
29	9050 MOU	*07:01		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
30	9021 RSH	*03:02		-	-	-	-	+	+	+	-	-	-	-	-	-	-	-	-	
31	9019 DUCAF	*03:01		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-	
32	9297 HAG	*13:03		-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	-	
33	9098 MT14B	*04:04		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	
34	9104 DHIF	*11:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	+	+	
35	9302 SSTO	*04:03		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	
36	9024 KT17	*04:03	*04:06	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	
37	9065 HHKB	*13:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	+	
38	9099 LZL	*14:02		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
39	9315 CML	*03:01	*04:01	-	-	-	-	+	+	-	+	-	-	-	-	-	-	-	-	
40	9134 WHONP199	*07:01	*09:01	-	-	-	-	-	-	-	-	-	+	-	+	-	-	-	-	
41	9055 H0301	*13:02		-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	+	
42	9066 TAB089	*08:03		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	
43	9076 T7526	*09:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	
44	9057 TEM	*14:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
45	9239 SHJO	*07:01		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
46	9013 SCHU	*15:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
47	9045 TUBO	*11:04	*12:01	-	-	-	-	-	-	-	-	-	-	-	-	+	+	+	+	
48	9303 TER-ND	*01:03		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

CELL LINE VALIDATION SHEET											
DR low resolution primer set ²											
				Well							
				89	90	91	92	93	94	95	
				Prod. No.:	201673417	201673418	201785419	201673420	201673429	201675122	201782731
IHCW cell line ¹		DRB1									
1	9001 SA	*01:01		-	-	-	-	-	-	-	-
2	9280 LK707	*15:02	*04:05	-	-	-	-	-	+	+	
3	9011 E4181324	*15:02		-	-	-	-	-	-	-	+
4	9275 GU373	*03:01		+	-	-	-	+	-	-	
5	9009 KAS011	*16:01		-	-	-	-	-	-	-	+
6	9353 SM	*04:07	*08:03	-	-	-	-	-	+	-	
7	9020 QBL	*03:01		+	-	-	-	+	-	-	
8	9025 DEU	*04:01		-	-	-	-	-	+	-	
9	9026 YAR	*04:02		-	-	-	-	-	+	-	
10	9107 LKT3	*04:05		-	-	-	-	-	+	-	
11	9051 PITOUT	*07:01		-	-	-	-	-	+	-	
12	9052 DBB	*07:01		-	-	-	-	-	+	-	
13	9004 JESTHOM	*01:01		-	-	-	-	-	-	-	
14	9071 OLGA	*08:02		-	-	-	-	-	-	-	
15	9075 DKB	*09:01		-	-	+	-	-	+	-	
16	9037 SWEIG007	*11:01		-	-	-	-	+	-	-	
17	9282 CTM3953540	*03:01	*13:01	+	-	-	-	+	-	-	
18	9257 32367	*09:01	*11:01	-	-	+	-	+	+	-	
19	9038 BM16	*12:01		-	-	-	-	+	-	-	
20	9059 SLE005	*13:02		+	-	-	-	+	-	-	
21	9064 AMALA	*14:02		+	-	+	-	+	-	-	
22	9056 KOSE	*13:02	*14:54	+	+	+	+	+	-	-	
23	9124 IHL	*08:03	*14:14	+	-	+	+	+	-	-	
24	9035 JBUSH	*11:01		-	-	-	-	+	-	-	
25	9049 IBW9	*07:01		-	-	-	-	-	+	-	
26	9285 WT49	*03:01		+	-	-	-	+	-	-	
27	9191 CH1007	*04:05	*10:01	-	-	-	-	-	+	-	
28	9320 BEL5GB	*04:16	*07:01	-	-	-	-	-	+	-	
29	9050 MOU	*07:01		-	-	-	-	-	+	-	
30	9021 RSH	*03:02		+	-	-	-	+	-	-	
31	9019 DUCAF	*03:01		+	-	-	-	+	-	-	
32	9297 HAG	*13:03		+	-	-	-	+	-	-	
33	9098 MT14B	*04:04		-	-	-	-	-	+	-	
34	9104 DHIF	*11:01		-	-	-	-	+	-	-	
35	9302 SSTO	*04:03		-	-	-	-	-	+	-	
36	9024 KT17	*04:03	*04:06	-	-	-	-	-	+	-	
37	9065 HHKB	*13:01		+	-	-	-	+	-	-	
38	9099 LZL	*14:02		+	-	+	-	+	-	-	
39	9315 CML	*03:01	*04:01	+	-	-	-	+	+	-	
40	9134 WHONP199	*07:01	*09:01	-	-	+	-	-	+	-	
41	9055 H0301	*13:02		+	-	-	-	+	-	-	
42	9066 TAB089	*08:03		-	-	-	-	-	-	-	
43	9076 T7526	*09:01		-	-	+	-	-	+	-	
44	9057 TEM	*14:01		+	+	+	+	+	-	-	
45	9239 SHJO	*07:01		-	-	-	-	-	+	-	
46	9013 SCHU	*15:01		-	-	-	-	-	-	+	
47	9045 TUBO	*11:04	*12:01	-	-	-	-	+	-	-	
48	9303 TER-ND	*01:03		-	-	-	-	-	-	-	

101.701-24/06 – including *Taq* pol., IFU-01
 101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
 “Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

Additional 5'- and 3'-primers in primer solutions 81, 83, 84, 87 and 90 were tested by separately adding one 3'-primer, respectively one or three 5'-primers.

Additional 5'-primers in primer solutions 78 and 88 were tested by separately adding one 3'-primer. Additional 3'-primers in primer solutions 73, 75, 77, 82 and 92 were tested by separately adding one or more 5'-primers.

One or more of the 5'-primers in primer solution 73, 75, 76, 80 to 82, 85 to 89, 92 and 94 were not possible to test. One or more of the 3'-primers in primer solution 73, 75 to 78, 81, 85, 87 and 94 were not possible to test.

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

101.701-24/06 – including *Taq* pol., IFU-01
101.701-24u/06u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **9F4**

Lot-specific information

ADDRESSES:

Manufacturer:

Olerup SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-717 88 27

Fax: +46-8-717 88 18

E-mail: olerup-se@caredx.com

Web page: <http://www.olerup.com>

Distributed by:

Olerup GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

Tel: +43-1-710 15 00

Fax: +43-1-710 15 00 10

E-mail: olerup-at@caredx.com

Web page: <http://www.olerup.com>

Olerup Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-OLERUP1

Fax: 610-344-7989

E-mail: olerup-us@caredx.com

Web page: <http://www.olerup.com>

For information on *Olerup* distributors worldwide, contact **Olerup GmbH**.