

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

Olerup SSP[®] HLA-A-B-DR SSP Combi Tray

Product number:	101.701-24/06 – including <i>Taq</i> pol. 101.701-24u/06u – without <i>Taq</i> pol.
Lot number:	6R6
Expiry date:	2027-01-01
Number of tests:	24 tests – Product No. 101.701-24/24u 6 tests – Product No. 101.701-06/06u
Number of wells per test:	95 +1
Storage - pre-aliquoted primers:	dark, between -15°C and -25°C
- PCR Master Mix:	between -15°C and -25°C
- Adhesive PCR seals	RT

This Product Description is only valid for Lot No. 6R6.

Complete product documentation consists of generic Instructions for Use (IFU), lot specific Product Insert, Worksheet and Certificate.

CHANGES COMPARED TO THE PREVIOUS OLERUP SSP[®] HLA-A-B-DR SSP COMBI TRAY LOT (0R8)

The **HLA-A low resolution** specificity and interpretation tables have been updated for the HLA-A alleles described since the previous *Olerup SSP[®]* HLA-A-B-DR SSP Combi Tray lot was made (**Lot No. 0R8**). The kit design is based on IMGT/HLA database 3.50.0.

The primers of the wells detailed below have been exchanged, modified or added compared to the previous lot.

Well	5'-primer	3'-primer	rationale
2	Added	-	5'-primer added for the A*02:01:207 allele.
6	-	Added	3'-primer added for the A*23:104 allele.
10	-	Added	3'-primer added for the A*26:01:63 allele.
16	-	Added	5'-primer added for the A*23:113N allele.

The **HLA-B low resolution** specificity and interpretation tables have been updated for the HLA-B alleles described since the previous *Olerup SSP[®]* HLA-A-B-DR SSP Combi Tray lot was made (**Lot No. 0R8**). The kit design is based on IMGT/HLA database 3.50.0.

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
35	Added	-	5'-primer added for the B*45:29 allele.
56	Added	-	5'-primer added for improved resolution of the B*47 alleles.



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

The **HLA-DR low resolution** specificity and interpretation tables have been updated for the HLA-DRB alleles described since the previous *Olerup* SSP[®] HLA-A-B-DR SSP Combi Tray lot was made (**Lot No. 0R8**). The kit design is based on IMGT/HLA database 3.50.0.

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
73	Added, removed		5'-primer added for the DRB1*01:131N and DRB1*01:134 alleles. Excess 5'-primer removed.
93	-	Added	3'-primer added for the DRB3*02:131 alleles.
95	-	Added	3'-primer added for DRB5*02:33 allele.



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

Well **96** contains Negative Control primer pairs, that will amplify the majority of the Olerup SSP[®] HLA Class I, DRB, DQB1, DPB1 and DQA1 amplicons as well as all the amplicons generated by the control primer pairs matching the human growth hormone gene.

HLA-specific PCR product sizes range from 75 to 200 base pairs.

The PCR product generated by the positive control primer pair is 200 base pairs.

Length of PCR product	105	200	105	80	75	80	85
5'-primer¹	164	340	440	45	45	43	36
	5'-CAC ^{3'}	5'-Agg ^{3'}	5'-TTA ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-TAC ^{3'}
							36
							5'-TAT ^{3'}
3'-primer²	231	2nd I	507	59	58	57	47
	5'-TgC ^{3'}	5'-AAA ^{3'}	5'-TTg ^{3'}	5'-CTC ^{3'}	5'-ggC ^{3'}	5'-CTC ^{3'}	5'-ACA ^{3'}
							48
							5'-gCA ^{3'}
							48
							5'-gCC ^{3'}
							52
							5'-TgT ^{3'}
A*	+	+	+				
B*	+	+	+				
C*	+	+	+				
DRB1				+	+		
DRB3				+	+		
DRB5				+			
DQB1					+		
DPB1						+	
DQA1							+

¹The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

²The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

PRODUCT DESCRIPTION

HLA-A-B-DR SSP Combi Tray

CONTENT

The primer set contains 5'- and 3'-primers for grouping the HLA-A*01:01 to A*80:09N alleles into the corresponding serological groups A1 to A80.

The primer set contains 5'- and 3'-primers for grouping the B*07:02 to B*83:01 alleles into the corresponding serological groups B7 to B81 as well as primer pairs for recognizing the Bw4 and Bw6 sequence motifs.

The primer set contains 5'- and 3'-primers for grouping the DRB1*01:01 to DRB1*10:45 alleles into the corresponding serological groups DR1 to DR18 as well as primer pairs for recognizing the DRB3, DRB4 and DRB5 groups of alleles.

PLATE LAYOUT

Each test consists of 96 PCR reactions in a 96 well PCR plate.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48
49	50	51	52	53	54	55	56
57	58	59	60	61	62	63	64
65	66	67	68	69	70	71	72
73	74	75	76	77	78	79	80
81	82	83	84	85	86	87	88
89	90	91	92	93	94	95	NC

The 96 well PCR plate is marked with 'A-B-DR' in silver/gray ink.

Well No. 1 is marked with the Lot No. '6R6'.

Wells 1 to 24 – HLA-A low resolution primers.

Wells 25 to 72 – HLA-B low resolution primers.

Wells 73 to 95 – HLA-DR low resolution primers.

Well 96 – Negative Control.

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays and should be disregarded.

The PCR plates are covered with a PCR-compatible foil.

INTERPRETATION

Only HLA-A alleles will be amplified by the 24 wells of the HLA-A low resolution primer set, **wells 1 to 24**, except that a few HLA-B and HLA-C alleles will be amplified by some primer mixes. For further details see Specificity Table.



0197

For *In Vitro* Diagnostic Use

MA123 v02 SSP PI Template

Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

Only HLA-B alleles will be amplified by the 48 wells of the HLA-B low resolution, primer set, **wells 25 to 72**, except that a few HLA-A and HLA-C alleles will be amplified by some primer mixes. For further details see Specificity Table.

Only HLA-DRB alleles will be amplified by the 23 wells of the DR low resolution primer set, **wells 73 to 95**. Thus, the interpretation of DR low resolution typings is not influenced by other HLA class II genes. For further details see Specificity Table.

UNIQUELY IDENTIFIED ALLELES

HLA-A

All alleles, i.e. **A*01:01 to A*80:09N**, recognized by the HLA Nomenclature Committee in October 2022^{1,2} will be amplified by the primers in the HLA-A low resolution primer set, **wells 1 to 24**. The HLA-A alleles will be grouped into their corresponding serological specificities³.

The HLA-A alleles will be grouped into their low-resolution groups, except that the following alleles give rise to identical amplification patterns:

Alleles	Alleles
A*01:26, 01:136, 01:192, 01:380:01:01-01:380:01:02, A*11:94, 11:112, 11:211, 11:290, 11:326	A*23:66, 23:99, A*24:14:01:01-24:15, 24:51-24:53, 24:57, 24:64, 24:94, 24:114, 24:138, 24:188, 24:222N, 24:228, 24:291, 24:296, 24:304, 24:316, 24:324, 24:412, 24:481, 24:515, 24:574, B*53:72, C*04:01:03
A*01:28, 01:229, 01:299, A*11:282	A*25:78, A*26:57, A*66:15
A*01:43, A*11:271	A*26:92, A*66:42
A*03:215, 03:225, 03:249, 03:322, 03:414, A*25:77	A*30:01:01:01-30:04:04, 30:06-30:07, 30:09:01:01-30:12:02, 30:14L-30:15, 30:17-30:20, 30:23-30:43, 30:45, 30:47-30:54, 30:56-30:88, 30:90-30:124, 30:126-30:139, 30:141-30:149, 30:151-30:162, 30:164-30:176, 30:178N-30:185, 30:187-30:204N, B*07:260, C*12:328
A*11:11, A*34:19	A*31:08, A*33:53
A*11:116, 11:140, 11:199:01-11:199:02, 11:222, A*66:23	A*31:109, A*33:125, 33:131
A*11:325, A*24:555	A*32:03:01:01-32:03:01:02, 32:38, 32:60, 32:71, 32:78, 32:80, 32:107, B*15:642
A*23:14:01-23:14:02, 23:104-23:105, A*24:24, 24:71, 24:315, 24:392, 24:527, 24:569	A*34:01:01:01-34:01:07, 34:05-34:06, 34:11-34:12, 34:14, 34:16-34:18, 34:23, 34:27, A*66:40

HLA-B

All alleles, i.e. **B*07:02 to B*83:01**, recognized by the HLA Nomenclature Committee in October 2022^{1,2} will be amplified by the primers in the HLA-B low resolution primer set, **wells 25 to 72**. Most HLA-B alleles will be also grouped into their corresponding serological specificities.

The HLA-B alleles will be grouped into the groups B*07:xx to B*83:xx, except that the following alleles give rise to identical amplification patterns:



0197

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

Alleles	Alleles
B*07:349, B*42:05:01-42:05:02 B*08:26:01-08:26:02, 08:50, 08:62, 08:85, 08:94, 08:242, 08:263, 08:265, 08:272, 08:299, B*42:07, 42:24, 42:26	B*40:496, B*41:62 B*53:30, B*58:111
B*13:35, 13:59, 13:71, 13:108, 13:140, B*44:135, 44:158, 44:184	B*53:39, B*57:128, B*58:100
B*14:08:01-14:08:02, 14:55, B*39:01:35, 39:30, 39:32-39:33, 39:47, 39:50, 39:74, 39:82, 39:102, 39:107, 39:112, 39:128, 39:136, 39:171, 39:174	B*55:04, 55:08, 55:49, 55:86, B*56:01:05, 56:01:13, 56:15, 56:19N, 56:22, 56:88
B*18:29, 18:72:01-18:72:03, 18:92, 18:175, B*35:01:55, 35:01:62, 35:32:01-35:32:03, 35:37, 35:53N, 35:64:01-35:64:02, 35:68:01-35:68:02, 35:99, 35:118-35:119, 35:174, 35:369, 35:409, 35:433, 35:523	B*57:67:01, B*58:36

Notes on Bw assignment

The Bw4 and Bw6 status of HLA-B alleles is based on their serological type:

	Serological Type
Bw4	B5, B5102, B5103, B13, B17, B27, B37, B38(16), B44(12), B47, B49(21), B51(5), B52(5), B53, B57(17), B58(17), B59, B63(15), B77(15)
Bw6	B7, B703, B8, B14, B18, B22, B2708, B35, B39(16), B3901, B3902, B40, B4005, B41, B42, B45(12), B46, B48, B50(21), B54(22), B55(22), B56(22), B60(40), B61(40), B62(15), B64(14), B65(14), B67, B70, B71(70), B72(70), B73, B75(15), B76(15), B78, B81, B82

However, according to the HLA Dictionary³ there are several exceptions to these associations. This product contains reactions specific to the sequences that code for the Bw4 and Bw6 epitopes. A small number of alleles classified as Bw4 or Bw6 in the HLA Dictionary are not amplified neither by the Bw4 nor the Bw6 primer mix:

Bw Reactivity	Allele
Bw6	B*15:377Q, B*35:01:39, B*35:02:09, B*39:01:15

HLA-DRB

All the HLA-DRB1, -DRB3, -DRB4 and -DRB5 alleles, i.e. **DRB1*01:01:01 to DRB1*10:45, DRB3*01:01 to DRB3*03:60, DRB4*01:01 to DRB4*03:01N, DRB5*01:01 to DRB5*02:36**, recognized by the HLA Nomenclature Committee in October 2022^{1,2} will be amplified by the primers in the DR low resolution primer set, **wells 73 to 95**. The HLA-DRB alleles will be grouped into their corresponding serological specificities³.

The HLA-DRB1, -DRB3, -DRB4 and -DRB5 alleles will be grouped into their low-resolution groups, except that the following alleles give rise to identical amplification patterns:



0197

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

Alleles	Alleles
DRB1*03:126, DRB1*13:193	DRB1*11:293, DRB1*13:45
DRB1*08:09:02, DRB1*14:15	DRB1*12:57, DRB3*02:111
DRB1*08:20, DRB1*13:18, 13:47, 13:55, 13:158:01-13:158:02, 13:164, 13:232, 13:286	DRB1*13:13, 13:119, 13:154, 13:156, DRB1*14:84, 14:116, 14:144
DRB1*08:31, DRB1*11:05, 11:209	DRB1*13:44, 13:86, 13:206, DRB1*14:183
DRB1*08:41, 08:75, DRB1*11:67, 11:193:01-11:193:02	DRB1*13:235, DRB1*14:98
DRB1*11:264, DRB1*13:05:01:01-13:05:03, 13:07:01-13:07:02, 13:11:01-13:11:02, 13:14:01-13:14:03, 13:21:01:01-13:21:02, 13:42:01-13:42:02, 13:46, 13:49:01-13:50:03, 13:62, 13:100, 13:108, 13:132, 13:136, 13:150, 13:169, 13:189, 13:199, 13:203, 13:229, 13:246, 13:284, 13:302	DRB1*13:314, DRB1*14:37

¹HLA-A, HLA-B and HLA-DRB alleles listed on the IMGT/HLA web page 2022-October-12, release 3.50.0, www.ebi.ac.uk/imgt/hla.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

³Holdsworth, R Hurley, C K Marsh, S G E Lau, M Noreen, H J Kempenich, J H Setterholm, M Maiers, M: 'The HLA dictionary 2008 : a summary of HLA-A , -B , -C , -DRB1 / 3 / 4 / 5 , and -DQB1 alleles and their association with serologically defined HLA-A, -B, -C, -DR and -DQ antigens', Tissue Antigens, 2009, 73(2), pp. 95-170.



0197

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information
SPECIFICITY TABLE

HLA-A low resolution primer set

Specificities and sizes of the PCR products of the 24 primer mixes used for
HLA-A low resolution SSP typing

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	HLA-A serology ³	Amplified HLA-A alleles ⁴
1⁵	120 bp, 145 bp, 225 bp	800 bp	A1, A36	*01:01:01:01-01:04:01:02N, 01:06-01:33, 01:35-01:424, 03:18, 03:135, 03:395, 11:94, 11:112, 11:211, 11:226, 11:271, 11:282, 11:290, 11:326, 25:59, 26:120, 31:200, 32:156, 36:01:01:01-36:12
2¹⁰	135 bp, 175 bp, 215 bp, 255 bp, 365 bp, 415 bp	800 bp	A2, A19, A28, A203, A210	*02:01:01:01-02:01:07, 02:01:09-02:01:15, 02:01:17-02:01:19, 02:01:21-02:01:81, 02:01:83-02:22:02, 02:24:01-02:35:01, 02:35:03-02:47, 02:49:01-02:77, 02:78 ^w , 02:79:01-02:97:02, 02:99:01-02:99:02, 02:101:01-02:128, 02:130-02:570:02, 02:572-02:643N, 02:645-02:828, 02:830:01-02:1079, 03:89:01, 30:13, 30:16, 30:22, 30:44, 30:46, B*35:01:60^w, B*58:01:24^w
3^{6,11}	205 bp, 235 bp	1070 bp	A1, A3, A11, A31, A32, A34, A36	*01:12, 01:19, 01:21, 01:126, 01:200, 01:244, 01:360, 03:01:01:01-03:17:02, 03:19-03:74, 03:76-03:94, 03:96-03:134, 03:136-03:176, 03:178N-03:186, 03:188-03:193, 03:195-03:199, 03:201-03:214, 03:216-03:224, 03:226-03:248, 03:250-03:259, 03:261-03:321, 03:323N-03:342N, 03:344-03:346, 03:348-03:376, 03:378-03:413, 03:415-03:453, 11:03:01:01-11:03:01:02, 11:20, 11:25:01-11:25:03, 11:60, 11:130, 11:158, 11:175, 11:183, 11:209, 11:348, 11:366, 11:381, 24:92, 24:458, 24:551, 24:557, 31:03 ^w -31:04:02 ^w , 31:176, 31:179, 32:04, 32:142, 33:49 ^w , 33:144 ^w , 33:179 ^w , 33:201 ^w , 34:02:01:01-34:04, 34:07-34:10N, 34:13, 34:15, 34:20-34:22, 34:24-34:26N, 34:28-34:30N, 36:02, 68:103:01-68:103:02, 74:23
4	190 bp	800 bp	A1, A3, A11, A30, A36, A68	*01:01:01:01-01:01:22, 01:01:24-01:01:47, 01:01:49-01:01:64, 01:01:67-01:01:78, 01:01:80-01:04:01:02N, 01:06-01:07, 01:09:01:01-01:11N, 01:13, 01:15N ^w , 01:16N-01:18N, 01:20-01:29, 01:31N-01:33, 01:35-01:78, 01:80-01:98, 01:100-01:144, 01:146, 01:148, 01:150-01:158, 01:160N-01:166, 01:168-01:177, 01:179N-01:199, 01:201-01:204, 01:206-01:207, 01:209-01:213, 01:215-01:227, 01:229-01:235, 01:237-01:243, 01:245-01:266:02, 01:268-01:270, 01:272-01:278, 01:280-01:313, 01:315-01:375, 01:377-01:390, 01:392-01:395, 01:397-01:402, 01:404-01:417N, 01:419-01:424, 02:78, 02:169, 02:829, 03:12, 03:18, 03:88, 03:135, 11:01:01:01-11:27, 11:29-11:52Q, 11:54:01-11:433N, 24:555, 25:43, 25:59, 26:19, 26:72, 26:184, 26:204:01:01-26:204:01:02, 29:67, 29:156,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				30:08:01-30:08:02, 30:150, 32:64, 34:19, 36:04, 66:23, 68:13:01, 68:66:01-68:66:02, 68:134, 68:178, 74:19, C*03:365, C*12:131
5	160 bp, 335 bp, 505 bp	1070 bp	A3, A9, A23, A24, A2403, A29, A31, A32	*01:301Q, 03:15, 03:19, 03:30, 03:152, 03:273, 03:347, 11:139, 11:325, 23:01:01:01-23:68, 23:70-23:126, 24:02:01:01-24:11N, 24:13:01-24:15, 24:17:01:01-24:64, 24:66-24:210, 24:212-24:470, 24:472-24:583, 25:81-25:82, 29:07, 29:49, 31:08, 31:29, 32:05, 32:79, 32:170, 33:19, 33:53, 33:187, 68:249, B*53:72, C*04:01:03
6	135 bp, 200 bp	800 bp	A9, A23, A24, A29, A80	*11:166, 23:01:01:01-23:56, 23:58-23:65, 23:67-23:68, 23:70-23:98, 23:100-23:112, 23:114-23:126, 24:24, 24:71, 24:315, 24:392, 24:527, 24:569, 29:07, 29:49, 31:29, 31:85, 32:72, 32:170, 33:187, 80:01:01:01-80:09N, B*18:27
7	175 bp, 205 bp	1070 bp	A2, A3, A9, A23, A24, A2403, A26	*01:301Q, 02:17:02:01 ^w -02:17:04 ^w , 02:804, 03:347, 11:139, 23:14:01-23:14:02, 23:104-23:105, 23:113N, 24:02:01:01-24:11N, 24:13:01-24:13:02, 24:17:01:01-24:50, 24:54-24:56:02, 24:58-24:63, 24:66-24:91, 24:93, 24:95-24:113, 24:115-24:137, 24:139-24:187, 24:189-24:210, 24:212-24:221, 24:223-24:227, 24:229-24:290, 24:292-24:295, 24:297-24:303N, 24:305-24:315, 24:317-24:323N, 24:325-24:411, 24:413-24:457, 24:459:01-24:470, 24:472-24:480, 24:482-24:514N, 24:516-24:550, 24:552-24:554, 24:556N, 24:558-24:573, 24:575-24:583, 26:16, 30:163, 33:19, 33:119, 68:45, 68:117
8	165 bp, 200 bp	800 bp	A2, A3, A10, A11, A25, A26, A28, A32, A34, A66, A68, A69	*01:51, 01:344, 01:381, 02:55, 02:527, 02:582, 02:644, 02:741, 02:815, 03:24, 03:50, 11:10, 11:183, 11:191, 11:269, 25:01:01:01-25:16, 25:18-25:33, 25:35-25:45, 25:47-25:53, 25:55-25:58, 25:60-25:70, 25:72-25:76, 25:78-25:80, 25:84, 26:01:01:01-26:01:47, 26:01:49-26:06, 26:08:01:01-26:15, 26:17-26:18, 26:20:01-26:43:02, 26:45-26:63, 26:65-26:71N, 26:73-26:88, 26:90-26:91, 26:93-26:176, 26:178-26:183, 26:185-26:206:02N, 26:208-26:224, 29:28, 32:15, 33:51, 34:01:01:01-34:18, 34:20-34:30N, 66:01:01:01-66:22, 66:24-66:41, 66:43-66:46, 68:01:01:01-68:177, 68:179-68:249, 68:251N-68:296, 69:01:01:01-69:09
9^{5,7}	80 bp	800 bp	A2, A3, A25, A32	*01:298, 02:81, 02:124, 02:829, 25:01:01:01-25:45, 25:47-25:58, 25:60-25:76, 25:79-25:80, 25:83-25:84, 32:01:01:01-32:02, 32:04, 32:06-32:37, 32:39-32:59, 32:61-32:70, 32:72-32:77, 32:81-32:106:01:02, 32:108-32:171, B*07:81, B*08:52, B*18:67, B*38:41, B*38:67, B*49:67, B*51:185, B*53:05:01-53:05:02, B*53:16, B*53:33, B*57:01:38, B*57:60
10^{5,7,9}	80 bp, 200 bp, 240 bp	800 bp	A10, A26, A43	*01:43, 01:51, 02:644, 11:17, 11:40, 11:223, 11:271, 25:59, 26:01:01:01-26:02:02, 26:04, 26:07:01-26:20:02, 26:22-26:29, 26:31-26:43:02, 26:45-26:56, 26:58-26:77, 26:79-26:91, 26:93-26:110, 26:112:01-26:176, 26:178-26:224, 29:105, 29:164, 33:13, 33:48, 43:01-43:03, 68:84, 68:227



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

Lot No.	Fragment Sizes	Reference Size	Alleles	SNPs
11 ⁵	80 bp, 175 bp, 500 bp	1070 bp	A1, A9, A10, A11, A26, A31, A34, A66	*01:13, 01:26, 01:28, 01:136, 01:176, 01:192, 01:229, 01:299, 01:380:01:01-01:380:01:02, 03:63, 03:88, 11:01:01:01-11:27, 11:29-11:52Q, 11:54:01- 11:270, 11:272Q-11:433N, 24:19, 24:44, 24:555, 26:03:01:01-26:03:01:02, 26:06, 26:21, 26:78, 26:92, 26:111, 26:152, 26:177, 26:193, 30:125, 30:140, 33:184, 34:01:01:01-34:08, 34:10N-34:30N, 66:01:01:01-66:01:07, 66:04-66:11, 66:13-66:14, 66:17-66:20, 66:22-66:24, 66:27N, 66:29-66:33, 66:35-66:38, 66:40-66:42, 66:44-66:46, 68:271, 69:02, 74:30, 80:02
12 ⁵	125 bp, 160 bp, 190 bp	800 bp	A1, A3, A10, A25, A26, A31, A34, A43, A66, A74	*01:01:56, 01:361N, 02:309, 02:454, 03:01:19, 03:103:02, 11:11, 25:05-25:06, 26:09, 26:54, 26:91, 26:217, 31:03-31:04:02, 31:123, 34:01:01:01- 34:30N, 43:01-43:03, 66:02:01:01-66:03:01:03, 66:16, 66:19, 66:21, 66:25-66:26Q, 66:28N, 66:34, 66:39N-66:40, 66:43, 68:130:01, 68:294, 74:01:03
13 ⁹	180 bp, 225 bp	1070 bp	A1, A2, A3, A10, A25, A26, A34, A43, A66	*01:13, 01:176, 01:194, 02:34-02:35:03, 02:56:01- 02:56:02, 02:62, 02:103, 02:135, 02:580, 02:741, 02:907, 03:01:01:01-03:01:22, 03:01:24-03:01:82, 03:01:84-03:01:93, 03:01:95-03:07:02, 03:09- 03:11N, 03:13-03:31, 03:33-03:35, 03:37-03:40, 03:42-03:56:01:02, 03:58, 03:60-03:71, 03:73-03:87, 03:90-03:106, 03:109-03:110, 03:112-03:141, 03:143-03:151, 03:153-03:171, 03:174-03:175, 03:177, 03:179-03:193, 03:195-03:197:02N, 03:201- 03:202, 03:204, 03:206-03:210, 03:212-03:218, 03:220-03:251, 03:253-03:259, 03:261-03:266N, 03:268-03:272, 03:274-03:329N, 03:331-03:381N, 03:383-03:397, 03:399-03:412, 03:414-03:428, 03:430-03:453, 11:116, 11:140, 11:199:01- 11:199:03, 11:222, 25:01:01:01-25:05, 25:07-25:35, 25:37-25:42N, 25:44-25:49N, 25:51, 25:53-25:58, 25:60-25:64, 25:66-25:84, 26:01:01:01-26:01:76, 26:02:01 ^w -26:02:02 ^w , 26:03:01:01-26:03:01:02, 26:05-26:08:03, 26:10-26:28, 26:29 ^w , 26:30- 26:43:02, 26:45-26:48, 26:49 ^w , 26:50-26:74, 26:76- 26:77, 26:79-26:90, 26:92-26:203, 26:206:01N- 26:224, 30:55, 30:186, 31:24-31:25, 32:26:01- 32:26:02, 32:152, 33:61, 33:190, 34:08, 43:01-43:03, 66:01:01:01-66:01:07, 66:04-66:09, 66:10 ^w , 66:11- 66:15, 66:17-66:20, 66:22-66:24, 66:27N, 66:29- 66:33, 66:35-66:38, 66:41-66:42, 66:44-66:46, 68:71, 69:07, 74:13, 74:37
14 ^{5,9}	80 bp, 115 bp, 200 bp, 240 bp, 460 bp	1070 bp	A26, A29, A31, A34, A36, A68	*02:237, 02:309, 02:454, 02:809, 02:927, 03:01:38, 03:95, 03:123:02, 03:171, 03:231:01, 03:343, 03:428, 11:130, 26:19, 26:22, 26:198, 29:01:01:01- 29:172, 31:03-31:04:02, 31:123, 32:42, 34:04, 36:02, 66:09, 68:19, 68:130:02
15 ^{5,10,11}	90 bp, 135 bp, 205 bp	1070 bp	A1, A30	*01:13, 01:28, 01:176, 01:194, 01:229, 01:299, 03:43, 03:82, 03:186, 03:395, 11:113, 11:162, 11:282, 30:01:01:01-30:04:04, 30:06-30:20, 30:22- 30:204N, 31:35, B*07:260, C*12:328



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

16	145 bp, 160 bp, 240 bp, 380 bp, 410 bp	1070 bp	A24, A31, A32	*02:237, 02:809, 02:927, 03:95, 03:343, 11:306, 23:11N, 23:113N, 29:14, 29:172, 31:01:02:01- 31:217, 32:05, 32:79, 32:128, 33:53, 33:125, 33:131, 74:13
17¹¹	140 bp, 180 bp, 240 bp, 260 bp	1070 bp	A1, A2, A3, A25, A32, A74	*01:07, 01:95, 01:289, 01:298, 03:32, 03:43, 03:82, 03:152, 03:186, 03:219, 03:273, 03:429, 23:64, 24:104, 24:243, 25:03, 25:30, 29:13, 29:39, 29:98, 30:89, 30:177, 31:21, 31:35, 32:01:01:01-32:171, 74:07, 74:33, 80:04, B*15:642
18	200 bp, 390 bp	1070 bp	A24, A31, A32, A33, A68, A74	*02:243:01-02:243:03, 02:929, 24:82, 26:177, 26:207, 29:48, 29:105, 29:164, 31:08, 31:109, 32:05, 32:15, 32:79, 32:93, 33:01:01:01-33:01:17, 33:03:01:01-33:37, 33:39-33:233, 68:29, 68:250, 74:04, 74:21, 74:43
19	370 bp	800 bp	A2, A19, A68, A74	*01:121, 02:65, 02:407, 02:449, 02:932, 03:246, 03:282, 11:07:02, 11:372, 68:25, 68:262, 74:01:01:01-74:43
20^{9,11,12}	210 bp, 240 bp	800 bp	A2, A25, A28, A68	*02:34-02:35:03, 02:46, 02:48, 02:56:01-02:56:02, 02:62, 02:70, 02:78, 02:103, 02:129, 02:571, 02:580, 02:651, 02:741, 02:907, 02:1050, 11:199:03, 25:05, 25:59, 26:54, 26:122, 68:01:01:01-68:217, 68:219- 68:257, 68:259-68:296
21^{5,12}	65 bp, 240 bp, 375 bp, 545 bp	800 bp	A2, A26, A28, A32, A66, A68, A69	*02:55, 02:149:01-02:149:02, 02:243:01-02:243:03, 02:309, 02:644, 02:741, 02:815, 02:842, 02:947, 23:57, 24:82, 25:30, 26:22, 26:177, 29:80, 29:155, 31:41, 32:06, 33:22, 66:06, 66:09, 68:08:02, 68:29, 68:105:01:01-68:105:01:02, 68:157, 68:250, 68:276Q, 69:01:01:01-69:09, 74:35
22^{5,8}	85 bp, 240 bp, 400 bp	800 bp	A2, A36	*02:34-02:35:03, 02:46, 02:48, 02:56:01-02:56:02, 02:62, 02:70, 02:78, 02:103, 02:129, 02:571, 02:576, 02:580, 02:682, 02:741, 02:907, 02:1050, 03:187, 11:155, 11:199:03, 11:226, 31:62, 34:22, 36:01:01:01-36:12, 68:11N, 68:41, 68:230, B*07:439, B*35:491, B*40:359, B*55:96, B*57:65, C*04:31, C*06:137, C*07:569, C*07:990, C*07:994
23^{5,9,10}	80 bp, 160 bp, 240 bp, 495 bp	800 bp	A2, A3, A24, A26, A28, A36, A68, A80	*02:55, 02:237, 02:671, 02:809, 02:927, 03:41, 03:63, 03:75, 03:88, 03:95, 03:177, 03:343, 03:347, 03:377, 03:382, 11:130, 24:18, 24:204, 24:213, 26:03:01:01-26:03:01:02, 26:05-26:06, 26:21, 26:30, 26:78, 26:111, 26:177, 33:24, 33:184, 36:02, 68:05, 68:15, 68:20, 68:109, 68:136, 68:255, 80:01:01:01- 80:01:02, 80:03-80:09N, C*06:187
24^{10,11}	360 bp, 445 bp	1070 bp	A10, A26, A31, A66	*02:135, 02:309, 02:454, 03:01:19, 03:103:02, 25:13, 25:46, 25:82-25:83, 26:07:01-26:07:02, 26:30, 26:65, 26:92, 31:04:01:01-31:04:02, 31:123, 34:09, 66:02:01:01-66:03:01:03, 66:12, 66:16, 66:21, 66:25- 66:26Q, 66:28N, 66:34, 66:39N, 66:42-66:43, 68:294, 74:01:03



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 6R6

Lot-specific information

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-A low resolution SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length. Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings. PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings. Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³The serological reactivity of all HLA-A alleles is not known. In this table we use the expert-assigned serological grouping in Tissue Antigens (2009) 73:95-170 and the serological grouping of the sequence-defined allele.

⁴For several HLA Class I alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁵HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁶Primer mix 3 may faintly amplify the A*30:04:01:01-30:04:04, 30:06, 30:09:01:01-30:09:01:02, 30:17, 30:29, 30:46, 30:77, 30:89-30:90, 30:99, 30:103, 30:105, 30:117, 30:139, 30:150, 30:155, 30:160 and 30:183 alleles.

⁷Primer mix 9 and 10 may weakly amplify the A*34 alleles.

⁸Primer mix 22 might faintly amplify most A*11 alleles.

⁹Primer mixes 10, 13, 14, 20 and 23 may have a tendency of giving rise to primer oligomer formation.

¹⁰Primer mixes 2, 15, 23 and 24 may give rise to a lower yield of HLA-specific PCR product than the other HLA-A low primer mixes.

¹¹Primer mixes 3, 15, 17, 20 and 24 may have tendencies of unspecific amplifications.

¹²Primer mixes 20 and 21 may generate a false positive band of about 560 base pairs. This band should be disregarded when interpreting HLA-A low resolution typings.

Abbreviations

w: might be weakly amplified.



0197

For *In Vitro* Diagnostic Use

MA123 v02 SSP PI Template

Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

SPECIFICITY TABLE

HLA-B low resolution primer set

Specificities and sizes of the PCR products of the 48 primer mixes used for HLA-B low resolution SSP typing

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	HLA-B serol. ³	Amplified HLA-B alleles ⁴
25^{5,9}	110 bp	800 bp	7, 703, 35, 40, 41, 42, 48, 61	*07:02:01:01-07:02:05, 07:02:07-07:02:84, 07:02:86-07:18:02, 07:20-07:32, 07:34-07:39, 07:41-07:47, 07:49N-07:52, 07:54-07:59, 07:61-07:99, 07:101-07:121, 07:123-07:138, 07:140-07:163, 07:166-07:208, 07:210-07:222, 07:224-07:225, 07:227-07:234, 07:236-07:300, 07:302-07:342, 07:344-07:345, 07:347-07:358, 07:360Q-07:362, 07:364-07:385, 07:387, 07:389-07:406, 07:408-07:410, 07:412-07:419, 07:421-07:441, 07:443-07:467Q, 07:469N, 08:79:01-08:79:02, 08:156, 08:187, 13:174, 15:138, 15:230, 15:241, 15:297, 35:66, 37:07, 40:15-40:16:01:02, 40:23, 40:32, 40:98, 40:128, 40:136, 40:158, 40:183, 40:198, 40:231, 40:279, 40:355, 40:466, 41:08, 42:05:01-42:05:02, 44:150, 46:87, 48:05, 48:08, 48:15, C*02:23, C*03:260, C*04:77
26	215 bp	1070 bp	8, 44	*08:01:01:01-08:05, 08:07-08:25, 08:27-08:49, 08:51-08:61, 08:63-08:64, 08:66-08:84, 08:86N-08:93, 08:95-08:122, 08:124-08:131, 08:133-08:225, 08:227, 08:229-08:238, 08:240-08:241, 08:243-08:262, 08:264, 08:266-08:271, 08:273-08:294, 08:296-08:298, 08:300Q-08:303, 13:93, 15:142, 15:180, 15:429, 35:218, 35:256, 37:62, 40:328, 44:49, 44:156, 49:25, 51:68, 51:176, 52:77, C*06:147
27⁹	140 bp, 235 bp	1070 bp	7, 8, 13, 15, 35, 4005, 44, 49, 61, 62, 77	*07:20, 07:24, 07:60, 07:100, 07:131, 07:220, 07:223, 07:358, 07:464, 08:21, 08:25, 08:175, 13:01:01:01-13:04, 13:06-13:08, 13:10-13:23, 13:25-13:38, 13:40-13:57, 13:59-13:66, 13:68-13:96, 13:99-13:174, 15:07:01:01-15:07:03, 15:36 ^w , 15:45, 15:55, 15:68, 15:89:01 ^w -15:89:02 ^w , 15:126, 15:197:02, 15:207, 15:324, 15:329, 15:331, 15:360, 15:383, 15:405, 15:431, 15:450, 15:524, 35:05:01:01-35:05:06, 35:16-35:17:05, 35:22:01:01-35:22:01:02, 35:30:01:01-35:30:01:02, 35:51:01-35:51:02, 35:58, 35:72, 35:89, 35:97, 35:113-35:114, 35:125, 35:164, 35:199, 35:232, 35:242, 35:318, 35:346, 35:368, 35:389, 35:416, 35:427N, 35:440, 35:459N, 35:472, 35:481, 35:498-35:499, 35:509, 35:519, 35:525, 35:541, 40:05:01:01-40:05:02, 40:71, 40:174, 40:496, 41:62, 44:08 ^w , 44:54, 44:57 ^w , 44:60, 44:106, 44:110, 44:135, 44:158, 44:184, 44:211, 44:213, 46:12, 51:64, 51:148, 53:14, 53:67, 55:51, 58:18, C*03:05:01:01-03:05:01:03, C*03:25, C*03:27, C*03:143, C*03:167, C*03:296:01-03:296:02, C*03:335, C*03:346, C*03:407, C*03:482, C*03:494, C*03:555



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 6R6

Lot-specific information

28 ^{5,6,7}	105 bp, 140 bp, 260 bp	800 bp	12, 13, 14, 15, 17, 18, 21, 27, 35, 39, 40, 41, 44, 45, 47, 48, 49, 50, 57, 58, 60, 61, 62, 64, 65	*07:68:01-07:68:03, 07:100, 07:218, 07:220, 07:237, 07:262:01-07:262:02, 07:424, 08:71, 08:118, 08:274, 08:296, 13:01:01:01-13:01:09, 13:01:11-13:04, 13:06-13:13:02, 13:15-13:23, 13:25-13:57, 13:59-13:82, 13:84-13:92, 13:94- 13:96, 13:99-13:101, 13:102 ^w , 13:103N-13:105, 13:107- 13:174, 14:01:01:01-14:02:08, 14:02:10-14:05, 14:07N, 14:09, 14:11-14:19, 14:21-14:36, 14:38-14:54, 14:57-14:95, 14:97-14:101N, 14:103-14:113N, 15:07:01:01-15:07:03, 15:45-15:46, 15:53, 15:55, 15:68, 15:106, 15:126, 15:197:02, 15:207, 15:212, 15:324, 15:329, 15:331, 15:336, 15:347, 15:360, 15:383, 15:399, 15:405, 15:431, 15:450, 15:510, 15:559, 15:588, 18:14, 18:44:01-18:44:02, 18:107, 27:32, 27:106, 27:138, 27:180, 27:224, 35:05:01:01-35:05:06, 35:22:01:01-35:22:01:02, 35:46-35:47, 35:51:01-35:51:02, 35:58, 35:63, 35:66, 35:72, 35:89, 35:97, 35:137, 35:147, 35:154, 35:182, 35:199, 35:207, 35:217, 35:232, 35:278, 35:317, 35:368, 35:389, 35:416, 35:427N, 35:440, 35:450, 35:459N, 35:472, 35:475, 35:481, 35:498-35:499, 35:509, 35:519, 35:540-35:541, 38:19, 39:03:01:01-39:03:02, 39:14:01:01-39:14:01:03, 39:24:01-39:24:03, 39:29, 39:37:01:01-39:37:01:02, 39:76, 39:79, 39:120, 39:144, 39:152, 39:182, 39:184-39:185, 40:01:01-40:01:39, 40:01:41-40:01:79, 40:07, 40:10:01:01-40:10:02, 40:14:01- 40:16:01:02, 40:22N-40:23, 40:25, 40:30-40:36, 40:38, 40:42-40:43, 40:45-40:49, 40:51-40:55, 40:58-40:63, 40:65- 40:67, 40:69, 40:71-40:73, 40:76-40:77, 40:79-40:81, 40:84, 40:87:01-40:88, 40:92, 40:100-40:102, 40:106, 40:108, 40:110, 40:112-40:114:02, 40:116-40:118N, 40:121, 40:123- 40:126, 40:128-40:130:02, 40:132, 40:134-40:136, 40:139- 40:141, 40:146-40:147, 40:150-40:156, 40:158, 40:160:01- 40:160:02, 40:163-40:164, 40:166, 40:168, 40:170-40:172, 40:174-40:175, 40:178-40:179, 40:182-40:183, 40:186:01- 40:188, 40:191-40:199, 40:204, 40:207-40:208, 40:212- 40:213, 40:215-40:219, 40:221-40:223, 40:227-40:228, 40:231, 40:233-40:242, 40:245:01-40:245:02, 40:247, 40:249-40:250, 40:252-40:253, 40:257-40:265N, 40:272- 40:273, 40:277-40:282, 40:285-40:286N, 40:288, 40:298:01- 40:301, 40:308, 40:310, 40:312, 40:315-40:316, 40:319, 40:321, 40:323:01-40:327, 40:329, 40:332-40:333, 40:336, 40:338N-40:339, 40:342, 40:344, 40:346, 40:352-40:353:02, 40:355, 40:357-40:358, 40:360, 40:367-40:368, 40:370, 40:372N, 40:377, 40:379:01:01-40:383, 40:385-40:386, 40:391-40:392, 40:394-40:398, 40:400-40:404, 40:406- 40:407, 40:416-40:422, 40:424, 40:426N, 40:428N, 40:431, 40:433, 40:436, 40:439, 40:443, 40:446-40:450, 40:452, 40:456, 40:461-40:463, 40:466, 40:469-40:470, 40:474- 40:475:02, 40:477, 40:479-40:480, 40:483N, 40:485-40:486, 40:490-40:491, 40:495-40:499, 40:501-40:503, 40:506N- 40:510N, 40:512-40:513, 40:516-40:517, 40:523-40:525, 41:01:01:01-41:58, 41:60-41:69, 41:71-41:79, 42:10, 42:17, 44:02:01:01-44:02:31, 44:02:33-44:02:39, 44:02:41- 44:02:56, 44:02:58-44:02:60, 44:02:61 ^w , 44:02:62-44:03:06, 44:03:08-44:03:18, 44:03:21-44:03:30, 44:03:32-44:03:59, 44:03:61-44:05:06, 44:07-44:48, 44:50:01-44:82,
---------------------	------------------------------	---------------	--	--



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				44:84:01-44:102, 44:104-44:129, 44:132-44:155, 44:157-44:193, 44:195N-44:211, 44:213-44:223, 44:226-44:230, 44:231 ^w , 44:232-44:245:02, 44:247-44:321, 44:322 ^w , 44:323-44:363, 44:437-44:444, 44:446-44:524, 44:526-44:527, 44:528 ^w , 44:529-44:546, 45:01:01:01-45:25, 45:27-45:30, 46:12, 46:84, 47:01:01:02-47:12, 48:12, 48:14, 48:29, 49:01:01:01-49:02, 49:04:01-49:14, 49:16-49:24, 49:26-49:30, 49:32-49:67, 49:69-49:76, 49:78, 50:01:01:01-50:02:01:02, 50:04:01-50:22, 50:31-50:33, 50:35-50:42, 50:44-50:70, 50:72N-50:79, 51:64, 51:203, 51:276, 52:106, 53:14, 53:30, 53:67, 54:38, 55:110, 56:67, 57:11, 57:45, 57:51, 57:69, 57:71, 57:88, 58:02:01:01-58:02:02, 58:06-58:07, 58:18, 58:25, 58:38, 58:43, 58:60, 58:104, 58:111, 58:132, C*01:30, C*01:32:01-01:32:02, C*06:20, C*06:207, C*12:87, C*14:76, C*14:82, C*16:85, C*16:98, C*16:102, C*16:110
29⁷	185 bp, 235 bp	800 bp	7, 16, 17, 27, 2708, 37, 38, 39, 3901, 3902, 64, 67	*07:02:32, 07:85:02, 07:197, 14:01:01:01-14:01:05, 14:01:07-14:01:13, 14:07N-14:08:02, 14:10, 14:12, 14:14, 14:19, 14:26, 14:32, 14:40, 14:46-14:47, 14:49, 14:53-14:55, 14:65, 14:71, 14:77-14:78, 14:80-14:83, 14:86-14:87, 14:95, 14:97, 14:102-14:105, 14:108Q, 14:112-14:113N, 18:01:37, 18:165, 27:01-27:05:15, 27:05:17-27:05:19, 27:05:21-27:05:37, 27:05:39-27:06:02, 27:07:04, 27:08-27:10, 27:12:01:01-27:13:02, 27:16-27:18, 27:20, 27:23, 27:26-27:27, 27:29, 27:31, 27:35-27:37, 27:39-27:42, 27:44-27:46, 27:48-27:61, 27:64N-27:69, 27:72-27:75, 27:77-27:80, 27:82-27:124, 27:126, 27:128-27:129, 27:131-27:137, 27:140-27:149, 27:151-27:152, 27:154-27:163, 27:165-27:179, 27:181-27:186, 27:188-27:212N, 27:214, 27:216-27:217, 27:220-27:232, 27:234-27:243N, 27:245-27:248, 27:250, 27:252-27:264, 37:02, 38:01:01:01-38:01:07, 38:01:09-38:01:18, 38:01:20-38:02:02:02, 38:02:04-38:03, 38:07-38:24, 38:26-38:32, 38:34N-38:83N, 38:84 ^w , 38:85-38:88, 38:90-38:102, 38:160-38:173Q, 39:01:01:01-39:01:01:29, 39:01:03:01-39:01:08, 39:01:10-39:01:25, 39:01:27-39:01:29, 39:01:31-39:02:01:02, 39:03:01:01-39:03:02, 39:05:01:01-39:06:02:08, 39:06:04-39:09:05, 39:11, 39:14:01:01-39:15, 39:18, 39:19:02, 39:22, 39:24:01-39:24:02, 39:25N-39:39:01, 39:40:01N-39:44, 39:46-39:48, 39:50-39:62, 39:64-39:71, 39:73-39:75, 39:77, 39:79-39:82, 39:83:02, 39:84-39:88, 39:90-39:91, 39:93-39:95N, 39:97N-39:98, 39:100, 39:102-39:104, 39:106-39:109, 39:111-39:125, 39:127-39:128, 39:130, 39:132, 39:134-39:142N, 39:144, 39:146N-39:161N, 39:163, 39:165-39:169, 39:171-39:174, 39:176-39:177, 39:178 ^w , 39:179, 39:182-39:183, 39:185, 39:187-39:193, 40:225:02, 40:429, 40:469, 47:04-47:05, 48:21, 48:26, 67:01:01, 67:03-67:04, 67:06-67:07
30	190 bp	800 bp	7, 14, 35, 38, 39, 44, 47, 57, 65	*07:28, 08:87, 14:02:01:01-14:02:02, 14:02:04-14:02:05, 14:02:07-14:06:02, 14:09, 14:11, 14:13, 14:15-14:18, 14:20, 14:22-14:23, 14:25, 14:27, 14:29-14:31, 14:33-14:39, 14:41N, 14:42:02-14:45, 14:48, 14:50-14:52, 14:56-14:64, 14:67-14:70Q, 14:72N-14:76N, 14:79N, 14:84-14:85N, 14:88-14:94, 14:96, 14:98-14:101N, 14:106-14:107, 14:109, 14:111, 15:77, 15:189:01-15:189:02, 15:233, 18:166, 27:227,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				35:01:64, 35:03:17, 35:26, 35:29:03, 38:05, 38:33, 39:04, 39:129, 40:414, 40:498, 44:16, 44:37:01-44:37:03, 44:64:01-44:64:02, 44:91, 44:132, 44:150, 44:190, 44:268, 44:275, 44:290, 44:347, 49:61, 51:01:68, 51:261, 51:274, 52:26, 57:04:01-57:04:03, 58:114, A*23:31, A*24:106, C*01:115, C*02:97, C*03:433, C*03:514, C*04:439, C*06:306, C*07:231, C*14:79, C*16:10
31	290 bp	1070 bp	15, 22, 46, 62, 63, 71, 72, 75, 76, 77	*15:01:01:01-15:01:03, 15:01:06-15:01:53, 15:01:55-15:02:12, 15:03:03, 15:04:01:01-15:08:02, 15:10:02, 15:11:01-15:11:03, 15:11:05-15:15:01:02, 15:17:01:01-15:17:07, 15:19-15:21:01:02, 15:24:01-15:28, 15:30:01:01-15:36, 15:38:01-15:40:02, 15:42-15:46, 15:48, 15:50, 15:55-15:58, 15:60, 15:63, 15:65-15:66, 15:70-15:71, 15:73:01:01-15:73:01:02, 15:75-15:79N, 15:81-15:89:02, 15:92, 15:94N, 15:96-15:97, 15:101-15:102, 15:104-15:107, 15:109-15:113, 15:116-15:118, 15:120-15:122, 15:125:01-15:126, 15:128-15:129, 15:135-15:150, 15:152, 15:154-15:155, 15:157, 15:159-15:160:01:02, 15:162-15:172, 15:174-15:175, 15:177-15:179:02, 15:181N-15:185, 15:187-15:196, 15:199, 15:201-15:209N, 15:211, 15:213-15:219, 15:223-15:225, 15:227-15:228, 15:230-15:234, 15:236-15:237, 15:239-15:241, 15:244-15:251:02, 15:256-15:262:02N, 15:264-15:265, 15:267-15:273, 15:276-15:280, 15:283-15:289, 15:291, 15:295-15:305, 15:308-15:310, 15:315-15:322, 15:324-15:328, 15:330-15:334, 15:339-15:346, 15:348-15:350, 15:352-15:353, 15:355-15:361, 15:363:01-15:368, 15:370-15:375N, 15:377Q-15:378, 15:381, 15:383-15:387, 15:390-15:396, 15:398-15:407, 15:409-15:413, 15:415-15:417, 15:419-15:425, 15:428-15:432, 15:437-15:438, 15:440-15:443, 15:446, 15:449-15:450, 15:452-15:454N, 15:456-15:474, 15:476-15:486, 15:488-15:497, 15:499-15:501, 15:504:02-15:505, 15:507-15:508, 15:510, 15:512-15:514, 15:516, 15:518-15:525, 15:528N-15:529, 15:531-15:538, 15:540N-15:545, 15:547, 15:549N, 15:551-15:552, 15:554-15:558, 15:561-15:562N, 15:564-15:569, 15:572, 15:574-15:581, 15:583-15:587, 15:589-15:601, 15:603-15:606, 15:608-15:609, 15:615, 15:618-15:619, 15:621-15:627, 15:629-15:633, 15:635-15:642, 15:644-15:647, 46:01:01:01-46:70, 46:72-46:76, 46:78-46:88, 46:90-46:98N, 51:367, C*01:73, C*03:02:23, C*05:01:28, C*06:02:93, C*07:02:30, C*08:16:02, C*14:02:27
32 ⁸	165 bp, 225 bp, 285 bp, 330 bp	1070 bp	5, 8, 12, 21, 22, 37, 41, 42, 44, 45, 48, 51, 56, 57, 60, 62, 70, 71, 72, 82	*07:349, 08:01:01:01-08:01:22, 08:01:24-08:01:62, 08:01:64-08:05, 08:08N-08:13, 08:15-08:27, 08:29-08:36, 08:38:01-08:48, 08:50-08:73, 08:75-08:76, 08:78, 08:80-08:86N, 08:88, 08:90-08:106, 08:109-08:110, 08:112-08:155, 08:157-08:163, 08:165-08:199, 08:201-08:226, 08:228-08:244, 08:246-08:250, 08:252N-08:254, 08:256-08:263, 08:265-08:289, 08:292-08:303, 13:46, 13:154, 15:03:01:01-15:03:13, 15:09:01-15:10:08, 15:18:01:01-15:18:10, 15:23, 15:29:01:01-15:29:01:02, 15:37, 15:46-15:47:02, 15:49, 15:51-15:54, 15:61-15:62, 15:64:01-15:64:02, 15:69, 15:72, 15:74, 15:80, 15:83:01:01-15:83:01:02, 15:90-15:91, 15:93, 15:98-15:99, 15:103, 15:106, 15:108, 15:114-15:115, 15:119, 15:123-15:124, 15:127, 15:131-15:134, 15:143, 15:151,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				15:153, 15:156, 15:158, 15:161, 15:173, 15:176, 15:186, 15:197:01-15:198, 15:200, 15:210, 15:212, 15:220:01:01-15:221, 15:226N, 15:229, 15:235, 15:238, 15:242:01-15:243, 15:252-15:253, 15:255, 15:259, 15:263, 15:266, 15:274-15:275:02, 15:281-15:282, 15:290, 15:292-15:294N, 15:306-15:307, 15:311-15:314, 15:323, 15:329, 15:335, 15:337-15:338, 15:351, 15:354, 15:369, 15:379-15:380N, 15:382, 15:388-15:389, 15:397, 15:414, 15:426-15:427, 15:433, 15:435, 15:439, 15:444-15:445, 15:447-15:448, 15:451, 15:455, 15:475, 15:487N, 15:498, 15:502-15:504:02, 15:509, 15:511, 15:517, 15:525-15:527, 15:530, 15:539, 15:548, 15:553, 15:560, 15:563, 15:570-15:571, 15:582, 15:588, 15:607, 15:610-15:612, 15:614:01:01-15:614:01:02, 15:616Q-15:617, 15:620, 15:628, 15:643, 27:195, 35:87, 37:12, 37:84, 38:30, 39:156, 40:12:01-40:12:02, 40:26, 40:40, 40:136, 40:231, 40:279, 40:516, 41:01:01:01-41:03:02, 41:05-41:09, 41:11-41:17, 41:19-41:48, 41:50, 41:53-41:61, 41:64-41:79, 42:01:01:01-42:02:02, 42:04-42:08, 42:10-42:15, 42:17-42:19, 42:21-42:32, 44:02:68, 44:14-44:15:01:02, 44:18:01:01-44:18:01:02, 44:20, 44:62, 44:100, 44:166, 44:213, 44:355, 45:01:01:01-45:25, 45:27-45:30, 51:08:01:01-51:08:06, 51:20, 51:36, 51:44N, 51:97, 51:141, 51:153:01-51:153:02, 51:263, 51:274-51:275, 51:342, 51:354, 51:362, 51:371, 52:19, 55:20, 55:56, 56:13, 57:09, 82:01:01:01-82:04, C*03:342, C*07:46, C*07:356, C*07:546, C*12:184, C*15:99, C*16:85
33^{7,8,10}	165 bp, 190 bp, 390 bp	1070 bp	5, 17, 21, 35, 51, 5102, 5103, 52, 56, 62, 63, 70, 78	*13:62, 15:01:02, 15:09:01-15:09:02, 15:16:01:01-15:17:07, 15:67, 15:95, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:222, 15:230, 15:243, 15:254, 15:268, 15:273, 15:277-15:278, 15:356, 15:361-15:362, 15:396, 15:403, 15:408, 15:411, 15:423-15:424, 15:435, 15:446, 15:462, 15:500, 15:516, 15:523, 15:532, 15:550, 15:555, 15:575N, 15:602-15:603, 15:613, 15:619, 15:642, 15:644, 18:205, 27:153, 35:01:10, 35:03:32, 35:04:02, 35:251:01, 40:26, 40:28, 44:03:40, 44:62, 44:459, 49:18:01, 50:01:21, 50:14, 51:01:01:01-51:01:03, 51:01:05-51:01:71, 51:01:73-51:02:03, 51:02:05-51:08:04, 51:08:06-51:09:03, 51:11N-51:12, 51:13:02-51:14, 51:16-51:23, 51:24:03-51:24:05, 51:26-51:33, 51:37-51:41N, 51:43-51:44N, 51:46, 51:48-51:61:02, 51:63:01, 51:64-51:80, 51:82-51:91, 51:92:02, 51:94-51:98N, 51:100-51:106:01, 51:107-51:130, 51:132-51:142, 51:145-51:147, 51:149N-51:156, 51:158:01-51:196, 51:198, 51:200-51:219, 51:221-51:224, 51:226-51:294, 51:296-51:299, 51:301-51:337, 51:339-51:351, 51:353-51:358, 51:360-51:375, 52:01:01:01-52:01:04, 52:01:06-52:01:14, 52:01:16-52:13, 52:14:02-52:24, 52:25:02-52:57, 52:59-52:66, 52:68-52:96N, 52:98-52:110N, 53:28, 56:05:01-56:06, 56:21, 57:01:32, 58:08:01-58:08:02, 78:01:01:02-78:10
34	180 bp	800 bp	13, 22, 55, 60	*07:78, 13:01:01:01-13:02:37, 13:07N-13:09, 13:11, 13:14-13:20, 13:22:01-13:23, 13:25, 13:27-13:34, 13:36-13:45, 13:47, 13:49N-13:50, 13:52, 13:54-13:58, 13:60-13:61, 13:63N-13:70, 13:73-13:78, 13:80-13:85, 13:87-13:107,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				13:109-13:119, 13:121-13:126, 13:128-13:130, 13:132-13:139N, 13:141-13:145, 13:147-13:159, 13:161N-13:163, 13:165-13:173, 40:48, 40:377, 45:10, 49:07, 50:31, 54:26, 55:09, 55:22, 55:24, 56:60, 59:09
35⁵	100 bp, 195 bp	1070 bp	8, 12, 27, 38, 39, 3902, 44, 45, 48, 60, 70, 71, 72, 82	*07:27, 07:50, 07:219, 07:453, 08:04:01-08:04:02, 08:17, 08:54:01-08:54:02, 08:110, 08:285, 13:97, 14:73, 15:03:01:01-15:03:09, 15:03:11, 15:03:13, 15:47:01-15:47:02, 15:54, 15:61-15:62, 15:64:01-15:64:02, 15:68-15:69, 15:91, 15:98, 15:103, 15:123, 15:127, 15:131-15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220:01:01-15:220:04, 15:235, 15:242:01-15:242:02, 15:253, 15:259, 15:266, 15:274, 15:281-15:282, 15:369, 15:379, 15:389, 15:397, 15:444, 15:448, 15:475, 15:502, 15:511, 15:517, 15:527, 15:539, 15:560, 15:571, 15:582, 15:588, 15:610, 15:614:01:01-15:614:01:02, 15:616Q-15:617, 15:628, 18:37:01-18:37:02, 18:125, 27:18, 27:29, 35:361, 35:545, 37:28, 38:03, 39:02:01:01-39:02:03, 39:08, 39:13:01-39:13:02, 39:23, 39:39:01-39:39:02, 39:49, 39:88, 39:101, 39:105, 39:110, 39:162, 39:164, 39:173, 39:184, 40:12:01-40:12:02, 40:46, 40:93, 40:430, 42:11, 44:10, 44:15:01:01-44:15:01:02, 44:18:01:01-44:18:01:02, 44:40, 44:44, 44:117, 44:130, 44:140, 44:156, 44:210:01-44:210:02, 44:251, 44:339-44:340, 44:355, 45:01:01:01-45:03, 45:05-45:07, 45:11-45:23, 45:25, 45:27-45:30, 48:01:01:01-48:05, 48:07-48:30, 48:32-48:40, 48:42-48:55, 49:20, 49:73, 50:02:01:01-50:02:01:02, 50:21, 50:56, 50:75, 52:16, 55:18, 55:122, 56:44, 56:74, 58:44, 82:01:01:01-82:03, C*04:345, C*05:102, C*05:184, C*05:262, C*08:93, C*15:25
36	280 bp	1070 bp	5, 7, 8, 13, 15, 16, 17, 18, 22, 35, 40, 42, 44, 46, 48, 49, 51, 53, 55, 56, 57, 58, 60, 61, 62, 63, 70, 71, 72, 75, 76, 77, 78	*07:09:01-07:09:02, 07:11, 07:17, 07:162, 07:235, 07:237, 07:255, 07:287, 08:28, 08:35, 08:37, 08:89, 08:107, 08:291, 13:04, 13:10, 13:26:01-13:26:02, 13:72, 13:127, 13:160, 14:56, 15:01:01:01-15:01:04, 15:01:06-15:02:11, 15:03:01:01-15:08:02, 15:11:01-15:16:05, 15:18:01:01-15:21:01:02, 15:23-15:29:01:02, 15:31:01:01-15:36, 15:38:01-15:40:02, 15:43-15:44, 15:46-15:47:02, 15:49-15:57, 15:60-15:62, 15:64:01-15:72, 15:74-15:76, 15:78:01:01-15:82, 15:84-15:85, 15:87-15:89:02, 15:91-15:98, 15:101-15:129, 15:131-15:132, 15:134-15:136, 15:138-15:149N, 15:151-15:161, 15:163-15:167, 15:169-15:176, 15:178-15:187, 15:189:01-15:195, 15:197:01-15:207, 15:209N-15:215, 15:217-15:223, 15:225-15:228, 15:231-15:232, 15:234-15:242:02, 15:244-15:247, 15:249-15:251:03, 15:254-15:268, 15:270-15:272N, 15:274-15:286, 15:288-15:302N, 15:304N-15:311, 15:314-15:322, 15:325-15:328, 15:330-15:333, 15:335-15:337, 15:340-15:341, 15:343-15:355, 15:357-15:360, 15:362-15:395, 15:397-15:399, 15:401-15:402, 15:404-15:410, 15:412-15:415, 15:417-15:419, 15:421-15:422, 15:425-15:426, 15:428-15:434, 15:436-15:438, 15:440-15:444, 15:447-15:454N, 15:456-15:459, 15:461, 15:463N-15:484, 15:487N, 15:489-15:499, 15:501-15:502, 15:504:01-15:505, 15:507-15:515, 15:517-15:522, 15:524-15:529, 15:531, 15:533-15:540N, 15:542-15:545, 15:547-15:551, 15:554, 15:556-15:562N, 15:564-15:574, 15:576-15:577, 15:579-15:584N,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

15:586-15:602, 15:604N-15:610, 15:612-15:618, 15:621-15:641, 15:643, 15:645, 15:647, 18:01:01:01-18:15, 18:17N-18:25, 18:27-18:40, 18:42-18:60, 18:62-18:121, 18:123-18:147, 18:149-18:164, 18:167, 18:169-18:188, 18:190, 18:192-18:204, 18:206-18:216, 18:218-18:226, 27:41, 27:107, 35:01:01:01-35:01:73, 35:05:01:01-35:05:06, 35:07-35:08:15, 35:10-35:11:03, 35:14:01-35:17:05, 35:19-35:21, 35:23-35:30:01:02, 35:32:01-35:32:03, 35:35, 35:37, 35:40N-35:43:01, 35:43:03-35:43:04, 35:45-35:54, 35:57-35:58, 35:61:01-35:64:02, 35:66-35:69, 35:71-35:72, 35:76-35:80, 35:82, 35:86, 35:89-35:94, 35:97, 35:99-35:105, 35:107-35:108:02, 35:110-35:126, 35:130N-35:135, 35:137-35:148, 35:158-35:159, 35:161, 35:165N-35:166, 35:168, 35:170-35:171, 35:173:01N-35:178, 35:180, 35:185-35:192, 35:194, 35:196-35:197, 35:200, 35:202-35:203, 35:206-35:210, 35:212, 35:215:01-35:219, 35:222, 35:224-35:229, 35:232, 35:238-35:241, 35:244-35:245, 35:247:01-35:250, 35:252-35:255, 35:259-35:260, 35:262-35:265, 35:268-35:269, 35:271-35:272, 35:275, 35:277, 35:280, 35:283-35:284, 35:286-35:289, 35:293-35:295, 35:297, 35:299-35:308, 35:310, 35:313-35:315, 35:317-35:319, 35:322, 35:324-35:325, 35:327-35:328, 35:330-35:334, 35:336-35:338, 35:340-35:343, 35:346-35:348, 35:351-35:353, 35:356, 35:359, 35:363, 35:365, 35:368-35:370, 35:373, 35:376, 35:380:01:01-35:380:01:02, 35:382-35:383, 35:385-35:386, 35:388-35:390N, 35:392-35:393, 35:400-35:402, 35:404-35:406, 35:408-35:409, 35:411-35:412, 35:414, 35:416-35:418, 35:421, 35:426-35:430N, 35:433, 35:437, 35:439:01:01-35:439:01:02, 35:441, 35:444-35:445, 35:447-35:449, 35:453N-35:454, 35:456N, 35:458-35:459N, 35:461N-35:462, 35:464-35:467, 35:470, 35:472-35:473, 35:475-35:476, 35:479-35:481, 35:484:01-35:485, 35:488, 35:490-35:491, 35:493-35:495, 35:497-35:500, 35:502-35:504, 35:506, 35:509, 35:512-35:520, 35:523-35:526, 35:528-35:530, 35:532-35:533, 35:536-35:538, 35:541-35:542N, 35:548-35:549, 35:551N-35:552, 35:556Q, 35:558-35:560, 35:564, 35:566, 37:35, 37:56, 37:94, 38:52, 38:88, 39:07, 39:43, 39:91, 39:134, 39:169, 40:03:01:01-40:03:02, 40:20:01:01-40:20:01:02, 40:38, 40:52, 40:59-40:60, 40:105, 40:158-40:159, 40:222, 40:267, 40:284, 40:352, 40:422, 41:24, 41:30, 41:51, 42:09, 42:32, 44:17, 44:43:01-44:43:02, 44:144, 44:230, 44:256, 44:276, 44:308, 45:09, 46:01:01:01-46:10, 46:12-46:17, 46:19-46:60, 46:62-46:69, 46:71-46:73, 46:75-46:80, 46:82-46:98N, 48:02:01-48:02:03, 48:14, 48:23, 48:25, 49:04:01-49:05:02, 50:33, 51:37, 51:45, 51:63:01-51:63:02, 51:97, 51:216, 51:220, 51:251, 52:39, 52:97, 53:01:01:01-53:03, 53:05:01-53:06, 53:08:01-53:18, 53:20-53:30, 53:32, 53:34-53:35, 53:37-53:41, 53:43:01-53:61, 53:63-53:72, 54:06, 54:09, 54:14, 55:14, 55:23, 55:32, 55:68, 56:03, 56:09, 56:18, 56:31-56:32, 56:50, 57:01:01:01-57:01:32, 57:01:34-57:01:46, 57:06, 57:08, 57:10-57:11, 57:13-57:16, 57:18-57:27, 57:29-57:31:02, 57:33-57:38, 57:40-57:41, 57:43-57:45, 57:47-57:56, 57:58-57:62, 57:64-57:65, 57:67:01-57:69, 57:71-57:75, 57:77-57:79N, 57:81,



0197

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				57:83, 57:85-57:93, 57:95, 57:97-57:100, 57:102-57:117, 57:119-57:123, 57:125-57:136, 57:138-57:147, 57:149-57:154, 57:156-57:164, 58:01:01:01-58:02:02, 58:04-58:07, 58:09-58:19, 58:21-58:26, 58:29, 58:31N-58:33, 58:35-58:64, 58:66-58:72N, 58:74-58:113, 58:115-58:123, 58:125-58:139
37	195 bp	1070 bp	27, 35, 37, 44, 58, 76	*07:27, 07:236, 07:240, 08:126, 08:131, 08:139, 15:12:01-15:12:03, 15:14, 15:19, 15:91, 15:131, 15:161, 15:270, 15:298, 15:304N, 15:518, 15:551, 15:592, 15:634, 18:54, 18:56, 18:108, 18:132, 18:169, 18:195, 35:45, 35:71, 37:01:01:01-37:01:11, 37:01:13-37:09, 37:12-37:13, 37:15-37:21, 37:23-37:28, 37:30N-37:33N, 37:35-37:36, 37:38-37:66, 37:68-37:94, 37:96-37:99, 37:101-37:107, 38:17, 44:17, 44:43:01-44:43:02, 44:144, 44:230, 44:256, 44:308, 45:09, 46:17, 53:22, 57:93, 57:119, 58:07, A*24:174, C*02:82, C*06:165, C*07:979, C*14:61
38^{5,8}	105 bp, 395 bp, 435 bp	800 bp	5, 7, 15, 41, 42, 62, 63, 70, 71, 75, 77	*07:04:01-07:04:02, 07:25, 07:146, 07:349, 08:187, 15:09:01-15:10:08, 15:13:01-15:13:03, 15:16:01:01-15:18:10, 15:21:01:01-15:21:01:02, 15:23-15:24:02, 15:37, 15:44, 15:51-15:52, 15:66-15:67, 15:72, 15:76, 15:80, 15:87, 15:90, 15:93, 15:95, 15:99, 15:101, 15:108, 15:114-15:115, 15:119, 15:124, 15:133-15:134, 15:153, 15:157, 15:161-15:162, 15:168, 15:176-15:177, 15:186, 15:189:01-15:189:02, 15:196-15:198, 15:200, 15:208, 15:216, 15:221-15:222, 15:226N, 15:229-15:230, 15:238, 15:243, 15:252, 15:254-15:255, 15:263, 15:268, 15:273, 15:275:01-15:275:02, 15:290, 15:292-15:294N, 15:306-15:307, 15:311-15:314, 15:323, 15:329, 15:335, 15:337-15:338, 15:351, 15:354, 15:356, 15:361-15:362, 15:376, 15:380N, 15:382, 15:388, 15:394, 15:396, 15:403-15:404, 15:408, 15:411, 15:414, 15:418, 15:420, 15:423-15:424, 15:426-15:427, 15:433-15:436, 15:439, 15:442, 15:445-15:447, 15:451, 15:455, 15:462, 15:487N, 15:498, 15:500, 15:503-15:504:02, 15:509, 15:516, 15:523, 15:526, 15:530, 15:532, 15:548, 15:550, 15:555, 15:563, 15:570, 15:573, 15:575N, 15:586, 15:602-15:603, 15:607, 15:611-15:613, 15:615, 15:619-15:620, 15:623, 15:642-15:644, 15:646, 40:136, 40:231, 40:279, 41:08, 42:05:01-42:05:02, 51:367, A*25:01:20, C*01:240, C*03:542, C*07:335, C*14:62
39⁵	115 bp, 150 bp	1070 bp	18	*08:176, 15:545, 18:01:01:01-18:15, 18:17N-18:28, 18:30-18:71, 18:73-18:91, 18:93-18:101, 18:103-18:106Q, 18:108-18:174, 18:176-18:226, 35:289, 35:487, 44:445, 48:28
40⁵	80 bp	1070 bp	7, 27, 2708, 44, 60	*07:73, 08:88, 08:102, 15:588, 27:01-27:05:22, 27:05:24-27:21:02, 27:23-27:51, 27:53-27:66N, 27:68-27:74, 27:76-27:139, 27:141-27:188, 27:190-27:197, 27:199-27:214, 27:216-27:235, 27:237-27:238, 27:240-27:264, 38:22, 40:46, 40:93, 44:40, 44:44, 44:130, 44:210:01-44:210:02, 49:26, 49:73, 50:75
41⁵	100 bp, 150 bp	800 bp	12, 13, 17, 18, 22, 27, 35, 37, 39, 44, 45, 47, 48, 51, 53, 56, 57, 58,	*07:137, 07:235, 07:359, 08:49, 08:60, 08:76, 08:129, 08:181, 08:255, 13:01:01:01-13:01:21, 13:06-13:07N, 13:12:01-13:13:02, 13:17, 13:20-13:23, 13:25-13:26:02, 13:28-13:29, 13:36, 13:39, 13:43, 13:50-13:52, 13:57, 13:60-13:61, 13:63N, 13:73, 13:76N-13:80, 13:83, 13:86-13:87, 13:90, 13:92, 13:95, 13:106-13:107, 13:109-13:110, 13:122,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

60, 61, 62, 13:129, 13:135, 13:138, 13:141-13:142, 13:144, 13:146,
70, 75, 77, 13:148, 13:161N-13:163, 13:168, 13:170-13:171, 13:173,
78 14:10, 14:102, 15:02:01:01-15:02:12, 15:13:01-15:13:03,
15:20-15:21:01:02, 15:25:01:01-15:25:04, 15:36, 15:44,
15:62, 15:77, 15:80, 15:85, 15:88-15:89:02, 15:106, 15:112,
15:121, 15:139, 15:144, 15:154, 15:165, 15:170, 15:194,
15:204, 15:213-15:214, 15:223, 15:240, 15:250, 15:265,
15:271, 15:283, 15:289, 15:291, 15:297, 15:301-15:302N,
15:308, 15:319, 15:325, 15:328, 15:330, 15:341, 15:345,
15:357-15:358, 15:366, 15:374, 15:378, 15:393-15:394,
15:399, 15:402, 15:404, 15:407, 15:418, 15:420, 15:425,
15:437-15:438, 15:454N, 15:458, 15:460, 15:463N, 15:468,
15:491, 15:495, 15:505, 15:537, 15:543-15:544N, 15:557,
15:564, 15:566-15:567, 15:569, 15:574, 15:584N, 15:587,
15:589, 15:592, 15:594, 15:596N, 15:599, 15:615, 15:623,
15:625, 15:627, 15:630, 15:645-15:646, 18:22, 18:69,
18:105, 18:149-18:150, 18:204, 27:19:01:01-27:19:01:02,
27:30, 27:127, 27:150, 27:188, 35:01:01:01-35:01:54,
35:01:56-35:01:61, 35:01:63-35:04:03, 35:06-35:17:04,
35:18-35:21, 35:23-35:30:01:02, 35:33-35:36, 35:38-
35:42:02, 35:45-35:50, 35:52, 35:54-35:57, 35:59:01-
35:59:02, 35:61:01-35:63, 35:65Q, 35:69-35:71, 35:74,
35:76-35:78, 35:80-35:85, 35:88, 35:90-35:96, 35:98,
35:100-35:101:02, 35:103-35:113, 35:115-35:116, 35:120-
35:126, 35:128-35:134N, 35:136-35:150:02, 35:152-
35:173:02N, 35:175-35:184, 35:186-35:198, 35:200-35:204,
35:206-35:231, 35:233, 35:236-35:250, 35:252-35:259,
35:261-35:264, 35:266-35:272, 35:274-35:275, 35:277-
35:291, 35:293-35:330, 35:332-35:366, 35:370-35:388,
35:390N-35:408, 35:410-35:415, 35:417-35:426, 35:428Q-
35:432, 35:434-35:439:01:02, 35:441-35:458, 35:460Q-
35:470, 35:473-35:474, 35:477-35:480, 35:482-35:497,
35:500-35:508N, 35:510-35:518, 35:520, 35:522, 35:524-
35:533, 35:536-35:540, 35:542N-35:548, 35:550-35:568,
37:01:01:01-37:01:07, 37:01:09, 37:01:11-37:01:29, 37:03N-
37:06:02, 37:08, 37:10-37:43, 37:45-37:50, 37:52-37:80,
37:82N-37:107, 38:20:01-38:20:02, 39:42, 40:04:01:01-
40:04:02, 40:28, 40:30, 40:34, 40:64:01:01-40:64:02, 40:68,
40:99, 40:129, 40:137, 40:160:01-40:160:02, 40:343, 40:458,
41:48, 42:13, 44:02:01:01-44:02:33, 44:02:35-44:02:57,
44:02:59-44:14, 44:16-44:17, 44:19N, 44:21-44:30, 44:32-
44:40, 44:42-44:46:02, 44:48-44:52N, 44:55-44:64:02, 44:66-
44:98, 44:101-44:105, 44:107-44:134, 44:136-44:137,
44:139-44:157, 44:159-44:165, 44:167-44:181, 44:183,
44:185-44:191, 44:194:01-44:196, 44:198N-44:212, 44:214-
44:229, 44:231-44:235, 44:237N-44:245:02, 44:247-44:304,
44:306N-44:338, 44:340-44:342, 44:344-44:354N, 44:356-
44:363, 44:437-44:439, 44:441-44:458, 44:460-44:490,
44:492-44:505, 44:507-44:517, 44:519-44:524, 44:526-
44:546, 46:33, 46:67, 48:02:01-48:02:03, 48:25, 49:61,
50:17, 51:04:01-51:04:02, 51:42, 51:46, 51:56:01-51:56:03,
51:139, 51:186, 51:227, 51:281, 51:310, 51:359, 52:88,
53:01:01:01-53:13, 53:15-53:43:02, 53:45-53:46, 53:48N-
53:66, 53:68-53:72, 55:14, 56:09, 56:11-56:12, 56:60, 56:80,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				57:01:01:01-57:01:04, 57:01:06-57:10, 57:12, 57:14:01-57:20, 57:22-57:30, 57:32-57:46, 57:48-57:81, 57:84-57:117, 57:119-57:144, 57:146-57:164, 58:01:01:01-58:01:02, 58:01:04-58:01:15, 58:01:17-58:01:42, 58:04-58:05, 58:09-58:15, 58:17N, 58:19, 58:21-58:24, 58:28:01-58:29, 58:31N-58:37, 58:39N-58:42, 58:44-58:45:02, 58:47-58:59:02, 58:62-58:103, 58:105-58:120, 58:123-58:126, 58:128N-58:131, 58:133N-58:139, 83:01, C*03:102, C*03:263:01-03:263:03, C*03:317, C*03:382, C*03:514, C*15:24, C*15:126
42^{7,8}	140 bp	1070 bp	8, 18, 22, 35, 39, 78	*07:65, 07:134, 07:400, 08:32, 08:133, 08:180, 08:239, 15:202, 15:239, 15:308, 15:376, 18:01:01:01-18:03:01:03, 18:04:01-18:08, 18:10-18:15, 18:17N-18:36, 18:38-18:47, 18:49-18:53, 18:55-18:65, 18:68-18:106Q, 18:108-18:124, 18:126-18:135:01:02, 18:137-18:177, 18:179-18:193, 18:196-18:226, 35:01:01:01-35:02:08, 35:02:10-35:18, 35:20:01-35:24:02, 35:27-35:43:01, 35:43:03-35:45, 35:48, 35:50-35:62, 35:64:01-35:72, 35:74-35:153, 35:155-35:197, 35:199-35:206, 35:208:01-35:209, 35:211-35:216N, 35:218-35:253, 35:255-35:270, 35:272-35:273, 35:275-35:310, 35:312-35:316, 35:318-35:360, 35:362-35:449, 35:451-35:470, 35:472-35:568, 37:05, 37:11, 37:14, 37:29, 37:37, 37:67, 39:19:01-39:19:02, 39:186, 40:185, 40:209, 40:308, 40:350, 40:365, 45:22, 46:40:01-46:40:02, 55:130, 56:06, 56:45, 58:64, 78:01:01:02-78:10, C*07:1011
43^{5,9}	55 bp, 245 bp, 390 bp, 415 bp	1070 bp	5, 15, 48, 60, 62, 63, 70, 71, 72, 75, 76, 77, 82	*15:01:01:01-15:01:04, 15:01:06-15:21:01:02, 15:23-15:40:02, 15:42-15:58, 15:60-15:67, 15:69-15:99, 15:102-15:129, 15:131-15:179:02, 15:181N-15:235, 15:237-15:254, 15:256-15:505, 15:507-15:514, 15:516-15:545, 15:547-15:593, 15:595-15:633, 15:635-15:647, 40:12:01-40:12:02, 40:26, 40:40, 40:516, 44:02:68, 46:50, 51:367, 82:01:01:01-82:03, C*03:179, C*04:345, C*06:147, C*12:301, C*14:61, C*16:85
44⁹	210 bp	800 bp	27, 37	*07:27, 07:236, 07:273, 08:49, 08:126, 18:54, 18:195, 37:01:01:01-37:107, 38:17, C*07:302, C*07:929, C*07:941
45	170 bp	800 bp	16, 35, 38, 39, 3901, 3902, 67, 72	*08:55, 15:69, 15:186, 15:449, 15:494, 18:168, 35:35, 35:276, 35:468, 35:471, 38:01:01:01-38:09, 38:11-38:21, 38:23-38:29, 38:31-38:35, 38:38-38:44, 38:46-38:49, 38:51, 38:53-38:62, 38:64, 38:66-38:69, 38:71-38:102, 38:160-38:167, 38:169, 38:171-38:173Q, 39:01:01:01-39:01:01:29, 39:01:03:01-39:01:34, 39:01:36-39:20, 39:22-39:29, 39:31:01:01-39:31:01:02, 39:35, 39:37:01:01-39:42, 39:44-39:46, 39:48-39:49, 39:51-39:73, 39:75-39:81, 39:83:01-39:101, 39:103-39:106, 39:108-39:111, 39:113-39:127, 39:129-39:135, 39:137-39:155, 39:157N-39:168, 39:170, 39:172-39:173, 39:175N-39:193, 40:184, 40:424, 50:49, 51:101, 52:71, 52:84, 55:123, 58:20, 58:84, 67:01:01-67:07, C*15:51
46⁵	110 bp	1070 bp	38	*14:93, 27:07:04, 38:01:01:01-38:01:07, 38:01:09-38:01:14, 38:01:16-38:02:02:02, 38:02:04-38:03, 38:05, 38:08-38:24, 38:26-38:34N, 38:36-38:80N, 38:82-38:102, 38:160-38:173Q, 44:03:59, 49:26, 51:01:68, 52:01:10, A*23:01:18, A*25:01:08



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

Lot No.	HLA-A	HLA-B	HLA-DR	Alleles
47	395 bp	1070 bp	14, 16, 38, 39, 3901, 61, 64, 65	*07:69, 07:85:01-07:85:02, 07:180, 07:306, 08:65, 08:226, 08:228, 08:295, 14:01:01:01-14:01:07, 14:01:09-14:20, 14:22-14:41N, 14:43-14:48, 14:50-14:56, 14:58-14:72N, 14:74-14:113N, 38:01:01:01-38:02:11, 38:04-38:05, 38:08-38:25, 38:27-38:53, 38:55Q-38:73, 38:75-38:102, 38:160-38:165N, 38:167-38:173Q, 39:01:01:01-39:01:01:29, 39:01:03:01-39:01:37, 39:03:01:01-39:07, 39:09:01:01-39:09:05, 39:11-39:12:01:02, 39:14:01:01-39:15, 39:18-39:19:02, 39:22, 39:24:01-39:38Q, 39:40:01N-39:48, 39:50-39:57, 39:59-39:62, 39:64-39:82, 39:84-39:87N, 39:90-39:95N, 39:97N-39:98, 39:100, 39:102-39:104, 39:106-39:109, 39:111-39:125, 39:127-39:130, 39:132-39:142N, 39:144, 39:146N-39:153, 39:155-39:161N, 39:163, 39:165-39:169, 39:171-39:174, 39:176-39:179, 39:181-39:183, 39:185, 39:187-39:188, 39:190-39:193, 40:75 [?] , 40:322 [?] , 40:517, 49:26
48	160 bp, 425 bp	1070 bp	7, 703, 8, 14, 15, 16, 21, 27, 39, 3901, 3902, 40, 41, 42, 45, 48, 50, 51, 54, 60, 61, 64, 65, 67, 70, 73, 81	*07:02:01:01-07:02:05, 07:02:07-07:02:66, 07:02:68-07:02:77, 07:02:79-07:02:81, 07:02:83-07:26, 07:28-07:35, 07:37:01-07:37:02, 07:39-07:64, 07:66-07:80, 07:82-07:148, 07:150-07:178, 07:181N-07:191, 07:193-07:194, 07:196-07:218, 07:220-07:235, 07:237-07:256, 07:258-07:270, 07:272N, 07:274-07:447, 07:449-07:452, 07:454-07:469N, 08:01:01:01-08:01:03, 08:01:05-08:01:68, 08:04:01-08:05, 08:07-08:35, 08:37-08:51, 08:53:01-08:77, 08:79:01-08:116, 08:118-08:125, 08:127-08:130, 08:132-08:135, 08:137-08:138, 08:140-08:228, 08:230-08:254, 08:256-08:266, 08:268-08:270Q, 08:272-08:288, 08:290-08:303, 14:01:01:01-14:01:04, 14:01:06-14:01:12, 14:02:01:01-14:02:04, 14:02:06-14:02:07, 14:02:09-14:29, 14:31-14:76N, 14:78-14:92, 14:94-14:113N, 15:68, 15:180, 15:634, 27:14, 27:81, 27:130, 35:516, 35:534, 39:01:01:01-39:01:01:29, 39:01:03:01-39:06:04, 39:06:06-39:20, 39:22-39:24:01, 39:24:03-39:103, 39:105-39:184, 39:186-39:188, 39:190-39:193, 40:01:01-40:01:04, 40:01:06-40:01:79, 40:06:01:01-40:06:04:03, 40:06:06-40:07, 40:09 [?] , 40:14:01 [?] -40:14:02 [?] , 40:14:03, 40:15 [?] , 40:16:01:01-40:16:01:02, 40:18 [?] , 40:21 [?] , 40:22N-40:23, 40:24 [?] -40:25 [?] , 40:27:02 [?] , 40:28 [?] -40:29 [?] , 40:30-40:33, 40:34 [?] -40:35:02 [?] , 40:36 [?] , 40:38, 40:42-40:44, 40:45 [?] -40:46 [?] , 40:48, 40:49 [?] , 40:51-40:55, 40:56 [?] -40:58 [?] , 40:59-40:60, 40:61 [?] -40:63 [?] , 40:65 [?] , 40:66, 40:67 [?] -40:69 [?] , 40:70:01-40:70:02, 40:71 [?] , 40:72:01-40:75, 40:77, 40:78 [?] , 40:79-40:80, 40:81 [?] -40:82 [?] , 40:83, 40:84 [?] -40:85 [?] , 40:86, 40:87:01 [?] -40:88 [?] , 40:91 [?] -40:92 [?] , 40:93, 40:95 [?] , 40:96, 40:98 [?] , 40:100 [?] -40:101 [?] , 40:102-40:103, 40:104 [?] -40:108 [?] , 40:109-40:110, 40:112 [?] , 40:113-40:114:01:02, 40:114:02 [?] -40:116 [?] , 40:118N [?] , 40:121 [?] , 40:123 [?] , 40:124:01, 40:124:02 [?] -40:125:01 [?] , 40:125:02-40:128, 40:130:01 [?] , 40:130:02-40:131, 40:132 [?] -40:134 [?] , 40:136 [?] -40:137 [?] , 40:138-40:140, 40:141 [?] , 40:145 [?] , 40:146-40:148, 40:149 [?] , 40:150, 40:151 [?] , 40:152-40:153, 40:154 [?] , 40:155:01N-40:155:02N, 40:156 [?] , 40:158 [?] , 40:159-40:160:01, 40:160:02 [?] , 40:161-40:162, 40:163 [?] -40:164 [?] , 40:165, 40:166 [?] , 40:167-40:168, 40:169 [?] , 40:170-40:172, 40:173 [?] -40:176 [?] , 40:177, 40:178 [?] -40:180 [?] , 40:182 [?] , 40:183, 40:184 [?] -40:185 [?] , 40:186:01,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				40:186:02 [?] -40:187 [?] , 40:189 [?] , 40:190-40:191, 40:192 [?] -40:200 [?] , 40:203 [?] -40:213 [?] , 40:215 [?] -40:225:01 [?] , 40:226 [?] -40:228 [?] , 40:230, 40:231 [?] -40:236 [?] , 40:237, 40:238 [?] -40:243 [?] , 40:244-40:245:02, 40:246 [?] -40:248 [?] , 40:249-40:253, 40:256N, 40:257 [?] , 40:258-40:265N, 40:268-40:270, 40:272-40:273, 40:274 [?] , 40:275, 40:276 [?] , 40:277-40:278, 40:279 [?] -40:280 [?] , 40:282 [?] -40:283 [?] , 40:285-40:286N, 40:287 [?] -40:289 [?] , 40:291N [?] , 40:293 [?] , 40:294-40:295, 40:297 [?] , 40:298:01-40:298:02, 40:299 [?] -40:300 [?] , 40:301, 40:306-40:307, 40:308 [?] , 40:310-40:311, 40:312 [?] , 40:313, 40:315 [?] , 40:316-40:318, 40:319 [?] , 40:321, 40:322 [?] , 40:323:01-40:323:02, 40:324 [?] -40:329 [?] , 40:330, 40:332, 40:333 [?] , 40:336 [?] -40:339 [?] , 40:340-40:342, 40:343 [?] , 40:346, 40:348-40:350, 40:351 [?] -40:352 [?] , 40:354-40:355, 40:357 [?] -40:360 [?] , 40:361N-40:364, 40:366-40:368, 40:370, 40:372N-40:374, 40:376-40:377, 40:379:01:01-40:383, 40:385-40:389, 40:391-40:393, 40:395, 40:396 [?] , 40:397-40:398, 40:400 [?] , 40:401-40:404, 40:406-40:408, 40:411, 40:416-40:419, 40:421Q-40:428N, 40:430-40:431, 40:433-40:434, 40:436, 40:439, 40:443, 40:445-40:454, 40:456, 40:459-40:461, 40:463-40:464, 40:466, 40:468-40:469, 40:471, 40:473-40:480, 40:482-40:488N, 40:490-40:492, 40:494-40:499, 40:501-40:503, 40:506N-40:509Q, 40:510N [?] , 40:512-40:513, 40:517, 40:519-40:520, 40:523 [?] , 40:524-40:526, 41:01:01:01-41:45N, 41:47-41:59, 41:62-41:79, 42:01:01:01-42:02:02, 42:04-42:21, 42:23-42:32, 45:01:01:01-45:01:09, 45:01:11-45:25, 45:27-45:30, 46:89, 48:01:01:01-48:17, 48:19-48:55, 49:01:14, 50:01:01:01-50:01:17, 50:01:19-50:02:01:02, 50:04:01-50:22, 50:31-50:79, 51:01:09, 51:01:21, 51:01:25, 51:10, 51:24:02-51:24:04, 51:197, 51:217, 51:279, 51:338, 52:01:50, 52:54, 54:11, 56:20:02, 56:68, 67:01:01-67:07, 73:01:01:01-73:03, 81:02:01-81:02:02, C*01:214, C*03:464, C*05:152, C*07:463, C*07:758, C*16:112
49¹¹	180 bp	1070 bp	15, 40, 44, 48, 60, 62, 71	*15:116, 15:124, 15:352, 15:412, 15:507, 27:125, 35:230, 35:331, 38:63, 40:01:01-40:01:06, 40:01:08-40:01:17, 40:01:19-40:01:73, 40:01:75-40:01:76, 40:01:78-40:01:79, 40:07, 40:10:01:01-40:10:02, 40:12:01-40:12:02, 40:21-40:23, 40:25, 40:30, 40:33-40:34, 40:36, 40:38, 40:42-40:43, 40:46-40:49, 40:51-40:52, 40:54-40:55, 40:59-40:63, 40:65-40:67, 40:69, 40:73-40:74, 40:76, 40:79, 40:81, 40:84, 40:87:01-40:88, 40:92, 40:100-40:102, 40:106, 40:108, 40:112-40:114:02, 40:116-40:118N, 40:123-40:126, 40:128, 40:130:01-40:130:02, 40:132, 40:134-40:135, 40:137-40:141, 40:146-40:147, 40:149-40:156, 40:160:01-40:160:02, 40:163, 40:166, 40:168, 40:170-40:172, 40:175, 40:178-40:179, 40:182-40:183, 40:185-40:188, 40:191, 40:193-40:196, 40:198-40:199, 40:204, 40:207-40:208, 40:210, 40:212-40:213, 40:215-40:218, 40:221-40:223, 40:227-40:228, 40:231, 40:234-40:236, 40:238-40:242, 40:245:01-40:245:02, 40:247, 40:249-40:253, 40:257-40:265N, 40:272-40:273, 40:277-40:282, 40:285-40:286N, 40:288, 40:299-40:301, 40:308, 40:310, 40:312, 40:315, 40:319, 40:321, 40:323:01-40:329, 40:332-40:333, 40:336,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

			40:338N-40:339, 40:344, 40:346, 40:352-40:353:02, 40:355, 40:357-40:358, 40:360, 40:365, 40:367-40:368, 40:370, 40:372N, 40:377, 40:379:01:01-40:379:01:02, 40:381:01-40:383, 40:385-40:386, 40:391-40:398, 40:400-40:402, 40:404, 40:406-40:407, 40:416-40:419, 40:421Q-40:422, 40:424, 40:426N, 40:428N, 40:430-40:431, 40:433, 40:436, 40:439, 40:443, 40:446-40:448, 40:450, 40:452, 40:456, 40:461-40:463, 40:466, 40:474-40:475:02, 40:479-40:480, 40:483N, 40:485-40:486, 40:490-40:491, 40:495, 40:497-40:499, 40:501, 40:503, 40:506N-40:510N, 40:512-40:513, 40:516-40:517, 40:523-40:525, 44:31, 44:169, 44:182, 46:06, 46:43, 48:03:01:01-48:03:02, 48:17, 48:23, 48:40, 57:59, C*03:129, C*03:234, C*04:366, C*04:368, C*05:136, C*17:16:01-17:16:02
50 ⁹	130 bp, 290 bp, 465 bp	1070 bp	13, 17, 18, 21, 35, 4005, 41, 47, 50, 57, 60, 61, 62 *13:01:18, 13:02:06, 13:26:01, 13:72, 13:127, 15:16:03, 15:20, 15:34:02, 15:85, 15:194, 15:393, 15:407, 15:587, 15:599, 15:623, 18:22, 18:48, 18:52, 18:105, 27:107, 35:46, 35:50, 35:84, 35:207, 35:231, 37:101, 40:02:01:01-40:06:30, 40:08:01:01-40:09, 40:10:02 [?] , 40:11:01:01-40:11:02, 40:13, 40:14:01 [?] -40:14:02 [?] , 40:15 [?] , 40:18-40:20:01:02, 40:24, 40:25 [?] , 40:26-40:29, 40:34 [?] , 40:35:01-40:35:03, 40:36 [?] , 40:37, 40:39, 40:44, 40:45 [?] , 40:47 [?] , 40:49 [?] , 40:50:01:01-40:50:01:02, 40:53 [?] , 40:56, 40:57 [?] -40:58 [?] , 40:61 [?] -40:63 [?] , 40:64:01:01-40:64:02, 40:65 [?] , 40:67 [?] , 40:68, 40:69 [?] , 40:70:01-40:71, 40:74-40:75, 40:78, 40:81 [?] , 40:82, 40:84 [?] , 40:85-40:86, 40:87:01 [?] -40:88 [?] , 40:89-40:91, 40:92 [?] , 40:94-40:99, 40:100 [?] -40:101 [?] , 40:103-40:105, 40:106 [?] , 40:107, 40:108 [?] , 40:109, 40:111, 40:112 [?] , 40:114:02 [?] , 40:115, 40:116 [?] , 40:118N [?] , 40:119-40:120, 40:121 [?] , 40:122, 40:123 [?] , 40:124:02 [?] -40:125:01 [?] , 40:127, 40:129, 40:130:01 [?] , 40:131, 40:132 [?] , 40:133Q, 40:134 [?] -40:136 [?] , 40:141 [?] , 40:142N-40:145, 40:148-40:149, 40:151 [?] , 40:154 [?] , 40:156 [?] , 40:157, 40:158 [?] , 40:159, 40:160:02 [?] , 40:161-40:162, 40:163 [?] -40:164 [?] , 40:165, 40:166 [?] , 40:167, 40:169, 40:173-40:174, 40:175 [?] , 40:176, 40:178 [?] -40:179 [?] , 40:180-40:181, 40:182 [?] , 40:184, 40:186:02 [?] -40:188 [?] , 40:189-40:190, 40:192 [?] -40:199 [?] , 40:200-40:203, 40:204 [?] , 40:205-40:206, 40:207 [?] -40:208 [?] , 40:210 [?] , 40:211, 40:212 [?] -40:213 [?] , 40:214, 40:215 [?] -40:219 [?] , 40:220, 40:221 [?] -40:223 [?] , 40:224-40:226, 40:227 [?] -40:228 [?] , 40:229-40:230, 40:231 [?] , 40:232, 40:233 [?] -40:236 [?] , 40:238 [?] -40:242 [?] , 40:243-40:244, 40:246, 40:247 [?] , 40:248, 40:254-40:256N, 40:257 [?] , 40:266-40:271, 40:274-40:276, 40:279 [?] -40:282 [?] , 40:283-40:284, 40:287, 40:288 [?] , 40:289-40:293, 40:295-40:297, 40:298:01 [?] -40:300 [?] , 40:302-40:307, 40:308 [?] , 40:309, 40:311, 40:312 [?] , 40:313-40:314, 40:315 [?] , 40:317, 40:318 [?] -40:319 [?] , 40:320, 40:322, 40:324 [?] -40:329 [?] , 40:330-40:331, 40:333 [?] , 40:334:01-40:335, 40:336 [?] , 40:337N, 40:338N [?] -40:339 [?] , 40:340-40:341, 40:342 [?] , 40:343, 40:344 [?] , 40:345N, 40:347-40:351, 40:354, 40:356, 40:357 [?] -40:358 [?] , 40:359, 40:360 [?] , 40:361N-40:364, 40:366, 40:369, 40:371, 40:374-40:376, 40:384, 40:387-40:390, 40:396 [?] , 40:399N, 40:400 [?] , 40:405, 40:408-40:415, 40:420, 40:423, 40:425, 40:427, 40:429, 40:432, 40:434-40:435, 40:437-40:438N, 40:440-40:442, 40:444-40:445, 40:449,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				40:451, 40:453-40:455, 40:458-40:460, 40:464-40:465, 40:467-40:468, 40:470-40:473, 40:476, 40:478, 40:481N-40:482, 40:484, 40:487N-40:489, 40:492-40:494, 40:500, 40:504-40:505, 40:510N ² , 40:511N, 40:514N-40:515, 40:518-40:522, 40:523 ² , 40:526, 44:02:59, 44:103, 44:131, 44:139 ² , 44:212, 46:01:07, 46:33, 47:01:01:02-47:01:02, 47:01:04-47:01:05, 47:02-47:03, 47:05-47:10, 57:05, 57:67:01, 58:36, C*01:73, C*03:504
51 ⁵	105 bp	1070 bp	7, 8, 22, 41, 42	*07:04:01-07:04:02, 07:25, 07:146, 07:349, 08:01:01:01-08:01:24, 08:01:26-08:01:61, 08:01:63-08:05, 08:07-08:12:03, 08:14-08:19N, 08:21-08:24, 08:26:01-08:30N, 08:32-08:39, 08:41-08:54:02, 08:56:01-08:69, 08:71-08:78, 08:80-08:102, 08:104-08:112, 08:114-08:133, 08:135-08:142, 08:144-08:155, 08:157-08:168, 08:170, 08:172-08:174, 08:176-08:204, 08:206-08:254, 08:256-08:273, 08:275-08:280, 08:282-08:303, 35:87, 40:136, 40:231, 40:279, 41:01:01:01-41:08, 41:10-41:28, 41:30-41:31, 41:33-41:61, 41:63-41:79, 42:01:01:01-42:02:02, 42:04-42:07, 42:09-42:18, 42:20-42:32, 44:166, 44:342, 55:20
52	325 bp	800 bp	7, 15, 22, 27, 2708, 35, 42, 45, 46, 54, 55, 56, 67, 73, 78, 81, 82	*07:02:01:01-07:02:05, 07:02:07-07:02:97, 07:04:01-07:07:03, 07:09:01-07:15, 07:17-07:26, 07:28-07:31:02, 07:33:01-07:36, 07:39-07:46, 07:47 ^w , 07:48-07:49N, 07:51-07:68:03, 07:70-07:84, 07:86-07:179, 07:181N-07:206, 07:208-07:218, 07:220-07:227, 07:229:01-07:231N, 07:233-07:305, 07:307-07:399, 07:401-07:452, 07:454-07:469N, 08:123, 08:132, 14:21, 14:57, 15:76, 15:101, 15:255, 15:509, 15:594, 27:01-27:11, 27:13:01-27:15, 27:17, 27:19:01:01-27:21:02, 27:24-27:28, 27:30-27:38, 27:40-27:58, 27:60-27:76, 27:78-27:91, 27:93-27:118, 27:120-27:128, 27:130-27:136, 27:138-27:152, 27:154-27:203, 27:205-27:214, 27:216-27:241, 27:243N-27:264, 35:76, 35:186, 35:229, 35:368, 35:561, 37:60, 38:26, 40:166, 40:513, 42:01:01:01-42:02:02, 42:04-42:06, 42:08-42:10, 42:12-42:23, 42:25, 42:27-42:32, 44:90, 44:97, 44:356, 45:06, 46:01:01:01-46:01:40, 46:02 ^w , 46:03-46:40:02, 46:42-46:48, 46:50-46:98N, 54:01:01:01-54:29, 54:31-54:46, 55:01:01:01-55:05, 55:07-55:17, 55:19-55:77, 55:79-55:132N, 56:01:01:01-56:16, 56:18-56:22, 56:24-56:32, 56:34-56:41, 56:43, 56:45-56:53, 56:55:01:02-56:66, 56:68-56:88, 67:01:01-67:07, 73:01:01:01-73:03, 81:01:01:01-81:10, 82:01:01:01-82:04, 83:01
53 ^{5,7}	80 bp, 115 bp, 160 bp, 195 bp, 225 bp, 260 bp	1070 bp	12, 35, 37, 44, 45, 47, 48, 57	*07:359, 08:49, 08:60, 08:76, 08:87, 08:129, 08:181, 08:213:01, 08:255, 15:578, 15:592, 27:127, 35:38, 35:45, 35:71, 35:115, 35:169, 35:260, 35:309, 35:452, 35:513N, 35:545, 37:01:01:01-37:01:07, 37:01:09, 37:01:11-37:01:29, 37:03N-37:06:02, 37:08, 37:10-37:38, 37:41-37:43, 37:45-37:50, 37:52-37:65, 37:67-37:82N, 37:84-37:93, 37:95-37:96, 37:98-37:107, 40:132, 40:178, 40:231, 40:279, 40:458, 40:523, 41:48, 42:13, 44:02:01:01-44:14, 44:16-44:17, 44:19N-44:36, 44:38-44:63, 44:65-44:134, 44:136-44:149N, 44:151-44:157, 44:159-44:165, 44:167-44:183, 44:185-44:212, 44:214-44:229, 44:231-44:235, 44:237N-44:245:02, 44:247-44:267N, 44:269-44:274, 44:277-44:289, 44:291-44:338, 44:340-44:346, 44:348-44:354N, 44:356-44:363,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				44:437-44:458, 44:460-44:524, 44:526-44:546, 48:11, 51:42, 51:274, 51:359, 53:22, 53:62, 57:07, 57:09, 57:24, 57:93, 57:119, 58:94N, 83:01, C*15:25
54	215 bp	800 bp	5, 8, 12, 17, 21, 22, 41, 42, 44, 45, 49, 51, 5102, 5103, 52, 53, 56, 62	*08:09, 08:84, 08:254, 13:46, 13:154, 15:83:01:01-15:83:01:02, 41:01:01:01-41:01:08, 41:05-41:07, 41:09, 41:12, 41:14, 41:16-41:17, 41:20-41:22, 41:25-41:26, 41:28-41:30, 41:32-41:35, 41:37, 41:44, 41:50, 41:53-41:54, 41:56Q-41:57, 41:59-41:61, 41:64, 41:66-41:67, 41:70, 41:73, 41:75-41:79, 42:04, 44:06, 44:15:01:01-44:15:01:02, 44:18:01:01-44:18:01:02, 44:20, 44:100, 44:355, 45:01:01:01-45:25, 45:27-45:30, 49:03, 51:01:01:01-51:01:01:108, 51:01:03-51:01:08, 51:01:10-51:01:39, 51:01:41-51:02:01:04, 51:02:03-51:04:02, 51:06:01:01-51:06:03, 51:07:01, 51:08:01:01-51:09:01, 51:09:03-51:14, 51:16-51:24:05, 51:26-51:46, 51:48-51:53, 51:55-51:77, 51:79-51:88, 51:90-51:92:02, 51:94-51:102, 51:104:01-51:115, 51:117-51:118N, 51:120-51:131, 51:133-51:143, 51:145-51:147, 51:149N-51:175, 51:177-51:178N, 51:181-51:217, 51:219-51:225, 51:227-51:237:01, 51:238-51:336, 51:338-51:340, 51:342-51:376, 52:01:02:01-52:01:02:09, 52:01:04, 52:01:09, 52:01:40, 52:01:43, 52:01:47, 52:01:49, 52:02:01, 52:03, 52:06:01-52:06:02, 52:09:01, 52:10:03, 52:11:02, 52:19, 52:21:01, 52:31:02, 52:85, 52:97, 52:102, 52:104, 53:01:13, 53:06, 55:20, 55:56, 56:13, 57:01:20, 58:08:01, A*23:01:14
55⁹	130 bp, 280 bp	800 bp	7, 17, 18, 22, 35, 37, 40, 46, 51, 53, 60, 62, 78	*07:12, 07:14:01-07:14:02, 07:18:01-07:18:02, 07:137, 07:162-07:163, 07:199, 07:225, 07:227, 07:248, 07:350, 14:10, 14:102, 15:57 ^w , 15:289, 15:297, 15:318, 15:579, 15:645, 18:22, 18:69, 18:105, 18:149-18:150, 18:204, 35:21, 35:24:01-35:24:02, 35:81, 35:96, 35:109, 35:157, 35:188, 35:190, 35:233, 35:287, 37:04:01-37:04:02, 40:28, 40:30, 40:34, 40:59, 40:137, 40:160:01-40:160:02, 44:150, 44:190, 44:203:01-44:203:02, 44:275, 44:290, 46:01:01:01-46:39, 46:41N-46:46, 46:48-46:71, 46:73-46:76, 46:78-46:85, 46:87-46:98N, 48:17, 51:04:01-51:04:02, 51:46, 51:56:01-51:56:03, 51:139, 51:227, 51:281, 51:310, 52:88, 53:02, 53:06, 53:28, 57:14:01-57:14:02, 58:09, 58:76, C*01:73, C*03:260, C*03:384, C*15:39, C*15:130, C*15:195
56⁵	90 bp, 410 bp	1070 bp	7, 15, 22, 27, 2708, 35, 47, 54, 55, 56, 78, 81, 82	*07:65, 07:197, 07:209, 07:271, 15:342, 18:165, 27:01 ^w , 27:02:01:01-27:21:02, 27:23-27:152, 27:154-27:214, 27:216-27:264, 35:76, 35:186, 35:229, 35:368, 35:561, 37:02, 37:60, 40:429, 44:90, 44:97, 44:356, 47:01:01:02-47:12, 48:26, 54:01:01:01-54:29, 54:31-54:46, 55:01:01:01-55:05, 55:07-55:17, 55:19-55:77, 55:79-55:132N, 56:01:01:01-56:16, 56:18-56:22, 56:24-56:32, 56:34-56:41, 56:43, 56:45-56:53, 56:55:01:02-56:66, 56:69-56:88, 81:01:01:01-81:01:05, 81:05, 81:09-81:10, 82:01:01:01-82:04, 83:01, C*04:439, C*14:79
57⁵	90 bp, 175 bp	1070 bp	15, 27, 35, 48, 60, 62, 75, 77, 81	*07:174, 07:202, 07:222, 07:279, 07:409 ^w , 08:213:01-08:213:02, 13:36, 15:02:01:01-15:02:12, 15:08:01:01-15:08:02, 15:11:01-15:11:07, 15:13:01-15:13:03, 15:15:01:01-15:15:01:02, 15:21:01:01-15:21:01:02, 15:31:01:01-15:31:01:03, 15:44, 15:55, 15:76 ^w ,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				15:88-15:89:02, 15:112, 15:121, 15:139, 15:144, 15:148, 15:170, 15:189:01-15:189:02, 15:191, 15:194, 15:209N, 15:213-15:215, 15:223, 15:265, 15:283-15:284, 15:291, 15:297, 15:301-15:303, 15:305, 15:319, 15:325, 15:328, 15:330, 15:341, 15:345, 15:350, 15:357-15:358, 15:367, 15:374, 15:378, 15:384:01-15:384:02, 15:394, 15:399, 15:404, 15:418, 15:420, 15:425, 15:437-15:438, 15:454N, 15:460, 15:479, 15:486, 15:491, 15:495, 15:505, 15:508, 15:513, 15:537, 15:557, 15:567, 15:569, 15:574, 15:592, 15:596N, 15:615, 15:622-15:623, 15:627, 15:629-15:630, 15:637, 15:646, 27:24, 35:46, 35:311, 40:31, 40:45, 40:80, 48:01:01:01-48:01:11, 48:04:01:01-48:04:02, 48:06-48:07, 48:09, 48:11, 48:15-48:16, 48:18-48:20, 48:22, 48:24, 48:27-48:39, 48:41-48:55, 51:104:02, 56:39 ^w , 81:01:01:01-81:10
58⁵	55 bp, 100 bp, 145 bp, 395 bp, 430 bp	1070 bp	15, 27, 44, 48, 49, 52, 59, 60, 61, 62, 63	*15:16:01:01-15:17:07, 15:24:01-15:24:02, 15:67, 15:87, 15:143, 15:152, 15:157, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:222, 15:230, 15:254, 15:268, 15:273, 15:356, 15:361-15:362, 15:396, 15:403, 15:408, 15:411, 15:423-15:424, 15:442, 15:446, 15:500, 15:516, 15:523, 15:532, 15:546Q, 15:550, 15:555, 15:575N, 15:586, 15:602-15:603, 15:613, 15:619, 15:642, 15:644, 27:02:01:01-27:02:04, 27:02:06, 27:30, 27:53, 27:57, 27:62, 27:65N, 27:75, 27:77, 27:83, 27:95, 27:102, 27:119, 27:126, 27:134, 27:163, 27:171-27:172, 27:176N, 27:181, 27:188, 27:197, 27:203, 27:213, 27:232, 27:236, 27:245-27:246N, 27:251, 37:34, 38:25, 40:13, 40:19, 40:61, 40:109, 40:117, 40:292, 40:340, 40:394, 41:46, 41:60, 44:18:01:01-44:18:01:02, 44:25, 44:50:01, 44:50:03, 44:95, 44:355, 48:18, 49:01:01:01-49:01:24, 49:03-49:14, 49:16-49:25, 49:27-49:69, 49:71-49:78, 51:112, 51:231, 51:310, 52:01:01:01-52:01:28, 52:01:30-52:19, 52:21:01-52:72, 52:74-52:110N, 53:17:01-53:17:02, 53:28, 53:30, 53:39, 54:12, 55:103, 56:21, 56:42, 56:54, 56:65, 57:45, 57:51, 57:69, 57:128, 58:100, 58:111, 59:01:01:01-59:07, 59:09-59:12, A*24:309, A*32:39
59^{5,8}	120 bp, 210 bp	800 bp	5, 13, 15, 17, 22, 35, 45, 49, 50, 51, 5102, 5103, 52, 55, 56, 62, 63, 78	*07:78, 07:84, 07:429, 13:16, 13:31 ^w , 13:48, 13:62, 13:128, 13:131 ^w , 13:145 ^w , 14:37, 15:04:01:01 ^w -15:04:04 ^w , 15:16:01:01-15:16:05, 15:42, 15:67, 15:95, 15:137 ^w , 15:254, 15:293, 15:303 ^w , 15:310, 15:362, 15:408, 15:517, 15:531 ^w , 15:550, 15:602, 15:609 ^w , 15:613, 27:153, 35:265 ^w , 40:06:30, 40:95, 40:148, 40:161, 44:459, 46:32:01 ^w -46:32:02 ^w , 49:01:01:01-49:01:01:18, 49:01:03-49:10, 49:12-49:14, 49:16-49:44, 49:46, 49:48-49:60N, 49:62-49:78, 50:01:01:01-50:02:01:02, 50:04:01-50:11, 50:13-50:16, 50:18-50:22, 50:31-50:46, 50:47 ^w , 50:48, 50:50, 50:52-50:65, 50:67-50:69, 50:71-50:79, 51:01:01:01-51:03, 51:05:01:01-51:05:01:02, 51:07:01-51:08:04, 51:08:06-51:09:03, 51:11N-51:14, 51:16-51:24:05, 51:26-51:41N, 51:43-51:44N, 51:48-51:55, 51:57-51:58, 51:60-51:61:02, 51:63:01-51:63:02, 51:65-51:80, 51:82-51:134, 51:136-51:138, 51:140-51:147, 51:149N-51:156, 51:158:01-51:193:02, 51:195-51:212, 51:214-51:219, 51:221-51:226, 51:228-51:280, 51:282-51:292, 51:294-51:309, 51:311-51:318N, 51:320-51:337, 51:339-51:340, 51:342-51:351, 51:354-51:358, 51:360-51:376,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				52:01:01:01-52:87, 52:89N-52:96N, 52:97 ^w , 52:98-52:110N, 54:20, 54:33 ^w , 55:01:01:01-55:01:06, 55:01:08-55:01:29, 55:02:11, 55:03, 55:05, 55:09, 55:11, 55:15, 55:17, 55:21 ^w , 55:24-55:25, 55:28:01-55:29, 55:31, 55:33, 55:36, 55:38, 55:40, 55:44-55:45, 55:52-55:55N, 55:58-55:60, 55:64, 55:66, 55:68, 55:73-55:76, 55:78-55:79, 55:84-55:85, 55:87, 55:90-55:95, 55:97N-55:106, 55:108-55:109, 55:111, 55:113-55:116, 55:118N-55:119, 55:121-55:122, 55:124-55:127, 56:05:01-56:06, 56:21, 56:25, 56:37, 58:08:01-58:08:02, 59:12, 78:01:01:02-78:10, C*03:464, C*05:152, C*12:183^w, C*14:92
60	390 bp, 430 bp	1070 bp	5, 17, 27, 35, 44, 51, 5102, 5103, 52, 53, 57, 58, 61, 78	*15:506, 18:33, 18:67, 18:136, 18:194:01:01-18:194:01:02, 27:02:01:01-27:02:06, 27:12:01:01-27:12:01:03, 27:16, 27:30, 27:53, 27:57, 27:62, 27:65N, 27:75, 27:77, 27:83, 27:92, 27:95, 27:102, 27:119, 27:126, 27:134, 27:153, 27:157, 27:163, 27:171-27:172, 27:176N, 27:181, 27:188, 27:197, 27:203, 27:213, 27:232, 27:236, 27:242, 27:245-27:246N, 27:251, 35:26, 35:82, 35:85, 35:135, 35:294, 35:471, 37:34, 37:57, 39:189, 40:13, 40:19, 40:75 [?] , 40:109 [?] , 40:292 [?] , 40:322 [?] , 40:340 [?] , 44:06, 44:25, 44:50:01, 44:50:03, 44:95, 44:263, 49:62, 51:01:01:01-51:01:39, 51:01:41-51:24:05, 51:26-51:46, 51:48-51:53, 51:55-51:77, 51:79-51:340, 51:342-51:366, 51:368-51:370, 51:372-51:376, 52:01:01:01-52:01:28, 52:01:30-52:19, 52:21:01-52:89N, 52:91-52:110N, 53:01:01:01-53:01:17, 53:01:19-53:02, 53:04-53:08:02, 53:10, 53:14-53:30, 53:32-53:35, 53:37, 53:39-53:54, 53:56-53:66, 53:68-53:72, 57:01:01:01-57:04:02, 57:05-57:11, 57:13-57:164, 58:01:01:01-58:02:02, 58:04-58:16:02, 58:18-58:29, 58:31N-58:63, 58:65-58:107, 58:109-58:136, 58:138-58:139, 78:03
61	145 bp, 210 bp	1070 bp	7, 12, 21, 27, 2708, 35, 4005, 41, 45, 46, 47, 50, 60, 61, 73	*07:133, 07:177, 08:102, 08:122, 08:155, 13:09, 13:39, 13:67, 14:48, 15:46, 15:53, 15:106, 15:143, 15:183, 15:212, 15:336, 15:345, 15:511, 15:525, 15:559, 18:48, 18:80, 18:102, 18:107, 27:08, 27:12:01:01-27:12:01:03, 27:18, 27:26, 27:33, 27:40, 27:42, 27:44, 27:89, 27:153, 27:165 ^w , 27:204:01:01-27:204:01:02, 27:242, 35:19, 35:47, 35:63, 35:84, 35:90, 35:150:01-35:150:02, 35:154, 35:184, 35:201, 35:217, 35:274, 35:317, 35:441, 35:450, 39:69, 40:01:01-40:11:02, 40:14:01-40:16:01:02, 40:18, 40:20:01:01-40:20:01:02, 40:22N-40:40, 40:42-40:46, 40:48-40:75, 40:77-40:95, 40:97-40:108, 40:111-40:116, 40:118N-40:136, 40:138-40:156, 40:158-40:187, 40:189-40:200, 40:202-40:291N, 40:293-40:339, 40:341-40:364, 40:366-40:393, 40:395-40:399N, 40:401-40:526, 41:01:01:01-41:04, 41:06-41:25, 41:27-41:45N, 41:47-41:59, 41:61-41:79, 44:09, 44:46:01-44:46:02, 44:75, 44:90, 44:129, 44:131, 44:254, 44:277, 44:356, 44:438N, 45:01:01:01-45:01:08, 45:01:10-45:15, 45:17-45:25, 45:27-45:30, 46:10, 47:02, 47:03 ^w , 48:41, 50:01:01:01-50:02:01:02, 50:04:01-50:22, 50:31-50:79, 56:36, 56:67, 73:01:01:01-73:03, C*03:504, C*07:828
62^s	300 bp	1070 bp	7, 12, 13, 15, 17, 18, 21, 27, 2708, 35,	*07:54, 07:123, 07:400, 07:453, 08:17, 08:38:01-08:38:02, 08:54:01-08:54:02, 08:101, 08:228, 08:285, 13:01:01:01-13:01:14, 13:01:16-13:04, 13:06, 13:08-13:23, 13:25-13:35, 13:37-13:174, 15:01:01:01-15:01:04, 15:01:06-15:01:18,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

37, 38, 39, 15:01:20-15:01:78, 15:03:01:01-15:03:04, 15:03:06-
3902, 40, 15:07:03, 15:12:01-15:12:03, 15:14, 15:19-15:20, 15:24:01-
4005, 41, 15:27:04, 15:30:01:01-15:30:01:03, 15:32:01-15:36,
44, 45, 47, 15:38:01-15:40:02, 15:42-15:43, 15:45-15:50, 15:53-15:54,
48, 49, 50, 15:56-15:58, 15:60-15:63, 15:65-15:66, 15:68-15:71,
52, 60, 61, 15:73:01:01-15:75, 15:77-15:79N, 15:81-15:87, 15:91-15:92,
62, 70, 72, 15:94N, 15:96-15:98, 15:101-15:107, 15:109-15:111N,
76 15:113, 15:116-15:118, 15:122-15:123, 15:125:01-15:129,
15:131-15:132, 15:135-15:138, 15:140-15:142, 15:145-
15:147, 15:150-15:152, 15:154-15:160:01:02, 15:163-
15:167, 15:169, 15:171-15:175, 15:178-15:179:02, 15:181N-
15:185, 15:187-15:188, 15:190N, 15:192-15:193, 15:195,
15:199, 15:201-15:207, 15:210-15:212, 15:217-15:220:04,
15:224-15:225, 15:227-15:228, 15:231-15:237, 15:239-
15:245:02Q, 15:247-15:251:03, 15:253, 15:255-15:262:02N,
15:264, 15:266-15:267, 15:269-15:272N, 15:274, 15:276-
15:282, 15:286-15:289, 15:295-15:296, 15:298-15:300,
15:304N, 15:309-15:310, 15:315-15:318, 15:320-15:322,
15:324, 15:326-15:327, 15:331-15:334, 15:336, 15:339-
15:340, 15:342-15:344, 15:346-15:349:02, 15:352-15:353,
15:355, 15:359-15:360, 15:363:01-15:366, 15:368-15:373,
15:375N-15:377Q, 15:379, 15:381, 15:383, 15:385-15:387,
15:389-15:391, 15:393, 15:395, 15:397-15:398, 15:400N-
15:402, 15:405-15:407, 15:409-15:410, 15:412-15:413,
15:415-15:417, 15:419, 15:421-15:422, 15:429-15:432,
15:440-15:444, 15:448-15:450, 15:452-15:453, 15:456-
15:459, 15:461, 15:463N-15:478, 15:480-15:485, 15:488-
15:490, 15:492-15:494, 15:496N-15:497, 15:499, 15:501-
15:502, 15:507, 15:510-15:512, 15:514-15:515, 15:517-
15:522, 15:524-15:525, 15:527-15:529, 15:531, 15:533-
15:536, 15:538-15:545, 15:547, 15:549N, 15:551-15:552,
15:554, 15:556, 15:558-15:562N, 15:564-15:566, 15:568,
15:571-15:572, 15:576-15:591, 15:593-15:595, 15:597-
15:601, 15:604N-15:606, 15:608-15:610, 15:614:01:01-
15:614:01:02, 15:616Q-15:618, 15:621, 15:624-15:626,
15:628, 15:631-15:636, 15:638-15:641, 15:645, 15:647,
18:12:01-18:12:02, 18:178, 18:195, 27:01-27:21:02, 27:24-
27:91, 27:93-27:108, 27:110-27:156, 27:158-27:214, 27:216-
27:241, 27:243N-27:245, 27:247-27:264, 35:10, 35:13,
35:16, 35:28:01:01-35:28:01:03, 35:63, 35:69, 35:80, 35:207,
35:226, 35:270, 35:317-35:318, 35:361, 35:477, 35:534,
35:541, 37:01:01:01-37:07, 37:09-37:66, 37:68-37:107,
38:03-38:04, 38:25, 39:02:01:01-39:02:03, 39:08, 39:13:01-
39:13:02, 39:22-39:23, 39:49, 39:101, 39:105, 39:110,
39:162, 39:164, 39:184, 40:01:01-40:07, 40:09-40:16:01:02,
40:18-40:24, 40:26-40:40, 40:42-40:67, 40:70:01-40:105,
40:107-40:165, 40:167-40:231, 40:233-40:312, 40:314-
40:356, 40:358-40:480, 40:482-40:526, 41:01:01:01-41:79,
42:11, 44:02:01:01-44:02:21, 44:02:23-44:02:64, 44:02:66-
44:03:35, 44:03:37-44:03:44, 44:03:46-44:05:06, 44:07-
44:11, 44:13-44:170, 44:172-44:231, 44:233-44:245:02,
44:247-44:262, 44:264-44:355, 44:357-44:363, 44:437-
44:524, 44:526-44:546, 45:01:01:01-45:21, 45:23-45:25,
45:27-45:30, 47:01:01:02-47:12, 48:01:01:01-48:05,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				48:07-48:15, 48:17-48:38, 48:40-48:55, 49:01:01:01-49:14, 49:16-49:25, 49:27-49:62, 49:64-49:66, 49:68-49:78, 50:01:01:01-50:02:01:02, 50:04:01-50:14, 50:16-50:22, 50:31-50:79, 51:310, 52:01:01:01-52:72, 52:74-52:110N, 53:17:01-53:17:02, 53:28, 53:38, 54:35, 55:18, 55:34, 55:40, 56:67, 78:05-78:06
63 ^{7,13,14}	160 bp, 370 bp	1070 bp	8, 22, 44, 45, 47, 51, 54, 57	*08:87, 08:255, 15:578, 27:127, 41:55, 44:02:01:01-44:02:78, 44:06, 44:08-44:09, 44:11-44:12, 44:16, 44:19N-44:24, 44:27:01:01-44:27:04, 44:33-44:34:02, 44:41:01-44:41:02, 44:44, 44:48-44:49, 44:52N-44:53:02, 44:55, 44:59:01-44:59:02, 44:63, 44:66-44:68, 44:71-44:74, 44:80, 44:83-44:84:02, 44:86-44:87, 44:89-44:91, 44:93, 44:95, 44:97, 44:99-44:102, 44:104, 44:106, 44:112-44:113, 44:116, 44:118-44:119, 44:121, 44:126:01-44:127, 44:131-44:132, 44:137-44:138Q, 44:142, 44:145, 44:148-44:149N, 44:151-44:152, 44:162, 44:168-44:173, 44:176-44:177, 44:179, 44:185, 44:187, 44:190-44:191, 44:195N-44:196, 44:200-44:201, 44:206, 44:208, 44:211-44:212, 44:214, 44:216, 44:218-44:221, 44:225, 44:229, 44:235, 44:238, 44:240-44:244, 44:249, 44:253-44:255, 44:257, 44:260-44:262, 44:265, 44:267N, 44:269-44:270:02, 44:273-44:274, 44:279, 44:287, 44:289, 44:291-44:293, 44:296-44:297, 44:301, 44:303N, 44:306N, 44:311-44:315, 44:320-44:322, 44:326-44:327, 44:329, 44:331, 44:334N-44:335, 44:337-44:338, 44:341N-44:344, 44:350-44:351, 44:353, 44:356-44:359, 44:361-44:363, 44:437, 44:439, 44:441-44:442, 44:445-44:446, 44:448N-44:449N, 44:451-44:453, 44:455, 44:457, 44:461-44:464:01:02, 44:467-44:468, 44:470, 44:473, 44:477-44:480Q, 44:482, 44:484, 44:492, 44:494-44:495, 44:497-44:499, 44:502-44:503, 44:505, 44:515, 44:517-44:518, 44:520-44:521, 44:523N, 44:526, 44:528, 44:532, 44:534-44:535, 44:537, 44:540, 44:543-44:544N, 51:274, 54:01:01:01-54:46, 83:01
64	180 bp, 210 bp, 240 bp	1070 bp	13, 15, 18, 35, 39, 40, 54, 55, 56, 59, 78	*07:65, 07:206, 08:09, 08:26:03, 08:56:03, 08:84, 08:181, 13:06, 13:46, 13:53, 13:131, 14:37, 15:01:76, 15:32:02, 15:42, 15:83:01:01-15:83:01:02, 15:86, 15:224, 15:485, 18:04:01, 18:98:02, 18:125, 18:191, 35:03:27, 35:08:05, 35:42:02, 35:60, 35:185, 35:276, 39:06:01:01-39:06:08, 39:17, 39:34, 39:57, 39:62-39:64, 39:83:01-39:83:03, 39:90, 39:109, 39:127, 39:129, 39:132, 39:140-39:141, 39:154, 39:178, 39:186-39:187, 40:01:11, 40:58, 41:66, 42:04, 45:08, 46:18, 46:40:01, 49:76, 50:45, 50:49, 51:21, 51:36, 51:101, 51:143, 51:153:01-51:153:02, 51:157, 51:197, 51:280, 51:295, 51:300, 51:340, 51:376, 52:37, 52:71, 52:84, 54:01:01:01-54:02, 54:04-54:05N, 54:07-54:08N, 54:10-54:25, 54:27-54:32, 54:34-54:39, 54:41-54:46, 55:01:01:01-55:03, 55:05, 55:07, 55:10-55:21, 55:23, 55:25-55:48, 55:50, 55:52-55:85, 55:87-55:109, 55:111-55:132N, 56:05:01, 56:06, 56:10, 56:23, 56:54, 56:66, 57:26, 59:01:01:01-59:03, 59:05-59:08, 59:10N-59:11, 78:01:01:02-78:01:02, 78:02:02-78:03, 78:07-78:10, C*01:90, C*06:72, C*07:536, C*07:758, C*08:207, C*12:183, C*14:92, C*15:91, C*15:128



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

65⁸	180 bp, 205 bp	1070 bp	7, 12, 15, 21, 22, 44, 45, 49, 50, 56, 61, 62, 82	*07:29, 07:186, 07:383, 08:123, 08:132, 13:03, 13:48, 13:120, 13:146, 13:164, 15:73:01:01-15:73:01:02, 15:303, 15:486, 15:552, 35:432, 35:446, 35:545, 40:71, 44:10, 44:15:01:01-44:15:01:02, 44:18:01:01-44:18:01:02, 44:140, 44:339-44:340, 44:355, 45:01:01:01-45:01:12, 45:04-45:07, 45:11-45:23, 45:24 ^w , 45:25, 45:27-45:28N, 45:30, 46:11, 49:01:01:01-49:01:07, 49:01:09-49:03, 49:06, 49:08-49:14, 49:16-49:17, 49:19N-49:49, 49:51-49:60N, 49:62-49:75, 49:77-49:78, 50:01:01:01-50:02:01:02, 50:04:01-50:08, 50:10-50:13, 50:15-50:19, 50:21-50:22, 50:32, 50:34-50:38, 50:40-50:44, 50:46-50:48, 50:50-50:79, 51:15, 51:62, 51:106:01-51:106:02, 51:223, 51:225, 52:25:01-52:25:02, 54:03, 54:33, 54:40, 56:01:01:01-56:01:04, 56:01:06- 56:01:12, 56:01:14-56:02:03, 56:04:01:01-56:04:04, 56:07- 56:08, 56:13-56:14, 56:16-56:17, 56:20:01-56:20:02, 56:24- 56:30, 56:33-56:49, 56:51-56:53, 56:55:01:02-56:59, 56:61- 56:65, 56:67-56:87, 56:89, 58:73, 59:04, 59:12, 82:01:01:01- 82:04, C*03:19, C*03:102
66^{5,8}	90 bp, 150 bp, 410 bp	800 bp	57	*15:462, 39:172, 57:01:01:01-57:01:22, 57:01:24-57:01:44, 57:01:46-57:04:01, 57:04:03-57:35, 57:37-57:88, 57:90- 57:108, 57:110-57:127, 57:129-57:164, 58:36, C*06:140
67^{5,12}	90 bp	1070 bp	5, 13, 17, 39, 55, 56, 58, 60, 61, 62, 71	*07:182N, 07:424, 13:11, 13:58, 15:18:03, 15:73:01:01- 15:73:01:02, 15:86, 15:224, 15:486, 18:191, 35:60, 35:291, 35:432, 35:446, 35:545, 39:17, 40:48, 40:71, 40:377, 40:522, 45:22, 49:03, 49:47:01-49:47:02, 50:51 ^w , 51:15, 51:62, 51:106:02, 51:157, 51:225, 51:300, 52:25:01, 55:130, 56:02:01:01-56:02:03, 56:04:01:01-56:04:04, 56:10, 56:45, 56:51, 56:56-56:57, 56:63, 56:66, 57:120, 58:01:01:01- 58:01:20, 58:01:21 ^w , 58:01:22-58:02:02, 58:04-58:29, 58:31N-58:35, 58:37-58:43, 58:45:01-58:99, 58:101-58:110, 58:112-58:139, 82:03, A*02:339, A*02:510, C*02:06:01- 02:06:02, C*02:47, C*05:132, C*06:168^w, C*06:224^w, C*07:536, C*07:725, C*08:207, C*12:03:66^w, C*12:15^w, C*12:208, C*15:58^w, C*15:72, C*15:225
68⁵	95 bp	1070 bp	8, 18, 22, 35, 37, 38, 39, 44, 51, 5102, 5103, 53, 78	*07:65 ^w , 07:134 ^w , 08:32, 08:133, 08:180, 08:239, 15:308, 18:01:01:01-18:03:01:03, 18:04:01-18:11, 18:13-18:15, 18:17N-18:36, 18:38-18:47, 18:49-18:65, 18:67-18:106Q, 18:108-18:124, 18:126-18:177, 18:179-18:194:01:02, 18:196-18:226, 35:01:01:01-35:01:55, 35:01:57-35:08:04, 35:08:06-35:09:04, 35:11:01-35:12:04, 35:14:01-35:15:02, 35:17:01:01-35:18, 35:20:01-35:24:02, 35:27, 35:29:01- 35:43:01, 35:43:03-35:45, 35:48, 35:50-35:62, 35:64:01- 35:68:02, 35:70-35:72, 35:74-35:75, 35:76 ^w , 35:77-35:79, 35:81-35:153, 35:155-35:185, 35:187-35:190, 35:192- 35:197, 35:199-35:206, 35:208:01-35:209, 35:211-35:216N, 35:218-35:225, 35:227-35:253, 35:255-35:269, 35:272- 35:273, 35:275-35:310, 35:312-35:316, 35:319-35:360, 35:362-35:449, 35:451-35:476, 35:478-35:533, 35:535- 35:540, 35:542N-35:560, 35:561 ^w , 35:562-35:568, 37:08, 37:67, 38:06-38:07, 39:19:01-39:19:02, 39:186, 44:06, 45:22, 51:01:01:01-51:24:05, 51:26-51:46, 51:48-51:103, 51:105- 51:111, 51:113-51:117, 51:119-51:146, 51:148- 51:230:01:03, 51:232:01-51:266, 51:268-51:283, 51:285- 51:309, 51:311-51:376, 53:01:01:01-53:16, 53:18-53:27,



0197

For *In Vitro* Diagnostic Use

MA123 v02 SSP PI Template

Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

			53:29, 53:31-53:37, 53:40-53:72, 55:130 ^w , 56:06 ^w , 56:45 ^w , 78:01:01:02-78:04, 78:07-78:10
69 ^{5,9}	115 bp, 150 bp	1070 bp	15, 18,22, 27, 35, 39, 46, 55, 62, 72, 73, 75, 76 *07:100, 07:428, 13:31, 13:41, 13:131, 13:145, 13:164, 15:01:01:01-15:01:04, 15:01:06-15:01:16, 15:01:18- 15:01:30, 15:01:32-15:01:78, 15:04:01:01-15:04:04, 15:07:01:01-15:08:02, 15:11:01-15:12:03, 15:14- 15:15:01:02, 15:19, 15:24:01-15:24:02, 15:26N-15:28, 15:30:01:01-15:30:01:03, 15:32:01-15:32:02, 15:34:01- 15:35:01:02, 15:38:01-15:38:02, 15:43, 15:45-15:46, 15:50, 15:53-15:54, 15:56-15:58, 15:60, 15:63, 15:66, 15:68, 15:70- 15:71, 15:73:01:01-15:73:01:02, 15:75-15:77, 15:79N, 15:81- 15:82, 15:85, 15:87, 15:92, 15:94N, 15:96-15:97, 15:101- 15:102, 15:104-15:105, 15:109-15:111N, 15:113, 15:117- 15:118, 15:120, 15:122, 15:125:01-15:126, 15:128-15:129, 15:135, 15:137, 15:140, 15:142-15:149N, 15:152, 15:154, 15:157, 15:159-15:160:01:02, 15:163-15:167, 15:169, 15:171-15:172, 15:174-15:175, 15:178, 15:180-15:184, 15:187, 15:189:01-15:193, 15:201-15:203, 15:205-15:207, 15:209N, 15:211-15:212, 15:215, 15:217, 15:225, 15:227- 15:228, 15:231-15:234, 15:236, 15:239, 15:241, 15:244- 15:247, 15:251:01-15:251:03, 15:256-15:262:02N, 15:264, 15:267, 15:269-15:270, 15:272N, 15:276-15:280, 15:284- 15:286, 15:296, 15:298-15:299, 15:303-15:305, 15:309, 15:315-15:318, 15:320-15:321Q, 15:324, 15:326-15:327, 15:331, 15:333-15:334, 15:336, 15:340, 15:342-15:344, 15:346-15:347, 15:349:01-15:350, 15:353, 15:355, 15:359, 15:363:01-15:363:02, 15:365, 15:367-15:368, 15:370- 15:373, 15:375N, 15:377Q, 15:381, 15:385-15:387, 15:390- 15:392, 15:398, 15:400N-15:401, 15:405-15:406, 15:409- 15:410, 15:412-15:413, 15:415-15:417, 15:419, 15:422, 15:428-15:432, 15:440-15:443, 15:450, 15:452-15:453, 15:456-15:457, 15:459, 15:461, 15:464-15:467, 15:469- 15:474, 15:476-15:484, 15:486, 15:488-15:490, 15:492- 15:493, 15:496N, 15:498-15:499, 15:501, 15:507, 15:510, 15:512-15:515, 15:520Q-15:522, 15:525, 15:528N-15:529, 15:531, 15:533-15:536, 15:538, 15:540N-15:542, 15:545, 15:547, 15:549N, 15:551, 15:554, 15:556, 15:558-15:559, 15:561-15:562N, 15:565, 15:568, 15:572, 15:576, 15:579- 15:581, 15:583, 15:585-15:586, 15:590-15:591, 15:593- 15:595, 15:597-15:598, 15:600-15:601, 15:604N-15:606, 15:608-15:609, 15:618, 15:621-15:622, 15:624, 15:631- 15:632, 15:634-15:635, 15:637-15:641, 15:647, 18:19, 18:164, 27:25, 27:75, 27:115, 27:161, 27:245, 35:14:01- 35:14:03, 35:43:01, 35:43:03-35:44, 35:62, 35:67, 35:79, 35:86, 35:102, 35:117, 35:135, 35:163, 35:185, 35:213, 35:265, 35:343, 35:390N, 35:413, 35:475-35:476, 35:512, 35:549, 37:60, 37:66, 39:18, 39:36, 44:146, 44:275, 44:299, 46:01:01:01-46:02, 46:04-46:05, 46:07N, 46:09-46:10, 46:12, 46:14-46:17, 46:20, 46:22-46:24, 46:26-46:32:02, 46:34- 46:42, 46:44, 46:46-46:53, 46:55N-46:57, 46:60-46:61, 46:63-46:82, 46:84-46:92, 46:94-46:98N, 50:47, 51:61:01- 51:61:02, 51:220, 52:21:01-52:21:02, 52:97, 53:70, 54:06, 54:33, 54:40, 55:21, 56:03, 56:43, 57:147, 57:161, 58:78,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

70 ^{8,15}	360 bp	1070 bp	Bw4	<p>73:01:01:01-73:03, A*25:66, A*26:68, A*26:100, A*26:200, A*68:56:01-68:56:02, C*02:56, C*06:20, C*12:50 *07:27, 07:36, 07:38, 07:81, 07:149, 07:180, 07:219, 07:236, 07:273, 07:453, 08:02-08:03, 08:52, 08:78, 08:117, 08:126, 08:255, 08:271, 13:01:01:01-13:04, 13:06-13:08, 13:10-13:23, 13:25-13:38, 13:40-13:66, 13:68-13:174, 14:93, 15:13:01-15:13:03, 15:16:01:01-15:17:07, 15:23-15:24:02, 15:36, 15:43, 15:67, 15:87, 15:89:01-15:89:02, 15:95, 15:115, 15:157, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:222, 15:230, 15:254, 15:256, 15:268, 15:273, 15:339, 15:356, 15:361-15:362, 15:396, 15:403, 15:408, 15:411, 15:418, 15:423-15:424, 15:442, 15:446, 15:462, 15:500, 15:505, 15:516, 15:523, 15:532, 15:546Q, 15:550, 15:555, 15:575N, 15:586, 15:594, 15:602-15:603, 15:613, 15:619, 15:642, 15:644, 18:09, 18:54, 18:67, 18:136, 18:194:01:01-18:195, 27:01-27:05:36, 27:05:38-27:07:06, 27:09-27:11, 27:13:01-27:17, 27:19:01:01-27:21:02, 27:23-27:25, 27:27-27:32, 27:34-27:39, 27:41, 27:43, 27:45-27:88, 27:90:01-27:152, 27:154-27:168, 27:170-27:203, 27:205-27:214, 27:216-27:241, 27:243N-27:259, 27:261N-27:264, 35:36, 35:329, 35:471, 37:01:01:01-37:04:02, 37:06:01-37:10, 37:12-37:36, 37:38-37:66, 37:68-37:95, 37:97-37:107, 38:01:01:01-38:102, 38:160-38:173Q, 40:13, 40:19, 40:47, 40:96, 40:109-40:110, 40:117, 40:157, 40:188, 40:201, 40:207, 40:246, 40:292, 40:340, 40:394, 41:46, 41:60, 44:02:01:01-44:08, 44:10-44:45, 44:47-44:74, 44:76-44:89, 44:91-44:128:02, 44:130, 44:132-44:245:02, 44:247-44:253, 44:255-44:276, 44:278-44:344, 44:346-44:355, 44:357-44:363, 44:437, 44:439-44:457, 44:459-44:524, 44:526-44:546, 47:01:01:02-47:01:06, 47:03-47:12, 48:18, 49:01:01:01-49:14, 49:16-49:20, 49:22-49:78, 51:01:01:01-51:24:05, 51:26-51:46, 51:48-51:376, 52:01:01:01-52:110N, 53:01:01:01-53:72, 54:12, 55:90, 55:103, 55:121, 56:07, 56:21, 56:46, 56:65, 57:01:01:01-57:11, 57:13-57:164, 58:01:01:01-58:02:02, 58:04-58:29, 58:31N-58:63, 58:65-58:127, 58:129-58:139, 59:01:01:01-59:12</p>
71	350 bp	1070 bp	Bw6	<p>*07:02:01:01-07:02:05, 07:02:07-07:02:77, 07:02:79-07:26, 07:28-07:35, 07:37:01-07:37:02, 07:39-07:80, 07:82-07:148, 07:150-07:179, 07:181N-07:191, 07:193-07:218, 07:220-07:235, 07:237-07:272N, 07:274-07:452, 07:454-07:469N, 08:01:01:01-08:01:68, 08:04:01-08:05, 08:07-08:51, 08:53:01-08:77, 08:79:01-08:116, 08:118-08:125, 08:127-08:130, 08:132-08:135, 08:137-08:138, 08:140-08:152, 08:154-08:156, 08:158-08:254, 08:256-08:270Q, 08:272-08:303, 13:09, 13:39, 13:67, 14:01:01:01-14:92, 14:94-14:113N, 15:01:01:01-15:01:04, 15:01:06-15:01:70, 15:01:72-15:12:03, 15:14-15:15:01:02, 15:18:01:01-15:21:01:02, 15:25:01:01-15:35:01:02, 15:37-15:40:02, 15:42, 15:44-15:58, 15:60-15:66, 15:68-15:86, 15:88, 15:90-15:94N, 15:96-15:99, 15:101-15:114, 15:116-15:129, 15:131-15:156, 15:158-15:161, 15:163-15:167, 15:169-15:176, 15:178-15:195, 15:197:01-15:207, 15:209N-15:215, 15:217-15:221, 15:223-15:229, 15:231-15:253, 15:255, 15:257-15:267, 15:269-15:272N, 15:274-15:276, 15:278-15:338,</p>



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				15:340-15:355, 15:357-15:360, 15:363:01-15:376, 15:378-15:395, 15:397-15:402, 15:404-15:407, 15:409-15:410, 15:412-15:417, 15:419-15:422, 15:425-15:429, 15:431-15:441, 15:443-15:445, 15:447-15:461, 15:463N-15:499, 15:501-15:504:02, 15:506-15:515, 15:517-15:522, 15:524-15:531, 15:533-15:545, 15:547-15:549N, 15:551-15:554, 15:556-15:574, 15:576-15:585, 15:587-15:593, 15:595-15:601, 15:604N-15:612, 15:614:01:01-15:618, 15:620-15:641, 15:643, 15:645-15:647, 18:01:01:01-18:08, 18:10-18:15, 18:17N-18:53, 18:55-18:66, 18:68-18:135:01:02, 18:137-18:178, 18:180-18:193, 18:196-18:226, 27:08, 27:12:01:01-27:12:01:03, 27:18, 27:26, 27:33, 27:40, 27:42, 27:44, 27:89, 27:153, 27:165 ^w , 27:204:01:01-27:204:01:02, 27:242, 35:01:01:01-35:01:38, 35:01:40-35:02:08, 35:02:10-35:43:01, 35:43:03-35:72, 35:74-35:221, 35:223-35:247:02, 35:249-35:292, 35:294-35:328, 35:330-35:423, 35:425-35:470, 35:472-35:521, 35:523-35:524, 35:526-35:568, 37:05, 37:11, 37:14 ^w , 37:37, 37:67, 39:01:01:01-39:01:01:29, 39:01:03:01-39:01:14, 39:01:16-39:20, 39:22-39:103, 39:105-39:193, 40:01:01-40:12:02, 40:14:01-40:16:01:02, 40:18, 40:20:01:01-40:40, 40:42-40:46, 40:48-40:75, 40:77-40:95, 40:97-40:108, 40:111-40:116, 40:118N-40:156, 40:158-40:187, 40:189-40:200, 40:202-40:291N, 40:293-40:339, 40:341-40:393, 40:395-40:399N, 40:401-40:526, 41:01:01:01-41:25, 41:27-41:45N, 41:47-41:59, 41:61-41:79, 42:01:01:01-42:02:02, 42:04-42:21, 42:23-42:32, 44:09, 44:46:01-44:46:02, 44:75, 44:90, 44:129, 44:131, 44:254, 44:277, 44:356, 44:438N, 45:01:01:01-45:01:08, 45:01:10-45:15, 45:17-45:25, 45:27-45:30, 46:01:01:01-46:98N, 47:02, 47:03 ^w , 48:01:01:01-48:17, 48:19-48:55, 50:01:01:01-50:02:01:02, 50:04:01-50:22, 50:31-50:79, 54:01:01:01-54:11, 54:13-54:46, 55:01:01:01-55:05, 55:07-55:89N, 55:91-55:102, 55:104-55:120, 55:122-55:132N, 56:01:01:01-56:06, 56:08-56:20:02, 56:22-56:45, 56:47-56:64, 56:66-56:89, 57:12, 58:64, 67:01:01-67:07, 73:01:01:01-73:03, 78:01:01:02-78:10, 81:01:01:01-81:10, 82:01:01:01-82:04, 83:01
72	285 bp	1070 bp	7, 703, 8, 35, 40, 41, 42, 48, 53, 60, 61, 81	*07:02:01:01-07:02:05, 07:02:07-07:02:38, 07:02:40-07:02:84, 07:02:86-07:06:05, 07:08:01-07:18:02, 07:20-07:32, 07:34-07:39, 07:41-07:52, 07:54-07:59, 07:61-07:76, 07:79-07:83, 07:85:01-07:99, 07:101-07:121, 07:123-07:138, 07:140-07:164, 07:166-07:192, 07:194-07:195, 07:199-07:213, 07:215-07:234, 07:236-07:238, 07:240-07:254, 07:256-07:342, 07:344-07:358, 07:360Q-07:385, 07:387-07:428, 07:430-07:441, 07:443-07:467Q, 07:469N, 08:01:01:01-08:01:20, 08:01:22-08:01:24, 08:01:26-08:05, 08:07-08:08N, 08:10-08:11, 08:13-08:15, 08:17-08:30N, 08:32-08:83, 08:85-08:88, 08:90-08:102, 08:104-08:106, 08:108-08:112, 08:114-08:115, 08:117-08:133, 08:135-08:168, 08:170-08:174, 08:176-08:183, 08:185-08:231, 08:233-08:253, 08:256-08:273, 08:275-08:290, 08:292-08:303, 13:174, 15:297, 27:164, 35:66, 35:87, 35:530, 37:07, 40:15-40:16:01:02, 40:30-40:32, 40:34, 40:45, 40:59, 40:80, 40:98, 40:137, 40:160:01-40:160:02, 41:02:01:01-41:02:11,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

41:04, 41:10-41:11, 41:13, 41:18-41:19, 41:23-41:24, 41:27,
41:31, 41:36, 41:38-41:43, 41:45N-41:49, 41:51-41:52,
41:55, 41:58, 41:65, 41:68-41:69, 41:71-41:72, 41:74,
42:01:01:01-42:02:02, 42:05:01-42:07, 42:09-42:13, 42:15-
42:32, 44:150, 44:166, 44:342, 48:01:01:01-48:01:11, 48:05-
48:12, 48:14-48:20, 48:22, 48:27-48:39, 48:41-48:50, 48:52,
48:54-48:55, 53:15, 81:01:01:01-81:10, **C*03:260**

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-B low resolution SSP typings. When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length. Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers. PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings. PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings. Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 25 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³The serological reactivity of all HLA-B alleles is not known. In this table we use the expert-assigned serological grouping in Tissue Antigens (2009) 73:95-170 and the serological grouping of the sequence-defined allele.

⁴For several HLA Class I alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁵HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁶Primer mix 28 may give rise to a lower yield of HLA-specific PCR product than the other HLA-B low resolution primer mixes in B*40, B*41, B*45, B*49 and B*50 alleles.

⁷Primer mixes 28, 29, 33, 42, 53 and 63 may give rise to a lower yield of HLA-specific PCR product than the other HLA-B low resolution primer mixes.

⁸Primer mixes 32, 33, 38, 42, 59, 62, 65, 66 and 70 have tendencies of unspecific amplifications, most pronounced in primer mix 33.

⁹Primer mixes 25, 27, 43, 44, 50, 55 and 69 have a tendency giving rise to primer oligomer formation.

¹⁰The B*57 and B*58 alleles might be faintly amplified by primer mix 33.

¹¹Many C*17 alleles may be faintly amplified by primer mix 49.

¹²Primer mix 67 may generate a false positive band of about 800 base pairs. This band should be disregarded when interpreting HLA-B low resolution typings.

¹³Primer mix 63 may give rise to a lower yield of B*54 alleles than the other B low primer mixes.

¹⁴Primer mix 63 may amplify B*44:05 samples.

¹⁵The Bw4-associated HLA-A specificities A9, A23, A24, A2403, A25 and A32 are not amplified by the primer pair in primer mix 70.

Abbreviations

w: might be weakly amplified.

?: nucleotide sequence information not available for the primer matching sequence.



0197

For *In Vitro* Diagnostic Use

MA123 v02 SSP PI Template

Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

SPECIFICITY TABLE

Specificities and sizes of the PCR products of the 23+1 primer mixes used for DR low resolution SSP typing

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	DR serology ³	Amplified HLA-DRB alleles ^{4,5}
73⁷	210 bp, 235 bp, 260 bp	515 bp	1	*01:01:01:01-01:02:13, 01:04-01:38, 01:40N-01:98, 01:100-01:101, 01:103-01:111, 01:113-01:126, 01:128-01:139
74	200 bp	430 bp	1, 103	*01:03:01:01-01:03:05, 01:39N, 01:42, 01:61, 01:89, 01:99, 01:102, 01:112, 01:127
75^{9,10}	200 bp, 230 bp	430 bp	2, 15	*15:01:01:01-15:211, 16:66
76	210 bp	430 bp	11, 16	*11:30, 11:288, 12:57, 16:01:01:01-16:05:02, 16:07-16:73, DRB3*02:111
77^{5,6}	75 bp, 205 bp	515 bp	3, 11, 13, 17, 18	*03:01:01:01-03:11:01, 03:11:03-03:125, 03:127-03:203, 11:07:01-11:07:02, 11:27:02-11:27:03, 11:84:02-11:84:03, 11:103:01-11:103:02, 11:105, 11:107, 11:124:02-11:125, 11:136, 11:173, 13:33:01, 13:61:02, 13:94:01, 13:96:01, 14:38:02, 14:104, 14:111, 14:171, 15:25, 15:37:01, 15:100, 15:104:01
78^{5,6}	75 bp, 210 bp	430 bp	3, 6, 11, 13, 14, 17	*03:01:01:01-03:01:23, 03:01:25-03:01:41, 03:04:01-03:06, 03:08-03:11:01, 03:11:03-03:16, 03:18-03:20, 03:22-03:23, 03:25:01-03:26, 03:28, 03:30-03:31, 03:33-03:34, 03:36-03:37, 03:43-03:48, 03:50-03:52:02, 03:54-03:68N, 03:70-03:73, 03:75-03:86, 03:89, 03:91-03:93, 03:95-03:96, 03:98-03:100:02, 03:104, 03:106-03:114, 03:116-03:118, 03:121-03:130, 03:132-03:149, 03:151-03:153, 03:156N-03:166, 03:168-03:170, 03:172, 03:174N-03:175, 03:177, 03:180-03:184, 03:186-03:187, 03:189N, 03:191-03:203, 08:40, 08:92, 11:02:01:01-11:03:04, 11:11:01, 11:11:03, 11:14:01:01-11:14:02, 11:16:01:01-11:16:01:02, 11:20-11:21, 11:36, 11:40-11:41, 11:48, 11:59, 11:63:01-11:63:02, 11:65:01-11:65:02, 11:68, 11:70, 11:73, 11:76, 11:79-11:80, 11:83, 11:85-11:87, 11:93, 11:118:01-11:118:02, 11:122, 11:124:01-11:124:02, 11:127, 11:131-11:132, 11:135, 11:138-11:139, 11:142, 11:151, 11:153, 11:161, 11:168, 11:171, 11:176, 11:182, 11:184, 11:213, 11:229, 11:231, 11:233, 11:237, 11:241-11:242, 11:251, 11:256, 11:260, 11:270, 11:276, 11:290, 11:293, 11:301N, 12:73, 13:01:01:01-13:02:01:14, 13:02:02 [?] , 13:02:03-13:04, 13:08, 13:10, 13:15-13:17, 13:19-13:20:01:03, 13:22:01:01-13:24:01:03, 13:27:01-13:29:01:02, 13:31-13:41, 13:43, 13:45, 13:48, 13:51-13:54, 13:57, 13:59, 13:61:01-13:61:02, 13:63-13:66:02, 13:68-13:76, 13:78-13:81, 13:83-13:85, 13:87-13:99, 13:101-13:102, 13:104-13:107, 13:109, 13:110 [?] , 13:111-13:117:01:02, 13:120-13:131, 13:133, 13:135, 13:137N-13:145, 13:147-13:149, 13:151-13:153, 13:155, 13:159, 13:162, 13:165-13:166, 13:167 [?] , 13:168, 13:170-13:180, 13:182-13:188, 13:190-13:191, 13:193-13:194, 13:196, 13:198, 13:200N-13:202, 13:204-13:205, 13:207-13:209, 13:212-13:219, 13:221, 13:223,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				13:225-13:227, 13:230, 13:233-13:234, 13:236-13:245, 13:247-13:283, 13:285, 13:287-13:288, 13:289N [?] , 13:290-13:294, 13:296-13:301, 13:303-13:313, 13:315-13:320, 13:322N-13:329N, 13:332, 14:16, 14:19, 14:21, 14:45, 14:82, 14:95, 14:109, 14:120, 14:122, 14:132, 14:178-14:179
79⁶	85 bp, 210 bp	430 bp	3, 6, 11, 13, 14, 1403, 18	*03:02:01:01-03:03, 03:27, 03:29, 03:38, 03:53, 03:74, 03:88, 03:90, 03:102-03:103, 03:115, 03:119, 03:126, 03:131, 03:154, 03:167, 03:176-03:177, 03:179, 11:13:01 ^w -11:13:02 ^w , 11:26, 11:34, 11:267, 13:15, 13:19, 13:26:01-13:26:02, 13:44, 13:53, 13:57, 13:85-13:86, 13:104, 13:193, 13:198, 13:206, 13:235, 14:02:01:01-14:03:02, 14:06:01-14:06:06, 14:09, 14:12:01-14:13, 14:17-14:21, 14:24, 14:27:01-14:27:02, 14:29-14:30, 14:32:01 ^w -14:32:02 ^w , 14:33, 14:40-14:41, 14:47-14:49, 14:51, 14:63, 14:65 ^w , 14:67, 14:77-14:78, 14:80-14:81, 14:83, 14:85, 14:89, 14:94, 14:98, 14:102, 14:106, 14:108-14:109, 14:115, 14:119, 14:121, 14:135, 14:146, 14:154, 14:156, 14:159-14:160, 14:165, 14:170, 14:174, 14:176-14:177, 14:179, 14:181, 14:183, 14:194-14:195N, 14:198, 14:200, 14:209, 14:212, 14:218, 14:221-14:223, 14:226, 14:239, 14:242
80^{6,7,9}	100 bp, 175 bp	430 bp	3, 4	*04:01:01:01-04:05:11, 04:05:13-04:353
81⁹	170 bp, 210 bp, 235 bp	430 bp	7, 13, 14	*07:01:01:01-07:01:01:02, 07:01:01:04-07:01:30, 07:03-07:141, 08:67, 08:92, 08:108, 10:11, 11:229, 11:270, 12:22, 12:73, 12:88, 13:17, 13:116, 13:175, 13:214, 14:50:01
82⁶	85 bp, 120 bp, 215 bp, 250 bp	515 bp	7, 8, 14	*07:06, 08:01:01:01-08:01:02, 08:01:04-08:19, 08:21-08:30:03, 08:32-08:56, 08:58-08:66, 08:68:01-08:79, 08:81-08:97, 08:99-08:117, 11:67, 11:193:01-11:193:02, 12:16:01-12:16:03, 12:22, 12:39:01-12:39:02, 13:171:01, 13:197, 14:15, 14:68:01-14:68:02, 14:93
83⁶	85 bp, 100 bp, 135 bp, 160 bp, 190 bp	430 bp	3, 9, 11	*03:08, 03:65, 03:140, 09:01:02:01-09:52N, 11:07:01-11:07:02, 11:103:01-11:103:02, 11:105, 11:107, 11:125, 11:173, DRB3*01:42, DRB5*02:31, DRB5*02:34
84^{6,9}	70 bp, 175 bp	430 bp	9, 10,11,13	*03:76, 03:125, 03:174N, 10:01:01:01, 10:01:01:03-10:45, 11:59, 11:80, 11:83, 11:87, 11:135, 11:142, 11:182, 13:27:01-13:27:02, 13:41, 13:71, 13:129, 13:176, 13:277, 13:326, 14:179
85⁶	100 bp, 175 bp	430 bp	3, 8, 11, 13, 14	*03:08, 03:65, 03:140, 08:31, 08:41, 08:63, 08:75, 11:01:01:01, 11:01:01:03-11:11:01, 11:11:03-11:70, 11:72-11:263, 11:265-11:292, 11:294N-11:307
86⁶	90 bp, 110 bp	430 bp	12, 14	*08:32, 08:53, 08:87, 12:01:01:01-12:56, 12:58-12:100, 14:52
87⁷	220 bp	430 bp	6, 8, 11, 13, 14, 1403	*03:76, 03:125, 03:174N, 08:20-08:21, 11:01:01:01, 11:01:01:03-11:01:27, 11:01:28 [?] , 11:01:29-11:04:22, 11:06:01-11:06:03, 11:08:01-11:11:01, 11:11:03-11:12:03, 11:14:01:01-11:16:01:02, 11:18-11:21, 11:23:01-11:25, 11:27:01-11:33, 11:35-11:51, 11:54:01-11:54:02, 11:56-11:66:02, 11:68, 11:70, 11:72-11:81, 11:83-11:88, 11:90-11:97, 11:99-11:102:02, 11:106, 11:108-11:124:02, 11:126-11:135, 11:137-11:142, 11:144-11:158, 11:160-11:172,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				11:174-11:192, 11:194:01-11:201, 11:203-11:208, 11:210-11:228, 11:230-11:260, 11:262-11:265, 11:268-11:269, 11:271-11:282, 11:284-11:307, 12:57, 13:01:01:01-13:02:01:14, 13:02:03-13:08, 13:10-13:16:01:02, 13:18-13:43, 13:45-13:85, 13:87-13:115, 13:117:01:01-13:128, 13:130-13:145, 13:147-13:174, 13:176-13:182, 13:184-13:191, 13:194-13:205, 13:207-13:209, 13:212-13:213, 13:215-13:221, 13:223-13:247, 13:249N-13:288, 13:290-13:298N, 13:300-13:313, 13:315-13:325, 13:327-13:330, 13:332, 14:03:01-14:03:02, 14:12:01-14:12:02, 14:16, 14:19, 14:21-14:22, 14:25:01-14:25:02, 14:27:01-14:27:02, 14:40, 14:53, 14:63, 14:67, 14:69, 14:74, 14:77-14:78, 14:84-14:85, 14:98, 14:102, 14:105, 14:109, 14:115-14:116, 14:128, 14:135, 14:144, 14:156, 14:173-14:174, 14:181, 14:218, 14:222N-14:223, 14:235, 14:239, DRB3*02:27, DRB3*02:111, DRB3*02:131, DRB3*02:160
88 ⁷	195 bp, 225 bp	430 bp	6, 8, 11, 12, 13, 14	*03:125, 08:01:01:01-08:01:02, 08:01:04-08:01:05, 08:01:07-08:02:05, 08:04:01:01-08:09:02, 08:11:01:01-08:11:01:02, 08:16-08:17, 08:20-08:22, 08:24:01-08:24:02, 08:26, 08:28, 08:31, 08:39, 08:41-08:44, 08:50, 08:52, 08:54-08:55, 08:57, 08:59, 08:64, 08:67, 08:70, 08:72-08:73, 08:75, 08:77-08:80, 08:82-08:90, 08:92-08:94, 08:96, 08:98-08:101, 08:103, 08:105-08:109, 08:111-08:114, 08:116-08:117, 11:01:01:01, 11:01:01:03-11:01:17, 11:01:20-11:04:20, 11:04:22-11:06:03, 11:09:01:01-11:11:01, 11:11:03-11:12:03, 11:14:01:01-11:16:01:02, 11:20-11:21, 11:23:01-11:25, 11:27:01-11:30, 11:32-11:33, 11:35-11:41, 11:43:01:01-11:44, 11:46:01-11:51, 11:54:01-11:56, 11:58:01-11:63:02, 11:65:01-11:70, 11:72, 11:74:01-11:78, 11:80-11:88, 11:90-11:97, 11:99-11:102:02, 11:106, 11:108-11:118:02, 11:120-11:124:02, 11:126-11:129, 11:133-11:135, 11:137-11:142, 11:144-11:152, 11:154-11:158, 11:160-11:172, 11:174-11:183, 11:185-11:190, 11:192-11:201, 11:203-11:247, 11:249-11:260, 11:262-11:265, 11:269-11:279, 11:281-11:282, 11:284-11:285, 11:287N-11:307, 12:02:01:01-12:02:13, 12:13, 12:15-12:16:03, 12:18-12:21, 12:23, 12:26-12:27, 12:31N-12:33, 12:37, 12:42-12:45, 12:49-12:52, 12:55-12:57, 12:60N-12:61, 12:64-12:65, 12:67, 12:69, 12:71, 12:73-12:76, 12:78, 12:80-12:81, 12:85, 12:92-12:94, 12:97-12:98N, 12:100, 13:01:01:01-13:02:01:14, 13:02:03-13:02:22, 13:04-13:05:03, 13:07:01-13:09, 13:11:01-13:11:02, 13:14:01-13:24:01:03, 13:26:01-13:29:01:02, 13:31-13:32, 13:34-13:36, 13:38-13:43, 13:45-13:55, 13:57, 13:59, 13:61:01-13:65, 13:67-13:76, 13:78-13:80, 13:83-13:84, 13:87, 13:91-13:93, 13:96:01-13:100, 13:102-13:109, 13:111-13:114, 13:116-13:117:01:02, 13:121, 13:123-13:132, 13:135-13:136, 13:138-13:150, 13:153, 13:155, 13:158:01-13:160, 13:162, 13:164-13:166, 13:168-13:169, 13:171:01-13:171:02, 13:173, 13:175, 13:177, 13:179, 13:182-13:187, 13:189-13:192, 13:195, 13:197-13:205, 13:207-13:215, 13:217-13:218, 13:221-13:223, 13:225-13:226, 13:229-13:230, 13:232-13:234, 13:236-13:252N, 13:254-13:288, 13:290-13:294, 13:296-13:298N, 13:300-13:316, 13:318, 13:320, 13:322N-13:328,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				13:331-13:332, 14:15-14:16, 14:22, 14:24-14:25:02, 14:27:01-14:27:02, 14:37, 14:53, 14:73, 14:105, 14:128, 14:226, DRB3*02:111, DRB3*02:160
89⁵	175 bp, 240 bp	430 bp	3, 4, 6, 8, 11, 13, 14, 1403,140 4, 17, 18	*03:01:01:01-03:01:05, 03:01:07-03:01:08, 03:01:10-03:07:02, 03:09, 03:11:01, 03:11:03-03:41:02, 03:43-03:45, 03:47-03:63, 03:66-03:68N, 03:70-03:86, 03:88-03:91, 03:93-03:110, 03:112-03:129, 03:131-03:139, 03:141-03:155, 03:157-03:161, 03:163-03:172, 03:174N-03:188, 03:191-03:203, 08:09:02, 08:20, 11:03:04, 11:04:17, 11:13:01-11:13:02, 11:17, 11:149, 11:202, 11:264, 12:20, 13:01:01:01-13:16:01:02, 13:18-13:42:02, 13:44, 13:46-13:66:02, 13:68-13:102, 13:104-13:115, 13:117:01:01-13:121, 13:123-13:158:02, 13:161-13:164, 13:166-13:170, 13:171:02-13:174, 13:176-13:178, 13:180-13:190, 13:192-13:194, 13:196, 13:198-13:213, 13:215-13:216, 13:218-13:227, 13:229-13:247, 13:249N-13:298N, 13:300-13:332, 14:01:01, 14:01:02 [?] -14:01:04 [?] , 14:01:05-14:04:01:03, 14:04:02 [?] -14:04:03 [?] , 14:04:04-14:04:06, 14:04:09-14:21, 14:22 [?] , 14:23:01:01-14:23:01:03, 14:23:02 [?] , 14:23:03-14:24, 14:25:01 [?] -14:26 [?] , 14:27:01-14:27:02, 14:28 [?] , 14:29-14:30, 14:31 [?] -14:32:03 [?] , 14:33, 14:34 [?] -14:35 [?] , 14:36-14:37, 14:38:01 [?] -14:39 [?] , 14:40-14:45, 14:47-14:48, 14:49 [?] -14:50:02 [?] , 14:51, 14:52 [?] , 14:54:01:01-14:54:06, 14:54:08-14:54:12, 14:55 [?] , 14:56-14:57, 14:58 [?] , 14:59, 14:60 [?] -14:62 [?] , 14:63-14:65, 14:67, 14:68:01 [?] -14:69 [?] , 14:70, 14:71 [?] -14:76 [?] , 14:77-14:78, 14:79 [?] , 14:80-14:85, 14:86 [?] -14:88 [?] , 14:89, 14:90 [?] , 14:91, 14:92N [?] -14:93 [?] , 14:94-14:96, 14:97 [?] , 14:98, 14:99 [?] , 14:100, 14:101 [?] , 14:102-14:103, 14:104 [?] -14:105 [?] , 14:106, 14:107 [?] , 14:108-14:109, 14:110 [?] , 14:111, 14:112 [?] -14:114 [?] , 14:115-14:116, 14:117 [?] -14:120 [?] , 14:121, 14:122 [?] , 14:123, 14:124 [?] -14:126:02 [?] , 14:127:01-14:127:02, 14:128 [?] -14:133 [?] , 14:134-14:137N, 14:138 [?] -14:140 [?] , 14:141, 14:142 [?] -14:143 [?] , 14:144, 14:145 [?] -14:151 [?] , 14:152N, 14:153 [?] , 14:154-14:156, 14:157 [?] -14:158 [?] , 14:159-14:162, 14:163 [?] -14:164 [?] , 14:165, 14:166N [?] , 14:167, 14:168 [?] , 14:169-14:171, 14:172 [?] -14:173 [?] , 14:174, 14:175 [?] , 14:176-14:177, 14:178 [?] , 14:179, 14:180 [?] , 14:181, 14:182 [?] , 14:183, 14:184 [?] -14:186 [?] , 14:187-14:188N, 14:189 [?] -14:190 [?] , 14:191-14:199, 14:201-14:202, 14:203 [?] , 14:204-14:207, 14:209-14:219, 14:220 [?] , 14:221, 14:222N [?] , 14:223-14:224, 14:225 [?] , 14:226-14:228, 14:229 [?] , 14:230-14:247
90⁶	75 bp, 100 bp, 140 bp, 185 bp, 240 bp	430 bp	4, 6, 8,13, 14, 1404	*04:62, 04:69, 04:73:01-04:73:02, 04:105:01-04:105:02, 04:122, 04:146, 04:162, 04:251, 04:284:01-04:284:02, 07:128, 08:08, 11:69, 11:82, 11:88, 11:293, 13:45, 13:197, 13:314, 14:01:01-14:01:02, 14:01:04-14:01:05, 14:04:01:01-14:04:09, 14:07:01-14:07:03, 14:10:01:01-14:10:01:02, 14:16, 14:22, 14:25:01-14:26, 14:28, 14:31-14:32:03, 14:35, 14:37-14:39, 14:49-14:50:02, 14:53-14:54:01:15, 14:54:03-14:55, 14:57-14:58, 14:60-14:62, 14:68:01-14:71, 14:73-14:76, 14:79, 14:82, 14:86-14:88, 14:90, 14:93, 14:99, 14:101, 14:104-14:105, 14:107, 14:110-14:114, 14:117-14:120, 14:122, 14:124-14:125, 14:128-14:129, 14:131, 14:137N-14:140, 14:142-14:143, 14:145-14:147,



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

				14:149-14:153, 14:157-14:158, 14:163-14:164, 14:166N, 14:168 ^w , 14:171-14:173, 14:178, 14:180, 14:182, 14:184, 14:186, 14:189-14:192, 14:196-14:197N, 14:199, 14:201-14:205:01, 14:207-14:208, 14:210Q, 14:213-14:217, 14:219-14:220, 14:224-14:225, 14:227-14:235, 14:237-14:238, 14:240-14:241, 14:243-14:245, 14:247, 15:167, DRB4*01:03:01:02N, DRB4*01:03:01:13N, DRB4*01:14N
91 ^{6,7,8}	110 bp, 145 bp, 170 bp	430 bp	3, 4, 6, 9, 11, 13, 14, 1404	*03:10, 03:126, 03:177, 04:211, 09:01:02:01-09:01:05, 09:01:07-09:02:02, 09:04-09:22, 09:24-09:28, 09:30-09:52N, 11:13:01-11:13:02, 11:17, 11:52, 11:202, 12:66, 12:96, 13:43, 13:159, 13:171:01, 13:179, 13:191, 13:193, 13:217, 14:01:01-14:02:09, 14:04:01:01-14:11:01:02, 14:13-14:14, 14:16-14:18, 14:19 ^w , 14:20, 14:21 ^w , 14:22-14:23:04, 14:26, 14:28-14:36, 14:38:01-14:39, 14:41, 14:43-14:52, 14:54:01:01-14:54:03, 14:54:04 ^w , 14:54:05-14:57, 14:59-14:62, 14:64-14:65, 14:68:01-14:68:02, 14:70-14:76, 14:79-14:83, 14:86-14:88, 14:90-14:97, 14:99-14:101, 14:103-14:108, 14:109 ^w , 14:110-14:114, 14:117-14:127:02, 14:129-14:134, 14:137N-14:140, 14:142-14:143, 14:145-14:155, 14:157-14:160, 14:162-14:173, 14:175-14:180, 14:182, 14:184-14:185, 14:187-14:207, 14:209-14:217, 14:219-14:221, 14:224-14:225, 14:227-14:238, 14:240-14:247, 15:27, 15:34, 15:54, 15:66:01-15:66:02, 15:122, 15:158, 15:187
92 ^{6,7,9}	110 bp, 150 bp, 180 bp, 220 bp	430 bp	2, 3, 6, 8, 11, 13, 14, 1403, 1404, 16	*03:10, 03:157, 08:03:11 ^w , 08:09:01-08:09:02, 08:20-08:21, 08:32, 08:35, 08:36:02, 08:53, 08:68:01 ^w , 08:70 ^w , 08:87 ^w , 11:13:01-11:13:02, 11:17, 11:23:01-11:23:02, 11:25, 11:31, 11:45, 11:52, 11:55, 11:64, 11:89, 11:96, 11:119, 11:148, 11:159, 11:202, 11:266, 11:283, 11:293, 12:96, 13:13, 13:18, 13:43, 13:45, 13:47, 13:55, 13:119, 13:144, 13:146, 13:154, 13:156, 13:158:01-13:159, 13:164, 13:171:01, 13:179, 13:191, 13:197, 13:217, 13:232, 13:286, 13:331, 14:01:01-14:01:05, 14:03:01-14:05:05, 14:07:01-14:08, 14:10:01:01-14:12:02, 14:14-14:16, 14:18, 14:22-14:23:04, 14:25:01-14:28, 14:31-14:32:03, 14:34-14:36, 14:38:01-14:40, 14:42-14:45, 14:49-14:50:02, 14:53-14:65, 14:67-14:79, 14:81-14:82, 14:84-14:93, 14:95-14:97, 14:99-14:105, 14:107, 14:110-14:120, 14:122-14:140, 14:142-14:158, 14:160-14:164, 14:166N-14:169, 14:171-14:175, 14:178, 14:180-14:182, 14:184-14:193, 14:196-14:197N, 14:199-14:208, 14:210Q-14:211, 14:213-14:220, 14:222N-14:225, 14:227-14:241, 14:243-14:247, 15:21 ^w , 15:158 ^w , 15:187 ^w , 16:04:01 ^w -16:04:02 ^w , 16:18 ^w , DRB3*02:131, DRB4*01:31
93 ^{5,6,7}	80 bp, 110 bp, 160 bp, 190 bp, 240 bp	515 bp	12, 52	*12:05, 12:15, 12:57, 14:141, DRB3*01:01:02:01-01:110, DRB3*02:01-02:185, DRB3*03:01:01:01-03:21, DRB3*03:23-03:60
94 ⁷	130 bp, 210 bp	430 bp	53	*01:139, 10:19, 14:217, 15:141, DRB4*01:01:01:01, DRB4*01:01:01:03-01:164, DRB4*02:01N, DRB4*03:01N
95	180 bp, 215 bp	430 bp	51	DRB5*01:01:01:01-01:129, DRB5*02:02:01-02:36
96 ¹¹	-	-	-	Negative control



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of DR low resolution SSP subtypings. When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length. Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers, e.g. the primers in wells 75 and 90 to 92. PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings. PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings. Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 73 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases. In the presence of a specific amplification the intensity of the control band often decreases.

³The serological split of all DRB1 alleles is not known. In this table we use the expert-assigned serological grouping in Tissue Antigens (2009) 73:95-170 and the serological grouping of the sequence-defined allele.

⁴For several DRB1 alleles 1st and/or 3rd exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁵Due to sharing of sequence motifs in codon 38 and 47, DRB3*01:14 will also be amplified in primer mixes 77, 78 and 89, and DRB3*01:23 and 02:32 in mix 77, in addition to primer mix 93.

⁶HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁷Primer mixes 73, 80, 87, 88 and 91 to 94 may have a tendency giving rise to primer oligomer formation.

⁸Primer mix 91 has a tendency of primer oligomer formation and also has an intense primer cloud due to the high number of primers present in the primer mix.

⁹Primer mixes 75, 80, 81, 84 and 92 may have tendencies of unspecific amplifications, most pronounced in primer mix 75.

¹⁰In primer mix 75 the positive control band may be weaker than for other A-B-DR primer mixes.

¹¹Primer mix 96 contains a negative control, which will amplify the majority of HLA amplicons as well as the amplicons generated by the control primer pairs matching the human growth hormone gene. HLA-specific PCR product sizes range from 75 to 200 base pairs and the PCR product generated by the HGH positive control primer pair is 200 base pairs.

Abbreviations

w: might be weakly amplified.

?: nucleotide sequence information not available for the primer matching sequence.



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 6R6

Lot-specific information

HLA-A LOW PRIMER SPECIFICATION

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec.	120	135	205	190	160	135	175	165	80	80	80	125
PCR product	145	175	235		335	200	205	200		200	175	160
		225	215		505					240	500	190
			255									
			365									
			415									
Length of int.	800	800	1070	800	1070	800	1070	800	800	800	1070	800
pos. control ¹												
5'-primer(s) ²	98	48	359	98	144	176	98	98	257	98	301	103
	5'-CTT 3'	5'-gCT 3'	5'-gCA 3'	5'-CTA 3'	5'-gCC 3'	5'-gCA 3'	5'-CTC 3'	5'-CTA 3'	5'-Cgg 3'	5'-CTA 3'	5'-Cgg 3'	5'-CCT 3'
	103	78	363	413	317	368	368	102	259	261	302	415
	5'-CCT 3'	5'-TCT 3'	5'-ATA 3'	5'-CCg 3'	5'-gCT 3'	5'-gTT 3'	5'-gTT 3'	5'-ACA 3'	5'-AgA 3'	5'-AAC 3'	5'-ggA 3'	5'-ggT 3'
	123	106	363					413	261		385	423
	5'-AgT 3'	5'-CCA 3'	5'-ATA 3'					5'-CCg 3'	5'-AAC 3'		5'-ggC 3'	5'-gCT 3'
	363	144							266			
	5'-ATA 3'	5'-gCA 3'							5'-ACg 3'			
	363											
	5'-ATg 3'											
3'-primer(s) ³	203	240	527	256	265	270	259	259	302	256	341	257
	5'-TCT 3'	5'-ggA 3'	5'-CCA 3'	5'-CTg 3'	5'-CCC 3'	5'-ACA 3'	5'-gTT 3'	5'-gTT 3'	5'-ggC 3'	5'-CCA 3'	5'-CgT 3'	5'-gCA 3'
	545	292	527	559	368	521	502	259	303	299	521	506
	5'-AgA 3'	5'-gTg 3'	5'-CCT 3'	5'-CCg 3'	5'-CAA 3'	5'-ggg 3'	5'-CTT 3'	5'-gTT 3'	5'-AgA 3'	5'-TCg 3'	5'-ggg 3'	5'-TgT 3'
			527			526	539	538		300		538
			5'-CCT 3'			5'-CCg 3'	5'-TCT 3'	5'-CCA 3'		5'-TTT 3'		5'-CCg 3'
			555			534				308		559
			5'-CCA 3'			5'-CgT 3'				5'-TCT 3'		5'-CTC 3'
			555									559
			5'-gCA 3'									5'-CgT 3'
												559
												5'-CgC 3'
												563
												5'-CgA 3'
Well No.	1	2	3	4	5	6	7	8	9	10	11	12



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

Well No.	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec.	180	80	90	145	140	200	370	210	65	85	80	360
PCR product	225	115	135	165	180	390		240	240	240	160	445
		200	205	240	240				375	400	240	
		240		380	260				545		495	
		460		410								
Length of int. pos. control ¹	1070	1070	1070	1070	1070	1070	800	800	800	800	800	1070
5'-primer(s) ²	98	98	203	41	98	41	302	78	28	47	176	257
	5'-CTT 3'	5'-CAC 3'	5'-gAA 3'	5'-CTT 3'	5'-CTT 3'	5'-CTT 3'	5'-ggA 3'	5'-TCT 3'	5'-TCg 3'	5'-g.T 3'	5'-gCA 3'	5'-Cgg 3'
	414	219	362	160	180	98	302	106	261	78	261	341
	5'-CAg 3'	5'-gCA 3'	5'-ggT 3'	5'-ACg 3'	5'-TTT 3'	5'-CAC 3'	5'-ggA 3'	5'-CCA 3'	5'-AAC 3'	5'-TCT 3'	5'-AAC 3'	5'-ggC 3'
	423	238	363	355	203		317	2 nd	502	521	341	
	5'-gCT 3'	5'-AgA 3'	5'-ATA 3'	5'-CCg 3'	5'-gAA 3'		5'-gCg 3'	5'-CCT 3'	5'-CCC 3'	5'-ggT 3'	5'-ggA 3'	
		355	363		418			2 nd		527	355	
		5'-CCg 3'	5'-ATA 3'		5'-AgC 3'			5'-ggC 3'		5'-TgC 3'	5'-CCC 3'	
		489	363								362	
		5'-gCA 3'	5'-ATA 3'								5'-gAg 3'	
			369								362	
			5'-TAC 3'								5'-gAg 3'	
3'-primer(s) ³	282	180	299	238	290	256	397	265	97	265	292	418
	5'-gAC 3'	5'-TCA 3'	5'-CCA 3'	5'-CCT 3'	5'-CAA 3'	5'-CCC 3'	5'-gAg 3'	5'-CCC 3'	5'-ggT 3'	5'-CCC 3'	5'-gTg 3'	5'-gTC 3'
	282	257	411	238	299	259		282	355	282	292	
	5'-gAC 3'	5'-gCA 3'	5'-TCA 3'	5'-CCT 3'	5'-TCT 3'	5'-gTT 3'		5'-gAC 3'	5'-gAC 3'	5'-gAC 3'	5'-gTT 3'	
	559	418	526	243	317	261		282	524	282	299	
	5'-CCC 3'	5'-gTC 3'	5'-CCA 3'	5'-TCA 3'	5'-ggA 3'	5'-gTg 3'		5'-gAC 3'	5'-CAT 3'	5'-gAC 3'	5'-TCT 3'	
	560	555		265	555	265		502		570	555	
	5'-ACC 3'	5'-CCA 3'		5'-CCC 3'	5'-CCA 3'	5'-CCC 3'		5'-CTT 3'		5'-CAC 3'	5'-CCA 3'	
				282		265		506				
				5'-gAC 3'		5'-CCC 3'		5'-TgT 3'				
				555		265						
				5'-CCA 3'		5'-CCA 3'						
Well No.	13	14	15	16	17	18	19	20	21	22	23	24

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 6R6

Lot-specific information

HLA-B LOW PRIMER SPECIFICATION

Well No.	25	26	27	28	29	30	31	32	33	34	35	36
Length of spec.	110	215	140	105	185	190	290	165	165	180	100	280
PCR product			235	140	235			225	190		195	
				260				285	390			
								330				
Length of int.	800	1070	1070	800	800	800	1070	1070	1070	800	1070	1070
pos. control ¹												
5'-primer(s) ²	527	97	209	103	103	103	45	45	45	420	206	419
	5'-TgA 3'	5'-TCg 3'	5'-ggC 3'	5'-CCg 3'	5'-CCT 3'	5'-CCg 3'	5'-ggA 3'	5'-ggA 3'	5'-ggA 3'	5'-TTA 3'	5'-AgA 3'	5'-gTC 3'
	539		363	103	103	418		357	435		412	
	5'-gCg 3'		5'-AgC 3'	5'-CCg 3'	5'-CCT 3'	5'-Agg 3'		5'-Tgg 3'	5'-AAA 3'		5'-ATA 3'	
	540		363	361	363			412			416	
	5'-gAC 3'		5'-AgC 3'	5'-AgT 3'	5'-AAT 3'			5'-ATA 3'			5'-CCg 3'	
				363							420	
				5'-AgC 3'							5'-TTA 3'	
											420	
											5'-TTC 3'	
3'-primer(s) ³	605	272	309	167	246	246	165	206	266	559	256	3 rd I
	5'-gCT 3'	5'-Tgg 3'	5'-gTg 3'	5'-ACA 3'	5'-TAT 3'	5'-TAT 3'	5'-Tgg 3'	5'-CCT 3'	5'-TCC 3'	5'-CTC 3'	5'-CCC 3'	5'-TAT 3'
		272	312	193	246	572		538	559		272	
		5'-TgA 3'	5'-gCC 3'	5'-CgT 3'	5'-TAT 3'	5'-gCC 3'		5'-gTC 3'	5'-CAg 3'		5'-Tgg 3'	
		272	559	193	559			603	583		572	
		5'-TgT 3'	5'-CAg 3'	5'-CgT 3'	5'-CTC 3'			5'-gTg 3'	5'-gTg 3'		5'-gCg 3'	
				463								
				5'-gCT 3'								
				583								
				5'-gTg 3'								
Well No.	25	26	27	28	29	30	31	32	33	34	35	36



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

Well No.	37	38	39	40	41	42	43	44	45	46	47	48
Length of spec.	195	105	115	80	100	140	55	210	170	110	395	160
PCR product		395	150		150		245					425
		435					390					
							415					
Length of int. pos. control ¹	1070	800	1070	1070	800	1070	1070	800	800	1070	1070	1070
5'-primer(s) ²	142	45	161	167	349	206	45	142	407	246	44	44
	5'-TCT 3'	5'-ggA 3'	5'-Cgg 3'	5'-gCT 3'	5'-CTT 3'	5'-gAC 3'	5'-ggA 3'	5'-TCT 3'	5'-gg 3'	5'-gAA 3'	5'-ggC 3'	5'-ggC 3'
	419	540			355		368	368	409			357
	5'-gTC 3'	5'-gAC 3'			5'-TCA 3'		5'-gTT 3'	5'-gTC 3'	5'-ggC 3'			5'-Tgg 3'
							557		409			
							5'-ggA 3'		5'-ggC 3'			
									420			
									5'-TTA 3'			
3'-primer(s) ³	301	272	234	204	412	301	259	311	544	317	272	302
	5'-gTC 3'	5'-TgC 3'	5'-TCT 3'	5'-TCT 3'	5'-gTT 3'	5'-gCT 3'	5'-gTT 3'	5'-ggg 3'	5'-ggT 3'	5'-ggA 3'	5'-TgC 3'	5'-ggC 3'
	301	272	272		463	302	259	311		320		302
	5'-gTC 3'	5'-TgT 3'	5'-TgA 3'		5'-gCT 3'	5'-ggC 3'	5'-gTT 3'	5'-ggg 3'		5'-AgT 3'		5'-ggT 3'
	570	309	272			311	266	538				477
	5'-CCg 3'	5'-ATC 3'	5'-Tgg 3'			5'-ggT 3'	5'-TCC 3'	5'-gTC 3'				5'-gCg 3'
	572	605				312	272					
	5'-gCg 3'	5'-gCT 3'				5'-AgT 3'	5'-Tgg 3'					
							292					
							5'-gTA 3'					
							572					
							5'-gCg 3'					
Well No.	37	38	39	40	41	42	43	44	45	46	47	48

Well No.	49	50	51	52	53	54	55	56	57	58	59	60
Length of spec.	180	130	105	325	80	215	130	90	90	55	120	390
PCR product		290			115		280	410	175	100	210	430
		465			160					145		
					195					395		
					225					430		
					260							
Length of int. pos. control ¹	1070	1070	1070	800	1070	800	800	1070	1070	1070	800	1070
5'-primer(s) ²	355	72	540	1 st I	355	141	209	41	209	48	357	49
	5'-TCA 3'	5'-gCT 3'	5'-gAC 3'	5'-CAg 3'	5'-TCA 3'	5'-ATT 3'	5'-ggC 3'	5'-CTg 3'	5'-ggC 3'	5'-gCC 3'	5'-Tgg 3'	5'-CAg 3'
	363	103			418	357	209	363	363	206		
	5'-Agg 3'	5'-CCT 3'			5'-Agg 3'	5'-Tgg 3'	5'-ggg 3'	5'-AAT 3'	5'-AgC 3'	5'-gAA 3'		
		209			499		355	368		259		
		5'-ggC 3'			5'-TCT 3'		5'-TCA 3'	5'-gTT 3'		5'-ggg 3'		
3'-primer(s) ³	499	193	603	282	538	309	299	282	259	272	435	272
	5'-ggA 3'	5'-CgT 3'	5'-gTg 3'	5'-gCC 3'	5'-gTC 3'	5'-ATC 3'	5'-TCA 3'	5'-gCC 3'	5'-gTT 3'	5'-TgT 3'	5'-TCT 3'	5'-TgC 3'
		193	603		572	538	583	418	499	309	527	309
		5'-TgT 3'	5'-gTg 3'		5'-gCA 3'	5'-gTC 3'	5'-gTg 3'	5'-gTC 3'	5'-ggA 3'	5'-ATC 3'	5'-CCT 3'	5'-ATC 3'
		387	605		572		605					
		5'-TCg 3'	5'-gCT 3'		5'-gCg 3'		5'-gCT 3'					
					573							
					5'-AgT 3'							
Well No.	49	50	51	52	53	54	55	56	57	58	59	60



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

Well No.	61	62	63	64	65	66	67	68	69	70	71	72
Length of spec.	145	300	160	180	180	90	90	95	115	360	350	285
PCR product	210		370	210	205	150			150			
				240		410						
Length of int. pos. control ¹	1070	1070	1070	1070	1070	800	1070	1070	1070	1070	1070	1070
5'-primer(s) ²	142	1 st I	15	141	97	15	207	206	165	1 st I	1 st I	355
	5'-TCA 3'	5'-CAG 3'	5'-gCA 3'	5'-ATT 3'	5'-TCg 3'	5'-gCA 3'	5'-ACC 3'	5'-gAC 3'	5'-ACC 3'	5'-CAg 3'	5'-CAg 3'	5'-TCA 3'
	206		418	362	420	209	208		463			363
	5'-gAA 3'		5'-Agg 3'	5'-gAC 3'	5'-TTA 3'	5'-ggC 3'	5'-CgT 3'		5'-TgA 3'			5'-AgC 3'
				362		209	209					
				5'-ggT 3'		5'-ggA 3'	5'-ggA 3'					
				420		449	369					
				5'-TTA 3'		5'-CCC 3'	5'-TAC 3'					
							369					
							5'-TAT 3'					
3'-primer(s) ³	311	259	207	311	259	256	256	259	272	317	311	603
	5'-ggT 3'	5'-CTC 3'	5'-TCC 3'	5'-ggT 3'	5'-gTT 3'	5'-CCC 3'	5'-CCC 3'	5'-gTT 3'	5'-TgC 3'	5'-ggA 3'	5'-ggT 3'	5'-gTg 3'
		259	226	559	559	259	256		538	320		605
		5'-CTC 3'	5'-CAC 3'	5'-CgT 3'	5'-CAG 3'	5'-CTT 3'	5'-CTC 3'		5'-CCA 3'	5'-AgT 3'		5'-gCT 3'
		262	538			263	420					
		5'-TgC 3'	5'-gTC 3'			5'-gTT 3'	5'-gCT 3'					
						559						
						5'-CgT 3'						
Well No.	61	62	63	64	65	66	67	68	69	70	71	72

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 25 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.



0197

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 6R6

Lot-specific information

DR LOW PRIMER SPECIFICATION

Well No.	73	74	75	76	77	78	79	80	81	82	83	84
Length of spec.	210	200	200	210	75	75	85	100	170	85	85	70
PCR product	235 260		230		205	210	210	175	210 235	120 215 250	100 135 160 190	175
Length of int.	515	430	430	430	515	430	430	430	430	515	430	430
pos. control ¹												
5'-primer(s) ²	11(119) 5'-gCT ^{3'}	14(129) 5'-gAA ^{3'}	12(121) 5'-CTg ^{3'}	11(119) 5'-gCT ^{3'}	13(125) 5'-gTC ^{3'}	12(122) 5'-TAC ^{3'}	13(125) 5'-gTC ^{3'}	13(124) 5'-AAC ^{3'}	11(118) 5'-AgA ^{3'}	16(133) 5'-gTT ^{3'}	26(165) 5'-TAT ^{3'}	26(164) 5'-gTA ^{3'}
	14(129) 5'-gAA ^{3'}		13(126) 5'-Agg ^{3'}	13(126) 5'-Agg ^{3'}	47(227) 5'-gTT ^{3'}	13(125) 5'-gTC ^{3'}		13(125) 5'-ATA ^{3'}	14(127) 5'-ATA ^{3'}	16(133) 5'-gTT ^{3'}	37(196) 5'-AgA ^{3'}	31(178) 5'-gCg ^{3'}
			13(126) 5'-AAg ^{3'}	13(126) 5'-AAg ^{3'}		16(133) 5'-gTT ^{3'}		13(125) 5'-gTC ^{3'}	14(127) 5'-ATA ^{3'}	16(133) 5'-gCT ^{3'}	58(261) 5'-gAg ^{3'}	
			13(126) 5'-AgA ^{3'}						14(127) 5'-gTA ^{3'}	60(265) 5'-AgT ^{3'}		
									16(133) 5'-gTT ^{3'}			
									29(174) 5'-AgA ^{3'}			
									31(178) 5'-gCg ^{3'}			
3'-primer(s) ³	67(286) 5'-gAg ^{3'}	67(286) 5'-gAT ^{3'}	67(286) 5'-gAT ^{3'}	67(286) 5'-gAA ^{3'}	73(305) 5'-ggC ^{3'}	26(164) 5'-ggT ^{3'}	28(171) 5'-CTC ^{3'}	33(184) 5'-gTg ^{3'}	71(298) 5'-CTC ^{3'}	74(307) 5'-CAg ^{3'}	57(257) 5'-CgA ^{3'}	38(200) 5'-gCg ^{3'}
	67(286) 5'-gAg ^{3'}	67(286) 5'-gAA ^{3'}	67(286) 5'-gAT ^{3'}	67(286) 5'-gAg ^{3'}	73(305) 5'-ggC ^{3'}	67(286) 5'-gAT ^{3'}	70(295) 5'-CTg ^{3'}	58(260) 5'-Cgg ^{3'}	73(305) 5'-ggC ^{3'}	86(344) 5'-CAC ^{3'}	73(305) 5'-ggC ^{3'}	72(303) 5'-gCA ^{3'}
	67(286) 5'-gAT ^{3'}	67(286) 5'-gAA ^{3'}	67(286) 5'-gAT ^{3'}	67(286) 5'-gAg ^{3'}	73(305) 5'-ggC ^{3'}	71(299) 5'-gCT ^{3'}			77(317) 5'-AAT ^{3'}		78(319) 5'-CAC ^{3'}	74(307) 5'-CgC ^{3'}
	70(297) 5'-CTg ^{3'}		70(295) 5'-CTg ^{3'}	70(296) 5'-TgT ^{3'}	73(305) 5'-ggC ^{3'}				78(319) 5'-CAC ^{3'}			
	71(299) 5'-gCg ^{3'}		70(295) 5'-Tg ^{3'}	72(302) 5'-Cgg ^{3'}	73(306) 5'-Tgg ^{3'}				78(319) 5'-CAA ^{3'}			
	73(305) 5'-ggC ^{3'}		70(296) 5'-TCC ^{3'}		74(308) 5'-CCC ^{3'}				78(319) 5'-CgC ^{3'}			
	77(317) 5'-AgT ^{3'}		71(299) 5'-gCg ^{3'}		75(310) 5'-CAA ^{3'}							
	86(344) 5'-CCA ^{3'}		71(299) 5'-gCT ^{3'}		77(317) 5'-AgT ^{3'}							
			72(301) 5'-Cg ^{3'}									
			73(305) 5'-ggC ^{3'}									
			77(317) 5'-AgT ^{3'}									
Well No.	73	74	75	76	77	78	79	80	81	82	83	84



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

Well No.	85	86	87	88	89	90	91	92	93	94	95
Length of spec.	100	90	220	195	175	75	110	110	80	130	180
PCR product	175	110		225	240	100	145	150	110	210	215
						140	170	180	160		
						185		220	190		
						240			240		
Length of int. pos. control ¹	430	430	430	430	430	430	430	430	515	430	430
5'-primer(s) ²	13(125) 5'-gTC 3'	13(124) 5'-Cgg 3'	10(116) 5'-gCT 3'	10(116) 5'-gCT 3'	13(125) 5'-gTC 3'	5(101) 5'-CAA 3'	26(164) 5'-gTA 3'	13(125) 5'-gTC 3'	10(116) 5'-gCT 3'	28(170) 5'-gAT 3'	9(112) 5'-TgC 3'
	16(133) 5'-gTC 3'	16(133) 5'-gTT 3'	12(122) 5'-TAT 3'	12(122) 5'-TAT 3'	13(125) 5'-gTg 3'	37(197) 5'-gTT 3'	34(189) 5'-CAg 3'	34(189) 5'-CAg 3'	10(116) 5'-gCT 3'	28(170) 5'-gAg 3'	13(125) 5'-gTA 3'
	38(200) 5'-CgT 3'		13(125) 5'-gTC 3'	12(123) 5'-ACA 3'	114(429) 5'-CTg 3'	37(197) 5'-gTA 3'	34(189) 5'-CAg 3'	37(196) 5'-AgC 3'	11(119) 5'-gCT 3'	105(401) 5'-AAA 3'	
				13(125) 5'-gTC 3'			34(189) 5'-CCg 3'		38(199) 5'-TCC 3'		
				13(125) 5'-gTg 3'							
				16(133) 5'-gTC 3'							
				16(133) 5'-gTT 3'							
3'-primer(s) ³	58(260) 5'-CCT 3'	28(171) 5'-CTC 3'	70(295) 5'-gTC 3'	67(286) 5'-gAA 3'	58(260) 5'-Cgg 3'	48(230) 5'-CCT 3'	57(257) 5'-CAg 3'	57(257) 5'-CAg 3'	51(239) 5'-CCC 3'	81(328) 5'-gTA 3'	57(256) 5'-gTC 3'
	58(260) 5'-CCT 3'	30(175) 5'-gTg 3'	71(299) 5'-gCT 3'	71(298) 5'-CTC 3'	58(260) 5'-CAg 3'	57(257) 5'-CAg 3'	70(295) 5'-CTg 3'	60(265) 5'-gTg 3'	57(256) 5'-gCT 3'	87(346) 5'-CTC 3'	57(256) 5'-gCT 3'
	58(260) 5'-CCT 3'	38(199) 5'-CAg 3'	71(299) 5'-ACT 3'	71(298) 5'-CgC 3'	181(630) 5'-CTT 3'	71(298) 5'-CgC 3'	70(296) 5'-TCC 3'	70(296) 5'-TCC 3'	60(266) 5'-Agg 3'	87(346) 5'-CTT 3'	57(258) 5'-gCg 3'
	58(261) 5'-TCT 3'		74(307) 5'-CgT 3'					74(307) 5'-CAg 3'	77(317) 5'-AAT 3'	135(490) 5'-gCT 3'	57(258) 5'-gCT 3'
			78(319) 5'-CAC 3'								58(260) 5'-CCT 3'
											70(296) 5'-TgT 3'
Well No.	85	86	87	88	89	90	91	92	93	94	95

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 73 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.



0197

101.701-24/06 – including *Taq* polymerase
 101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

CELL LINE VALIDATION SHEET																				
HLA-A low resolution primer set ²																				
				Well																
				1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
				Prod. No.:	202135301	202243602	202135303	202135304	202135305	202244206	202135307	202135308	202135309	202244210	202135311	202135312	202135313	202135314	202135315	202244216
	IHWC cell line ¹	A*	A*																	
1	9001 SA	*24:02		-	-	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-
2	9280 LK707	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011 E4181324	*01:01		+	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275 GU373	*30:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
5	9009 KAS011	*01:01		+	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM	*02:01	*26:03	-	+	-	-	-	-	-	+	-	-	+	-	+	-	-	-	-
7	9020 QBL	*26:01		-	-	-	-	-	-	-	+	-	+	-	-	+	-	-	-	-
8	9025 DEU	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
9	9026 YAR	*26:01		-	-	-	-	-	-	-	+	-	+	-	-	+	-	-	-	-
10	9107 LKT3	*24:02		-	-	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT	*29:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-
12	9052 DBB	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
15	9075 DKB	*24:02		-	-	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007	*29:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-
17	9282 CTM3953540	*03:01	*80:01	-	-	+	-	-	+	-	-	-	-	-	-	+	-	-	-	-
18	9257 32367	*33:03	*74:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA	*02:17		-	+	-	-	-	-	W	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL	*02:01	*34:01	-	+	-	-	-	-	-	+	-	-	+	+	-	-	-	-	-
24	9035 JBUSH	*32:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
25	9049 IBW9	*33:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49	*02:05		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007	*24:10	*29:01	-	-	-	-	+	-	+	-	-	-	-	-	-	+	-	-	-
28	9320 BEL5GB	*02:01	*29:02	-	+	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
29	9050 MOU	*29:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
30	9021 RSH	*30:01	*68:02	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	+	-
31	9019 DUCAF	*30:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
32	9297 HAG	*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
34	9104 DHIF	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
35	9302 SSTO	*32:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
36	9024 KT17	*02:06	*11:01	-	+	-	+	-	-	-	-	-	-	+	-	-	-	-	-	-
37	9065 HHKB	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	+	-	-	-	-
38	9099 LZL	*02:17		-	+	-	-	-	-	W	-	-	-	-	-	-	-	-	-	-
39	9315 CML	*01:01	*03:01	+	-	+	+	-	-	-	-	-	-	-	-	+	-	-	-	-
40	9134 WHONP199	*02:07	*30:01	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
41	9055 H0301	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	+	-	-	-	-
42	9066 TAB089	*02:07		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526	*02:06	*02:07	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM	*66:01		-	-	-	-	-	-	-	+	-	-	+	-	+	-	-	-	-
45	9239 SHJO	*23:01	*24:02	-	-	-	-	+	+	+	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	+	-	-	-	-
47	9045 TUBO	*02:16	*03:01	-	+	+	-	-	-	-	-	-	-	-	-	+	-	-	-	-
48	9303 TER-ND	*02:01	*11:01	-	+	-	+	-	-	-	-	-	-	+	-	-	-	-	-	-



0197

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

CELL LINE VALIDATION SHEET												
HLA-A low resolution primer set ²												
				Well								
				17	18	19	20	21	22	23	24	
				Prod. No.:	202135317	202135318	202135319	202135320	202135321	202136922	202135323	202135324
	IHWC cell line ¹	A*	A*									
1	9001 SA	*24:02		-	-	-	-	-	-	-	-	-
2	9280 LK707	*02:01		-	-	-	-	-	-	-	-	-
3	9011 E4181324	*01:01		-	-	-	-	-	-	-	-	-
4	9275 GU373	*30:01		-	-	-	-	-	-	-	-	-
5	9009 KAS011	*01:01		-	-	-	-	-	-	-	-	-
6	9353 SM	*02:01	*26:03	-	-	-	-	-	-	-	+	-
7	9020 QBL	*26:01		-	-	-	-	-	-	-	-	-
8	9025 DEU	*31:01		-	-	-	-	-	-	-	-	-
9	9026 YAR	*26:01		-	-	-	-	-	-	-	-	-
10	9107 LKT3	*24:02		-	-	-	-	-	-	-	-	-
11	9051 PITOUT	*29:02		-	-	-	-	-	-	-	-	-
12	9052 DBB	*02:01		-	-	-	-	-	-	-	-	-
13	9004 JESTHOM	*02:01		-	-	-	-	-	-	-	-	-
14	9071 OLGA	*31:01		-	-	-	-	-	-	-	-	-
15	9075 DKB	*24:02		-	-	-	-	-	-	-	-	-
16	9037 SWEIG007	*29:02		-	-	-	-	-	-	-	-	-
17	9282 CTM3953540	*03:01	*80:01	-	-	-	-	-	-	-	+	-
18	9257 32367	*33:03	*74:01	-	+	+	-	-	-	-	-	-
19	9038 BM16	*02:01		-	-	-	-	-	-	-	-	-
20	9059 SLE005	*02:01		-	-	-	-	-	-	-	-	-
21	9064 AMALA	*02:17		-	-	-	-	-	-	-	-	-
22	9056 KOSE	*02:01		-	-	-	-	-	-	-	-	-
23	9124 IHL	*02:01	*34:01	-	-	-	-	-	-	-	-	-
24	9035 JBUSH	*32:01		+	-	-	-	-	-	-	-	-
25	9049 IBW9	*33:01		-	+	-	-	-	-	-	-	-
26	9285 WT49	*02:05		-	-	-	-	-	-	-	-	-
27	9191 CH1007	*24:10	*29:01	-	-	-	-	-	-	-	-	-
28	9320 BEL5GB	*02:01	*29:02	-	-	-	-	-	-	-	-	-
29	9050 MOU	*29:02		-	-	-	-	-	-	-	-	-
30	9021 RSH	*30:01	*68:02	-	-	-	+	-	-	-	-	-
31	9019 DUCAF	*30:02		-	-	-	-	-	-	-	-	-
32	9297 HAG	*02:01		-	-	-	-	-	-	-	-	-
33	9098 MT14B	*31:01		-	-	-	-	-	-	-	-	-
34	9104 DHIF	*31:01		-	-	-	-	-	-	-	-	-
35	9302 SSTO	*32:01		+	-	-	-	-	-	-	-	-
36	9024 KT17	*02:06	*11:01	-	-	-	-	-	-	-	-	-
37	9065 HHKB	*03:01		-	-	-	-	-	-	-	-	-
38	9099 LZL	*02:17		-	-	-	-	-	-	-	-	-
39	9315 CML	*01:01	*03:01	-	-	-	-	-	-	-	-	-
40	9134 WHONP199	*02:07	*30:01	-	-	-	-	-	-	-	-	-
41	9055 H0301	*03:01		-	-	-	-	-	-	-	-	-
42	9066 TAB089	*02:07		-	-	-	-	-	-	-	-	-
43	9076 T7526	*02:06	*02:07	-	-	-	-	-	-	-	-	-
44	9057 TEM	*66:01		-	-	-	-	-	-	-	-	-
45	9239 SHJO	*23:01	*24:02	-	-	-	-	-	-	-	-	-
46	9013 SCHU	*03:01		-	-	-	-	-	-	-	-	-
47	9045 TUBO	*02:16	*03:01	-	-	-	-	-	-	-	-	-
48	9303 TER-ND	*02:01	*11:01	-	-	-	-	-	-	-	-	-



0197

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

One or more 5'- and 3'-primers in primer solutions 4, 5, 7, 8, 12, 14 to 17 and 20 to 23 were tested by separately adding one or more 3'-primers, and one or more 5'-primers respectively.

One or more 5'-primers in primer solutions 1, 2, 9, 13 and 24 were tested by separately adding one or two 3'-primers. One or more 3'-primers in primer solutions 3, 6, 10 and 18 were tested by separately adding one or two 5'-primers. In primer solutions 2, 3, 9, 11, 12, 15, 16, 20 and 22 one or two 5'-primers were not possible to be tested, and in primer solutions 3, 6, 8 to 10, 12, 13, 16 and 18 one or two 3'-primers were not possible to be tested.



0197

For *In Vitro* Diagnostic Use

MA123 v02 SSP PI Template

Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

CELL LINE VALIDATION SHEET																				
HLA-B low resolution SSP kit ²																				
				Well																
				25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	
				Prod. No.:	202241501	202240902	202240903	202240904	202240905	202240906	202240907	202240908	202240909	202240910	202243711	202241512	202240913	202348014	202240915	202240916
IHCW cell line ¹		B*																		
1	9001 SA	*07:02		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707	*52:01	*73:01	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
3	9011 E4181324	*52:01		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
4	9275 GU373	*15:10	*53:01	-	-	-	-	-	-	-	-	+	-	-	-	+	-	+	-	-
5	9009 KAS011	*37:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
6	9353 SM	*39:01	*51:01	-	-	-	-	+	-	-	-	-	+	-	-	-	-	-	-	-
7	9020 QBL	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-
8	9025 DEU	*35:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
9	9026 YAR	*38:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3	*54:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT	*44:03		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB	*57:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
13	9004 JESTHOM	*27:05		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	+
14	9071 OLGA	*15:01	*15:20	-	-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	-
15	9075 DKB	*40:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007	*40:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282 CTM3953540	*08:01	*55:01	-	+	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
18	9257 32367	*14:01	*56:01	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-
20	9059 SLE005	*40:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA	*15:01		-	-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	-
22	9056 KOSE	*35:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL	*40:02	*56:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035 JBUSH	*38:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9	*14:02		-	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49	*58:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
27	9191 CH1007	*07:05	*51:01	+	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
28	9320 BEL5GB	*44:02	*44:03	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU	*44:03		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH	*42:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
31	9019 DUCAF	*18:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-
32	9297 HAG	*41:02		-	-	-	+	-	-	-	+	-	-	-	-	-	-	-	-	-
33	9098 MT14B	*40:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF	*38:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO	*44:02		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17	*15:01	*35:01	-	-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	-
37	9065 HHKB	*07:02		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL	*15:01		-	-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	-
39	9315 CML	*08:01	*27:05	-	+	-	-	+	-	-	+	-	-	-	-	-	-	-	-	+
40	9134 WHONP199	*13:02	*46:01	-	-	+	+	-	-	+	-	-	+	-	+	-	-	-	-	-
41	9055 H0301	*14:02		-	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089	*46:01		-	-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	-
43	9076 T7526	*46:01		-	-	-	-	-	-	-	+	-	-	-	-	+	-	-	-	-
44	9057 TEM	*38:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO	*42:01	*50:01	-	-	-	+	-	-	-	+	-	-	-	-	-	-	-	-	-
46	9013 SCHU	*07:02		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO	*51:01		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
48	9303 TER-ND	*35:01	*44:03	-	-	-	+	-	-	-	-	-	-	-	-	+	-	-	-	-



0197

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

CELL LINE VALIDATION SHEET																			
HLA-B low resolution SSP kit ²																			
			Prod. No.:	Well															
				41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56
				202240917	202240918	202240919	202240920	202240921	202240922	202240923	202240924	202240925	202241526	202240927	202240928	202240929	202241530	202240931	202243732
	IHWC cell line ¹	B*																	
1	9001 SA	*07:02		-	-	-	-	-	-	-	+	-	-	-	+	-	-	-	-
2	9280 LK707	*52:01	*73:01	-	-	-	-	-	-	-	+	-	-	-	+	-	-	-	-
3	9011 E4181324	*52:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275 GU373	*15:10	*53:01	+	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011	*37:01		+	-	-	+	-	-	-	-	-	-	-	-	+	-	-	-
6	9353 SM	*39:01	*51:01	-	-	-	-	+	-	+	+	-	-	-	-	-	+	-	-
7	9020 QBL	*18:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU	*35:01		+	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026 YAR	*38:01		-	-	-	-	+	+	+	-	-	-	-	-	-	-	-	-
10	9107 LKT3	*54:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	+
11	9051 PITOUT	*44:03		+	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
12	9052 DBB	*57:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM	*27:05		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	+
14	9071 OLGA	*15:01	*15:20	+	-	+	-	-	-	-	-	-	+	-	-	-	-	-	-
15	9075 DKB	*40:01		-	-	-	-	-	-	-	+	+	-	-	-	-	-	-	-
16	9037 SWEIG007	*40:02		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
17	9282 CTM3953540	*08:01	*55:01	-	-	-	-	-	-	-	+	-	-	+	+	-	-	-	+
18	9257 32367	*14:01	*56:01	-	-	-	-	-	-	+	+	-	-	+	-	-	-	-	+
19	9038 BM16	*18:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005	*40:01		-	-	-	-	-	-	-	+	+	-	-	-	-	-	-	-
21	9064 AMALA	*15:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE	*35:03		+	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL	*40:02	*56:02	-	-	-	-	-	-	-	-	-	+	-	+	-	-	-	+
24	9035 JBUSH	*38:01		-	-	-	-	+	+	+	-	-	-	-	-	-	-	-	-
25	9049 IBW9	*14:02		-	-	-	-	-	-	+	+	-	-	-	-	-	-	-	-
26	9285 WT49	*58:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007	*07:05	*51:01	-	-	-	-	-	-	-	+	-	-	-	+	-	+	-	-
28	9320 BEL5GB	*44:02	*44:03	+	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
29	9050 MOU	*44:03		+	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
30	9021 RSH	*42:01		-	-	-	-	-	-	-	+	-	-	+	+	-	-	-	-
31	9019 DUCAF	*18:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG	*41:02		-	-	-	-	-	-	-	+	-	-	+	-	-	-	-	-
33	9098 MT14B	*40:01		-	-	-	-	-	-	-	+	+	-	-	-	-	-	-	-
34	9104 DHIF	*38:01		-	-	-	-	+	+	+	-	-	-	-	-	-	-	-	-
35	9302 SSTO	*44:02		+	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
36	9024 KT17	*15:01	*35:01	+	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB	*07:02		-	-	-	-	-	-	-	+	-	-	-	+	-	-	-	-
38	9099 LZL	*15:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML	*08:01	*27:05	-	-	-	-	-	-	-	+	-	-	+	+	-	-	-	+
40	9134 WHONP199	*13:02	*46:01	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	+
41	9055 H0301	*14:02		-	-	-	-	-	-	+	+	-	-	-	-	-	-	-	-
42	9066 TAB089	*46:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-
43	9076 T7526	*46:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-
44	9057 TEM	*38:01		-	-	-	-	+	+	+	-	-	-	-	-	-	-	-	-
45	9239 SHJO	*42:01	*50:01	-	-	-	-	-	-	-	+	-	-	+	+	-	-	-	-
46	9013 SCHU	*07:02		-	-	-	-	-	-	-	+	-	-	-	+	-	-	-	-
47	9045 TUBO	*51:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-
48	9303 TER-ND	*35:01	*44:03	+	+	-	-	-	-	-	-	-	-	-	-	+	-	-	-



0197

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

CELL LINE VALIDATION SHEET																			
HLA-B low resolution SSP kit ²																			
			Prod. No.:	Well															
				57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72
				202240933	202240934	202240935	202240936	202240937	202240938	202240939	202240940	202348041	202240942	202348043	202240944	202240945	202240946	202240947	202240948
	IHWC cell line ¹	B*																	
1	9001 SA	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	+	
2	9280 LK707	*52:01	*73:01	-	+	+	+	-	+	-	-	-	-	-	-	+	+	-	
3	9011 E4181324	*52:01		-	+	+	+	-	+	-	-	-	-	-	-	+	-	-	
4	9275 GU373	*15:10	*53:01	-	-	-	+	-	-	-	-	-	-	-	+	-	+	-	
5	9009 KAS011	*37:01		-	-	-	-	-	+	-	-	-	-	-	-	+	-	-	
6	9353 SM	*39:01	*51:01	-	-	+	+	-	-	-	-	-	-	-	+	-	+	-	
7	9020 QBL	*18:01		-	-	-	-	-	-	-	-	-	-	-	+	-	+	-	
8	9025 DEU	*35:01		-	-	-	-	-	-	-	-	-	-	-	+	-	+	-	
9	9026 YAR	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	
10	9107 LKT3	*54:01		-	-	-	-	-	-	+	+	-	-	-	-	-	+	-	
11	9051 PITOUT	*44:03		-	-	-	-	-	+	-	-	-	-	-	-	+	-	-	
12	9052 DBB	*57:01		-	-	-	+	-	-	-	-	-	+	-	-	-	+	-	
13	9004 JESTHOM	*27:05		-	-	-	-	-	+	-	-	-	-	-	-	+	-	-	
14	9071 OLGA	*15:01	*15:20	-	-	-	-	-	+	-	-	-	-	-	+	-	+	-	
15	9075 DKB	*40:01		-	-	-	-	+	+	-	-	-	-	-	-	-	+	-	
16	9037 SWEIG007	*40:02		-	-	-	-	+	+	-	-	-	-	-	-	-	+	-	
17	9282 CTM3953540	*08:01	*55:01	-	-	+	-	-	-	-	+	-	-	-	-	-	+	+	
18	9257 32367	*14:01	*56:01	-	-	-	-	-	-	-	-	+	-	-	-	-	+	-	
19	9038 BM16	*18:01		-	-	-	-	-	-	-	-	-	+	-	-	-	+	-	
20	9059 SLE005	*40:01		-	-	-	-	+	+	-	-	-	-	-	-	-	+	-	
21	9064 AMALA	*15:01		-	-	-	-	-	+	-	-	-	-	-	+	-	+	-	
22	9056 KOSE	*35:03		-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
23	9124 IHL	*40:02	*56:02	-	-	-	-	+	+	-	-	+	-	+	-	-	+	-	
24	9035 JBUSH	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	
25	9049 IBW9	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
26	9285 WT49	*58:01		-	-	-	+	-	-	-	-	-	+	-	-	+	-	-	
27	9191 CH1007	*07:05	*51:01	-	-	+	+	-	-	-	-	-	-	-	+	-	+	+	
28	9320 BEL5GB	*44:02	*44:03	-	-	-	-	-	+	+	-	-	-	-	-	+	-	-	
29	9050 MOU	*44:03		-	-	-	-	-	+	-	-	-	-	-	-	+	-	-	
30	9021 RSH	*42:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	+	
31	9019 DUCAF	*18:01		-	-	-	-	-	-	-	-	-	-	-	+	-	+	-	
32	9297 HAG	*41:02		-	-	-	+	+	-	-	-	-	-	-	-	-	+	+	
33	9098 MT14B	*40:01		-	-	-	-	+	+	-	-	-	-	-	-	-	+	-	
34	9104 DHIF	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	
35	9302 SSTO	*44:02		-	-	-	-	-	+	+	-	-	-	-	-	+	-	-	
36	9024 KT17	*15:01	*35:01	-	-	-	-	-	+	-	-	-	-	+	+	-	+	-	
37	9065 HHKB	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	+	
38	9099 LZL	*15:01		-	-	-	-	+	-	-	-	-	-	-	+	-	+	-	
39	9315 CML	*08:01	*27:05	-	-	-	-	-	+	-	-	-	-	-	-	+	+	+	
40	9134 WHONP199	*13:02	*46:01	-	-	-	-	-	+	-	-	-	-	-	+	+	+	-	
41	9055 H0301	*14:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
42	9066 TAB089	*46:01		-	-	-	-	-	-	-	-	-	-	-	-	+	+	-	
43	9076 T7526	*46:01		-	-	-	-	-	-	-	-	-	-	-	-	+	+	-	
44	9057 TEM	*38:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	
45	9239 SHJO	*42:01	*50:01	-	-	+	-	+	+	-	-	+	-	-	-	+	+	+	
46	9013 SCHU	*07:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	+	
47	9045 TUBO	*51:01		-	-	+	+	-	-	-	-	-	-	-	+	-	-	-	
48	9303 TER-ND	*35:01	*44:03	-	-	-	-	-	+	-	-	-	-	-	+	-	+	-	



0197

For *In Vitro* Diagnostic Use
MA123 v02 SSP PI Template
Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

One or more 5'-primers and 3'-primers in primer solutions 27, 30, 35, 38, 43, 53, 55 and 64 to 66 were tested by separately adding one or two additional 3'-primers, and one or two additional 5'-primers respectively.

One or more 5'-primers in primer solutions 25, 45, 49, 50, 67 and 72 were tested by separately adding one or two additional 3'-primers. One or more 3'-primers in primer solutions 26, 28, 32, 37, 39, 42, 58 and 60 were tested by separately adding one, two or three additional 5'-primers.

In primer mixes 35, 41, 45, 55, 66 and 67 one or two 5'-primers were not possible to be tested, and in primer mixes 27, 46, 51, 53, 62, 63, 66 and 70 one or two 3'-primers were not possible to be tested.



0197

For *In Vitro* Diagnostic Use

MA123 v02 SSP PI Template

Date: February 2023, Rev. No: 00

101.701-24/06 – including *Taq* polymerase
 101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

CELL LINE VALIDATION SHEET																				
DR low resolution primer set ²																				
				Well																
				73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	
				Prod. No.:	202243801	202243802	202243803	202243804	202243805	202243806	202243807	202243808	202135509	202243810	202243811	202243812	202243813	202243814	202243815	202243816
IHCW cell line ¹		DRB1																		
1	9001 SA	*01:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
2	9280 LK707	*15:02	*04:05	-	-	+	-	-	-	-	+	-	-	-	-	-	-	-	-	
3	9011 E4181324	*15:02		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
4	9275 GU373	*03:01		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-	
5	9009 KAS011	*16:01		-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
6	9353 SM	*04:07	*08:03	-	-	-	-	-	-	-	+	-	+	-	-	-	-	-	-	
7	9020 QBL	*03:01		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-	
8	9025 DEU	*04:01		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	
9	9026 YAR	*04:02		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	
10	9107 LKT3	*04:05		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	
11	9051 PITOUT	*07:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	
12	9052 DBB	*07:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	
13	9004 JESTHOM	*01:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
14	9071 OLGA	*08:02		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	+	
15	9075 DKB	*09:01		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	
16	9037 SWEIG007	*11:01		-	-	-	-	-	-	-	-	-	-	-	+	-	+	+	+	
17	9282 CTM3953540	*03:01	*13:01	-	-	-	-	+	+	-	-	-	-	-	-	-	-	+	+	
18	9257 32367	*09:01	*11:01	-	-	-	-	-	-	-	-	-	-	+	-	+	-	+	+	
19	9038 BM16	*12:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	
20	9059 SLE005	*13:02		-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	+	
21	9064 AMALA	*14:02		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
22	9056 KOSE	*13:02	*14:54	-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	+	
23	9124 IHL	*08:03	*14:14	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
24	9035 JBUSH	*11:01		-	-	-	-	-	-	-	-	-	-	-	+	-	+	+	+	
25	9049 IBW9	*07:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	
26	9285 WT49	*03:01		-	-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	
27	9191 CH1007	*04:05	*10:01	-	-	-	-	-	-	-	+	-	-	-	+	-	-	-	-	
28	9320 BEL5GB	*04:16	*07:01	-	-	-	-	-	-	-	+	+	-	-	-	-	-	-	-	
29	9050 MOU	*07:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	
30	9021 RSH	*03:02		-	-	-	-	+	-	+	-	-	-	-	-	-	-	-	-	
31	9019 DUCAF	*03:01		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-	-	
32	9297 HAG	*13:03		-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	-	
33	9098 MT14B	*04:04		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	
34	9104 DHIF	*11:01		-	-	-	-	-	-	-	-	-	-	-	+	-	+	+	+	
35	9302 SSTO	*04:03		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	
36	9024 KT17	*04:03	*04:06	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	
37	9065 HHKB	*13:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	+	
38	9099 LZL	*14:02		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
39	9315 CML	*03:01	*04:01	-	-	-	-	+	+	-	+	-	-	-	-	-	-	-	-	
40	9134 WHONP199	*07:01	*09:01	-	-	-	-	-	-	-	-	+	-	+	-	-	-	-	-	
41	9055 H0301	*13:02		-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	+	
42	9066 TAB089	*08:03		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
43	9076 T7526	*09:01		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	
44	9057 TEM	*14:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
45	9239 SHJO	*07:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	
46	9013 SCHU	*15:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
47	9045 TUBO	*11:04	*12:01	-	-	-	-	-	-	-	-	-	-	-	-	+	+	+	+	
48	9303 TER-ND	*01:03		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	



0197

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

CELL LINE VALIDATION SHEET											
DR low resolution primer set ²											
				Well ³							
				89	90	91	92	93	94	95	
				Prod. No.:	202128517	202243818	202243819	202243820	202348129	202243830	202243831
IHCW cell line ¹		DRB1									
1	9001 SA	*01:01		-	-	-	-	-	-	-	-
2	9280 LK707	*15:02	*04:05	-	-	-	-	-	-	+	+
3	9011 E4181324	*15:02		-	-	-	-	-	-	-	+
4	9275 GU373	*03:01		+	-	-	-	+	-	-	-
5	9009 KAS011	*16:01		-	-	-	-	-	-	-	+
6	9353 SM	*04:07	*08:03	-	-	-	-	-	-	+	-
7	9020 QBL	*03:01		+	-	-	-	+	-	-	-
8	9025 DEU	*04:01		-	-	-	-	-	-	+	-
9	9026 YAR	*04:02		-	-	-	-	-	-	+	-
10	9107 LKT3	*04:05		-	-	-	-	-	-	+	-
11	9051 PITOUT	*07:01		-	-	-	-	-	-	+	-
12	9052 DBB	*07:01		-	+	-	-	-	-	+	-
13	9004 JESTHOM	*01:01		-	-	-	-	-	-	-	-
14	9071 OLGA	*08:02		-	-	-	-	-	-	-	-
15	9075 DKB	*09:01		-	-	+	-	-	+	-	-
16	9037 SWEIG007	*11:01		-	-	-	-	+	-	-	-
17	9282 CTM3953540	*03:01	*13:01	+	-	-	-	+	-	-	-
18	9257 32367	*09:01	*11:01	-	-	+	-	+	+	-	-
19	9038 BM16	*12:01		-	-	-	-	+	-	-	-
20	9059 SLE005	*13:02		+	-	-	-	+	-	-	-
21	9064 AMALA	*14:02		+	-	+	-	+	-	-	-
22	9056 KOSE	*13:02	*14:54	+	+	+	+	+	-	-	-
23	9124 IHL	*08:03	*14:14	+	-	+	+	+	-	-	-
24	9035 JBUSH	*11:01		-	-	-	-	+	-	-	-
25	9049 IBW9	*07:01		-	-	-	-	-	+	-	-
26	9285 WT49	*03:01		+	-	-	-	+	-	-	-
27	9191 CH1007	*04:05	*10:01	-	-	-	-	-	-	+	-
28	9320 BEL5GB	*04:16	*07:01	-	-	-	-	-	-	+	-
29	9050 MOU	*07:01		-	-	-	-	-	-	+	-
30	9021 RSH	*03:02		+	-	-	-	+	-	-	-
31	9019 DUCAF	*03:01		+	-	-	-	+	-	-	-
32	9297 HAG	*13:03		+	-	-	-	+	-	-	-
33	9098 MT14B	*04:04		-	-	-	-	-	+	-	-
34	9104 DHIF	*11:01		-	-	-	-	+	-	-	-
35	9302 SSTO	*04:03		-	-	-	-	-	-	+	-
36	9024 KT17	*04:03	*04:06	-	-	-	-	-	-	+	-
37	9065 HHKB	*13:01		+	-	-	-	+	-	-	-
38	9099 LZL	*14:02		+	-	+	-	+	-	-	-
39	9315 CML	*03:01	*04:01	+	-	-	-	+	+	-	-
40	9134 WHONP199	*07:01	*09:01	-	-	+	-	-	+	-	-
41	9055 H0301	*13:02		+	-	-	-	+	-	-	-
42	9066 TAB089	*08:03		-	-	-	-	-	-	-	-
43	9076 T7526	*09:01		-	-	+	-	-	+	-	-
44	9057 TEM	*14:01		+	+	+	+	+	-	-	-
45	9239 SHJO	*07:01		-	-	-	-	-	-	+	-
46	9013 SCHU	*15:01		-	-	-	-	-	-	-	+
47	9045 TUBO	*11:04	*12:01	-	-	-	-	+	-	-	-
48	9303 TER-ND	*01:03		-	-	-	-	-	-	-	-



0197

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

³The DRB4*01:03:01:02N allele is amplified by primer mix 90 in the 9052 (DBB) cell line. One or more 5'- and 3'-primers in primer solutions 81, 83, 87 and 90 were tested by separately adding one 3'-primer, and one or three 5'-primers respectively.

One or more 5'-primers in primer solutions 76, 78, 84 and 88 were tested by separately adding one 3'-primer. One or more 3'-primers in primer solutions 73 to 75, 77, 92, 93 and 95 were tested by separately adding one or more 5'-primers.

One or more 5'-primers in primer solutions 73, 75, 76, 80 to 82, 85 to 89, 91, 92 and 94 were not possible to be tested. One or more 3'-primers in primer solutions 73 to 78, 81, 85, 87, 94 and 95 were not possible to be tested.



0197

101.701-24/06 – including *Taq* polymerase
101.701-24u/06u – without *Taq* polymerase

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **6R6**

Lot-specific information

ADDRESSES:

Manufacturer:

CareDx AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-508 939 00

Fax: +46-8-717 88 18

E-mail: orders-se@caredx.com

Web page: www.caredx.com

CareDx Lab Solutions Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-653-78171

Fax: 610-344-7989

E-mail: orders-us@caredx.com

Web page: www.caredx.com

For information on CareDx distributors worldwide, contact **CareDx AB**.



0197

For *In Vitro* Diagnostic Use

MA123 v02 SSP PI Template

Date: February 2023, Rev. No: 00