

101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information
Olerup SSP® DPB1

Product number:	101.311-24/04 – including <i>Taq</i> pol. 101.311-24u/04u – without <i>Taq</i> pol.
Lot number:	8L1
Expiry date:	2025-08-01
Number of tests:	24 tests – Product No. 101.311-24/24u 4 tests – Product No. 101.311-04/04u
Number of wells per test:	95+1
Storage - pre-aliquoted primers:	dark, between -15°C and -25°C
- PCR Master Mix:	between -15°C and -25°C
- Adhesive PCR seals	RT

This Product Description is only valid for Lot No. 8L1.

Complete product documentation consists of generic Instructions for Use (IFU), lot specific Product Insert, Worksheet and Certificate.

CHANGES COMPARED TO THE PREVIOUS OLERUP SSP® DPB1 Lot (5L8)

- The product documentation has been updated for new alleles of IMGT 3.44.0.
- The kit resolution focuses on common and well documented (CWD) alleles¹.

¹As described in section Uniquely Identified Alleles.

The DPB1 primer set, specificity and interpretation tables have been updated for the DPB1 alleles described since the previous *Olerup SSP® DPB1* lot was made (Lot No. 5L8).

¹S. J. Mack, P. Cano, J. A. Hollenbach et al.
Common and well-documented HLA alleles: 2012 update to
the CWD catalogue. Tissue Antigens, 2013, 81, 194–203



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

The primers of the wells detailed below have been exchanged, added, or modified compared to the previous lot (**Lot No. 5L8**).

Well	5'-primer	3'-primer	rationale
15	Modified	-	5'-primer modified for improving HLA-specific amplification.
16	Modified	-	5'-primer modified for improving HLA-specific amplification.
25	Modified	-	5'-primer modified for improving HLA-specific amplification.
45	Added	-	5'-primer modified for improved yield of the DPB1*115:01 allele.
51	Modified	-	5'-primer modified for improving HLA-specific amplification.
55	Added	-	5'-primer added for improved yield of the DPB1*115:01 allele.
56	-	Added	3'-primer added from mix 85 for improving HLA-specific amplification.
63	-	Added	3'-primer added from mix 65 for reducing tendency of primer oligomer formation.
65	-	Moved	3'-primer moved to mix 63 for reducing tendency of primer oligomer formation.
70	Modified	-	5'-primer modified for improving HLA-specific amplification.
72	Added	-	5'-primer added for improved yield.
84	-	Added	3'-primer added from mix 85 for reducing tendency of primer oligomer formation.
85	Removed	Moved	5'-primer removed, 3'-primer moved to mixes 46, 56 and 94 for reducing tendency of primer oligomer formation.
91	-	-	Exchanged positive control primer pair for reducing tendency of primer oligomer formation.
94	-	Added	3'-primer added from mix 85 for improving HLA-specific amplification.

Changes in revision R01 compared to R00:

1. In primer mix 48 the specific PCR product of 280 base pairs may give rise to a lower yield of HLA-specific PCR product than the other DPB1 primer mixes. Optimized gel electrophoresis run time may need to be considered to ensure separation from the control band.



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

Well **96** contains Negative Control primer pairs, that will amplify the majority of the *Olerup* SSP® HLA Class I, DRB, DQB1, DPB1 and DQA1 amplicons as well as all the amplicons generated by the control primer pairs matching the human growth hormone gene.

HLA-specific PCR product sizes range from 75 to 200 base pairs.

The PCR product generated by the positive control primer pair is 200 base pairs.

Length of PCR product	105	200	105	80	75	80	85
5'-primer¹	164	340	440	45	45	43	36
	5'-CAC ^{3'}	5'-Agg ^{3'}	5'-TTA ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-TAC ^{3'}
							36
							5'-TAT ^{3'}
3'-primer²	231	2nd I	507	59	58	57	47
	5'-TgC ^{3'}	5'-AAA ^{3'}	5'-TTg ^{3'}	5'-CTC ^{3'}	5'-ggC ^{3'}	5'-CTC ^{3'}	5'-ACA ^{3'}
							48
							5'-gCA ^{3'}
							48
							5'-gCC ^{3'}
							52
							5'-TgT ^{3'}
A*	+	+	+				
B*	+	+	+				
C*	+	+	+				
DRB1				+	+		
DRB3				+	+		
DRB5				+			
DQB1					+		
DPB1						+	
DQA1							+

¹The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

²The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

PRODUCT DESCRIPTION

DPB1 SSP typing

CONTENT

The primer set contains 5'- and 3'-primers for identifying the DPB1*01:01 to DPB1*1195:01 alleles.

PLATE LAYOUT

Each test consists of 96 PCR reactions in a 96 well cut PCR plate.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48
49	50	51	52	53	54	55	56
57	58	59	60	61	62	63	64
65	66	67	68	69	70	71	72
73	74	75	76	77	78	79	80
81	82	83	84	85	86	87	88
89	90	91	92	93	94	95	NC

The 96 well cut PCR plate is marked with ‘DPB1’ in silver/gray ink.

Well No. 1 is marked with the Lot No. ‘8L1’.

Wells 1 to 95 – DPB1 high resolution primers.

Well 96 – Negative Control (NC).

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays and should be disregarded.

The PCR plates are heat-sealed with a PCR-compatible foil.

INTERPRETATION

Only DPB1 alleles will be amplified by the DPB1 typing kit. Thus, the interpretation of DPB1 typings is not influenced by the DPB2 gene or other HLA class II genes.

For further details see Specificity Table.

UNIQUELY IDENTIFIED ALLELES

All the DPB1 alleles, i.e. **DPB1*01:01 to DPB1*1195:01 alleles**, recognized by the HLA Nomenclature Committee in April 2021^{1,2} will be amplified by the primers in the DPB1 subtyping kit.



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

The DPB1 kit enables separation of the confirmed DPB1 alleles as listed in the IMGT/HLA database 3.33.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources. Current allele confirmation status for DPB1 alleles is listed below.

The DPB1 kit also enables identification of many null and alternatively expressed alleles.

The following DPB1 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix	Alleles	Primer mix
DPB1*32:01, 339:01	42	DPB1*250:01, 268:01	71, 81
DPB1*120:01N, 453:01, 521:01	44	DPB1*281:01, 591:01	81
DPB1*149:01, 840:01	49	DPB1*498:01, 572:01	85
DPB1*153:01, 321:01, 595:01	53		

¹DPB1 alleles listed on the IMGT/HLA web page 2021-April-20, release 3.44.0, www.ebi.ac.uk/imgt/hla.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

ALLELE CONFIRMATION STATUS

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DPB1*01:01:01:01	Confirmed	DPB1*02:01:03	Confirmed	DPB1*03:01:08	Unconfirmed	DPB1*04:01:22	Confirmed
DPB1*01:01:01:02	Confirmed	DPB1*02:01:04	Confirmed	DPB1*04:01:01:01	Confirmed	DPB1*04:01:23	Unconfirmed
DPB1*01:01:01:03	Confirmed	DPB1*02:01:05	Unconfirmed	DPB1*04:01:01:02	Confirmed	DPB1*04:01:24	Confirmed
DPB1*01:01:01:04	Unconfirmed	DPB1*02:01:06	Unconfirmed	DPB1*04:01:01:03	Confirmed	DPB1*04:01:25	Unconfirmed
DPB1*01:01:02:01	Confirmed	DPB1*02:01:07	Unconfirmed	DPB1*04:01:01:04	Confirmed	DPB1*04:01:26	Unconfirmed
DPB1*01:01:02:02	Unconfirmed	DPB1*02:01:08	Confirmed	DPB1*04:01:01:05	Confirmed	DPB1*04:01:27	Unconfirmed
DPB1*01:01:03	Unconfirmed	DPB1*02:01:09	Confirmed	DPB1*04:01:01:06	Unconfirmed	DPB1*04:01:28	Unconfirmed
DPB1*01:01:04	Unconfirmed	DPB1*02:01:10	Unconfirmed	DPB1*04:01:01:07	Unconfirmed	DPB1*04:01:29	Confirmed
DPB1*01:01:05	Unconfirmed	DPB1*02:01:11	Confirmed	DPB1*04:01:01:08	Unconfirmed	DPB1*04:01:30	Unconfirmed
DPB1*01:01:06	Unconfirmed	DPB1*02:01:12	Confirmed	DPB1*04:01:01:09	Unconfirmed	DPB1*04:01:31	Unconfirmed
DPB1*02:01:02:01	Confirmed	DPB1*02:01:13	Confirmed	DPB1*04:01:01:10	Unconfirmed	DPB1*04:01:32	Confirmed
DPB1*02:01:02:02	Confirmed	DPB1*02:01:14	Unconfirmed	DPB1*04:01:01:11	Unconfirmed	DPB1*04:01:33	Unconfirmed
DPB1*02:01:02:03	Confirmed	DPB1*02:01:15	Confirmed	DPB1*04:01:01:12	Unconfirmed	DPB1*04:01:34	Unconfirmed
DPB1*02:01:02:04	Confirmed	DPB1*02:01:16	Unconfirmed	DPB1*04:01:01:13	Unconfirmed	DPB1*04:01:35	Unconfirmed
DPB1*02:01:02:05	Confirmed	DPB1*02:01:17	Unconfirmed	DPB1*04:01:01:14	Unconfirmed	DPB1*04:01:36	Unconfirmed
DPB1*02:01:02:06	Confirmed	DPB1*02:01:18	Confirmed	DPB1*04:01:01:15	Unconfirmed	DPB1*04:01:37	Unconfirmed
DPB1*02:01:02:07	Confirmed	DPB1*02:01:19	Unconfirmed	DPB1*04:01:01:16	Unconfirmed	DPB1*04:02:01:01	Confirmed
DPB1*02:01:02:08	Confirmed	DPB1*02:01:20	Confirmed	DPB1*04:01:01:17	Unconfirmed	DPB1*04:02:01:02	Confirmed
DPB1*02:01:02:09	Confirmed	DPB1*02:01:21	Unconfirmed	DPB1*04:01:01:18	Confirmed	DPB1*04:02:01:03	Confirmed
DPB1*02:01:02:10	Confirmed	DPB1*02:01:22	Confirmed	DPB1*04:01:01:19	Unconfirmed	DPB1*04:02:01:04	Unconfirmed
DPB1*02:01:02:11	Confirmed	DPB1*02:01:23	Confirmed	DPB1*04:01:01:20	Confirmed	DPB1*04:02:01:05	Unconfirmed
DPB1*02:01:02:12	Confirmed	DPB1*02:01:24	Confirmed	DPB1*04:01:01:21	Confirmed	DPB1*04:02:01:06	Unconfirmed
DPB1*02:01:02:13	Unconfirmed	DPB1*02:01:25	Confirmed	DPB1*04:01:01:22	Confirmed	DPB1*04:02:01:07	Unconfirmed
DPB1*02:01:02:14	Unconfirmed	DPB1*02:01:26	Unconfirmed	DPB1*04:01:01:23	Unconfirmed	DPB1*04:02:01:08	Unconfirmed
DPB1*02:01:02:15	Unconfirmed	DPB1*02:01:27	Unconfirmed	DPB1*04:01:01:24N	Unconfirmed	DPB1*04:02:01:09	Unconfirmed
DPB1*02:01:02:16	Unconfirmed	DPB1*02:01:28	Unconfirmed	DPB1*04:01:01:25	Unconfirmed	DPB1*04:02:01:10	Unconfirmed
DPB1*02:01:02:17	Unconfirmed	DPB1*02:01:29	Unconfirmed	DPB1*04:01:01:26	Confirmed	DPB1*04:02:01:11	Unconfirmed
DPB1*02:01:02:18	Unconfirmed	DPB1*02:02:01:01	Confirmed	DPB1*04:01:01:27	Unconfirmed	DPB1*04:02:02	Unconfirmed
DPB1*02:01:02:19	Unconfirmed	DPB1*02:02:01:02	Confirmed	DPB1*04:01:01:28	Confirmed	DPB1*04:02:03	Confirmed
DPB1*02:01:02:20	Unconfirmed	DPB1*02:02:01:03	Unconfirmed	DPB1*04:01:02	Unconfirmed	DPB1*04:02:04	Unconfirmed
DPB1*02:01:02:21	Unconfirmed	DPB1*02:02:01:04	Unconfirmed	DPB1*04:01:03	Confirmed	DPB1*04:02:05	Confirmed
DPB1*02:01:02:22	Unconfirmed	DPB1*02:02:01:05	Unconfirmed	DPB1*04:01:04:01	Confirmed	DPB1*04:02:06	Unconfirmed
DPB1*02:01:02:23	Confirmed	DPB1*02:02:01:06	Unconfirmed	DPB1*04:01:04:02	Unconfirmed	DPB1*04:02:07	Confirmed
DPB1*02:01:02:24	Confirmed	DPB1*02:02:01:07	Unconfirmed	DPB1*04:01:05	Unconfirmed	DPB1*04:02:08	Unconfirmed
DPB1*02:01:02:25	Confirmed	DPB1*03:01:01:01	Confirmed	DPB1*04:01:06	Confirmed	DPB1*04:02:09	Unconfirmed
DPB1*02:01:02:26	Unconfirmed	DPB1*03:01:01:02	Confirmed	DPB1*04:01:07	Confirmed	DPB1*04:02:10	Confirmed
DPB1*02:01:02:27	Unconfirmed	DPB1*03:01:01:03	Confirmed	DPB1*04:01:08	Confirmed	DPB1*04:02:11	Unconfirmed
DPB1*02:01:02:28	Unconfirmed	DPB1*03:01:01:04	Unconfirmed	DPB1*04:01:09	Confirmed	DPB1*04:02:12	Confirmed
DPB1*02:01:02:29	Unconfirmed	DPB1*03:01:01:05	Unconfirmed	DPB1*04:01:10	Confirmed	DPB1*04:02:13	Confirmed
DPB1*02:01:02:30	Unconfirmed	DPB1*03:01:01:06	Unconfirmed	DPB1*04:01:11	Confirmed	DPB1*05:01:01:01	Confirmed
DPB1*02:01:02:31	Unconfirmed	DPB1*03:01:01:07	Unconfirmed	DPB1*04:01:12	Confirmed	DPB1*05:01:01:02	Unconfirmed
DPB1*02:01:02:32	Unconfirmed	DPB1*03:01:01:08	Unconfirmed	DPB1*04:01:13	Confirmed	DPB1*05:01:01:03	Unconfirmed
DPB1*02:01:02:33	Unconfirmed	DPB1*03:01:01:09	Unconfirmed	DPB1*04:01:14	Confirmed	DPB1*05:01:01:04	Confirmed
DPB1*02:01:02:34	Unconfirmed	DPB1*03:01:01:10	Unconfirmed	DPB1*04:01:15	Confirmed	DPB1*05:01:01:05	Unconfirmed
DPB1*02:01:02:35	Unconfirmed	DPB1*03:01:02	Unconfirmed	DPB1*04:01:16	Confirmed	DPB1*05:01:01:06	Unconfirmed
DPB1*02:01:02:36	Unconfirmed	DPB1*03:01:03	Unconfirmed	DPB1*04:01:17	Unconfirmed	DPB1*05:01:01:07	Unconfirmed
DPB1*02:01:02:37	Unconfirmed	DPB1*03:01:04	Confirmed	DPB1*04:01:18	Unconfirmed	DPB1*05:01:01:08	Confirmed
DPB1*02:01:02:38	Unconfirmed	DPB1*03:01:05	Unconfirmed	DPB1*04:01:19	Unconfirmed	DPB1*05:01:01:09	Unconfirmed
DPB1*02:01:02:39	Unconfirmed	DPB1*03:01:06	Unconfirmed	DPB1*04:01:20	Unconfirmed	DPB1*05:01:01:10	Unconfirmed
DPB1*02:01:02:40	Unconfirmed	DPB1*03:01:07	Unconfirmed	DPB1*04:01:21	Unconfirmed	DPB1*05:01:02	Unconfirmed



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

Allele	Status¹	Allele	Status¹	Allele	Status¹	Allele	Status¹
DPB1*05:01:03	Confirmed	DPB1*20:01:01:02	Unconfirmed	DPB1*49:01:01:01	Confirmed	DPB1*88:01	Confirmed
DPB1*05:01:04	Unconfirmed	DPB1*20:01:02	Unconfirmed	DPB1*49:01:01:02	Unconfirmed	DPB1*89:01	Confirmed
DPB1*05:01:05	Unconfirmed	DPB1*20:01:03	Unconfirmed	DPB1*49:01:01:03	Confirmed	DPB1*90:01:01	Confirmed
DPB1*05:01:06	Unconfirmed	DPB1*20:01:04	Unconfirmed	DPB1*50:01	Confirmed	DPB1*90:01:02	Confirmed
DPB1*05:01:07	Unconfirmed	DPB1*21:01	Confirmed	DPB1*51:01	Confirmed	DPB1*91:01:01:01	Confirmed
DPB1*06:01:01:01	Confirmed	DPB1*22:01	Confirmed	DPB1*52:01	Unconfirmed	DPB1*91:01:01:02	Unconfirmed
DPB1*06:01:01:02	Confirmed	DPB1*23:01:01	Confirmed	DPB1*53:01	Unconfirmed	DPB1*92:01	Unconfirmed
DPB1*06:01:01:03	Unconfirmed	DPB1*23:01:02	Unconfirmed	DPB1*54:01	Unconfirmed	DPB1*93:01	Confirmed
DPB1*06:01:02	Unconfirmed	DPB1*24:01	Confirmed	DPB1*55:01:01:01	Confirmed	DPB1*94:01	Unconfirmed
DPB1*06:01:03	Unconfirmed	DPB1*25:01	Confirmed	DPB1*55:01:01:02	Confirmed	DPB1*95:01	Unconfirmed
DPB1*06:01:04	Unconfirmed	DPB1*26:01:01	Unconfirmed	DPB1*56:01	Unconfirmed	DPB1*96:01	Unconfirmed
DPB1*06:01:05	Unconfirmed	DPB1*26:01:02	Confirmed	DPB1*57:01	Confirmed	DPB1*97:01	Unconfirmed
DPB1*08:01	Unconfirmed	DPB1*26:01:03	Unconfirmed	DPB1*58:01	Confirmed	DPB1*98:01	Confirmed
DPB1*09:01:01	Confirmed	DPB1*27:01	Confirmed	DPB1*59:01	Confirmed	DPB1*99:01	Unconfirmed
DPB1*09:01:02	Confirmed	DPB1*28:01	Confirmed	DPB1*60:01	Confirmed	DPB1*100:01	Unconfirmed
DPB1*09:01:03	Unconfirmed	DPB1*29:01	Confirmed	DPB1*61:01N	Confirmed	DPB1*101:01	Unconfirmed
DPB1*09:01:04	Confirmed	DPB1*30:01:01:01	Confirmed	DPB1*62:01	Confirmed	DPB1*102:01	Unconfirmed
DPB1*10:01:01:01	Confirmed	DPB1*30:01:01:02	Unconfirmed	DPB1*63:01	Confirmed	DPB1*103:01	Unconfirmed
DPB1*10:01:01:02	Unconfirmed	DPB1*31:01	Confirmed	DPB1*64:01N	Unconfirmed	DPB1*104:01:01:01	Confirmed
DPB1*10:01:02	Unconfirmed	DPB1*32:01	Confirmed	DPB1*65:01:01	Unconfirmed	DPB1*104:01:01:02	Confirmed
DPB1*10:01:03	Unconfirmed	DPB1*33:01:01:01	Confirmed	DPB1*65:01:02	Unconfirmed	DPB1*104:01:01:03	Confirmed
DPB1*11:01:01:01	Confirmed	DPB1*33:01:01:02	Confirmed	DPB1*66:01	Unconfirmed	DPB1*104:01:01:04	Unconfirmed
DPB1*11:01:01:02	Unconfirmed	DPB1*33:01:01:03	Unconfirmed	DPB1*67:01	Unconfirmed	DPB1*104:01:01:05	Unconfirmed
DPB1*11:01:02	Unconfirmed	DPB1*33:01:01:04	Unconfirmed	DPB1*68:01	Confirmed	DPB1*105:01:01:01	Confirmed
DPB1*11:01:03	Confirmed	DPB1*33:01:01:05	Unconfirmed	DPB1*69:01:01:01	Confirmed	DPB1*105:01:01:02	Confirmed
DPB1*11:01:04	Unconfirmed	DPB1*34:01:01:01	Confirmed	DPB1*69:01:01:02	Unconfirmed	DPB1*105:01:01:03	Confirmed
DPB1*13:01:01:01	Confirmed	DPB1*34:01:01:02	Unconfirmed	DPB1*70:01	Unconfirmed	DPB1*105:01:01:04	Unconfirmed
DPB1*13:01:01:02	Unconfirmed	DPB1*34:01:02	Unconfirmed	DPB1*71:01:01	Confirmed	DPB1*105:01:01:05	Unconfirmed
DPB1*13:01:01:03	Unconfirmed	DPB1*35:01:01	Confirmed	DPB1*71:01:02	Unconfirmed	DPB1*105:01:01:06	Unconfirmed
DPB1*13:01:01:04	Unconfirmed	DPB1*36:01	Confirmed	DPB1*72:01:01:01	Confirmed	DPB1*105:01:01:07	Unconfirmed
DPB1*13:01:01:05	Unconfirmed	DPB1*37:01	Unconfirmed	DPB1*72:01:01:02	Confirmed	DPB1*105:01:01:08	Unconfirmed
DPB1*13:01:01:06	Unconfirmed	DPB1*38:01	Confirmed	DPB1*72:01:01:03	Unconfirmed	DPB1*105:01:01:09	Unconfirmed
DPB1*13:01:02	Confirmed	DPB1*39:01:01:01	Confirmed	DPB1*73:01	Unconfirmed	DPB1*106:01	Confirmed
DPB1*14:01:01:01	Confirmed	DPB1*39:01:01:02	Unconfirmed	DPB1*74:01	Unconfirmed	DPB1*107:01	Unconfirmed
DPB1*14:01:01:02	Unconfirmed	DPB1*39:01:01:03	Unconfirmed	DPB1*75:01	Unconfirmed	DPB1*108:01	Unconfirmed
DPB1*14:01:01:03	Unconfirmed	DPB1*39:01:01:04	Unconfirmed	DPB1*76:01	Confirmed	DPB1*109:01	Confirmed
DPB1*14:01:02	Unconfirmed	DPB1*39:01:02	Unconfirmed	DPB1*77:01	Confirmed	DPB1*110:01	Unconfirmed
DPB1*14:01:03	Unconfirmed	DPB1*40:01:01:01	Confirmed	DPB1*78:01	Confirmed	DPB1*111:01	Unconfirmed
DPB1*14:01:04	Unconfirmed	DPB1*40:01:01:02	Unconfirmed	DPB1*79:01	Unconfirmed	DPB1*112:01	Unconfirmed
DPB1*14:01:05	Unconfirmed	DPB1*41:01:01:01	Confirmed	DPB1*80:01	Confirmed	DPB1*113:01	Unconfirmed
DPB1*15:01:01	Confirmed	DPB1*41:01:01:02	Unconfirmed	DPB1*81:01:01:01	Confirmed	DPB1*114:01	Unconfirmed
DPB1*15:01:02	Unconfirmed	DPB1*41:01:02	Unconfirmed	DPB1*81:01:01:02	Unconfirmed	DPB1*115:01	Confirmed
DPB1*16:01:01	Confirmed	DPB1*44:01	Confirmed	DPB1*82:01	Unconfirmed	DPB1*116:01	Unconfirmed
DPB1*16:01:02	Unconfirmed	DPB1*45:01	Confirmed	DPB1*83:01	Unconfirmed	DPB1*117:01	Unconfirmed
DPB1*17:01:01:01	Confirmed	DPB1*46:01:01	Confirmed	DPB1*84:01	Unconfirmed	DPB1*118:01	Unconfirmed
DPB1*17:01:01:02	Unconfirmed	DPB1*46:01:02	Unconfirmed	DPB1*85:01:01:01	Confirmed	DPB1*119:01	Unconfirmed
DPB1*18:01	Confirmed	DPB1*47:01:01:01	Confirmed	DPB1*85:01:01:02	Unconfirmed	DPB1*120:01N	Unconfirmed
DPB1*19:01:01:01	Confirmed	DPB1*47:01:01:02	Unconfirmed	DPB1*85:01:02	Unconfirmed	DPB1*121:01	Unconfirmed
DPB1*19:01:01:02	Unconfirmed	DPB1*47:01:01:03	Unconfirmed	DPB1*86:01	Confirmed	DPB1*122:01	Unconfirmed
DPB1*20:01:01:01	Confirmed	DPB1*48:01	Confirmed	DPB1*87:01	Unconfirmed	DPB1*123:01	Unconfirmed



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DPB1*124:01:01:01	Confirmed	DPB1*169:01	Unconfirmed	DPB1*218:01N	Confirmed	DPB1*267:01	Unconfirmed
DPB1*124:01:01:02	Unconfirmed	DPB1*170:01	Confirmed	DPB1*219:01	Unconfirmed	DPB1*268:01	Unconfirmed
DPB1*125:01	Unconfirmed	DPB1*171:01	Unconfirmed	DPB1*220:01	Confirmed	DPB1*269:01	Unconfirmed
DPB1*126:01:01:01	Confirmed	DPB1*172:01	Confirmed	DPB1*221:01	Confirmed	DPB1*270:01	Unconfirmed
DPB1*126:01:01:02	Unconfirmed	DPB1*173:01	Confirmed	DPB1*222:01	Confirmed	DPB1*271:01	Unconfirmed
DPB1*127:01	Confirmed	DPB1*174:01	Unconfirmed	DPB1*223:01	Unconfirmed	DPB1*272:01	Unconfirmed
DPB1*128:01	Confirmed	DPB1*175:01	Unconfirmed	DPB1*224:01	Confirmed	DPB1*273:01	Unconfirmed
DPB1*129:01	Unconfirmed	DPB1*176:01	Unconfirmed	DPB1*225:01	Confirmed	DPB1*274:01	Unconfirmed
DPB1*130:01	Confirmed	DPB1*177:01	Confirmed	DPB1*226:01	Unconfirmed	DPB1*275:01	Unconfirmed
DPB1*131:01:01	Confirmed	DPB1*178:01	Confirmed	DPB1*227:01	Unconfirmed	DPB1*276:01	Unconfirmed
DPB1*131:01:02	Confirmed	DPB1*179:01	Unconfirmed	DPB1*228:01	Unconfirmed	DPB1*277:01	Unconfirmed
DPB1*132:01	Unconfirmed	DPB1*180:01	Confirmed	DPB1*229:01	Unconfirmed	DPB1*278:01	Unconfirmed
DPB1*133:01	Confirmed	DPB1*181:01	Confirmed	DPB1*230:01	Unconfirmed	DPB1*279:01:01	Unconfirmed
DPB1*134:01	Unconfirmed	DPB1*182:01	Confirmed	DPB1*231:01	Unconfirmed	DPB1*279:01:02	Unconfirmed
DPB1*135:01	Confirmed	DPB1*183:01	Unconfirmed	DPB1*232:01	Unconfirmed	DPB1*280:01	Unconfirmed
DPB1*136:01	Confirmed	DPB1*184:01	Confirmed	DPB1*233:01	Confirmed	DPB1*281:01	Unconfirmed
DPB1*137:01	Unconfirmed	DPB1*185:01	Confirmed	DPB1*234:01	Unconfirmed	DPB1*282:01	Confirmed
DPB1*138:01	Confirmed	DPB1*186:01	Confirmed	DPB1*235:01	Unconfirmed	DPB1*283:01	Unconfirmed
DPB1*139:01	Unconfirmed	DPB1*187:01	Confirmed	DPB1*236:01:01	Confirmed	DPB1*284:01	Unconfirmed
DPB1*140:01	Unconfirmed	DPB1*188:01	Unconfirmed	DPB1*236:01:02	Confirmed	DPB1*285:01	Unconfirmed
DPB1*141:01	Confirmed	DPB1*189:01	Confirmed	DPB1*237:01	Unconfirmed	DPB1*286:01	Unconfirmed
DPB1*142:01	Unconfirmed	DPB1*190:01	Confirmed	DPB1*238:01	Unconfirmed	DPB1*287:01	Unconfirmed
DPB1*143:01	Confirmed	DPB1*191:01	Confirmed	DPB1*239:01	Unconfirmed	DPB1*288:01	Unconfirmed
DPB1*144:01	Confirmed	DPB1*192:01	Confirmed	DPB1*240:01	Unconfirmed	DPB1*289:01	Unconfirmed
DPB1*145:01	Confirmed	DPB1*193:01	Confirmed	DPB1*241:01	Unconfirmed	DPB1*290:01	Unconfirmed
DPB1*146:01	Confirmed	DPB1*194:01	Confirmed	DPB1*242:01	Unconfirmed	DPB1*291:01	Unconfirmed
DPB1*147:01	Unconfirmed	DPB1*195:01	Confirmed	DPB1*243:01	Unconfirmed	DPB1*292:01	Unconfirmed
DPB1*148:01	Unconfirmed	DPB1*196:01	Unconfirmed	DPB1*244:01	Confirmed	DPB1*293:01	Unconfirmed
DPB1*149:01	Confirmed	DPB1*197:01	Unconfirmed	DPB1*245:01	Confirmed	DPB1*294:01	Unconfirmed
DPB1*150:01	Confirmed	DPB1*198:01	Confirmed	DPB1*246:01	Unconfirmed	DPB1*295:01	Unconfirmed
DPB1*151:01	Confirmed	DPB1*199:01	Confirmed	DPB1*247:01	Unconfirmed	DPB1*296:01	Confirmed
DPB1*152:01	Confirmed	DPB1*200:01	Confirmed	DPB1*248:01	Unconfirmed	DPB1*297:01	Unconfirmed
DPB1*153:01	Confirmed	DPB1*201:01	Confirmed	DPB1*249:01	Unconfirmed	DPB1*298:01	Unconfirmed
DPB1*154:01N	Confirmed	DPB1*202:01	Unconfirmed	DPB1*250:01	Unconfirmed	DPB1*299:01	Unconfirmed
DPB1*155:01:01	Confirmed	DPB1*203:01:01	Confirmed	DPB1*251:01	Unconfirmed	DPB1*300:01	Unconfirmed
DPB1*155:01:02	Confirmed	DPB1*203:01:02	Unconfirmed	DPB1*252:01	Unconfirmed	DPB1*301:01	Confirmed
DPB1*156:01	Confirmed	DPB1*204:01	Confirmed	DPB1*253:01	Confirmed	DPB1*302:01	Unconfirmed
DPB1*157:01	Confirmed	DPB1*205:01	Unconfirmed	DPB1*254:01	Unconfirmed	DPB1*303:01	Confirmed
DPB1*158:01	Confirmed	DPB1*206:01	Confirmed	DPB1*255:01	Unconfirmed	DPB1*304:01	Confirmed
DPB1*159:01N	Unconfirmed	DPB1*207:01	Confirmed	DPB1*256:01	Unconfirmed	DPB1*305:01	Unconfirmed
DPB1*160:01	Unconfirmed	DPB1*208:01	Confirmed	DPB1*257:01	Confirmed	DPB1*306:01	Confirmed
DPB1*161:01N	Confirmed	DPB1*209:01	Confirmed	DPB1*258:01	Unconfirmed	DPB1*307:01	Unconfirmed
DPB1*162:01:01	Unconfirmed	DPB1*210:01	Confirmed	DPB1*259:01	Unconfirmed	DPB1*308:01	Unconfirmed
DPB1*162:01:02	Confirmed	DPB1*211:01	Confirmed	DPB1*260:01	Confirmed	DPB1*309:01	Confirmed
DPB1*163:01	Unconfirmed	DPB1*212:01	Confirmed	DPB1*261:01	Unconfirmed	DPB1*310:01	Unconfirmed
DPB1*164:01	Unconfirmed	DPB1*213:01	Unconfirmed	DPB1*262:01	Unconfirmed	DPB1*311:01	Confirmed
DPB1*165:01	Confirmed	DPB1*214:01	Confirmed	DPB1*263:01	Unconfirmed	DPB1*312:01	Confirmed
DPB1*166:01	Confirmed	DPB1*215:01	Confirmed	DPB1*264:01	Unconfirmed	DPB1*313:01	Confirmed
DPB1*167:01	Unconfirmed	DPB1*216:01N	Unconfirmed	DPB1*265:01	Confirmed	DPB1*314:01	Unconfirmed
DPB1*168:01	Unconfirmed	DPB1*217:01	Unconfirmed	DPB1*266:01	Unconfirmed	DPB1*315:01	Confirmed



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

Allele	Status¹	Allele	Status¹	Allele	Status¹	Allele	Status¹
DPB1*316:01	Confirmed	DPB1*364:01	Unconfirmed	DPB1*414:01:01:01	Confirmed	DPB1*461:01	Unconfirmed
DPB1*317:01	Unconfirmed	DPB1*365:01	Confirmed	DPB1*414:01:01:02	Unconfirmed	DPB1*462:01	Unconfirmed
DPB1*318:01	Confirmed	DPB1*366:01	Unconfirmed	DPB1*415:01	Confirmed	DPB1*463:01:01:01	Unconfirmed
DPB1*319:01	Confirmed	DPB1*367:01	Unconfirmed	DPB1*416:01:01:01	Confirmed	DPB1*463:01:01:02	Unconfirmed
DPB1*320:01	Unconfirmed	DPB1*368:01	Unconfirmed	DPB1*416:01:01:02	Confirmed	DPB1*463:01:01:03	Unconfirmed
DPB1*321:01	Confirmed	DPB1*369:01	Confirmed	DPB1*417:01	Unconfirmed	DPB1*464:01	Unconfirmed
DPB1*322:01	Unconfirmed	DPB1*370:01	Confirmed	DPB1*418:01	Confirmed	DPB1*465:01	Unconfirmed
DPB1*323:01	Unconfirmed	DPB1*371:01	Unconfirmed	DPB1*419:01	Unconfirmed	DPB1*466:01	Unconfirmed
DPB1*324:01	Unconfirmed	DPB1*372:01	Unconfirmed	DPB1*420:01	Unconfirmed	DPB1*467:01	Confirmed
DPB1*325:01	Unconfirmed	DPB1*373:01	Unconfirmed	DPB1*421:01	Unconfirmed	DPB1*468:01	Unconfirmed
DPB1*326:01	Unconfirmed	DPB1*374:01	Confirmed	DPB1*422:01	Confirmed	DPB1*469:01	Unconfirmed
DPB1*327:01	Confirmed	DPB1*375:01	Confirmed	DPB1*423:01:01	Confirmed	DPB1*470:01	Confirmed
DPB1*328:01N	Confirmed	DPB1*376:01	Unconfirmed	DPB1*423:01:02	Unconfirmed	DPB1*471:01	Unconfirmed
DPB1*329:01	Confirmed	DPB1*377:01	Confirmed	DPB1*424:01	Confirmed	DPB1*472:01	Unconfirmed
DPB1*330:01	Unconfirmed	DPB1*378:01	Unconfirmed	DPB1*425:01	Unconfirmed	DPB1*473:01	Unconfirmed
DPB1*331:01	Unconfirmed	DPB1*379:01	Unconfirmed	DPB1*426:01	Unconfirmed	DPB1*474:01	Unconfirmed
DPB1*332:01	Unconfirmed	DPB1*380:01	Confirmed	DPB1*427:01	Unconfirmed	DPB1*475:01	Confirmed
DPB1*333:01	Confirmed	DPB1*381:01	Unconfirmed	DPB1*428:01	Unconfirmed	DPB1*476:01	Confirmed
DPB1*334:01	Confirmed	DPB1*382:01N	Confirmed	DPB1*429:01	Unconfirmed	DPB1*477:01	Unconfirmed
DPB1*335:01	Unconfirmed	DPB1*383:01	Confirmed	DPB1*430:01	Unconfirmed	DPB1*478:01	Unconfirmed
DPB1*336:01	Unconfirmed	DPB1*384:01	Confirmed	DPB1*431:01	Unconfirmed	DPB1*479:01	Unconfirmed
DPB1*337:01	Unconfirmed	DPB1*385:01	Confirmed	DPB1*432:01	Unconfirmed	DPB1*480:01	Unconfirmed
DPB1*338:01	Confirmed	DPB1*386:01	Unconfirmed	DPB1*433:01	Confirmed	DPB1*481:01	Unconfirmed
DPB1*339:01	Unconfirmed	DPB1*387:01	Unconfirmed	DPB1*434:01	Unconfirmed	DPB1*482:01	Confirmed
DPB1*340:01	Confirmed	DPB1*388:01	Unconfirmed	DPB1*435:01	Unconfirmed	DPB1*483:01	Unconfirmed
DPB1*341:01	Unconfirmed	DPB1*389:01	Unconfirmed	DPB1*436:01	Unconfirmed	DPB1*484:01	Confirmed
DPB1*342:01	Unconfirmed	DPB1*390:01	Confirmed	DPB1*437:01	Confirmed	DPB1*485:01	Confirmed
DPB1*343:01	Unconfirmed	DPB1*391:01	Confirmed	DPB1*438:01	Unconfirmed	DPB1*486:01	Confirmed
DPB1*344:01	Unconfirmed	DPB1*392:01	Confirmed	DPB1*439:01	Unconfirmed	DPB1*487:01	Unconfirmed
DPB1*345:01	Unconfirmed	DPB1*393:01	Confirmed	DPB1*440:01	Confirmed	DPB1*488:01	Confirmed
DPB1*346:01	Unconfirmed	DPB1*394:01	Confirmed	DPB1*441:01	Confirmed	DPB1*489:01	Confirmed
DPB1*347:01	Unconfirmed	DPB1*395:01	Unconfirmed	DPB1*442:01	Unconfirmed	DPB1*490:01	Confirmed
DPB1*348:01:01	Unconfirmed	DPB1*396:01	Confirmed	DPB1*443:01	Unconfirmed	DPB1*491:01	Unconfirmed
DPB1*348:01:02	Unconfirmed	DPB1*397:01	Unconfirmed	DPB1*444:01	Unconfirmed	DPB1*492:01	Unconfirmed
DPB1*349:01	Confirmed	DPB1*398:01	Confirmed	DPB1*445:01	Confirmed	DPB1*493:01	Confirmed
DPB1*350:01	Confirmed	DPB1*399:01	Unconfirmed	DPB1*446:01	Unconfirmed	DPB1*494:01	Unconfirmed
DPB1*351:01	Unconfirmed	DPB1*400:01	Unconfirmed	DPB1*447:01	Unconfirmed	DPB1*495:01	Confirmed
DPB1*352:01	Unconfirmed	DPB1*401:01N	Unconfirmed	DPB1*448:01	Confirmed	DPB1*496:01	Unconfirmed
DPB1*353:01	Unconfirmed	DPB1*402:01	Unconfirmed	DPB1*449:01	Unconfirmed	DPB1*497:01	Unconfirmed
DPB1*354:01:01	Unconfirmed	DPB1*403:01N	Unconfirmed	DPB1*450:01N	Unconfirmed	DPB1*498:01	Confirmed
DPB1*354:01:02	Confirmed	DPB1*404:01	Confirmed	DPB1*451:01	Unconfirmed	DPB1*499:01	Unconfirmed
DPB1*355:01	Confirmed	DPB1*405:01	Unconfirmed	DPB1*452:01	Confirmed	DPB1*500:01	Unconfirmed
DPB1*356:01	Unconfirmed	DPB1*406:01	Unconfirmed	DPB1*453:01	Unconfirmed	DPB1*501:01	Confirmed
DPB1*357:01N	Confirmed	DPB1*407:01	Unconfirmed	DPB1*454:01	Confirmed	DPB1*502:01	Unconfirmed
DPB1*358:01	Unconfirmed	DPB1*408:01	Unconfirmed	DPB1*455:01N	Unconfirmed	DPB1*503:01	Confirmed
DPB1*359:01	Confirmed	DPB1*409:01	Unconfirmed	DPB1*456:01	Unconfirmed	DPB1*504:01	Confirmed
DPB1*360:01	Confirmed	DPB1*410:01	Unconfirmed	DPB1*457:01	Unconfirmed	DPB1*505:01	Confirmed
DPB1*361:01	Unconfirmed	DPB1*411:01	Unconfirmed	DPB1*458:01	Confirmed	DPB1*506:01	Unconfirmed
DPB1*362:01	Unconfirmed	DPB1*412:01	Confirmed	DPB1*459:01	Unconfirmed	DPB1*507:01N	Confirmed
DPB1*363:01	Unconfirmed	DPB1*413:01	Confirmed	DPB1*460:01	Unconfirmed	DPB1*508:01	Confirmed



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DPB1*509:01	Unconfirmed	DPB1*559:01	Unconfirmed	DPB1*609:01	Unconfirmed	DPB1*658:01	Confirmed
DPB1*510:01	Unconfirmed	DPB1*560:01	Unconfirmed	DPB1*610:01	Unconfirmed	DPB1*659:01	Unconfirmed
DPB1*511:01	Confirmed	DPB1*561:01	Unconfirmed	DPB1*611:01	Unconfirmed	DPB1*660:01	Unconfirmed
DPB1*512:01	Confirmed	DPB1*562:01	Unconfirmed	DPB1*612:01	Unconfirmed	DPB1*661:01N	Unconfirmed
DPB1*513:01	Unconfirmed	DPB1*563:01	Unconfirmed	DPB1*613:01	Unconfirmed	DPB1*662:01	Confirmed
DPB1*514:01	Confirmed	DPB1*564:01	Unconfirmed	DPB1*614:01	Unconfirmed	DPB1*663:01	Unconfirmed
DPB1*515:01	Unconfirmed	DPB1*565:01	Unconfirmed	DPB1*615:01	Unconfirmed	DPB1*664:01	Unconfirmed
DPB1*516:01	Unconfirmed	DPB1*566:01	Unconfirmed	DPB1*616:01	Unconfirmed	DPB1*665:01	Unconfirmed
DPB1*517:01	Unconfirmed	DPB1*567:01	Unconfirmed	DPB1*617:01	Unconfirmed	DPB1*666:01	Unconfirmed
DPB1*518:01	Unconfirmed	DPB1*568:01	Unconfirmed	DPB1*618:01	Unconfirmed	DPB1*667:01	Unconfirmed
DPB1*519:01	Confirmed	DPB1*569:01	Unconfirmed	DPB1*619:01	Unconfirmed	DPB1*668:01:01:01	Unconfirmed
DPB1*520:01	Confirmed	DPB1*570:01N	Confirmed	DPB1*620:01	Unconfirmed	DPB1*668:01:01:02	Unconfirmed
DPB1*521:01	Unconfirmed	DPB1*571:01	Unconfirmed	DPB1*621:01	Unconfirmed	DPB1*669:01	Unconfirmed
DPB1*522:01	Confirmed	DPB1*572:01	Confirmed	DPB1*622:01	Unconfirmed	DPB1*670:01	Unconfirmed
DPB1*523:01	Confirmed	DPB1*573:01	Confirmed	DPB1*623:01	Unconfirmed	DPB1*671:01	Unconfirmed
DPB1*524:01	Unconfirmed	DPB1*574:01	Unconfirmed	DPB1*624:01	Unconfirmed	DPB1*672:01	Unconfirmed
DPB1*525:01	Unconfirmed	DPB1*575:01	Unconfirmed	DPB1*625:01	Unconfirmed	DPB1*673:01	Unconfirmed
DPB1*526:01	Confirmed	DPB1*576:01	Confirmed	DPB1*626:01	Unconfirmed	DPB1*674:01	Confirmed
DPB1*527:01	Unconfirmed	DPB1*577:01	Confirmed	DPB1*627:01	Unconfirmed	DPB1*675:01	Unconfirmed
DPB1*528:01	Unconfirmed	DPB1*578:01	Unconfirmed	DPB1*628:01	Unconfirmed	DPB1*676:01	Unconfirmed
DPB1*529:01	Unconfirmed	DPB1*579:01	Unconfirmed	DPB1*629:01	Unconfirmed	DPB1*677:01	Unconfirmed
DPB1*530:01	Confirmed	DPB1*580:01	Unconfirmed	DPB1*630:01	Unconfirmed	DPB1*678:01	Unconfirmed
DPB1*531:01	Confirmed	DPB1*581:01	Unconfirmed	DPB1*631:01	Unconfirmed	DPB1*679:01	Confirmed
DPB1*532:01	Unconfirmed	DPB1*582:01	Unconfirmed	DPB1*632:01	Unconfirmed	DPB1*680:01	Unconfirmed
DPB1*533:01	Unconfirmed	DPB1*583:01	Unconfirmed	DPB1*633:01	Unconfirmed	DPB1*681:01	Confirmed
DPB1*534:01	Unconfirmed	DPB1*584:01	Confirmed	DPB1*634:01	Confirmed	DPB1*682:01	Confirmed
DPB1*535:01	Confirmed	DPB1*585:01	Unconfirmed	DPB1*635:01	Unconfirmed	DPB1*683:01	Confirmed
DPB1*536:01	Unconfirmed	DPB1*586:01	Unconfirmed	DPB1*636:01	Unconfirmed	DPB1*684:01	Unconfirmed
DPB1*537:01	Unconfirmed	DPB1*587:01	Confirmed	DPB1*637:01	Unconfirmed	DPB1*685:01	Unconfirmed
DPB1*538:01	Unconfirmed	DPB1*588:01	Unconfirmed	DPB1*638:01	Unconfirmed	DPB1*686:01	Unconfirmed
DPB1*539:01	Confirmed	DPB1*589:01	Unconfirmed	DPB1*639:01	Unconfirmed	DPB1*687:01	Unconfirmed
DPB1*540:01	Unconfirmed	DPB1*590:01	Unconfirmed	DPB1*640:01	Unconfirmed	DPB1*688:01	Unconfirmed
DPB1*541:01	Unconfirmed	DPB1*591:01	Unconfirmed	DPB1*641:01	Unconfirmed	DPB1*689:01	Unconfirmed
DPB1*542:01	Unconfirmed	DPB1*592:01	Unconfirmed	DPB1*642:01	Unconfirmed	DPB1*690:01	Unconfirmed
DPB1*543:01	Unconfirmed	DPB1*593:01	Unconfirmed	DPB1*643:01	Unconfirmed	DPB1*691:01N	Unconfirmed
DPB1*544:01	Unconfirmed	DPB1*594:01	Unconfirmed	DPB1*644:01	Unconfirmed	DPB1*692:01	Unconfirmed
DPB1*545:01	Unconfirmed	DPB1*595:01	Unconfirmed	DPB1*645:01	Unconfirmed	DPB1*693:01N	Unconfirmed
DPB1*546:01	Unconfirmed	DPB1*596:01	Unconfirmed	DPB1*646:01	Unconfirmed	DPB1*694:01	Unconfirmed
DPB1*547:01	Unconfirmed	DPB1*597:01	Unconfirmed	DPB1*647:01	Unconfirmed	DPB1*695:01	Unconfirmed
DPB1*548:01	Unconfirmed	DPB1*598:01N	Unconfirmed	DPB1*648:01:01:01	Unconfirmed	DPB1*696:01N	Unconfirmed
DPB1*549:01	Unconfirmed	DPB1*599:01	Confirmed	DPB1*648:01:01:02	Unconfirmed	DPB1*697:01Q	Unconfirmed
DPB1*550:01	Unconfirmed	DPB1*600:01	Unconfirmed	DPB1*649:01	Unconfirmed	DPB1*698:01	Confirmed
DPB1*551:01N	Unconfirmed	DPB1*601:01	Unconfirmed	DPB1*650:01	Unconfirmed	DPB1*699:01	Confirmed
DPB1*552:01	Unconfirmed	DPB1*602:01	Unconfirmed	DPB1*651:01	Confirmed	DPB1*700:01N	Confirmed
DPB1*553:01	Unconfirmed	DPB1*603:01	Unconfirmed	DPB1*652:01	Confirmed	DPB1*701:01	Unconfirmed
DPB1*554:01	Unconfirmed	DPB1*604:01	Confirmed	DPB1*653:01	Unconfirmed	DPB1*702:01	Unconfirmed
DPB1*555:01	Unconfirmed	DPB1*605:01	Unconfirmed	DPB1*654:01	Unconfirmed	DPB1*703:01	Unconfirmed
DPB1*556:01	Unconfirmed	DPB1*606:01	Unconfirmed	DPB1*655:01	Unconfirmed	DPB1*704:01	Unconfirmed
DPB1*557:01	Unconfirmed	DPB1*607:01	Unconfirmed	DPB1*656:01	Confirmed	DPB1*705:01	Unconfirmed
DPB1*558:01	Confirmed	DPB1*608:01	Unconfirmed	DPB1*657:01N	Unconfirmed	DPB1*706:01	Unconfirmed



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

Allele	Status ¹
DPB1*707:01	Unconfirmed
DPB1*708:01	Unconfirmed
DPB1*709:01	Unconfirmed
DPB1*710:01	Unconfirmed
DPB1*711:01	Unconfirmed
DPB1*712:01N	Unconfirmed
DPB1*713:01	Unconfirmed
DPB1*714:01	Unconfirmed
DPB1*715:01	Unconfirmed
DPB1*716:01	Unconfirmed
DPB1*717:01	Unconfirmed
DPB1*718:01	Unconfirmed
DPB1*719:01	Unconfirmed
DPB1*720:01	Unconfirmed

¹Allele status “confirmed” or “unconfirmed” as listed on the IMGT/HLA web page 2018-July-11, release 3.33.0, www.ebi.ac.uk/imgt/hla.

RESOLUTION IN HOMO- AND HETEROZYGOTES

Results file with resolution in DPB1 homo- and heterozygotes is available upon request.



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

SPECIFICITY TABLE

DPB1 SSP typing

Specificities and sizes of the PCR products of the 95+1 primer mixes used for DPB1 SSP typing³

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified DPB1 alleles ³
1	245 bp	515 bp	*01:01:01:01-01:01:04, 01:01:06-01:01:07, 03:01:01:01-03:01:14, 13:01:01:01-13:01:04, 25:01, 26:01:01-26:01:03, 29:01, 37:01, 44:01, 50:01, 52:01, 56:01, 61:01N, 78:01, 79:01, 88:01, 90:01:01-90:01:02, 92:01, 104:01:01:01-104:01:04, 107:01, 110:01, 111:01, 114:01, 118:01, 122:01, 124:01:01:01-124:01:02:03, 127:01, 132:01, 133:01, 142:01, 159:01N, 160:01, 162:01:01-162:01:02, 201:01, 203:01:01-203:01:02, 209:01, 220:01, 221:01, 222:01, 246:01, 248:01, 249:01, 267:01, 269:01, 270:01, 291:01, 293:01, 295:01, 315:01, 316:01, 326:01, 327:01, 329:01, 351:01, 362:01, 363:01, 389:01, 390:01, 393:01, 394:01, 403:01N, 404:01, 405:01, 411:01, 412:01, 417:01:01-417:01:02, 436:01, 439:01, 440:01, 448:01, 449:01, 462:01, 472:01, 481:01, 483:01, 506:01, 507:01N, 508:01, 509:01, 516:01, 518:01, 519:01, 536:01, 542:01, 543:01, 544:01, 545:01, 548:01, 562:01, 564:01, 565:01, 567:01, 568:01, 597:01, 599:01, 600:01, 609:01, 610:01, 612:01, 616:01, 631:01, 632:01, 634:01, 636:01, 645:01, 661:01N, 664:01, 667:01, 669:01, 675:01, 676:01, 688:01, 689:01, 691:01N, 703:01, 704:01, 706:01, 708:01, 709:01, 710:01, 711:01, 718:01, 727:01, 728:01, 733:01, 746:01, 760:01, 792:01N, 800:01N, 801:01, 807:01, 810:01, 822:01, 825:01, 829:01, 839:01, 853:01, 855:01, 862:01N, 866:01N, 870:01N, 873:01N, 876:01N, 879:01:01:01-879:01:01:03, 883:01, 888:01, 893:01, 909:01, 917:01N, 920:01, 924:01, 931:01, 935:01Q, 938:01, 941:01N, 946:01, 947:01, 948:01, 950:01N, 952:01, 953:01, 960:01N, 969:01, 970:01, 979:01, 983:01, 984:01N, 991:01, 996:01, 998:01, 999:01, 1000:01, 1006:01, 1008:01, 1012:01, 1014:01, 1021:01, 1024:01, 1027:01, 1030:01, 1038:01Q, 1041:01N, 1045:01N, 1049:01, 1050:01, 1057:01, 1065:01, 1068:01, 1069:01, 1071:01, 1073:01, 1076:01, 1088:01, 1090:01, 1093:01, 1096:01, 1098:01N, 1103:01, 1105:01, 1114:01, 1116:01, 1123:01, 1130:01, 1131:01, 1133:01, 1134:01, 1137:01, 1141:01, 1147:01, 1151:01, 1157:01, 1162:01, 1166:01, 1169:01N, 1183:01, 1185:01
2	220 bp	515 bp	*01:01:01:01-01:01:07, 03:01:01:01-03:01:14, 18:01:01:01-18:01:03, 20:01:01:01-20:01:04, 25:01, 26:01:01-26:01:03, 27:01, 36:01, 50:01, 52:01, 56:01, 61:01N, 78:01, 79:01, 85:01:01:01-85:01:02, 87:01, 89:01, 90:01:01-90:01:02, 92:01, 104:01:01:01-104:01:04, 111:01, 114:01, 118:01, 122:01, 124:01:01:01-124:01:02:03, 127:01, 130:01, 132:01, 142:01, 158:01, 159:01N, 162:01:01-162:01:02, 164:01, 199:01, 201:01, 203:01:01-203:01:02, 209:01, 220:01, 221:01, 222:01, 248:01, 267:01, 270:01, 271:01, 291:01, 292:01, 293:01, 294:01, 315:01, 316:01, 326:01, 327:01, 329:01, 346:01, 351:01, 363:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			369:01, 387:01, 390:01, 391:01, 392:01, 394:01, 403:01N, 404:01, 405:01, 411:01, 412:01, 417:01:01-417:01:02, 436:01, 439:01, 440:01, 462:01, 472:01, 481:01, 483:01, 506:01, 507:01N, 508:01, 509:01, 516:01, 532:01, 536:01, 542:01, 543:01, 545:01, 548:01, 562:01, 568:01, 599:01, 610:01, 613:01, 616:01, 623:01, 632:01, 635:01, 636:01, 645:01, 661:01N, 664:01, 667:01, 669:01, 675:01, 676:01, 688:01, 689:01, 691:01N, 704:01, 706:01, 708:01, 709:01, 710:01, 712:01N, 713:01, 718:01, 727:01, 728:01, 733:01, 749:01, 760:01, 792:01N, 800:01N, 801:01, 807:01, 810:01, 825:01, 829:01, 839:01, 853:01, 855:01, 862:01N, 866:01N, 870:01N, 876:01N, 879:01:01:01-879:01:01:03, 883:01, 897:01, 900:01, 901:01, 905:01, 917:01N, 931:01, 935:01Q, 938:01, 941:01N, 942:01, 943:01, 946:01, 948:01, 950:01N, 952:01, 953:01, 961:01, 969:01, 970:01, 979:01, 980:01, 983:01, 984:01N, 998:01, 999:01, 1000:01, 1012:01, 1014:01, 1021:01, 1024:01, 1027:01, 1030:01, 1034:01, 1038:01Q, 1041:01N, 1049:01, 1050:01, 1068:01, 1071:01, 1073:01, 1076:01, 1080:01, 1088:01, 1090:01, 1093:01, 1096:01, 1098:01N, 1103:01, 1114:01, 1116:01, 1130:01, 1133:01, 1134:01, 1137:01, 1141:01, 1147:01, 1151:01, 1157:01, 1162:01, 1165:01, 1166:01, 1169:01N, 1174:01, 1183:01
3⁶	140 bp	430 bp	*01:01:01:01-01:01:07, 04:01:01:01-04:01:52, 11:01:01:01-11:01:05, 13:01:01:01-13:01:04, 15:01:01:01-15:01:04, 23:01:01:01-23:01:02, 26:01:01-26:01:03, 27:01, 31:01:01:01-31:01:01:02, 33:01:01:01-33:01:01:05, 34:01:01:01-34:01:02, 39:01:01:01-39:01:03, 40:01:01:01-40:01:02, 52:01, 55:01:01:01-55:01:03, 56:01, 58:01, 62:01, 63:01, 65:01:01-65:01:02, 66:01, 67:01, 71:01:01-71:01:02, 72:01:01:01-72:01:01:03, 74:01, 85:01:01:01-85:01:02, 87:01, 89:01, 90:01:01-90:01:02, 95:01, 96:01, 99:01, 102:01, 103:01, 107:01, 110:01, 112:01, 117:01, 118:01, 120:01N, 121:01, 125:01, 126:01:01:01-126:01:01:03, 127:01, 128:01, 133:01, 134:01, 138:01, 142:01, 147:01, 149:01, 150:01, 158:01, 159:01N, 160:01, 162:01:01-162:01:02, 169:01 ^w , 173:01, 174:01, 175:01, 176:01, 177:01, 178:01, 179:01, 180:01, 181:01, 192:01, 193:01, 194:01, 195:01, 199:01, 201:01, 202:01, 206:01, 207:01, 209:01, 213:01:01-213:01:02, 215:01, 216:01N, 220:01, 224:01, 225:01, 227:01:01-227:01:03, 228:01, 230:01, 231:01, 232:01, 240:01, 244:01, 246:01, 247:01, 250:01, 253:01, 255:01, 262:01, 264:01, 267:01, 268:01, 272:01, 275:01, 276:01, 278:01, 279:01:01-279:01:02, 280:01, 281:01, 282:01, 283:01, 290:01, 294:01, 295:01, 298:01, 299:01, 303:01, 304:01, 305:01, 306:01, 314:01, 318:01, 319:01, 320:01, 322:01, 325:01, 326:01, 327:01, 328:01N, 333:01, 334:01:01-334:01:02, 335:01, 336:01, 340:01, 341:01, 345:01, 346:01, 348:01:01-348:01:02, 350:01, 353:01, 354:01:01-354:01:02, 356:01, 357:01N, 360:01, 362:01, 370:01, 371:01, 375:01, 376:01, 377:01, 378:01, 387:01, 388:01, 389:01, 392:01, 393:01, 396:01, 397:01, 398:01, 399:01, 403:01N, 411:01, 412:01, 415:01, 417:01:01-417:01:02, 418:01, 425:01, 426:01, 428:01, 435:01, 436:01, 437:01, 438:01, 440:01, 449:01, 450:01N, 451:01, 453:01, 454:01, 455:01N, 456:01, 458:01:01-458:01:02, 459:01, 462:01, 464:01, 465:01, 468:01, 471:01, 474:01, 475:01, 476:01, 480:01, 481:01, 482:01, 483:01, 485:01, 486:01, 487:01, 490:01, 493:01, 497:01, 500:01, 503:01, 512:01, 516:01, 517:01, 518:01, 519:01, 520:01, 521:01, 522:01, 523:01:01-523:01:02, 524:01, 529:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			531:01, 534:01, 538:01, 542:01, 543:01, 544:01, 553:01, 554:01, 556:01, 559:01, 561:01, 562:01, 563:01, 564:01, 565:01, 569:01, 575:01, 576:01, 578:01, 580:01, 583:01, 584:01:01:01-584:01:02:02, 585:01:01:01-585:01:01:02, 591:01, 592:01, 593:01, 597:01, 599:01, 600:01, 609:01, 612:01, 614:01, 615:01, 616:01, 618:01, 623:01, 625:01, 626:01, 631:01, 632:01, 634:01, 635:01, 643:01, 644:01, 649:01, 654:01, 658:01, 661:01N, 666:01, 667:01, 670:01, 672:01, 677:01, 679:01, 682:01, 683:01, 686:01, 687:01, 694:01, 695:01, 696:01N, 699:01, 702:01, 703:01, 707:01, 708:01, 709:01, 713:01, 716:01, 722:01, 732:01N, 733:01, 739:01, 742:01, 743:01N, 745:01, 747:01, 748:01N, 749:01, 750:01, 753:01, 754:01N, 755:01, 756:01N, 757:01, 758:01, 761:01, 765:01, 767:01, 768:01, 769:01, 772:01, 773:01, 784:01, 787:01, 788:01, 789:01, 795:01, 800:01N, 803:01, 804:01, 806:01:01:01-806:01:01:02, 807:01, 808:01, 810:01, 811:01, 812:01, 813:01, 814:01, 820:01, 821:01N, 822:01, 824:01:01-824:01:02, 826:01, 828:01, 830:01, 835:01:01:01-835:01:01:03, 837:01, 838:01N, 840:01, 842:01, 844:01N, 849:01, 850:01, 852:01, 853:01, 856:01, 859:01, 862:01N, 866:01N, 867:01N, 868:01N, 870:01N, 871:01N, 873:01N, 874:01N, 875:01N, 878:01N, 879:01:01:01-879:01:01:03, 880:01, 882:01, 888:01, 893:01, 895:01, 896:01, 901:01, 904:01, 907:01, 908:01, 910:01, 911:01N, 912:01, 913:01, 915:01, 916:01, 921:01, 922:01, 924:01, 926:01, 930:01, 931:01, 932:01, 934:01Q, 937:01, 939:01, 945:01, 947:01, 953:01, 957:01, 960:01N, 966:01, 969:01, 972:01, 974:01N, 976:01, 978:01, 979:01, 984:01N, 988:01, 989:01, 991:01, 992:01, 993:01, 996:01, 997:01, 998:01, 999:01, 1001:01, 1002:01, 1003:01, 1004:01, 1010:01, 1011:01, 1012:01, 1016:01, 1023:01, 1024:01, 1033:01, 1034:01, 1038:01Q, 1040:01, 1042:01, 1044:01N, 1045:01N, 1046:01, 1048:01, 1050:01, 1054:01, 1057:01, 1060:01, 1062:01, 1063:01, 1064:01, 1065:01, 1066:01, 1068:01, 1069:01, 1073:01, 1074:01, 1076:01, 1078:01, 1079:01N, 1080:01, 1081:01, 1086:01, 1088:01, 1091:01, 1096:01, 1097:01, 1100:01, 1105:01, 1108:01, 1109:01, 1112:01N, 1113:01, 1121:01N, 1122:01, 1123:01, 1129:01, 1131:01, 1132:01, 1133:01, 1135:01N, 1137:01, 1138:01, 1139:01, 1141:01, 1144:01, 1145:01, 1146:01, 1147:01, 1148:01Q, 1151:01, 1152:01, 1155:01, 1161:01, 1162:01, 1164:01, 1166:01, 1167:01, 1170:01, 1173:01, 1177:01, 1181:01, 1183:01, 1184:01, 1185:01, 1192:01, 1193:01N, 1195:01
4	140 bp	430 bp	*01:01:01:01-01:01:07, 04:01:01:01-04:01:01:68, 04:01:03-04:01:15, 04:01:17-04:01:52, 24:01, 26:01:01-26:01:03, 27:01, 28:01, 31:01:01:01-31:01:01:02, 39:01:01:01-39:01:03, 40:01:01:01-40:01:02, 49:01:01:01-49:01:01:03, 51:01:01:01-51:01:01:02, 53:01, 56:01, 65:01:01-65:01:02, 66:01, 72:01:01-72:01:01:03, 76:01, 85:01:01:01-85:01:02, 89:01, 90:01:01-90:01:02, 92:01, 96:01, 99:01, 102:01, 118:01, 120:01N, 121:01, 122:01, 125:01, 126:01:01:01-126:01:01:03, 127:01, 128:01, 134:01, 139:01, 142:01, 149:01, 150:01, 158:01, 159:01N, 162:01:01-162:01:02, 169:01, 174:01, 175:01, 176:01, 177:01, 178:01, 180:01, 181:01, 192:01, 193:01, 194:01, 195:01, 199:01, 203:01:01-203:01:02, 212:01, 213:01:01-213:01:02, 214:01, 215:01, 216:01N, 220:01, 224:01, 225:01, 227:01:03, 228:01, 230:01, 231:01, 232:01, 240:01, 253:01, 254:01, 255:01, 262:01, 267:01, 272:01, 275:01, 276:01, 278:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			279:01:01-279:01:02, 280:01, 281:01, 282:01, 283:01, 291:01, 294:01, 296:01, 298:01, 306:01, 315:01, 319:01, 320:01 ^w , 323:01, 326:01, 327:01, 334:01:01-334:01:02, 335:01, 336:01, 337:01, 341:01, 346:01, 350:01, 354:01:01-354:01:02, 355:01, 356:01, 357:01N, 360:01, 363:01, 372:01, 375:01, 376:01, 377:01, 378:01, 385:01, 387:01, 396:01, 398:01, 399:01, 403:01N, 411:01, 412:01, 415:01, 417:01:01-417:01:02, 418:01, 425:01, 426:01, 427:01, 428:01, 434:01, 436:01, 440:01, 442:01, 453:01, 454:01, 455:01N, 456:01, 459:01, 462:01, 464:01, 465:01, 468:01, 475:01, 476:01 ^w , 479:01, 481:01, 483:01, 485:01, 486:01, 487:01, 490:01, 497:01, 499:01, 500:01, 512:01, 516:01, 520:01, 521:01, 522:01, 523:01:01-523:01:02, 524:01, 529:01, 534:01, 536:01, 538:01, 542:01, 543:01, 553:01, 554:01, 556:01, 559:01, 562:01, 569:01, 574:01, 578:01, 583:01, 584:01:01:01-584:01:02:02, 591:01, 592:01, 593:01, 607:01, 614:01, 615:01, 616:01, 618:01, 623:01, 625:01, 632:01, 636:01, 642:01, 643:01, 658:01, 661:01N, 666:01, 667:01, 670:01, 677:01, 679:01, 682:01, 683:01, 686:01, 687:01, 696:01N, 699:01, 702:01, 708:01, 709:01, 710:01, 713:01, 722:01, 733:01, 736:01, 739:01 ^w , 742:01, 743:01N, 745:01, 747:01, 748:01N, 749:01, 750:01, 752:01, 753:01, 754:01N, 755:01, 756:01N, 757:01, 758:01, 761:01, 765:01, 767:01, 769:01, 772:01, 773:01, 784:01, 788:01, 795:01, 800:01N, 803:01, 804:01, 806:01:01:01-806:01:01:02, 807:01, 808:01, 810:01, 811:01, 812:01, 813:01, 814:01, 820:01, 824:01:01-824:01:02, 825:01, 826:01, 828:01, 837:01, 838:01N, 840:01, 842:01, 844:01N, 849:01, 850:01, 852:01, 853:01, 859:01, 862:01N, 866:01N, 867:01N, 868:01N, 870:01N, 871:01N, 874:01N, 875:01N, 878:01N, 879:01:01:01-879:01:01:03, 880:01, 882:01, 895:01, 900:01, 901:01, 911:01N, 915:01, 916:01, 918:01, 921:01, 922:01, 926:01, 931:01, 932:01, 939:01, 943:01, 945:01, 953:01, 957:01, 966:01, 969:01, 974:01N, 976:01, 978:01, 979:01, 980:01, 984:01N, 988:01, 989:01, 992:01, 993:01, 997:01, 998:01, 999:01, 1001:01, 1002:01, 1003:01, 1004:01, 1010:01, 1011:01, 1012:01, 1016:01, 1023:01, 1024:01, 1033:01, 1034:01, 1038:01Q, 1039:01, 1040:01, 1042:01, 1044:01N, 1046:01, 1048:01, 1050:01, 1060:01, 1062:01, 1064:01, 1066:01, 1068:01, 1074:01, 1076:01, 1079:01N, 1080:01, 1081:01, 1086:01, 1088:01, 1091:01, 1097:01, 1100:01, 1108:01, 1112:01N, 1113:01, 1121:01N, 1129:01, 1132:01, 1135:01N, 1137:01, 1138:01, 1139:01, 1141:01, 1144:01, 1146:01, 1147:01, 1148:01Q, 1150:01, 1151:01, 1152:01, 1155:01, 1161:01, 1162:01, 1164:01, 1166:01, 1167:01, 1169:01N, 1173:01, 1177:01, 1181:01, 1183:01, 1184:01, 1193:01N, 1195:01
5	235 bp	515 bp	*01:01:01:01-01:01:07, 08:01, 19:01:01:01-19:01:01:03, 50:01, 57:01, 65:01:01-65:01:02, 68:01, 73:01, 75:01, 84:01, 90:01:01-90:01:02, 106:01, 118:01, 122:01, 127:01, 150:01, 159:01N, 162:01:01-162:01:02, 166:01, 173:01, 188:01, 201:01, 220:01, 234:01, 236:01:01-236:01:02, 246:01, 248:01, 249:01, 267:01, 278:01, 285:01, 291:01, 304:01, 306:01, 315:01, 316:01, 326:01, 327:01, 340:01, 365:01, 377:01, 379:01, 400:01, 403:01N, 411:01, 417:01:01-417:01:02, 422:01, 436:01, 440:01, 442:01, 462:01, 466:01, 481:01, 483:01, 495:01, 499:01, 533:01, 535:01, 536:01, 538:01, 542:01, 543:01, 558:01, 599:01, 616:01, 620:01, 626:01, 631:01, 632:01, 648:01:01:01-648:01:01:02, 661:01N, 667:01, 684:01, 710:01, 711:01, 720:01, 733:01, 785:01, 800:01N, 807:01, 810:01, 853:01, 862:01N,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			866:01N, 870:01N, 879:01:01:01-879:01:01:03, 904:01, 916:01, 931:01, 953:01, 960:01N, 965:01:01:01-965:01:01:02, 969:01, 976:01, 979:01, 982:01, 984:01N, 998:01, 999:01, 1006:01, 1012:01, 1024:01, 1038:01Q, 1050:01, 1068:01, 1073:01, 1076:01, 1101:01, 1137:01, 1140:01, 1141:01, 1150:01, 1151:01, 1162:01, 1169:01N, 1183:01, 1189:01
6	165 bp	430 bp	*01:01:01:01-01:01:07, 13:01:01:01-13:01:04, 26:01:01-26:01:03, 65:01:01-65:01:02, 107:01, 118:01, 122:01, 127:01, 133:01, 142:01, 159:01N, 160:01, 162:01:01-162:01:02, 173:01, 203:01:01-203:01:02, 220:01, 246:01, 267:01, 269:01, 291:01, 315:01, 326:01, 327:01, 362:01, 363:01, 385:01, 389:01, 393:01, 403:01N, 411:01, 412:01, 417:01:01-417:01:02, 436:01, 440:01, 449:01, 462:01, 481:01, 483:01, 490:01, 503:01, 514:01, 516:01, 518:01, 519:01, 536:01, 542:01, 543:01, 544:01, 562:01, 564:01, 565:01, 597:01, 600:01, 609:01, 612:01, 616:01, 631:01, 632:01, 634:01, 636:01, 661:01N, 667:01, 703:01, 708:01, 710:01, 711:01, 733:01, 800:01N, 807:01, 810:01, 822:01, 825:01, 853:01, 862:01N, 866:01N, 870:01N, 873:01N, 879:01:01:01-879:01:01:03, 888:01, 893:01, 916:01, 924:01, 931:01, 947:01, 953:01, 960:01N, 969:01, 979:01, 984:01N, 991:01, 996:01, 998:01, 999:01, 1008:01, 1024:01, 1038:01Q, 1045:01N, 1050:01, 1057:01, 1065:01, 1068:01, 1069:01, 1076:01, 1088:01, 1105:01, 1123:01, 1131:01, 1137:01, 1141:01, 1147:01, 1150:01, 1151:01, 1162:01, 1166:01, 1169:01N, 1183:01, 1185:01
7	160 bp	430 bp	*01:01:01:01-01:01:07, 13:01:01:01-13:01:04, 26:01:01-26:01:03, 56:01, 65:01:01-65:01:02, 76:01, 90:01:01-90:01:02, 92:01, 107:01, 110:01, 118:01, 122:01, 125:01, 127:01, 133:01, 142:01, 150:01, 159:01N, 160:01, 162:01:01-162:01:02, 173:01, 203:01:01-203:01:02, 220:01, 246:01, 267:01, 269:01, 278:01, 288:01, 291:01, 304:01, 306:01, 315:01, 326:01, 327:01, 362:01, 363:01, 377:01, 385:01, 389:01, 393:01, 403:01N, 411:01, 412:01, 417:01:01-417:01:02, 436:01, 440:01, 442:01, 449:01, 462:01, 470:01, 481:01, 483:01, 490:01, 499:01, 503:01, 514:01, 516:01, 518:01, 519:01, 536:01, 538:01, 542:01, 543:01, 544:01, 562:01, 564:01, 565:01, 597:01, 600:01, 609:01, 612:01, 616:01, 626:01, 631:01, 632:01, 636:01, 661:01N, 667:01, 703:01, 708:01, 709:01, 710:01, 711:01, 720:01, 733:01, 800:01N, 807:01, 810:01, 822:01, 825:01, 853:01, 862:01N, 866:01N, 870:01N, 873:01N, 879:01:01:01-879:01:01:03, 888:01, 893:01, 904:01, 916:01, 924:01, 931:01, 947:01, 953:01, 960:01N, 969:01, 976:01, 979:01, 984:01N, 991:01, 996:01, 998:01, 999:01, 1008:01, 1012:01, 1024:01, 1038:01Q, 1045:01N, 1050:01, 1057:01, 1065:01, 1068:01, 1069:01, 1076:01, 1088:01, 1105:01, 1123:01, 1131:01, 1137:01, 1140:01, 1141:01, 1147:01, 1150:01, 1151:01, 1162:01, 1166:01, 1169:01N, 1183:01, 1185:01
8	210 bp	515 bp	*02:01:02:01-02:02:04, 04:01:01:01-04:02:17, 05:01:01:01-05:01:14, 08:01, 16:01:01:01-16:01:03, 19:01:01:01-19:01:01:03, 22:01:01:01-22:01:01:02, 23:01:01:01-23:01:02, 24:01, 32:01, 33:01:01:01-33:01:01:05, 38:01, 39:01:01:01-39:01:03, 40:01:01:01-40:01:02, 46:01:01:01-46:01:02, 47:01:01:01-47:01:01:03, 48:01, 49:01:01:01-49:01:01:03, 51:01:01:01-51:01:01:02, 53:01, 60:01, 62:01, 63:01, 65:01:01-65:01:02, 68:01, 71:01:01-71:01:02, 75:01, 77:01, 80:01, 81:01:01:01-81:01:02, 82:01, 84:01, 94:01, 96:01, 97:01, 99:01, 100:01, 101:01, 103:01, 105:01:01:01-105:01:01:12, 106:01, 112:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

115:01, 116:01, 120:01N, 121:01, 123:01, 126:01:01:01-126:01:01:03,
128:01, 134:01, 135:01, 138:01, 140:01, 141:01, 143:01, 144:01,
145:01, 146:01, 148:01, 149:01, 150:01, 151:01, 153:01, 154:01N,
155:01:01-155:01:02, 163:01, 165:01, 169:01, 170:01, 171:01, 172:01,
173:01, 174:01, 175:01, 176:01, 177:01, 178:01, 180:01, 181:01,
182:01, 184:01, 185:01, 186:01, 187:01, 188:01, 190:01, 191:01,
192:01, 193:01, 194:01, 195:01, 196:01, 210:01, 211:01, 212:01,
213:01:01-213:01:02, 214:01, 215:01, 216:01N, 217:01, 218:01N,
219:01, 224:01, 225:01, 227:01:01-227:01:03, 228:01, 230:01, 231:01,
232:01, 233:01, 234:01, 236:01:01-236:01:02, 237:01, 238:01, 239:01,
240:01, 252:01, 253:01, 254:01, 255:01, 256:01, 257:01, 258:01,
260:01, 261:01, 262:01, 263:01, 272:01, 273:01, 274:01, 275:01,
276:01, 277:01, 278:01, 279:01:01-279:01:02, 280:01, 281:01, 282:01,
283:01, 284:01, 286:01, 297:01, 298:01, 300:01, 301:01, 302:01,
303:01, 304:01, 305:01, 307:01, 308:01, 310:01, 311:01, 317:01,
318:01, 320:01, 321:01, 323:01, 330:01, 331:01, 332:01, 334:01:01-
334:01:02, 335:01, 336:01, 337:01, 338:01, 339:01, 340:01, 341:01,
342:01, 350:01, 352:01:01-352:01:02:02, 354:01:01-354:01:02,
356:01, 357:01N, 358:01, 359:01, 360:01, 364:01, 365:01, 366:01,
367:01, 368:01, 372:01, 373:01, 374:01, 375:01, 376:01, 377:01,
378:01, 380:01, 381:01, 382:01N, 395:01, 396:01, 397:01, 398:01,
400:01, 401:01N, 402:01, 406:01, 414:01:01:01-414:01:01:02, 415:01,
416:01:01:01-416:01:02, 418:01, 419:01, 420:01, 421:01, 422:01,
423:01:01-423:01:02, 424:01, 425:01, 426:01, 427:01, 428:01, 429:01,
430:01, 431:01, 432:01, 433:01, 434:01, 435:01, 441:01, 443:01,
444:01, 453:01, 454:01, 455:01N, 456:01, 457:01, 459:01, 461:01,
463:01:01:01-463:01:01:03, 464:01, 465:01, 466:01, 467:01, 468:01,
469:01, 473:01, 474:01, 475:01, 476:01, 477:01, 478:01, 479:01,
485:01, 486:01, 487:01, 488:01, 489:01, 494:01, 495:01, 496:01,
497:01, 500:01, 501:01, 502:01, 510:01, 511:01, 512:01, 513:01,
520:01, 521:01, 522:01, 523:01:01-523:01:02, 524:01, 525:01, 526:01,
527:01, 528:01, 529:01, 533:01, 534:01, 535:01, 537:01, 538:01,
539:01, 547:01, 549:01, 551:01N, 553:01, 554:01, 556:01, 557:01,
558:01, 559:01, 560:01, 569:01, 570:01N, 571:01, 573:01, 575:01,
577:01, 578:01, 579:01, 582:01, 583:01, 584:01:01:01-584:01:02:02,
586:01, 587:01, 588:01, 589:01, 590:01, 591:01, 592:01, 593:01,
594:01, 595:01, 596:01, 598:01N, 602:01, 603:01, 604:01, 605:01,
606:01, 607:01, 611:01, 614:01, 615:01, 617:01, 618:01, 619:01,
622:01, 624:01, 625:01, 626:01, 627:01, 628:01, 637:01, 639:01,
640:01, 641:01, 642:01, 643:01, 646:01, 647:01, 652:01, 653:01,
655:01, 656:01, 657:01N, 658:01, 659:01, 660:01, 665:01, 666:01,
668:01:01:01-668:01:01:02, 670:01, 674:01, 677:01, 678:01, 680:01,
681:01, 682:01, 683:01, 684:01, 685:01, 686:01, 687:01, 690:01,
692:01, 693:01N, 694:01, 696:01N, 697:01Q, 699:01, 700:01N,
701:01, 702:01, 720:01, 721:01, 722:01, 723:01, 724:01N, 725:01,
726:01, 729:01, 730:01, 731:01, 734:01, 735:01, 736:01, 739:01,
740:01, 741:01, 742:01, 743:01N, 744:01, 745:01, 747:01, 748:01N,
750:01, 751:01, 753:01, 754:01N, 755:01, 756:01N, 757:01, 758:01,
759:01, 761:01, 762:01, 763:01, 764:01, 765:01, 766:01, 767:01,
769:01, 770:01, 771:01, 772:01, 773:01, 774:01, 775:01, 776:01,
777:01N, 778:01, 779:01, 780:01, 781:01, 783:01, 784:01, 785:01,
786:01:01N-786:01:02N, 787:01, 788:01, 789:01, 790:01, 795:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			796:01, 798:01, 799:01, 802:01, 803:01, 804:01, 805:01, 806:01:01:01-806:01:01:02, 808:01, 809:01, 811:01, 812:01, 813:01, 814:01, 816:01, 817:01, 818:01, 819:01, 820:01, 823:01, 824:01:01-824:01:02, 826:01, 827:01, 828:01, 830:01, 831:01N, 836:01, 837:01, 838:01N, 840:01, 841:01, 842:01, 843:01, 844:01N, 845:01, 847:01, 848:01, 849:01, 850:01, 851:01, 852:01, 857:01, 858:01, 859:01, 860:01, 861:01, 863:01, 864:01, 865:01N, 867:01N, 868:01N, 869:01N, 871:01N, 872:01N, 874:01N, 875:01N, 877:01N, 878:01N, 880:01, 881:01, 882:01, 884:01, 885:01, 886:01, 887:01, 889:01, 890:01, 894:01N, 895:01, 898:01, 904:01, 915:01, 916:01, 918:01, 919:01N, 921:01, 922:01, 923:01, 925:01N, 926:01, 927:01, 928:01, 929:01, 932:01, 933:01, 934:01Q, 936:01Q, 939:01, 940:01, 951:01, 954:01, 955:01, 957:01, 958:01, 959:01N, 962:01, 963:01, 964:01, 965:01:01:01-965:01:01:02, 966:01, 967:01, 968:01, 971:01, 972:01, 973:01, 974:01N, 975:01, 978:01, 981:01, 982:01, 985:01N, 988:01, 989:01, 992:01, 993:01, 994:01, 995:01N, 997:01, 1001:01, 1002:01, 1003:01, 1004:01, 1005:01, 1010:01, 1011:01, 1013:01, 1015:01, 1016:01, 1018:01, 1020:01, 1023:01, 1025:01, 1026:01, 1028:01, 1031:01, 1033:01, 1035:01, 1036:01, 1037:01, 1039:01, 1040:01, 1042:01, 1044:01N, 1046:01, 1048:01, 1051:01, 1053:01, 1055:01, 1056:01, 1058:01, 1059:01, 1060:01, 1062:01, 1064:01, 1066:01, 1070:01N, 1072:01, 1074:01, 1075:01, 1077:01, 1078:01, 1079:01N, 1082:01, 1083:01, 1084:01N, 1085:01, 1086:01, 1089:01, 1091:01, 1092:01Q, 1094:01, 1095:01, 1097:01, 1099:01, 1100:01, 1101:01, 1102:01, 1106:01, 1107:01, 1108:01, 1110:01, 1112:01N, 1113:01, 1115:01, 1117:01, 1118:01, 1119:01, 1120:01, 1121:01N, 1122:01, 1124:01, 1129:01, 1132:01, 1135:01N, 1136:01, 1139:01, 1140:01, 1143:01, 1144:01, 1146:01, 1148:01Q, 1150:01, 1152:01, 1153:01, 1154:01N, 1155:01, 1156:01, 1159:01, 1160:01, 1161:01, 1163:01, 1164:01, 1167:01, 1171:01, 1172:01, 1173:01, 1175:01, 1176:01, 1177:01, 1179:01, 1180:01, 1181:01, 1184:01, 1188:01, 1189:01, 1191:01N, 1193:01N, 1194:01, 1195:01
9⁵	165 bp	430 bp	*02:01:02:01-02:01:04, 02:01:05 ^w , 02:01:07-02:01:43, 02:01:45-02:02:04, 04:02:01:01-04:02:17, 05:01:01:01-05:01:09, 05:01:11-05:01:14, 06:01:01:01-06:01:05, 16:01:01:01-16:01:03, 17:01:01:01-17:01:04, 18:01:01:01-18:01:03, 19:01:01:01-19:01:03, 20:01:01:01-20:01:04, 21:01, 22:01:01:01-22:01:01:02, 23:01:01:01-23:01:02, 30:01:01:01-30:01:01:02, 32:01, 34:01:01:01-34:01:02, 36:01, 38:01, 41:01:01:01-41:01:02, 46:01:01:01-46:01:02, 47:01:01:01-47:01:01:03, 48:01, 55:01:01:01-55:01:03, 58:01, 59:01, 60:01, 62:01, 63:01, 64:01N, 69:01:01:01-69:01:01:02, 71:01:01-71:01:02, 77:01, 80:01, 82:01, 83:01, 86:01, 87:01, 91:01:01:01-91:01:02, 93:01, 94:01, 95:01, 97:01, 98:01, 100:01, 101:01, 103:01, 105:01:01:01-105:01:01:12, 106:01, 108:01, 109:01, 112:01, 113:01, 115:01, 116:01, 123:01, 129:01, 130:01, 131:01:01:01-131:01:03, 135:01, 136:01, 137:01, 138:01, 140:01, 141:01, 143:01, 144:01, 145:01, 146:01, 148:01, 151:01, 152:01, 153:01, 154:01N, 155:01:01-155:01:02, 161:01N, 163:01, 164:01, 165:01, 166:01, 168:01, 170:01, 171:01, 182:01, 183:01, 184:01, 185:01, 186:01, 189:01, 190:01, 191:01, 196:01, 197:01, 200:01, 208:01, 210:01, 211:01, 217:01, 218:01N, 219:01, 223:01, 226:01, 227:01:03, 233:01, 235:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

236:01:01-236:01:02, 237:01, 238:01, 239:01, 241:01, 244:01, 245:01,
251:01, 252:01, 256:01, 257:01, 258:01, 259:01, 260:01, 261:01,
263:01, 271:01, 273:01, 274:01^w, 284:01, 286:01, 292:01, 295:01,
297:01, 300:01, 301:01, 307:01, 308:01, 309:01, 310:01, 311:01,
317:01, 321:01, 322:01, 325:01, 330:01, 331:01, 332:01, 333:01,
338:01, 339:01, 340:01, 342:01, 343:01^w, 347:01, 349:01, 352:01:01-
352:01:02:02, 358:01, 359:01, 361:01, 364:01, 365:01, 366:01,
367:01, 368:01, 369:01, 371:01, 373:01, 374:01^w, 379:01, 380:01,
381:01, 382:01N, 384:01, 392:01, 395:01, 400:01, 402:01, 406:01,
407:01, 408:01, 413:01, 414:01:01:01-414:01:01:02, 416:01:01:01-
416:01:02, 419:01, 420:01, 421:01, 422:01, 424:01, 429:01, 430:01,
431:01, 432:01, 433:01, 435:01, 443:01, 444:01, 447:01, 448:01,
451:01, 452:01, 457:01, 458:01:01-458:01:02, 460:01, 461:01,
463:01:01:01-463:01:01:03, 466:01, 467:01, 469:01, 473:01, 477:01,
478:01, 484:01, 488:01, 489:01, 494:01, 495:01, 496:01, 501:01,
502:01, 510:01, 511:01, 513:01, 525:01, 526:01, 528:01, 530:01,
532:01, 533:01, 535:01, 537:01, 539:01, 540:01, 546:01, 547:01,
549:01, 550:01, 551:01N, 557:01, 558:01, 560:01, 563:01, 567:01,
570:01N, 571:01, 573:01, 576:01, 577:01, 579:01, 580:01, 581:01,
582:01, 586:01, 587:01, 588:01, 589:01, 590:01, 594:01, 595:01,
596:01, 598:01N, 601:01, 602:01, 604:01, 605:01, 606:01, 608:01,
611:01, 613:01, 617:01, 619:01, 624:01, 627:01, 628:01, 635:01,
637:01, 638:01, 639:01, 640:01, 646:01, 647:01, 652:01, 653:01,
655:01, 656:01, 657:01N, 659:01, 660:01, 662:01, 665:01,
668:01:01:01-668:01:01:02, 674:01, 678:01, 680:01, 681:01, 684:01,
685:01, 689:01, 690:01, 692:01, 693:01N, 697:01Q, 698:01, 700:01N,
701:01, 717:01, 719:01, 721:01, 723:01, 724:01N, 725:01, 726:01,
729:01, 730:01, 731:01, 734:01, 735:01, 737:01, 738:01N, 740:01,
741:01, 744:01, 746:01, 751:01, 759:01, 762:01, 763:01, 764:01,
766:01, 771:01, 774:01, 775:01, 776:01, 777:01N, 778:01, 779:01,
780:01, 781:01, 782:01, 783:01, 785:01, 786:01:01N-786:01:02N,
787:01, 789:01, 790:01, 791:01, 796:01, 798:01, 799:01, 802:01,
805:01, 809:01, 817:01, 818:01, 819:01, 823:01, 827:01, 830:01,
831:01N, 835:01:01:01-835:01:01:03, 836:01, 841:01, 843:01, 845:01,
846:01, 847:01, 848:01, 851:01, 857:01, 858:01, 860:01, 861:01,
863:01, 864:01, 865:01N, 869:01N, 872:01N, 877:01N, 881:01,
884:01, 885:01, 886:01, 887:01, 889:01^w, 890:01, 891:01, 894:01N,
897:01, 898:01, 905:01, 906:01, 908:01, 913:01, 914:01, 919:01N,
920:01, 923:01, 925:01N, 927:01, 928:01, 929:01, 933:01, 936:01Q,
940:01, 942:01, 944:01, 951:01, 954:01, 955:01, 956:01, 958:01,
959:01N, 961:01, 962:01, 963:01, 964:01, 965:01:01:01-965:01:01:02,
967:01, 968:01, 971:01, 972:01, 973:01, 975:01, 981:01, 982:01,
985:01N, 994:01, 995:01N, 1005:01, 1013:01, 1015:01, 1017:01,
1018:01, 1019:01, 1020:01, 1022:01, 1025:01, 1026:01, 1028:01,
1031:01, 1032:01, 1035:01, 1036:01, 1037:01, 1047:01, 1051:01,
1052:01, 1053:01, 1055:01, 1056:01, 1058:01, 1059:01, 1061:01,
1067:01, 1070:01N, 1072:01, 1075:01, 1077:01, 1078:01, 1082:01,
1083:01, 1084:01N, 1085:01, 1087:01, 1089:01, 1092:01Q, 1094:01,
1095:01, 1099:01, 1101:01, 1102:01, 1104:01, 1106:01, 1107:01,
1111:01, 1115:01, 1117:01, 1118:01, 1119:01, 1120:01, 1124:01,
1127:01, 1142:01, 1143:01, 1153:01, 1154:01N, 1156:01, 1158:01,
1159:01, 1160:01, 1163:01, 1165:01, 1170:01, 1171:01, 1172:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			1174:01, 1175:01, 1176:01, 1179:01, 1180:01, 1186:01, 1188:01, 1189:01, 1190:01, 1191:01N, 1194:01
10	220 bp	515 bp	*02:01:02:01-02:01:50, 02:01:52-02:02:04, 08:01, 16:01:01:01-16:01:03, 19:01:01:01-19:01:03, 22:01:01:01-22:01:02, 32:01, 33:01:01:01-33:01:05, 41:01:01:01-41:01:02, 46:01:01:01-46:01:02, 47:01:01:01-47:01:03, 48:01, 71:01:01-71:01:02, 81:01:01:01-81:01:02, 95:01, 101:01, 103:01, 106:01, 115:01, 123:01, 141:01, 145:01, 146:01, 148:01, 151:01, 152:01, 154:01N, 155:01:01-155:01:02, 163:01, 165:01, 166:01, 172:01, 173:01, 182:01, 186:01, 187:01, 188:01, 189:01, 191:01, 196:01, 211:01, 217:01, 218:01N, 229:01, 234:01, 235:01, 236:01:01-236:01:02, 237:01, 238:01, 239:01, 258:01, 261:01, 273:01, 277:01, 297:01, 302:01, 304:01, 305:01, 310:01, 311:01, 332:01, 339:01, 340:01, 352:01:01-352:01:02:02, 358:01, 359:01, 366:01, 367:01, 381:01, 382:01N, 400:01, 401:01N, 406:01, 414:01:01:01-414:01:02, 416:01:01:01-416:01:02, 419:01, 422:01, 423:01:01-423:01:02, 432:01, 433:01, 441:01, 444:01, 457:01, 461:01, 467:01, 469:01, 474:01, 478:01, 489:01, 496:01, 501:01, 502:01, 510:01, 511:01, 513:01, 526:01, 528:01, 533:01, 535:01, 537:01, 547:01, 551:01N, 555:01, 557:01, 558:01, 570:01N, 576:01, 579:01, 581:01, 582:01, 590:01, 602:01, 603:01, 604:01, 617:01, 622:01, 626:01, 628:01, 637:01, 639:01, 640:01, 641:01, 652:01, 653:01, 655:01, 656:01, 660:01, 678:01, 680:01, 681:01, 684:01, 685:01, 692:01, 700:01N, 720:01, 721:01, 723:01, 740:01, 751:01, 766:01, 770:01, 771:01, 779:01, 780:01, 781:01, 783:01, 785:01, 786:01:01N-786:01:02N, 789:01, 796:01, 799:01, 816:01, 819:01, 827:01, 830:01, 831:01N, 836:01, 845:01, 857:01, 861:01, 863:01, 864:01, 865:01N, 869:01N, 872:01N, 877:01N, 884:01, 886:01, 889:01, 890:01, 894:01N, 904:01, 919:01N, 928:01, 940:01, 944:01, 955:01, 963:01, 964:01, 965:01:01:01-965:01:01:02, 967:01, 968:01, 973:01, 975:01, 982:01, 1026:01, 1028:01, 1036:01, 1051:01, 1053:01, 1055:01, 1056:01, 1059:01, 1077:01, 1082:01, 1089:01, 1094:01, 1101:01, 1102:01, 1107:01, 1115:01, 1136:01, 1140:01, 1154:01N, 1159:01, 1160:01, 1163:01, 1175:01, 1176:01, 1179:01, 1188:01, 1191:01N
11	145 bp 185 bp	515 bp	*69:01:01:01-69:01:02, 136:01, 179:01, 207:01, 303:01, 309:01, 322:01, 458:01:01-458:01:02, 644:01, 659:01, 930:01, 1109:01 *02:01:02:01-02:01:42, 02:01:44-02:01:57, 04:01:01:01-04:01:06, 04:01:08-04:01:10, 04:01:12-04:02:17, 23:01:01:01-23:01:02, 24:01, 32:01, 33:01:01:01-33:01:05, 41:01:01:01-41:01:02, 46:01:01:01-46:01:02, 47:01:01:01-47:01:03, 51:01:01:01-51:01:02, 59:01, 60:01, 66:01, 71:01:01-71:01:02, 72:01:01:01-72:01:03, 73:01, 75:01, 77:01, 80:01, 81:01:01:01-81:01:02, 82:01, 83:01, 86:01, 94:01, 99:01, 105:01:01:01-105:01:01:12, 108:01, 109:01, 113:01, 115:01, 119:01, 120:01N, 121:01, 123:01, 126:01:01:01-126:01:01:03, 128:01, 129:01, 134:01, 138:01, 141:01, 143:01, 144:01, 145:01, 146:01, 149:01, 151:01, 153:01, 154:01N, 155:01:01-155:01:02, 163:01, 164:01, 169:01, 174:01, 175:01, 176:01, 177:01, 178:01, 179:01, 180:01, 181:01, 183:01, 185:01, 186:01, 187:01, 188:01, 189:01, 190:01, 191:01, 192:01, 193:01, 194:01, 195:01, 196:01, 199:01, 200:01, 210:01, 211:01, 212:01, 213:01:01-213:01:02, 214:01, 215:01, 216:01N, 217:01, 218:01N, 219:01, 224:01, 225:01, 227:01:01-227:01:02, 228:01, 229:01, 231:01, 232:01, 233:01, 235:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

236:01:01-236:01:02, 237:01, 238:01, 239:01, 240:01, 252:01, 253:01, 254:01, 255:01, 256:01, 257:01, 258:01, 260:01, 261:01, 262:01, 263:01, 271:01, 272:01, 273:01, 275:01, 276:01, 278:01, 281:01, 282:01, 283:01, 286:01, 294:01, 297:01, 298:01, 302:01, 303:01, 306:01, 307:01, 308:01, 309:01, 310:01, 311:01, 318:01, 319:01, 320:01, 321:01, 322:01, 323:01, 334:01:01-334:01:02, 335:01, 336:01, 338:01, 339:01, 341:01, 342:01, 344:01, 350:01, 352:01:01-352:01:02:02, 355:01, 356:01, 357:01N, 359:01, 360:01, 364:01, 365:01, 366:01, 367:01, 368:01, 369:01, 372:01, 373:01, 374:01, 375:01, 376:01, 377:01, 378:01, 380:01, 381:01, 382:01N, 392:01, 396:01, 397:01, 399:01, 401:01N, 402:01, 414:01:01:01-414:01:01:02, 415:01, 416:01:01:01-416:01:02, 418:01, 419:01, 420:01, 421:01, 423:01:01-423:01:02, 424:01, 425:01, 426:01, 427:01, 428:01, 429:01, 430:01, 432:01, 433:01, 434:01, 435:01, 443:01, 444:01, 451:01, 452:01, 453:01, 455:01N, 456:01, 457:01, 459:01, 461:01, 463:01:01:01-463:01:01:03, 464:01, 465:01, 468:01, 469:01, 475:01, 476:01, 477:01, 479:01, 485:01, 486:01, 487:01, 488:01, 489:01, 494:01, 497:01, 500:01, 501:01, 502:01, 508:01, 510:01, 511:01, 520:01, 521:01, 522:01, 523:01:01-523:01:02, 524:01, 525:01, 528:01, 529:01, 534:01, 537:01, 539:01, 540:01, 549:01, 551:01N, 553:01, 554:01, 555:01, 556:01, 557:01, 559:01, 569:01, 570:01N, 571:01, 574:01, 575:01, 576:01, 577:01, 578:01, 579:01, 581:01, 582:01, 583:01, 586:01, 591:01, 592:01, 593:01, 594:01, 595:01, 596:01, 601:01, 602:01, 603:01, 604:01, 605:01, 606:01, 607:01, 614:01, 615:01, 617:01, 618:01, 620:01, 622:01, 624:01, 625:01, 626:01, 627:01, 628:01, 637:01, 639:01, 640:01, 641:01, 642:01, 643:01, 646:01, 647:01, 655:01, 656:01, 658:01, 659:01, 660:01, 665:01, 666:01, 670:01, 674:01, 677:01, 678:01, 679:01, 680:01, 681:01, 682:01, 683:01, 685:01, 687:01, 690:01, 692:01, 693:01N, 694:01, 696:01N, 699:01, 700:01N, 701:01, 702:01, 722:01, 723:01, 724:01N, 725:01, 726:01, 730:01, 731:01, 734:01, 735:01, 736:01, 739:01, 740:01, 741:01, 742:01, 743:01N, 747:01, 748:01N, 750:01, 751:01, 753:01, 754:01N, 755:01, 756:01N, 757:01, 758:01, 759:01, 761:01, 762:01, 763:01, 765:01, 767:01, 769:01, 770:01, 771:01, 772:01, 773:01, 774:01, 775:01, 776:01, 777:01N, 780:01, 781:01, 782:01, 783:01, 784:01, 786:01:01N-786:01:02N, 787:01, 788:01, 795:01, 796:01, 799:01, 803:01, 804:01, 805:01, 806:01:01:01-806:01:01:02, 808:01, 809:01, 811:01, 812:01, 813:01, 814:01, 817:01, 818:01, 819:01, 820:01, 823:01, 824:01:01-824:01:02, 826:01, 827:01, 828:01, 831:01N, 836:01, 837:01, 838:01N, 840:01, 841:01, 842:01, 843:01, 844:01N, 845:01, 849:01, 850:01, 852:01, 857:01, 858:01, 859:01, 861:01, 863:01, 867:01N, 868:01N, 869:01N, 871:01N, 872:01N, 874:01N, 875:01N, 877:01N, 878:01N, 880:01, 881:01, 882:01, 884:01, 885:01, 887:01, 889:01, 890:01, 894:01N, 895:01, 911:01N, 915:01, 919:01N, 921:01, 925:01N, 926:01, 927:01, 932:01, 933:01, 936:01Q, 939:01, 954:01, 955:01, 957:01, 958:01, 959:01N, 963:01, 964:01, 966:01, 967:01, 972:01, 973:01, 974:01N, 975:01, 978:01, 981:01, 985:01N, 987:01, 988:01, 989:01, 992:01, 993:01, 997:01, 1001:01, 1002:01, 1003:01, 1004:01, 1005:01, 1010:01, 1011:01, 1013:01, 1016:01, 1020:01, 1023:01, 1025:01, 1028:01, 1031:01, 1033:01, 1035:01, 1036:01, 1037:01, 1040:01, 1042:01, 1044:01N, 1048:01, 1051:01, 1053:01, 1055:01, 1056:01, 1059:01, 1060:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			1062:01, 1064:01, 1066:01, 1072:01, 1074:01, 1075:01, 1077:01, 1078:01, 1079:01N, 1081:01, 1082:01, 1083:01, 1084:01N, 1085:01, 1086:01, 1089:01, 1091:01, 1094:01, 1097:01, 1100:01, 1102:01, 1104:01, 1106:01, 1107:01, 1108:01, 1110:01, 1112:01N, 1113:01, 1115:01, 1121:01N, 1124:01, 1129:01, 1132:01, 1135:01N, 1138:01, 1139:01, 1144:01, 1146:01, 1148:01Q, 1152:01, 1153:01, 1154:01N, 1155:01, 1159:01, 1160:01, 1161:01, 1163:01, 1164:01, 1167:01, 1171:01, 1173:01, 1174:01, 1175:01, 1176:01, 1177:01, 1179:01, 1181:01, 1184:01, 1191:01N, 1194:01, 1195:01
12^{4,6}	80 bp	515 bp	*02:01:02:01-02:01:50, 02:01:52-02:01:57, 06:01:01:01-06:01:05, 08:01, 09:01:01-09:01:05, 10:01:01:01-10:01:05, 16:01:01:01-16:01:03, 17:01:01:01-17:01:04, 29:01, 32:01, 37:01, 41:01:01:01-41:01:02, 44:01, 46:01:01:01-46:01:02, 48:01, 64:01N, 81:01:01:01-81:01:02, 86:01, 88:01, 93:01, 109:01, 113:01, 115:01, 123:01, 131:01:01:01-131:01:03, 137:01, 141:01, 145:01, 146:01, 151:01, 152:01, 154:01N, 155:01:01-155:01:02, 156:01, 157:01, 161:01N, 163:01, 165:01, 166:01, 167:01, 168:01, 172:01, 182:01, 186:01, 187:01, 188:01, 189:01, 191:01, 196:01, 198:01, 200:01, 208:01, 211:01, 217:01, 218:01N, 229:01, 234:01, 235:01, 236:01:01-236:01:02, 237:01, 238:01, 239:01, 241:01, 245:01, 249:01, 251:01, 261:01, 265:01, 266:01, 269:01, 273:01, 277:01, 287:01, 288:01, 289:01, 297:01, 310:01, 311:01, 313:01, 324:01, 332:01, 339:01, 343:01, 344:01, 347:01, 349:01, 352:01:01-352:01:02:02, 359:01, 361:01, 366:01, 367:01, 381:01, 382:01N, 384:01, 386:01, 401:01N, 407:01, 410:01, 413:01, 414:01:01:01-414:01:01:02, 416:01:01:01-416:01:02, 419:01, 422:01, 423:01:01-423:01:02, 432:01, 433:01, 441:01, 444:01, 447:01, 448:01, 452:01, 457:01, 460:01, 461:01, 467:01, 469:01, 470:01, 489:01, 501:01, 502:01, 504:01, 505:01, 510:01, 511:01, 513:01, 514:01, 515:01, 526:01, 528:01, 530:01, 537:01, 546:01, 551:01N, 552:01, 555:01, 557:01, 566:01, 567:01, 570:01N, 579:01, 581:01, 582:01, 601:01, 602:01, 603:01, 604:01, 608:01, 617:01, 622:01, 628:01, 630:01, 633:01, 637:01, 639:01, 640:01, 641:01, 650:01, 652:01, 653:01, 655:01, 656:01, 660:01, 662:01, 663:01, 673:01, 678:01, 680:01, 681:01, 684:01, 685:01, 692:01, 698:01, 700:01N, 711:01, 714:01, 719:01, 723:01, 737:01, 738:01N, 740:01, 751:01, 770:01, 771:01, 780:01, 781:01, 783:01, 786:01:01N-786:01:02N, 791:01, 797:01, 799:01, 815:01, 816:01, 819:01, 827:01, 831:01N, 833:01, 836:01, 845:01, 846:01, 857:01, 861:01, 863:01, 864:01, 865:01N, 869:01N, 872:01N, 877:01N, 884:01, 886:01, 889:01, 890:01, 894:01N, 899:01, 902:01, 906:01, 909:01, 914:01, 919:01N, 920:01, 940:01, 944:01, 955:01, 956:01, 963:01, 967:01, 968:01, 973:01, 975:01, 977:01, 1006:01, 1007:01, 1009:01, 1017:01, 1022:01, 1028:01, 1032:01, 1036:01, 1043:01, 1047:01, 1051:01, 1052:01, 1053:01, 1055:01, 1056:01, 1059:01, 1067:01, 1077:01, 1082:01, 1087:01, 1089:01, 1094:01, 1102:01, 1107:01, 1111:01, 1115:01, 1126:01, 1127:01, 1128:01, 1136:01, 1140:01, 1149:01, 1158:01, 1159:01, 1160:01, 1163:01, 1168:01, 1175:01, 1176:01, 1178:01, 1179:01, 1182:01, 1191:01N
13	180 bp	430 bp	*02:02:01:01-02:02:04, 05:01:01:01-05:01:02, 05:01:04-05:01:14, 19:01:01:01-19:01:01:03, 22:01:01:01-22:01:01:02, 24:01, 38:01, 47:01:01:01-47:01:01:03, 97:01, 100:01, 101:01, 106:01, 135:01, 139:01, 140:01, 170:01, 171:01, 226:01, 233:01, 284:01, 300:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			301:01, 302:01, 317:01, 330:01, 331:01, 337:01, 358:01, 395:01, 400:01, 406:01, 473:01, 478:01, 495:01, 496:01, 527:01, 533:01, 535:01, 547:01, 558:01, 560:01, 573:01, 587:01, 588:01, 589:01, 590:01, 598:01N, 611:01, 619:01, 638:01, 657:01N, 668:01:01:01-668:01:01:02, 697:01Q, 720:01, 721:01, 729:01, 744:01, 764:01, 766:01, 778:01, 779:01, 785:01, 790:01, 798:01, 802:01, 847:01, 848:01, 851:01, 860:01, 923:01, 928:01, 929:01, 951:01, 962:01, 965:01:01:01-965:01:01:02, 971:01, 982:01, 1015:01, 1018:01, 1026:01, 1058:01, 1061:01, 1070:01N, 1095:01, 1099:01, 1101:01, 1117:01, 1118:01, 1119:01, 1120:01, 1143:01, 1156:01, 1172:01, 1180:01, 1188:01, 1189:01
14⁴	105 bp	430 bp	*02:02:01:01-02:02:04, 05:01:01:01-05:01:02, 05:01:04-05:01:14, 21:01, 22:01:01:01-22:01:01:02, 36:01, 38:01, 97:01, 100:01, 101:01, 114:01, 135:01, 140:01, 170:01, 171:01, 223:01, 226:01, 300:01, 301:01, 317:01, 330:01, 331:01, 395:01, 406:01, 473:01, 478:01, 495:01, 496:01, 547:01, 550:01, 560:01, 573:01, 587:01, 588:01, 589:01, 590:01, 598:01N, 611:01, 619:01, 638:01, 657:01N, 668:01:01:01-668:01:01:02, 697:01Q, 721:01, 729:01, 744:01, 746:01, 764:01, 766:01, 778:01, 779:01, 790:01, 798:01, 802:01, 847:01, 848:01, 851:01, 860:01, 923:01, 928:01, 929:01, 951:01, 961:01, 962:01, 971:01, 1015:01, 1018:01, 1019:01, 1026:01, 1058:01, 1061:01, 1070:01N, 1095:01, 1099:01, 1117:01, 1118:01, 1119:01, 1120:01, 1143:01, 1156:01, 1172:01, 1180:01, 1186:01, 1188:01, 1190:01
	175 bp		*21:01, 30:01:01:01-30:01:01:02, 36:01, 54:01, 114:01, 390:01, 408:01, 550:01, 715:01, 717:01, 718:01, 746:01, 961:01, 980:01, 1008:01, 1019:01, 1186:01, 1190:01
15	200 bp	430 bp	*02:02:01:01-02:02:03, 48:01, 95:01, 99:01, 100:01, 101:01, 112:01, 116:01, 148:01, 227:01:03, 332:01, 406:01, 478:01, 496:01, 547:01, 590:01, 656:01, 721:01, 766:01, 779:01, 830:01, 928:01, 961:01, 962:01, 1138:01, 1188:01
16⁴	80 bp	430 bp	*03:01:01:01-03:01:14, 04:02:01:01-04:02:01:20, 04:02:02, 04:02:03-04:02:17, 14:01:01:01-14:01:12, 18:01:01:01-18:01:03, 20:01:01:01-20:01:04, 25:01, 28:01, 35:01:01, 45:01, 49:01:01:01-49:01:01:03, 50:01, 51:01:01:01-51:01:01:02, 53:01, 57:01, 59:01, 60:01, 61:01N, 68:01, 70:01, 73:01, 75:01, 76:01, 77:01, 78:01, 79:01, 80:01, 82:01, 83:01, 91:01:01:01-91:01:01:02, 92:01, 98:01, 104:01:01:01-104:01:04, 105:01:01:01-105:01:01:12, 108:01, 111:01, 116:01, 119:01, 122:01, 124:01:01:01-124:01:02:03, 129:01, 130:01, 132:01, 143:01, 144:01, 153:01, 164:01, 183:01, 184:01, 185:01, 197:01, 203:01:01-203:01:02, 210:01, 214:01, 219:01, 221:01, 222:01, 242:01, 243:01, 248:01, 252:01, 256:01, 257:01, 259:01, 260:01, 263:01, 270:01, 271:01, 274:01, 285:01, 286:01, 292:01, 293:01, 296:01, 307:01, 308:01, 312:01, 316:01, 321:01, 329:01, 338:01, 342:01, 351:01, 355:01, 363:01, 364:01, 365:01, 368:01, 369:01, 373:01, 374:01, 379:01, 380:01, 383:01, 385:01, 391:01, 394:01, 402:01, 404:01, 405:01, 409:01, 420:01, 424:01, 429:01, 431:01, 439:01, 442:01, 443:01, 445:01, 446:01, 463:01:01:01-463:01:01:03, 466:01, 472:01, 477:01, 484:01 ^w , 488:01, 491:01, 492:01, 494:01, 498:01, 499:01, 506:01, 507:01N, 508:01, 509:01, 525:01, 532:01, 536:01, 539:01, 540:01, 545:01, 548:01, 549:01 ^w , 568:01, 571:01, 572:01, 574:01, 577:01, 586:01, 594:01, 595:01, 596:01, 605:01, 606:01, 610:01, 613:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			620:01, 621:01, 624:01, 627:01, 629:01, 645:01, 646:01, 647:01, 648:01:01:01-648:01:01:02, 651:01, 664:01, 665:01, 669:01, 671:01, 674:01, 675:01, 676:01, 688:01, 689:01, 690:01, 691:01N, 693:01N, 701:01, 704:01, 705:01, 706:01, 712:01N, 724:01N, 725:01, 726:01, 727:01, 728:01, 730:01, 731:01, 734:01, 735:01, 736:01, 741:01, 752:01, 759:01, 760:01, 762:01, 763:01, 775:01, 776:01, 777:01N, 782:01, 792:01N, 794:01N, 801:01, 805:01, 809:01, 817:01, 818:01, 823:01, 829:01, 832:01, 834:01, 839:01, 841:01, 854:01, 855:01, 858:01, 876:01N, 881:01, 883:01, 885:01, 887:01, 891:01, 892:01, 897:01, 898:01, 900:01, 905:01, 917:01N, 918:01, 925:01N, 927:01, 933:01, 935:01Q, 936:01Q, 938:01, 941:01N, 942:01, 943:01, 946:01, 948:01, 949:01, 950:01N, 952:01, 954:01, 958:01, 959:01N, 970:01, 981:01, 983:01, 985:01N, 986:01N, 987:01, 990:01, 994:01, 995:01N, 1000:01, 1005:01, 1013:01, 1014:01, 1020:01, 1021:01, 1025:01, 1027:01, 1030:01, 1031:01, 1035:01, 1037:01, 1039:01, 1041:01N, 1049:01, 1071:01, 1072:01, 1075:01, 1083:01, 1084:01N, 1085:01, 1090:01, 1092:01Q, 1093:01, 1098:01N, 1103:01, 1104:01, 1106:01, 1110:01, 1114:01, 1116:01, 1124:01, 1125:01, 1130:01, 1134:01, 1142:01, 1150:01, 1153:01, 1157:01, 1165:01, 1171:01, 1174:01, 1187:01, 1194:01
	165 bp		*99:01, 220:01, 1138:01
	210 bp		*127:01
17	260 bp	515 bp	*03:01:01:01-03:01:14, 06:01:01:01-06:01:05, 09:01:01-09:01:05, 10:01:01:01-10:01:05, 11:01:01:01-11:01:05, 13:01:01:01-13:01:04, 14:01:01:01-14:01:12, 17:01:01:01-17:01:04, 20:01:01:01-20:01:04, 21:01, 25:01, 26:01:01-26:01:03, 27:01, 29:01, 30:01:01:01-30:01:01:02, 35:01:01, 36:01, 37:01, 44:01, 45:01, 52:01, 54:01, 55:01:01:01-55:01:03, 56:01, 58:01, 61:01N, 64:01N, 67:01, 69:01:01:01-69:01:01:02, 70:01, 76:01, 78:01, 79:01, 85:01:01:01-85:01:02, 87:01, 88:01, 91:01:01:01-91:01:02, 92:01, 93:01, 98:01, 104:01:01:01-104:01:04, 107:01, 110:01, 111:01, 114:01, 119:01, 124:01:01:01-124:01:02:03, 125:01, 130:01, 131:01:01:01-131:01:03, 132:01, 133:01, 136:01, 137:01, 142:01, 156:01, 157:01, 167:01, 168:01, 197:01, 202:01, 203:01:01-203:01:02, 204:01, 205:01, 206:01, 207:01, 208:01, 209:01, 221:01, 222:01, 241:01, 243:01, 244:01, 245:01, 250:01, 251:01, 264:01, 265:01, 266:01, 268:01, 269:01, 270:01, 287:01, 288:01, 289:01, 293:01, 295:01, 312:01, 313:01, 314:01, 324:01, 325:01, 328:01N, 329:01, 343:01, 348:01:01-348:01:02, 349:01, 351:01, 353:01, 361:01, 362:01, 363:01, 370:01, 371:01, 383:01, 384:01, 385:01, 386:01, 388:01, 389:01, 390:01, 391:01, 393:01, 394:01, 404:01, 405:01, 407:01, 408:01, 409:01, 410:01, 412:01, 413:01, 437:01, 438:01, 439:01, 445:01, 446:01, 447:01, 448:01, 449:01, 450:01N, 458:01:01-458:01:02, 460:01, 470:01, 472:01, 490:01, 491:01, 492:01, 498:01, 503:01, 504:01, 505:01, 506:01, 509:01, 514:01, 515:01, 517:01, 518:01, 519:01, 530:01, 532:01, 541:01, 544:01, 545:01, 548:01, 550:01, 552:01, 562:01, 563:01, 564:01, 565:01, 566:01, 567:01, 568:01, 572:01, 597:01, 600:01, 609:01, 610:01, 612:01, 613:01, 621:01, 623:01, 629:01, 630:01, 633:01, 634:01, 635:01, 636:01, 645:01, 649:01, 650:01, 651:01, 654:01, 662:01, 664:01, 669:01, 671:01, 672:01, 673:01, 675:01, 676:01, 688:01, 689:01, 691:01N, 698:01, 703:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			704:01, 705:01, 706:01, 707:01, 708:01, 709:01, 712:01N, 713:01, 714:01, 715:01, 716:01, 717:01, 718:01, 719:01, 727:01, 728:01, 732:01N, 737:01, 738:01N, 746:01, 749:01, 760:01, 791:01, 792:01N, 794:01N, 797:01, 801:01, 815:01, 822:01, 825:01, 829:01, 832:01, 833:01, 834:01, 839:01, 846:01, 854:01, 855:01, 856:01, 873:01N, 876:01N, 883:01, 888:01, 891:01, 892:01, 893:01, 899:01, 901:01, 902:01, 905:01, 906:01, 907:01, 908:01, 909:01, 912:01, 914:01, 917:01N, 920:01, 924:01, 930:01, 935:01Q, 937:01, 938:01, 941:01N, 946:01, 947:01, 948:01, 949:01, 950:01N, 952:01, 956:01, 970:01, 977:01, 980:01, 983:01, 986:01N, 990:01, 991:01, 996:01, 1000:01, 1007:01, 1008:01, 1009:01, 1014:01, 1017:01, 1019:01, 1021:01, 1022:01, 1027:01, 1029:01N, 1030:01, 1032:01, 1034:01, 1041:01N, 1043:01, 1049:01, 1052:01, 1057:01, 1063:01, 1065:01, 1067:01, 1069:01, 1071:01, 1087:01, 1088:01, 1090:01, 1093:01, 1096:01, 1098:01N, 1103:01, 1105:01, 1111:01, 1114:01, 1116:01, 1123:01, 1125:01, 1126:01, 1127:01, 1128:01, 1130:01, 1131:01, 1133:01, 1134:01, 1145:01, 1147:01, 1149:01, 1157:01, 1158:01, 1166:01, 1168:01, 1170:01, 1178:01, 1182:01, 1185:01, 1186:01, 1187:01, 1190:01
18	180 bp	430 bp	*03:01:01:01-03:01:14, 06:01:01:01-06:01:05, 09:01:01:09:01:05, 14:01:01:01-14:01:12, 17:01:01:01-17:01:04, 20:01:01:01-20:01:04, 29:01, 35:01:01, 44:01, 61:01N, 64:01N, 69:01:01:01-69:01:01:02, 70:01, 76:01, 78:01, 86:01, 88:01, 91:01:01:01-91:01:02, 92:01, 98:01, 104:01:01:01-104:01:04, 111:01, 119:01, 124:01:01:01-124:01:02:03, 130:01, 131:01:01:01-131:01:03, 132:01, 156:01, 157:01, 161:01N, 164:01, 168:01, 197:01, 203:01:01-203:01:02, 205:01, 208:01, 221:01, 222:01, 241:01, 242:01, 243:01, 245:01, 251:01, 266:01, 270:01, 287:01, 288:01, 289:01, 293:01, 329:01, 343:01, 351:01, 361:01, 363:01, 383:01, 384:01, 385:01, 386:01, 391:01, 394:01, 404:01, 405:01, 407:01, 409:01, 413:01, 439:01, 445:01, 446:01, 447:01, 460:01, 472:01, 491:01, 492:01, 498:01, 504:01, 505:01, 506:01, 507:01N, 508:01, 509:01, 530:01, 540:01, 541:01, 545:01, 546:01, 548:01, 566:01, 567:01, 568:01, 572:01, 601:01, 610:01, 613:01, 621:01, 629:01, 630:01, 645:01, 651:01, 662:01, 664:01, 669:01, 671:01, 675:01, 676:01, 688:01, 689:01, 691:01N, 698:01, 704:01, 705:01, 706:01, 712:01N, 714:01, 719:01, 727:01, 728:01, 737:01, 738:01N, 760:01, 792:01N, 794:01N, 797:01, 801:01, 815:01, 829:01, 833:01, 834:01, 839:01, 846:01, 854:01, 855:01, 876:01N, 883:01, 899:01, 905:01, 906:01, 909:01, 914:01, 917:01N, 920:01, 935:01Q, 938:01, 941:01N, 946:01, 948:01, 949:01, 950:01N, 952:01, 956:01, 970:01, 977:01, 983:01, 986:01N, 987:01, 990:01, 1000:01, 1009:01, 1014:01, 1017:01, 1021:01, 1022:01, 1027:01, 1029:01N, 1030:01, 1032:01, 1041:01N, 1043:01, 1049:01, 1052:01, 1067:01, 1071:01, 1087:01, 1090:01, 1093:01, 1096:01, 1098:01N, 1103:01, 1104:01, 1111:01, 1114:01, 1116:01, 1125:01, 1127:01, 1128:01, 1130:01, 1133:01, 1134:01, 1149:01, 1157:01, 1158:01, 1174:01, 1178:01, 1182:01, 1187:01
19	200 bp	430 bp	*03:01:01:01-03:01:14, 06:01:01:01-06:01:05, 11:01:01:01-11:01:05, 14:01:01:01-14:01:12, 20:01:01:01-20:01:04, 25:01, 29:01, 44:01, 45:01, 52:01, 56:01, 61:01N, 64:01N, 67:01, 69:01:01:01-69:01:01:02, 70:01, 74:01, 76:01, 78:01, 87:01, 91:01:01:01-91:01:02, 92:01, 104:01:01:01-104:01:04, 111:01, 114:01, 119:01,



101.311-24/04 – including Taq pol., IFU-01
101.311-24u/04u – without Taq pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			124:01:01:01-124:01:02:03, 130:01, 136:01, 142:01, 161:01N, 164:01, 167:01, 197:01, 202:01, 203:01:01-203:01:02, 204:01, 205:01, 207:01, 208:01, 221:01, 222:01, 242:01, 243:01, 250:01, 264:01, 266:01, 268:01, 270:01, 293:01, 328:01N, 329:01, 348:01:01-348:01:02, 349:01, 351:01, 353:01, 363:01, 371:01, 383:01, 385:01, 388:01, 391:01, 394:01, 404:01, 405:01, 409:01, 437:01, 439:01, 445:01, 446:01, 450:01N, 458:01:01-458:01:02, 472:01, 491:01, 492:01, 498:01, 506:01, 507:01N, 508:01, 509:01, 517:01, 531:01, 546:01, 548:01, 566:01, 567:01, 568:01, 572:01, 601:01, 610:01, 629:01, 645:01, 649:01, 651:01, 654:01, 662:01, 664:01, 669:01, 671:01, 672:01, 675:01, 676:01, 688:01, 689:01, 691:01N, 704:01, 705:01, 706:01, 707:01, 712:01N, 715:01, 718:01, 727:01, 728:01, 732:01N, 737:01, 749:01, 760:01, 792:01N, 794:01N, 801:01, 829:01, 832:01, 834:01, 839:01, 854:01, 855:01, 876:01N, 883:01, 892:01, 905:01, 906:01, 907:01, 909:01, 914:01, 917:01N, 930:01, 935:01Q, 937:01, 938:01, 941:01N, 946:01, 948:01, 949:01, 950:01N, 952:01, 970:01, 983:01, 986:01N, 987:01, 990:01, 1000:01, 1014:01, 1021:01, 1022:01, 1029:01N, 1030:01, 1041:01N, 1049:01, 1063:01, 1071:01, 1087:01, 1090:01, 1093:01, 1096:01, 1098:01N, 1103:01, 1104:01, 1111:01, 1114:01, 1116:01, 1125:01, 1127:01, 1130:01, 1133:01, 1134:01, 1157:01, 1187:01
20	235 bp	515 bp	*03:01:01:01-03:01:14, 09:01:01-09:01:05, 10:01:01:01-10:01:05, 13:01:01:01-13:01:04, 14:01:01:01-14:01:12, 25:01, 26:01:01-26:01:03, 29:01, 35:01:01, 37:01, 44:01, 45:01, 52:01, 54:01, 56:01, 61:01N, 67:01, 70:01, 76:01, 78:01, 79:01, 88:01, 92:01, 104:01:01:01-104:01:04, 107:01, 110:01, 111:01, 114:01, 119:01, 124:01:01:01-124:01:02:03, 125:01, 132:01, 133:01, 142:01, 157:01, 160:01, 167:01, 197:01, 203:01:01-203:01:02, 209:01, 221:01, 222:01, 242:01, 265:01, 266:01, 269:01, 270:01, 288:01, 289:01, 293:01, 295:01, 312:01, 313:01, 314:01, 324:01, 329:01, 344:01, 351:01, 362:01, 363:01, 383:01, 385:01, 386:01, 389:01, 390:01, 393:01, 394:01, 404:01, 405:01, 409:01, 410:01, 412:01, 439:01, 445:01, 446:01, 448:01, 449:01, 470:01, 472:01, 490:01, 491:01, 492:01, 498:01, 503:01, 504:01, 506:01, 507:01N, 508:01, 509:01, 514:01, 515:01, 518:01, 519:01, 530:01, 541:01, 544:01, 545:01, 548:01, 552:01, 562:01, 564:01, 565:01, 567:01, 568:01, 572:01, 597:01, 600:01, 609:01, 610:01, 612:01, 621:01, 629:01, 630:01, 634:01, 636:01, 645:01, 650:01, 651:01, 664:01, 669:01, 671:01, 673:01, 675:01, 676:01, 688:01, 689:01, 691:01N, 703:01, 704:01, 705:01, 706:01, 708:01, 709:01, 714:01, 715:01, 718:01, 727:01, 728:01, 746:01, 760:01, 792:01N, 794:01N, 797:01, 801:01, 815:01, 822:01, 825:01, 829:01, 832:01, 833:01, 834:01, 839:01, 854:01, 855:01, 873:01N, 876:01N, 883:01, 888:01, 893:01, 899:01, 902:01, 908:01, 909:01, 912:01, 917:01N, 920:01, 924:01, 935:01Q, 938:01, 941:01N, 946:01, 947:01, 948:01, 949:01, 950:01N, 952:01, 970:01, 977:01, 983:01, 986:01N, 987:01, 990:01, 991:01, 996:01, 1000:01, 1007:01, 1008:01, 1009:01, 1014:01, 1021:01, 1027:01, 1029:01N, 1030:01, 1041:01N, 1043:01, 1045:01N, 1049:01, 1057:01, 1065:01, 1069:01, 1071:01, 1088:01, 1090:01, 1093:01, 1096:01, 1098:01N, 1103:01, 1105:01, 1114:01, 1116:01, 1123:01, 1125:01, 1126:01, 1128:01, 1130:01, 1131:01, 1133:01, 1134:01, 1147:01, 1149:01, 1157:01, 1166:01, 1168:01, 1178:01, 1182:01, 1185:01, 1187:01



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

21⁴	100 bp	430 bp	*03:01:01:01-03:01:14, 09:01:01-09:01:05, 14:01:01:01-14:01:12, 29:01, 35:01:01, 44:01, 50:01, 57:01, 61:01N, 70:01, 76:01, 78:01, 88:01, 92:01, 104:01:01:01-104:01:04, 111:01, 119:01, 124:01:01:01-124:01:02:03, 132:01, 157:01, 166:01, 197:01, 203:01:01-203:01:02, 221:01, 222:01, 234:01, 242:01, 248:01, 249:01, 266:01, 270:01, 288:01, 289:01, 293:01, 329:01, 351:01, 363:01, 379:01, 383:01, 385:01, 386:01, 394:01, 404:01, 405:01, 409:01, 439:01, 442:01, 445:01, 446:01, 472:01, 491:01, 492:01, 498:01, 504:01, 506:01, 507:01N, 508:01, 509:01, 530:01, 536:01, 541:01, 545:01, 548:01, 567:01, 568:01, 572:01, 610:01, 620:01, 621:01, 629:01, 630:01, 645:01, 648:01:01:01-648:01:01:02, 651:01, 664:01, 669:01, 671:01, 675:01, 676:01, 684:01, 688:01, 689:01, 691:01N, 704:01, 705:01, 706:01, 714:01, 727:01, 728:01, 760:01, 792:01N, 794:01N, 797:01, 801:01, 815:01, 829:01, 833:01, 834:01, 839:01, 854:01, 855:01, 876:01N, 883:01, 899:01, 909:01, 917:01N, 920:01, 935:01Q, 938:01, 941:01N, 946:01, 948:01, 949:01, 950:01N, 952:01, 970:01, 977:01, 983:01, 986:01N, 987:01, 990:01, 1000:01, 1009:01, 1014:01, 1021:01, 1027:01, 1030:01, 1041:01N, 1043:01, 1049:01, 1071:01, 1090:01, 1093:01, 1098:01N, 1103:01, 1114:01, 1116:01, 1125:01, 1128:01, 1130:01, 1134:01, 1149:01, 1157:01, 1178:01, 1182:01, 1187:01
22	185 bp	430 bp	*04:01:01:01-04:01:06, 04:01:08-04:01:10, 04:01:12-04:01:52, 24:01, 33:01:01:01-33:01:01:05, 39:01:01:01-39:01:03, 49:01:01:01-49:01:01:03, 51:01:01:01-51:01:01:02, 66:01, 72:01:01:01-72:01:01:03, 81:01:01:01-81:01:02, 96:01, 99:01, 117:01, 120:01N, 121:01, 126:01:01:01-126:01:01:03, 128:01, 134:01, 149:01, 158:01, 169:01, 172:01, 174:01, 175:01, 176:01, 177:01, 178:01, 179:01, 180:01, 181:01, 192:01, 193:01, 194:01, 195:01, 199:01, 212:01, 213:01:01-213:01:02, 214:01, 215:01, 216:01N, 224:01, 225:01, 227:01:03, 228:01, 229:01, 231:01, 232:01, 240:01, 253:01, 254:01, 255:01, 262:01, 272:01, 275:01, 276:01, 278:01, 281:01, 282:01, 283:01, 294:01, 298:01, 302:01, 303:01, 306:01, 319:01, 320:01, 323:01, 334:01:01-334:01:02, 335:01, 336:01, 341:01, 350:01, 354:01:01-354:01:02, 355:01, 356:01, 357:01N, 360:01, 372:01, 375:01, 376:01, 377:01, 378:01, 396:01, 397:01, 399:01, 415:01, 418:01, 425:01, 426:01, 427:01, 428:01, 434:01, 440:01, 453:01, 455:01N, 456:01, 459:01, 464:01, 465:01, 468:01, 474:01, 475:01, 476:01, 479:01, 480:01, 485:01, 486:01, 487:01, 497:01, 500:01, 520:01, 521:01, 522:01, 523:01:01-523:01:02, 524:01, 529:01, 531:01, 534:01, 553:01, 554:01, 555:01, 556:01, 559:01, 561:01, 569:01, 574:01, 575:01, 578:01, 583:01, 584:01:01:01-584:01:02:02, 591:01, 592:01, 593:01, 607:01, 614:01, 615:01, 618:01, 625:01, 626:01, 642:01, 643:01, 658:01, 663:01, 666:01, 670:01, 677:01, 679:01, 682:01, 683:01, 686:01, 687:01, 694:01, 696:01N, 699:01, 702:01, 722:01, 736:01, 739:01, 742:01, 743:01N, 747:01, 748:01N, 750:01, 753:01, 754:01N, 755:01, 756:01N, 757:01, 758:01, 761:01, 765:01, 767:01, 769:01, 772:01, 773:01, 784:01, 788:01, 795:01, 803:01, 804:01, 806:01:01:01-806:01:01:02, 808:01, 811:01, 812:01, 813:01, 814:01, 820:01, 824:01:01-824:01:02, 826:01, 828:01, 837:01, 838:01N, 840:01, 842:01, 844:01N, 849:01, 850:01, 852:01, 859:01, 867:01N, 868:01N, 871:01N, 874:01N, 875:01N, 878:01N, 880:01, 882:01, 895:01, 911:01N, 915:01, 916:01, 921:01, 926:01, 932:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			939:01, 943:01, 957:01, 966:01, 974:01N, 978:01, 988:01, 989:01, 992:01, 993:01, 997:01, 1001:01, 1002:01, 1003:01, 1004:01, 1010:01, 1011:01, 1016:01, 1023:01, 1033:01, 1039:01, 1040:01, 1042:01, 1044:01N, 1046:01, 1048:01, 1060:01, 1062:01, 1064:01, 1066:01, 1074:01, 1079:01N, 1080:01, 1081:01, 1086:01, 1091:01, 1097:01, 1100:01, 1108:01, 1112:01N, 1113:01, 1121:01N, 1129:01, 1132:01, 1135:01N, 1138:01, 1139:01, 1144:01, 1146:01, 1148:01Q, 1152:01, 1155:01, 1161:01, 1164:01, 1167:01, 1173:01, 1177:01, 1181:01, 1184:01, 1195:01
23⁷	145 bp	430 bp	*05:01:01:01-05:01:14, 34:01:01:01-34:01:02, 36:01, 38:01, 62:01, 63:01, 97:01, 100:01, 112:01, 114:01, 116:01, 135:01, 140:01, 170:01, 171:01, 223:01, 226:01, 227:01:03, 300:01, 301:01, 317:01, 330:01, 331:01, 333:01, 394:01, 395:01, 473:01, 484:01, 495:01, 560:01, 573:01, 580:01, 587:01, 589:01, 598:01N, 599:01, 611:01, 619:01, 638:01, 657:01N, 668:01:01:01-668:01:01:02, 697:01Q, 729:01, 744:01, 764:01, 778:01, 790:01, 798:01, 802:01, 835:01:01:01-835:01:01:03, 847:01, 848:01, 851:01, 860:01, 898:01, 913:01, 923:01, 929:01, 951:01, 961:01, 962:01, 971:01, 994:01, 1015:01, 1018:01, 1058:01, 1061:01, 1070:01N, 1095:01, 1099:01, 1117:01, 1118:01, 1119:01, 1120:01, 1143:01, 1156:01, 1172:01, 1180:01
	210 bp		*127:01, 383:01
24	195 bp 265 bp	515 bp	*101:01 *05:01:01:01-05:01:14, 08:01, 16:01:01:01-16:01:03, 19:01:01:01-19:01:01:03, 22:01:01:01-22:01:01:02, 31:01:01:01-31:01:01:02, 38:01, 57:01, 63:01, 65:01:01-65:01:02, 68:01, 84:01, 97:01, 102:01, 103:01, 106:01, 135:01, 140:01, 147:01, 150:01, 152:01, 165:01, 166:01, 170:01, 171:01, 173:01, 182:01, 184:01, 226:01, 234:01, 259:01, 277:01, 284:01, 285:01, 300:01, 301:01, 304:01, 305:01, 317:01, 331:01, 337:01, 340:01, 358:01, 379:01, 395:01, 398:01, 400:01, 422:01, 442:01, 466:01, 473:01, 495:01, 527:01, 533:01, 535:01, 538:01, 558:01, 560:01, 573:01, 587:01, 588:01, 598:01N, 611:01, 619:01, 638:01, 648:01:01:01-648:01:01:02, 652:01, 653:01, 657:01N, 668:01:01:01-668:01:01:02, 684:01, 697:01Q, 720:01, 729:01, 740:01, 744:01, 764:01, 778:01, 785:01, 789:01, 790:01, 798:01, 802:01, 847:01, 848:01, 851:01, 860:01, 864:01, 886:01, 898:01, 904:01, 922:01, 923:01, 929:01, 940:01, 944:01, 945:01, 951:01, 965:01:01:01-965:01:01:02, 968:01, 971:01, 976:01, 982:01, 994:01, 1015:01, 1018:01, 1026:01, 1058:01, 1061:01, 1095:01, 1099:01, 1101:01, 1117:01, 1118:01, 1119:01, 1120:01, 1140:01, 1143:01, 1150:01, 1156:01, 1172:01, 1180:01, 1189:01
25⁷	200 bp	515 bp	*05:01:01:01-05:01:14, 21:01, 22:01:01:01-22:01:01:02, 36:01, 38:01, 44:01, 58:01, 63:01, 97:01, 114:01, 135:01, 140:01, 170:01, 171:01, 220:01, 223:01, 226:01, 265:01, 300:01, 301:01, 317:01, 331:01, 384:01, 394:01, 395:01, 473:01, 495:01, 550:01, 560:01, 563:01, 573:01, 587:01, 588:01, 598:01N, 599:01, 611:01, 619:01, 630:01, 638:01, 657:01N, 668:01:01:01-668:01:01:02, 697:01Q, 729:01, 744:01, 746:01, 764:01, 778:01, 789:01, 790:01, 798:01, 802:01, 847:01, 848:01, 851:01, 860:01, 898:01, 923:01, 929:01, 951:01, 971:01, 994:01, 1015:01, 1018:01, 1019:01, 1026:01, 1058:01, 1061:01, 1095:01, 1099:01, 1117:01, 1118:01, 1119:01, 1120:01, 1143:01, 1156:01, 1172:01, 1180:01, 1186:01, 1190:01



101.311-24/04 – including *Taq pol.*, IFU-01
 101.311-24u/04u – without *Taq pol.*, IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

	265 bp		*70:01
26 ^{4,6}	75 bp	515 bp	*06:01:01:01-06:01:05, 09:01:01-09:01:05, 17:01:01:01-17:01:04, 29:01, 44:01, 46:01:01:01-46:01:02, 64:01N, 86:01, 88:01, 101:01, 131:01:01:01-131:01:03, 152:01, 156:01, 157:01, 161:01N, 166:01, 168:01, 182:01, 208:01, 234:01, 235:01, 241:01, 245:01, 249:01, 251:01, 266:01, 287:01, 288:01, 289:01, 332:01, 343:01, 361:01, 384:01, 386:01, 407:01, 413:01, 447:01, 460:01, 504:01, 505:01, 530:01, 546:01, 555:01, 566:01, 567:01, 581:01, 601:01, 630:01, 662:01, 684:01, 698:01, 714:01, 719:01, 737:01, 738:01N, 797:01, 815:01, 833:01, 846:01, 899:01, 906:01, 909:01, 914:01, 920:01, 944:01, 956:01, 977:01, 1009:01, 1017:01, 1022:01, 1032:01, 1043:01, 1047:01, 1052:01, 1067:01, 1087:01, 1111:01, 1127:01, 1128:01, 1149:01, 1158:01, 1178:01, 1182:01
	205 bp		*123:01
27 ^{4,5}	95 bp	430 bp	*06:01:01:01-06:01:05, 17:01:01:01-17:01:04, 20:01:01:01-20:01:04, 46:01:01:01-46:01:02, 64:01N, 69:01:01:01-69:01:01:02, 80:01, 86:01, 91:01:01:01-91:01:02, 98:01, 108:01, 130:01, 131:01:01:01-131:01:03, 152:01, 156:01, 161:01N, 164:01, 166:01, 168:01, 182:01, 197:01, 205:01, 208:01, 214:01, 214:01, 235:01, 241:01, 245:01, 251:01, 259:01, 287:01, 292:01, 332:01, 343:01, 355:01, 361:01, 379:01, 384:01, 407:01, 413:01, 447:01, 460:01, 484:01, 530:01, 540:01, 546:01, 555:01, 566:01, 567:01, 581:01, 601:01, 613:01, 662:01, 684:01, 689:01, 698:01, 719:01, 737:01, 738:01N, 762:01, 846:01, 905:01, 906:01, 914:01, 920:01, 944:01, 956:01, 994:01, 1017:01, 1022:01, 1032:01, 1047:01, 1052:01, 1067:01, 1087:01, 1104:01, 1111:01, 1127:01, 1158:01, 1174:01
	215 bp		*169:01, 252:01, 582:01, 772:01
	240 bp		*15:01:01:01-15:01:01:04, 15:01:03-15:01:04, 18:01:01:01-18:01:03, 89:01, 158:01, 198:01, 199:01, 200:01, 219:01, 246:01, 247:01, 290:01, 292:01, 345:01, 347:01, 369:01, 387:01, 471:01, 480:01, 482:01, 493:01, 561:01, 585:01:01:01-585:01:01:02, 608:01, 644:01, 747:01, 768:01, 821:01N, 896:01, 897:01, 900:01, 903:01, 910:01, 942:01, 1047:01, 1054:01, 1165:01, 1192:01
28	200 bp	430 bp	*09:01:01-09:01:05, 10:01:01:01-10:01:05, 13:01:01:01-13:01:04, 17:01:01:01-17:01:04, 21:01, 26:01:01-26:01:03, 27:01, 30:01:01:01-30:01:01:02, 35:01:01, 36:01, 37:01, 54:01, 55:01:01:01-55:01:03, 58:01, 66:01, 79:01, 85:01:01:01-85:01:02, 86:01, 88:01, 93:01, 98:01, 107:01, 109:01, 110:01, 113:01, 117:01, 125:01, 129:01, 131:01:01:01-131:01:03, 132:01, 133:01, 137:01, 156:01, 157:01, 160:01, 168:01, 206:01, 209:01, 241:01, 244:01, 245:01, 251:01, 265:01, 269:01, 271:01, 287:01, 288:01, 289:01, 294:01, 295:01, 312:01, 313:01, 314:01, 324:01, 325:01, 343:01, 344:01, 361:01, 362:01, 370:01, 384:01, 386:01, 389:01, 390:01, 392:01, 393:01, 407:01, 408:01, 410:01, 412:01, 413:01, 438:01, 447:01, 448:01, 449:01, 460:01, 470:01, 490:01, 503:01, 504:01, 505:01, 514:01, 515:01, 518:01, 519:01, 530:01, 532:01, 540:01, 541:01, 544:01, 550:01, 552:01, 562:01, 563:01, 564:01, 565:01, 597:01, 600:01, 609:01, 612:01, 613:01, 621:01, 623:01, 630:01, 633:01, 634:01, 635:01, 636:01, 650:01, 663:01, 673:01, 698:01, 703:01, 708:01, 709:01, 713:01, 714:01, 716:01, 717:01, 719:01, 738:01N, 746:01, 791:01, 797:01, 815:01, 822:01, 825:01, 846:01, 873:01N, 888:01, 891:01, 893:01, 899:01, 901:01, 902:01, 908:01, 912:01, 920:01,



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

			924:01, 943:01, 947:01, 956:01, 961:01, 977:01, 980:01, 991:01, 996:01, 1007:01, 1008:01, 1009:01, 1017:01, 1019:01, 1027:01, 1032:01, 1034:01, 1043:01, 1045:01N, 1052:01, 1057:01, 1065:01, 1067:01, 1069:01, 1080:01, 1088:01, 1105:01, 1123:01, 1126:01, 1128:01, 1131:01, 1145:01, 1147:01, 1149:01, 1158:01, 1166:01, 1168:01, 1170:01, 1174:01, 1178:01, 1182:01, 1185:01, 1186:01, 1190:01
	260 bp		*119:01
29⁴	105 bp 160 bp 220 bp	515 bp	*141:01, 763:01, 826:01 *617:01, 699:01 *09:01:01-09:01:05, 10:01:01:01-10:01:05, 17:01:01:01-17:01:04, 30:01:01:01-30:01:01:02, 54:01, 55:01:01:01-55:01:03, 58:01, 86:01, 113:01, 131:01:01:01-131:01:03, 137:01, 156:01, 157:01, 167:01, 168:01, 241:01, 245:01, 265:01, 266:01, 287:01, 288:01, 289:01, 313:01, 314:01, 324:01, 325:01, 343:01, 344:01, 361:01, 384:01, 386:01, 407:01, 408:01, 410:01, 447:01, 460:01, 470:01, 503:01, 504:01, 505:01, 514:01, 515:01, 530:01, 630:01, 650:01, 673:01, 698:01, 714:01, 716:01, 719:01, 738:01N, 791:01, 797:01, 815:01, 833:01, 846:01, 899:01, 902:01, 908:01, 956:01, 977:01, 1007:01, 1009:01, 1017:01, 1032:01, 1043:01, 1052:01, 1067:01, 1126:01, 1128:01, 1149:01, 1158:01, 1168:01, 1178:01, 1182:01
30	185 bp	430 bp	*09:01:01-09:01:05, 14:01:01:01-14:01:11, 17:01:01:01-17:01:04, 35:01:01, 76:01, 86:01, 91:01:01:01-91:01:02, 98:01, 119:01, 131:01:01:01-131:01:03, 156:01, 157:01, 168:01, 197:01, 241:01, 242:01, 243:01, 245:01, 266:01, 287:01, 288:01, 289:01, 343:01, 361:01, 383:01, 384:01, 385:01, 386:01, 407:01, 409:01, 445:01, 446:01, 447:01, 460:01, 491:01, 492:01, 498:01, 504:01, 505:01, 530:01, 540:01, 541:01, 572:01, 629:01, 630:01, 651:01, 671:01, 698:01, 705:01, 714:01, 719:01, 738:01N, 794:01N, 797:01, 815:01, 833:01, 834:01, 846:01, 854:01, 899:01, 949:01, 956:01, 977:01, 986:01N, 987:01, 990:01, 1009:01, 1017:01, 1029:01N, 1032:01, 1043:01, 1052:01, 1067:01, 1104:01, 1128:01, 1149:01, 1158:01, 1178:01, 1182:01, 1187:01
31⁴	120 bp	515 bp	*01:01:01:01-01:01:07, 05:01:01:01-05:01:02, 05:01:04-05:01:12, 05:01:14, 08:01, 10:01:01:01-10:01:02, 10:01:04-10:01:05, 11:01:01:01-11:01:05, 13:01:01:01-13:01:04, 16:01:01:01-16:01:03, 19:01:01:01-19:01:01:03, 21:01, 22:01:01:01-22:01:01:02, 25:01, 26:01:01-26:01:03, 27:01, 30:01:01:01-30:01:01:02, 31:01:01:01-31:01:01:02, 36:01, 37:01, 38:01, 45:01, 52:01, 54:01, 55:01:01:01-55:01:03, 56:01, 58:01, 63:01, 65:01:01-65:01:02, 67:01, 68:01, 79:01, 84:01, 85:01:01:01-85:01:02, 87:01, 89:01, 90:01:01-90:01:02, 93:01, 97:01, 102:01, 103:01, 106:01, 107:01, 110:01, 114:01, 118:01, 122:01, 125:01, 127:01, 133:01, 135:01, 136:01, 137:01, 140:01, 142:01, 147:01, 150:01, 159:01N, 162:01:01-162:01:02, 165:01, 167:01, 170:01, 171:01, 173:01, 184:01, 201:01, 202:01, 204:01, 206:01, 207:01, 209:01, 220:01, 223:01, 226:01, 244:01, 246:01, 247:01, 250:01, 264:01, 265:01, 267:01, 268:01, 269:01, 277:01, 284:01, 285:01, 291:01, 295:01, 300:01, 301:01, 304:01, 305:01, 312:01, 313:01, 314:01, 315:01, 316:01, 317:01, 324:01, 325:01, 326:01, 327:01, 328:01N, 331:01, 337:01, 340:01, 346:01, 348:01:01-348:01:02, 349:01, 353:01, 358:01, 362:01, 370:01, 371:01, 388:01, 389:01, 390:01, 393:01, 395:01, 398:01, 400:01, 403:01N, 408:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			410:01, 411:01, 412:01, 417:01:01-417:01:02, 422:01, 436:01, 437:01, 438:01, 448:01, 449:01, 450:01N, 458:01:01-458:01:02, 462:01, 466:01, 470:01, 473:01, 481:01, 483:01, 490:01, 495:01, 503:01, 514:01, 515:01, 516:01, 517:01, 518:01, 519:01, 527:01, 532:01, 533:01, 535:01, 538:01, 542:01, 543:01, 544:01, 550:01, 552:01, 558:01, 560:01, 562:01, 563:01, 564:01, 565:01, 573:01, 587:01, 588:01, 597:01, 598:01N, 599:01, 600:01, 608:01, 609:01, 611:01, 612:01, 616:01, 619:01, 623:01, 631:01, 632:01, 633:01, 634:01, 635:01, 636:01, 638:01, 649:01, 650:01, 652:01, 653:01, 654:01, 657:01N, 661:01N, 667:01, 668:01:01:01-668:01:01:02, 672:01, 673:01, 697:01Q, 703:01, 707:01, 708:01, 709:01, 710:01, 711:01, 713:01, 715:01, 716:01, 717:01, 718:01, 720:01, 729:01, 732:01N, 733:01, 740:01, 744:01, 746:01, 749:01, 764:01, 778:01, 785:01, 789:01, 790:01, 791:01, 798:01, 800:01N, 802:01, 807:01, 810:01, 822:01, 825:01, 832:01, 847:01, 848:01, 851:01, 853:01, 856:01, 860:01, 862:01N, 864:01, 866:01N, 870:01N, 873:01N, 879:01:01:01-879:01:01:03, 886:01, 888:01, 891:01, 892:01, 893:01, 898:01, 901:01, 902:01, 904:01, 907:01, 908:01, 912:01, 922:01, 923:01, 924:01, 929:01, 930:01, 931:01, 937:01, 940:01, 945:01, 947:01, 951:01, 953:01, 960:01N, 965:01:01:01-965:01:01:02, 968:01, 969:01, 971:01, 976:01, 979:01, 980:01, 982:01, 984:01N, 991:01, 996:01, 998:01, 999:01, 1006:01, 1007:01, 1008:01, 1012:01, 1015:01, 1018:01, 1019:01, 1024:01, 1026:01, 1034:01, 1038:01Q, 1050:01, 1057:01, 1061:01, 1063:01, 1065:01, 1068:01, 1069:01, 1073:01, 1076:01, 1088:01, 1095:01, 1099:01, 1101:01, 1105:01, 1117:01, 1118:01, 1119:01, 1120:01, 1123:01, 1126:01, 1131:01, 1137:01, 1140:01, 1141:01, 1143:01, 1145:01, 1147:01, 1150:01, 1151:01, 1156:01, 1162:01, 1166:01, 1168:01, 1169:01N, 1170:01, 1172:01, 1180:01, 1183:01, 1185:01, 1186:01, 1189:01, 1190:01
	215 bp		*222:01
32^{4,7}	75 bp 200 bp	430 bp	*78:01 *11:01:01:01-11:01:01:02, 11:01:02-11:01:05, 130:01, 147:01, 202:01, 203:01:01-203:01:02, 204:01, 205:01, 250:01, 264:01, 268:01, 328:01N, 343:01, 353:01, 370:01, 388:01, 389:01, 405:01, 450:01N, 517:01, 623:01, 649:01, 654:01, 672:01, 707:01, 732:01N, 778:01, 856:01, 907:01, 930:01, 937:01, 1063:01, 1147:01
	260 bp		*223:01, 293:01, 1017:01
33⁴	105 bp 145 bp	430 bp	*76:01, 92:01, 156:01, 203:01:01-203:01:02, 205:01, 214:01, 215:01, 287:01, 288:01, 355:01, 363:01, 385:01, 442:01, 536:01, 555:01, 566:01 *13:01:01:01-13:01:04, 33:01:01:01-33:01:01:05, 81:01:01:01-81:01:01:02, 107:01, 110:01, 117:01, 133:01, 156:01, 160:01, 172:01, 173:01, 198:01, 206:01, 229:01, 246:01, 247:01, 269:01, 277:01, 287:01, 288:01, 302:01, 304:01, 305:01, 362:01, 389:01, 393:01, 438:01, 449:01, 470:01, 474:01, 503:01, 514:01, 518:01, 519:01, 544:01, 555:01, 561:01, 564:01, 565:01, 566:01, 597:01, 600:01, 609:01, 612:01, 626:01, 631:01, 633:01, 663:01, 703:01, 711:01, 716:01, 720:01, 822:01, 873:01N, 888:01, 893:01, 904:01, 924:01, 947:01, 991:01, 996:01, 1008:01, 1045:01N, 1057:01, 1065:01, 1069:01, 1105:01, 1123:01, 1131:01, 1136:01, 1140:01, 1145:01, 1185:01



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

Lot No.	Fragment Size(s)	Reference Size	Sequence Information
34	195 bp		*123:01, 178:01, 218:01N
	210 bp	430 bp	*14:01:01:01-14:01:11, 45:01, 67:01, 76:01, 91:01:01:01-91:01:02, 119:01, 167:01, 197:01, 242:01, 243:01, 264:01, 266:01, 383:01, 385:01, 409:01, 445:01, 446:01, 491:01, 492:01, 498:01, 572:01, 629:01, 651:01, 671:01, 705:01, 715:01, 794:01N, 832:01, 834:01, 854:01, 892:01, 949:01, 986:01N, 987:01, 990:01, 1029:01N, 1104:01, 1187:01
35	195 bp	430 bp	*15:01:01:01-15:01:04, 74:01, 198:01, 290:01, 299:01, 345:01, 387:01, 471:01, 493:01, 585:01:01:01-585:01:01:02, 695:01, 768:01, 821:01N, 896:01, 903:01, 910:01, 1054:01, 1192:01
36	260 bp	515 bp	*15:01:01:01-15:01:04, 18:01:01:01-18:01:03, 28:01, 34:01:01:01-34:01:02, 40:01:01:01-40:01:02, 53:01, 62:01, 139:01, 198:01, 290:01, 292:01, 296:01, 299:01, 333:01, 345:01, 347:01, 387:01, 471:01, 482:01, 484:01, 493:01, 499:01, 512:01, 526:01, 580:01, 585:01:01:01-585:01:01:02, 644:01, 695:01, 745:01, 768:01, 821:01N, 835:01:01:01-835:01:01:03, 896:01, 897:01, 900:01, 903:01, 910:01, 913:01, 918:01, 942:01, 1054:01, 1109:01, 1136:01, 1142:01, 1165:01, 1192:01
37 ⁶	160 bp	430 bp	*13:01:01:01-13:01:04, 107:01, 110:01, 133:01, 160:01, 173:01, 246:01, 269:01, 278:01, 304:01, 362:01, 389:01, 393:01, 412:01, 449:01, 503:01, 518:01, 519:01, 538:01, 544:01, 564:01, 565:01, 597:01, 600:01, 612:01, 703:01, 720:01, 822:01, 873:01N, 888:01, 893:01, 924:01, 947:01, 991:01, 996:01, 1008:01, 1045:01N, 1057:01, 1065:01, 1069:01, 1105:01, 1123:01, 1131:01, 1185:01
	185 bp		*15:01:01:01-15:01:04, 28:01, 40:01:01:01-40:01:02, 53:01, 74:01, 139:01, 198:01, 290:01, 296:01, 299:01, 345:01, 387:01, 471:01, 482:01, 493:01, 499:01, 512:01, 585:01:01:01-585:01:01:02, 644:01, 683:01, 695:01, 745:01, 752:01, 768:01, 821:01N, 896:01, 900:01, 903:01, 910:01, 918:01, 1054:01, 1109:01, 1136:01, 1192:01
38 ^{4,7}	95 bp	515 bp	*23:01:01:01-23:01:02, 34:01:01:01-34:01:02, 52:01, 55:01:01:01-55:01:03, 58:01, 62:01, 63:01, 67:01, 71:01:01-71:01:02, 87:01, 95:01, 103:01, 112:01, 121:01, 138:01, 201:01, 209:01, 227:01:03, 244:01, 295:01, 314:01, 322:01, 325:01, 333:01, 340:01, 371:01, 392:01, 435:01, 451:01, 458:01:01-458:01:02, 563:01, 576:01, 580:01, 599:01, 635:01, 787:01, 789:01, 830:01, 835:01:01:01-835:01:01:03, 908:01, 912:01, 913:01, 972:01, 1078:01, 1096:01, 1133:01, 1170:01
	170 bp 200 bp		*186:01, 243:01, 657:01N, 787:01, 1090:01 *57:01 [?] , 185:01, 611:01, 674:01, 1180:01
39	210 bp	515 bp	*28:01, 31:01:01:01-31:01:01:02, 34:01:01:01-34:01:02, 57:01, 59:01, 72:01:01:01-72:01:01:03, 73:01, 95:01, 102:01, 108:01, 139:01, 147:01, 152:01, 166:01, 179:01, 189:01, 235:01, 259:01, 285:01, 296:01, 299:01, 306:01, 309:01, 322:01, 333:01, 355:01, 379:01, 442:01, 451:01, 484:01, 499:01, 555:01, 574:01, 576:01, 580:01, 620:01, 648:01:01:01-648:01:01:02, 752:01, 782:01, 835:01:01:01-835:01:01:03, 913:01, 944:01, 945:01, 976:01, 1061:01, 1109:01, 1142:01
40	135 bp	430 bp	*275:01
	180 bp		*30:01:01:01-30:01:01:02, 54:01, 408:01, 715:01
	230 bp		*14:01:01:01-14:01:11, 35:01:01, 45:01, 66:01, 67:01, 76:01, 91:01:01:01-91:01:02, 98:01, 119:01, 125:01, 129:01, 197:01, 242:01, 243:01, 244:01, 312:01, 383:01, 385:01, 409:01, 445:01, 446:01, 490:01, 491:01, 492:01, 498:01, 503:01, 530:01, 540:01, 572:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

		629:01, 651:01, 671:01, 705:01, 715:01, 794:01N, 832:01, 834:01, 854:01, 891:01, 892:01, 908:01, 912:01, 949:01, 986:01N, 987:01, 990:01, 1029:01N, 1104:01, 1187:01	
41^{4,6,7}	75 bp 225 bp 250 bp 280 bp	515 bp	*77:01, 245:01, 559:01, 1095:01 *154:01N *31:01:01:01-31:01:01:02, 34:01:01:01-34:01:02, 94:01, 176:01, 484:01, 493:01, 575:01, 580:01, 805:01, 835:01:01:01-835:01:01:03, 843:01, 913:01, 945:01 *134:01, 177:01
42⁷	140 bp 185 bp 210 bp 245 bp 285 bp	430 bp	*97:01, 1108:01 *32:01 *176:01, 180:01, 805:01 *128:01 *82:01, 301:01, 339:01, 400:01
43⁵	200 bp	430 bp	*34:01:01:01-34:01:02, 62:01, 301:01, 333:01, 484:01, 580:01, 835:01:01:01-835:01:01:03, 913:01
44⁴	60 bp 115 bp 150 bp 185 bp 240 bp 270 bp	430 bp	*282:01, 619:01 *99:01, 118:01, 220:01, 437:01, 453:01, 521:01, 1138:01 *38:01, 619:01, 698:01, 815:01, 1099:01 *111:01, 118:01, 120:01N, 220:01, 437:01, 1182:01 *38:01, 241:01, 1099:01, 1130:01 *64:01N
45^{4,7}	120 bp 160 bp 190 bp 240 bp	430 bp	*46:01:01:01-46:01:02, 80:01, 86:01, 108:01, 164:01, 214:01, 235:01, *254:01, 332:01, 355:01, 508:01, 540:01, 555:01, 581:01, 601:01, 620:01, 762:01, 987:01, 1104:01, 1174:01 *253:01 *39:01:01:01-39:01:03, 49:01:01:01-49:01:01:03, 96:01, 117:01, 158:01, 172:01, 354:01:01-354:01:02, 440:01, 441:01, 474:01, 480:01, 531:01, 561:01, 584:01:01:01-584:01:02:02, 663:01, 916:01, 943:01, 1039:01, 1080:01 *115:01, 240:01, 962:01
46⁴	75 bp 210 bp	430 bp	*41:01:01:01-41:01:02, 60:01, 69:01:01:01-69:01:01:02, 83:01, 136:01, 204:01, 205:01, 229:01, 309:01, 581:01, 659:01, 903:01 *41:01:01:01-41:01:02, 60:01, 83:01, 96:01, 147:01, 175:01, 179:01, 229:01, 299:01, 303:01, 309:01, 322:01, 419:01, 429:01, 581:01, 659:01, 1109:01
47⁵	165 bp 215 bp	430 bp	*44:01, 114:01, 265:01, 394:01, 495:01, 599:01, 630:01, 746:01 *140:01
48^{4,5,10}	115 bp 210 bp 280 bp	430 bp	*56:01, 66:01, 76:01, 92:01, 110:01, 125:01, 156:01, 207:01, 288:01, 294:01, 470:01, 566:01, 709:01, 749:01, 930:01 *61:01N *85:01:01:01-85:01:02, 515:01, 713:01, 901:01, 1034:01
49^{4,7}	95 bp 130 bp 205 bp	430 bp	*149:01 *198:01, 290:01, 343:01, 354:01:01-354:01:02, 370:01, 375:01, 387:01, 389:01, 522:01, 623:01, 663:01, 903:01, 1039:01, 1147:01 *143:01, 251:01, 552:01, 840:01
50⁷	160 bp 210 bp	430 bp	*144:01, 146:01, 542:01, 637:01 *01:01:01:01-01:01:02:02, 01:01:04-01:01:07, 18:01:01:01-18:01:03, 89:01, 90:01:01-90:01:02, 118:01, 122:01, 127:01, 158:01, 159:01N, 162:01:01-162:01:02, 198:01, 199:01, 200:01, 201:01, 220:01, 246:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			247:01, 248:01, 249:01, 267:01, 290:01, 291:01, 292:01, 315:01, 316:01, 326:01, 327:01, 346:01, 347:01, 369:01, 387:01, 403:01N, 417:01:01-417:01:02, 436:01, 440:01, 462:01, 481:01, 483:01, 536:01, 542:01, 543:01, 561:01, 599:01, 608:01, 616:01, 631:01, 632:01, 661:01N, 667:01, 710:01, 711:01, 733:01, 747:01, 800:01N, 807:01, 810:01, 853:01, 862:01N, 866:01N, 870:01N, 879:01:01:01-879:01:01:03, 897:01, 900:01, 903:01, 931:01, 942:01, 953:01, 960:01N, 969:01, 979:01, 984:01N, 998:01, 999:01, 1006:01, 1012:01, 1024:01, 1038:01Q, 1047:01, 1050:01, 1068:01, 1073:01, 1076:01, 1137:01, 1151:01, 1162:01, 1165:01, 1169:01N, 1183:01
51	170 bp 225 bp	430 bp	*145:01, 382:01N *452:01, 570:01N
52	130 bp 160 bp 215 bp	430 bp	*151:01 *511:01 *165:01, 289:01, 452:01, 570:01N
53^{4,8,9}	90 bp 145 bp 175 bp 200 bp	430 bp	*157:01, 321:01, 801:01 *221:01, 501:01, 573:01, 595:01 *391:01, 467:01, 712:01N, 770:01 *153:01, 163:01, 242:01
54⁴	115 bp 160 bp 200 bp	430 bp	*161:01N, 505:01, 507:01N, 546:01 *683:01 *04:01:01:24N, 159:01N, 160:01, 328:01N, 454:01, 683:01, 1045:01N, 1121:01N
55	135 bp 165 bp 190 bp	430 bp	*155:01:01-155:01:02 *382:01N, 511:01 *115:01, 218:01N, 393:01, 1007:01
56	145 bp 180 bp 230 bp 290 bp	430 bp	*174:01, 658:01 *84:01, 148:01 *170:01, 239:01, 260:01, 487:01 *08:01 [?] , 60:01 [?] , 62:01 [?] , 73:01 [?] , 75:01 [?] , 102:01 [?] , 237:01, 722:01
57	255 bp	430 bp	*04:02:01:01-04:02:17, 05:01:03 [?] -05:01:06 [?] , 06:01:02 [?] -06:01:03 [?] , 08:01 [?] , 09:01:03 [?] -09:01:04 [?] , 10:01:02 [?] -10:01:03 [?] , 11:01:02 [?] -11:01:03 [?] , 20:01:02 [?] -20:01:03 [?] , 26:01:01 [?] , 26:01:03 [?] , 32:01 [?] , 37:01 [?] , 41:01:02 [?] , 44:01 [?] , 52:01 [?] , 56:01 [?] , 58:01 [?] , 59:01, 60:01 [?] , 61:01N [?] , 62:01 [?] , 65:01:01 [?] -65:01:02 [?] , 66:01 [?] , 67:01 [?] , 68:01 [?] , 70:01 [?] , 71:01:02 [?] , 73:01 [?] , 74:01 [?] , 75:01 [?] , 77:01, 79:01 [?] , 80:01, 82:01 [?] , 83:01 [?] , 84:01 [?] , 85:01:02 [?] , 86:01 [?] , 87:01 [?] , 89:01 [?] , 90:01:02 [?] , 92:01 [?] , 94:01 [?] , 95:01 [?] , 96:01 [?] , 97:01 [?] , 98:01 [?] , 99:01 [?] , 100:01 [?] , 101:01 [?] , 102:01 [?] , 103:01 [?] , 108:01 [?] , 110:01 [?] , 111:01 [?] , 112:01 [?] , 113:01 [?] , 114:01 [?] , 116:01 [?] , 117:01 [?] , 118:01 [?] , 119:01 [?] , 120:01N [?] , 122:01 [?] , 125:01 [?] , 126:01:01:01-126:01:01:03, 129:01 [?] , 138:01, 139:01 [?] , 140:01 [?] , 142:01 [?] , 143:01 [?] , 144:01 [?] , 145:01 [?] , 146:01 [?] , 147:01 [?] , 148:01 [?] , 149:01 [?] , 150:01 [?] , 151:01 [?] , 153:01 [?] , 154:01N [?] , 155:01:01 [?] -155:01:02 [?] , 156:01 [?] , 157:01 [?] , 158:01 [?] , 159:01N [?] , 160:01 [?] , 161:01N [?] , 163:01 [?] , 164:01 [?] , 166:01 [?] , 167:01 [?] , 169:01 [?] , 170:01 [?] , 171:01 [?] , 172:01 [?] , 173:01 [?] , 174:01 [?] , 175:01 [?] , 176:01 [?] , 177:01 [?] , 179:01 [?] , 180:01 [?] , 181:01 [?] , 182:01 [?] , 183:01 [?] , 184:01, 185:01 [?] , 186:01 [?] , 188:01 [?] , 189:01 [?] , 190:01, 192:01 [?] , 193:01 [?] , 194:01 [?] , 195:01 [?] , 196:01 [?] , 197:01 [?] , 198:01 [?] , 199:01 [?] , 200:01 [?] , 201:01 [?] , 202:01 [?] , 203:01:01 [?] -203:01:02 [?] , 204:01 [?] , 205:01 [?] , 206:01 [?] , 207:01 [?] , 209:01 [?] , 210:01 [?] , 211:01 [?] , 212:01 [?] , 213:01:01 [?] , 214:01 [?] , 215:01 [?] , 216:01N [?] , 217:01 [?] , 218:01N [?] , 219:01 [?] , 220:01 [?] ,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

221:01[?], 222:01[?], 223:01[?], 224:01[?], 226:01[?], 227:01:01[?], 228:01[?],
229:01[?], 230:01[?], 231:01[?], 232:01[?], 233:01[?], 234:01[?], 235:01[?],
236:01:02[?], 237:01[?], 238:01[?], 239:01[?], 240:01[?], 241:01[?], 242:01[?],
243:01[?], 245:01[?], 246:01[?], 247:01[?], 248:01[?], 249:01[?], 250:01[?],
251:01[?], 252:01[?], 253:01[?], 254:01[?], 255:01[?], 256:01[?], 257:01, 258:01[?],
259:01[?], 260:01, 261:01[?], 262:01[?], 263:01[?], 264:01[?], 265:01[?], 266:01[?],
267:01[?], 268:01[?], 269:01[?], 270:01[?], 271:01[?], 272:01[?], 273:01[?],
274:01[?], 275:01[?], 276:01[?], 277:01[?], 278:01[?], 279:01:01[?]-279:01:02[?],
280:01[?], 281:01[?], 282:01[?], 283:01[?], 284:01[?], 285:01[?], 286:01[?],
287:01[?], 288:01[?], 289:01[?], 290:01[?], 291:01[?], 292:01[?], 293:01[?],
294:01[?], 295:01[?], 297:01[?], 298:01[?], 299:01[?], 300:01[?], 301:01[?],
302:01[?], 303:01[?], 304:01[?], 305:01[?], 306:01[?], 307:01[?], 308:01[?],
309:01[?], 310:01[?], 311:01[?], 312:01[?], 313:01[?], 314:01[?], 315:01[?],
316:01[?], 317:01[?], 318:01[?], 319:01[?], 320:01[?], 321:01[?], 322:01[?],
323:01[?], 324:01[?], 325:01[?], 326:01[?], 327:01[?], 329:01[?], 330:01[?],
331:01[?], 332:01[?], 333:01[?], 334:01:01[?], 335:01[?], 336:01[?], 337:01[?],
338:01[?], 339:01[?], 340:01[?], 341:01[?], 342:01[?], 343:01[?], 344:01[?],
345:01[?], 346:01[?], 347:01[?], 348:01:01[?]-348:01:02[?], 349:01[?], 351:01,
354:01:01[?]-354:01:02[?], 355:01[?], 356:01[?], 357:01N[?], 358:01[?], 359:01[?],
360:01[?], 361:01[?], 362:01[?], 363:01[?], 364:01[?], 365:01[?], 366:01[?],
367:01[?], 368:01[?], 369:01[?], 370:01[?], 371:01[?], 372:01[?], 373:01[?],
374:01[?], 375:01[?], 376:01[?], 377:01[?], 378:01[?], 379:01[?], 380:01[?],
381:01[?], 382:01N[?], 383:01[?], 384:01[?], 385:01[?], 386:01[?], 387:01[?],
388:01[?], 389:01[?], 390:01[?], 391:01[?], 392:01[?], 393:01[?], 394:01[?],
395:01[?], 396:01[?], 397:01[?], 399:01[?], 400:01[?], 401:01N[?], 402:01[?],
403:01N[?], 404:01[?], 405:01[?], 406:01[?], 407:01[?], 408:01[?], 409:01[?],
410:01[?], 411:01[?], 412:01[?], 413:01[?], 416:01:01:01-416:01:02, 418:01[?],
419:01[?], 420:01[?], 421:01[?], 422:01[?], 423:01:01[?]-423:01:02[?], 424:01[?],
425:01[?], 426:01[?], 427:01[?], 428:01[?], 429:01[?], 430:01[?], 431:01[?],
432:01[?], 433:01[?], 434:01[?], 435:01[?], 436:01[?], 437:01[?], 438:01[?],
439:01[?], 441:01[?], 442:01[?], 443:01[?], 444:01[?], 445:01[?], 446:01[?],
447:01[?], 448:01[?], 449:01[?], 450:01N[?], 451:01[?], 452:01[?], 453:01[?],
454:01[?], 455:01N[?], 456:01[?], 457:01[?], 458:01:01[?], 460:01, 462:01,
465:01[?], 466:01[?], 467:01[?], 468:01[?], 469:01[?], 470:01[?], 471:01[?],
472:01[?], 473:01[?], 474:01[?], 475:01[?], 476:01[?], 477:01[?], 478:01[?],
479:01[?], 480:01[?], 481:01[?], 482:01[?], 483:01[?], 484:01[?], 485:01[?],
486:01[?], 487:01[?], 488:01[?], 489:01[?], 490:01[?], 491:01[?], 492:01[?],
493:01[?], 494:01[?], 495:01[?], 496:01[?], 497:01[?], 499:01[?], 500:01[?],
501:01[?], 502:01[?], 503:01[?], 504:01[?], 505:01[?], 506:01[?], 507:01N[?],
508:01[?], 509:01[?], 510:01[?], 511:01[?], 512:01[?], 513:01[?], 514:01[?],
515:01[?], 516:01[?], 517:01[?], 520:01[?], 521:01[?], 522:01[?], 523:01:01[?],
524:01[?], 525:01[?], 526:01[?], 527:01[?], 528:01[?], 529:01[?], 530:01[?],
531:01[?], 532:01[?], 537:01[?], 538:01[?], 539:01[?], 540:01[?], 541:01[?],
542:01[?], 543:01[?], 544:01[?], 545:01[?], 546:01[?], 548:01[?], 549:01[?],
550:01[?], 551:01N[?], 552:01[?], 553:01[?], 554:01[?], 555:01[?], 556:01[?],
557:01[?], 559:01[?], 560:01[?], 561:01[?], 562:01[?], 563:01[?], 564:01[?],
565:01[?], 566:01[?], 567:01[?], 568:01[?], 569:01[?], 570:01N[?], 573:01[?],
574:01[?], 575:01[?], 576:01[?], 577:01[?], 578:01[?], 579:01[?], 580:01[?],
581:01[?], 582:01[?], 583:01[?], 586:01[?], 587:01[?], 588:01[?], 589:01[?],
590:01[?], 591:01[?], 592:01[?], 593:01[?], 594:01[?], 595:01[?], 596:01[?],
597:01[?], 598:01N[?], 599:01[?], 600:01[?], 601:01[?], 602:01[?], 603:01[?],
604:01[?], 605:01[?], 606:01[?], 607:01[?], 608:01[?], 609:01[?], 610:01[?],



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			611:01 [?] , 612:01 [?] , 614:01 [?] , 619:01 [?] , 620:01 [?] , 622:01 [?] , 623:01 [?] , 624:01 [?] , 625:01 [?] , 626:01 [?] , 627:01 [?] , 628:01 [?] , 629:01 [?] , 630:01 [?] , 631:01 [?] , 632:01 [?] , 633:01 [?] , 634:01 [?] , 635:01 [?] , 636:01 [?] , 637:01 [?] , 638:01 [?] , 639:01 [?] , 641:01 [?] , 642:01 [?] , 643:01 [?] , 644:01 [?] , 645:01 [?] , 646:01 [?] , 647:01, 655:01 [?] , 658:01 [?] , 659:01 [?] , 660:01 [?] , 662:01 [?] , 663:01 [?] , 664:01 [?] , 674:01, 680:01 [?] , 684:01 [?] , 685:01 [?] , 686:01 [?] , 688:01 [?] , 689:01 [?] , 690:01 [?] , 691:01N [?] , 692:01 [?] , 693:01N [?] , 694:01 [?] , 697:01Q [?] , 701:01, 712:01N [?] , 716:01 [?] , 717:01 [?] , 718:01 [?] , 720:01 [?] , 721:01, 724:01N, 725:01, 726:01, 730:01, 731:01, 734:01, 735:01, 736:01, 741:01, 749:01 [?] , 753:01 [?] , 758:01 [?] , 759:01, 760:01 [?] , 775:01, 776:01, 777:01N, 780:01 [?] , 785:01, 788:01 [?] , 789:01 [?] , 790:01, 796:01 [?] , 805:01, 809:01, 817:01, 818:01, 823:01, 837:01 [?] , 840:01 [?] , 841:01, 854:01, 858:01, 868:01N [?] , 869:01N [?] , 870:01N [?] , 872:01N [?] , 873:01N [?] , 881:01, 885:01, 887:01, 891:01 [?] , 893:01 [?] , 895:01 [?] , 903:01 [?] , 906:01, 911:01N [?] , 925:01N, 927:01, 933:01, 936:01Q, 937:01, 954:01, 958:01, 959:01N, 960:01N [?] , 961:01, 968:01, 981:01, 985:01N, 995:01N, 996:01, 1005:01, 1006:01 [?] , 1020:01, 1025:01, 1035:01, 1041:01N [?] , 1075:01, 1078:01 [?] , 1079:01N [?] , 1080:01 [?] , 1083:01, 1084:01N, 1085:01, 1092:01Q, 1110:01, 1124:01, 1138:01 [?] , 1141:01 [?] , 1144:01, 1150:01 [?] , 1153:01, 1154:01N [?] , 1158:01 [?] , 1159:01 [?] , 1163:01 [?] , 1166:01 [?] , 1174:01 [?] , 1176:01 [?] , 1178:01 [?] , 1182:01 [?] , 1194:01, 1195:01 [?]
58^{4,5}	95 bp	430 bp	*02:01:24 [?] , 02:01:26 [?] -02:01:33 [?] , 02:01:35 [?] -02:01:36 [?] , 02:01:38 [?] - 02:01:41 [?] , 02:01:43 [?] -02:01:57 [?] , 02:02:02 [?] , 02:02:04 [?] , 03:01:09 [?] , 03:01:13 [?] -03:01:14 [?] , 04:01:41 [?] , 04:01:44 [?] -04:01:45 [?] , 04:01:47 [?] , 04:02:15 [?] -04:02:17 [?] , 05:01:08 [?] -05:01:11 [?] , 05:01:13 [?] , 06:01:05 [?] , 08:01 [?] , 09:01:03 [?] -09:01:05 [?] , 10:01:02 [?] -10:01:05 [?] , 11:01:02 [?] - 11:01:04 [?] , 14:01:02 [?] -14:01:04 [?] , 14:01:06 [?] -14:01:12 [?] , 17:01:02 [?] - 17:01:04 [?] , 18:01:02 [?] -18:01:03 [?] , 20:01:02 [?] -20:01:03 [?] , 23:01:02 [?] , 26:01:01 [?] , 26:01:03 [?] , 32:01 [?] , 34:01:02 [?] , 37:01 [?] , 40:01:02 [?] , 41:01:02 [?] , 44:01 [?] , 46:01:02 [?] , 52:01 [?] , 55:01:02 [?] -55:01:03 [?] , 56:01 [?] , 58:01 [?] , 60:01 [?] , 61:01N [?] , 62:01 [?] , 65:01:01 [?] -65:01:02 [?] , 66:01 [?] , 67:01 [?] , 68:01 [?] , 70:01 [?] , 71:01:02 [?] , 73:01 [?] , 74:01 [?] , 75:01 [?] , 79:01 [?] , 81:01:02 [?] , 82:01 [?] , 83:01 [?] , 84:01 [?] , 85:01:02 [?] , 86:01 [?] , 87:01 [?] , 89:01 [?] , 90:01:02 [?] , 91:01:02 [?] , 92:01 [?] , 94:01 [?] , 95:01 [?] , 96:01 [?] , 97:01 [?] , 98:01 [?] , 99:01 [?] , 100:01 [?] , 101:01 [?] , 102:01 [?] , 103:01 [?] , 104:01:03 [?] -104:01:04 [?] , 107:01, 108:01 [?] , 110:01 [?] , 111:01 [?] , 112:01 [?] , 113:01 [?] , 114:01 [?] , 116:01 [?] , 117:01 [?] , 118:01 [?] , 119:01 [?] , 120:01N [?] , 121:01 [?] , 122:01 [?] , 123:01 [?] , 125:01 [?] , 129:01 [?] , 131:01:02 [?] -131:01:03 [?] , 132:01 [?] , 134:01 [?] , 137:01 [?] , 139:01 [?] , 140:01 [?] , 142:01 [?] , 143:01 [?] , 144:01 [?] , 145:01 [?] , 146:01 [?] , 147:01 [?] , 148:01 [?] , 149:01 [?] , 150:01 [?] , 151:01 [?] , 152:01 [?] , 153:01 [?] , 154:01N [?] , 155:01:01 [?] -155:01:02 [?] , 156:01 [?] , 157:01 [?] , 158:01 [?] , 159:01N [?] , 160:01 [?] , 161:01N [?] , 162:01:01 [?] , 163:01 [?] , 164:01 [?] , 165:01 [?] , 166:01 [?] , 167:01 [?] , 168:01 [?] , 169:01 [?] , 170:01 [?] , 171:01 [?] , 172:01 [?] , 173:01 [?] , 174:01 [?] , 175:01 [?] , 176:01 [?] , 177:01 [?] , 179:01 [?] , 180:01 [?] , 181:01 [?] , 182:01 [?] , 183:01 [?] , 185:01 [?] , 186:01 [?] , 188:01 [?] , 189:01 [?] , 192:01 [?] , 193:01 [?] , 194:01 [?] , 195:01 [?] , 196:01 [?] , 197:01 [?] , 198:01 [?] , 199:01 [?] , 200:01 [?] , 201:01 [?] , 202:01 [?] , 203:01:01 [?] -203:01:02 [?] , 204:01 [?] , 205:01 [?] , 206:01 [?] , 207:01 [?] , 209:01 [?] , 210:01 [?] , 211:01 [?] , 212:01 [?] , 213:01:01 [?] -213:01:02 [?] , 214:01 [?] , 215:01 [?] , 216:01N [?] , 217:01 [?] , 218:01N [?] , 219:01 [?] , 220:01 [?] , 221:01 [?] , 222:01 [?] , 223:01 [?] , 224:01 [?] , 225:01 [?] , 226:01 [?] , 227:01:01 [?] - 227:01:03 [?] , 228:01 [?] , 229:01 [?] , 230:01 [?] , 231:01 [?] , 232:01 [?] , 233:01 [?] , 234:01 [?] , 235:01 [?] , 236:01:02 [?] , 237:01 [?] , 238:01 [?] , 239:01 [?] , 240:01 [?] ,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

241:01[?], 242:01[?], 243:01[?], 245:01[?], 246:01[?], 247:01[?], 248:01[?],
249:01[?], 250:01[?], 251:01[?], 252:01[?], 253:01[?], 254:01[?], 255:01[?],
256:01[?], 257:01[?], 258:01[?], 259:01[?], 261:01[?], 262:01[?], 263:01[?],
264:01[?], 265:01[?], 266:01[?], 267:01[?], 268:01[?], 269:01[?], 270:01[?],
271:01[?], 272:01[?], 273:01[?], 274:01[?], 275:01[?], 276:01[?], 277:01[?],
278:01[?], 279:01:01[?]-279:01:02[?], 280:01[?], 281:01[?], 282:01[?], 283:01[?],
284:01[?], 285:01[?], 286:01[?], 287:01[?], 288:01[?], 289:01[?], 290:01[?],
291:01[?], 292:01[?], 293:01[?], 294:01[?], 295:01[?], 296:01, 297:01[?], 298:01[?],
299:01[?], 300:01[?], 301:01[?], 302:01[?], 303:01[?], 304:01[?], 305:01[?],
306:01[?], 307:01[?], 308:01[?], 309:01[?], 310:01[?], 311:01[?], 312:01[?],
313:01[?], 314:01[?], 315:01[?], 316:01[?], 317:01[?], 318:01[?], 319:01[?],
320:01[?], 321:01[?], 322:01[?], 323:01[?], 324:01[?], 325:01[?], 326:01[?],
327:01[?], 329:01[?], 330:01[?], 331:01[?], 332:01[?], 333:01[?], 334:01:01[?]-
334:01:02[?], 335:01[?], 336:01[?], 337:01[?], 338:01[?], 339:01[?], 340:01[?],
341:01[?], 342:01[?], 343:01[?], 344:01[?], 345:01[?], 346:01[?], 347:01[?],
348:01:01[?]-348:01:02[?], 349:01[?], 350:01[?], 351:01[?], 352:01:01[?],
353:01[?], 354:01:01[?]-354:01:02[?], 355:01[?], 356:01[?], 357:01N[?], 358:01[?],
359:01[?], 360:01[?], 361:01[?], 362:01[?], 363:01[?], 364:01[?], 365:01[?],
366:01[?], 367:01[?], 368:01[?], 369:01[?], 370:01[?], 371:01[?], 372:01[?],
373:01[?], 374:01[?], 375:01[?], 376:01[?], 377:01[?], 378:01[?], 379:01[?],
380:01[?], 381:01[?], 382:01N[?], 383:01[?], 384:01[?], 385:01[?], 386:01[?],
387:01[?], 388:01[?], 389:01[?], 390:01[?], 391:01[?], 392:01[?], 393:01[?],
394:01[?], 395:01[?], 396:01[?], 397:01[?], 399:01[?], 400:01[?], 401:01N[?],
402:01[?], 403:01N[?], 404:01[?], 405:01[?], 406:01[?], 407:01[?], 408:01[?],
409:01[?], 410:01[?], 411:01[?], 412:01[?], 413:01[?], 417:01:01[?]-417:01:02[?],
418:01[?], 419:01[?], 420:01[?], 421:01[?], 422:01[?], 423:01:01[?]-423:01:02[?],
424:01[?], 425:01[?], 426:01[?], 427:01[?], 428:01[?], 429:01[?], 430:01[?],
431:01[?], 432:01[?], 433:01[?], 434:01[?], 435:01[?], 436:01[?], 437:01[?],
438:01[?], 439:01[?], 440:01[?], 441:01[?], 442:01[?], 443:01[?], 444:01[?],
445:01[?], 446:01[?], 447:01[?], 448:01[?], 449:01[?], 450:01N[?], 451:01[?],
452:01[?], 453:01[?], 454:01[?], 455:01N[?], 456:01[?], 457:01[?], 458:01:01[?]-
458:01:02[?], 459:01[?], 461:01[?], 462:01[?], 465:01[?], 466:01[?], 467:01[?],
468:01[?], 469:01[?], 470:01[?], 471:01[?], 472:01[?], 473:01[?], 474:01[?],
475:01[?], 476:01[?], 477:01[?], 478:01[?], 479:01[?], 480:01[?], 481:01[?],
482:01[?], 483:01[?], 484:01[?], 485:01[?], 486:01[?], 487:01[?], 488:01[?],
489:01[?], 490:01[?], 491:01[?], 492:01[?], 493:01[?], 494:01[?], 495:01[?],
496:01[?], 497:01[?], 499:01[?], 500:01[?], 501:01[?], 502:01[?], 503:01[?],
504:01[?], 505:01[?], 506:01[?], 507:01N[?], 508:01[?], 509:01[?], 510:01[?],
511:01[?], 512:01[?], 513:01[?], 514:01[?], 515:01[?], 516:01[?], 517:01[?],
520:01[?], 521:01[?], 522:01[?], 523:01:01[?]-523:01:02[?], 524:01[?], 525:01[?],
526:01[?], 527:01[?], 528:01[?], 529:01[?], 530:01[?], 531:01[?], 532:01[?],
533:01[?], 534:01[?], 536:01[?], 537:01[?], 538:01[?], 539:01[?], 540:01[?],
541:01[?], 542:01[?], 543:01[?], 544:01[?], 545:01[?], 546:01[?], 548:01[?],
549:01[?], 550:01[?], 551:01N[?], 552:01[?], 553:01[?], 554:01[?], 555:01[?],
556:01[?], 557:01[?], 558:01[?], 559:01[?], 560:01[?], 561:01[?], 562:01[?],
563:01[?], 564:01[?], 565:01[?], 566:01[?], 567:01[?], 568:01[?], 569:01[?],
570:01N[?], 571:01[?], 573:01[?], 574:01[?], 575:01[?], 576:01[?], 577:01[?],
578:01[?], 579:01[?], 580:01[?], 581:01[?], 582:01[?], 583:01[?], 586:01[?],
587:01[?], 588:01[?], 589:01[?], 590:01[?], 591:01[?], 592:01[?], 593:01[?],
594:01[?], 595:01[?], 596:01[?], 597:01[?], 598:01N[?], 599:01[?], 600:01[?],
601:01[?], 602:01[?], 603:01[?], 604:01[?], 605:01[?], 606:01[?], 607:01[?],
608:01[?], 609:01[?], 610:01[?], 611:01[?], 612:01[?], 613:01[?], 614:01[?],



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

615:01[?], 616:01[?], 617:01[?], 618:01[?], 619:01[?], 620:01[?], 621:01[?],
622:01[?], 623:01[?], 624:01[?], 625:01[?], 626:01[?], 627:01[?], 628:01[?],
629:01[?], 630:01[?], 631:01[?], 632:01[?], 633:01[?], 634:01[?], 635:01[?],
636:01[?], 637:01[?], 638:01[?], 639:01[?], 640:01[?], 641:01[?], 642:01[?],
643:01[?], 644:01[?], 645:01[?], 646:01[?], 647:01[?], 654:01[?], 655:01[?],
656:01[?], 657:01N[?], 658:01[?], 659:01[?], 660:01[?], 661:01N[?], 662:01[?],
663:01[?], 664:01[?], 665:01[?], 679:01[?], 680:01[?], 681:01[?], 682:01[?], 683:01[?],
684:01[?], 685:01[?], 686:01[?], 687:01[?], 688:01[?], 689:01[?], 690:01[?],
691:01N[?], 692:01[?], 693:01N[?], 694:01[?], 695:01[?], 696:01N[?], 697:01Q[?],
698:01[?], 702:01[?], 703:01[?], 707:01[?], 708:01[?], 709:01[?], 710:01[?],
711:01[?], 712:01N[?], 713:01[?], 714:01[?], 715:01[?], 716:01[?], 717:01[?],
718:01[?], 719:01[?], 720:01[?], 721:01[?], 722:01[?], 723:01[?], 725:01[?],
727:01[?], 729:01[?], 730:01[?], 731:01[?], 732:01N[?], 739:01[?], 740:01[?],
741:01[?], 742:01[?], 743:01N[?], 744:01[?], 745:01[?], 746:01[?], 747:01[?],
748:01N[?], 749:01[?], 750:01[?], 751:01[?], 752:01[?], 753:01[?], 754:01N[?],
755:01[?], 756:01N[?], 757:01[?], 758:01[?], 759:01[?], 760:01[?], 761:01[?],
763:01[?], 765:01[?], 766:01[?], 767:01[?], 768:01[?], 769:01[?], 770:01[?],
771:01[?], 772:01[?], 773:01[?], 774:01[?], 775:01[?], 776:01[?], 777:01N[?],
778:01[?], 779:01[?], 780:01[?], 781:01[?], 782:01[?], 783:01[?], 785:01[?],
786:01:02N[?], 787:01[?], 788:01[?], 789:01[?], 790:01[?], 791:01[?], 792:01N[?],
794:01N[?], 795:01[?], 796:01[?], 797:01[?], 798:01[?], 799:01[?], 800:01N[?],
801:01[?], 802:01[?], 803:01[?], 805:01[?], 807:01[?], 808:01[?], 810:01[?],
811:01[?], 812:01[?], 813:01[?], 814:01[?], 815:01[?], 816:01[?], 817:01[?],
818:01[?], 820:01[?], 821:01N[?], 822:01[?], 824:01:01[?]-824:01:02[?], 825:01[?],
826:01[?], 827:01[?], 828:01[?], 829:01[?], 830:01[?], 831:01N[?], 833:01[?],
836:01[?], 837:01[?], 838:01N[?], 839:01[?], 840:01[?], 841:01[?], 842:01[?],
843:01[?], 844:01N[?], 845:01[?], 846:01[?], 847:01[?], 848:01[?], 850:01[?],
851:01[?], 852:01[?], 853:01[?], 854:01[?], 855:01[?], 856:01[?], 857:01[?],
858:01[?], 859:01[?], 860:01[?], 861:01[?], 862:01N[?], 863:01[?], 864:01[?],
865:01N[?], 866:01N[?], 867:01N[?], 868:01N[?], 869:01N[?], 870:01N[?],
871:01N[?], 872:01N[?], 873:01N[?], 874:01N[?], 875:01N[?], 876:01N[?],
877:01N[?], 883:01[?], 884:01[?], 885:01[?], 886:01[?], 887:01[?], 888:01[?],
889:01[?], 890:01[?], 891:01[?], 892:01[?], 893:01[?], 894:01N[?], 895:01[?],
896:01[?], 897:01[?], 898:01[?], 899:01[?], 900:01[?], 901:01[?], 902:01[?],
903:01[?], 904:01[?], 906:01[?], 907:01[?], 908:01[?], 910:01[?], 911:01N[?],
912:01[?], 913:01[?], 914:01[?], 915:01[?], 916:01[?], 917:01N[?], 918:01[?],
919:01N[?], 920:01[?], 921:01[?], 922:01[?], 923:01[?], 924:01[?], 925:01N[?],
926:01[?], 927:01[?], 928:01[?], 929:01[?], 930:01[?], 931:01[?], 932:01[?],
933:01[?], 934:01Q[?], 935:01Q[?], 936:01Q[?], 937:01[?], 938:01[?], 939:01[?],
940:01[?], 941:01N[?], 942:01[?], 943:01[?], 944:01[?], 945:01[?], 958:01[?],
959:01N[?], 960:01N[?], 961:01[?], 962:01[?], 963:01[?], 964:01[?], 969:01[?],
970:01[?], 971:01[?], 972:01[?], 973:01[?], 974:01N[?], 975:01[?], 976:01[?],
977:01[?], 978:01[?], 979:01[?], 980:01[?], 981:01[?], 982:01[?], 983:01[?],
984:01N[?], 985:01N[?], 986:01N[?], 987:01[?], 988:01[?], 989:01[?], 990:01[?],
991:01[?], 992:01[?], 993:01[?], 994:01[?], 995:01N[?], 996:01[?], 997:01[?],
998:01[?], 999:01[?], 1000:01[?], 1001:01[?], 1002:01[?], 1003:01[?], 1004:01[?],
1005:01[?], 1006:01[?], 1007:01[?], 1008:01[?], 1009:01[?], 1010:01[?],
1011:01[?], 1012:01[?], 1013:01[?], 1014:01[?], 1015:01[?], 1016:01[?],
1017:01[?], 1018:01[?], 1019:01[?], 1020:01[?], 1021:01[?], 1022:01[?],
1023:01[?], 1024:01[?], 1025:01[?], 1026:01[?], 1027:01[?], 1028:01[?],
1029:01N[?], 1030:01[?], 1031:01[?], 1032:01[?], 1033:01[?], 1034:01[?],
1035:01[?], 1036:01[?], 1038:01Q[?], 1039:01[?], 1040:01[?], 1041:01N[?],



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			1042:01 [?] , 1043:01 [?] , 1044:01N [?] , 1045:01N [?] , 1046:01 [?] , 1047:01 [?] , 1048:01 [?] , 1049:01 [?] , 1050:01 [?] , 1051:01 [?] , 1052:01 [?] , 1053:01 [?] , 1054:01 [?] , 1055:01 [?] , 1056:01 [?] , 1057:01 [?] , 1058:01 [?] , 1059:01 [?] , 1060:01 [?] , 1061:01 [?] , 1062:01 [?] , 1063:01 [?] , 1064:01 [?] , 1065:01 [?] , 1066:01 [?] , 1067:01 [?] , 1068:01 [?] , 1069:01 [?] , 1070:01N [?] , 1071:01 [?] , 1073:01 [?] , 1074:01 [?] , 1075:01 [?] , 1076:01 [?] , 1078:01 [?] , 1079:01N [?] , 1080:01 [?] , 1081:01 [?] , 1084:01N [?] , 1085:01 [?] , 1086:01 [?] , 1087:01 [?] , 1088:01 [?] , 1089:01 [?] , 1090:01 [?] , 1091:01 [?] , 1092:01Q [?] , 1093:01 [?] , 1094:01 [?] , 1095:01 [?] , 1096:01 [?] , 1097:01 [?] , 1098:01N [?] , 1101:01 [?] , 1102:01 [?] , 1103:01 [?] , 1104:01 [?] , 1105:01 [?] , 1106:01 [?] , 1107:01 [?] , 1108:01 [?] , 1110:01 [?] , 1111:01 [?] , 1112:01N [?] , 1113:01 [?] , 1114:01 [?] , 1115:01 [?] , 1116:01 [?] , 1117:01 [?] , 1118:01 [?] , 1119:01 [?] , 1120:01 [?] , 1121:01N [?] , 1122:01 [?] , 1123:01 [?] , 1124:01 [?] , 1135:01N [?] , 1136:01 [?] , 1137:01 [?] , 1138:01 [?] , 1140:01 [?] , 1141:01 [?] , 1142:01 [?] , 1143:01 [?] , 1146:01 [?] , 1148:01Q [?] , 1149:01 [?] , 1150:01 [?] , 1151:01 [?] , 1152:01 [?] , 1154:01N [?] , 1155:01 [?] , 1156:01 [?] , 1157:01 [?] , 1158:01 [?] , 1159:01 [?] , 1160:01 [?] , 1161:01 [?] , 1162:01 [?] , 1163:01 [?] , 1164:01 [?] , 1165:01 [?] , 1166:01 [?] , 1167:01 [?] , 1168:01 [?] , 1169:01N [?] , 1170:01 [?] , 1171:01 [?] , 1172:01 [?] , 1173:01 [?] , 1174:01 [?] , 1175:01 [?] , 1176:01 [?] , 1177:01 [?] , 1178:01 [?] , 1179:01 [?] , 1180:01 [?] , 1181:01 [?] , 1182:01 [?] , 1183:01 [?] , 1184:01 [?] , 1185:01 [?] , 1186:01 [?] , 1187:01 [?] , 1189:01 [?] , 1190:01 [?] , 1191:01N [?] , 1192:01 [?] , 1193:01N [?] , 1194:01 [?] , 1195:01 [?]
59⁴	80 bp	430 bp	*02:01:24 [?] , 02:01:26 [?] -02:01:33 [?] , 02:01:35 [?] -02:01:36 [?] , 02:01:38 [?] - 02:01:41 [?] , 02:01:43 [?] -02:01:49 [?] , 02:01:51 [?] , 02:02:02 [?] , 02:02:04 [?] , 05:01:01:01-05:01:14, 06:01:05 [?] , 08:01 [?] , 09:01:03 [?] -09:01:05 [?] , 10:01:02 [?] -10:01:03 [?] , 11:01:02 [?] -11:01:04 [?] , 14:01:02 [?] -14:01:03 [?] , 14:01:06 [?] -14:01:09 [?] , 15:01:02 [?] , 17:01:02 [?] -17:01:03 [?] , 18:01:02 [?] - 18:01:03 [?] , 19:01:01:01-19:01:01:03, 20:01:02 [?] -20:01:03 [?] , 22:01:01:01-22:01:01:02, 26:01:01 [?] , 26:01:03 [?] , 29:01, 32:01 [?] , 34:01:01:01-34:01:01:03, 34:01:02 [?] , 37:01 [?] , 38:01, 40:01:02 [?] , 41:01:02 [?] , 44:01 [?] , 46:01:02 [?] , 52:01 [?] , 55:01:02 [?] -55:01:03 [?] , 56:01 [?] , 58:01 [?] , 60:01 [?] , 61:01N [?] , 62:01 [?] , 63:01, 65:01:01 [?] -65:01:02 [?] , 66:01 [?] , 67:01 [?] , 68:01 [?] , 70:01 [?] , 71:01:02 [?] , 73:01 [?] , 74:01 [?] , 75:01 [?] , 78:01, 79:01 [?] , 81:01:02 [?] , 82:01 [?] , 83:01 [?] , 84:01 [?] , 85:01:02 [?] , 86:01 [?] , 87:01 [?] , 89:01 [?] , 90:01:02 [?] , 91:01:02 [?] , 92:01 [?] , 94:01 [?] , 95:01 [?] , 96:01 [?] , 97:01 [?] , 98:01 [?] , 99:01 [?] , 100:01 [?] , 101:01 [?] , 102:01 [?] , 103:01 [?] , 104:01:01:01- 104:01:04, 108:01 [?] , 110:01 [?] , 111:01 [?] , 112:01 [?] , 113:01 [?] , 114:01 [?] , 116:01 [?] , 117:01 [?] , 118:01 [?] , 119:01 [?] , 120:01N [?] , 121:01 [?] , 122:01 [?] , 123:01 [?] , 125:01 [?] , 129:01 [?] , 131:01:02 [?] -131:01:03 [?] , 132:01 [?] , 134:01 [?] , 137:01 [?] , 139:01 [?] , 140:01 [?] , 142:01 [?] , 143:01 [?] , 144:01 [?] , 145:01 [?] , 146:01 [?] , 147:01 [?] , 148:01 [?] , 149:01 [?] , 150:01 [?] , 151:01 [?] , 152:01 [?] , 153:01 [?] , 154:01N [?] , 155:01:01 [?] -155:01:02 [?] , 156:01 [?] , 157:01 [?] , 158:01 [?] , 159:01N [?] , 160:01 [?] , 161:01N [?] , 163:01 [?] , 164:01 [?] , 166:01 [?] , 167:01 [?] , 168:01 [?] , 169:01 [?] , 170:01 [?] , 171:01 [?] , 172:01 [?] , 173:01 [?] , 174:01 [?] , 175:01 [?] , 176:01 [?] , 177:01 [?] , 179:01 [?] , 180:01 [?] , 181:01 [?] , 182:01 [?] , 183:01 [?] , 185:01 [?] , 186:01 [?] , 188:01 [?] , 189:01 [?] , 192:01 [?] , 193:01 [?] , 194:01 [?] , 195:01 [?] , 196:01 [?] , 197:01 [?] , 198:01 [?] , 199:01 [?] , 200:01 [?] , 201:01 [?] , 202:01 [?] , 203:01:01 [?] -203:01:02 [?] , 204:01 [?] , 205:01 [?] , 206:01 [?] , 207:01 [?] , 209:01 [?] , 210:01 [?] , 211:01 [?] , 212:01 [?] , 213:01:01 [?] -213:01:02 [?] , 214:01 [?] , 215:01 [?] , 216:01N [?] , 217:01 [?] , 218:01N [?] , 219:01 [?] , 220:01 [?] , 221:01 [?] , 222:01 [?] , 223:01 [?] , 224:01 [?] , 225:01 [?] , 226:01 [?] , 227:01:01 [?] - 227:01:02 [?] , 228:01 [?] , 229:01 [?] , 230:01 [?] , 231:01 [?] , 232:01 [?] , 233:01 [?] ,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

234:01[?], 235:01[?], 236:01:02[?], 237:01[?], 238:01[?], 239:01[?], 240:01[?],
241:01[?], 242:01[?], 243:01[?], 245:01[?], 246:01[?], 247:01[?], 248:01[?],
249:01[?], 250:01[?], 251:01[?], 252:01[?], 253:01[?], 254:01[?], 255:01[?],
256:01[?], 258:01[?], 259:01[?], 261:01[?], 262:01[?], 263:01[?], 264:01[?],
265:01[?], 266:01[?], 267:01[?], 268:01[?], 269:01[?], 270:01[?], 271:01[?],
272:01[?], 273:01[?], 274:01[?], 275:01[?], 276:01[?], 277:01[?], 278:01[?],
279:01:01[?]-279:01:02[?], 280:01[?], 281:01[?], 282:01[?], 283:01[?], 284:01[?],
285:01[?], 286:01[?], 287:01[?], 288:01[?], 289:01[?], 290:01[?], 291:01[?],
292:01[?], 293:01[?], 294:01[?], 295:01[?], 297:01[?], 298:01[?], 299:01[?],
300:01[?], 301:01[?], 302:01[?], 303:01[?], 304:01[?], 305:01[?], 306:01[?],
307:01[?], 308:01[?], 309:01[?], 310:01[?], 311:01[?], 312:01[?], 313:01[?],
314:01[?], 315:01[?], 316:01[?], 317:01[?], 318:01[?], 319:01[?], 320:01[?],
321:01[?], 322:01[?], 323:01[?], 324:01[?], 325:01[?], 326:01[?], 327:01[?],
329:01[?], 330:01[?], 331:01[?], 332:01[?], 333:01[?], 334:01:01[?]-334:01:02[?],
335:01[?], 336:01[?], 337:01[?], 338:01[?], 339:01[?], 340:01[?], 341:01[?],
342:01[?], 343:01[?], 344:01[?], 345:01[?], 346:01[?], 347:01[?], 348:01:01[?]-
348:01:02[?], 349:01[?], 354:01:01[?]-354:01:02[?], 355:01[?], 356:01[?],
357:01N[?], 358:01[?], 359:01[?], 360:01[?], 361:01[?], 362:01[?], 363:01[?],
364:01[?], 365:01[?], 366:01[?], 367:01[?], 368:01[?], 369:01[?], 370:01[?],
371:01[?], 372:01[?], 373:01[?], 374:01[?], 375:01[?], 376:01[?], 377:01[?],
378:01[?], 379:01[?], 380:01[?], 381:01[?], 382:01N[?], 383:01[?], 384:01[?],
385:01[?], 386:01[?], 387:01[?], 388:01[?], 389:01[?], 390:01[?], 391:01[?],
392:01[?], 393:01[?], 394:01[?], 395:01[?], 396:01[?], 397:01[?], 399:01[?],
400:01[?], 401:01N[?], 402:01[?], 403:01N[?], 404:01[?], 405:01[?], 406:01[?],
407:01[?], 408:01[?], 409:01[?], 410:01[?], 411:01[?], 412:01[?], 413:01[?],
414:01:01:01-414:01:01:02, 418:01[?], 419:01[?], 420:01[?], 421:01[?],
422:01[?], 423:01:01[?]-423:01:02[?], 424:01[?], 425:01[?], 426:01[?], 427:01[?],
428:01[?], 429:01[?], 430:01[?], 431:01[?], 432:01[?], 433:01[?], 434:01[?],
435:01[?], 436:01[?], 437:01[?], 438:01[?], 439:01[?], 441:01[?], 442:01[?],
443:01[?], 444:01[?], 445:01[?], 446:01[?], 447:01[?], 448:01[?], 449:01[?],
450:01N[?], 451:01[?], 452:01[?], 453:01[?], 454:01[?], 455:01N[?], 456:01[?],
457:01[?], 458:01:01[?]-458:01:02[?], 465:01[?], 466:01[?], 467:01[?], 468:01[?],
469:01[?], 470:01[?], 471:01[?], 472:01[?], 473:01[?], 474:01[?], 475:01[?],
476:01[?], 477:01[?], 478:01[?], 479:01[?], 480:01[?], 481:01[?], 482:01[?],
483:01[?], 484:01[?], 485:01[?], 486:01[?], 487:01[?], 488:01[?], 489:01[?],
490:01[?], 491:01[?], 492:01[?], 493:01[?], 494:01[?], 495:01[?], 496:01[?],
497:01[?], 499:01[?], 500:01[?], 501:01[?], 502:01[?], 503:01[?], 504:01[?],
505:01[?], 506:01[?], 507:01N[?], 508:01[?], 509:01[?], 510:01[?], 511:01[?],
512:01[?], 513:01[?], 514:01[?], 515:01[?], 516:01[?], 517:01[?], 519:01, 520:01[?],
521:01[?], 522:01[?], 523:01:01[?]-523:01:02[?], 524:01[?], 525:01[?], 526:01[?],
527:01[?], 528:01[?], 529:01[?], 530:01[?], 531:01[?], 532:01[?], 533:01[?],
534:01[?], 536:01[?], 537:01[?], 538:01[?], 539:01[?], 540:01[?], 541:01[?],
542:01[?], 543:01[?], 544:01[?], 545:01[?], 546:01[?], 547:01, 548:01[?], 549:01[?],
550:01[?], 551:01N[?], 552:01[?], 553:01[?], 554:01[?], 555:01[?], 556:01[?],
557:01[?], 559:01[?], 560:01[?], 561:01[?], 562:01[?], 563:01[?], 564:01[?],
565:01[?], 566:01[?], 567:01[?], 568:01[?], 569:01[?], 570:01N[?], 573:01[?],
574:01[?], 575:01[?], 576:01[?], 577:01[?], 578:01[?], 579:01[?], 580:01[?],
581:01[?], 582:01[?], 583:01[?], 586:01[?], 587:01[?], 588:01[?], 589:01[?],
590:01[?], 591:01[?], 592:01[?], 593:01[?], 594:01[?], 595:01[?], 596:01[?],
597:01[?], 598:01N[?], 599:01[?], 600:01[?], 601:01[?], 602:01[?], 603:01[?],
604:01[?], 605:01[?], 606:01[?], 607:01[?], 608:01[?], 609:01[?], 610:01[?],
611:01[?], 612:01[?], 613:01[?], 614:01[?], 616:01[?], 617:01[?], 618:01[?],



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

619:01[?], 620:01[?], 621:01[?], 622:01[?], 623:01[?], 624:01[?], 625:01[?],
626:01[?], 627:01[?], 628:01[?], 629:01[?], 630:01[?], 631:01[?], 632:01[?],
633:01[?], 634:01[?], 635:01[?], 636:01[?], 637:01[?], 638:01[?], 639:01[?],
641:01[?], 642:01[?], 643:01[?], 644:01[?], 645:01[?], 646:01[?], 648:01:01:01-
648:01:01:02, 650:01, 652:01, 655:01[?], 656:01[?], 657:01N[?], 658:01[?],
659:01[?], 660:01[?], 661:01N[?], 662:01[?], 663:01[?], 664:01[?], 668:01:01:01-
668:01:01:02, 679:01[?], 680:01[?], 681:01[?], 682:01[?], 683:01[?], 684:01[?],
685:01[?], 686:01[?], 688:01[?], 689:01[?], 690:01[?], 691:01N[?], 692:01[?],
693:01N[?], 694:01[?], 695:01[?], 696:01N[?], 697:01Q[?], 698:01[?], 702:01,
703:01[?], 707:01[?], 708:01[?], 709:01[?], 710:01[?], 711:01[?], 712:01N[?],
713:01[?], 714:01[?], 715:01[?], 716:01[?], 717:01[?], 718:01[?], 720:01[?],
721:01[?], 722:01[?], 723:01[?], 725:01[?], 727:01[?], 729:01[?], 730:01[?],
731:01[?], 732:01N[?], 739:01[?], 740:01[?], 741:01[?], 742:01[?], 743:01N[?],
744:01[?], 745:01[?], 746:01[?], 747:01[?], 748:01N[?], 749:01[?], 750:01[?],
751:01[?], 752:01[?], 753:01[?], 754:01N[?], 755:01[?], 756:01N[?], 757:01[?],
758:01[?], 759:01[?], 760:01[?], 761:01[?], 764:01, 765:01[?], 766:01[?], 767:01[?],
768:01[?], 769:01[?], 770:01[?], 771:01[?], 772:01[?], 773:01[?], 774:01[?],
775:01[?], 776:01[?], 777:01N[?], 778:01[?], 779:01[?], 780:01[?], 781:01[?],
782:01[?], 783:01[?], 785:01[?], 786:01:02N[?], 787:01[?], 788:01[?], 789:01[?],
790:01[?], 791:01[?], 792:01N[?], 794:01N[?], 795:01[?], 796:01[?], 797:01[?],
798:01[?], 799:01[?], 800:01N[?], 801:01[?], 802:01[?], 803:01[?], 805:01[?],
807:01[?], 808:01[?], 810:01[?], 811:01[?], 812:01[?], 813:01[?], 814:01[?],
815:01[?], 816:01[?], 817:01[?], 818:01[?], 820:01[?], 821:01N[?], 822:01[?],
824:01:01[?]-824:01:02[?], 825:01[?], 826:01[?], 827:01[?], 828:01[?], 829:01[?],
830:01[?], 831:01N[?], 833:01[?], 834:01, 836:01[?], 837:01[?], 838:01N[?],
839:01[?], 840:01[?], 841:01[?], 842:01[?], 843:01[?], 844:01N[?], 845:01[?],
846:01[?], 847:01[?], 848:01[?], 850:01[?], 851:01[?], 852:01[?], 853:01[?],
854:01[?], 855:01[?], 856:01[?], 857:01[?], 858:01[?], 859:01[?], 860:01[?],
861:01[?], 862:01N[?], 863:01[?], 864:01[?], 865:01N[?], 866:01N[?], 867:01N[?],
868:01N[?], 869:01N[?], 870:01N[?], 871:01N[?], 872:01N[?], 873:01N[?],
874:01N[?], 875:01N[?], 876:01N[?], 877:01N[?], 883:01[?], 884:01[?], 885:01[?],
886:01[?], 887:01[?], 888:01[?], 889:01[?], 890:01[?], 891:01[?], 892:01[?],
893:01[?], 894:01N[?], 895:01[?], 897:01[?], 898:01[?], 899:01[?], 900:01[?],
901:01[?], 902:01[?], 903:01[?], 904:01[?], 908:01[?], 910:01[?], 911:01N[?],
912:01[?], 913:01[?], 914:01[?], 915:01[?], 916:01[?], 917:01N[?], 918:01[?],
919:01N[?], 920:01[?], 921:01[?], 922:01[?], 923:01[?], 924:01[?], 925:01N[?],
926:01[?], 927:01[?], 928:01[?], 929:01[?], 930:01[?], 931:01[?], 932:01[?],
933:01[?], 934:01Q[?], 935:01Q[?], 936:01Q[?], 937:01[?], 940:01[?], 941:01N[?],
942:01[?], 943:01[?], 944:01[?], 945:01[?], 946:01, 951:01, 960:01N[?],
961:01[?], 962:01[?], 963:01[?], 964:01[?], 965:01:01:01-965:01:01:02,
969:01[?], 970:01[?], 971:01[?], 972:01[?], 973:01[?], 974:01N[?], 976:01[?],
977:01[?], 980:01[?], 982:01[?], 984:01N[?], 985:01N[?], 986:01N[?], 987:01[?],
988:01[?], 989:01[?], 990:01[?], 991:01[?], 992:01[?], 993:01[?], 994:01[?],
995:01N[?], 996:01[?], 997:01[?], 998:01[?], 999:01[?], 1000:01[?], 1001:01[?],
1002:01[?], 1003:01[?], 1004:01[?], 1005:01[?], 1006:01[?], 1007:01[?],
1008:01[?], 1009:01[?], 1010:01[?], 1012:01[?], 1013:01[?], 1014:01[?],
1015:01[?], 1016:01[?], 1017:01[?], 1018:01[?], 1019:01[?], 1020:01[?],
1021:01[?], 1022:01[?], 1023:01[?], 1024:01[?], 1025:01[?], 1026:01[?],
1027:01[?], 1028:01[?], 1029:01N[?], 1030:01[?], 1031:01[?], 1034:01,
1037:01, 1038:01Q[?], 1039:01[?], 1041:01N[?], 1042:01[?], 1043:01[?],
1045:01N[?], 1046:01[?], 1047:01[?], 1048:01[?], 1050:01, 1056:01[?],
1058:01, 1059:01[?], 1060:01[?], 1061:01, 1062:01[?], 1063:01[?], 1064:01[?],



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			1065:01 [?] , 1068:01 [?] , 1070:01N [?] , 1073:01 [?] , 1074:01 [?] , 1075:01 [?] , 1076:01 [?] , 1078:01 [?] , 1079:01N [?] , 1080:01 [?] , 1092:01Q [?] , 1093:01 [?] , 1094:01 [?] , 1095:01, 1099:01, 1101:01, 1112:01N [?] , 1114:01 [?] , 1115:01 [?] , 1117:01, 1118:01, 1119:01, 1120:01, 1121:01N [?] , 1133:01, 1138:01 [?] , 1141:01 [?] , 1143:01 [?] , 1149:01, 1150:01 [?] , 1154:01N [?] , 1156:01, 1158:01 [?] , 1159:01 [?] , 1160:01 [?] , 1161:01 [?] , 1162:01 [?] , 1163:01 [?] , 1164:01 [?] , 1166:01 [?] , 1167:01 [?] , 1169:01N [?] , 1172:01 [?] , 1174:01 [?] , 1175:01 [?] , 1176:01 [?] , 1178:01 [?] , 1180:01 [?] , 1181:01 [?] , 1182:01 [?] , 1183:01 [?] , 1184:01 [?] , 1186:01 [?] , 1187:01 [?] , 1189:01, 1190:01 [?] , 1194:01, 1195:01 [?]
60	255 bp	430 bp	*02:01:02:01-02:02:04, 04:01:01:01-04:01:52, 08:01 [?] , 09:01:03 [?] - 09:01:04 [?] , 10:01:02 [?] -10:01:03 [?] , 11:01:02 [?] -11:01:04 [?] , 14:01:02 [?] - 14:01:03 [?] , 15:01:02 [?] , 17:01:01:01-17:01:04, 23:01:01:01-23:01:01:02, 23:01:02 [?] , 24:01, 26:01:01 [?] , 26:01:03 [?] , 28:01, 30:01:01:01- 30:01:01:02, 31:01:01:01-31:01:01:02, 32:01 [?] , 33:01:01:01- 33:01:01:05, 37:01 [?] , 39:01:01:01-39:01:03, 40:01:01:01-40:01:02, 41:01:01:01-41:01:01:02, 41:01:02 [?] , 44:01 [?] , 46:01:01:01-46:01:01:02, 46:01:02 [?] , 47:01:01:01-47:01:01:03, 48:01, 49:01:01:01-49:01:01:03, 51:01:01:01-51:01:01:02, 52:01 [?] , 53:01, 55:01:01:01-55:01:03, 56:01 [?] , 58:01 [?] , 60:01 [?] , 61:01N [?] , 62:01 [?] , 65:01:01 [?] -65:01:02 [?] , 66:01 [?] , 67:01 [?] , 68:01 [?] , 70:01 [?] , 71:01:01, 71:01:02 [?] , 72:01:01:01-72:01:01:03, 73:01 [?] , 74:01 [?] , 75:01 [?] , 79:01 [?] , 81:01:01:01-81:01:02, 82:01 [?] , 83:01 [?] , 84:01 [?] , 85:01:02 [?] , 86:01 [?] , 87:01 [?] , 89:01 [?] , 90:01:02 [?] , 92:01 [?] , 94:01 [?] , 95:01 [?] , 96:01 [?] , 97:01 [?] , 98:01 [?] , 99:01 [?] , 100:01 [?] , 101:01 [?] , 102:01 [?] , 103:01 [?] , 105:01:01:01-105:01:01:12, 106:01, 108:01 [?] , 109:01, 110:01 [?] , 111:01 [?] , 112:01 [?] , 113:01 [?] , 114:01 [?] , 115:01, 116:01 [?] , 117:01 [?] , 118:01 [?] , 119:01 [?] , 120:01N [?] , 121:01, 122:01 [?] , 123:01, 124:01:01:01- 124:01:02:03, 125:01 [?] , 128:01, 129:01 [?] , 133:01, 134:01, 137:01, 139:01 [?] , 140:01 [?] , 141:01, 142:01 [?] , 143:01 [?] , 144:01 [?] , 145:01 [?] , 146:01 [?] , 147:01 [?] , 148:01 [?] , 149:01 [?] , 150:01 [?] , 151:01 [?] , 153:01 [?] , 154:01N [?] , 155:01:01 [?] -155:01:02 [?] , 156:01 [?] , 157:01 [?] , 158:01 [?] , 159:01N [?] , 160:01 [?] , 161:01N [?] , 162:01:01-162:01:02, 163:01 [?] , 164:01 [?] , 166:01 [?] , 167:01 [?] , 168:01, 169:01 [?] , 170:01 [?] , 171:01 [?] , 172:01 [?] , 173:01 [?] , 174:01 [?] , 175:01 [?] , 176:01 [?] , 177:01 [?] , 178:01, 179:01 [?] , 180:01 [?] , 181:01 [?] , 182:01 [?] , 183:01 [?] , 185:01 [?] , 186:01 [?] , 187:01, 188:01 [?] , 189:01 [?] , 191:01, 192:01 [?] , 193:01 [?] , 194:01 [?] , 195:01 [?] , 196:01 [?] , 197:01 [?] , 198:01 [?] , 199:01 [?] , 200:01 [?] , 201:01 [?] , 202:01 [?] , 203:01:01 [?] -203:01:02 [?] , 204:01 [?] , 205:01 [?] , 206:01 [?] , 207:01 [?] , 209:01 [?] , 210:01 [?] , 211:01 [?] , 212:01 [?] , 213:01:01 [?] , 213:01:02, 214:01 [?] , 215:01 [?] , 216:01N [?] , 217:01 [?] , 218:01N [?] , 219:01 [?] , 220:01 [?] , 221:01 [?] , 222:01 [?] , 223:01 [?] , 224:01 [?] , 225:01, 226:01 [?] , 227:01:01 [?] , 227:01:02-227:01:03, 228:01 [?] , 229:01 [?] , 230:01 [?] , 231:01 [?] , 232:01 [?] , 233:01 [?] , 234:01 [?] , 235:01 [?] , 236:01:01, 236:01:02 [?] , 237:01 [?] , 238:01 [?] , 239:01 [?] , 240:01 [?] , 241:01 [?] , 242:01 [?] , 243:01 [?] , 245:01 [?] , 246:01 [?] , 247:01 [?] , 248:01 [?] , 249:01 [?] , 250:01 [?] , 251:01 [?] , 252:01 [?] , 253:01 [?] , 254:01 [?] , 255:01 [?] , 256:01 [?] , 258:01 [?] , 259:01 [?] , 261:01 [?] , 262:01 [?] , 263:01 [?] , 264:01 [?] , 265:01 [?] , 266:01 [?] , 267:01 [?] , 268:01 [?] , 269:01 [?] , 270:01 [?] , 271:01 [?] , 272:01 [?] , 273:01 [?] , 274:01 [?] , 275:01 [?] , 276:01 [?] , 277:01 [?] , 278:01 [?] , 279:01:01 [?] -279:01:02 [?] , 280:01 [?] , 281:01 [?] , 282:01 [?] , 283:01 [?] , 284:01 [?] , 285:01 [?] , 286:01 [?] , 287:01 [?] , 288:01 [?] , 289:01 [?] , 290:01 [?] , 291:01 [?] , 292:01 [?] , 293:01 [?] , 294:01 [?] , 295:01 [?] , 297:01 [?] , 298:01 [?] , 299:01 [?] , 300:01 [?] , 301:01 [?] , 302:01 [?] , 303:01 [?] , 304:01 [?] , 305:01 [?] , 306:01 [?] , 307:01 [?] , 308:01 [?] , 309:01 [?] , 310:01 [?] , 311:01 [?] ,

101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

312:01[?], 313:01[?], 314:01[?], 315:01[?], 316:01[?], 317:01[?], 318:01[?],
319:01[?], 320:01[?], 321:01[?], 322:01[?], 323:01[?], 324:01[?], 325:01[?],
326:01[?], 327:01[?], 329:01[?], 330:01[?], 331:01[?], 332:01[?], 333:01[?],
334:01:01[?], 334:01:02, 335:01[?], 336:01[?], 337:01[?], 338:01[?], 339:01[?],
340:01[?], 341:01[?], 342:01[?], 343:01[?], 344:01[?], 345:01[?], 346:01[?],
347:01[?], 348:01:01[?]-348:01:02[?], 349:01[?], 354:01:01[?]-354:01:02[?],
355:01[?], 356:01[?], 357:01N[?], 358:01[?], 359:01[?], 360:01[?], 361:01[?],
362:01[?], 363:01[?], 364:01[?], 365:01[?], 366:01[?], 367:01[?], 368:01[?],
369:01[?], 370:01[?], 371:01[?], 372:01[?], 373:01[?], 374:01[?], 375:01[?],
376:01[?], 377:01[?], 378:01[?], 379:01[?], 380:01[?], 381:01[?], 382:01N[?],
383:01[?], 384:01[?], 385:01[?], 386:01[?], 387:01[?], 388:01[?], 389:01[?],
390:01[?], 391:01[?], 392:01[?], 393:01[?], 394:01[?], 395:01[?], 396:01[?],
397:01[?], 398:01, 399:01[?], 400:01[?], 401:01N[?], 402:01[?], 403:01N[?],
404:01[?], 405:01[?], 406:01[?], 407:01[?], 408:01[?], 409:01[?], 410:01[?],
411:01[?], 412:01[?], 413:01[?], 415:01, 418:01[?], 419:01[?], 420:01[?], 421:01[?],
422:01[?], 423:01:01[?]-423:01:02[?], 424:01[?], 425:01[?], 426:01[?], 427:01[?],
428:01[?], 429:01[?], 430:01[?], 431:01[?], 432:01[?], 433:01[?], 434:01[?],
435:01[?], 436:01[?], 437:01[?], 438:01[?], 439:01[?], 441:01[?], 442:01[?],
443:01[?], 444:01[?], 445:01[?], 446:01[?], 447:01[?], 448:01[?], 449:01[?],
450:01N[?], 451:01[?], 452:01[?], 453:01[?], 454:01[?], 455:01N[?], 456:01[?],
457:01[?], 458:01:01[?], 459:01, 464:01, 465:01[?], 466:01[?], 467:01[?],
468:01[?], 469:01[?], 470:01[?], 471:01[?], 472:01[?], 473:01[?], 474:01[?],
475:01[?], 476:01[?], 477:01[?], 478:01[?], 479:01[?], 480:01[?], 481:01[?],
482:01[?], 483:01[?], 484:01[?], 485:01[?], 486:01[?], 487:01[?], 488:01[?],
489:01[?], 490:01[?], 491:01[?], 492:01[?], 493:01[?], 494:01[?], 495:01[?],
496:01[?], 497:01[?], 499:01[?], 500:01[?], 501:01[?], 502:01[?], 503:01[?],
504:01[?], 505:01[?], 506:01[?], 507:01N[?], 508:01[?], 509:01[?], 510:01[?],
511:01[?], 512:01[?], 513:01[?], 514:01[?], 515:01[?], 516:01[?], 517:01[?],
520:01[?], 521:01[?], 522:01[?], 523:01:01[?], 523:01:02, 524:01[?], 525:01[?],
526:01[?], 527:01[?], 528:01[?], 529:01[?], 530:01[?], 531:01[?], 532:01[?], 534:01,
535:01, 537:01[?], 538:01[?], 539:01[?], 540:01[?], 541:01[?], 542:01[?], 543:01[?],
544:01[?], 545:01[?], 546:01[?], 548:01[?], 549:01[?], 550:01[?], 551:01N[?],
552:01[?], 553:01[?], 554:01[?], 555:01[?], 556:01[?], 557:01[?], 558:01, 559:01[?],
560:01[?], 561:01[?], 562:01[?], 563:01[?], 564:01[?], 565:01[?], 566:01[?],
567:01[?], 568:01[?], 569:01[?], 570:01N[?], 571:01, 573:01[?], 574:01[?],
575:01[?], 576:01[?], 577:01[?], 578:01[?], 579:01[?], 580:01[?], 581:01[?],
582:01[?], 583:01[?], 586:01[?], 587:01[?], 588:01[?], 589:01[?], 590:01[?],
591:01[?], 592:01[?], 593:01[?], 594:01[?], 595:01[?], 596:01[?], 597:01[?],
598:01N[?], 599:01[?], 600:01[?], 601:01[?], 602:01[?], 603:01[?], 604:01[?],
605:01[?], 606:01[?], 607:01[?], 608:01[?], 609:01[?], 610:01[?], 611:01[?],
612:01[?], 614:01[?], 615:01, 617:01, 618:01, 619:01[?], 620:01[?], 622:01[?],
623:01[?], 624:01[?], 625:01[?], 626:01[?], 627:01[?], 628:01[?], 629:01[?],
630:01[?], 631:01[?], 632:01[?], 633:01[?], 634:01[?], 635:01[?], 636:01[?],
637:01[?], 638:01[?], 639:01[?], 640:01, 641:01[?], 642:01[?], 643:01[?], 644:01[?],
645:01[?], 646:01[?], 651:01, 655:01[?], 656:01, 658:01[?], 659:01[?], 660:01[?],
662:01[?], 663:01[?], 664:01[?], 665:01, 666:01, 670:01, 677:01, 678:01,
679:01, 680:01[?], 681:01, 682:01, 683:01, 684:01[?], 685:01[?], 686:01[?],
687:01, 688:01[?], 689:01[?], 690:01[?], 691:01N[?], 692:01[?], 693:01N[?],
694:01[?], 696:01N, 697:01Q[?], 698:01, 699:01, 700:01N, 712:01N[?],
716:01[?], 717:01[?], 718:01[?], 719:01, 720:01[?], 722:01, 723:01, 728:01,
738:01N, 739:01, 742:01, 743:01N, 744:01, 747:01, 748:01N, 749:01[?],
750:01, 751:01, 753:01[?], 754:01N, 755:01, 756:01N, 757:01, 758:01[?],



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			760:01 [?] , 761:01, 762:01, 763:01, 765:01, 766:01, 767:01, 769:01, 770:01, 771:01, 772:01, 773:01, 774:01, 779:01, 780:01 [?] , 781:01, 782:01, 783:01, 784:01, 786:01:01N-786:01:02N, 787:01, 788:01 [?] , 789:01 [?] , 795:01, 796:01 [?] , 799:01, 803:01, 804:01, 806:01:01:01-806:01:01:02, 808:01, 811:01, 812:01, 813:01, 814:01, 816:01, 819:01, 820:01, 824:01:01-824:01:02, 826:01, 827:01, 828:01, 830:01, 831:01N, 835:01:01:01-835:01:01:03, 836:01, 837:01 [?] , 838:01N, 840:01 [?] , 842:01, 843:01, 844:01N, 845:01, 846:01, 849:01, 850:01, 852:01, 857:01, 859:01, 861:01, 863:01, 865:01N, 867:01N, 868:01N [?] , 869:01N [?] , 870:01N [?] , 871:01N, 872:01N [?] , 873:01N [?] , 874:01N, 875:01N, 877:01N, 878:01N, 880:01, 882:01, 884:01, 889:01, 890:01, 891:01 [?] , 893:01 [?] , 894:01N, 895:01 [?] , 899:01, 901:01, 902:01, 903:01 [?] , 904:01, 905:01, 907:01, 910:01, 911:01N [?] , 913:01, 914:01, 915:01, 916:01, 919:01N, 921:01, 928:01, 932:01, 934:01Q, 939:01, 940:01, 942:01, 943:01, 944:01, 955:01, 956:01, 957:01, 960:01N [?] , 962:01, 963:01, 964:01, 966:01, 967:01, 972:01, 973:01, 974:01N, 975:01, 978:01, 988:01, 989:01, 992:01, 993:01, 997:01, 1001:01, 1002:01, 1003:01, 1004:01, 1006:01 [?] , 1010:01, 1011:01, 1012:01, 1013:01, 1016:01, 1017:01, 1023:01, 1026:01, 1027:01, 1028:01, 1031:01, 1032:01, 1033:01, 1036:01, 1039:01, 1040:01, 1041:01N [?] , 1042:01, 1044:01N, 1046:01, 1047:01, 1048:01, 1051:01, 1052:01, 1053:01, 1055:01, 1056:01, 1059:01, 1060:01, 1062:01, 1064:01, 1066:01, 1072:01, 1074:01, 1077:01, 1078:01 [?] , 1079:01N [?] , 1080:01 [?] , 1081:01, 1082:01, 1086:01, 1089:01, 1091:01, 1094:01, 1097:01, 1100:01, 1102:01, 1106:01, 1107:01, 1108:01, 1109:01, 1112:01N, 1113:01, 1115:01, 1121:01N, 1122:01, 1129:01, 1132:01, 1135:01N, 1138:01 [?] , 1139:01, 1140:01, 1141:01 [?] , 1146:01, 1150:01 [?] , 1152:01, 1154:01N [?] , 1155:01, 1158:01 [?] , 1159:01 [?] , 1160:01, 1161:01, 1163:01 [?] , 1164:01, 1166:01 [?] , 1167:01, 1171:01, 1173:01, 1174:01 [?] , 1175:01, 1176:01 [?] , 1177:01, 1178:01 [?] , 1179:01, 1181:01, 1182:01 [?] , 1184:01, 1188:01, 1190:01, 1191:01N, 1193:01N, 1195:01 [?]
61	145 bp 225 bp	430 bp	*168:01 *459:01
62⁷	260 bp	515 bp	*01:01:01:01-01:01:07, 03:01:01:01-03:01:14, 05:01:01:01-05:01:14, 06:01:01:01-06:01:01:05, 06:01:02 [?] -06:01:03 [?] , 06:01:04-06:01:05, 08:01 [?] , 09:01:01-09:01:05, 10:01:01:01-10:01:01:02, 10:01:02 [?] -10:01:03 [?] , 10:01:05, 11:01:01:01-11:01:05, 13:01:01:01-13:01:04, 14:01:01:01-14:01:01:05, 14:01:02 [?] -14:01:03 [?] , 14:01:04-14:01:12, 15:01:01:01 ^w -15:01:01:04 ^w , 15:01:02 [?] , 15:01:03 ^w -15:01:04 ^w , 16:01:01:01-16:01:03, 18:01:01:01-18:01:02, 18:01:03 ^w , 19:01:01:01-19:01:01:03, 20:01:01:01-20:01:02, 20:01:03 [?] , 21:01, 22:01:01:01-22:01:01:02, 25:01, 26:01:01 [?] , 26:01:02, 26:01:03 [?] , 27:01, 29:01, 32:01 [?] , 34:01:01:01-34:01:02, 35:01:01, 36:01, 37:01 [?] , 38:01, 41:01:02 [?] , 44:01 [?] , 45:01, 46:01:02 [?] , 50:01, 52:01 [?] , 54:01, 56:01 [?] , 57:01, 58:01 [?] , 60:01 [?] , 61:01N [?] , 62:01 [?] , 63:01, 64:01N, 65:01:01 [?] -65:01:02 [?] , 66:01 [?] , 67:01 [?] , 68:01 [?] , 69:01:01:01-69:01:01:02, 70:01 [?] , 71:01:02 [?] , 73:01 [?] , 74:01 [?] , 75:01 [?] , 76:01, 78:01, 79:01 [?] , 82:01 [?] , 83:01 [?] , 84:01 [?] , 85:01:01:01-85:01:01:02, 85:01:02 [?] , 86:01 [?] , 87:01 [?] , 88:01, 89:01 [?] , 90:01:01, 90:01:02 [?] , 91:01:01:01-91:01:02, 92:01 [?] , 93:01, 94:01 [?] , 95:01 [?] , 96:01 [?] , 97:01 [?] , 98:01 [?] , 99:01 [?] , 100:01 [?] , 101:01 [?] , 102:01 [?] , 103:01 [?] , 104:01:01:01-104:01:04, 107:01, 108:01 [?] , 110:01 [?] , 111:01 [?] , 112:01 [?] , 113:01 [?] , 114:01 [?] , 116:01 [?] , 117:01 [?] , 118:01 [?] ,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

119:01[?], 120:01N[?], 122:01[?], 125:01[?], 127:01, 129:01[?], 130:01,
131:01:01:01-131:01:01:03, 131:01:02^w, 131:01:03, 132:01, 135:01,
136:01, 139:01[?], 140:01[?], 142:01[?], 143:01[?], 144:01[?], 145:01[?], 146:01[?],
147:01[?], 148:01[?], 149:01[?], 150:01[?], 151:01[?], 152:01, 153:01[?],
154:01N[?], 155:01:01[?]-155:01:02[?], 156:01[?], 157:01[?], 158:01[?],
159:01N[?], 160:01[?], 161:01N[?], 163:01[?], 164:01[?], 165:01, 166:01[?],
167:01[?], 169:01[?], 170:01[?], 171:01[?], 172:01[?], 173:01[?], 174:01[?],
175:01[?], 176:01[?], 177:01[?], 179:01[?], 180:01[?], 181:01[?], 182:01[?],
183:01[?], 185:01[?], 186:01[?], 188:01[?], 189:01[?], 192:01[?], 193:01[?],
194:01[?], 195:01[?], 196:01[?], 197:01[?], 198:01[?], 199:01[?], 200:01[?],
201:01[?], 202:01[?], 203:01:01[?]-203:01:02[?], 204:01[?], 205:01[?], 206:01[?],
207:01[?], 208:01, 209:01[?], 210:01[?], 211:01[?], 212:01[?], 213:01:01[?],
214:01[?], 215:01[?], 216:01N[?], 217:01[?], 218:01N[?], 219:01[?], 220:01[?],
221:01[?], 222:01[?], 223:01[?], 224:01[?], 226:01[?], 227:01:01[?], 228:01[?],
229:01[?], 230:01[?], 231:01[?], 232:01[?], 233:01[?], 234:01[?], 235:01[?],
236:01:02[?], 237:01[?], 238:01[?], 239:01[?], 240:01[?], 241:01[?], 242:01[?],
243:01[?], 244:01, 245:01[?], 246:01[?], 247:01[?], 248:01[?], 249:01[?], 250:01[?],
251:01[?], 252:01[?], 253:01[?], 254:01[?], 255:01[?], 256:01[?], 258:01[?],
259:01[?], 261:01[?], 262:01[?], 263:01[?], 264:01[?], 265:01[?], 266:01[?],
267:01[?], 268:01[?], 269:01[?], 270:01[?], 271:01[?], 272:01[?], 273:01[?],
274:01[?], 275:01[?], 276:01[?], 277:01[?], 278:01[?], 279:01:01[?]-279:01:02[?],
280:01[?], 281:01[?], 282:01[?], 283:01[?], 284:01[?], 285:01[?], 286:01[?],
287:01[?], 288:01[?], 289:01[?], 290:01[?], 291:01[?], 292:01[?], 293:01[?],
294:01[?], 295:01[?], 296:01, 297:01[?], 298:01[?], 299:01[?], 300:01[?], 301:01[?],
302:01[?], 303:01[?], 304:01[?], 305:01[?], 306:01[?], 307:01[?], 308:01[?],
309:01[?], 310:01[?], 311:01[?], 312:01[?], 313:01[?], 314:01[?], 315:01[?],
316:01[?], 317:01[?], 318:01[?], 319:01[?], 320:01[?], 321:01[?], 322:01[?],
323:01[?], 324:01[?], 325:01[?], 326:01[?], 327:01[?], 328:01N, 329:01[?],
330:01[?], 331:01[?], 332:01[?], 333:01[?], 334:01:01[?], 335:01[?], 336:01[?],
337:01[?], 338:01[?], 339:01[?], 340:01[?], 341:01[?], 342:01[?], 343:01[?],
344:01[?], 345:01[?], 346:01[?], 347:01[?], 348:01:01[?]-348:01:02[?], 349:01[?],
350:01, 352:01:01-352:01:02:02, 353:01, 354:01:01[?]-354:01:02[?],
355:01[?], 356:01[?], 357:01N[?], 358:01[?], 359:01[?], 360:01[?], 361:01[?],
362:01[?], 363:01[?], 364:01[?], 365:01[?], 366:01[?], 367:01[?], 368:01[?],
369:01[?], 370:01[?], 371:01[?], 372:01[?], 373:01[?], 374:01[?], 375:01[?],
376:01[?], 377:01[?], 378:01[?], 379:01[?], 380:01[?], 381:01[?], 382:01N[?],
383:01[?], 384:01[?], 385:01[?], 386:01[?], 387:01[?], 388:01[?], 389:01[?],
390:01[?], 391:01[?], 392:01[?], 393:01[?], 394:01[?], 395:01[?], 396:01[?],
397:01[?], 399:01[?], 400:01[?], 401:01N[?], 402:01[?], 403:01N[?], 404:01[?],
405:01[?], 406:01[?], 407:01[?], 408:01[?], 409:01[?], 410:01[?], 411:01[?],
412:01[?], 413:01[?], 414:01:01:01-414:01:01:02, 417:01:01-417:01:02,
418:01[?], 419:01[?], 420:01[?], 421:01[?], 422:01[?], 423:01:01[?]-423:01:02[?],
424:01[?], 425:01[?], 426:01[?], 427:01[?], 428:01[?], 429:01[?], 430:01[?],
431:01[?], 432:01[?], 433:01[?], 434:01[?], 435:01[?], 436:01[?], 437:01[?],
438:01[?], 439:01[?], 440:01, 441:01[?], 442:01[?], 443:01[?], 444:01[?], 445:01[?],
446:01[?], 447:01[?], 448:01[?], 449:01[?], 450:01N[?], 451:01[?], 452:01[?],
453:01[?], 454:01[?], 455:01N[?], 456:01[?], 457:01[?], 458:01:01[?], 458:01:02,
461:01, 463:01:01:01-463:01:01:03, 465:01[?], 466:01[?], 467:01[?],
468:01[?], 469:01[?], 470:01[?], 471:01[?], 472:01[?], 473:01[?], 474:01[?],
475:01[?], 476:01[?], 477:01[?], 478:01[?], 479:01[?], 480:01[?], 481:01[?],
482:01[?], 483:01[?], 484:01[?], 485:01[?], 486:01[?], 487:01[?], 488:01[?],
489:01[?], 490:01[?], 491:01[?], 492:01[?], 493:01[?], 494:01[?], 495:01[?],



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

496:01[?], 497:01[?], 498:01, 499:01[?], 500:01[?], 501:01[?], 502:01[?], 503:01[?],
504:01[?], 505:01[?], 506:01[?], 507:01N[?], 508:01[?], 509:01[?], 510:01[?],
511:01[?], 512:01[?], 513:01[?], 514:01[?], 515:01[?], 516:01[?], 517:01[?], 518:01,
519:01, 520:01[?], 521:01[?], 522:01[?], 523:01:01[?], 524:01[?], 525:01[?],
526:01[?], 527:01[?], 528:01[?], 529:01[?], 530:01[?], 531:01[?], 532:01[?], 533:01,
536:01, 537:01[?], 538:01[?], 539:01[?], 540:01[?], 541:01[?], 542:01[?], 543:01[?],
544:01[?], 545:01[?], 546:01[?], 547:01, 548:01[?], 549:01[?], 550:01[?],
551:01N[?], 552:01[?], 553:01[?], 554:01[?], 555:01[?], 556:01[?], 557:01[?],
559:01[?], 560:01[?], 561:01[?], 562:01[?], 563:01[?], 564:01[?], 565:01[?],
566:01[?], 567:01[?], 568:01[?], 569:01[?], 570:01N[?], 572:01, 573:01[?],
574:01[?], 575:01[?], 576:01[?], 577:01[?], 578:01[?], 579:01[?], 580:01[?],
581:01[?], 582:01[?], 583:01[?], 584:01:01:01-584:01:02:02, 585:01:01:01^w-
585:01:01:02^w, 586:01[?], 587:01[?], 588:01[?], 589:01[?], 590:01[?], 591:01[?],
592:01[?], 593:01[?], 594:01[?], 595:01[?], 596:01[?], 597:01[?], 598:01N[?],
599:01[?], 600:01[?], 601:01[?], 602:01[?], 603:01[?], 604:01[?], 605:01[?],
606:01[?], 607:01[?], 608:01[?], 609:01[?], 610:01[?], 611:01[?], 612:01[?], 613:01,
614:01[?], 616:01, 619:01[?], 620:01[?], 621:01, 622:01[?], 623:01[?], 624:01[?],
625:01[?], 626:01[?], 627:01[?], 628:01[?], 629:01[?], 630:01[?], 631:01[?],
632:01[?], 633:01[?], 634:01[?], 635:01[?], 636:01[?], 637:01[?], 638:01[?],
639:01[?], 641:01[?], 642:01[?], 643:01[?], 644:01[?], 645:01[?], 646:01[?],
648:01:01:01-648:01:01:02, 649:01, 650:01, 652:01, 653:01, 654:01,
655:01[?], 657:01N, 658:01[?], 659:01[?], 660:01[?], 661:01N, 662:01[?],
663:01[?], 664:01[?], 667:01, 668:01:01:01-668:01:01:02, 669:01, 671:01,
672:01, 673:01, 675:01, 676:01, 680:01[?], 684:01[?], 685:01[?], 686:01[?],
688:01[?], 689:01[?], 690:01[?], 691:01N[?], 692:01[?], 693:01N[?], 694:01[?],
695:01^w, 697:01Q[?], 702:01, 703:01, 704:01, 705:01, 706:01, 707:01,
708:01, 709:01, 710:01, 711:01, 712:01N[?], 713:01, 714:01, 715:01,
716:01[?], 717:01[?], 718:01[?], 720:01[?], 727:01, 729:01, 732:01N, 733:01,
737:01, 740:01, 745:01, 746:01, 749:01[?], 752:01, 753:01[?], 758:01[?],
760:01[?], 764:01, 768:01^w, 778:01, 780:01[?], 788:01[?], 789:01[?], 791:01,
792:01N, 794:01N, 796:01[?], 797:01, 798:01, 800:01N, 801:01, 802:01,
807:01, 810:01, 815:01, 821:01N^w, 822:01, 825:01, 829:01, 832:01,
833:01, 834:01, 837:01[?], 839:01, 840:01[?], 847:01, 848:01, 851:01,
855:01, 856:01, 860:01, 862:01N, 864:01, 866:01N, 868:01N[?],
869:01N[?], 870:01N[?], 872:01N[?], 873:01N[?], 876:01N, 879:01:01:01-
879:01:01:03, 883:01, 886:01, 888:01, 891:01[?], 892:01, 893:01[?],
895:01[?], 896:01, 897:01, 898:01, 900:01, 903:01[?], 908:01, 909:01,
911:01N[?], 912:01, 917:01N, 918:01, 920:01, 922:01, 923:01, 924:01,
929:01, 930:01, 931:01, 935:01Q, 938:01, 941:01N, 945:01, 946:01,
947:01, 948:01, 949:01, 950:01N, 951:01, 952:01, 953:01, 960:01N[?],
965:01:01:01-965:01:01:02, 969:01, 970:01, 971:01, 976:01, 977:01,
979:01, 980:01, 982:01, 983:01, 984:01N, 986:01N, 987:01, 990:01,
991:01, 994:01, 998:01, 999:01, 1000:01, 1006:01[?], 1007:01, 1008:01,
1009:01, 1014:01, 1015:01, 1018:01, 1019:01, 1021:01, 1022:01,
1024:01, 1029:01N, 1030:01, 1034:01, 1037:01, 1038:01Q,
1041:01N[?], 1043:01, 1045:01N, 1049:01, 1050:01, 1054:01^w, 1057:01,
1058:01, 1061:01, 1063:01, 1065:01, 1067:01, 1068:01, 1069:01,
1070:01N, 1071:01, 1073:01, 1076:01, 1078:01[?], 1079:01N[?],
1080:01[?], 1087:01, 1088:01, 1090:01, 1093:01, 1095:01, 1096:01,
1098:01N, 1099:01, 1101:01, 1103:01, 1104:01, 1105:01, 1111:01,
1114:01, 1116:01, 1117:01, 1118:01, 1119:01, 1120:01, 1123:01,
1125:01, 1126:01, 1127:01, 1128:01, 1130:01, 1131:01, 1133:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			1134:01, 1136:01, 1137:01, 1138:01 [?] , 1141:01 [?] , 1142:01, 1145:01, 1147:01, 1149:01, 1150:01 [?] , 1151:01, 1154:01N [?] , 1156:01, 1157:01, 1158:01 [?] , 1159:01 [?] , 1162:01, 1163:01 [?] , 1166:01 [?] , 1168:01, 1169:01N, 1170:01, 1172:01, 1174:01 [?] , 1176:01 [?] , 1178:01 [?] , 1180:01, 1182:01 [?] , 1183:01, 1185:01, 1186:01, 1187:01, 1189:01, 1192:01 ^w , 1195:01 [?]
63	170 bp 210 bp 245 bp 275 bp	515 bp	*475:01 *183:01, 319:01 *216:01N *57:01 [?] , 280:01, 356:01, 539:01
64	150 bp 180 bp 215 bp	430 bp	*217:01, 256:01, 658:01, 1040:01 *232:01, 475:01 *226:01
65	170 bp 205 bp 240 bp	430 bp	*181:01, 257:01, 258:01 *604:01 *238:01, 1066:01
66	185 bp 235 bp 280 bp	515 bp	*276:01, 418:01, 432:01, 697:01Q *239:01, 260:01, 487:01 *57:01 [?] , 356:01, 455:01N, 539:01
67⁴	90 bp 155 bp 180 bp 210 bp	430 bp	*504:01, 781:01 *501:01, 712:01N *261:01, 380:01, 657:01N *674:01
68⁷	175 bp 200 bp 265 bp	430 bp	*228:01, 334:01:01, 522:01, 594:01 *191:01, 192:01, 310:01, 342:01 *109:01, 117:01, 158:01, 164:01, 199:01, 200:01, 271:01, 294:01, 369:01, 392:01, 440:01, 480:01, 508:01, 531:01, 561:01, 601:01, 663:01, 679:01, 943:01, 961:01, 1080:01, 1174:01
69	180 bp 215 bp 245 bp	515 bp	*227:01:01, 228:01, 646:01, 682:01 *520:01 *210:01, 211:01, 240:01, 272:01
70⁴	75 bp 115 bp 170 bp 215 bp 260 bp	430 bp	*476:01, 751:01 *551:01N *253:01, 286:01, 334:01:01, 522:01, 594:01, 660:01 *262:01, 263:01, 681:01 *196:01, 596:01
71	140 bp 180 bp 260 bp	430 bp	*208:01, 250:01 *315:01, 825:01 *268:01, 1103:01
72	145 bp 175 bp 250 bp	430 bp	*327:01 *270:01 *267:01, 293:01, 393:01, 662:01, 1007:01, 1017:01
73	170 bp 190 bp	430 bp	*231:01 *18:01:01:01-18:01:03, 28:01, 139:01, 178:01, 230:01, 279:01:01, 280:01, 292:01, 296:01, 347:01, 431:01, 489:01, 499:01, 512:01, 526:01, 644:01, 752:01, 897:01, 918:01, 942:01, 1109:01, 1136:01, 1142:01, 1165:01
74^{4,7}	115 bp 160 bp 215 bp	430 bp	*181:01, 406:01 *255:01, 274:01, 374:01, 889:01 *252:01, 273:01, 582:01, 624:01, 772:01



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

75	140 bp	430 bp	*03:01:01:01-03:01:14, 04:02:01:01-04:02:01:20, 04:02:02 ^w , 04:02:03-04:02:17, 05:01:01:01-05:01:14, 14:01:01:01-14:01:01:05, 14:01:03-14:01:12, 18:01:01:01-18:01:03, 20:01:01:01-20:01:04, 23:01:01:01-23:01:02, 25:01, 34:01:01:01 ^w -34:01:02 ^w , 35:01:01, 36:01, 38:01, 45:01, 50:01, 52:01, 57:01, 59:01, 60:01, 61:01N, 62:01, 63:01, 67:01, 68:01, 70:01, 73:01, 75:01, 77:01, 78:01, 79:01, 80:01, 82:01, 83:01, 84:01, 87:01, 91:01:01:01-91:01:02, 97:01, 98:01, 100:01, 104:01:01:01-104:01:04, 105:01:01:01-105:01:01:12, 108:01, 111:01, 112:01, 114:01, 116:01, 119:01, 124:01:01:01-124:01:02:03, 129:01, 130:01, 132:01, 135:01, 138:01, 140:01, 143:01, 144:01, 153:01, 164:01, 170:01 ^w , 171:01, 183:01, 184:01, 185:01, 190:01, 197:01, 201:01, 209:01, 210:01, 219:01, 221:01, 222:01, 223:01, 226:01, 233:01, 242:01, 243:01, 244:01, 248:01, 252:01, 256:01, 257:01, 259:01, 260:01, 263:01, 270:01, 271:01, 274:01 ^w , 284:01, 285:01, 286:01, 292:01, 293:01, 300:01, 301:01, 307:01, 308:01, 312:01, 316:01, 317:01, 321:01, 329:01, 330:01, 331:01, 333:01, 338:01, 342:01, 351:01, 364:01, 365:01, 368:01, 369:01, 373:01, 374:01 ^w , 379:01, 380:01, 383:01, 390:01, 391:01, 392:01, 394:01, 395:01, 402:01, 404:01, 405:01, 409:01, 420:01, 424:01, 429:01, 431:01, 435:01, 439:01, 443:01 ^w , 445:01, 446:01, 451:01, 463:01:01:01-463:01:01:03, 466:01, 472:01, 473:01, 477:01, 484:01 ^w , 488:01, 491:01 ^w , 492:01, 494:01, 495:01, 498:01, 507:01N, 508:01, 509:01, 525:01, 532:01, 539:01, 540:01, 545:01, 548:01, 549:01 ^w , 560:01, 568:01, 571:01, 572:01, 573:01, 577:01, 580:01 ^w , 586:01, 587:01, 589:01, 594:01, 595:01, 596:01, 598:01N, 599:01, 605:01, 606:01, 610:01, 611:01, 613:01, 619:01, 620:01, 621:01, 624:01, 627:01, 629:01, 635:01, 638:01, 645:01, 646:01, 647:01, 648:01:01:01-648:01:01:02, 651:01, 657:01N, 664:01, 665:01, 668:01:01:01-668:01:01:02, 669:01, 671:01, 674:01, 675:01, 676:01, 688:01, 689:01, 690:01, 691:01N, 693:01N, 697:01Q, 701:01, 704:01, 705:01, 706:01, 712:01N, 715:01, 718:01, 724:01N, 725:01, 726:01, 727:01, 728:01, 729:01, 730:01, 731:01, 734:01, 735:01, 741:01, 744:01, 759:01, 760:01, 762:01, 763:01, 764:01, 775:01, 776:01, 777:01N, 782:01, 787:01, 790:01, 792:01N, 794:01N, 798:01, 801:01, 802:01, 805:01, 809:01, 817:01, 818:01, 823:01, 829:01, 832:01, 834:01, 835:01:01:01 ^w -835:01:01:03 ^w , 839:01, 841:01, 847:01, 848:01, 851:01, 854:01, 855:01, 858:01, 860:01, 876:01N, 881:01, 883:01, 885:01, 887:01, 891:01, 892:01, 897:01, 898:01, 905:01, 912:01, 913:01 ^w , 923:01, 925:01N, 927:01, 929:01, 933:01, 935:01Q, 936:01Q, 938:01, 941:01N, 942:01, 946:01, 948:01, 949:01, 950:01N, 951:01, 952:01, 954:01, 958:01, 959:01N, 961:01, 962:01, 970:01, 971:01, 972:01, 981:01, 983:01, 985:01N, 986:01N, 987:01, 990:01, 994:01, 995:01N, 1000:01, 1005:01, 1013:01, 1014:01, 1015:01, 1018:01, 1020:01, 1021:01, 1025:01, 1027:01, 1029:01N, 1030:01, 1031:01, 1035:01, 1037:01, 1041:01N, 1049:01, 1058:01, 1061:01, 1070:01N, 1071:01, 1072:01, 1075:01, 1078:01, 1083:01, 1084:01N, 1085:01, 1090:01, 1092:01Q, 1093:01, 1095:01, 1096:01, 1098:01N, 1099:01, 1103:01, 1104:01, 1106:01, 1110:01, 1114:01, 1116:01, 1117:01, 1118:01, 1119:01, 1120:01, 1124:01, 1125:01, 1130:01, 1133:01, 1134:01, 1142:01, 1143:01, 1153:01, 1156:01, 1157:01, 1165:01, 1171:01, 1172:01, 1174:01, 1180:01, 1187:01, 1189:01, 1194:01
----	--------	--------	--



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

76 ⁴	120 bp	430 bp	*215:01
	180 bp		*283:01, 359:01, 423:01:01-423:01:02, 523:01:01, 577:01, 928:01
	220 bp		*193:01, 681:01
	250 bp		*225:01
77 ⁴	120 bp	430 bp	*215:01
	170 bp		*488:01, 557:01
	220 bp		*194:01, 360:01
	250 bp		*311:01, 368:01, 687:01, 696:01N
78	145 bp	430 bp	*282:01, 1056:01
	180 bp		*318:01, 485:01, 523:01:01, 641:01
	220 bp		*195:01, 424:01, 510:01, 837:01
	250 bp		*224:01, 696:01N, 1159:01
79 ⁴	75 bp	515 bp	*476:01, 751:01
	115 bp		*551:01N
	170 bp		*275:01, 283:01, 359:01, 396:01, 577:01
	245 bp		*212:01, 224:01, 1159:01
	270 bp		*679:01
80 ⁶	160 bp	430 bp	*213:01:01, 282:01, 286:01, 396:01, 433:01, 488:01, 523:01:02, 557:01, 646:01, 660:01, 682:01, 1056:01
	250 bp		*225:01
81 ⁴	65 bp	430 bp	*591:01
	95 bp		*276:01
	145 bp		*250:01
	180 bp		*202:01, 230:01, 313:01, 1103:01
	205 bp		*281:01, 472:01
	270 bp		*268:01, 472:01, 1103:01
82 ⁴	80 bp	430 bp	*01:01:01:01-01:01:07, 02:01:02:01-02:02:03, 02:02:04 [?] , 03:01:01:01-03:01:14, 04:01:01:01-04:02:17, 06:01:01:01-06:01:01:05, 06:01:02 [?] -06:01:03 [?] , 06:01:04, 06:01:05 [?] , 08:01 [?] , 09:01:01-09:01:02, 09:01:03 [?] -09:01:05 [?] , 10:01:01:01-10:01:01:02, 10:01:02 [?] -10:01:03 [?] , 10:01:04-10:01:05, 11:01:01:01-11:01:01:02, 11:01:02 [?] -11:01:04 [?] , 11:01:05, 13:01:01:01-13:01:04, 14:01:01:01-14:01:01:05, 14:01:02 [?] -14:01:03 [?] , 14:01:04-14:01:05, 14:01:06 [?] -14:01:09 [?] , 14:01:10-14:01:12, 15:01:01:01-15:01:01:04, 15:01:02 [?] , 15:01:03-15:01:04, 16:01:01:01-16:01:01:02, 16:01:02 [?] , 16:01:03, 17:01:01:01-17:01:01:05, 17:01:02 [?] -17:01:03 [?] , 17:01:04, 18:01:01:01-18:01:01:04, 18:01:02 [?] -18:01:03 [?] , 20:01:01:01-20:01:01:02, 20:01:02 [?] -20:01:03 [?] , 20:01:04, 21:01, 23:01:01:01-23:01:01:02, 23:01:02 [?] , 24:01, 25:01, 26:01:01 [?] , 26:01:02, 26:01:03 [?] , 27:01, 28:01, 30:01:01:01-30:01:01:02, 31:01:01:01-31:01:01:02, 32:01 [?] , 33:01:01:01-33:01:01:05, 34:01:02 [?] , 35:01:01, 36:01, 37:01 [?] , 39:01:01:01-39:01:03, 40:01:01:01-40:01:01:03, 40:01:02 [?] , 41:01:01:01-41:01:01:02, 41:01:02 [?] , 44:01 [?] , 45:01, 46:01:01:01-46:01:01:02, 46:01:02 [?] , 47:01:01:01-47:01:01:03, 48:01, 49:01:01:01-49:01:01:03, 50:01, 51:01:01:01-51:01:01:02, 52:01 [?] , 53:01, 54:01, 55:01:01:01-55:01:01:05, 55:01:02 [?] -55:01:03 [?] , 56:01 [?] , 57:01, 58:01 [?] , 59:01, 60:01 [?] , 61:01N [?] , 62:01 [?] , 64:01N, 65:01:01 [?] -65:01:02 [?] , 66:01 [?] , 67:01 [?] , 68:01 [?] , 69:01:01:01-69:01:01:02, 70:01 [?] , 71:01:01, 71:01:02 [?] , 72:01:01:01-72:01:01:03, 73:01 [?] , 74:01 [?] , 75:01 [?] , 76:01, 77:01, 79:01 [?] , 80:01, 81:01:01:01-81:01:01:02, 81:01:02 [?] , 82:01 [?] , 83:01 [?] , 84:01 [?] , 85:01:01:01-85:01:01:02, 85:01:02 [?] , 86:01 [?] , 87:01 [?] , 88:01, 89:01 [?] , 90:01:01, 90:01:02 [?] ,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

91:01:01:01-91:01:01:02, 91:01:02[?], 92:01[?], 93:01, 94:01[?], 95:01[?],
96:01[?], 97:01[?], 98:01[?], 99:01[?], 100:01[?], 101:01[?], 102:01[?], 103:01[?],
105:01:01:01-105:01:01:12, 106:01, 107:01, 108:01[?], 109:01, 110:01[?],
111:01[?], 112:01[?], 113:01[?], 114:01[?], 115:01, 116:01[?], 117:01[?], 118:01[?],
119:01[?], 120:01N[?], 121:01[?], 122:01[?], 123:01[?], 124:01:01:01-
124:01:02:03, 125:01[?], 126:01:01:01-126:01:01:03, 127:01, 128:01,
129:01[?], 130:01, 131:01:01:01-131:01:01:03, 131:01:02[?]-131:01:03[?],
132:01[?], 133:01, 134:01[?], 135:01, 136:01, 137:01[?], 138:01, 139:01[?],
140:01[?], 141:01, 142:01[?], 143:01[?], 144:01[?], 145:01[?], 146:01[?], 147:01[?],
148:01[?], 149:01[?], 150:01[?], 151:01[?], 152:01[?], 153:01[?], 154:01N[?],
155:01:01[?]-155:01:02[?], 156:01[?], 157:01[?], 158:01[?], 159:01N[?], 160:01[?],
161:01N[?], 162:01:01-162:01:02, 163:01[?], 164:01[?], 165:01, 166:01[?],
167:01[?], 168:01[?], 169:01[?], 170:01[?], 171:01[?], 172:01[?], 173:01[?],
174:01[?], 175:01[?], 176:01[?], 177:01[?], 178:01, 179:01[?], 180:01[?], 181:01[?],
182:01[?], 183:01[?], 184:01, 185:01[?], 186:01[?], 187:01, 188:01[?], 189:01[?],
190:01, 191:01, 192:01[?], 193:01[?], 194:01[?], 195:01[?], 196:01[?], 197:01[?],
198:01[?], 199:01[?], 200:01[?], 201:01[?], 202:01[?], 203:01:01[?]-203:01:02[?],
204:01[?], 205:01[?], 206:01[?], 207:01[?], 208:01, 209:01[?], 210:01[?], 211:01[?],
212:01[?], 213:01:01[?]-213:01:02[?], 214:01[?], 215:01[?], 216:01N[?], 217:01[?],
218:01N[?], 219:01[?], 220:01[?], 221:01[?], 222:01[?], 223:01[?], 224:01[?],
225:01[?], 226:01[?], 227:01:01[?]-227:01:02[?], 227:01:03, 228:01[?], 229:01[?],
230:01[?], 231:01[?], 232:01[?], 233:01[?], 234:01[?], 235:01[?], 236:01:01,
236:01:02[?], 237:01[?], 238:01[?], 239:01[?], 240:01[?], 241:01[?], 242:01[?],
243:01[?], 244:01, 245:01[?], 246:01[?], 247:01[?], 248:01[?], 249:01[?], 250:01[?],
251:01[?], 252:01[?], 253:01[?], 254:01[?], 255:01[?], 256:01[?], 257:01, 258:01[?],
259:01[?], 260:01, 261:01[?], 262:01[?], 263:01[?], 264:01[?], 265:01[?], 266:01[?],
267:01[?], 268:01[?], 269:01[?], 270:01[?], 271:01[?], 272:01[?], 273:01[?],
274:01[?], 275:01[?], 276:01[?], 277:01[?], 278:01[?], 279:01:01[?]-279:01:02[?],
280:01[?], 281:01[?], 282:01[?], 283:01[?], 284:01[?], 285:01[?], 286:01[?],
287:01[?], 288:01[?], 289:01[?], 290:01[?], 291:01[?], 292:01[?], 293:01[?],
294:01[?], 295:01[?], 296:01, 297:01[?], 298:01[?], 299:01[?], 300:01[?], 301:01[?],
302:01[?], 303:01[?], 304:01[?], 305:01[?], 306:01[?], 307:01[?], 308:01[?],
309:01[?], 310:01[?], 311:01[?], 312:01[?], 313:01[?], 314:01[?], 315:01[?],
316:01[?], 317:01[?], 318:01[?], 319:01[?], 320:01[?], 321:01[?], 322:01[?],
323:01[?], 324:01[?], 325:01[?], 326:01[?], 327:01[?], 328:01N, 329:01[?],
330:01[?], 331:01[?], 332:01[?], 333:01[?], 334:01:01[?]-334:01:02[?], 335:01[?],
336:01[?], 337:01[?], 338:01[?], 339:01[?], 340:01[?], 341:01[?], 342:01[?],
343:01[?], 344:01[?], 345:01[?], 346:01[?], 347:01[?], 348:01:01[?]-348:01:02[?],
349:01[?], 350:01, 351:01, 352:01:01-352:01:02:02, 353:01, 354:01:01[?]-
354:01:02[?], 355:01[?], 356:01[?], 357:01N[?], 358:01[?], 359:01[?], 360:01[?],
361:01[?], 362:01[?], 363:01[?], 364:01[?], 365:01[?], 366:01[?], 367:01[?],
368:01[?], 369:01[?], 370:01[?], 371:01[?], 372:01[?], 373:01[?], 374:01[?],
375:01[?], 376:01[?], 377:01[?], 378:01[?], 379:01[?], 380:01[?], 381:01[?],
382:01N[?], 383:01[?], 384:01[?], 385:01[?], 386:01[?], 387:01[?], 388:01[?],
389:01[?], 390:01[?], 391:01[?], 392:01[?], 393:01[?], 394:01[?], 395:01[?],
396:01[?], 397:01[?], 398:01, 399:01[?], 400:01[?], 401:01N[?], 402:01[?],
403:01N[?], 404:01[?], 405:01[?], 406:01[?], 407:01[?], 408:01[?], 409:01[?],
410:01[?], 411:01[?], 412:01[?], 413:01[?], 415:01, 416:01:01:01-416:01:02,
417:01:01-417:01:02, 418:01[?], 419:01[?], 420:01[?], 421:01[?], 422:01[?],
423:01:01[?]-423:01:02[?], 424:01[?], 425:01[?], 426:01[?], 427:01[?], 428:01[?],
429:01[?], 430:01[?], 431:01[?], 432:01[?], 433:01[?], 434:01[?], 435:01[?],
436:01[?], 437:01[?], 438:01[?], 439:01[?], 440:01, 441:01[?], 442:01[?], 443:01[?],



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

444:01[?], 445:01[?], 446:01[?], 447:01[?], 448:01[?], 449:01[?], 450:01N[?],
451:01[?], 452:01[?], 453:01[?], 454:01[?], 455:01N[?], 456:01[?], 457:01[?],
458:01:01[?]-458:01:02[?], 459:01, 460:01, 461:01, 462:01, 463:01:01:01-
463:01:01:03, 464:01, 465:01[?], 466:01[?], 467:01[?], 468:01[?], 469:01[?],
470:01[?], 471:01[?], 472:01[?], 473:01[?], 474:01[?], 475:01[?], 476:01[?],
477:01[?], 478:01[?], 479:01[?], 480:01[?], 481:01[?], 482:01[?], 483:01[?],
484:01[?], 485:01[?], 486:01[?], 487:01[?], 488:01[?], 489:01[?], 490:01[?],
491:01[?], 492:01[?], 493:01[?], 494:01[?], 495:01[?], 496:01[?], 497:01[?], 498:01,
499:01[?], 500:01[?], 501:01[?], 502:01[?], 503:01[?], 504:01[?], 505:01[?],
506:01[?], 507:01N[?], 508:01[?], 509:01[?], 510:01[?], 511:01[?], 512:01[?],
513:01[?], 514:01[?], 515:01[?], 516:01[?], 517:01[?], 518:01, 520:01[?], 521:01[?],
522:01[?], 523:01:01[?]-523:01:02[?], 524:01[?], 525:01[?], 526:01[?], 527:01[?],
528:01[?], 529:01[?], 530:01[?], 531:01[?], 532:01[?], 533:01[?], 534:01[?], 535:01,
536:01[?], 537:01[?], 538:01[?], 539:01[?], 540:01[?], 541:01[?], 542:01[?],
543:01[?], 544:01[?], 545:01[?], 546:01[?], 548:01[?], 549:01[?], 550:01[?],
551:01N[?], 552:01[?], 553:01[?], 554:01[?], 555:01[?], 556:01[?], 557:01[?],
558:01, 559:01[?], 560:01[?], 561:01[?], 562:01[?], 563:01[?], 564:01[?], 565:01[?],
566:01[?], 567:01[?], 568:01[?], 569:01[?], 570:01N[?], 571:01, 572:01,
573:01[?], 574:01[?], 575:01[?], 576:01[?], 577:01[?], 578:01[?], 579:01[?],
580:01[?], 581:01[?], 582:01[?], 583:01[?], 584:01:01:01-584:01:02:02,
585:01:01:01-585:01:01:02, 586:01[?], 587:01[?], 588:01[?], 589:01[?],
590:01[?], 591:01[?], 592:01[?], 593:01[?], 594:01[?], 595:01[?], 596:01[?],
597:01[?], 598:01N[?], 599:01[?], 600:01[?], 601:01[?], 602:01[?], 603:01[?],
604:01[?], 605:01[?], 606:01[?], 607:01[?], 608:01[?], 609:01[?], 610:01[?],
611:01[?], 612:01[?], 613:01[?], 614:01[?], 615:01, 616:01[?], 617:01[?], 618:01[?],
619:01[?], 620:01[?], 621:01[?], 622:01[?], 623:01[?], 624:01[?], 625:01[?],
626:01[?], 627:01[?], 628:01[?], 629:01[?], 630:01[?], 631:01[?], 632:01[?],
633:01[?], 634:01[?], 635:01[?], 636:01[?], 637:01[?], 638:01[?], 639:01[?], 640:01,
641:01[?], 642:01[?], 643:01[?], 644:01[?], 645:01[?], 646:01[?], 647:01, 649:01,
651:01, 653:01, 654:01, 655:01[?], 656:01[?], 657:01N[?], 658:01[?], 659:01[?],
660:01[?], 661:01N[?], 662:01[?], 663:01[?], 664:01[?], 665:01, 666:01, 667:01,
669:01, 670:01, 671:01, 672:01, 673:01, 674:01, 675:01, 676:01,
677:01, 678:01, 679:01[?], 680:01[?], 681:01[?], 682:01[?], 683:01[?], 684:01[?],
685:01[?], 686:01[?], 687:01, 688:01[?], 689:01[?], 690:01[?], 691:01N[?],
692:01[?], 693:01N[?], 694:01[?], 695:01[?], 696:01N[?], 697:01Q[?], 698:01[?],
699:01, 700:01N, 701:01, 703:01[?], 704:01, 705:01, 706:01, 707:01[?],
708:01[?], 709:01[?], 710:01[?], 711:01[?], 712:01N[?], 713:01[?], 714:01[?],
715:01[?], 716:01[?], 717:01[?], 718:01[?], 719:01, 720:01[?], 721:01[?], 722:01[?],
723:01[?], 724:01N, 725:01[?], 726:01, 727:01[?], 728:01, 729:01[?], 730:01[?],
731:01[?], 732:01N[?], 733:01, 734:01, 735:01, 736:01, 737:01, 738:01N,
739:01[?], 740:01[?], 741:01[?], 742:01[?], 743:01N[?], 744:01[?], 745:01[?],
746:01[?], 747:01[?], 748:01N[?], 749:01[?], 750:01[?], 751:01[?], 752:01[?],
753:01[?], 754:01N[?], 755:01[?], 756:01N[?], 757:01[?], 758:01[?], 759:01[?],
760:01[?], 761:01[?], 762:01, 763:01, 765:01[?], 766:01[?], 767:01[?], 768:01[?],
769:01[?], 770:01[?], 771:01[?], 772:01[?], 773:01[?], 774:01[?], 775:01[?],
776:01[?], 777:01N[?], 778:01[?], 779:01[?], 780:01[?], 781:01[?], 782:01[?],
783:01[?], 784:01, 785:01[?], 786:01:01N, 786:01:02N[?], 787:01[?], 788:01[?],
789:01[?], 790:01[?], 791:01[?], 792:01N[?], 794:01N[?], 795:01[?], 796:01[?],
797:01[?], 798:01[?], 799:01[?], 800:01N[?], 801:01[?], 802:01[?], 803:01[?],
804:01, 805:01[?], 806:01:01:01-806:01:01:02, 807:01[?], 808:01[?],
809:01, 810:01[?], 811:01[?], 812:01[?], 813:01[?], 814:01[?], 815:01[?], 816:01[?],
817:01[?], 818:01[?], 819:01, 820:01[?], 821:01N[?], 822:01[?], 823:01,



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

824:01:01[?]-824:01:02[?], 825:01[?], 826:01[?], 827:01[?], 828:01[?], 829:01[?],
830:01[?], 831:01N[?], 832:01, 833:01[?], 835:01:01:01-835:01:01:03,
836:01[?], 837:01[?], 838:01N[?], 839:01[?], 840:01[?], 841:01[?], 842:01[?],
843:01[?], 844:01N[?], 845:01[?], 846:01[?], 847:01[?], 848:01[?], 849:01,
850:01[?], 851:01[?], 852:01[?], 853:01[?], 854:01[?], 855:01[?], 856:01[?],
857:01[?], 858:01[?], 859:01[?], 860:01[?], 861:01[?], 862:01N[?], 863:01[?],
864:01[?], 865:01N[?], 866:01N[?], 867:01N[?], 868:01N[?], 869:01N[?],
870:01N[?], 871:01N[?], 872:01N[?], 873:01N[?], 874:01N[?], 875:01N[?],
876:01N[?], 877:01N[?], 878:01N, 879:01:01:01-879:01:01:03, 880:01,
881:01, 882:01, 883:01[?], 884:01[?], 885:01[?], 886:01[?], 887:01[?], 888:01[?],
889:01[?], 890:01[?], 891:01[?], 892:01[?], 893:01[?], 894:01N[?], 895:01[?],
896:01, 897:01[?], 898:01[?], 899:01[?], 900:01[?], 901:01[?], 902:01[?], 903:01[?],
904:01[?], 905:01, 906:01, 907:01, 908:01[?], 909:01, 910:01[?], 911:01N[?],
912:01[?], 913:01[?], 914:01[?], 915:01[?], 916:01[?], 917:01N[?], 918:01[?],
919:01N[?], 920:01[?], 921:01[?], 922:01[?], 923:01[?], 924:01[?], 925:01N[?],
926:01[?], 927:01[?], 928:01[?], 929:01[?], 930:01[?], 931:01[?], 932:01[?],
933:01[?], 934:01Q[?], 935:01Q[?], 936:01Q[?], 937:01[?], 938:01, 939:01,
940:01[?], 941:01N[?], 942:01[?], 943:01[?], 944:01[?], 945:01[?], 947:01,
948:01, 949:01, 950:01N, 952:01, 953:01, 954:01, 955:01, 956:01,
957:01, 958:01, 959:01N, 960:01N[?], 961:01[?], 962:01[?], 963:01[?],
964:01[?], 966:01, 967:01, 968:01, 969:01[?], 970:01[?], 971:01[?], 972:01[?],
973:01[?], 974:01N[?], 975:01, 976:01[?], 977:01[?], 978:01, 979:01, 980:01[?],
981:01, 982:01[?], 983:01, 984:01N[?], 985:01N[?], 986:01N[?], 987:01[?],
988:01[?], 989:01[?], 990:01[?], 991:01[?], 992:01[?], 993:01[?], 994:01[?],
995:01N[?], 996:01[?], 997:01[?], 998:01[?], 999:01[?], 1000:01[?], 1001:01[?],
1002:01[?], 1003:01[?], 1004:01[?], 1005:01[?], 1006:01[?], 1007:01[?],
1008:01[?], 1009:01[?], 1010:01[?], 1011:01, 1012:01[?], 1013:01[?], 1014:01[?],
1015:01[?], 1016:01[?], 1017:01[?], 1018:01[?], 1019:01[?], 1020:01[?],
1021:01[?], 1022:01[?], 1023:01[?], 1024:01[?], 1025:01[?], 1026:01[?],
1027:01[?], 1028:01[?], 1029:01N[?], 1030:01[?], 1031:01[?], 1032:01,
1033:01, 1035:01, 1036:01, 1038:01Q[?], 1039:01[?], 1040:01,
1041:01N[?], 1042:01[?], 1043:01[?], 1044:01N, 1045:01N[?], 1046:01[?],
1047:01[?], 1048:01[?], 1049:01, 1051:01, 1052:01, 1053:01, 1054:01,
1055:01, 1056:01[?], 1057:01, 1059:01[?], 1060:01[?], 1062:01[?], 1063:01[?],
1064:01[?], 1065:01[?], 1066:01, 1067:01, 1068:01[?], 1069:01, 1070:01N[?],
1071:01, 1072:01, 1073:01[?], 1074:01[?], 1075:01[?], 1076:01[?], 1077:01,
1078:01[?], 1079:01N[?], 1080:01[?], 1081:01, 1082:01, 1083:01,
1084:01N, 1085:01, 1086:01, 1087:01, 1088:01, 1089:01, 1090:01,
1091:01, 1092:01Q[?], 1093:01[?], 1094:01[?], 1096:01, 1097:01,
1098:01N, 1100:01, 1102:01, 1103:01, 1104:01, 1105:01, 1106:01,
1107:01, 1108:01, 1109:01, 1110:01, 1111:01, 1112:01N[?], 1113:01,
1114:01[?], 1115:01[?], 1116:01, 1121:01N[?], 1122:01, 1123:01, 1124:01,
1125:01, 1126:01, 1127:01, 1128:01, 1129:01, 1130:01, 1131:01,
1132:01, 1134:01, 1135:01N, 1136:01, 1137:01, 1138:01[?], 1139:01,
1140:01, 1141:01[?], 1142:01, 1143:01[?], 1144:01, 1145:01, 1146:01,
1147:01, 1148:01Q, 1150:01[?], 1151:01, 1152:01, 1153:01, 1154:01N[?],
1155:01, 1157:01, 1158:01[?], 1159:01[?], 1160:01[?], 1161:01[?], 1162:01[?],
1163:01[?], 1164:01[?], 1165:01, 1166:01[?], 1167:01[?], 1168:01,
1169:01N[?], 1170:01, 1171:01, 1172:01[?], 1173:01, 1174:01[?], 1175:01[?],
1176:01[?], 1177:01, 1178:01[?], 1179:01, 1180:01[?], 1181:01[?], 1182:01[?],
1183:01[?], 1184:01[?], 1185:01, 1186:01[?], 1187:01[?], 1188:01, 1190:01[?],
1191:01N, 1192:01, 1193:01N, 1195:01[?]



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

83	160 bp	430 bp	*04:01:01:01-04:01:52, 11:01:01:01-11:01:05, 15:01:01:01-15:01:04, 24:01, 27:01, 28:01, 31:01:01:01-31:01:01:02, 33:01:01:01-33:01:01:05, 39:01:01:01-39:01:03, 40:01:01:01-40:01:02, 49:01:01:01-49:01:01:03, 51:01:01:01-51:01:01:02, 53:01, 66:01, 72:01:01:01-72:01:01:03, 74:01, 81:01:01:01-81:01:02, 85:01:01:01-85:01:02, 89:01, 96:01, 99:01, 102:01, 117:01, 120:01N, 121:01, 126:01:01:01-126:01:01:03, 139:01, 147:01, 149:01, 156:01, 158:01, 169:01, 172:01, 174:01, 175:01, 176:01, 177:01, 178:01, 179:01, 180:01, 181:01, 192:01, 193:01, 194:01, 195:01, 198:01, 199:01, 202:01, 204:01, 205:01, 206:01, 207:01, 212:01, 213:01:01-213:01:02, 214:01, 215:01, 224:01, 225:01, 227:01:03, 228:01, 229:01, 230:01, 232:01, 240:01, 247:01, 250:01, 253:01, 254:01, 255:01, 262:01, 264:01, 268:01, 272:01, 275:01, 276:01, 277:01, 279:01:01-279:01:02, 280:01, 281:01, 282:01, 283:01, 287:01, 290:01, 294:01, 296:01, 298:01, 299:01, 302:01, 303:01, 305:01, 319:01, 320:01, 323:01, 328:01N, 334:01:01-334:01:02, 335:01, 336:01, 337:01, 341:01, 345:01, 348:01:01-348:01:02, 350:01, 354:01:01-354:01:02, 355:01, 356:01, 357:01N, 360:01, 370:01, 372:01, 375:01, 376:01, 378:01, 387:01, 388:01, 396:01, 397:01, 398:01, 399:01, 415:01, 418:01, 425:01, 426:01, 427:01, 428:01, 434:01, 437:01, 450:01N, 453:01, 454:01, 455:01N, 456:01, 459:01, 464:01, 465:01, 468:01, 471:01, 474:01, 475:01, 476:01, 479:01, 480:01, 482:01, 485:01, 486:01, 487:01, 493:01, 512:01, 517:01, 520:01, 521:01, 522:01, 523:01:01-523:01:02, 524:01, 529:01, 531:01, 534:01, 553:01, 554:01, 555:01, 556:01, 559:01, 561:01, 566:01, 569:01, 574:01, 575:01, 578:01, 583:01, 584:01:01:01-584:01:02:02, 585:01:01:01-585:01:01:02, 591:01, 592:01, 593:01, 607:01, 614:01, 615:01, 618:01, 623:01, 625:01, 633:01, 642:01, 643:01, 644:01, 649:01, 654:01, 658:01, 663:01, 666:01, 670:01, 672:01, 677:01, 679:01, 682:01, 683:01, 686:01, 687:01, 694:01, 696:01N, 699:01, 702:01, 707:01, 713:01, 716:01, 722:01, 732:01N, 736:01, 739:01, 742:01, 743:01N, 745:01, 747:01, 748:01N, 749:01, 750:01, 752:01, 753:01, 754:01N, 755:01, 756:01N, 757:01, 758:01, 761:01, 765:01, 767:01, 768:01, 769:01, 772:01, 773:01, 784:01, 788:01, 795:01, 803:01, 804:01, 806:01:01:01-806:01:01:02, 808:01, 811:01, 812:01, 813:01, 814:01, 820:01, 821:01N, 824:01:01-824:01:02, 826:01, 828:01, 837:01, 838:01N, 840:01, 842:01, 844:01N, 849:01, 850:01, 852:01, 856:01, 859:01, 867:01N, 868:01N, 871:01N, 874:01N, 875:01N, 878:01N, 880:01, 882:01, 895:01, 896:01, 900:01, 901:01, 903:01, 907:01, 910:01, 911:01N, 915:01, 918:01, 921:01, 922:01, 926:01, 930:01, 932:01, 937:01, 939:01, 943:01, 945:01, 957:01, 966:01, 974:01N, 978:01, 980:01, 988:01, 989:01, 992:01, 993:01, 997:01, 1001:01, 1002:01, 1003:01, 1004:01, 1010:01, 1011:01, 1016:01, 1023:01, 1033:01, 1034:01, 1039:01, 1040:01, 1042:01, 1046:01, 1048:01, 1054:01, 1060:01, 1062:01, 1063:01, 1064:01, 1066:01, 1074:01, 1079:01N, 1080:01, 1081:01, 1086:01, 1091:01, 1097:01, 1100:01, 1108:01, 1109:01, 1112:01N, 1113:01, 1121:01N, 1129:01, 1132:01, 1135:01N, 1136:01, 1138:01, 1139:01, 1144:01, 1146:01, 1148:01Q, 1152:01, 1155:01, 1161:01, 1164:01, 1167:01, 1173:01, 1177:01, 1181:01, 1184:01, 1192:01, 1193:01N, 1195:01
	215 bp		*169:01
84⁷	200 bp	430 bp	*180:01, 418:01, 432:01, 604:01



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

	260 bp		*171:01, 177:01, 186:01, 338:01, 787:01
85^{4,7}	110 bp	515 bp	*572:01, 923:01, 1021:01
	195 bp		*498:01
	280 bp		*585:01:01:01-585:01:01:02, 1119:01
86^{4,7}	125 bp	430 bp	*190:01, 796:01
	170 bp		*485:01
	195 bp		*354:01:01-354:01:02, 375:01, 480:01, 531:01, 663:01, 693:01N, 1039:01
	220 bp		*360:01, 520:01
87	145 bp	430 bp	*327:01, 619:01, 698:01, 815:01
	230 bp		*329:01, 587:01, 929:01
88⁴	80 bp	430 bp	*15:01:01:01-15:01:04, 18:01:01:01-18:01:03, 28:01, 40:01:01:01-40:01:02, 53:01, 62:01, 74:01, 139:01, 198:01, 290:01, 292:01, 296:01, 299:01, 333:01, 347:01, 387:01, 471:01, 482:01, 499:01, 512:01, 526:01, 585:01:01:01-585:01:01:02, 644:01, 695:01, 745:01, 752:01, 821:01N, 896:01, 897:01, 900:01, 903:01, 910:01, 918:01, 942:01, 1054:01, 1109:01, 1136:01, 1142:01, 1165:01, 1192:01
89⁷	245 bp	430 bp	*187:01, 357:01N, 401:01N
90^{4,6}	80 bp	515 bp	*571:01
	125 bp		*700:01N
	215 bp		*415:01, 534:01
91	155 bp	515 bp	*464:01, 535:01, 617:01, 699:01
	220 bp		*534:01
92⁴	115 bp	430 bp	*403:01N
	165 bp		*445:01
	195 bp		*691:01N
	230 bp		*533:01
93⁷	185 bp	430 bp	*634:01
	225 bp		*404:01, 598:01N, 661:01N
	255 bp		*450:01N
94	160 bp	430 bp	*01:01:01:01-01:01:04, 01:01:06-01:01:07, 26:01:01-26:01:03, 56:01, 65:01:01-65:01:02, 76:01, 90:01:01-90:01:02, 92:01, 118:01, 122:01, 125:01, 127:01, 142:01, 150:01, 159:01N, 162:01:01-162:01:02, 203:01:01-203:01:02, 220:01, 267:01, 288:01, 291:01, 306:01, 315:01, 326:01, 327:01, 346:01, 363:01, 377:01, 385:01, 403:01N, 411:01, 417:01:01-417:01:02, 436:01, 440:01, 442:01, 462:01, 470:01, 481:01, 483:01, 490:01, 499:01, 514:01, 516:01, 536:01, 542:01, 543:01, 562:01, 609:01, 616:01, 626:01, 631:01, 632:01, 636:01, 661:01N, 667:01, 708:01, 709:01, 710:01, 711:01, 733:01, 800:01N, 807:01, 810:01, 825:01, 853:01, 862:01N, 866:01N, 870:01N, 879:01:01:01-879:01:01:03, 904:01, 916:01, 931:01, 953:01, 960:01N, 969:01, 976:01, 979:01, 984:01N, 998:01, 999:01, 1012:01, 1024:01, 1038:01Q, 1050:01, 1068:01, 1076:01, 1088:01, 1137:01, 1140:01, 1141:01, 1147:01, 1150:01, 1151:01, 1162:01, 1166:01, 1169:01N, 1183:01
	190 bp		*11:01:02 [?] , 53:01 [?] , 56:01 [?] , 74:01 [?] , 102:01 [?] , 454:01, 486:01, 722:01
	240 bp		*08:01, 57:01, 65:01:01-65:01:02, 68:01, 73:01, 75:01, 84:01, 150:01, 188:01, 234:01, 285:01, 306:01, 377:01, 442:01, 499:01, 620:01, 626:01, 648:01:01:01-648:01:01:02, 904:01, 916:01, 976:01, 1140:01, 1150:01
	270 bp		



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

			*08:01 [?] , 53:01 [?] , 60:01 [?] , 62:01 [?] , 73:01 [?] , 75:01 [?] , 102:01 [?] , 237:01, 454:01, 486:01, 722:01, 775:01
95	140 bp	430 bp	*13:01:01:01-13:01:04, 107:01, 117:01, 133:01, 160:01, 172:01, 173:01, 198:01 ^w , 206:01, 246:01, 247:01, 269:01, 287:01, 362:01, 389:01 ^w , 393:01, 438:01, 441:01, 449:01, 474:01, 503:01, 514:01, 518:01, 519:01, 544:01, 561:01, 564:01, 565:01, 597:01, 600:01, 609:01, 612:01, 631:01, 633:01, 634:01, 663:01 ^w , 703:01, 711:01, 716:01, 822:01, 873:01N, 888:01, 893:01, 924:01, 947:01, 991:01, 996:01, 1008:01, 1045:01N, 1057:01, 1065:01, 1069:01, 1105:01, 1123:01, 1131:01, 1145:01, 1185:01
	220 bp		*06:01:01:01-06:01:05, 13:01:01:01-13:01:04, 21:01, 29:01, 37:01, 44:01, 88:01, 93:01, 107:01, 109:01, 110:01, 117:01, 133:01, 160:01, 161:01N, 198:01, 200:01, 206:01, 208:01, 246:01, 247:01, 249:01, 269:01, 295:01, 347:01, 349:01, 362:01, 371:01, 389:01, 393:01, 413:01, 438:01, 448:01, 449:01, 518:01, 519:01, 544:01, 546:01, 550:01, 561:01, 563:01, 564:01, 565:01, 566:01, 567:01, 597:01, 600:01, 601:01, 608:01, 609:01, 612:01, 631:01, 633:01, 634:01, 662:01, 663:01, 703:01, 711:01, 717:01, 737:01, 746:01, 822:01, 873:01N, 888:01, 893:01, 906:01, 909:01, 914:01, 920:01, 924:01, 947:01, 991:01, 996:01, 1006:01, 1008:01, 1019:01, 1022:01, 1045:01N, 1047:01, 1057:01, 1065:01, 1069:01, 1087:01, 1105:01, 1111:01, 1123:01, 1127:01, 1131:01, 1145:01, 1170:01, 1185:01, 1186:01, 1190:01
96¹¹	-	-	Negative Control

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of DPB1 SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³For several DPB1 alleles 1st and/or 3rd exon(s) and beyond, as well as intron nucleotide sequences, are not available and thus the primer pairs in primer mixes 57, 58, 59, 60, 62 and 82 may amplify other DPB1 alleles. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁴HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: 8L1

Lot-specific information

⁵Primer mixes 9, 27, 43, 47, 48 and 58 may give rise to a lower yield of HLA-specific PCR product than the other DPB1 primer mixes.

⁶Primer mixes 3, 12, 26, 37, 41, 80 and 90 have a tendency to give rise to primer oligomer formation.

⁷Primer mixes 23, 25, 32, 38, 41, 42, 45, 49, 50, 62, 68, 74, 84, 85, 86, 89 and 93 may have tendencies of unspecific amplifications.

⁸Primer mix 53 has a tendency of primer oligomer formation and additionally an intense primer cloud due to the high number of primers present in the primer mix.

⁹In primer mix 53 the positive control band may be weaker than for other DPB1 primer mixes.

Primer mix 53 may weakly amplify DPB1*104:01 alleles.

¹⁰In primer mix 48 the specific PCR product of 280 base pairs may give rise to a lower yield of HLA-specific PCR product than the other DPB1 primer mixes. Optimized gel electrophoresis run time may need to be considered to ensure separation from the control band.

¹¹Primer mix 96 contains a negative control, which will amplify the majority of the HLA amplicons as well as the amplicons generated by the control primer pairs matching the human growth hormone gene. HLA-specific PCR product sizes range from 75 to 200 base pairs and the PCR product generated by the HGH positive control primer pair is 200 base pairs.

Abbreviations

w: might be weakly amplified.

?: nucleotide sequence information not available for the primer matching sequence.



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information
PRIMER SPECIFICATION

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec. PCR product	245	220	140	140	235	165	160	210	165	220	145	80
Length of int. pos. control ¹	515	515	430	430	515	430	430	515	430	515	515	515
5'-primer(s) ²	9(112) 5'-TgT 3' 5'	9(112) 5'-TgT 3' 5'	21(150) 5'-ACA 3' 5'	36(194) 5'-CgC 3' 5'	11(119) 5'-ggg 3' 5'	35(191) 5'-gTA 3' 5'	36(194) 5'-CgC 3' 5'	9(113) 5'-TTT 3' 5'	35(190) 5'-AgC 3' 5'	9(113) 5'-TTT 3' 5'	35(191) 5'-gTT 3' 5'	55(252) 5'-gAT 3'
	9(112) 5'-TgT 3' 5'	9(112) 5'-TgT 3' 5'		36(194) 5'-CgC 3' 5'	11(119) 5'-ggg 3' 5'	35(191) 5'-gTA 3' 5'	36(194) 5'-CgC 3' 5'		36(194) 5'-CgT 3'			55(252) 5'-gAT 3'
												55(252) 5'-gAT 3'
3'-primer(s) ³	76(315) 5'-CAT 3' 5'	69(292) 5'-CTT 3' 5'	55(251) 5'-CAg 3' 5'	69(292) 5'-CTT 3' 5'	76(313) 5'-CAC 3' 5'	76(313) 5'-CAC 3' 5'	76(313) 5'-CAC 3' 5'	65(280) 5'-TAT 3' 5'	76(313) 5'-CAT 3' 5'	69(292) 5'-CTC 3' 5'	69(293) 5'-gCC 3' 5'	69(292) 5'-CTC 3'
				69(292) 5'-CTT 3' 5'	76(315) 5'-CAT 3' 5'	76(315) 5'-CAT 3' 5'	76(315) 5'-CAT 3' 5'	65(280) 5'-gAT 3' 5'	76(313) 5'-TAT 3'		84(337) 5'-gTT 3'	
					76(315) 5'-CAg 3'			65(280) 5'-AAT 3' 5'	76(315) 5'-CAg 3'		84(338) 5'-CgC 3'	
								65(281) 5'-ggT 3'				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Well No.	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec. PCR product	180	105	200	80	260	180	200	235	100	185	145	195
		175		165							210	265
				210								
Length of int. pos. control ¹	430	430	430	430	515	430	430	515	430	430	430	515
5'-primer(s) ²	9(113) 5'-TTT 3' 5'	11(119) 5'-gTT 3' 5'	28(169) 5'-gAg 3' 5'	12(121) 5'-gAT 3' 5'	11(119) 5'-gTT 3' 5'	11(119) 5'-gTT 3' 5'	11(119) 5'-gTT 3' 5'	11(119) 5'-gTT 3' 5'	57(258) 5'-gAC 3' 5'	36(194) 5'-CgC 3' 5'	12(121) 5'-gAT 3' 5'	9(113) 5'-TTT 3' 5'
		11(119) 5'-gTT 3' 5'	35(190) 5'-AgC 3' 5'	28(169) 5'-gAg 3' 5'	11(119) 5'-gTT 3' 5'	11(119) 5'-gTT 3' 5'	11(119) 5'-gTT 3' 5'	11(119) 5'-gTT 3' 5'		36(194) 5'-CgC 3' 5'	12(122) 5'-ACA 3'	
			35(190) 5'-AgC 3'	56(254) 5'-TgA 3'							35(190) 5'-AgC 3'	
3'-primer(s) ³	55(252) 5'-gCC 3' 5'	55(252) 5'-gCC 3' 5'	84(338) 5'-CgC 3' 5'	69(292) 5'-CTT 3' 5'	84(338) 5'-CgT 3' 5'	57(258) 5'-TAg 3' 5'	65(280) 5'-gAg 3' 5'	76(315) 5'-CAT 3' 5'	76(315) 5'-CAT 3' 5'	84(338) 5'-CgC 3' 5'	69(292) 5'-CTT 3' 5'	60(265) 5'-gTg 3'
				69(292) 5'-CTT 3'							69(292) 5'-CTT 3' 5'	84(338) 5'-CgT 3'
Well No.	13	14	15	16	17	18	19	20	21	22	23	24



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

Well No.	25	26	27	28	29	30	31	32	33	34	35	36
Length of spec.	200	75	95	200	105	185	120	75	105	210	195	260
PCR product	265	205	215	260	160		215	200	145			
			240		220			260	195			
Length of int. pos. control ¹	515	515	430	430	515	430	515	430	430	430	430	515
5'-primer(s) ²	9(112)	15(130)	10(116)	11(119)	9(112)	9(112)	26(164)	9(113)	15(130)	9(112)	33(184)	11(119)
	5' -Tgg 3'	5' -AAG 3'	5' -CCg 3'	5' -gTT 3'	5' -TgC 3'	5' -TgC 3'	5' -ggT 3'	5' -TTC 3'	5' -AAG 3'	5' -TgC 3'	5' -ggC 3'	5' -ggg 3'
	28(169)	57(258)	11(119)	11(119)	107(406)		57(258)	14(127)	22(151)		33(184)	11(119)
	5' -gAg 3'	5' -gAC 3'	5' -ggg 3'	5' -gTT 3'	5' -CCT 3'		5' -gAg 3'	5' -AgC 3'	5' -CAT 3'		5' -ggC 3'	5' -ggg 3'
	35(190)	60(265)	16(133)					32(181)	36(194)			11(119)
	5' -AgC 3'	5' -ggC 3'	5' -gCC 3'					5' -ACT 3'	5' -CgC 3'			5' -ggg 3'
			18(139)					33(184)	36(194)			
			5' -Cgg 3'					5' -ggC 3'	5' -CgC 3'			
			57(258)					33(184)				
			5' -gAC 3'					5' -ggC 3'				
								73(305)				
								5' -gCT 3'				
3'-primer(s) ³	84(338)	69(292)	76(313)	65(280)	69(292)	57(258)	84(338)	84(338)	57(258)	65(280)	84(338)	84(338)
	5' -CgT 3'	5' -CTC 3'	5' -CAT 3'	5' -gAT 3'	5' -CTC 3'	5' -TAG 3'	5' -CgT 3'	5' -CgT 3'	5' -TAG 3'	5' -gAg 3'	5' -CgA 3'	5' -CgA 3'
			76(313)	84(337)	128(470)		84(338)		69(292)			
			5' -TAT 3'	5' -gTT 3'	5' -ATT 3'		5' -CgT 3'		5' -CTC 3'			
					147(527)				85(341)			
					5' -TAT 3'				5' -gCC 3'			
Well No.	25	26	27	28	29	30	31	32	33	34	35	36



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

Well No.	37	38	39	40	41	42	43	44	45	46	47	48
Length of spec.	160	95	210	135	75	140	200	60	120	75	165	115
PCR product	185	170 200		180 230	225 250 280	185 210 245 285		115 150 185 240 270	160 190	210	215	210 280
Length of int. pos. control ¹	430	515	515	430	515	430	430	430	430	430	430	430
5'-primer(s) ²	36(194) 5'-CgC 3'	35(190) 5'-AgC 3'	9(113) 5'-TTT 3' 5'	9(112) 5'-TgC 3' 5'	1 st I 5'-CAT 3' 5'	9(113) 5'-TTT 3' 5'	35(190) 5'-AgC 3' 5'	7(108) 5'-TAA 3' 5'	17(136) 5'-ACA 3' 5'	9(113) 5'-TTT 3' 5'	35(190) 5'-AgC 3' 5'	11(119) 5'-gTT 3'
	36(194) 5'-CgC 3'	36(194) 5'-CgT 3'		37(197) 5'-gCC 3'				10(115) 5'-TCT 3' 5'	17(137) 5'-CgT 3' 5'	56(254) 5'-TgA 3'		11(119) 5'-gTT 3'
		39(202) 5'-TCT 3'						17(136) 5'-ACC 3' 5'	35(191) 5'-gTA 3'			
								18(140) 5'-gTg 3' 5'	35(191) 5'-gTA 3'			
								28(169) 5'-gAg 3' 5'	45(220) 5'-Agg 3'			
								32(182) 5'-CCC 3' 5'	56(254) 5'-TgT 3'			
								32(182) 5'-CCA 3' 5'	57(258) 5'-gAC 3'			
								49(233) 5'-gAT 3'				
3'-primer(s) ³	76(313) 5'-TAT 3'	55(251) 5'-CAg 3' 5'	65(280) 5'-gAg 3' 5'	55(252) 5'-gCC 3' 5'	12(122) 5'-gCA 3' 5'	43(214) 5'-CCA 3' 5'	84(338) 5'-CgA 3' 5'	55(251) 5'-CAg 3' 5'	84(338) 5'-CgC 3' 5'	59(264) 5'-TTg 3' 5'	76(315) 5'-CAT 3' 5'	35(191) 5'-CgA 3'
	81(328) 5'-gTC 3'	79(324) 5'-TTC 3'		69(292) 5'-CTT 3' 5'	62(271) 5'-CTA 3' 5'	57(257) 5'-ACA 3' 5'	91(359) 5'-ggT 3' 5'	55(252) 5'-gCC 3'		64(277) 5'-gTg 3' 5'	92(362) 5'-CTg 3' 5'	67(286) 5'-CTA 3'
	84(338) 5'-CgA 3'	80(326) 5'-AgC 3'		76(313) 5'-TAT 3' 5'	67(286) 5'-CTT 3' 5'	63(275) 5'-CCg 3'		84(338) 5'-CgT 3'		65(280) 5'-gAA 3'		91(359) 5'-ggT 3'
		81(330) 5'-TCT 3'			70(295) 5'-CCA 3' 5'	67(286) 5'-CTT 3'				65(281) 5'-ggT 3'		
		88(350) 5'-ggA 3'			72(301) 5'-CAA 3' 5'	77(316) 5'-gCg 3'				69(293) 5'-gCC 3'		
		93(364) 5'-T 3'			78(320) 5'-gTg 3' 5'	91(359) 5'-ggT 3'						
					84(337) 5'-gCT 3'							
Well No.	37	38	39	40	41	42	43	44	45	46	47	48



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

Well No.	49	50	51	52	53	54	55	56	57	58	59	60
Length of spec.	95	160	170	130	90	115	135	145	255	95	80	255
PCR product	130	210	225	160	145	160	165	180				
	205			215	175	200	190	230				
					200			290				
Length of int. pos. control ¹	430	430	430	430	430	430	430	430	430	430	430	430
5'-primer(s) ²	8(109)	11(119)	8(110)	7(107)	36(194)	36(194)	17(136)	9(113)	107(406)	-22(24)	191(660)	107(406)
	5'-ACA 3'	5'-ggg 3'	5'-CCg 3'	5'-TTg 3'	5'-CgT 3'	5'-CgC 3'	5'-ACA 3'	5'-TTT 3'	5'-CCT 3'	5'-gCA 3'	5'-gAT 3'	5'-CCT 3'
	9(114)	23(155)	10(115)	8(110)		36(194)	22(151)					
	5'-TTA 3'	5'-gCA 3'	5'-TCT 3'	5'-CCg 3'		5'-CgC 3'	5'-CAT 3'					
	33(184)	29(173)	26(163)	10(115)		57(258)	27(166)					
	5'-ggC 3'	5'-CAG 3'	5'-TgC 3'	5'-TCT 3'		5'-gAC 3'	5'-AgT 3'					
	37(197)		27(166)	14(129)			30(175)					
	5'-gCA 3'		5'-AgT 3'	5'-gAC 3'			5'-TCC 3'					
	47(226)			30(175)			37(197)					
	5'-ggA 3'			5'-TCC 3'			5'-gCT 3'					
				39(204)								
				5'-gAA 3'								
3'-primer(s) ³	65(280)	65(280)	69(292)	69(292)	53(245)	79(323)	69(292)	44(217)	178(619)	-5(76)	205(700)	178(619)
	5'-gAT 3'	5'-gAT 3'	5'-CTC 3'	5'-CTC 3'	5'-gCT 3'	5'-TgC 3'	5'-CTC 3'	5'-CTT 3'	5'-CAT 3'	5'-CAC 3'	5'-CAT 3'	5'-CAg 3'
		69(292)			68(289)	81(328)		45(220)				178(619)
		5'-CTC 3'			5'-CTT 3'	5'-gTC 3'		5'-gAT 3'				5'-CAA 3'
					74(308)	85(340)		55(252)				
					5'-TgC 3'	5'-CTg 3'		5'-TCC 3'				
					75(310)	85(340)		71(299)				
					5'-CCC 3'	5'-CTA 3'		5'-CTT 3'				
					78(319)	85(341)		73(305)				
					5'-TCA 3'	5'-gCT 3'		5'-CCC 3'				
					79(322)	92(361)		92(362)				
					5'-gTA 3'	5'-TCA 3'		5'-CTT 3'				
					85(340)	2 nd I						
					5'-CCT 3'	5'-C 3'						
					85(341)							
					5'-CCC 3'							
					90(355)							
					5'-CTC 3'							
					91(358)							
					5'-gCC 3'							
Well No.	49	50	51	52	53	54	55	56	57	58	59	60



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

Well No.	61	62	63	64	65	66	67	68	69	70	71	72
Length of spec.	145	260	170	150	170	185	90	175	180	75	140	145
PCR product	225		210	180	205	235	155	200	215	115	180	175
			245	215	240	280	180	265	245	170	260	250
			275				210			215		
										260		
Length of int. pos. control ¹	430	515	515	430	430	515	430	430	515	430	430	430
5'-primer(s) ²	116(433) 5'-gCT 3'	107(406) 5'-CCC 3'	9(113) 5'-TTT 3'	9(113) 5'-TTT 3'	9(113) 5'-TTT 3'	9(113) 5'-TTT 3'	36(194) 5'-CgT 3'	8(110) 5'-CCC 3'	12(122) 5'-ACC 3'	12(121) 5'-gAT 3'	9(112) 5'-TgT 3'	12(122) 5'-ACC 3'
	144(517) 5'-CCT 3'							9(112) 5'-TgT 3'	17(137) 5'-CgT 3'	23(155) 5'-gCC 3'		14(127) 5'-AgC 3'
								31(179) 5'-CAG 3'	26(165) 5'-gAC 3'	29(172) 5'-ACg 3'		15(130) 5'-AAC 3'
								32(181) 5'-ACT 3'	35(192) 5'-TTA 3'	37(197) 5'-gCA 3'		17(136) 5'-ACA 3'
								37(197) 5'-gCA 3'	39(202) 5'-TCC 3'	41(210) 5'-gAA 3'		39(202) 5'-TCC 3'
								41(208) 5'-gCC 3'	41(208) 5'-gCC 3'	45(220) 5'-Agg 3'		49(232) 5'-Tgg 3'
								41(210) 5'-gAA 3'		46(223) 5'-TCT 3'		
										59(264) 5'-TgA 3'		
										73(305) 5'-gCT 3'		
3'-primer(s) ³	178(619) 5'-CAg 3'	179(624) 5'-CTg 3'	52(241) 5'-CCT 3'	44(217) 5'-CTT 3'	48(230) 5'-TCg 3'	53(245) 5'-gCT 3'	52(241) 5'-CCT 3'	84(338) 5'-CgC 3'	84(338) 5'-CgC 3'	84(338) 5'-CgC 3'	42(212) 5'-CCC 3'	84(338) 5'-CgT 3'
			66(284) 5'-CCg 3'	47(227) 5'-CCT 3'	51(238) 5'-CAC 3'	60(265) 5'-gTC 3'	74(308) 5'-TgC 3'				44(217) 5'-CTg 3'	
			66(284) 5'-CCT 3'	52(241) 5'-CCT 3'	53(245) 5'-gCg 3'	62(271) 5'-CTg 3'	78(319) 5'-TCA 3'				55(251) 5'-CAA 3'	
			78(319) 5'-TCA 3'	60(266) 5'-TgC 3'	64(278) 5'-TgC 3'	73(305) 5'-CCC 3'	79(323) 5'-TgA 3'				82(331) 5'-CTg 3'	
			85(341) 5'-gCg 3'	67(287) 5'-CCg 3'	75(312) 5'-ATA 3'	88(350) 5'-ggT 3'	81(330) 5'-TCT 3'					
			88(350) 5'-ggT 3'			90(355) 5'-CTA 3'	83(335) 5'-CCg 3'					
							93(364) 5'-T 3'					
Well No.	61	62	63	64	65	66	67	68	69	70	71	72



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

Well No.	73	74	75	76	77	78	79	80	81	82	83	84
Length of spec.	170	115	140	120	120	145	75	160	65	80	160	200
PCR product	190	160		180	170	180	115	250	95		215	260
		215		220	220	220	170		145			
				250	250	250	245		180			
							270		205			
									270			
Length of int.	430	430	430	430	430	430	515	430	430	430	430	430
pos. control ¹												
5'-primer(s) ²	35(191) 5'-gTT 3'	16(133) 5'-gCC 3'	36(194) 5'-CgT 3'	14(127) 5'-AgA 3'	11(119) 5'-gTT 3'	12(122) 5'-ACA 3'	8(110) 5'-CCC 3'	14(127) 5'-AgA 3'	9(112) 5'-TgT 3'	191(660) 5'-gAT 3'	18(139) 5'-Cgg 3'	9(113) 5'-TTT 3'
		20(146) 5'-TgC 3'	36(194) 5'-CgT 3'	23(155) 5'-gCC 3'	16(135) 5'-AC 3'	16(135) 5'-AC 3'	12(122) 5'-ACA 3'	39(202) 5'-TCC 3'	35(191) 5'-gTT 3'		36(194) 5'-CgC 3'	
		33(184) 5'-ggA 3'	36(194) 5'-CgT 3'	25(160) 5'-TCg 3'	22(152) 5'-ACT 3'	22(152) 5'-ACg 3'	20(145) 5'-ATA 3'	40(207) 5'-Agg 3'			36(194) 5'-CgC 3'	
		34(187) 5'-AgA 3'		36(193) 5'-TCT 3'	24(159) 5'-TTg 3'	26(163) 5'-TgA 3'	37(197) 5'-gCC 3'	42(211) 5'-ACA 3'				
		37(197) 5'-gCT 3'		36(193) 5'-TCC 3'	42(211) 5'-ACA 3'	36(193) 5'-TCA 3'	41(208) 5'-gCT 3'	42(212) 5'-CgC 3'				
		51(238) 5'-Agg 3'		40(207) 5'-AgA 3'	57(258) 5'-gAC 3'	40(207) 5'-AgA 3'	42(212) 5'-CgC 3'	45(220) 5'-AgC 3'				
				41(208) 5'-gCT 3'		41(208) 5'-gCA 3'	59(264) 5'-TgA 3'	46(223) 5'-TCT 3'				
				57(258) 5'-gAC 3'		49(233) 5'-gAT 3'	73(305) 5'-gCT 3'	49(233) 5'-gAT 3'				
3'-primer(s) ³	78(321) 5'-TgA 3'	76(313) 5'-CAT 3'	69(292) 5'-CTT 3'	84(338) 5'-CgC 3'	84(338) 5'-CgC 3'	84(338) 5'-CgC 3'	84(338) 5'-CgC 3'	84(338) 5'-CgC 3'	44(217) 5'-CTg 3'	205(700) 5'-CAC 3'	76(313) 5'-CAT 3'	60(265) 5'-gTC 3'
	84(338) 5'-CgA 3'								53(245) 5'-gCT 3'			63(275) 5'-CCg 3'
	85(340) 5'-CCg 3'								82(331) 5'-CTg 3'			64(278) 5'-TgC 3'
	85(340) 5'-CCA 3'								85(340) 5'-CCA 3'			80(326) 5'-AgC 3'
	85(341) 5'-gCC 3'								90(355) 5'-CTT 3'			82(331) 5'-CTT 3'
	85(341) 5'-gCg 3'											83(334) 5'-CAC 3'
												84(337) 5'-gCT 3'
Well No.	73	74	75	76	77	78	79	80	81	82	83	84



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

Well No.	85	86	87	88	89	90	91	92	93	94	95
Length of spec.	110	125	145	80	245	80	155	115	185	160	140
PCR product	195	170	230			125	220	165	225	190	220
	280	195				215		195	255	240	
		220						230		270	
Length of int.	515	430	430	430	430	515	515	430	430	430	430
pos. control ¹											
5'-primer(s) ²	107(406) 5'-CCC 3'	22(152) 5'-ACT 3'	21(149) 5'-gAT 3'	72(301) 5'-CAg 3'	11(119) 5'-ggg 3'	107(406) 5'-CCT 3'	107(406) 5'-CCT 3'	33(184) 5'-ggT 3'	13(124) 5'-ggT 3'	9(113) 5'-TTT 3'	9(112) 5'-TgT 3'
		26(165) 5'-gAC 3'	25(161) 5'-CCC 3'		11(119) 5'-ggg 3'			43(214) 5'-TgC 3'	22(151) 5'-CAT 3'	36(194) 5'-CgC 3'	35(191) 5'-gTA 3'
		33(184) 5'-ggC 3'	49(232) 5'-Tgg 3'					59(264) 5'-TgA 3'	25(161) 5'-CCA 3'	36(194) 5'-CgC 3'	
		33(186) 5'-Ag 3'	49(233) 5'-gAT 3'					117(436) 5'-ACA 3'	27(166) 5'-AAT 3'		
		41(208) 5'-gCA 3'							36(194) 5'-CgT 3'		
		56(254) 5'-TgA 3'									
3'-primer(s) ³	129(474) 5'-AAg 3'	84(338) 5'-CgC 3'	84(338) 5'-CgT 3'	84(338) 5'-CgA 3'	77(317) 5'-TgA 3'	119(443) 5'-AAC 3'	141(508) 5'-CAT 3'	84(338) 5'-CgT 3'	84(338) 5'-CgT 3'	76(313) 5'-TAC 3'	69(292) 5'-CTC 3'
	158(561) 5'-AgT 3'				81(330) 5'-TCC 3'	135(490) 5'-CTA 3'	141(509) 5'-ACC 3'	179(624) 5'-CTg 3'		83(334) 5'-CAC 3'	
	187(646) 5'-CTC 3'					162(572) 5'-ggT 3'	147(527) 5'-TAT 3'			85(341) 5'-gCT 3'	
						167(586) 5'-ATT 3'	167(586) 5'-ATT 3'			92(362) 5'-CTT 3'	
Well No.	85	86	87	88	89	90	91	92	93	94	95

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

CELL LINE VALIDATION SHEET																				
DPB1 SSP kit ²																				
				Well ³																
				1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
				Prod. No.:	202127001	202127002	202127003	202127004	202127005	202127006	202127007	202127008	202127009	202127010	202127011	202127012	202127013	202127014	202127015	202127016
IHCW cell line ¹		DPB1																		
1	9001 SA	*04:02		-	-	-	-	-	-	-	-	+	+	-	+	-	-	-	-	+
2	9280 LK707	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
3	9011 E4181324	*02:01	*04:01	-	-	+	+	-	-	-	-	+	-	+	+	-	-	-	-	-
4	9275 GU373	*04:02	*17:01	-	-	-	-	-	-	-	-	+	-	+	+	-	-	-	-	+
5	9009 KAS011	*04:01	*14:01	-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	+
6	9353 SM	*05:01		-	-	-	-	-	-	-	-	+	+	-	-	+	+	-	-	-
7	9020 QBL	*02:02		-	-	-	-	-	-	-	-	+	+	+	-	+	+	+	-	-
8	9025 DEU	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
9	9026 YAR	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
10	9107 LKT3	*05:01		-	-	-	-	-	-	-	-	+	+	-	-	+	+	-	-	-
11	9051 PTOU	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
12	9052 DBB	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
13	9004 JESTHOM	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
14	9071 OLGA	*03:01	*04:02	+	+	-	-	-	-	-	-	+	+	-	+	-	-	-	-	+
15	9075 DKB	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
16	9037 SWEIG007	*04:02		-	-	-	-	-	-	-	-	+	+	-	+	-	-	-	-	+
17	9282 CTM3953540	*01:01	*02:01	+	+	+	+	+	+	+	+	+	+	+	+	-	-	-	-	-
18	9257 32367	*18:01	*105:01	-	+	-	-	-	-	-	-	+	-	+	-	-	-	-	-	+
19	9038 BM16	*02:01		-	-	-	-	-	-	-	-	+	+	+	+	-	-	-	-	-
20	9059 SLE005	*03:01		+	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
21	9064 AMALA	*04:02		-	-	-	-	-	-	-	-	+	+	-	+	-	-	-	-	+
22	9056 KOSE	*02:01	*13:01	+	-	+	-	-	+	+	+	+	+	+	+	-	-	-	-	-
23	9124 IHL	*02:01		-	-	-	-	-	-	-	-	+	+	+	+	-	-	-	-	-
24	9035 JBUSH	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
25	9049 IBW9	*01:01		+	+	+	+	+	+	+	+	-	-	-	-	-	-	-	-	-
26	9285 WT49	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
27	9191 CH1007	*03:01	*13:01	+	+	+	-	-	+	+	-	-	-	-	-	-	-	-	-	+
28	9320 BEL5GB	*02:01	*04:01	-	-	+	+	-	-	-	-	+	+	+	+	-	-	-	-	-
29	9050 MOU	*02:01		-	-	-	-	-	-	-	-	+	+	+	+	-	-	-	-	-
30	9021 RSH	*01:01	*105:01	+	+	+	+	+	+	+	+	+	+	-	+	-	-	-	-	+
31	9019 DUCAF	*02:02		-	-	-	-	-	-	-	-	+	+	+	-	+	+	+	-	-
32	9297 HAG	*02:01	*03:01	+	+	-	-	-	-	-	-	+	+	+	+	-	-	-	-	+
33	9098 MT14B	*04:02		-	-	-	-	-	-	-	-	+	-	+	-	-	-	-	-	+
34	9104 DHIF	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
35	9302 SSTO	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
36	9024 KT17	*05:01		-	-	-	-	-	-	-	-	+	-	-	-	+	+	-	-	-
37	9065 HHKB	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
38	9099 LZL	*04:02		-	-	-	-	-	-	-	-	+	+	-	+	-	-	-	-	+
39	9315 CML	*04:01	*14:01	-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	+
40	9134 WHONP199	*05:01		-	-	-	-	-	-	-	-	+	+	-	-	+	+	-	-	-
41	9055 H0301	*05:01		-	-	-	-	-	-	-	-	+	+	-	-	+	+	-	-	-
42	9066 TAB089	*02:02		-	-	-	-	-	-	-	-	+	+	+	-	+	+	+	-	-
43	9076 T7526	*107:01		+	-	+	-	-	+	+	-	-	-	-	-	-	-	-	-	-
44	9057 TEM	*04:01		-	-	+	+	-	-	-	-	+	-	-	+	-	-	-	-	-
45	9239 SHJO	*04:02	*104:01	+	+	-	-	-	-	-	-	+	+	-	+	-	-	-	-	+
46	9013 SCHU	*04:02		-	-	-	-	-	-	-	-	+	-	+	-	-	-	-	-	+
47	9045 TUBO	*02:01	*04:02	-	-	-	-	-	-	-	-	+	+	+	+	-	-	-	-	+
48	9303 TER-ND	*02:01	*04:01	-	-	+	+	-	-	-	-	+	+	+	+	-	-	-	-	-



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

CELL LINE VALIDATION SHEET																				
DPB1 SSP kit ²																				
				Well ³																
				17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
				Prod. No.:	202127017	202127018	202127019	202127020	202127021	202127022	202127023	202127024	202127025	202127026	202127027	202127028	202127029	202127030	202127031	202127032
	IHWC cell line ¹	DPB1																		
1	9001 SA	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
2	9280 LK707	*04:01		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
3	9011 E4181324	*02:01	*04:01	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
4	9275 GU373	*04:02	*17:01	+	+	-	-	-	-	-	-	-	+	+	+	+	-	-	-	
5	9009 KAS011	*04:01	*14:01	+	+	+	+	+	+	-	-	-	-	-	-	-	+	-	-	
6	9353 SM	*05:01		-	-	-	-	-	-	-	+	+	+	-	-	-	-	+	-	
7	9020 QBL	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
8	9025 DEU	*04:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
9	9026 YAR	*04:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
10	9107 LKT3	*05:01		-	-	-	-	-	-	-	+	+	+	-	-	-	-	+	-	
11	9051 PITOUT	*04:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
12	9052 DBB	*04:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
13	9004 JESTHOM	*04:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
14	9071 OLGA	*03:01	*04:02	+	+	+	+	+	-	-	-	-	-	-	-	-	-	-	-	
15	9075 DKB	*04:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
16	9037 SWEIG007	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
17	9282 CTM3953540	*01:01	*02:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
18	9257 32367	*18:01	*105:01	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	
19	9038 BM16	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
20	9059 SLE005	*03:01		+	+	+	+	+	-	-	-	-	-	-	-	-	-	-	-	
21	9064 AMALA	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
22	9056 KOSE	*02:01	*13:01	+	-	-	+	-	-	-	-	-	-	-	+	-	-	+	-	
23	9124 IHL	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
24	9035 JBUSH	*04:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
25	9049 IBW9	*01:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
26	9285 WT49	*04:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
27	9191 CH1007	*03:01	*13:01	+	+	+	+	+	-	-	-	-	-	-	+	-	-	+	-	
28	9320 BEL5GB	*02:01	*04:01	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
29	9050 MOU	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
30	9021 RSH	*01:01	*105:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
31	9019 DUCAF	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
32	9297 HAG	*02:01	*03:01	+	+	+	+	+	-	-	-	-	-	-	-	-	-	-	-	
33	9098 MT14B	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
34	9104 DHIF	*04:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
35	9302 SSTO	*04:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
36	9024 KT17	*05:01		-	-	-	-	-	-	-	+	+	+	-	-	-	-	+	-	
37	9065 HHKB	*04:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
38	9099 LZL	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
39	9315 CML	*04:01	*14:01	+	+	+	+	+	+	-	-	-	-	-	-	-	+	-	-	
40	9134 WHONP199	*05:01		-	-	-	-	-	-	-	+	+	+	-	-	-	-	+	-	
41	9055 H0301	*05:01		-	-	-	-	-	-	-	+	+	+	-	-	-	-	+	-	
42	9066 TAB089	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
43	9076 T7526	*107:01		+	-	-	+	-	-	-	-	-	-	-	+	-	-	+	-	
44	9057 TEM	*04:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	
45	9239 SHJO	*04:02	*104:01	+	+	+	+	+	-	-	-	-	-	-	-	-	-	-	-	
46	9013 SCHU	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
47	9045 TUBO	*02:01	*04:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
48	9303 TER-ND	*02:01	*04:01	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	



101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

CELL LINE VALIDATION SHEET																			
DPB1 SSP kit ²																			
				Well ³															
				33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48
				Prod. No.:															
				202127033	202127034	202127035	202127036	202127037	202127038	202127039	202127040	202127041	202127042	202127043	202127044	202127045	202127046	202127047	202127048
	IHWC cell line ¹		DPB1																
1	9001 SA		*04:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707		*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011 E4181324		*02:01	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275 GU373		*04:02	*17:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011		*04:01	*14:01	-	+	-	-	-	-	-	+	-	-	-	-	-	-	-
6	9353 SM		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
7	9020 QBL		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026 YAR		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA		*03:01	*04:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075 DKB		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007		*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282 CTM3953540		*01:01	*02:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257 32367		*18:01	*105:01	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA		*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE		*02:01	*13:01	+	-	-	-	+	-	-	-	-	-	-	-	-	-	-
23	9124 IHL		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035 JBUSH		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9		*01:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007		*03:01	*13:01	+	-	-	-	+	-	-	-	-	-	-	-	-	-	-
28	9320 BEL5GB		*02:01	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH		*01:01	*105:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019 DUCAF		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG		*02:01	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B		*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL		*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML		*04:01	*14:01	-	+	-	-	-	-	-	+	-	-	-	-	-	-	-
40	9134 WHONP199		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055 H0301		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526		*107:01		+	-	-	-	+	-	-	-	-	-	-	-	-	-	-
44	9057 TEM		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO		*04:02	*104:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU		*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO		*02:01	*04:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303 TER-ND		*02:01	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

101.311-24/04 – including *Taq* pol., IFU-01
 101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
 “Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

CELL LINE VALIDATION SHEET																				
DPB1 SSP kit ²																				
				Well ³																
				49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
				Prod. No.:	202127049	202127050	202127051	202127052	202127053	202127054	202127055	202127056	202127057	202127058	202127059	202127060	202127061	202127062	202127063	202127064
IHC cell line ¹		DPB1																		
1	9001 SA	*04:02		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
2	9280 LK707	*04:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	
3	9011 E4181324	*02:01	*04:01	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	
4	9275 GU373	*04:02	*17:01	-	-	-	-	-	-	-	-	-	-	+	-	+	-	+	-	
5	9009 KAS011	*04:01	*14:01	-	-	-	-	-	-	-	-	-	-	-	-	+	-	+	-	
6	9353 SM	*05:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
7	9020 QBL	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
8	9025 DEJ	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
9	9026 YAR	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
10	9107 LKT3	*05:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
11	9051 PITOUT	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
12	9052 DBB	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
13	9004 JESTHOM	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
14	9071 OLGA	*03:01	*04:02	-	-	-	-	-	-	-	-	-	+	-	-	-	-	+	-	
15	9075 DKB	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
16	9037 SWEIG007	*04:02		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
17	9282 CTM3953540	*01:01	*02:01	-	+	-	-	-	-	-	-	-	-	-	-	+	-	+	-	
18	9257 32367	*18:01	*105:01	-	+	-	-	-	-	-	-	-	-	-	-	+	-	+	-	
19	9038 BM16	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
20	9059 SLE005	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
21	9064 AMALA	*04:02		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
22	9056 KOSE	*02:01	*13:01	-	-	-	-	-	-	-	-	-	-	-	-	+	-	+	-	
23	9124 IHL	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
24	9035 JBUSH	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
25	9049 IBW9	*01:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	+	-	
26	9285 WT49	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
27	9191 CH1007	*03:01	*13:01	-	-	-	-	-	-	-	-	-	-	+	+	-	-	+	-	
28	9320 BEL5GB	*02:01	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
29	9050 MOU	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
30	9021 RSH	*01:01	*105:01	-	+	-	-	-	-	-	-	-	-	-	-	+	-	+	-	
31	9019 DUCAF	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
32	9297 HAG	*02:01	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	+	-	+	-	
33	9098 MT14B	*04:02		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
34	9104 DHIF	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
35	9302 SSTO	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
36	9024 KT17	*05:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
37	9065 HHKB	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
38	9099 LZL	*04:02		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
39	9315 CML	*04:01	*14:01	-	-	-	-	-	-	-	-	-	-	-	-	+	-	+	-	
40	9134 WHONP199	*05:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
41	9055 H0301	*05:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	
42	9066 TAB089	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
43	9076 T7526	*107:01		-	-	-	-	-	-	-	-	-	-	+	-	-	-	+	-	
44	9057 TEM	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
45	9239 SHJO	*04:02	*104:01	-	-	-	-	-	-	-	-	-	-	-	+	+	-	+	-	
46	9013 SCHU	*04:02		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	
47	9045 TUBO	*02:01	*04:02	-	-	-	-	-	-	-	-	-	+	-	-	+	-	-	-	
48	9303 TER-ND	*02:01	*04:01	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

CELL LINE VALIDATION SHEET																					
DPB1 SSP kit ²																					
					Well ³																
					65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	
					Prod. No.:	202127065	202127066	202127067	202127068	202127069	202127070	202127071	202127072	202127073	202127074	202127075	202127076	202127077	202127078	202127079	202127080
IHCW cell line ¹		DPB1																			
1	9001 SA	*04:02			-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
2	9280 LK707	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011 E4181324	*02:01	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275 GU373	*04:02	*17:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
5	9009 KAS011	*04:01	*14:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
6	9353 SM	*05:01			-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
7	9020 QBL	*02:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026 YAR	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3	*05:01			-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
11	9051 PITOUT	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA	*03:01	*04:02		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
15	9075 DKB	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007	*04:02			-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
17	9282 CTM3953540	*01:01	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257 32367	*18:01	*105:01		-	-	-	-	-	-	-	-	-	+	-	+	-	-	-	-	-
19	9038 BM16	*02:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005	*03:01			-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
21	9064 AMALA	*04:02			-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
22	9056 KOSE	*02:01	*13:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL	*02:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035 JBUSH	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9	*01:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007	*03:01	*13:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
28	9320 BEL5GB	*02:01	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU	*02:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH	*01:01	*105:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
31	9019 DUCAF	*02:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG	*02:01	*03:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
33	9098 MT14B	*04:02			-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
34	9104 DHIF	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17	*05:01			-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
37	9065 HHKB	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL	*04:02			-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
39	9315 CML	*04:01	*14:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
40	9134 WHONP199	*05:01			-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
41	9055 H0301	*05:01			-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
42	9066 TAB089	*02:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526	*107:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM	*04:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO	*04:02	*104:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
46	9013 SCHU	*04:02			-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
47	9045 TUBO	*02:01	*04:02		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
48	9303 TER-ND	*02:01	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

CELL LINE VALIDATION SHEET																			
DPB1 SSP kit ²																			
				Well ³															
				81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	
				202127081	202127082	202127083	202127084	202127085	202127086	202127087	202127088	202127089	202127090	202127091	202127092	202127093	202127094	202127095	
			Prod. No.:																
	IHWC cell line ¹		DPB1																
1	9001 SA		*04:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
2	9280 LK707		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
3	9011 E4181324		*02:01	*04:01	-	+	+	-	-	-	-	-	-	-	-	-	-	-	
4	9275 GU373		*04:02	*17:01	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
5	9009 KAS011		*04:01	*14:01	-	+	+	-	-	-	-	-	-	-	-	-	-	-	
6	9353 SM		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	
7	9020 QBL		*02:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
8	9025 DEU		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
9	9026 YAR		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
10	9107 LKT3		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	
11	9051 PITOUT		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
12	9052 DBB		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
13	9004 JESTHOM		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
14	9071 OLGA		*03:01	*04:02	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
15	9075 DKB		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
16	9037 SWEIG007		*04:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
17	9282 CTM3953540		*01:01	*02:01	-	+	-	-	-	-	-	-	-	-	-	-	+	-	
18	9257 32367		*18:01	*105:01	-	+	-	-	-	-	+	-	-	-	-	-	-	-	
19	9038 BM16		*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
20	9059 SLE005		*03:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
21	9064 AMALA		*04:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
22	9056 KOSE		*02:01	*13:01	-	+	-	-	-	-	-	-	-	-	-	-	-	+	
23	9124 IHL		*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
24	9035 JBUSH		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
25	9049 IBW9		*01:01		-	+	-	-	-	-	-	-	-	-	-	-	+	-	
26	9285 WT49		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
27	9191 CH1007		*03:01	*13:01	-	+	-	-	-	-	-	-	-	-	-	-	-	+	
28	9320 BEL5GB		*02:01	*04:01	-	+	+	-	-	-	-	-	-	-	-	-	-	-	
29	9050 MOU		*02:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
30	9021 RSH		*01:01	*105:01	-	+	-	-	-	-	-	-	-	-	-	-	+	-	
31	9019 DUCAF		*02:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
32	9297 HAG		*02:01	*03:01	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
33	9098 MT14B		*04:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
34	9104 DHIF		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
35	9302 SSTO		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
36	9024 KT17		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	
37	9065 HHKB		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
38	9099 LZL		*04:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
39	9315 CML		*04:01	*14:01	-	+	+	-	-	-	-	-	-	-	-	-	-	-	
40	9134 WHONP199		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	
41	9055 H0301		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	
42	9066 TAB089		*02:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
43	9076 T7526		*107:01		-	+	-	-	-	-	-	-	-	-	-	-	-	+	
44	9057 TEM		*04:01		-	+	+	-	-	-	-	-	-	-	-	-	-	-	
45	9239 SHJO		*04:02	*104:01	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
46	9013 SCHU		*04:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	
47	9045 TUBO		*02:01	*04:02	-	+	-	-	-	-	-	-	-	-	-	-	-	-	
48	9303 TER-ND		*02:01	*04:01	-	+	+	-	-	-	-	-	-	-	-	-	-	-	



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 41 to 44, 47, 49, 51 to 54, 56, 61, 63 to 72, 74, 76 to 81, 84 to 87 and 89 to 93 were available.

The specificities of the primers in primer solutions 42, 43, 47, 49, 68 and 86 were tested by separately adding one or two 5'-primers, and one 3'-primer respectively. In primer solutions 41, 53, 54, 56, 63 to 67, 71, 81, 84, 85 and 89 to 91 it was only possible to test the 5'-primers, the 3'-primers were not possible to be tested. In primer solutions 44, 51, 52, 61, 69, 70, 72, 74, 76 to 80, 87, 92 and 93 it was only possible to test the 3'-primers, the 5'-primers were not possible to be tested. In primer solutions 1, 2, 12, 15, 16, 23, 25 to 27, 31 to 33, 36, 38, 40, 45, 49, 50, 55, 68, 75, 83 and 86 one or more 5'-primers were not possible to be tested. In primer solutions 5 to 9, 11, 16, 24, 28, 29, 31, 33, 37, 38, 42, 43, 46 to 48, 60, 73 and 94 one or more 3'-primers were not possible to be tested.

In addition, one or more primers in primer solutions 14, 29, 33, 38, 40, 45, 46, 48, 50 and 94 were tested by separately adding one or more 5'-primers or 3'-primers accordingly.

³The DPB1*04:02 allele is not amplified by primer mix 57 in the 9275 (GU373) and 9239 (SHJO) cell lines. For further details see the specificity table.

The DPB1*04:02 allele is not amplified by primer mix 58 in the 9275 (GU373) cell line. For further details see the specificity table.

The DPB1*04:02 allele is not amplified by primer mix 60 in the 9239 (SHJO) cell line. For further details see the specificity table.



101.311-24/04 – including *Taq* pol., IFU-01
101.311-24u/04u – without *Taq* pol., IFU-02

Visit www.caredx.com for
“Instructions for Use” (IFU)

Lot No.: **8L1**

Lot-specific information

ADDRESSES:

Manufacturer:

CareDx AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-508 939 00

Fax: +46-8-717 88 18

E-mail: orders-se@caredx.com

Web page: <https://labproducts.caredx.com/>

CareDx Lab Solutions Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-653-78171

Fax: 610-344-7989

E-mail: orders-us@caredx.com

Web page: <https://labproducts.caredx.com/>

For information on CareDx distributors worldwide, contact **CareDx AB**.

