

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

Olerup SSP® DQB1*03

Product number:	101.214-24/06 – including <i>Taq</i> pol. 101.214-24u/06u – without <i>Taq</i> pol.
Lot number:	6F6
Expiry date:	2020-06-01
Number of tests:	24 test – Product No. 101.214-24/24u 6 tests – Product No. 101.214-06/06u
Number of wells per test:	63+1
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
- Adhesive PCR seals	RT
- Product Insert	RT

This Product Description is only valid for Lot No. 6F6.

Complete product documentation consists of generic Instructions for Use (IFU), lot specific Product Insert, Worksheet and Certificate

**CHANGES COMPARED TO THE PREVIOUS OLERUP SSP®
DQB1*03 LOT (0F9)**

The DQB1*03 kit is updated for new alleles to enable separation of:

- Confirmed¹ alleles as listed in the IMGT/HLA database.
- Null and Alternatively expressed alleles.

The format of the Worksheet has been changed.

One well has been added to DQB1*03, well **64**.

¹As described in section Uniquely Identified Alleles.

The DQB1*03 primer set, specificity and interpretation tables have been updated for the HLA-DQB1 alleles described since the previous *Olerup SSP®* DQB1*03 lot was made (**Lot No. 0F9**). The kit design is based on IMGT/HLA database 3.29.0.



101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
13	-	Added	3'-primers added for the DQB1*03:01:36 and DQB1*03:255 alleles.
16	Added	-	5'-primer added for the DQB1*03:03:15 allele.
19	-	Modified, added	3'-primer modified for increased yield, 3'-primer added for the DQB1*03:01:36 allele.
20	Added	-	5'-primer added for the DQB1*03:246 allele.
29	Added	-	5'-primer added for the DQB1*03:245 allele.
30	-	-	Strength of control band has been optimized.
31	-	Added	3'-primer added for the DQB1*03:255 allele.
33	Added	-	5'-primer added from well 45 for improved HLA-specific amplification.
34	-	-	Strength of control band has been optimized.
36	Added	-	5'-primer added for the DQB1*03:269N allele.
38	-	Added	3'-primer added for the DQB1*03:194 allele.
41	Added	-	5'-primer added for the DQB1*03:269N allele.
43	-	Added	3'-primer added for the DQB1*03:247 allele.
45	Moved	-	5'-primer moved to well 33 for improved HLA-specific amplification.
48	Exchanged	-	5'-primer exchanged for decreased tendency of primer oligomer formation.
50	-	-	Exchange of positive control primer pair for improved HLA-specific amplification.
51	-	Added	3'-primer added for the DQB1*03:239 allele.
56	Moved	Moved	Primer pair moved to well 63 for decreased tendency of primer oligomer formation.
57	Added	-	5'-primer added for the DQB1*03:246 allele.
62	-	Added	3'-primer added for the DQB1*03:247 allele.
63	Added	Added	Negative control moved to well 63, primer pair added from well 56 for decreased tendency of primer oligomer formation.
64	-	-	Negative control added from well 63.

Change in revision R01 compared to R00:

1. Primer mix 43 does not amplify the DQB1*03:49 and the DQB1*02:10, 02:81 and 05:152 alleles. Thus, this lot of the DQB1*03 kit cannot distinguish the DQB1*03:49 and the DQB1*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:02, 03:01:04-03:01:07, 03:01:09-03:01:12, 03:01:14, 03:01:16-03:01:24, 03:01:26-03:01:33, 03:01:35-03:01:36, 03:01:38-03:01:39, 03:73, 03:92-03:93, 03:241-03:242, 03:253-03:254, 03:264 and 03:266 alleles. This has been corrected in the Interpretation and Specificity Tables.

Changes in revision R02 compared to R01:

1. Primer mix 28 does not amplify the DQB1*03:38 allele. This correction has been implemented in the specificity and interpretation tables.

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

Well **64** contains Negative Control primer pairs, that will amplify more than 95% of the *Olerup* SSP[®] HLA Class I, DRB, DQB1, DPB1 and DQA1 amplicons as well as all the amplicons generated by the control primer pairs matching the human growth hormone gene.

HLA-specific PCR product sizes range from 75 to 200 base pairs.
The PCR product generated by the positive control primer pair is 430 base pairs.

Length of PCR product	105	200	105	80	75	80	85
5'-primer¹	164	340	440	45	45	43	36
	5'-CAC ^{3'}	5'-Agg ^{3'}	5'-TTA ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-TAC ^{3'}
							36
							5'-TAT ^{3'}
3'-primer²	231	2nd I	507	59	58	57	47
	5'-TgC ^{3'}	5'-AAA ^{3'}	5'-TTg ^{3'}	5'-CTC ^{3'}	5'-ggC ^{3'}	5'-CTC ^{3'}	5'-ACA ^{3'}
							48
							5'-gCA ^{3'}
							48
							5'-gCC ^{3'}
							52
							5'-TgT ^{3'}
A*	+	+	+				
B*	+	+	+				
C*	+	+	+				
DRB1				+	+		
DRB3				+	+		
DRB5				+			
DQB1					+		
DPB1						+	
DQA1							+

¹The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

²The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

PRODUCT DESCRIPTION

DQB1*03 SSP subtyping

CONTENT

The primer set contains 5'- and 3'-primers for identifying the DQB1*03:01 to DQB1*03:269N alleles.

Please note that DQB1 amplifications usually are somewhat less pronounced than e.g. DRB and DQA1 amplifications even when using the same DNA preparation and exactly the same experimental procedures.

PLATE LAYOUT

Each test consists of 64 PCR reactions in a 64 well PCR plate.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48
49	50	51	52	53	54	55	56
57	58	59	60	61	62	63	NC

The 64 well cut PCR plate is marked with 'DQB1*03' in silver/gray ink.

Well No. 1 is marked with the Lot No. '6F6'.

Wells 1 to 63 – DQB1*03 high resolution primers.

Well 64 – Negative Control (NC).

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded.

The PCR plates are covered with a PCR-compatible foil.

INTERPRETATION

Due to the sharing of sequence motifs, non-DQB1*03 alleles will be amplified by some primer mixes. The interpretation of DQB1*03 subtypings is not influenced by the DQB2 and DQB3 genes. For further details see Specificity Table.

UNIQUELY IDENTIFIED ALLELES

All the DQB1*03 alleles, i.e. **DQB1*03:01 to DQB1*03:269N**, recognized by the HLA Nomenclature Committee in October 2016^{1,2} will be amplified by the primers in the DQB1*03 subtyping kit^{3,4}.

The DQB1*03 kit enables separation of the confirmed DQB1*03 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources. Current allele confirmation status for DQB1*03 alleles is listed below.

The DQB1*03 also enables identification of null and alternatively expressed alleles.

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

The following DQB1*03 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix
DQB1*03:40, 03:137	9
DQB1*03:169, 03:196	57

¹HLA-DQB1 alleles listed on the IMGT/HLA web page 2017-August-10, release 3.29.0, www.ebi.ac.uk/imgt/hla.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

³The HLA-DQB1*03 primer set cannot separate the DQB1*03:25:01 from the DQB1*04:03:03 alleles. These alleles can be distinguished by the DQ low resolution kit and/or the HLA-DQB1*04 high resolution kit.

⁴This lot of the DQB1*03 kit cannot distinguish the DQB1*03:49 and the DQB1*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:02, 03:01:04-03:01:07, 03:01:09-03:01:12, 03:01:14, 03:01:16-03:01:24, 03:01:26-03:01:33, 03:01:35-03:01:36, 03:01:38-03:01:39, 03:73, 03:92-03:93, 03:241-03:242, 03:253-03:254, 03:264 and 03:266 alleles.

RESOLUTION IN HOMO- AND HETEROZYGOTES

Results file with resolution in DQB1*03 homo- and heterozygotes is available upon request.

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

ALLELE CONFIRMATION STATUS

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DQB1*03:01:01:01	Confirmed	DQB1*03:01:35	Unconfirmed	DQB1*03:04:01	Confirmed	DQB1*03:43	Unconfirmed
DQB1*03:01:01:02	Unconfirmed	DQB1*03:01:36	Unconfirmed	DQB1*03:04:02	Unconfirmed	DQB1*03:44	Confirmed
DQB1*03:01:01:03	Confirmed	DQB1*03:01:37	Unconfirmed	DQB1*03:04:03	Unconfirmed	DQB1*03:45	Confirmed
DQB1*03:01:01:04	Unconfirmed	DQB1*03:01:38	Unconfirmed	DQB1*03:05:01	Confirmed	DQB1*03:46	Unconfirmed
DQB1*03:01:01:05	Confirmed	DQB1*03:01:39	Unconfirmed	DQB1*03:05:02	Unconfirmed	DQB1*03:47	Confirmed
DQB1*03:01:01:06	Unconfirmed	DQB1*03:02:01:01	Confirmed	DQB1*03:05:03	Confirmed	DQB1*03:48	Unconfirmed
DQB1*03:01:01:07	Confirmed	DQB1*03:02:01:02	Confirmed	DQB1*03:05:04	Unconfirmed	DQB1*03:49	Unconfirmed
DQB1*03:01:01:08	Confirmed	DQB1*03:02:01:03	Unconfirmed	DQB1*03:06	Unconfirmed	DQB1*03:50	Unconfirmed
DQB1*03:01:01:09	Unconfirmed	DQB1*03:02:01:04	Unconfirmed	DQB1*03:07	Unconfirmed	DQB1*03:51	Unconfirmed
DQB1*03:01:01:10	Unconfirmed	DQB1*03:02:02	Unconfirmed	DQB1*03:08	Unconfirmed	DQB1*03:52	Confirmed
DQB1*03:01:01:11	Unconfirmed	DQB1*03:02:03	Unconfirmed	DQB1*03:09	Confirmed	DQB1*03:53	Confirmed
DQB1*03:01:01:12	Unconfirmed	DQB1*03:02:04	Unconfirmed	DQB1*03:10:01	Unconfirmed	DQB1*03:54	Confirmed
DQB1*03:01:01:14	Unconfirmed	DQB1*03:02:05	Unconfirmed	DQB1*03:10:02	Confirmed	DQB1*03:55	Confirmed
DQB1*03:01:01:15	Unconfirmed	DQB1*03:02:06	Confirmed	DQB1*03:11	Confirmed	DQB1*03:56	Confirmed
DQB1*03:01:01:16	Unconfirmed	DQB1*03:02:07	Confirmed	DQB1*03:12	Unconfirmed	DQB1*03:57	Unconfirmed
DQB1*03:01:01:17	Unconfirmed	DQB1*03:02:08	Unconfirmed	DQB1*03:13	Confirmed	DQB1*03:58	Confirmed
DQB1*03:01:01:18	Unconfirmed	DQB1*03:02:09	Confirmed	DQB1*03:14:01	Unconfirmed	DQB1*03:59	Confirmed
DQB1*03:01:02	Unconfirmed	DQB1*03:02:10	Confirmed	DQB1*03:14:02	Confirmed	DQB1*03:60	Unconfirmed
DQB1*03:01:03	Confirmed	DQB1*03:02:11	Confirmed	DQB1*03:15	Unconfirmed	DQB1*03:61	Unconfirmed
DQB1*03:01:04	Confirmed	DQB1*03:02:12	Confirmed	DQB1*03:16	Confirmed	DQB1*03:62	Confirmed
DQB1*03:01:05	Unconfirmed	DQB1*03:02:13	Confirmed	DQB1*03:17:01	Confirmed	DQB1*03:63	Confirmed
DQB1*03:01:06	Unconfirmed	DQB1*03:02:14	Unconfirmed	DQB1*03:17:02	Confirmed	DQB1*03:64	Confirmed
DQB1*03:01:07	Unconfirmed	DQB1*03:02:15	Confirmed	DQB1*03:18	Unconfirmed	DQB1*03:65	Confirmed
DQB1*03:01:08	Confirmed	DQB1*03:02:16	Confirmed	DQB1*03:19:01	Confirmed	DQB1*03:66N	Unconfirmed
DQB1*03:01:09	Unconfirmed	DQB1*03:02:17	Unconfirmed	DQB1*03:19:02	Unconfirmed	DQB1*03:67	Confirmed
DQB1*03:01:10	Unconfirmed	DQB1*03:02:18	Unconfirmed	DQB1*03:20	Unconfirmed	DQB1*03:68	Confirmed
DQB1*03:01:11	Unconfirmed	DQB1*03:02:19	Confirmed	DQB1*03:21	Unconfirmed	DQB1*03:69	Unconfirmed
DQB1*03:01:12	Unconfirmed	DQB1*03:02:20	Unconfirmed	DQB1*03:22	Unconfirmed	DQB1*03:70	Confirmed
DQB1*03:01:13	Unconfirmed	DQB1*03:02:21	Confirmed	DQB1*03:23:01	Unconfirmed	DQB1*03:71	Unconfirmed
DQB1*03:01:14	Confirmed	DQB1*03:02:22	Unconfirmed	DQB1*03:23:02	Unconfirmed	DQB1*03:72	Confirmed
DQB1*03:01:15	Unconfirmed	DQB1*03:02:23	Unconfirmed	DQB1*03:24	Unconfirmed	DQB1*03:73	Unconfirmed
DQB1*03:01:16	Unconfirmed	DQB1*03:02:24	Unconfirmed	DQB1*03:25:01	Confirmed	DQB1*03:74	Unconfirmed
DQB1*03:01:17	Confirmed	DQB1*03:02:25	Unconfirmed	DQB1*03:25:02	Unconfirmed	DQB1*03:75	Confirmed
DQB1*03:01:18	Unconfirmed	DQB1*03:03:02:01	Confirmed	DQB1*03:26	Confirmed	DQB1*03:76	Unconfirmed
DQB1*03:01:19	Unconfirmed	DQB1*03:03:02:02	Unconfirmed	DQB1*03:27	Confirmed	DQB1*03:77	Unconfirmed
DQB1*03:01:20	Unconfirmed	DQB1*03:03:02:03	Unconfirmed	DQB1*03:28	Unconfirmed	DQB1*03:78	Unconfirmed
DQB1*03:01:21	Confirmed	DQB1*03:03:02:04	Confirmed	DQB1*03:29	Unconfirmed	DQB1*03:79	Unconfirmed
DQB1*03:01:22	Confirmed	DQB1*03:03:03	Unconfirmed	DQB1*03:30	Unconfirmed	DQB1*03:80	Unconfirmed
DQB1*03:01:23	Unconfirmed	DQB1*03:03:04	Confirmed	DQB1*03:31	Unconfirmed	DQB1*03:81	Confirmed
DQB1*03:01:24	Unconfirmed	DQB1*03:03:05	Confirmed	DQB1*03:32	Unconfirmed	DQB1*03:82	Unconfirmed
DQB1*03:01:25	Unconfirmed	DQB1*03:03:06	Unconfirmed	DQB1*03:33	Unconfirmed	DQB1*03:83	Unconfirmed
DQB1*03:01:26	Unconfirmed	DQB1*03:03:07	Confirmed	DQB1*03:34	Unconfirmed	DQB1*03:84N	Unconfirmed
DQB1*03:01:27	Confirmed	DQB1*03:03:08	Confirmed	DQB1*03:35	Unconfirmed	DQB1*03:85	Unconfirmed
DQB1*03:01:28	Unconfirmed	DQB1*03:03:09	Unconfirmed	DQB1*03:36	Unconfirmed	DQB1*03:86	Unconfirmed
DQB1*03:01:29	Unconfirmed	DQB1*03:03:10	Unconfirmed	DQB1*03:37	Confirmed	DQB1*03:87	Unconfirmed
DQB1*03:01:30	Unconfirmed	DQB1*03:03:11	Confirmed	DQB1*03:38	Unconfirmed	DQB1*03:88	Unconfirmed
DQB1*03:01:31	Unconfirmed	DQB1*03:03:12	Unconfirmed	DQB1*03:39	Unconfirmed	DQB1*03:89	Unconfirmed
DQB1*03:01:32	Unconfirmed	DQB1*03:03:13	Unconfirmed	DQB1*03:40	Confirmed	DQB1*03:90N	Unconfirmed
DQB1*03:01:33	Unconfirmed	DQB1*03:03:14	Unconfirmed	DQB1*03:41	Unconfirmed	DQB1*03:91Q	Unconfirmed
DQB1*03:01:34	Unconfirmed	DQB1*03:03:15	Unconfirmed	DQB1*03:42	Unconfirmed	DQB1*03:92	Unconfirmed

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DQB1*03:93	Unconfirmed	DQB1*03:143	Unconfirmed	DQB1*03:193	Unconfirmed	DQB1*03:243	Confirmed
DQB1*03:94	Unconfirmed	DQB1*03:144	Unconfirmed	DQB1*03:194	Confirmed	DQB1*03:244	Unconfirmed
DQB1*03:95N	Unconfirmed	DQB1*03:145	Unconfirmed	DQB1*03:195	Confirmed	DQB1*03:245	Confirmed
DQB1*03:96	Unconfirmed	DQB1*03:146	Unconfirmed	DQB1*03:196	Confirmed	DQB1*03:246	Confirmed
DQB1*03:97	Unconfirmed	DQB1*03:147	Unconfirmed	DQB1*03:197Q	Confirmed	DQB1*03:247	Confirmed
DQB1*03:98	Unconfirmed	DQB1*03:148	Unconfirmed	DQB1*03:198	Unconfirmed	DQB1*03:248	Unconfirmed
DQB1*03:99Q	Unconfirmed	DQB1*03:149	Unconfirmed	DQB1*03:199	Confirmed	DQB1*03:249	Unconfirmed
DQB1*03:100	Unconfirmed	DQB1*03:150	Confirmed	DQB1*03:200	Confirmed	DQB1*03:250	Unconfirmed
DQB1*03:101	Unconfirmed	DQB1*03:151	Unconfirmed	DQB1*03:201	Confirmed	DQB1*03:251	Unconfirmed
DQB1*03:102	Confirmed	DQB1*03:152	Unconfirmed	DQB1*03:202	Unconfirmed	DQB1*03:252	Unconfirmed
DQB1*03:103	Unconfirmed	DQB1*03:153	Unconfirmed	DQB1*03:203	Unconfirmed	DQB1*03:253	Unconfirmed
DQB1*03:104	Unconfirmed	DQB1*03:154	Unconfirmed	DQB1*03:204	Confirmed	DQB1*03:254	Unconfirmed
DQB1*03:105	Confirmed	DQB1*03:155	Unconfirmed	DQB1*03:205	Confirmed	DQB1*03:255	Confirmed
DQB1*03:106	Confirmed	DQB1*03:156	Unconfirmed	DQB1*03:206	Unconfirmed	DQB1*03:256	Unconfirmed
DQB1*03:107	Unconfirmed	DQB1*03:157	Unconfirmed	DQB1*03:207	Confirmed	DQB1*03:257	Unconfirmed
DQB1*03:108	Unconfirmed	DQB1*03:158	Confirmed	DQB1*03:208	Unconfirmed	DQB1*03:258	Unconfirmed
DQB1*03:109	Confirmed	DQB1*03:159	Unconfirmed	DQB1*03:209	Unconfirmed	DQB1*03:259	Unconfirmed
DQB1*03:110	Unconfirmed	DQB1*03:160	Unconfirmed	DQB1*03:210	Unconfirmed	DQB1*03:260	Unconfirmed
DQB1*03:111	Unconfirmed	DQB1*03:161	Unconfirmed	DQB1*03:211	Confirmed	DQB1*03:261	Unconfirmed
DQB1*03:112	Unconfirmed	DQB1*03:162	Confirmed	DQB1*03:212	Confirmed	DQB1*03:262	Unconfirmed
DQB1*03:113	Confirmed	DQB1*03:163	Unconfirmed	DQB1*03:213N	Unconfirmed	DQB1*03:263	Unconfirmed
DQB1*03:114	Confirmed	DQB1*03:164	Unconfirmed	DQB1*03:214	Unconfirmed	DQB1*03:264	Unconfirmed
DQB1*03:115	Unconfirmed	DQB1*03:165	Unconfirmed	DQB1*03:215	Unconfirmed	DQB1*03:265	Unconfirmed
DQB1*03:116	Unconfirmed	DQB1*03:166	Unconfirmed	DQB1*03:216	Unconfirmed	DQB1*03:266	Unconfirmed
DQB1*03:117	Unconfirmed	DQB1*03:167	Unconfirmed	DQB1*03:217	Unconfirmed	DQB1*03:267	Unconfirmed
DQB1*03:118N	Confirmed	DQB1*03:168	Unconfirmed	DQB1*03:218	Unconfirmed	DQB1*03:268	Unconfirmed
DQB1*03:119	Confirmed	DQB1*03:169	Unconfirmed	DQB1*03:219	Unconfirmed	DQB1*03:269N	Unconfirmed
DQB1*03:120	Unconfirmed	DQB1*03:170	Confirmed	DQB1*03:220	Unconfirmed		
DQB1*03:121	Unconfirmed	DQB1*03:171	Unconfirmed	DQB1*03:221	Confirmed		
DQB1*03:122	Unconfirmed	DQB1*03:172	Unconfirmed	DQB1*03:222	Unconfirmed		
DQB1*03:123	Unconfirmed	DQB1*03:173	Unconfirmed	DQB1*03:223	Unconfirmed		
DQB1*03:124	Unconfirmed	DQB1*03:174	Unconfirmed	DQB1*03:224	Unconfirmed		
DQB1*03:125	Unconfirmed	DQB1*03:175	Unconfirmed	DQB1*03:225	Confirmed		
DQB1*03:126	Unconfirmed	DQB1*03:176	Unconfirmed	DQB1*03:226	Unconfirmed		
DQB1*03:127	Unconfirmed	DQB1*03:177	Confirmed	DQB1*03:227	Unconfirmed		
DQB1*03:128	Unconfirmed	DQB1*03:178	Unconfirmed	DQB1*03:228	Confirmed		
DQB1*03:129	Unconfirmed	DQB1*03:179	Unconfirmed	DQB1*03:229	Unconfirmed		
DQB1*03:130	Confirmed	DQB1*03:180	Unconfirmed	DQB1*03:230	Unconfirmed		
DQB1*03:131	Unconfirmed	DQB1*03:181	Unconfirmed	DQB1*03:231	Unconfirmed		
DQB1*03:132	Confirmed	DQB1*03:182	Unconfirmed	DQB1*03:232	Confirmed		
DQB1*03:133	Unconfirmed	DQB1*03:183	Unconfirmed	DQB1*03:233	Confirmed		
DQB1*03:134	Unconfirmed	DQB1*03:184	Confirmed	DQB1*03:234	Unconfirmed		
DQB1*03:135	Unconfirmed	DQB1*03:185	Unconfirmed	DQB1*03:235	Unconfirmed		
DQB1*03:136	Unconfirmed	DQB1*03:186	Unconfirmed	DQB1*03:236	Unconfirmed		
DQB1*03:137	Unconfirmed	DQB1*03:187	Unconfirmed	DQB1*03:237N	Unconfirmed		
DQB1*03:138	Confirmed	DQB1*03:188	Confirmed	DQB1*03:238	Unconfirmed		
DQB1*03:139	Confirmed	DQB1*03:189	Unconfirmed	DQB1*03:239	Confirmed		
DQB1*03:140	Unconfirmed	DQB1*03:190	Unconfirmed	DQB1*03:240	Confirmed		
DQB1*03:141	Unconfirmed	DQB1*03:191	Confirmed	DQB1*03:241	Unconfirmed		
DQB1*03:142	Unconfirmed	DQB1*03:192	Unconfirmed	DQB1*03:242	Unconfirmed		

¹Allele status “confirmed” or “unconfirmed” as listed on the IMGT/HLA web page 2017-August-10, release 3.29.0, www.ebi.ac.uk/imgt/hla.

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information
SPECIFICITY TABLE

DQB1*03 SSP subtyping

Specificities and sizes of the PCR products of the 63+1 primer mixes used for DQB1*03 SSP subtyping

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified DQB1*03 alleles ³	Other amplified DQB1 alleles
1 ⁷	170 bp	515 bp	*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:39, 03:04:01-03:04:03, 03:09, 03:13, 03:16, 03:19:01-03:19:02, 03:21-03:22, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:60, 03:69, 03:71, 03:73-03:78, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108, 03:113-03:116, 03:118N-03:122, 03:127-03:131, 03:134-03:135, 03:139-03:140, 03:142-03:144, 03:147-03:148, 03:150-03:152, 03:157-03:160, 03:162-03:167, 03:169-03:173, 03:182, 03:184, 03:188, 03:191-03:194, 03:196-03:198, 03:201-03:202, 03:206-03:208, 03:216, 03:218-03:219, 03:231-03:232, 03:235-03:236, 03:241-03:243, 03:246, 03:252-03:257, 03:260, 03:264, 03:266-03:268	*04:09
2 ⁷	220 bp	515 bp	*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:07, 03:01:09-03:01:39, 03:04:01-03:04:03, 03:09-03:10:02, 03:13-03:14:02, 03:16, 03:19:01-03:19:02, 03:21-03:22, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:60, 03:69, 03:71, 03:73, 03:75-03:77, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108, 03:114-03:116, 03:118N-03:122, 03:127-03:131, 03:133-03:135, 03:138-03:140, 03:142-03:144, 03:147-03:148, 03:150, 03:152, 03:157-03:160, 03:162-03:167, 03:169-03:173, 03:180, 03:182-03:183, 03:186-03:188, 03:191-03:198, 03:201-03:202, 03:206-03:208, 03:216, 03:218-03:219, 03:231-03:232, 03:235-03:236, 03:241-03:243, 03:246, 03:252-03:257, 03:260, 03:264, 03:266-03:268	
3 ⁵	135 bp	515 bp	*03:02:01:01-03:02:25, 03:07-03:08, 03:11, 03:18, 03:32, 03:37, 03:45, 03:62-03:64, 03:66N-03:68, 03:70, 03:81, 03:85, 03:106-03:107, 03:125, 03:146, 03:153, 03:161, 03:174-03:175, 03:178-03:179, 03:184-03:185, 03:189-03:190, 03:199, 03:203-03:205, 03:210-03:211, 03:213N-03:215, 03:220-03:221, 03:223-03:225, 03:228-03:229, 03:233, 03:237N, 03:240, 03:245, 03:247, 03:251, 03:261, 03:263, 03:265, 03:269N	*02:01:01-02:01:20, 02:01:22-02:02:04, 02:04-02:16, 02:18N-02:36, 02:38-02:76, 02:78-02:94, 02:96N-02:101, 06:29, 06:63, 06:123, 06:139

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

4	220 bp	515 bp	*03:02:01:01-03:03:15, 03:06-03:08, 03:11-03:12, 03:15, 03:18, 03:20, 03:23:01-03:23:02, 03:25:01-03:26, 03:30-03:34, 03:37-03:41, 03:43, 03:45, 03:62-03:68, 03:70, 03:74, 03:79, 03:81, 03:85-03:91Q, 03:95N-03:99Q, 03:106- 03:107, 03:110-03:113, 03:117, 03:123- 03:126, 03:136-03:137, 03:141, 03:145- 03:146, 03:149, 03:153, 03:155-03:156, 03:161, 03:168, 03:174-03:179, 03:184- 03:185, 03:190, 03:199-03:200, 03:203- 03:205, 03:209-03:215, 03:217, 03:220- 03:225, 03:227-03:230, 03:233-03:234, 03:237N-03:240, 03:245, 03:247-03:249, 03:251, 03:259, 03:261, 03:263, 03:265, 03:269N	*02:01:01-02:36, 02:38- 02:94, 02:96N-02:101, 04:03:01-04:03:03
5⁵	135 bp	430 bp	*03:04:01-03:04:03, 03:14:01-03:14:02, 03:80, 03:138	
6^{4,6}	95 bp 130 bp	430 bp	*03:20 *03:05:01-03:05:04, 03:17:01-03:17:02, 03:61, 03:72, 03:100, 03:132, 03:181, 03:226, 03:244, 03:250, 03:262	*04:09, 05:131
7⁴	115 bp	430 bp	*03:03:06, 03:06, 03:25:01-03:25:02, 03:42, 03:88 *03:84N	*04:01:03, 04:03:03
8⁴	150 bp 110 bp 145 bp	430 bp	*03:15 *03:07, 03:16	
9⁴	65 bp 145 bp	430 bp	*03:40 *03:08, 03:119, 03:137, 03:194	*06:02:02, 06:03:02
10	135 bp 260 bp	430 bp	*03:09 *03:11, 03:26, 03:244	
11⁶	135 bp	430 bp	*03:10:01-03:10:02, 03:12, 03:14:01- 03:14:02, 03:70, 03:179, 03:183, 03:195	*06:01:01-06:01:06, 06:01:08-06:01:15, 06:06 [?] , 06:43, 06:54N- 06:58, 06:98, 06:99:02- 06:105, 06:108, 06:120, 06:132, 06:140, 06:153, 06:157, 06:167-06:168, 06:177, 06:181, 06:194, 06:205, 06:209, 06:214, 06:229, 06:239
12	255 bp 260 bp	515 bp	*03:47 *03:01:01:01-03:01:01:12, 03:01:01:14- 03:01:07, 03:01:09-03:01:39, 03:04:01- 03:04:03, 03:09-03:10:02, 03:12- 03:14:02, 03:16, 03:19:01-03:19:02, 03:21-03:22, 03:24, 03:27-03:29, 03:35- 03:36, 03:42, 03:44, 03:46-03:60, 03:70- 03:71, 03:73-03:78, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108, 03:114-03:116, 03:118N-03:122, 03:127- 03:131, 03:133-03:135, 03:139-03:140, 03:143-03:144, 03:147-03:148, 03:150- 03:152, 03:157-03:160, 03:162-03:167, 03:169-03:173, 03:179, 03:182-03:183,	

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

Lot No.: 6F6			Lot-specific information	
			03:186-03:188, 03:191-03:198, 03:201-03:202, 03:206, 03:208, 03:216, 03:218-03:219, 03:231-03:232, 03:235-03:236, 03:241-03:243, 03:246, 03:252-03:257, 03:260, 03:264, 03:266-03:268	
13⁴	115 bp 170 bp 205 bp	430 bp	*03:255 *03:13, 03:232 *03:48	
14⁶	130 bp	430 bp	*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:12, 03:01:14, 03:01:16-03:01:24, 03:01:26-03:01:39, 03:04:01-03:04:03, 03:09-03:10:02, 03:13-03:14:02, 03:19:01-03:19:02, 03:21-03:22, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:60, 03:69, 03:71, 03:73, 03:75-03:77, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108-03:109, 03:114-03:116, 03:118N-03:122, 03:127-03:129, 03:131, 03:133-03:135, 03:138-03:140, 03:142-03:144, 03:147-03:148, 03:150, 03:152, 03:154, 03:157-03:160, 03:162-03:167, 03:169-03:171, 03:173, 03:180, 03:182-03:183, 03:186-03:188, 03:191-03:196, 03:198, 03:201, 03:206-03:208, 03:216, 03:218-03:219, 03:231-03:232, 03:235-03:236, 03:241-03:243, 03:246, 03:252-03:257, 03:260, 03:264, 03:266-03:268	*06:209
15	135 bp	515 bp	*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:07, 03:01:09-03:01:33, 03:01:35-03:01:36, 03:01:38-03:15, 03:17:01-03:22, 03:24, 03:26-03:103, 03:106-03:108, 03:110-03:152, 03:155-03:171, 03:173-03:188, 03:190-03:204, 03:206-03:213N, 03:215-03:216, 03:218-03:233, 03:235-03:257, 03:260, 03:262-03:269N	*04:09
16⁶	130 bp	430 bp	*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:02, 03:01:04-03:01:07, 03:01:09-03:01:12, 03:01:14, 03:01:16-03:01:24, 03:01:26-03:01:39, 03:03:02:01-03:03:15, 03:06, 03:09-03:10:02, 03:12-03:13, 03:15-03:17:02, 03:19:01-03:24, 03:26-03:31, 03:33-03:36, 03:38-03:44, 03:46-03:60, 03:65, 03:69, 03:71-03:79, 03:82-03:84N, 03:86-03:103, 03:108, 03:111-03:124, 03:126-03:129, 03:131, 03:133-03:137, 03:139-03:145, 03:147-03:152, 03:155-03:160, 03:162-03:171, 03:173, 03:176-03:177, 03:180, 03:182-03:183, 03:186-03:188, 03:191-03:198, 03:200-03:201, 03:206-03:209, 03:212, 03:216-03:219, 03:222, 03:227, 03:230-03:232, 03:235-03:236, 03:238, 03:241-03:244, 03:246, 03:248-03:249, 03:252-03:257, 03:260, 03:264, 03:266-03:268	*02:03, 02:77, 04:09-04:10

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

17	165 bp 225 bp	430 bp	*03:43, 03:87 *03:18, 03:85	
18⁷	175 bp	515 bp	*03:02:01:01-03:02:24, 03:02:25 [?] , 03:03:02:01-03:03:13, 03:03:14 [?] - 03:03:15 [?] , 03:05:01-03:05:04, 03:06 [?] - 03:08 [?] , 03:11 [?] -03:15 [?] , 03:17:01 [?] -03:18 [?] , 03:19:01-03:19:02, 03:20 [?] , 03:23:01 [?] - 03:23:02 [?] , 03:25:01, 03:25:02 [?] -03:26 [?] , 03:30-03:32, 03:33 [?] -03:34 [?] , 03:37 [?] , 03:38-03:39, 03:40 [?] , 03:41, 03:43, 03:45, 03:48 [?] , 03:52 [?] -03:71 [?] , 03:72, 03:74 [?] - 03:78 [?] , 03:79, 03:81 [?] -03:82 [?] , 03:85- 03:91Q, 03:95N-03:100, 03:101 [?] -03:104 [?] , 03:105, 03:106 [?] -03:112 [?] , 03:117, 03:118N [?] , 03:120 [?] -03:131 [?] , 03:132, 03:133 [?] -03:137 [?] , 03:138, 03:140 [?] - 03:149 [?] , 03:150, 03:151 [?] -03:163 [?] , 03:165 [?] -03:167 [?] , 03:168, 03:170 [?] - 03:179 [?] , 03:180-03:181, 03:183 [?] -03:185 [?] , 03:187 [?] -03:189 [?] , 03:190, 03:192 [?] - 03:194 [?] , 03:195, 03:199-03:200, 03:201 [?] - 03:205 [?] , 03:207 [?] -03:210 [?] , 03:211, 03:212 [?] -03:221 [?] , 03:222, 03:223 [?] - 03:230 [?] , 03:232 [?] -03:238 [?] , 03:239, 03:240 [?] , 03:244 [?] , 03:245, 03:247-03:248, 03:250-03:251, 03:255 [?] -03:262 [?] , 03:263, 03:265 [?] , 03:267 [?] -03:269N [?]	*04:01:01:01- 04:02:01:01, 04:02:02- 04:41N
19	175 bp	430 bp	*03:01:01:01-03:01:01:12, 03:01:01:14- 03:01:39, 03:04:01-03:04:03, 03:06 [?] - 03:08 [?] , 03:09-03:10:02, 03:11 [?] -03:15 [?] , 03:16, 03:17:01 [?] -03:18 [?] , 03:20 [?] , 03:21- 03:22, 03:23:01 [?] -03:23:02 [?] , 03:24, 03:25:02 [?] -03:26 [?] , 03:27-03:29, 03:33 [?] - 03:34 [?] , 03:35-03:36, 03:37 [?] , 03:40 [?] , 03:42, 03:44, 03:46-03:47, 03:48 [?] , 03:49- 03:51, 03:52 [?] -03:71 [?] , 03:73, 03:74 [?] - 03:78 [?] , 03:80, 03:81 [?] -03:82 [?] , 03:83- 03:84N, 03:92-03:94, 03:101 [?] -03:104 [?] , 03:106 [?] -03:112 [?] , 03:113-03:116, 03:118N [?] , 03:119, 03:120 [?] -03:131 [?] , 03:133 [?] -03:137 [?] , 03:139, 03:140 [?] - 03:149 [?] , 03:151 [?] -03:163 [?] , 03:164, 03:165 [?] -03:167 [?] , 03:169, 03:170 [?] - 03:179 [?] , 03:182, 03:183 [?] -03:185 [?] , 03:186, 03:187 [?] -03:189 [?] , 03:191, 03:192 [?] -03:194 [?] , 03:196-03:198, 03:201 [?] - 03:205 [?] , 03:206, 03:207 [?] -03:210 [?] , 03:212 [?] -03:221 [?] , 03:223 [?] -03:238 [?] , 03:240 [?] , 03:241-03:243, 03:244 [?] , 03:246, 03:252-03:254, 03:255 [?] -03:262 [?] , 03:264, 03:266, 03:267 [?] -03:269N [?]	
20	150 bp 240 bp	430 bp	*03:21 *03:246	*02:64, 02:79, 05:35
21	160 bp	515 bp	*03:22, 03:96, 03:251	*04:23, 05:103 ^w , 06:04:01 ^w -06:04:10 ^w , 06:04:11 [?] , 06:05:01 ^w - 06:05:02 ^w , 06:06 [?] ,

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

				06:07:01 ^w -06:07:02 ^w , 06:08:01 [?] -06:08:03 [?] , 06:09:01:01 ^w -06:09:06 ^w , 06:09:07 [?] , 06:11:01 [?] - 06:11:03 [?] , 06:12 ^w , 06:13:01 [?] -06:13:02 [?] , 06:15:01 [?] -06:17 [?] , 06:18:01 ^w -06:18:02 ^w , 06:19:01 [?] -06:20 [?] , 06:21 ^w , 06:22:01 [?] - 06:32:02 [?] , 06:34 ^w , 06:35 [?] , 06:36 ^w , 06:38 ^w , 06:40 [?] , 06:41 ^w -06:42 ^w , 06:45 [?] -06:46 [?] , 06:49 [?] - 06:50 [?] , 06:51:02 [?] , 06:52 ^w , 06:53:01 [?] - 06:54N [?] , 06:56 [?] -06:57 [?] , 06:58 ^w , 06:59 [?] -06:61 [?] , 06:63 [?] -06:76 [?] , 06:79:02 [?] -06:83 [?] , 06:84 ^w , 06:85 [?] , 06:89 [?] , 06:91 [?] -06:97 [?] , 06:118:01 ^w -06:118:03 ^w , 06:119 [?] -06:124 [?] , 06:126 [?] , 06:128 [?] , 06:129 ^w , 06:131 [?] - 06:147 [?] , 06:149 [?] - 06:184 [?] , 06:186 ^w , 06:189 ^w , 06:190:01 [?] - 06:199 [?] , 06:201 [?] - 06:204 [?] , 06:206:01 [?] , 06:206:02 ^w , 06:207 [?] - 06:208 [?] , 06:210 [?] - 06:215 [?] , 06:217 ^w , 06:229 [?] -06:236 [?] , 06:238 [?] , 06:241 [?]
22	160 bp	515 bp	*03:23:01	*05:03:12, 05:43:02, 06:03:01:01-06:03:03, 06:03:05-06:03:07, 06:03:11-06:03:25, 06:03:27, 06:04:02, 06:07:01, 06:08:01, 06:09:02, 06:11:01- 06:11:03, 06:26N, 06:28, 06:30-06:32:01, 06:40-06:41, 06:44, 06:59-06:62, 06:64- 06:65, 06:67, 06:90- 06:92, 06:94, 06:110, 06:118:01, 06:118:03, 06:128, 06:133-06:134, 06:141, 06:143-06:145, 06:148, 06:154, 06:165, 06:169-06:170, 06:184- 06:185, 06:187, 06:191, 06:195-06:196, 06:199, 06:203, 06:210, 06:218,

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
 101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

				06:221-06:223, 06:230, 06:234, 06:238
23^{4,5}	95 bp 125 bp 150 bp 185 bp	430 bp	*03:41, 03:224 *03:24-03:25:02, 03:79, 03:89, 03:234, 03:239 *03:84N *03:45	*04:03:01-04:03:03
24^{4,6}	120 bp 210 bp	430 bp	*03:39 *03:27, 03:48, 03:56	
25^{4,5}	105 bp 130 bp 155 bp 180 bp 240 bp	430 bp	*03:206 *03:24, 03:79 *03:80, 03:257 *03:29, 03:118N, 03:237N *03:34, 02:14:01	
26	180 bp	430 bp	*03:30, 03:52, 03:72, 03:100, 03:132, 03:179, 03:215	*04:09, 06:29, 06:96, 06:139, 06:209
27⁷	170 bp 225 bp 250 bp	430 bp	*03:32, 03:87 *03:85 *03:31	
28	200 bp 230 bp	430 bp	*03:37, 03:46, 03:252 *03:225	
29⁴	110 bp 165 bp 255 bp	515 bp	*03:33, 03:245 *03:35 *03:243	*06:228
30⁵	135 bp 175 bp	430 bp	*03:36, 03:149 *03:150	
31⁴	115 bp 195 bp 225 bp	430 bp	*03:255 *03:28 *03:56, 03:83	
32⁴	100 bp	430 bp	*03:44	
33⁴	85 bp 160 bp	430 bp	*03:59 *03:204	*06:37
34⁵	190 bp	430 bp	*03:58, 03:82, 03:198, 03:222	
35	135 bp 190 bp 235 bp	430 bp	*03:190, 03:205, 03:214 *03:51, 03:54 *03:117	*04:15, 05:102, 06:223 *04:35
36	145 bp 175 bp 200 bp	430 bp	*03:205, 03:214 *03:55, 03:60, 03:250 *03:269N	*04:15
37	180 bp 230 bp	430 bp	*03:60, 03:77 *03:53, 03:62, 03:225	
38	170 bp 225 bp	430 bp	*03:58, 03:194, 03:232 *03:57, 03:83	*06:82, 06:233
39⁴	105 bp 140 bp 190 bp 235 bp	430 bp	*03:240 *03:25:01-03:25:02, 03:66N, 03:234, 03:239 *03:68 *03:86	*04:03:01-04:03:03, 06:193N
40	180 bp 230 bp	430 bp	*03:67-03:68, 03:228 *03:212	*06:17, 06:24, 06:30, 06:42 *02:47
41	175 bp	430 bp	*03:63, 03:75, 03:168	

101.214-24/06 – including *Taq* polymerase, IFU-01
101.214-24u/06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

	210 bp		*03:201, 03:269N	
42⁴	80 bp	430 bp	*03:188	
	165 bp		*03:64, 03:76-03:77	
	220 bp		*03:201	
43⁵	175 bp	430 bp	*03:198, 03:222, 03:247	
	250 bp		*03:97	*06:86, 06:104
44	160 bp	430 bp	*03:199, 03:211	
	220 bp		*03:78	
	270 bp		*03:71	
45⁴	90 bp	430 bp	*03:65	
	170 bp		*03:133, 03:200	*02:84, 04:14
	230 bp		*03:81, 03:159	
46^{4,5}	115 bp	430 bp	*03:88-03:89	
	175 bp		*03:191	
47	175 bp	430 bp	*03:06 [?] -03:08 [?] , 03:11 [?] -03:15 [?] , 03:17:01 [?] - 03:18 [?] , 03:19:01-03:19:02, 03:20 [?] , 03:23:01 [?] -03:23:02 [?] , 03:25:02 [?] -03:26 [?] , 03:37 [?] , 03:40 [?] , 03:48 [?] , 03:52 [?] -03:71 [?] , 03:74 [?] -03:78 [?] , 03:81 [?] -03:82 [?] , 03:101 [?] - 03:104 [?] , 03:106 [?] -03:112 [?] , 03:118N [?] , 03:120 [?] -03:131 [?] , 03:133 [?] -03:137 [?] , 03:140 [?] -03:149 [?] , 03:150, 03:151 [?] - 03:163 [?] , 03:165 [?] -03:167 [?] , 03:170 [?] - 03:179 [?] , 03:183 [?] -03:185 [?] , 03:187 [?] - 03:189 [?] , 03:192 [?] -03:194 [?] , 03:201 [?] - 03:205 [?] , 03:207 [?] -03:210 [?] , 03:212 [?] - 03:221 [?] , 03:223 [?] -03:230 [?] , 03:232 [?] - 03:238 [?] , 03:240 [?] , 03:244 [?] , 03:255 [?] - 03:262 [?] , 03:267 [?] -03:269N [?]	
48	145 bp	430 bp	*03:105-03:106, 03:109, 03:120	
49	150 bp	430 bp	*03:94, 03:164, 03:182	
	175 bp		*03:191	
50^{4,6}	105 bp	515 bp	*03:206	
	150 bp		*03:98, 03:165, 03:182	
51⁴	105 bp	430 bp	*03:90N, 03:240	
	140 bp		*03:213N, 03:239	
	165 bp		*03:115, 03:228	*06:17, 06:24, 06:30, 06:42, 06:109
52⁴	80 bp	430 bp	*03:188	
	185 bp		*03:99Q, 03:118N, 03:133, 03:162, 03:183, 03:200, 03:223, 03:233, 03:237N	*02:84, 04:14
53⁴	125 bp	430 bp	*03:130	
	155 bp		*03:139, 03:199	
	190 bp		*03:162, 03:183, 03:204, 03:223, 03:233	
54⁴	95 bp	430 bp	*03:91Q	
	170 bp		*03:146, 03:158, 03:170, 03:207	
	220 bp		*03:181	
55⁸	155 bp	430 bp	*03:102	
	195 bp		*03:114	
56	155 bp	430 bp	*03:95N, 03:210	
	200 bp		*03:150, 03:170	
57⁴	100 bp	515 bp	*03:196	*05:45
	215 bp		*03:169	
	240 bp		*03:246	

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

58⁸	155 bp 180 bp	430 bp	*03:102 *03:197Q	
59⁴	100 bp	430 bp	*03:10:01-03:10:02, 03:14:01-03:14:02, 03:138, 03:180, 03:183, 03:195	*06:01:01-06:01:15, 06:43, 06:45, 06:54N- 06:57, 06:98, 06:99:02- 06:105, 06:108, 06:120, 06:132, 06:140, 06:153, 06:157, 06:177, 06:181, 06:194, 06:205, 06:209, 06:214, 06:229, 06:239
60	185 bp	430 bp	*03:132, 03:215	*06:29, 06:123, 06:139
61	155 bp	430 bp	*03:221	*06:51:01-06:51:02
62	175 bp	430 bp	*03:50, 03:247	
63⁴	115 bp 255 bp	430 bp	*03:116 *03:243	
64⁹	-	-	Negative Control	

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of DQB1*03 SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³For several DQB1 alleles 1st and/or 3rd exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁴HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁵Primer mixes 3, 5, 23, 25, 30, 43 and 46 may have a tendency to giving rise to primer oligomer formation.

⁶Primer mixes 6, 11, 14, 16, 24 and 50 may have tendencies of unspecific amplifications.

⁷Primer mixes 1, 2, 18 and 27 may give rise to a lower yield of HLA-specific PCR product than the other DQB1*03 primer mixes.

⁸Primer mixes 55 and 58 may give rise to a long unspecific amplification product. This should be disregarded when interpreting the DQB1*03 typings.

⁹Primer mix 64 contains a negative control, which will amplify more than 95% of HLA amplicons as well as the amplicons generated by the control primer pairs matching the human growth hormone gene. HLA-specific PCR product sizes range from 75 to 200 base pairs and the PCR product generated by the HGH positive control primer pair is 430 base pairs.

'w', might be weakly amplified.

'?', nucleotide sequence information not available for the primer matching region.

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

PRIMER SPECIFICATION

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec. PCR product	170	220	135	220	135	95	115	110	65	135	135	260
Length of int. pos. control ¹	515	515	515	515	430	430	430	430	430	430	430	515
5'-primer(s) ²	45(230) 5'-ggA 3'	26(173) 5'-TTA 3'	26(173) 5'-TCT 3'	26(173) 5'-TCT 3'	26(173) 5'-TTA 3'	26(173) 5'-ggg 3'	38(210) 5'-gCA 3'	49(242) 5'-ggT 3'	38(210) 5'-gCA 3'	14(136) 5'-gCC 3'	13(134) 5'-ggC 3'	13(134) 5'-ggC 3'
						38(208) 5'-ACA 3'	140(514) 5'-CAA 3'	55(260) 5'-gCA 3'		135(500) 5'-TgA 3'		
								63(285) 5'-Agg 3'				
3'-primer(s) ³	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	57(266) 5'-Cgg 3'	86(353) 5'-gCT 3'	57(266) 5'-Cgg 3'	55(260) 5'-gCg 3'	62(282) 5'-CTA 3'	86(353) 5'-gCT 3'	46(232) 5'-CAA 3'	86(353) 5'-gCT 3'	45(230) 5'-CCC 3'	86(353) 5'-gCT 3'
	86(353) 5'-gTT 3'	87(357) 5'-CgT 3'	57(266) 5'-CAg 3'				163(583) 5'-CAC 3'		70(304) 5'-CCC 3'	167(596) 5'-CAT 3'	85(349) 5'-CAT 3'	87(357) 5'-CgT 3'
	87(357) 5'-CgT 3'		57(266) 5'-Cgg 3'				165(589) 5'-ggT 3'		75(320) 5'-CCC 3'			
							176(622) 5'-CTA 3'					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Well No.	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec. PCR product	115	130	135	130	165	175	175	150	160	160	95	120
	170				225			240			125	210
	205										150	
											185	
Length of int. pos. control ¹	430	430	515	430	430	515	430	430	515	515	430	430
5'-primer(s) ²	26(173) 5'-TTA 3'	26(173) 5'-TTA 3'	55(260) 5'-gCC 3'	57(266) 5'-TgA 3'	27(175) 5'-TTC 3'	140(514) 5'-CAA 3'	140(514) 5'-CAA 3'	104(406) 5'-CAC 3'	130(485) 5'-CCA 3'	9(122) 5'-gTA 3'	26(173) 5'-TCT 3'	26(173) 5'-TTA 3'
				57(266) 5'-TgA 3'	123(464) 5'-CTT 3'			133(493) 5'-TTT 3'			140(514) 5'-CAA 3'	102(400) 5'-TCT 3'
					142(520) 5'-gCC 3'							
					146(533) 5'-CCA 3'							
3'-primer(s) ³	51(248) 5'-gCA 3'	55(260) 5'-gCg 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	185(650) 5'-CgA 3'	185(650) 5'-Cgg 3'	170(604) 5'-gAC 3'	170(604) 5'-gAC 3'	48(240) 5'-gCg 3'	45(229) 5'-CCA 3'	81(338) 5'-TgC 3'
	67(296) 5'-ggT 3'	55(260) 5'-gCg 3'		87(357) 5'-CgT 3'	185(650) 5'-CgA 3'		185(650) 5'-CCg 3'				56(263) 5'-CAA 3'	83(343) 5'-gTC 3'
	71(308) 5'-ggA 3'	55(260) 5'-ggg 3'									74(318) 5'-AAg 3'	83(343) 5'-gTg 3'
	81(338) 5'-TgC 3'										167(595) 5'-ACA 3'	127(477) 5'-TTC 3'
											170(604) 5'-gAT 3'	
											176(622) 5'-CTA 3'	
Well No.	13	14	15	16	17	18	19	20	21	22	23	24

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

Well No.	25	26	27	28	29	30	31	32	33	34	35	36
Length of spec.	105	180	170	200	110	135	115	100	85	190	135	145
PCR product	130		225	230	165	175	195		160		190	175
	155		250		255		225				235	200
	180											
	240											
Length of int.	430	430	430	430	515	430	430	430	430	430	430	430
pos. control ¹												
5'-primer(s) ²	20(154) 5'-ACA 3'	9(121) 5'-TgC 3'	116(442) 5'-TgA 3'	24(166) 5'-gCC 3'	99(391) 5'-CAA 3'	10(125) 5'-CCg 3'	26(173) 5'-TTA 3'	13(134) 5'-ggC 3'	37(205) 5'-AgA 3'	26(173) 5'-TTA 3'	34(197) 5'-CCA 3'	34(196) 5'-ACT 3'
	37(207) 5'-TAA 3'	9(122) 5'-gTT 3'	123(464) 5'-CTT 3'	30(184) 5'-gAC 3'	129(481) 5'-AAA 3'	24(166) 5'-gCA 3'			62(282) 5'-AAg 3'	102(400) 5'-TCT 3'	52(251) 5'-gCT 3'	39(211) 5'-CAT 3'
	48(239) 5'-CCA 3'		141(517) 5'-CTA 3'	31(187) 5'-ACg 3'	145(530) 5'-CAT 3'						102(400) 5'-TCT 3'	41(217) 5'-TCC 3'
	140(514) 5'-CAA 3'		146(533) 5'-CCA 3'	36(204) 5'-gAC 3'	147(535) 5'-CCT 3'							52(251) 5'-gCT 3'
3'-primer(s) ³	86(353) 5'-gCT 3'	55(260) 5'-gCg 3'	185(650) 5'-CgA 3'	86(353) 5'-gCT 3'	170(604) 5'-gAC 3'	55(260) 5'-gCg 3'	51(248) 5'-gCA 3'	32(191) 5'-TAC 3'	77(326) 5'-CCg 3'	71(308) 5'-ggC 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'
	161(577) 5'-CAT 3'	55(260) 5'-gCg 3'				55(260) 5'-gCg 3'	77(325) 5'-CgC 3'			79(332) 5'-TgT 3'	130(484) 5'-CCA 3'	
	170(604) 5'-gAT 3'						83(343) 5'-gTC 3'			150(546) 5'-CCC 3'	150(544) 5'-gTC 3'	
	182(641) 5'-ggC 3'						88(359) 5'-TgA 3'				167(596) 5'-CAA 3'	
Well No.	25	26	27	28	29	30	31	32	33	34	35	36

Well No.	37	38	39	40	41	42	43	44	45	46	47	48
Length of spec.	180	170	105	180	175	80	175	160	90	115	175	145
PCR product	230	225	140	230	210	165	250	220	170	175		
			190			220		270	230			
			235									
Length of int.	430	430	430	430	430	430	430	430	430	430	430	430
pos. control ¹												
5'-primer(s) ²	23(164) 5'-gCA 3'	26(173) 5'-TTA 3'	26(173) 5'-TCT 3'	26(173) 5'-TCT 3'	27(175) 5'-ATA 3'	27(175) 5'-ATA 3'	102(400) 5'-TCT 3'	10(126) 5'-CAT 3'	23(163) 5'-AgA 3'	140(514) 5'-CAA 3'	140(516) 5'-ACC 3'	55(260) 5'-gCC 3'
	24(166) 5'-gCC 3'				34(196) 5'-ACT 3'	41(219) 5'-gAg 3'		26(173) 5'-TTT 3'	44(226) 5'-ACA 3'			
	39(211) 5'-CAT 3'				38(209) 5'-CgA 3'	47(235) 5'-TgA 3'		47(235) 5'-TgC 3'	70(306) 5'-AgC 3'			
	41(219) 5'-gAg 3'				44(227) 5'-CgC 3'	48(239) 5'-TCC 3'		48(239) 5'-TCA 3'				
						74(316) 5'-CgA 3'						
3'-primer(s) ³	86(353) 5'-gCT 3'	70(304) 5'-CCC 3'	48(238) 5'-CCA 3'	71(307) 5'-ggC 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	142(520) 5'-AAT 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	165(589) 5'-ggT 3'	185(650) 5'-CgA 3'	87(355) 5'-gAC 3'
		71(308) 5'-ggC 3'	56(263) 5'-CAA 3'	73(314) 5'-CCT 3'			146(532) 5'-ggC 3'			167(595) 5'-ACA 3'		88(358) 5'-gCC 3'
		71(308) 5'-ggA 3'	61(278) 5'-TCT 3'	76(324) 5'-gTT 3'			150(546) 5'-CCC 3'			184(647) 5'-Tgg 3'		93(374) 5'-gCA 3'
		88(359) 5'-TgA 3'	76(324) 5'-gTT 3'	89(361) 5'-CgC 3'			172(611) 5'-AgA 3'					
		89(362) 5'-TCA 3'	90(365) 5'-AgT 3'									
Well No.	37	38	39	40	41	42	43	44	45	46	47	48

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

Well No.	49	50	51	52	53	54	55	56	57	58	59	60
Length of spec.	150	105	105	80	125	95	155	155	100	155	100	185
PCR product	175	150	140	185	155	170	195	200	215	180		
			165		190	220			240			
Length of int. pos. control ¹	430	515	430	430	430	430	430	430	515	430	430	430
5'-primer(s) ²	140(514) 5'-CAA 3'	140(514) 5'-CAA 3'	26(173) 5'-TCT 3'	32(192) 5'-..T 3'	37(205) 5'-AgA 3'	8(119) 5'-CgC 3'	13(134) 5'-ggC 3'	10(125) 5'-CCg 3'	104(406) 5'-CAC 3'	13(134) 5'-ggC 3'	26(173) 5'-TTA 3'	9(122) 5'-gTT 3'
			135(500) 5'-TgA 3'	37(205) 5'-Agg 3'	37(205) 5'-Agg 3'	13(133) 5'-AgA 3'		16(142) 5'-gCC 3'	112(430) 5'-ACT 3'			
				37(207) 5'-TAA 3'	48(239) 5'-CCC 3'	16(142) 5'-gCC 3'		29(182) 5'-CAA 3'	151(547) 5'-ACA 3'			
				44(226) 5'-ACA 3'	48(239) 5'-TCA 3'	125(470) 5'-Agg 3'						
				74(316) 5'-CgA 3'	58(268) 5'-ACT 3'	167(595) 5'-AgA 3'						
3'-primer(s) ³	170(604) 5'-gAA 3'	161(577) 5'-CAT 3'	48(238) 5'-Cg 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	55(260) 5'-gCg 3'	52(250) 5'-CgC 3'	66(294) 5'-ACT 3'	170(604) 5'-gAC 3'	52(250) 5'-CgC 3'	45(230) 5'-CCC 3'	57(266) 5'-Cgg 3'
	174(617) 5'-CgC 3'	172(611) 5'-AgC 3'	48(238) 5'-CCA 3'			185(650) 5'-CgA 3'	64(287) 5'-TCC 3'			60(274) 5'-gTA 3'		
	177(625) 5'-gTT 3'	177(625) 5'-gTT 3'	55(260) 5'-gCg 3'									
	184(647) 5'-Tgg 3'	179(632) 5'-ggT 3'	64(286) 5'-CTA 3'									
			71(307) 5'-ggC 3'									
			175(619) 5'-CAT 3'									
Well No.	49	50	51	52	53	54	55	56	57	58	59	60

Well No.	61	62	63
Length of spec.	155	175	115
PCR product			255
Length of int. pos. control ¹	430	430	430
5'-primer(s) ²	28(179) 5'-gAC 3'	102(400) 5'-TCT 3'	99(391) 5'-CAA 3'
		145(530) 5'-CAA 3'	
3'-primer(s) ³	66(294) 5'-ATg 3'	146(532) 5'-ggT 3'	170(604) 5'-gAC 3'
		146(532) 5'-ggC 3'	
Well No.	61	62	63

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.



101.214-24/06 – including Taq polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without Taq polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

CELL LINE VALIDATION SHEET																				
DQB1*03 SSP subtyping kit ²																				
				Well																
				1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
				Production No.	201784901	201784902	201784903	201784904	201784905	201784906	201784907	201784908	201784909	201784910	201784911	201784912	201784913	201784914	201784915	201784916
	IHWC cell line ¹	DQB1																		
1	9001 SA	*05:01		- -																
2	9280 LK707	*06:01	*02:02	- - + + -																
3	9011 E4181324	*06:01		- -																
4	9275 GU373	*02:01		- - + + -																
5	9009 KAS011	*05:02		- -																
6	9353 SM	*03:02	*06:01	- - + + -																
7	9020 QBL	*02:01		- - + + -																
8	9025 DEU	*03:01		+ + -																
9	9026 YAR	*03:02		- - + + -																
10	9107 LKT3	*04:01		- -																
11	9051 PITOUT	*02:02		- - + + -																
12	9052 DBB	*03:03		- - - - + -																
13	9004 JESTHOM	*05:01		- -																
14	9071 OLGA	*04:02		- -																
15	9075 DKB	*03:03		- - - - + -																
16	9037 SWEIG007	*03:01		+ + -																
17	9282 CTM3953540	*02:01	*06:03	- - + + -																
18	9257 32367	*06:02	*02:02	- - + + -																
19	9038 BM16	*03:01		+ + -																
20	9059 SLE005	*06:04		- -																
21	9064 AMALA	*03:01		+ + -																
22	9056 KOSE	*05:03	*06:04	- -																
23	9124 IHL	*05:03	*06:01	- -																
24	9035 JBUSH	*03:01		+ + -																
25	9049 IBW9	*02:02		- - + + -																
26	9285 WT49	*02:01		- - + + -																
27	9191 CH1007	*04:01	*05:01	- -																
28	9320 BEL5GB	*02:02	*03:01	+ + + + -																
29	9050 MOU	*02:02		- - + + -																
30	9021 RSH	*04:02		- -																
31	9019 DUCAF	*02:01		- - + + -																
32	9297 HAG	*03:01		+ + -																
33	9098 MT14B	*03:02		- - + + -																
34	9104 DHIF	*03:01		+ + -																
35	9302 SSTO	*03:05		- - - - - - - - + - - - - - - - - - - - - - - - - - - -																
36	9024 KT17	*03:02		- - + + -																
37	9065 HHKB	*06:03		- -																
38	9099 LZL	*03:01		+ + -																
39	9315 CML	*02:01	*03:01	+ + + + -																
40	9134 WHONP199	*02:02	*03:03	- - + + -																
41	9055 H0301	*06:09		- -																
42	9066 TAB089	*06:01		- -																
43	9076 T7526	*03:03		- - - - + -																
44	9057 TEM	*05:03		- -																
45	9239 SHJO	*02:02		- - + + -																
46	9013 SCHU	*06:02		- -																
47	9045 TUBO	*03:01		+ + -																
48	9303 TER-ND	*05:01		- -																



101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

CELL LINE VALIDATION SHEET																				
DQB1*03 SSP subtyping kit ²																				
				Well																
				17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
				Production No.	201784917	201784918	201784919	201784920	201784921	201784922	201784923	201784924	201784925	201784926	201784927	201784928	201784929	201784930	201784931	201784932
IHCW cell line ¹		DQB1																		
1	9001 SA	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
2	9280 LK707	*06:01	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
3	9011 E4181324	*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
4	9275 GU373	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
5	9009 KAS011	*05:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
6	9353 SM	*03:02	*06:01	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
7	9020 QBL	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
8	9025 DEU	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
9	9026 YAR	*03:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
10	9107 LKT3	*04:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
11	9051 PITOUT	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
12	9052 DBB	*03:03		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
13	9004 JESTHOM	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
14	9071 OLGA	*04:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
15	9075 DKB	*03:03		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
16	9037 SWEIG007	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
17	9282 CTM3953540	*02:01	*06:03	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
18	9257 32367	*06:02	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
19	9038 BM16	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
20	9059 SLE005	*06:04		-	-	-	-	W	-	-	-	-	-	-	-	-	-	-	-	
21	9064 AMALA	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
22	9056 KOSE	*05:03	*06:04	-	-	-	-	W	-	-	-	-	-	-	-	-	-	-	-	
23	9124 IHL	*05:03	*06:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
24	9035 JBUSH	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
25	9049 IBW9	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
26	9285 WT49	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
27	9191 CH1007	*04:01	*05:01	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
28	9320 BEL5GB	*02:02	*03:01	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
29	9050 MOU	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
30	9021 RSH	*04:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
31	9019 DUCAF	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
32	9297 HAG	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
33	9098 MT14B	*03:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
34	9104 DHIF	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
35	9302 SSTO	*03:05		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
36	9024 KT17	*03:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
37	9065 HHKB	*06:03		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	
38	9099 LZL	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
39	9315 CML	*02:01	*03:01	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
40	9134 WHONP199	*02:02	*03:03	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
41	9055 H0301	*06:09		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
42	9066 TAB089	*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
43	9076 T7526	*03:03		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
44	9057 TEM	*05:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
45	9239 SHJO	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
46	9013 SCHU	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
47	9045 TUBO	*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
48	9303 TER-ND	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
 101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

CELL LINE VALIDATION SHEET																				
DQB1*03 SSP subtyping kit ²																				
				Well																
				33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	
				Production No.	201784933	201784934	201784935	201784936	201784937	201784938	201784939	201784940	201784941	201784942	201784943	201784944	201784945	201784946	201784947	201784948
IHWC cell line ¹		DQB1																		
1	9001 SA	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707	*06:01	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011 E4181324	*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275 GU373	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011	*05:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM	*03:02	*06:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
7	9020 QBL	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026 YAR	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075 DKB	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282 CTM3953540	*02:01	*06:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257 32367	*06:02	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005	*06:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE	*05:03	*06:04	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL	*05:03	*06:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035 JBUSH	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007	*04:01	*05:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320 BEL5GB	*02:02	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019 DUCAF	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO	*03:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB	*06:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML	*02:01	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134 WHONP199	*02:02	*03:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055 H0301	*06:09		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089	*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM	*05:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303 TER-ND	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-



101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
 101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

CELL LINE VALIDATION SHEET																			
DQB1*03 SSP subtyping kit ²																			
				Well															
				49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	
				Production No.	201784949	201784950	201784951	201784952	201784953	201784954	201784955	201784956	201784957	201784958	201784959	201784960	201784961	201784962	201784963
IHCW cell line ¹		DQB1																	
1	9001 SA	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
2	9280 LK707	*06:01	*02:02	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
3	9011 E4181324	*06:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
4	9275 GU373	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
5	9009 KAS011	*05:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
6	9353 SM	*03:02	*06:01	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
7	9020 QBL	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
8	9025 DEU	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
9	9026 YAR	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
10	9107 LKT3	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
11	9051 PITOUT	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
12	9052 DBB	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
13	9004 JESTHOM	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
14	9071 OLGA	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
15	9075 DKB	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
16	9037 SWEIG007	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
17	9282 CTM3953540	*02:01	*06:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
18	9257 32367	*06:02	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
19	9038 BM16	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
20	9059 SLE005	*06:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
21	9064 AMALA	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
22	9056 KOSE	*05:03	*06:04	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
23	9124 IHL	*05:03	*06:01	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
24	9035 JBUSH	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
25	9049 IBW9	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
26	9285 WT49	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
27	9191 CH1007	*04:01	*05:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
28	9320 BEL5GB	*02:02	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
29	9050 MOU	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
30	9021 RSH	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
31	9019 DUCAF	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
32	9297 HAG	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
33	9098 MT14B	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
34	9104 DHIF	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
35	9302 SSTO	*03:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
36	9024 KT17	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
37	9065 HHKB	*06:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
38	9099 LZL	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
39	9315 CML	*02:01	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
40	9134 WHONP199	*02:02	*03:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
41	9055 H0301	*06:09		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
42	9066 TAB089	*06:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	
43	9076 T7526	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
44	9057 TEM	*05:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
45	9239 SHJO	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
46	9013 SCHU	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
47	9045 TUBO	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
48	9303 TER-ND	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
 101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution of the kit has been tested against 48 well characterized IHWC cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 7 to 10, 13, 17, 20, 23 to 46, 48 to 58 and 60 to 63 were available.

The specificities of the primers in primer solutions 7, 9, 10, 23, 25 to 28, 35, 38, 39, 52 to 54, 60 and 61 were tested by separately adding additional 5'-primers respectively 3'-primers.

In primer solutions 8, 17, 20, 29, 30, 33, 36, 37, 41, 42, 44, 45, 56, 57 and 63 it was only possible to test the 3'-primers, the 5'-primers were not possible to test. In primer solutions 13, 24, 31, 32, 34, 40, 43, 46, 48 to 51, 55, 58 and 62 it was only possible to test the 5'-primers, the 3'-primers were not possible to test.

In primer solutions 1 to 3, 7, 9 to 12, 14, 16, 19, 23, 25, 35, 38 and 39 one or more of the 3'-primers were not possible to test, and in primer solution 6, 16, 25 to 28, 35 and 52 to 54 one or more of the 5'-primers were not possible to test.

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
 101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **6F6**

Lot-specific information

ADDRESSES:

Manufacturer:

Caredx AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-508 939 00

Fax: +46-8-717 88 18

E-mail: orders-se@caredx.com

Web page: <https://labproducts.caredx.com/>

Distributed by:

CareDx GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

Tel: +43-1-710 15 00

Fax: +43-1-710 15 00 10

E-mail: orders-at@caredx.com

Web page: <https://labproducts.caredx.com/>

CareDx Lab Solutions Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-653-78171

Fax: 610-344-7989

E-mail: orders-us@caredx.com

Web page: <https://labproducts.caredx.com/>

For information on CareDx distributors worldwide, contact **CareDx GmbH**.