

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

Olerup SSP[®] DQB1*03

Product number:	101.214-24/06 – including <i>Taq</i> pol. 101.214-24u/06u – without <i>Taq</i> pol.
Lot number:	9K1
Expiry date:	2024-03-01
Number of tests:	24 test – Product No. 101.214-24/24u 6 tests – Product No. 101.214-06/06u
Number of wells per test:	63+1
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
- Adhesive PCR seals	RT
- Product Insert	RT

This Product Description is only valid for Lot No. 9K1.

Complete product documentation consists of generic Instructions for Use (IFU), lot specific Product Insert, Worksheet and Certificate

**CHANGES COMPARED TO THE PREVIOUS OLERUP SSP[®]
DQB1*03 LOT (7H0)**

- The product documentation has been updated for new alleles of IMGT 3.38.0.
- The kit resolution focuses on common and well documented (CWD) alleles¹.

- ¹As described in section Uniquely Identified Alleles

The DQB1*03 primer set, specificity and interpretation tables have been updated for the HLA-DQB1 alleles described since the previous *Olerup SSP[®] DQB1*03* lot was made (**Lot No. 7H0**).

¹S. J. Mack, P. Cano, J. A. Hollenbach et al.
Common and well-documented HLA alleles: 2012 update to the CWD catalogue. Tissue Antigens, 2013, 81, 194–203

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
28	Added	-	5'-primer added for improved yield of the DQB1*03:38:01-03:38:02.

Change in revision R01 compared to R00:

1. Primer mixes 23 and 25 do not amplify the DQB1*03:24 and 03:79 alleles. The corrections above have been implemented in the Specificity and Interpretation tables.

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

Well **64** contains Negative Control primer pairs, that will amplify a majority of the Olerup SSP[®] HLA Class I, DRB, DQB1, DPB1 and DQA1 amplicons as well as all the amplicons generated by the control primer pairs matching the human growth hormone gene.

HLA-specific PCR product sizes range from 75 to 200 base pairs.
The PCR product generated by the positive control primer pair is 200 base pairs.

Length of PCR product	105	200	105	80	75	80	85
5'-primer¹	164	340	440	45	45	43	36
	5'-CAC ^{3'}	5'-Agg ^{3'}	5'-TTA ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-TAC ^{3'}
							36
							5'-TAT ^{3'}
3'-primer²	231	2nd I	507	59	58	57	47
	5'-TgC ^{3'}	5'-AAA ^{3'}	5'-TTg ^{3'}	5'-CTC ^{3'}	5'-ggC ^{3'}	5'-CTC ^{3'}	5'-ACA ^{3'}
							48
							5'-gCA ^{3'}
							48
							5'-gCC ^{3'}
							52
							5'-TgT ^{3'}
A*	+	+	+				
B*	+	+	+				
C*	+	+	+				
DRB1				+	+		
DRB3				+	+		
DRB5				+			
DQB1					+		
DPB1						+	
DQA1							+

¹The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

²The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

PRODUCT DESCRIPTION

DQB1*03 SSP subtyping

CONTENT

The primer set contains 5'- and 3'-primers for identifying the DQB1*03:01 to DQB1*03:412 alleles.

Please note that DQB1 amplifications usually are somewhat less pronounced than e.g. DRB and DQA1 amplifications even when using the same DNA preparation and exactly the same experimental procedures.

PLATE LAYOUT

Each test consists of 64 PCR reactions in a 64 well PCR plate.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48
49	50	51	52	53	54	55	56
57	58	59	60	61	62	63	NC

The 64 well cut PCR plate is marked with ‘DQB1*03’ in silver/gray ink.

Well No. 1 is marked with the Lot No. ‘9K1’.

Wells 1 to 63 – DQB1*03 high resolution primers.

Well 64 – Negative Control (NC).

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded. The PCR plates are covered with a PCR-compatible foil.

INTERPRETATION

Due to the sharing of sequence motifs, non-DQB1*03 alleles will be amplified by some primer mixes. The interpretation of DQB1*03 subtypings is not influenced by the DQB2 and DQB3 genes. For further details see Specificity Table.

UNIQUELY IDENTIFIED ALLELES

All the DQB1*03 alleles, i.e. **DQB1*03:01 to DQB1*03:412**, recognized by the HLA Nomenclature Committee in October 2019^{1,2} will be amplified by the primers in the DQB1*03 subtyping kit³.

The DQB1*03 kit enables separation of the confirmed DQB1*03 alleles as listed in the IMGT/HLA database 3.29.0. An HLA allele is listed as confirmed by IMGT/HLA

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: 9K1

Lot-specific information

if it has been sequenced by more than a single laboratory or from multiple sources. Current allele confirmation status for DQB1*03 alleles is listed below.

The DQB1*03 kit also enables identification of many null and alternatively expressed alleles.

The following DQB1*03 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix
DQB1*03:40, 03:137	9
DQB1*03:169, 03:196	57
DQB1*03:190, 03:287	35

¹HLA-DQB1 alleles listed on the IMGT/HLA web page 2019-October-17, release 3.38.0, www.ebi.ac.uk/imgt/hla.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

³The following alleles cannot be separated by the HLA-DQB1*03 primer set. These alleles can be distinguished by the DQ low resolution kit and/or the HLA-DQB1*04 high resolution kit.

Alleles

DQB1*03:25:01, DQB1*04:03:03
DQB1*03:132, DQB1*04:75

RESOLUTION IN HOMO- AND HETEROZYGOTES

Results file with resolution in DQB1*03 homo- and heterozygotes is available upon request.

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

ALLELE CONFIRMATION STATUS

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DQB1*03:01:01:01	Confirmed	DQB1*03:01:35	Unconfirmed	DQB1*03:04:01	Confirmed	DQB1*03:43	Unconfirmed
DQB1*03:01:01:02	Unconfirmed	DQB1*03:01:36	Unconfirmed	DQB1*03:04:02	Unconfirmed	DQB1*03:44	Confirmed
DQB1*03:01:01:03	Confirmed	DQB1*03:01:37	Unconfirmed	DQB1*03:04:03	Unconfirmed	DQB1*03:45	Confirmed
DQB1*03:01:01:04	Unconfirmed	DQB1*03:01:38	Unconfirmed	DQB1*03:05:01	Confirmed	DQB1*03:46	Unconfirmed
DQB1*03:01:01:05	Confirmed	DQB1*03:01:39	Unconfirmed	DQB1*03:05:02	Unconfirmed	DQB1*03:47	Confirmed
DQB1*03:01:01:06	Unconfirmed	DQB1*03:02:01:01	Confirmed	DQB1*03:05:03	Confirmed	DQB1*03:48	Unconfirmed
DQB1*03:01:01:07	Confirmed	DQB1*03:02:01:02	Confirmed	DQB1*03:05:04	Unconfirmed	DQB1*03:49	Unconfirmed
DQB1*03:01:01:08	Confirmed	DQB1*03:02:01:03	Unconfirmed	DQB1*03:06	Unconfirmed	DQB1*03:50	Unconfirmed
DQB1*03:01:01:09	Unconfirmed	DQB1*03:02:01:04	Unconfirmed	DQB1*03:07	Unconfirmed	DQB1*03:51	Unconfirmed
DQB1*03:01:01:10	Unconfirmed	DQB1*03:02:02	Unconfirmed	DQB1*03:08	Unconfirmed	DQB1*03:52	Confirmed
DQB1*03:01:01:11	Unconfirmed	DQB1*03:02:03	Unconfirmed	DQB1*03:09	Confirmed	DQB1*03:53	Confirmed
DQB1*03:01:01:12	Unconfirmed	DQB1*03:02:04	Unconfirmed	DQB1*03:10:01	Unconfirmed	DQB1*03:54	Confirmed
DQB1*03:01:01:14	Unconfirmed	DQB1*03:02:05	Unconfirmed	DQB1*03:10:02	Confirmed	DQB1*03:55	Confirmed
DQB1*03:01:01:15	Unconfirmed	DQB1*03:02:06	Confirmed	DQB1*03:11	Confirmed	DQB1*03:56	Confirmed
DQB1*03:01:01:16	Unconfirmed	DQB1*03:02:07	Confirmed	DQB1*03:12	Unconfirmed	DQB1*03:57	Unconfirmed
DQB1*03:01:01:17	Unconfirmed	DQB1*03:02:08	Unconfirmed	DQB1*03:13	Confirmed	DQB1*03:58	Confirmed
DQB1*03:01:01:18	Unconfirmed	DQB1*03:02:09	Confirmed	DQB1*03:14:01	Unconfirmed	DQB1*03:59	Confirmed
DQB1*03:01:02	Unconfirmed	DQB1*03:02:10	Confirmed	DQB1*03:14:02	Confirmed	DQB1*03:60	Unconfirmed
DQB1*03:01:03	Confirmed	DQB1*03:02:11	Confirmed	DQB1*03:15	Unconfirmed	DQB1*03:61	Unconfirmed
DQB1*03:01:04	Confirmed	DQB1*03:02:12	Confirmed	DQB1*03:16	Confirmed	DQB1*03:62	Confirmed
DQB1*03:01:05	Unconfirmed	DQB1*03:02:13	Confirmed	DQB1*03:17:01	Confirmed	DQB1*03:63	Confirmed
DQB1*03:01:06	Unconfirmed	DQB1*03:02:14	Unconfirmed	DQB1*03:17:02	Confirmed	DQB1*03:64	Confirmed
DQB1*03:01:07	Unconfirmed	DQB1*03:02:15	Confirmed	DQB1*03:18	Unconfirmed	DQB1*03:65	Confirmed
DQB1*03:01:08	Confirmed	DQB1*03:02:16	Confirmed	DQB1*03:19:01	Confirmed	DQB1*03:66N	Unconfirmed
DQB1*03:01:09	Unconfirmed	DQB1*03:02:17	Unconfirmed	DQB1*03:19:02	Unconfirmed	DQB1*03:67	Confirmed
DQB1*03:01:10	Unconfirmed	DQB1*03:02:18	Unconfirmed	DQB1*03:20	Unconfirmed	DQB1*03:68	Confirmed
DQB1*03:01:11	Unconfirmed	DQB1*03:02:19	Confirmed	DQB1*03:21	Unconfirmed	DQB1*03:69	Unconfirmed
DQB1*03:01:12	Unconfirmed	DQB1*03:02:20	Unconfirmed	DQB1*03:22	Unconfirmed	DQB1*03:70	Confirmed
DQB1*03:01:13	Unconfirmed	DQB1*03:02:21	Confirmed	DQB1*03:23:01	Unconfirmed	DQB1*03:71	Unconfirmed
DQB1*03:01:14	Confirmed	DQB1*03:02:22	Unconfirmed	DQB1*03:23:02	Unconfirmed	DQB1*03:72	Confirmed
DQB1*03:01:15	Unconfirmed	DQB1*03:02:23	Unconfirmed	DQB1*03:24	Unconfirmed	DQB1*03:73	Unconfirmed
DQB1*03:01:16	Unconfirmed	DQB1*03:02:24	Unconfirmed	DQB1*03:25:01	Confirmed	DQB1*03:74	Unconfirmed
DQB1*03:01:17	Confirmed	DQB1*03:02:25	Unconfirmed	DQB1*03:25:02	Unconfirmed	DQB1*03:75	Confirmed
DQB1*03:01:18	Unconfirmed	DQB1*03:03:02:01	Confirmed	DQB1*03:26	Confirmed	DQB1*03:76	Unconfirmed
DQB1*03:01:19	Unconfirmed	DQB1*03:03:02:02	Unconfirmed	DQB1*03:27	Confirmed	DQB1*03:77	Unconfirmed
DQB1*03:01:20	Unconfirmed	DQB1*03:03:02:03	Unconfirmed	DQB1*03:28	Unconfirmed	DQB1*03:78	Unconfirmed
DQB1*03:01:21	Confirmed	DQB1*03:03:02:04	Confirmed	DQB1*03:29	Unconfirmed	DQB1*03:79	Unconfirmed
DQB1*03:01:22	Confirmed	DQB1*03:03:03	Unconfirmed	DQB1*03:30	Unconfirmed	DQB1*03:80	Unconfirmed
DQB1*03:01:23	Unconfirmed	DQB1*03:03:04	Confirmed	DQB1*03:31	Unconfirmed	DQB1*03:81	Confirmed
DQB1*03:01:24	Unconfirmed	DQB1*03:03:05	Confirmed	DQB1*03:32	Unconfirmed	DQB1*03:82	Unconfirmed
DQB1*03:01:25	Unconfirmed	DQB1*03:03:06	Unconfirmed	DQB1*03:33	Unconfirmed	DQB1*03:83	Unconfirmed
DQB1*03:01:26	Unconfirmed	DQB1*03:03:07	Confirmed	DQB1*03:34	Unconfirmed	DQB1*03:84N	Unconfirmed
DQB1*03:01:27	Confirmed	DQB1*03:03:08	Confirmed	DQB1*03:35	Unconfirmed	DQB1*03:85	Unconfirmed
DQB1*03:01:28	Unconfirmed	DQB1*03:03:09	Unconfirmed	DQB1*03:36	Unconfirmed	DQB1*03:86	Unconfirmed
DQB1*03:01:29	Unconfirmed	DQB1*03:03:10	Unconfirmed	DQB1*03:37	Confirmed	DQB1*03:87	Unconfirmed
DQB1*03:01:30	Unconfirmed	DQB1*03:03:11	Confirmed	DQB1*03:38	Unconfirmed	DQB1*03:88	Unconfirmed
DQB1*03:01:31	Unconfirmed	DQB1*03:03:12	Unconfirmed	DQB1*03:39	Unconfirmed	DQB1*03:89	Unconfirmed
DQB1*03:01:32	Unconfirmed	DQB1*03:03:13	Unconfirmed	DQB1*03:40	Confirmed	DQB1*03:90N	Unconfirmed
DQB1*03:01:33	Unconfirmed	DQB1*03:03:14	Unconfirmed	DQB1*03:41	Unconfirmed	DQB1*03:91Q	Unconfirmed
DQB1*03:01:34	Unconfirmed	DQB1*03:03:15	Unconfirmed	DQB1*03:42	Unconfirmed	DQB1*03:92	Unconfirmed

101.214-24/06 – including *Taq* polymerase, IFU-01
101.214-24u/06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DQB1*03:93	Unconfirmed	DQB1*03:143	Unconfirmed	DQB1*03:193	Unconfirmed	DQB1*03:243	Confirmed
DQB1*03:94	Unconfirmed	DQB1*03:144	Unconfirmed	DQB1*03:194	Confirmed	DQB1*03:244	Unconfirmed
DQB1*03:95N	Unconfirmed	DQB1*03:145	Unconfirmed	DQB1*03:195	Confirmed	DQB1*03:245	Confirmed
DQB1*03:96	Unconfirmed	DQB1*03:146	Unconfirmed	DQB1*03:196	Confirmed	DQB1*03:246	Confirmed
DQB1*03:97	Unconfirmed	DQB1*03:147	Unconfirmed	DQB1*03:197Q	Confirmed	DQB1*03:247	Confirmed
DQB1*03:98	Unconfirmed	DQB1*03:148	Unconfirmed	DQB1*03:198	Unconfirmed	DQB1*03:248	Unconfirmed
DQB1*03:99Q	Unconfirmed	DQB1*03:149	Unconfirmed	DQB1*03:199	Confirmed	DQB1*03:249	Unconfirmed
DQB1*03:100	Unconfirmed	DQB1*03:150	Confirmed	DQB1*03:200	Confirmed	DQB1*03:250	Unconfirmed
DQB1*03:101	Unconfirmed	DQB1*03:151	Unconfirmed	DQB1*03:201	Confirmed	DQB1*03:251	Unconfirmed
DQB1*03:102	Confirmed	DQB1*03:152	Unconfirmed	DQB1*03:202	Unconfirmed	DQB1*03:252	Unconfirmed
DQB1*03:103	Unconfirmed	DQB1*03:153	Unconfirmed	DQB1*03:203	Unconfirmed	DQB1*03:253	Unconfirmed
DQB1*03:104	Unconfirmed	DQB1*03:154	Unconfirmed	DQB1*03:204	Confirmed	DQB1*03:254	Unconfirmed
DQB1*03:105	Confirmed	DQB1*03:155	Unconfirmed	DQB1*03:205	Confirmed	DQB1*03:255	Confirmed
DQB1*03:106	Confirmed	DQB1*03:156	Unconfirmed	DQB1*03:206	Unconfirmed	DQB1*03:256	Unconfirmed
DQB1*03:107	Unconfirmed	DQB1*03:157	Unconfirmed	DQB1*03:207	Confirmed	DQB1*03:257	Unconfirmed
DQB1*03:108	Unconfirmed	DQB1*03:158	Confirmed	DQB1*03:208	Unconfirmed	DQB1*03:258	Unconfirmed
DQB1*03:109	Confirmed	DQB1*03:159	Unconfirmed	DQB1*03:209	Unconfirmed	DQB1*03:259	Unconfirmed
DQB1*03:110	Unconfirmed	DQB1*03:160	Unconfirmed	DQB1*03:210	Unconfirmed	DQB1*03:260	Unconfirmed
DQB1*03:111	Unconfirmed	DQB1*03:161	Unconfirmed	DQB1*03:211	Confirmed	DQB1*03:261	Unconfirmed
DQB1*03:112	Unconfirmed	DQB1*03:162	Confirmed	DQB1*03:212	Confirmed	DQB1*03:262	Unconfirmed
DQB1*03:113	Confirmed	DQB1*03:163	Unconfirmed	DQB1*03:213N	Unconfirmed	DQB1*03:263	Unconfirmed
DQB1*03:114	Confirmed	DQB1*03:164	Unconfirmed	DQB1*03:214	Unconfirmed	DQB1*03:264	Unconfirmed
DQB1*03:115	Unconfirmed	DQB1*03:165	Unconfirmed	DQB1*03:215	Unconfirmed	DQB1*03:265	Unconfirmed
DQB1*03:116	Unconfirmed	DQB1*03:166	Unconfirmed	DQB1*03:216	Unconfirmed	DQB1*03:266	Unconfirmed
DQB1*03:117	Unconfirmed	DQB1*03:167	Unconfirmed	DQB1*03:217	Unconfirmed	DQB1*03:267	Unconfirmed
DQB1*03:118N	Confirmed	DQB1*03:168	Unconfirmed	DQB1*03:218	Unconfirmed	DQB1*03:268	Unconfirmed
DQB1*03:119	Confirmed	DQB1*03:169	Unconfirmed	DQB1*03:219	Unconfirmed	DQB1*03:269N	Unconfirmed
DQB1*03:120	Unconfirmed	DQB1*03:170	Confirmed	DQB1*03:220	Unconfirmed		
DQB1*03:121	Unconfirmed	DQB1*03:171	Unconfirmed	DQB1*03:221	Confirmed		
DQB1*03:122	Unconfirmed	DQB1*03:172	Unconfirmed	DQB1*03:222	Unconfirmed		
DQB1*03:123	Unconfirmed	DQB1*03:173	Unconfirmed	DQB1*03:223	Unconfirmed		
DQB1*03:124	Unconfirmed	DQB1*03:174	Unconfirmed	DQB1*03:224	Unconfirmed		
DQB1*03:125	Unconfirmed	DQB1*03:175	Unconfirmed	DQB1*03:225	Confirmed		
DQB1*03:126	Unconfirmed	DQB1*03:176	Unconfirmed	DQB1*03:226	Unconfirmed		
DQB1*03:127	Unconfirmed	DQB1*03:177	Confirmed	DQB1*03:227	Unconfirmed		
DQB1*03:128	Unconfirmed	DQB1*03:178	Unconfirmed	DQB1*03:228	Confirmed		
DQB1*03:129	Unconfirmed	DQB1*03:179	Unconfirmed	DQB1*03:229	Unconfirmed		
DQB1*03:130	Confirmed	DQB1*03:180	Unconfirmed	DQB1*03:230	Unconfirmed		
DQB1*03:131	Unconfirmed	DQB1*03:181	Unconfirmed	DQB1*03:231	Unconfirmed		
DQB1*03:132	Confirmed	DQB1*03:182	Unconfirmed	DQB1*03:232	Confirmed		
DQB1*03:133	Unconfirmed	DQB1*03:183	Unconfirmed	DQB1*03:233	Confirmed		
DQB1*03:134	Unconfirmed	DQB1*03:184	Confirmed	DQB1*03:234	Unconfirmed		
DQB1*03:135	Unconfirmed	DQB1*03:185	Unconfirmed	DQB1*03:235	Unconfirmed		
DQB1*03:136	Unconfirmed	DQB1*03:186	Unconfirmed	DQB1*03:236	Unconfirmed		
DQB1*03:137	Unconfirmed	DQB1*03:187	Unconfirmed	DQB1*03:237N	Unconfirmed		
DQB1*03:138	Confirmed	DQB1*03:188	Confirmed	DQB1*03:238	Unconfirmed		
DQB1*03:139	Confirmed	DQB1*03:189	Unconfirmed	DQB1*03:239	Confirmed		
DQB1*03:140	Unconfirmed	DQB1*03:190	Unconfirmed	DQB1*03:240	Confirmed		
DQB1*03:141	Unconfirmed	DQB1*03:191	Confirmed	DQB1*03:241	Unconfirmed		
DQB1*03:142	Unconfirmed	DQB1*03:192	Unconfirmed	DQB1*03:242	Unconfirmed		

¹Allele status “confirmed” or “unconfirmed” as listed on the IMGT/HLA web page 2017-August-10, release 3.29.0, www.ebi.ac.uk/imgt/hla.



101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information
SPECIFICITY TABLE

DQB1*03 SSP subtyping

Specificities and sizes of the PCR products of the 63+1 primer mixes used for DQB1*03 SSP subtyping

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified DQB1*03 alleles ³	Other amplified DQB1 alleles
1 ⁷	170 bp	515 bp	*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:48, 03:04:01-03:04:04, 03:09, 03:13, 03:16, 03:19:01:01-03:19:04, 03:21-03:22:02, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:60, 03:69, 03:71, 03:73-03:78, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108, 03:113-03:116, 03:118N-03:122, 03:127-03:131, 03:134-03:135, 03:139-03:140, 03:142-03:144, 03:147-03:148, 03:150-03:152, 03:157-03:160, 03:162-03:167, 03:169-03:173, 03:182, 03:184, 03:188, 03:191-03:194, 03:196-03:198:02, 03:201-03:202, 03:206-03:208, 03:216, 03:218-03:219, 03:231-03:232, 03:235-03:236, 03:241-03:243, 03:246, 03:252-03:257, 03:260, 03:264, 03:266-03:268, 03:271, 03:275-03:276N, 03:281, 03:284-03:286, 03:288, 03:290-03:294, 03:297, 03:302-03:303N, 03:305-03:307, 03:309, 03:311-03:312, 03:314, 03:317:01-03:318, 03:326-03:331, 03:338N, 03:340N-03:342, 03:347, 03:350, 03:353, 03:355, 03:358N, 03:360-03:361, 03:366, 03:370, 03:372-03:373, 03:376N-03:378, 03:380-03:382, 03:385N, 03:387, 03:389-03:391, 03:394, 03:396, 03:399N-03:401, 03:404, 03:407N-03:408	*02:152, 04:09
2 ⁷	220 bp	515 bp	*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:07, 03:01:09-03:01:48, 03:04:01-03:04:04, 03:09-03:10:03, 03:13-03:14:02, 03:16, 03:19:01:01-03:19:04, 03:21-03:22:02, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:60, 03:69, 03:71, 03:73, 03:75-03:77, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108, 03:114-03:116, 03:118N-03:122, 03:127-03:131, 03:133-03:135, 03:138-03:140, 03:142-03:144, 03:147-03:148, 03:150, 03:152, 03:157-03:160, 03:162-03:167, 03:169-03:173, 03:180, 03:182-03:183, 03:186-03:188,	

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

			03:191-03:198:02, 03:201-03:202, 03:206-03:208, 03:216, 03:218-03:219, 03:231-03:232, 03:235-03:236, 03:241- 03:243, 03:246, 03:252-03:257, 03:260, 03:264, 03:266-03:268, 03:271, 03:275- 03:276N, 03:281, 03:284-03:286, 03:288, 03:290-03:294, 03:297, 03:302-03:303N, 03:305-03:307, 03:309, 03:311-03:312, 03:317:01-03:318, 03:326-03:331, 03:338N, 03:340N-03:342, 03:347, 03:350, 03:353-03:355, 03:358N, 03:360- 03:361, 03:366, 03:370, 03:372-03:373, 03:376N-03:378, 03:380-03:381, 03:385N, 03:387, 03:389-03:391, 03:394, 03:396, 03:399N-03:401, 03:404, 03:407N-03:408
3⁵	135 bp	515 bp	*03:02:01:01-03:02:31, 03:07-03:08, 03:11, 03:18, 03:32, 03:37, 03:45:01- 03:45:02, 03:62-03:64, 03:66N-03:68, 03:70, 03:81, 03:85, 03:106-03:107, 03:125, 03:146, 03:153, 03:161, 03:174- 03:175, 03:178-03:179, 03:184-03:185, 03:189-03:190, 03:199, 03:203-03:205, 03:210-03:211, 03:213N-03:215, 03:220- 03:221, 03:223-03:225, 03:228-03:229, 03:233, 03:237N, 03:240, 03:245, 03:247, 03:251, 03:261, 03:263, 03:265, 03:269N, 03:273-03:274, 03:277-03:279, 03:287, 03:289, 03:295-03:296, 03:298-03:301, 03:308, 03:310N, 03:315, 03:320-03:324, 03:333-03:334N, 03:339N, 03:343- 03:345, 03:348-03:349, 03:352, 03:362, 03:364, 03:367-03:369, 03:371, 03:379, 03:383, 03:386, 03:388, 03:392, 03:403N, 03:409-03:410, 03:412
4	220 bp	515 bp	*03:02:01:01-03:03:26, 03:06-03:08, 03:11-03:12, 03:15, 03:18, 03:20, 03:23:01-03:23:03, 03:25:01-03:26, 03:30-03:34, 03:37-03:41, 03:43, 03:45:01-03:45:02, 03:62-03:68, 03:70, 03:74, 03:79, 03:81, 03:85-03:91Q, 03:95N-03:99Q, 03:104, 03:106-03:107, 03:110-03:113, 03:117, 03:123-03:126, 03:136-03:137, 03:141, 03:145-03:146, 03:149, 03:153, 03:155-03:156, 03:161, 03:168, 03:174-03:179, 03:184-03:185, 03:190, 03:199-03:200, 03:203-03:205, 03:209-03:215, 03:217, 03:220-03:225, 03:227-03:230, 03:233-03:234, 03:237N- 03:240, 03:245, 03:247-03:249, 03:251, 03:259, 03:261, 03:263, 03:265, 03:269N- 03:270, 03:273-03:274, 03:277-03:280, 03:282N-03:283, 03:287, 03:289, 03:295- 03:296, 03:298-03:301, 03:304, 03:308, 03:310N, 03:313, 03:315-03:316,

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

			03:319-03:324, 03:332-03:334N, 03:336-03:337, 03:339N, 03:343-03:345, 03:348-03:349, 03:351-03:352, 03:356N-03:357N, 03:359, 03:362-03:365, 03:367-03:369, 03:371, 03:374-03:375N, 03:379, 03:382-03:384, 03:386, 03:388, 03:392-03:393, 03:395, 03:397-03:398, 03:402-03:403N, 03:405-03:406, 03:409-03:412	
5⁵	135 bp	430 bp	*03:04:01-03:04:04, 03:14:01-03:14:02, 03:80, 03:138, 03:318, 03:327	*06:246
6^{4,6}	95 bp 130 bp	430 bp	*03:20 *03:05:01-03:05:05, 03:17:01-03:17:02, 03:61, 03:72, 03:100, 03:132, 03:181, 03:226, 03:244, 03:250, 03:262, 03:272, 03:346	*04:09, 04:63, 04:75, 05:131
7⁴	115 bp	430 bp	*03:03:06, 03:06, 03:25:01-03:25:02, 03:42, 03:88	*04:01:03, 04:02:16, 04:03:03
	150 bp		*03:84N	
8⁴	110 bp 145 bp	430 bp	*03:15 *03:07, 03:16	
9⁴	65 bp 145 bp	430 bp	*03:40, 03:353 *03:08, 03:119, 03:137, 03:194	*06:02:02, 06:03:02, 06:03:34
10	135 bp 260 bp	430 bp	*03:09 *03:11, 03:26, 03:244	
11⁶	135 bp	430 bp	*03:10:01-03:10:02:02, 03:12, 03:14:01-03:14:02, 03:70, 03:179, 03:183, 03:195	*06:01:01:01-06:01:06, 06:01:08-06:01:26, 06:06 [?] , 06:43, 06:54N-06:58, 06:98, 06:99:02-06:105, 06:108, 06:120, 06:132, 06:140, 06:153:01-06:153:02, 06:157, 06:167-06:168, 06:177, 06:181, 06:194, 06:205, 06:209, 06:214, 06:229, 06:239, 06:243, 06:245-06:247, 06:251, 06:257-06:258, 06:260, 06:263, 06:268, 06:274, 06:277, 06:285, 06:302, 06:305, 06:307, 06:309-06:310, 06:312, 06:321, 06:323, 06:330N, 06:342, 06:359
	255 bp		*03:47	
12	260 bp	515 bp	*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:07, 03:01:09-03:01:48, 03:04:01-03:04:04, 03:09-03:10:03, 03:12-03:14:02, 03:16, 03:19:01:01-03:19:04, 03:21-03:22:02, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:60, 03:70-03:71, 03:73-03:78, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108, 03:114-03:116, 03:118N-03:122,	

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

			03:127-03:131, 03:133-03:135, 03:139-03:140, 03:143-03:144, 03:147-03:148, 03:150-03:152, 03:157-03:160, 03:162-03:167, 03:169-03:173, 03:179, 03:182-03:183, 03:186-03:188, 03:191-03:198:02, 03:201-03:202, 03:206, 03:208, 03:216, 03:218-03:219, 03:231-03:232, 03:235-03:236, 03:241-03:243, 03:246, 03:252-03:257, 03:260, 03:264, 03:266-03:268, 03:271, 03:275-03:276N, 03:281, 03:284-03:286, 03:288, 03:290-03:294, 03:297, 03:302-03:303N, 03:305-03:307, 03:309, 03:311-03:312, 03:314, 03:317:01-03:318, 03:326-03:331, 03:338N, 03:340N-03:342, 03:347, 03:350, 03:353-03:355, 03:358N, 03:360-03:361, 03:366, 03:370, 03:372-03:373, 03:376N-03:378, 03:380-03:381, 03:385N, 03:387, 03:389-03:391, 03:394, 03:396, 03:399N-03:401, 03:404, 03:407N-03:408
13⁴	115 bp 170 bp 205 bp	430 bp	*03:255 *03:13, 03:232 *03:48
14⁶	130 bp	430 bp	*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:12, 03:01:14, 03:01:16-03:01:24, 03:01:26-03:01:48, 03:04:01-03:04:04, 03:09-03:10:03, 03:13-03:14:02, 03:19:01:01-03:19:04, 03:21-03:22:02, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:60, 03:69, 03:71, 03:73, 03:75-03:77, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108-03:109, 03:114-03:116, 03:118N-03:122, 03:127-03:129, 03:131, 03:133-03:135, 03:138-03:140, 03:142-03:144, 03:147-03:148, 03:150, 03:152, 03:154, 03:157-03:160, 03:162-03:167, 03:169-03:171, 03:173, 03:180, 03:182-03:183, 03:186-03:188, 03:191-03:196, 03:198:01-03:198:02, 03:201, 03:206-03:208, 03:216, 03:218-03:219, 03:231-03:232, 03:235-03:236, 03:241-03:243, 03:246, 03:252-03:257, 03:260, 03:264, 03:266-03:268, 03:271, 03:275-03:276N, 03:281, 03:284-03:286, 03:288, 03:290-03:294, 03:297, 03:302-03:303N, 03:305-03:307, 03:309, 03:311-03:312, 03:317:01-03:318, 03:326-03:331, 03:335, 03:338N, 03:340N-03:342, 03:347, 03:350, 03:353-03:354N, 03:358N, 03:360-03:361, 03:366, 03:370, 03:372-03:373, 03:376N-03:378, 03:380-03:381, 03:385N, 03:387, 03:389-03:391, 03:394, 03:396, 03:399N-03:401, 03:404, 03:407N-03:408

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

Lot No.	135 bp	515 bp	Lot-specific information	
15	135 bp	515 bp	*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:07, 03:01:09-03:01:33, 03:01:35-03:01:36, 03:01:38-03:03:22, 03:03:24-03:15, 03:17:01-03:22:02, 03:24, 03:26-03:103, 03:106-03:108, 03:110-03:152, 03:155-03:171, 03:173-03:188, 03:190-03:204, 03:206-03:213N, 03:215-03:216, 03:218-03:233, 03:235-03:257, 03:260, 03:262-03:281, 03:283-03:324, 03:326-03:332, 03:334N, 03:336-03:338N, 03:340N-03:354N, 03:356N-03:386, 03:388-03:412	*04:09, 04:63, 04:75
16 ⁶	130 bp	430 bp	*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:02, 03:01:04-03:01:07, 03:01:09-03:01:12, 03:01:14, 03:01:16-03:01:24, 03:01:26-03:01:48, 03:03:02:01-03:03:26, 03:06, 03:09-03:10:03, 03:12-03:13, 03:15-03:17:02, 03:19:01:01-03:24, 03:26-03:31, 03:33-03:36, 03:38:01-03:44, 03:46-03:60, 03:65, 03:69, 03:71-03:79, 03:82-03:84N, 03:86-03:103, 03:108, 03:111-03:124, 03:126-03:129, 03:131, 03:133-03:137, 03:139-03:145, 03:147-03:152, 03:155-03:160, 03:162-03:171, 03:173, 03:176-03:177, 03:180, 03:182-03:183, 03:186-03:188, 03:191-03:198:02, 03:200-03:201, 03:206-03:209, 03:212, 03:216-03:219, 03:222, 03:227, 03:230-03:232, 03:235-03:236, 03:238, 03:241-03:244, 03:246, 03:248-03:249, 03:252-03:257, 03:260, 03:264, 03:266-03:268, 03:270-03:272, 03:275-03:276N, 03:280-03:286, 03:288, 03:290-03:294, 03:297, 03:302-03:307, 03:309, 03:311-03:314, 03:316-03:317:02, 03:319, 03:326, 03:328-03:332, 03:336-03:338N, 03:340N-03:342, 03:347, 03:350-03:351, 03:353-03:361, 03:363, 03:365-03:366, 03:370, 03:372-03:378, 03:380-03:382, 03:384-03:385N, 03:387, 03:389-03:391, 03:393-03:402, 03:404-03:408, 03:411N	*02:03:01-02:03:02, 02:77, 04:09-04:10, 04:50, 04:63
17	165 bp 225 bp	430 bp	*03:43, 03:87 *03:18, 03:85	
18 ⁷	175 bp	515 bp	*03:02:01:01-03:03:26, 03:05:01-03:05:05, 03:06 [?] -03:07 [?] , 03:08, 03:11 [?] , 03:14:01 [?] -03:15 [?] , 03:17:01 [?] -03:18 [?] , 03:19:01:01-03:19:04, 03:20 [?] , 03:23:01 [?] , 03:23:02, 03:23:03 [?] , 03:25:01, 03:25:02 [?] -03:26 [?] , 03:30-03:32, 03:33 [?] -03:34 [?] , 03:37 [?] , 03:38:01-03:39, 03:40 [?] , 03:41, 03:43, 03:45:01-03:45:02, 03:48 [?] , 03:52 [?] -03:54 [?] , 03:56 [?] -03:71 [?] , 03:72, 03:74 [?] -03:78 [?] , 03:79, 03:81 [?] -03:82 [?] ,	*04:01:01:01-04:02:01:01, 04:02:01:04-04:44, 04:46N-04:77

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

			03:85-03:91Q, 03:95N-03:100, 03:101 [?] -03:104 [?] , 03:105, 03:106 [?] -03:112 [?] , 03:117, 03:118N [?] , 03:120 [?] -03:131 [?] , 03:132, 03:133 [?] -03:135 [?] , 03:136, 03:137 [?] , 03:138, 03:140 [?] -03:149 [?] , 03:150, 03:151 [?] -03:163 [?] , 03:165 [?] -03:167 [?] , 03:168, 03:170 [?] -03:179 [?] , 03:180-03:181, 03:183 [?] , 03:184, 03:185 [?] , 03:187 [?] -03:189 [?] , 03:190, 03:192 [?] -03:194 [?] , 03:195, 03:199-03:200, 03:201 [?] -03:205 [?] , 03:207 [?] -03:210 [?] , 03:211, 03:212 [?] -03:221 [?] , 03:222, 03:223 [?] -03:230 [?] , 03:232 [?] -03:238 [?] , 03:239, 03:240 [?] , 03:244 [?] , 03:245, 03:247-03:248, 03:250-03:251, 03:255 [?] -03:262 [?] , 03:263, 03:265, 03:267 [?] -03:270 [?] , 03:272, 03:273 [?] -03:275 [?] , 03:277, 03:278 [?] , 03:279-03:280, 03:282N-03:283, 03:286-03:287, 03:289, 03:295-03:296, 03:298-03:301, 03:304, 03:308, 03:310N, 03:313, 03:315-03:316, 03:319-03:324, 03:332-03:337, 03:339N, 03:343-03:346, 03:348-03:349, 03:351 [?] -03:353 [?] , 03:356N [?] , 03:357N, 03:359 [?] , 03:362-03:365, 03:367-03:369, 03:371, 03:374-03:376N, 03:379, 03:381-03:384, 03:386, 03:387 [?] , 03:388, 03:392-03:395, 03:398, 03:401-03:403N, 03:405-03:406, 03:409-03:412
19 ⁶	175 bp	515 bp	*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:48, 03:04:01-03:04:04, 03:06 [?] -03:07 [?] , 03:09-03:10:03, 03:11 [?] , 03:12-03:13, 03:14:01 [?] -03:15 [?] , 03:16, 03:17:01 [?] -03:18 [?] , 03:20 [?] , 03:21-03:22:02, 03:23:01 [?] , 03:23:03 [?] , 03:24, 03:25:02 [?] -03:26 [?] , 03:27-03:29, 03:33 [?] -03:34 [?] , 03:35-03:36, 03:37 [?] , 03:40 [?] , 03:42, 03:44, 03:46-03:47, 03:48 [?] , 03:49-03:51, 03:52 [?] -03:54 [?] , 03:55, 03:56 [?] -03:71 [?] , 03:73, 03:74 [?] -03:78 [?] , 03:80, 03:81 [?] -03:82 [?] , 03:83-03:84N, 03:92-03:94, 03:101 [?] -03:104 [?] , 03:106 [?] -03:112 [?] , 03:113-03:116, 03:118N [?] , 03:119, 03:120 [?] -03:131 [?] , 03:133 [?] -03:135 [?] , 03:137 [?] , 03:139, 03:140 [?] -03:149 [?] , 03:151 [?] -03:163 [?] , 03:164, 03:165 [?] -03:167 [?] , 03:169, 03:170 [?] -03:179 [?] , 03:182, 03:183 [?] , 03:185 [?] , 03:186, 03:187 [?] -03:189 [?] , 03:191, 03:192 [?] -03:194 [?] , 03:196-03:198:02, 03:201 [?] -03:205 [?] , 03:206, 03:207 [?] -03:210 [?] , 03:212 [?] -03:221 [?] , 03:223 [?] -03:230 [?] , 03:231, 03:232 [?] -03:238 [?] , 03:240 [?] , 03:241-03:243, 03:244 [?] , 03:246, 03:252-03:254, 03:255 [?] -03:262 [?] , 03:264, 03:266, 03:267 [?] -03:270 [?] , 03:271, 03:273 [?] -03:275 [?] ,

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

			03:276N, 03:278 [?] , 03:281, 03:284- 03:285, 03:288, 03:290-03:294, 03:297, 03:302-03:303N, 03:305-03:307, 03:309, 03:311, 03:314, 03:317:01-03:318, 03:326-03:331, 03:338N, 03:340N- 03:342, 03:347, 03:350, 03:351 [?] -03:353 [?] , 03:354N-03:355, 03:356N [?] , 03:359 [?] , 03:361, 03:366, 03:372-03:373, 03:377- 03:378, 03:380, 03:385N, 03:387 [?] , 03:389-03:391, 03:396, 03:399N- 03:400N, 03:404, 03:407N-03:408	
20	150 bp	430 bp	*03:21	*02:64, 02:79, 05:35, 05:182, 06:243, 06:255 *06:310
	240 bp		*03:246	
21	160 bp	515 bp	*03:22:01-03:22:02, 03:96, 03:251	*02:108, 02:144, 04:23, 05:103 ^w , 05:226, 05:231 ^w , 06:04:01:01 ^w - 06:05:02 ^w , 06:06 [?] , 06:07:01 ^w -06:07:02 ^w , 06:08:02 [?] -06:08:03 [?] , 06:09:01:01 ^w -06:09:11 ^w , 06:09:12 [?] , 06:11:02 [?] - 06:11:03 [?] , 06:12 ^w , 06:13:02 [?] -06:13:03 [?] , 06:15:01 [?] -06:15:02 [?] , 06:17 [?] , 06:18:01 ^w - 06:18:02 ^w , 06:19:01 [?] , 06:21 ^w , 06:22:01 [?] - 06:25 [?] , 06:27:01 [?] - 06:32:02 [?] , 06:34 ^w , 06:35 [?] , 06:36 ^w , 06:38 ^w , 06:40 [?] , 06:41 ^w -06:42 ^w , 06:45 [?] , 06:50 [?] , 06:51:02 [?] , 06:52 ^w , 06:53:01 [?] -06:54N [?] , 06:56 [?] -06:57 [?] , 06:58 ^w , 06:59 [?] -06:61 [?] , 06:63 [?] - 06:67 [?] , 06:69:01 [?] - 06:76 [?] , 06:79:02 [?] - 06:83 [?] , 06:84 ^w , 06:85 [?] , 06:89 [?] , 06:91 [?] - 06:92:01 [?] , 06:93 [?] - 06:97 [?] , 06:118:01 ^w - 06:118:03 ^w , 06:119 [?] - 06:124 [?] , 06:126 [?] , 06:128 [?] , 06:129 ^w , 06:131 [?] -06:147 [?] , 06:149 [?] -06:153:01 [?] , 06:154 [?] -06:184 [?] , 06:186 ^w , 06:189 ^w , 06:190:01 [?] -06:199 [?] , 06:201 [?] -06:204 [?] , 06:206:01 [?] , 06:206:02 ^w , 06:207 [?] -06:208 [?] , 06:210 [?] -06:215 [?] ,

101.214-24/06 – including *Taq* polymerase, IFU-01
101.214-24u/06u – without *Taq* polymerase, IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

				06:217 ^w , 06:229 [?] - 06:236 [?] , 06:238 [?] , 06:241 [?] , 06:244 [?] , 06:246 [?] -06:247 [?] , 06:249 [?] , 06:252N [?] , 06:254 ^w , 06:260 [?] , 06:261 ^w , 06:263 [?] , 06:265 ^w -06:267 ^w , 06:268 [?] , 06:275 ^w , 06:277 [?] , 06:280 ^w - 06:283 ^w , 06:287 ^w - 06:288 ^w , 06:291 ^w - 06:292 ^w , 06:299 ^w , 06:301 ^w , 06:303N ^w , 06:313 ^w , 06:318 ^w - 06:319 ^w , 06:320 [?] , 06:325 ^w , 06:332 ^w , 06:339 ^w , 06:343 ^w , 06:348 ^w -06:349 ^w , 06:351 ^w , 06:353 [?] , 06:358 [?]
22	160 bp	515 bp	*03:23:01	*05:03:12, 05:43:02, 06:03:01:01-06:03:03, 06:03:05-06:03:07, 06:03:11-06:03:25, 06:03:27-06:03:33, 06:03:35-06:03:38, 06:04:02, 06:07:01, 06:08:01, 06:09:02, 06:11:01-06:11:03, 06:26N, 06:28, 06:30- 06:32:01, 06:40-06:41, 06:44, 06:59-06:62, 06:64-06:65, 06:67, 06:90-06:92:02, 06:94, 06:110, 06:118:01, 06:118:03, 06:128, 06:133-06:134, 06:141, 06:143-06:145, 06:148, 06:154, 06:165, 06:169- 06:170, 06:184-06:185, 06:187, 06:191, 06:195- 06:196, 06:199, 06:203, 06:210, 06:218, 06:221- 06:223, 06:230, 06:234, 06:238, 06:248, 06:250, 06:253, 06:259, 06:272, 06:276, 06:279, 06:316, 06:327-06:329, 06:331, 06:334, 06:336, 06:345N-06:346, 06:350, 06:352
23 ^{4,5}	95 bp 125 bp 150 bp	430 bp	*03:41, 03:224 *03:25:02, 03:89, 03:234, 03:239 *03:84N	*04:03:01-04:03:03

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

24 ^{4,6}	185 bp		*03:45:01	
	220 bp		*03:258	
	120 bp	430 bp	*03:39	
	210 bp		*03:27, 03:48, 03:56	
25 ^{4,5}	105 bp	430 bp	*03:206	
	155 bp		*03:80, 03:257	
	180 bp		*03:29, 03:118N, 03:237N	
	240 bp		*03:34	*02:14:01
26	180 bp	430 bp	*03:30, 03:52, 03:72, 03:100, 03:132, 03:179, 03:215, 03:270	*04:09, 04:63, 04:75, 06:29, 06:96:01- 06:96:02, 06:139, 06:209, 06:322, 06:337
27 ⁷	170 bp	430 bp	*03:32, 03:87	
	225 bp		*03:85	
	250 bp		*03:31	
28	200 bp	430 bp	*03:37-03:38:02, 03:46, 03:252, 03:277	
	230 bp		*03:225	
29 ⁴	110 bp	515 bp	*03:33, 03:245, 03:328	*06:228, 06:304N
	165 bp		*03:35, 03:315	
	255 bp		*03:243	*02:126
30 ⁵	135 bp	515 bp	*03:36, 03:149	
	175 bp		*03:150	
31 ⁴	115 bp	430 bp	*03:255	
	195 bp		*03:28	
	225 bp		*03:56, 03:83	
32 ⁴	100 bp	430 bp	*03:44	
33 ^{4,7}	85 bp	430 bp	*03:59	*06:37
	160 bp		*03:204	
34 ⁵	190 bp	430 bp	*03:58, 03:82, 03:198:01, 03:222	
35	135 bp	430 bp	*03:190, 03:205, 03:214	*02:157, 04:15, 04:42, 05:102, 06:223
	190 bp		*03:51, 03:54, 03:287	*04:35
36	235 bp		*03:117	
	145 bp	430 bp	*03:205, 03:214	*04:15
	175 bp		*03:55, 03:60, 03:250	
37	200 bp		*03:269N	
	180 bp	430 bp	*03:60, 03:77	
38	230 bp		*03:53, 03:62, 03:225	
	170 bp	430 bp	*03:58, 03:194, 03:232	*06:82, 06:233, 06:302
39 ⁴	225 bp		*03:57, 03:83	
	105 bp	430 bp	*03:240	
40	140 bp		*03:25:01-03:25:02, 03:66N, 03:234, 03:239	*04:03:01-04:03:03, 06:193N
	190 bp		*03:68	
	235 bp		*03:86	*02:151
	180 bp	430 bp	*03:67-03:68, 03:228	*06:17, 06:24, 06:30, 06:42, 06:318
41	230 bp		*03:104, 03:212	*02:47
	135 bp	430 bp	*03:282N	
41	175 bp		*03:63, 03:75, 03:168	
	210 bp		*03:201, 03:269N	

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

42⁴	80 bp 165 bp 220 bp	430 bp	*03:188 *03:64, 03:76-03:77 *03:201	
43⁵	175 bp 250 bp	430 bp	*03:49, 03:198:01, 03:222, 03:247, 03:285-03:286, 03:402 *03:97, 03:329, 03:349, 03:381	*02:10, 02:81, 05:152, 06:359 *06:86, 06:104
44	160 bp 220 bp 270 bp	430 bp	*03:199, 03:211 *03:78 *03:71	
45^{4,5}	90 bp 135 bp 170 bp 230 bp	430 bp	*03:65 *03:282N *03:133, 03:200 *03:81, 03:159	*02:84, 04:14 *02:116
46^{4,5}	115 bp 175 bp	515 bp	*03:88-03:89 *03:191	
47	175 bp	430 bp	*03:06 [?] -03:07 [?] , 03:11 [?] , 03:14:01 [?] -03:15 [?] , 03:17:01 [?] -03:18 [?] , 03:19:01:01-03:19:04, 03:20 [?] , 03:23:01 [?] , 03:23:03 [?] , 03:25:02 [?] - 03:26 [?] , 03:37 [?] , 03:40 [?] , 03:48 [?] , 03:52 [?] - 03:54 [?] , 03:56 [?] -03:71 [?] , 03:74 [?] -03:78 [?] , 03:81 [?] -03:82 [?] , 03:101 [?] -03:104 [?] , 03:106 [?] - 03:112 [?] , 03:118N [?] , 03:120 [?] -03:131 [?] , 03:133 [?] -03:135 [?] , 03:137 [?] , 03:140 [?] - 03:149 [?] , 03:150, 03:151 [?] -03:163 [?] , 03:165 [?] -03:167 [?] , 03:170 [?] -03:179 [?] , 03:183 [?] , 03:185 [?] , 03:187 [?] -03:189 [?] , 03:192 [?] -03:194 [?] , 03:201 [?] -03:205 [?] , 03:207 [?] -03:210 [?] , 03:212 [?] -03:221 [?] , 03:223 [?] -03:230 [?] , 03:232 [?] -03:238 [?] , 03:240 [?] , 03:244 [?] , 03:255 [?] -03:262 [?] , 03:267 [?] -03:270 [?] , 03:273 [?] -03:275 [?] , 03:278 [?] , 03:286, 03:335, 03:351 [?] - 03:353 [?] , 03:356N [?] , 03:359 [?] , 03:376N, 03:381, 03:387 [?] , 03:394, 03:401	
48⁵	145 bp	430 bp	*03:105-03:106, 03:109, 03:120	
49	150 bp 175 bp	430 bp	*03:94, 03:164, 03:182 *03:191	*04:44
50^{4,6}	105 bp 150 bp	515 bp	*03:206 *03:98, 03:165, 03:182	
51⁴	105 bp 140 bp 165 bp 220 bp	430 bp	*03:90N, 03:240 *03:239 *03:115, 03:213N, 03:228, 03:368 *03:258	*02:106, 04:43, 06:17, 06:24, 06:30, 06:42, 06:109
52⁴	80 bp 185 bp	430 bp	*03:188 *03:99Q, 03:118N, 03:133, 03:162, 03:183, 03:200, 03:223, 03:233, 03:237N, 03:351	*02:84, 04:14
53⁴	125 bp 155 bp 190 bp	430 bp	*03:130, 03:271 *03:139, 03:199 *03:162, 03:183, 03:204, 03:223, 03:233, 03:351	

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

54⁴	95 bp 170 bp 220 bp	430 bp	*03:91Q *03:146, 03:158, 03:170, 03:207 *03:181	
55⁸	155 bp 195 bp	430 bp	*03:102 *03:114	
56	155 bp 200 bp	430 bp	*03:95N, 03:210 *03:150, 03:170	
57⁴	100 bp 215 bp 240 bp	515 bp	*03:196 *03:169 *03:246	*05:45, 05:201, 06:323 *06:327 *06:310
58⁸	155 bp 180 bp	430 bp	*03:102 *03:197Q	
59⁴	100 bp	430 bp	*03:10:01-03:10:02:02, 03:14:01- 03:14:02, 03:138, 03:180, 03:183, 03:195	*06:01:01:01-06:01:26, 06:43, 06:45, 06:54N- 06:57, 06:98, 06:99:02- 06:105, 06:108, 06:120, 06:132, 06:140, 06:153:01-06:153:02, 06:157, 06:177, 06:181, 06:194, 06:205, 06:209, 06:214, 06:229, 06:239, 06:243, 06:245-06:247, 06:251, 06:257-06:258, 06:263, 06:268, 06:274, 06:277, 06:285, 06:302, 06:305, 06:307, 06:309- 06:310, 06:312, 06:321, 06:323, 06:330N, 06:359
60	185 bp	430 bp	*03:132, 03:215	*04:75, 06:29, 06:123, 06:139, 06:337
61^{4,5}	115 bp 155 bp 255 bp	430 bp	*03:116 *03:221 *03:243	*06:51:01-06:51:02 *02:126
62	175 bp	430 bp	*03:50, 03:247, 03:285	
63	160 bp	430 bp	*03:276N	
64⁹	-	-	Negative Control	

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of DQB1*03 SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³For several DQB1 alleles 1st and/or 3rd exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁴HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁵Primer mixes 3, 5, 23, 25, 30, 43, 45, 46, 48 and 61 may have a tendency to giving rise to primer oligomer formation.

⁶Primer mixes 6, 11, 14, 16, 19, 24 and 50 may have tendencies of unspecific amplifications.

⁷Primer mixes 1, 2, 18, 27 och 33 may give rise to a lower yield of HLA-specific PCR product than the other DQB1*03 primer mixes.

⁸Primer mixes 55 and 58 may give rise to a long unspecific amplification product at approximately 1000 bp. This should be disregarded when interpreting the DQB1*03 typings.

⁹Primer mix 64 contains a negative control, which will amplify the majority of the HLA amplicons as well as the amplicons generated by the control primer pairs matching the human growth hormone gene. HLA-specific PCR product sizes range from 75 to 200 base pairs and the PCR product generated by the HGH positive control primer pair is 200 base pairs.

Abbreviations

‘w’, might be weakly amplified.

‘?’, nucleotide sequence information not available for the primer matching region.

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

PRIMER SPECIFICATION

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec.	170	220	135	220	135	95	115	110	65	135	135	260
PCR product						130	150	145	145	260	255	
Length of int. pos. control ¹	515	515	515	515	430	430	430	430	430	430	430	515
5'-primer(s) ²	45(230) 5'-ggA 3'	26(173) 5'-TTA 3'	26(173) 5'-TCT 3'	26(173) 5'-TCT 3'	26(173) 5'-TTA 3'	26(173) 5'-ggg 3'	38(210) 5'-gCA 3'	49(242) 5'-ggT 3'	38(210) 5'-gCA 3'	14(136) 5'-gCC 3'	13(134) 5'-ggC 3'	13(134) 5'-ggC 3'
						38(208) 5'-ACA 3'	140(514) 5'-CAA 3'	55(260) 5'-gCA 3'		135(500) 5'-TgA 3'		
								63(285) 5'-Agg 3'				
3'-primer(s) ³	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	57(266) 5'-Cgg 3'	86(353) 5'-gCT 3'	57(266) 5'-Cgg 3'	55(260) 5'-gCg 3'	62(282) 5'-CTA 3'	86(353) 5'-gCT 3'	46(232) 5'-CAA 3'	86(353) 5'-gCT 3'	45(230) 5'-CCC 3'	86(353) 5'-gCT 3'
	86(353) 5'-gTT 3'	87(357) 5'-CgT 3'	57(266) 5'-CAG 3'	87(356) 5'-ggg 3'			163(583) 5'-CAC 3'		70(304) 5'-CCC 3'	167(596) 5'-CAT 3'	85(349) 5'-CAT 3'	87(357) 5'-CgT 3'
	87(357) 5'-CgT 3'		57(266) 5'-Cgg 3'	87(357) 5'-CgA 3'			165(589) 5'-ggT 3'		75(320) 5'-CCC 3'			
							176(622) 5'-CTA 3'					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Well No.	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec.	115	130	135	130	165	175	175	150	160	160	95	120
PCR product	170				225			240			125	210
	205										150	
											185	
											220	
Length of int. pos. control ¹	430	430	515	430	430	515	515	430	515	515	430	430
5'-primer(s) ²	26(173) 5'-TTA 3'	26(173) 5'-TTA 3'	55(260) 5'-gCC 3'	57(266) 5'-TgA 3'	27(175) 5'-TTC 3'	140(514) 5'-CAA 3'	140(514) 5'-CAA 3'	104(406) 5'-CAC 3'	130(485) 5'-CCA 3'	9(122) 5'-gTA 3'	26(173) 5'-TCT 3'	26(173) 5'-TTA 3'
				57(266) 5'-TgA 3'	123(464) 5'-CTT 3'			133(493) 5'-TTT 3'			140(514) 5'-CAA 3'	102(400) 5'-TCT 3'
					142(520) 5'-gCC 3'							
					146(533) 5'-CCA 3'							
3'-primer(s) ³	51(248) 5'-gCA 3'	55(260) 5'-gCg 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	185(650) 5'-CgA 3'	185(650) 5'-Cgg 3'	170(604) 5'-gAC 3'	170(604) 5'-gAC 3'	48(240) 5'-gCg 3'	45(229) 5'-CCA 3'	81(338) 5'-TgC 3'
	67(296) 5'-ggT 3'	55(260) 5'-gCg 3'	87(357) 5'-CgA 3'	87(357) 5'-CgT 3'	185(650) 5'-CgA 3'		185(650) 5'-CCg 3'				56(263) 5'-CAA 3'	83(343) 5'-gTC 3'
	71(308) 5'-ggA 3'	55(260) 5'-ggg 3'		87(357) 5'-CgA 3'							74(318) 5'-AAg 3'	83(343) 5'-gTg 3'
	81(338) 5'-TgC 3'										86(354) 5'-Agg 3'	127(477) 5'-TTC 3'
											167(595) 5'-ACA 3'	
											170(604) 5'-gAT 3'	
											176(622) 5'-CTA 3'	
Well No.	13	14	15	16	17	18	19	20	21	22	23	24

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

Well No.	25	26	27	28	29	30	31	32	33	34	35	36
Length of spec.	105	180	170	200	110	135	115	100	85	190	135	145
PCR product	155		225	230	165	175	195		160		190	175
	180		250		255		225				235	200
	240											
Length of int.	430	430	430	430	515	515	430	430	430	430	430	430
pos. control ¹												
5'-primer(s) ²	20(154) 5'-ACA 3'	9(121) 5'-TgC 3'	116(442) 5'-TgA 3'	24(166) 5'-gCC 3'	99(391) 5'-CAA 3'	10(125) 5'-CCg 3'	26(173) 5'-TTA 3'	13(134) 5'-ggC 3'	37(205) 5'-AgA 3'	26(173) 5'-TTA 3'	34(197) 5'-CCA 3'	34(196) 5'-ACT 3'
	37(207) 5'-TAA 3'	9(122) 5'-gTT 3'	123(464) 5'-CTT 3'	30(184) 5'-gAC 3'	129(481) 5'-AAA 3'	24(166) 5'-gCA 3'			62(282) 5'-AAg 3'	102(400) 5'-TCT 3'	52(251) 5'-gCT 3'	39(211) 5'-CAT 3'
	48(239) 5'-CCA 3'		141(517) 5'-CTA 3'	30(184) 5'-gAC 3'	145(530) 5'-CAT 3'						102(400) 5'-TCT 3'	41(217) 5'-TCC 3'
	140(514) 5'-CAA 3'		146(533) 5'-CCA 3'	31(187) 5'-ACg 3'	147(535) 5'-CCT 3'							52(251) 5'-gCT 3'
				36(204) 5'-gAC 3'								
3'-primer(s) ³	86(353) 5'-gCT 3'	55(260) 5'-gCg 3'	185(650) 5'-CgA 3'	86(353) 5'-gCT 3'	170(604) 5'-gAC 3'	55(260) 5'-gCg 3'	51(248) 5'-gCA 3'	32(191) 5'-TAC 3'	77(326) 5'-CCg 3'	71(308) 5'-ggC 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'
	161(577) 5'-CAT 3'	55(260) 5'-gCg 3'				55(260) 5'-gCg 3'	77(325) 5'-CgC 3'			79(332) 5'-TgT 3'	130(484) 5'-CCA 3'	
	170(604) 5'-gAT 3'						83(343) 5'-gTC 3'			150(546) 5'-CCC 3'	150(544) 5'-gTC 3'	
	182(641) 5'-ggC 3'						88(359) 5'-TgA 3'				167(596) 5'-CAA 3'	
Well No.	25	26	27	28	29	30	31	32	33	34	35	36

Well No.	37	38	39	40	41	42	43	44	45	46	47	48
Length of spec.	180	170	105	180	135	80	175	160	90	115	175	145
PCR product	230	225	140	230	175	165	250	220	135	175		
			190		210	220		270	170			
			235						230			
Length of int.	430	430	430	430	430	430	430	430	430	515	430	430
pos. control ¹												
5'-primer(s) ²	23(164) 5'-gCA 3'	26(173) 5'-TTA 3'	26(173) 5'-TCT 3'	26(173) 5'-TCT 3'	27(175) 5'-ATA 3'	27(175) 5'-ATA 3'	102(400) 5'-TCT 3'	10(126) 5'-CAT 3'	23(163) 5'-AgA 3'	140(514) 5'-CAA 3'	140(516) 5'-ACC 3'	55(260) 5'-gCC 3'
	24(166) 5'-gCC 3'				34(196) 5'-ACT 3'	41(219) 5'-gAg 3'		26(173) 5'-TTT 3'	44(226) 5'-ACA 3'			
	39(211) 5'-CAT 3'				38(209) 5'-CgA 3'	47(235) 5'-TgA 3'		47(235) 5'-TgC 3'	55(259) 5'-g.C 3'			
	41(219) 5'-gAg 3'				44(227) 5'-CgC 3'	48(239) 5'-TCC 3'		48(239) 5'-TCA 3'	70(306) 5'-AgC 3'			
					55(259) 5'-g.C 3'	74(316) 5'-CgA 3'						
3'-primer(s) ³	86(353) 5'-gCT 3'	70(304) 5'-CCC 3'	48(238) 5'-CCA 3'	71(307) 5'-ggC 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	142(520) 5'-AAT 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	165(589) 5'-ggT 3'	185(650) 5'-CgA 3'	87(355) 5'-gAC 3'
		71(308) 5'-ggC 3'	56(263) 5'-CAA 3'	73(314) 5'-CCT 3'			146(532) 5'-ggC 3'			167(595) 5'-ACA 3'		88(358) 5'-gCC 3'
		71(308) 5'-ggA 3'	61(278) 5'-TCT 3'	76(324) 5'-gTT 3'			150(546) 5'-CCC 3'			184(647) 5'-Tgg 3'		93(374) 5'-gCA 3'
		88(359) 5'-TgA 3'	76(324) 5'-gTT 3'	87(356) 5'-ggg 3'			172(611) 5'-AgA 3'					
		89(362) 5'-TCA 3'	90(365) 5'-AgT 3'	89(361) 5'-CgC 3'								
Well No.	37	38	39	40	41	42	43	44	45	46	47	48

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

Well No.	49	50	51	52	53	54	55	56	57	58	59	60
Length of spec.	150	105	105	80	125	95	155	155	100	155	100	185
PCR product	175	150	140	185	155	170	195	200	215	180		
			165		190	220			240			
			220									
Length of int. pos. control ¹	430	515	430	430	430	430	430	430	515	430	430	430
5'-primer(s) ²	140(514) 5'-CAA 3'	140(514) 5'-CAA 3'	26(173) 5'-TCT 3'	32(192) 5'-..T 3'	37(205) 5'-AgA 3'	8(119) 5'-CgC 3'	13(134) 5'-ggC 3'	10(125) 5'-CCg 3'	104(406) 5'-CAC 3'	13(134) 5'-ggC 3'	26(173) 5'-TTA 3'	9(122) 5'-gTT 3'
			135(500) 5'-TgA 3'	37(205) 5'-Agg 3'	37(205) 5'-Agg 3'	13(133) 5'-AgA 3'		16(142) 5'-gCC 3'	112(430) 5'-ACT 3'			
				37(207) 5'-TAA 3'	48(239) 5'-CCC 3'	16(142) 5'-gCC 3'		29(182) 5'-CAA 3'	151(547) 5'-ACA 3'			
				44(226) 5'-ACA 3'	48(239) 5'-TCA 3'	125(470) 5'-Agg 3'						
				74(316) 5'-CgA 3'	58(268) 5'-ACT 3'	167(595) 5'-AgA 3'						
3'-primer(s) ³	170(604) 5'-gAA 3'	161(577) 5'-CAT 3'	48(238) 5'-.Cg 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	55(260) 5'-gCg 3'	52(250) 5'-CgC 3'	66(294) 5'-ACT 3'	170(604) 5'-gAC 3'	52(250) 5'-CgC 3'	45(230) 5'-CCC 3'	57(266) 5'-Cgg 3'
	174(617) 5'-CgC 3'	172(611) 5'-AgC 3'	48(238) 5'-CCA 3'			185(650) 5'-CgA 3'	64(287) 5'-TCC 3'			60(274) 5'-gTA 3'		
	177(625) 5'-gTT 3'	177(625) 5'-gTT 3'	55(260) 5'-gCg 3'									
	184(647) 5'-Tgg 3'	179(632) 5'-ggT 3'	64(286) 5'-CTA 3'									
			71(307) 5'-ggC 3'									
			86(354) 5'-Agg 3'									
			175(619) 5'-CAT 3'									
Well No.	49	50	51	52	53	54	55	56	57	58	59	60

Well No.	61	62	63
Length of spec.	115	175	160
PCR product	155		
	255		
Length of int. pos. control ¹	430	430	430
5'-primer(s) ²	28(179) 5'-gAC 3'	102(400) 5'-TCT 3'	1 st I 5'-... 3'
	99(391) 5'-CAA 3'		
	145(530) 5'-CAA 3'		
3'-primer(s) ³	66(294) 5'-ATg 3'	146(532) 5'-ggT 3'	1 st I 5'-CCT 3'
	170(604) 5'-gAC 3'	146(532) 5'-ggC 3'	
Well No.	61	62	63

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 **“Instructions for Use” (IFU)**

Lot No.: 9K1

Lot-specific information

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.214-24/06 – including Taq polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without Taq polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

CELL LINE VALIDATION SHEET																					
DQB1*03 SSP subtyping kit ²																					
					Production No.	Well															
						1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
						201784901	201784902	201784903	201905104	201784905	201784906	201784907	201784908	201784909	201784910	201784911	201784912	201905113	201784914	201895915	201895916
	IHWC cell line ¹	DQB1																			
1	9001 SA	*05:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707	*06:01	*02:02			-	-	+	+	-	-	-	-	-	-	+	-	-	-	-	-
3	9011 E4181324	*06:01				-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
4	9275 GU373	*02:01				-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011	*05:02				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM	*03:02	*06:01			-	-	+	+	-	-	-	-	-	-	+	-	-	-	+	-
7	9020 QBL	*02:01				-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU	*03:01				+	+	-	-	-	-	-	-	-	-	+	-	-	+	+	+
9	9026 YAR	*03:02				-	-	+	+	-	-	-	-	-	-	-	-	-	-	+	-
10	9107 LKT3	*04:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT	*02:02				-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB	*03:03				-	-	-	+	-	-	-	-	-	-	-	-	-	-	+	+
13	9004 JESTHOM	*05:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA	*04:02				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075 DKB	*03:03				-	-	-	+	-	-	-	-	-	-	-	-	-	-	+	+
16	9037 SWEIG007	*03:01				+	+	-	-	-	-	-	-	-	-	+	-	-	+	+	+
17	9282 CTM3953540	*02:01	*06:03			-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
18	9257 32367	*06:02	*02:02			-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16	*03:01				+	+	-	-	-	-	-	-	-	-	+	-	-	+	+	+
20	9059 SLE005	*06:04				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA	*03:01				+	+	-	-	-	-	-	-	-	-	-	+	-	+	+	+
22	9056 KOSE	*05:03	*06:04			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL	*05:03	*06:01			-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
24	9035 JBUSH	*03:01				+	+	-	-	-	-	-	-	-	-	+	-	-	+	+	+
25	9049 IBW9	*02:02				-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49	*02:01				-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007	*04:01	*05:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320 BEL5GB	*02:02	*03:01			+	+	+	+	-	-	-	-	-	-	+	-	-	+	+	+
29	9050 MOU	*02:02				-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH	*04:02				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019 DUCAF	*02:01				-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG	*03:01				+	+	-	-	-	-	-	-	-	-	+	-	-	+	+	+
33	9098 MT14B	*03:02				-	-	+	+	-	-	-	-	-	-	-	-	-	-	+	-
34	9104 DHIF	*03:01				+	+	-	-	-	-	-	-	-	-	+	-	-	+	+	+
35	9302 SSTO	*03:05				-	-	-	-	-	+	-	-	-	-	-	-	-	-	+	-
36	9024 KT17	*03:02				-	-	+	+	-	-	-	-	-	-	-	-	-	-	+	-
37	9065 HHKB	*06:03				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL	*03:01				+	+	-	-	-	-	-	-	-	-	+	-	-	+	+	+
39	9315 CML	*02:01	*03:01			+	+	+	+	-	-	-	-	-	-	+	-	-	+	+	+
40	9134 WHONP199	*02:02	*03:03			-	-	+	+	-	-	-	-	-	-	-	-	-	-	+	+
41	9055 H0301	*06:09				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089	*06:01				-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
43	9076 T7526	*03:03				-	-	-	+	-	-	-	-	-	-	-	-	-	-	+	+
44	9057 TEM	*05:03				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO	*02:02				-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU	*06:02				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO	*03:01				+	+	-	-	-	-	-	-	-	-	+	-	-	+	+	+
48	9303 TER-ND	*05:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-



101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

CELL LINE VALIDATION SHEET																			
DQB1*03 SSP subtyping kit ²																			
				Well															
				17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
				201784917	201784918	201905119	201784920	201784921	201784922	201905123	201784924	201784925	201784926	201784927	202016828	201784929	201905130	201784931	201784932
	IHWC cell line ¹		DQB1	Production No.															
1	9001 SA		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707		*06:01	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011 E4181324		*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275 GU373		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011		*05:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM		*03:02	*06:01	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
7	9020 QBL		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU		*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
9	9026 YAR		*03:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3		*04:01		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB		*03:03		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA		*04:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075 DKB		*03:03		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007		*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
17	9282 CTM3953540		*02:01	*06:03	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-
18	9257 32367		*06:02	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16		*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005		*06:04		-	-	-	-	W	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA		*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE		*05:03	*06:04	-	-	-	-	W	-	-	-	-	-	-	-	-	-	-
23	9124 IHL		*05:03	*06:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035 JBUSH		*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007		*04:01	*05:01	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320 BEL5GB		*02:02	*03:01	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH		*04:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019 DUCAF		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG		*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B		*03:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF		*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO		*03:05		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17		*03:02		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB		*06:03		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-
38	9099 LZL		*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML		*02:01	*03:01	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
40	9134 WHONP199		*02:02	*03:03	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055 H0301		*06:09		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089		*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526		*03:03		-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM		*05:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU		*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO		*03:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
48	9303 TER-ND		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

101.214-24/06 – including Taq polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without Taq polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

CELL LINE VALIDATION SHEET																			
DQB1*03 SSP subtyping kit ²																			
				Well															
				33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48
				201784933	201784934	201784935	201784936	201784937	201784938	201784939	201905140	201895941	201784942	201905143	201784944	201905145	202016846	201784947	201905148
	IHWC cell line ¹	DQB1	Production No.																
1	9001 SA	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707	*06:01	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011 E4181324	*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275 GU373	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011	*05:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM	*03:02	*06:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
7	9020 QBL	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026 YAR	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075 DKB	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282 CTM3953540	*02:01	*06:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257 32367	*06:02	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005	*06:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE	*05:03	*06:04	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL	*05:03	*06:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035 JBUSH	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007	*04:01	*05:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320 BEL5GB	*02:02	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019 DUCAF	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO	*03:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB	*06:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML	*02:01	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134 WHONP199	*02:02	*03:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055 H0301	*06:09		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089	*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM	*05:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303 TER-ND	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-



101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: 9K1

Lot-specific information

CELL LINE VALIDATION SHEET																			
DQB1*03 SSP subtyping kit ²																			
				Production No.	Well														
					49	50	51	52	53	54	55	56	57	58	59	60	61	62	63
					201784949	201905150	201905151	201784952	201784953	201784954	201784955	201905156	201905157	201784958	201784959	201784960	201895961	201784962	201905163
IHWC cell line ¹			DQB1																
1	9001 SA		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707		*06:01	*02:02	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
3	9011 E4181324		*06:01		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
4	9275 GU373		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011		*05:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM		*03:02	*06:01	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
7	9020 QBL		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026 YAR		*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB		*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA		*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075 DKB		*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282 CTM3953540		*02:01	*06:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257 32367		*06:02	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005		*06:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE		*05:03	*06:04	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL		*05:03	*06:01	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
24	9035 JBUSH		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007		*04:01	*05:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320 BEL5GB		*02:02	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH		*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019 DUCAF		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B		*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO		*03:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17		*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB		*06:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML		*02:01	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134 WHONP199		*02:02	*03:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055 H0301		*06:09		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089		*06:01		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
43	9076 T7526		*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM		*05:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU		*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303 TER-ND		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution of the kit has been tested against 48 well characterized IHWC cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 7 to 10, 13, 17, 20, 23 to 46, 48 to 58 and 60 to 63 were available. The specificities of the primers in primer solutions 7, 9, 10, 23, 25 to 28, 35, 38, 39, 52 to 54, 60 and 61 were tested by separately adding additional 5'-primers respectively 3'-primers. In primer solutions 8, 17, 20, 29, 30, 33, 36, 37, 41, 42, 44, 45, 56, 57 and 63 it was only possible to test the 3'-primers, the 5'-primers were not possible to test. In primer solutions 13, 24, 31, 32, 34, 40, 43, 46, 48 to 51, 55, 58 and 62 it was only possible to test the 5'-primers, the 3'-primers were not possible to test. In primer solutions 1 to 4, 7, 9 to 12, 14 to 16, 19, 23, 25, 35, 38 and 39 one or more of the 3'-primers were not possible to test, and in primer solution 6, 16, 25 to 28, 35, 52 to 54 and 61 one or more of the 5'-primers were not possible to test.

101.214-24/06 – including *Taq* polymerase, IFU-01 Visit <https://labproducts.caredx.com> for
101.214-24u/06u – without *Taq* polymerase, IFU-02 “Instructions for Use” (IFU)

Lot No.: **9K1**

Lot-specific information

ADDRESSES:

Manufacturer:

CareDx AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-508 939 00

Fax: +46-8-717 88 18

E-mail: orders-se@caredx.com

Web page: <https://labproducts.caredx.com/>

Distributed by:

CareDx GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

Tel: +43-1-710 15 00

Fax: +43-1-710 15 00 10

E-mail: orders-at@caredx.com

Web page: <https://labproducts.caredx.com/>

CareDx Lab Solutions Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-653-78171

Fax: 610-344-7989

E-mail: orders-us@caredx.com

Web page: <https://labproducts.caredx.com/>

For information on CareDx distributors worldwide, contact **CareDx GmbH**.