

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Olerup SSP[®] DQB1 high resolution for frequent alleles

Product number:	101.221-12 – including <i>Taq</i> polymerase 101.221-12u – without <i>Taq</i> polymerase
Lot number:	2F4
Expiry date:	2021-07-01
Number of tests:	12 tests
Number of wells per test:	95+ 1
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
- Adhesive PCR seals	RT
- Product Insert	RT

This Product Description is only valid for Lot No. 2F4.

Complete product documentation consists of generic Instructions for Use (IFU), lot specific Product Insert, Worksheet and Certificate.

**CHANGES COMPARED TO THE PREVIOUS OLERUP SSP[®]
DQB1 HIGH RESOLUTION LOT (9D8)**

The DQB1*high kit is updated for new alleles to enable separation of:

- Confirmed¹ alleles as listed in the IMGT/HLA database.
- Polymorphisms in exons outside of the region encoding the peptide binding domain.
- Null and Alternatively expressed alleles.

The format of the Worksheet has been changed.

¹As described in section Uniquely Identified Alleles.

The DQB1 high resolution specificity and interpretation tables have been updated with the DQB1 alleles described since the previous *Olerup SSP[®]* DQB1 high resolution lot was made (**Lot No. 9D8**). The kit design is based on IMGT/HLA database 3.27.0.

Product Insert

101.221-12 – including *Taq* pol., IFU-01

Visit <https://labproducts.caredx.com> for

101.221-12u – without *Taq* pol., IFU-02

“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
5	Added	-	5'-primer added for the DQB1*05:106 allele.
10	Modified	Modified	5'-and 3'-primers modified for the increased yield and improved HLA-specific amplification.
12	Added	Added	5'-primer added for the DQB1*06:210 allele, 3'-primer added for the DQB1*03:239 allele.
13	Added	-	5'-primer added for the DQB1*02:62 allele.
15	-	Added	3'-primer added for the DQB1*02:78 allele.
17	-	Added	3'-primers added for the DQB1*03:232 and DQB1*04:32 alleles.
18	Added	-	5'-primer added for the DQB1*03:237N allele.
25	-	Removed, added	3'-primer removed for decreased tendency of primer oligomer formation, 3'-primers added for the DQB1*05:03:14 and DQB1*05:128N allele.
30	-	Added	3'-primer added for the DQB1*02:67N allele.
33	Added	-	5'-primer added for the DQB1*03:199 allele.
34	Added	-	5'-primer added for the DQB1*05:112 allele.
35	Added	Added	5'-primer added for the DQB1*06:212 allele, 3'-primer added for the DQB1*05:103 allele.
38	Added	Added	5'-primer added for the DQB1*06:194 allele, 3'-primer added for the DQB1*04:33 and DQB1*06:188 alleles.
39	Added	Added, modified	5'-primers added for the DQB1*05:111 and DQB1*06:202 alleles, 3'-primer added for the DQB1*03:191 allele, 3'-primer modified for improved HLA-specific amplification.
40	Added	-	5'-primers added for the DQB1*03:245 and DQB1*05:133 alleles.
41	Removed	Added	5'-primer removed for improved HLA-specific amplification, 3'-primer added for the DQB1*03:212 allele.
42	-	Added	3'-primer added for the DQB1*03:239 allele.
44	-	Modified, added	3'-primer modified for decreased tendency of primer oligomer formation, 3'-primer added for the DQB1*03:191 allele.
45	Added	Added	5'-primer added for the DQB1*06:212 allele, 3'-primers added for the DQB1*06:216N and DQB1*06:224 alleles.
49	Added	-	5'-primer added for the DQB1*03:246 allele.
50	Added	-	5'-primer added for the DQB1*03:201 allele.
51	Modified, moved	Moved	5'-primer modified, 5'-primer moved to well 91 and 3'-primer moved to well 90 for decreased tendency of primer oligomer formation.
54	-	Added	3'-primer added for the DQB1*04:36N allele.
56	-	Added	3'-primer added for the DQB1*06:224 allele.
58	-	Added	3'-primer added for the DQB1*03:228 allele.
61	-	Added	3'-primer added for the DQB1*05:128N allele.
62	Added	-	5'-primer added for the DQB1*03:205 allele.
63	-	Added	3'-primer added for the DQB1*03:240 allele.
66	Added	-	5'-primers added for the DQB1*03:243 and DQB1*05:133 alleles.
67	Added	-	5'-primer added for the DQB1*02:26 allele.
69	Added, modified	-	5'-primer added for the DQB1*03:225 allele, 5'-primers modified for improved HLA-specific amplification.
71	Added, modified	Exchanged	5'-primer added for the DQB1*03:207 allele, 5'-primer modified and 3'-primer exchanged for improved HLA-specific

Product Insert

101.221-12 – including *Taq* pol., IFU-01

Visit <https://labproducts.caredx.com> for

101.221-12u – without *Taq* pol., IFU-02

“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

			amplification and decreased tendency of primer oligomer formation.
72	Added	Added	5'-primers added for the DQB1*03:150 and DQB1*03:201 alleles, 3'-primer added for the DQB1*06:218 allele.
73	Added	-	5'-primer added for the DQB1*03:150 allele.
74	-	Added	3'-primer added for the DQB1*05:97 allele.
75	Added, modified	Added	5'-primer added for the DQB1*06:202 allele, 5'-primer modified for improved HLA-specific amplification, 3'-primer added for the DQB1*06:122 allele.
76	Modified	Modified, added	5'-and 3'-primers modified for improved HLA-specific amplification, 3'-primer added for the DQB1*06:53:02 allele.
77	-	Added, modified	3'-primer added for the DQB1*06:190:01 allele, 3'-primer modified for improved HLA-specific amplification.
78	Added	Added	5'-primer added for the DQB1*02:232 allele, 3'-primer added for the DQB1*06:190:01 allele.
79	Added	-	5'-primer added for the DQB1*06:187 allele.
84	Added	Added	Primer pair added for the DQB1*06:217 allele.
85	Added	-	5'-primer added for the DQB1*06:184 allele.
86	Added, modified	-	5'-primer added for the DQB1*06:210 allele, 5'-primer modified for improved HLA-specific amplification.
87	Added	-	5'-primer added for the DQB1*06:184 allele.
90	-	Added	3'-primer added for the DQB1*03:247 allele, 3'-primer added from well 51 for decreased tendency of primer oligomer formation.
91	Added	-	5'-primer added from well 51 for decreased tendency of primer oligomer formation.
92	-	Added	3'-primer added for the DQB1*02:78 allele.
94	Modified	-	5'-primer modified for increased yield.

Changes in revision R01 compared to R00:

1. Primer mixes 17 and 51 do not amplify the DQB1*02:10, 02:81 and 03:49 alleles. Thus, this lot of the DQB1*high kit cannot distinguish the DQB1*02:10 and the DQB1*02:02:01:01-02:02:04 and 02:65 alleles, the DQB1*02:81 and the DQB1*02:01:01-02:01:24, 02:34-02:35, 02:38-02:39, 02:43-02:45, 02:49, 02:52, 02:55-02:56, 02:59-02:61, 02:63, 02:66, 02:68-02:71, 02:73-02:76, 02:82-02:83, and 02:85-02:88 alleles or the DQB1*03:49 and the DQB1*03:01:01-03:01:02, 03:01:05-03:01:07, 03:01:09-03:01:12, 03:01:14, 03:01:16-03:01:24, 03:01:26-03:01:33, 03:01:35, 03:92-03:93, 03:169, 03:182, 03:198, 03:206, 03:241-03:242 and 03:253-03:254 alleles. This has been corrected in the Interpretation and Specificity Tables.

Changes in revision R02 compared to R01:

1. The expiration date has been altered due to extension of shelf-life.

Changes in revision R04 compared to R02:

1. In primer mix 25 the DQB1*04:28 allele does not get amplified at 125 bp. This has been corrected in the specificity table.

Changes in revision R06 compared to R04:

1. Primer mix 18 does not amplify the DQB1*03:24 and 03:79 alleles.
2. Primer mix 51 does not amplify the DQB1*02:06, 02:48, 03:24 and 03:79 alleles.

The corrections above have been implemented in the Specificity and Interpretation tables.

Changes in revision R07 compared to R06:

1. Primer mix 49 does not amplify the DQB1*05:84 and DQB1*06:103 alleles. This correction has been implemented in the Specificity and Interpretation tables.

Product Insert

101.221-12 – including *Taq* pol., IFU-01

Visit <https://labproducts.caredx.com> for

101.221-12u – without *Taq* pol., IFU-02

“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Well **96** contains Negative Control primer pairs, that will amplify more than 95% of the *Olerup* SSP[®] HLA Class I, DRB, DQB1, DPB1 and DQA1 amplicons as well as all the amplicons generated by the control primer pairs matching the human growth hormone gene.

HLA-specific PCR product sizes range from 75 to 200 base pairs.

The PCR product generated by the positive control primer pair is 430 base pairs.

Length of PCR product	105	200	105	80	75	80	85
5'-primer¹	164	340	440	45	45	43	36
	5'-CAC ^{3'}	5'-Agg ^{3'}	5'-TTA ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-TAC ^{3'}
							36
							5'-TAT ^{3'}
3'-primer²	231	2nd I	507	59	58	57	47
	5'-TgC ^{3'}	5'-AAA ^{3'}	5'-TTg ^{3'}	5'-CTC ^{3'}	5'-ggC ^{3'}	5'-CTC ^{3'}	5'-ACA ^{3'}
							48
							5'-gCA ^{3'}
							48
							5'-gCC ^{3'}
							52
							5'-TgT ^{3'}
A*	+	+	+				
B*	+	+	+				
C*	+	+	+				
DRB1				+	+		
DRB3				+	+		
DRB5				+			
DQB1					+		
DPB1						+	
DQA1							+

¹The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

²The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

PRODUCT DESCRIPTION

DQB1 SSP typing

CONTENT

The primer set contains 5'- and 3'-primers for identifying the DQB1*02:01:01 to DQB1*06:230 alleles.

Please note that DQB1 amplifications usually are somewhat less pronounced than e.g. DRB and DQA1 amplifications even when using the same DNA preparation and exactly the same experimental procedures.

PLATE LAYOUT

Each test consists of 96 PCR reactions in a 96 well cut PCR plate.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48
49	50	51	52	53	54	55	56
57	58	59	60	61	62	63	64
65	66	67	68	69	70	71	72
73	74	75	76	77	78	79	80
81	82	83	84	85	86	87	88
89	90	91	92	93	94	95	NC

The 96 well cut PCR plate is marked with ‘DQB1 “HIGH”’ in silver/gray ink.

Well No. 1 is marked with the Lot No. ‘2F4’.

Wells 1 to 95 – DQB1 high resolution primers.

Well 96 – Negative Control.

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded. The PCR plates are covered with a PCR-compatible foil.

INTERPRETATION

Only DQB1 alleles will be amplified by the DQB1 high resolution for frequent alleles typing kit. Thus, the interpretation of DQB1 high resolution for frequent alleles typings is not influenced by the DQB2 or DQB3 genes or other HLA class II genes.

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

UNIQUELY IDENTIFIED ALLELES

All the phenotypically different DQB1 alleles, i.e. **DQB1*02:01:01 to 02:88, DQB1*03:01:01:01 to 03:256, DQB1*04:01:01 to 04:39, DQB1*05:01:01:01 to 05:143 and DQB1*06:01:01 to 06:230**, recognized by the HLA Nomenclature Committee in January 2017^{1,2} will give rise to unique amplification patterns by the primers in the DQB1 high resolution for frequent alleles typing kit^{3,4}.

The following DQB1 alleles can be distinguished by the different sizes of the specific PCR product:

Alleles	Primer mix	Alleles	Primer mix
DQB1*02:15, 02:30	67	DQB1*05:13, 05:24	59
DQB1*02:16, 02:20N	64	DQB1*05:17, 05:47	54
DQB1*02:23, 02:28	68	DQB1*05:25, 05:112	34
DQB1*02:29, 02:31	69	DQB1*06:02:04, 06:02:16	47
DQB1*03:31, 03:87	35	DQB1*06:02:19, 06:48	62
DQB1*03:35, 03:246	49	DQB1*06:31, 06:145	35
DQB1*03:59, 03:71	43	DQB1*06:35, 06:45	35
DQB1*03:66N, 03:224	42	DQB1*06:70, 06:76	80
DQB1*03:116, 03:243	66	DQB1*06:75N, 06:80	84
DQB1*04:05, 04:15	30	DQB1*06:102N, 06:104	90
DQB1*05:10, 05:23	54		

¹DQB1 alleles listed on the IMGT/HLA web page 2017-January-20, release 3.27.0, www.ebi.ac.uk/imgt/hla.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

³The DQB1 high primer set cannot separate the DQB1*05:03:14 from the DQB1*05:43:01 alleles. These alleles can be distinguished by the DQB1*05 kit.

⁴This lot of the DQB1*high kit cannot distinguish the DQB1*02:10 and the DQB1*02:02:01:01-02:02:04 and 02:65 alleles, the DQB1*02:81 and the DQB1*02:01:01-02:01:24, 02:34-02:35, 02:38-02:39, 02:43-02:45, 02:49, 02:52, 02:55-02:56, 02:59-02:61, 02:63, 02:66, 02:68-02:71, 02:73-02:76, 02:82-02:83 and 02:85-02:88 alleles or the DQB1*03:49 and the DQB1*03:01:01:01-03:01:02, 03:01:05-03:01:07, 03:01:09-03:01:12, 03:01:14, 03:01:16-03:01:24, 03:01:26-03:01:33, 03:01:35, 03:92-03:93, 03:169, 03:182, 03:198, 03:206, 03:241-03:242 and 03:253-03:254 alleles.

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

ALLELE CONFIRMATION STATUS

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DQB1*02:01:01	Confirmed	DQB1*02:04	Unconfirmed	DQB1*02:32	Unconfirmed	DQB1*02:62	Confirmed
DQB1*02:01:02	Unconfirmed	DQB1*02:05	Confirmed	DQB1*02:33	Confirmed	DQB1*02:63	Unconfirmed
DQB1*02:01:03	Unconfirmed	DQB1*02:06	Unconfirmed	DQB1*02:34	Unconfirmed	DQB1*02:64	Unconfirmed
DQB1*02:01:04	Confirmed	DQB1*02:07:01	Confirmed	DQB1*02:35	Unconfirmed	DQB1*02:65	Unconfirmed
DQB1*02:01:05	Confirmed	DQB1*02:07:02	Unconfirmed	DQB1*02:36	Unconfirmed	DQB1*02:66	Unconfirmed
DQB1*02:01:06	Unconfirmed	DQB1*02:08	Unconfirmed	DQB1*02:37	Unconfirmed	DQB1*02:67N	Unconfirmed
DQB1*02:01:07	Confirmed	DQB1*02:09	Unconfirmed	DQB1*02:38	Unconfirmed	DQB1*02:68	Unconfirmed
DQB1*02:01:08	Confirmed	DQB1*02:10	Unconfirmed	DQB1*02:39	Unconfirmed	DQB1*02:69	Unconfirmed
DQB1*02:01:09	Unconfirmed	DQB1*02:11	Unconfirmed	DQB1*02:40	Unconfirmed	DQB1*02:70	Unconfirmed
DQB1*02:01:10	Confirmed	DQB1*02:12	Confirmed	DQB1*02:41	Confirmed	DQB1*02:71	Unconfirmed
DQB1*02:01:11	Confirmed	DQB1*02:13	Unconfirmed	DQB1*02:42	Confirmed	DQB1*02:72	Unconfirmed
DQB1*02:01:12	Unconfirmed	DQB1*02:14:01	Confirmed	DQB1*02:43	Unconfirmed	DQB1*02:73	Unconfirmed
DQB1*02:01:13	Unconfirmed	DQB1*02:14:02	Unconfirmed	DQB1*02:44	Unconfirmed	DQB1*02:74	Unconfirmed
DQB1*02:01:14	Confirmed	DQB1*02:15	Unconfirmed	DQB1*02:45	Unconfirmed	DQB1*02:75	Unconfirmed
DQB1*02:01:15	Unconfirmed	DQB1*02:16	Confirmed	DQB1*02:46	Confirmed	DQB1*02:76	Unconfirmed
DQB1*02:01:16	Unconfirmed	DQB1*02:17	Confirmed	DQB1*02:47	Unconfirmed	DQB1*02:77	Unconfirmed
DQB1*02:01:17	Unconfirmed	DQB1*02:18N	Confirmed	DQB1*02:48	Unconfirmed	DQB1*02:78	Confirmed
DQB1*02:01:18	Unconfirmed	DQB1*02:19	Confirmed	DQB1*02:49	Unconfirmed	DQB1*02:79	Unconfirmed
DQB1*02:01:19	Confirmed	DQB1*02:20N	Confirmed	DQB1*02:50	Confirmed	DQB1*02:80	Unconfirmed
DQB1*02:01:20	Confirmed	DQB1*02:21	Unconfirmed	DQB1*02:51	Confirmed	DQB1*02:81	Unconfirmed
DQB1*02:01:21	Unconfirmed	DQB1*02:22	Unconfirmed	DQB1*02:52	Unconfirmed	DQB1*02:82	Unconfirmed
DQB1*02:01:22	Unconfirmed	DQB1*02:23	Confirmed	DQB1*02:53Q	Confirmed	DQB1*02:83	Unconfirmed
DQB1*02:01:23	Confirmed	DQB1*02:24	Unconfirmed	DQB1*02:54	Confirmed	DQB1*02:84	Unconfirmed
DQB1*02:01:24	Unconfirmed	DQB1*02:25	Unconfirmed	DQB1*02:55	Unconfirmed	DQB1*02:85	Unconfirmed
DQB1*02:02:01:01	Confirmed	DQB1*02:26	Confirmed	DQB1*02:56	Unconfirmed	DQB1*02:86	Unconfirmed
DQB1*02:02:01:02	Unconfirmed	DQB1*02:27	Confirmed	DQB1*02:57	Unconfirmed	DQB1*02:87	Unconfirmed
DQB1*02:02:02	Unconfirmed	DQB1*02:28	Confirmed	DQB1*02:58N	Unconfirmed	DQB1*02:88	Unconfirmed
DQB1*02:02:03	Unconfirmed	DQB1*02:29	Confirmed	DQB1*02:59	Unconfirmed		
DQB1*02:02:04	Unconfirmed	DQB1*02:30	Confirmed	DQB1*02:60	Unconfirmed		
DQB1*02:03	Confirmed	DQB1*02:31	Unconfirmed	DQB1*02:61	Unconfirmed		



Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DQB1*03:01:01:01	Confirmed	DQB1*03:02:09	Confirmed	DQB1*03:16	Confirmed	DQB1*03:62	Confirmed
DQB1*03:01:01:02	Unconfirmed	DQB1*03:02:10	Confirmed	DQB1*03:17:01	Confirmed	DQB1*03:63	Confirmed
DQB1*03:01:01:03	Confirmed	DQB1*03:02:11	Confirmed	DQB1*03:17:02	Confirmed	DQB1*03:64	Confirmed
DQB1*03:01:01:04	Unconfirmed	DQB1*03:02:12	Confirmed	DQB1*03:18	Unconfirmed	DQB1*03:65	Confirmed
DQB1*03:01:01:05	Confirmed	DQB1*03:02:13	Confirmed	DQB1*03:19:01	Confirmed	DQB1*03:66N	Unconfirmed
DQB1*03:01:02	Unconfirmed	DQB1*03:02:14	Unconfirmed	DQB1*03:19:02	Unconfirmed	DQB1*03:67	Confirmed
DQB1*03:01:03	Confirmed	DQB1*03:02:15	Confirmed	DQB1*03:20	Unconfirmed	DQB1*03:68	Confirmed
DQB1*03:01:04	Confirmed	DQB1*03:02:16	Confirmed	DQB1*03:21	Unconfirmed	DQB1*03:69	Unconfirmed
DQB1*03:01:05	Unconfirmed	DQB1*03:02:17	Unconfirmed	DQB1*03:22	Unconfirmed	DQB1*03:70	Confirmed
DQB1*03:01:06	Unconfirmed	DQB1*03:02:18	Unconfirmed	DQB1*03:23:01	Unconfirmed	DQB1*03:71	Unconfirmed
DQB1*03:01:07	Unconfirmed	DQB1*03:02:19	Confirmed	DQB1*03:23:02	Unconfirmed	DQB1*03:72	Confirmed
DQB1*03:01:08	Confirmed	DQB1*03:02:20	Unconfirmed	DQB1*03:24	Unconfirmed	DQB1*03:73	Unconfirmed
DQB1*03:01:09	Unconfirmed	DQB1*03:02:21	Confirmed	DQB1*03:25:01	Confirmed	DQB1*03:74	Unconfirmed
DQB1*03:01:10	Unconfirmed	DQB1*03:02:22	Unconfirmed	DQB1*03:25:02	Unconfirmed	DQB1*03:75	Confirmed
DQB1*03:01:11	Unconfirmed	DQB1*03:02:23	Unconfirmed	DQB1*03:26	Confirmed	DQB1*03:76	Unconfirmed
DQB1*03:01:12	Unconfirmed	DQB1*03:02:24	Unconfirmed	DQB1*03:27	Confirmed	DQB1*03:77	Unconfirmed
DQB1*03:01:13	Unconfirmed	DQB1*03:03:02:01	Confirmed	DQB1*03:28	Unconfirmed	DQB1*03:78	Unconfirmed
DQB1*03:01:14	Confirmed	DQB1*03:03:02:02	Unconfirmed	DQB1*03:29	Unconfirmed	DQB1*03:79	Unconfirmed
DQB1*03:01:15	Unconfirmed	DQB1*03:03:02:03	Unconfirmed	DQB1*03:30	Unconfirmed	DQB1*03:80	Unconfirmed
DQB1*03:01:16	Unconfirmed	DQB1*03:03:02:04	Unconfirmed	DQB1*03:31	Unconfirmed	DQB1*03:81	Confirmed
DQB1*03:01:17	Confirmed	DQB1*03:03:03	Unconfirmed	DQB1*03:32	Unconfirmed	DQB1*03:82	Unconfirmed
DQB1*03:01:18	Unconfirmed	DQB1*03:03:04	Confirmed	DQB1*03:33	Unconfirmed	DQB1*03:83	Unconfirmed
DQB1*03:01:19	Unconfirmed	DQB1*03:03:05	Confirmed	DQB1*03:34	Unconfirmed	DQB1*03:84N	Unconfirmed
DQB1*03:01:20	Unconfirmed	DQB1*03:03:06	Unconfirmed	DQB1*03:35	Unconfirmed	DQB1*03:85	Unconfirmed
DQB1*03:01:21	Confirmed	DQB1*03:03:07	Confirmed	DQB1*03:36	Unconfirmed	DQB1*03:86	Unconfirmed
DQB1*03:01:22	Confirmed	DQB1*03:03:08	Confirmed	DQB1*03:37	Confirmed	DQB1*03:87	Unconfirmed
DQB1*03:01:23	Unconfirmed	DQB1*03:03:09	Unconfirmed	DQB1*03:38	Unconfirmed	DQB1*03:88	Unconfirmed
DQB1*03:01:24	Unconfirmed	DQB1*03:03:10	Unconfirmed	DQB1*03:39	Unconfirmed	DQB1*03:89	Unconfirmed
DQB1*03:01:25	Unconfirmed	DQB1*03:03:11	Confirmed	DQB1*03:40	Confirmed	DQB1*03:90N	Unconfirmed
DQB1*03:01:26	Unconfirmed	DQB1*03:03:12	Unconfirmed	DQB1*03:41	Unconfirmed	DQB1*03:91Q	Unconfirmed
DQB1*03:01:27	Confirmed	DQB1*03:03:13	Unconfirmed	DQB1*03:42	Unconfirmed	DQB1*03:92	Unconfirmed
DQB1*03:01:28	Unconfirmed	DQB1*03:04:01	Confirmed	DQB1*03:43	Unconfirmed	DQB1*03:93	Unconfirmed
DQB1*03:01:29	Unconfirmed	DQB1*03:04:02	Unconfirmed	DQB1*03:44	Confirmed	DQB1*03:94	Unconfirmed
DQB1*03:01:30	Unconfirmed	DQB1*03:04:03	Unconfirmed	DQB1*03:45	Confirmed	DQB1*03:95N	Unconfirmed
DQB1*03:01:31	Unconfirmed	DQB1*03:05:01	Confirmed	DQB1*03:46	Unconfirmed	DQB1*03:96	Unconfirmed
DQB1*03:01:32	Unconfirmed	DQB1*03:05:02	Unconfirmed	DQB1*03:47	Confirmed	DQB1*03:97	Unconfirmed
DQB1*03:01:33	Unconfirmed	DQB1*03:05:03	Confirmed	DQB1*03:48	Unconfirmed	DQB1*03:98	Unconfirmed
DQB1*03:01:34	Unconfirmed	DQB1*03:05:04	Unconfirmed	DQB1*03:49	Unconfirmed	DQB1*03:99Q	Unconfirmed
DQB1*03:01:35	Unconfirmed	DQB1*03:06	Unconfirmed	DQB1*03:50	Unconfirmed	DQB1*03:100	Unconfirmed
DQB1*03:01:36	Unconfirmed	DQB1*03:07	Unconfirmed	DQB1*03:51	Unconfirmed	DQB1*03:101	Unconfirmed
DQB1*03:01:37	Unconfirmed	DQB1*03:08	Unconfirmed	DQB1*03:52	Confirmed	DQB1*03:102	Confirmed
DQB1*03:02:01:01	Confirmed	DQB1*03:09	Confirmed	DQB1*03:53	Confirmed	DQB1*03:103	Unconfirmed
DQB1*03:02:01:02	Confirmed	DQB1*03:10:01	Unconfirmed	DQB1*03:54	Confirmed	DQB1*03:104	Unconfirmed
DQB1*03:02:02	Unconfirmed	DQB1*03:10:02	Confirmed	DQB1*03:55	Confirmed	DQB1*03:105	Confirmed
DQB1*03:02:03	Unconfirmed	DQB1*03:11	Confirmed	DQB1*03:56	Confirmed	DQB1*03:106	Confirmed
DQB1*03:02:04	Unconfirmed	DQB1*03:12	Unconfirmed	DQB1*03:57	Unconfirmed	DQB1*03:107	Unconfirmed
DQB1*03:02:05	Unconfirmed	DQB1*03:13	Confirmed	DQB1*03:58	Confirmed	DQB1*03:108	Unconfirmed
DQB1*03:02:06	Confirmed	DQB1*03:14:01	Unconfirmed	DQB1*03:59	Confirmed	DQB1*03:109	Confirmed
DQB1*03:02:07	Confirmed	DQB1*03:14:02	Confirmed	DQB1*03:60	Unconfirmed	DQB1*03:110	Unconfirmed
DQB1*03:02:08	Unconfirmed	DQB1*03:15	Unconfirmed	DQB1*03:61	Unconfirmed	DQB1*03:111	Unconfirmed



Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DQB1*03:112	Unconfirmed	DQB1*03:152	Unconfirmed	DQB1*03:192	Unconfirmed	DQB1*03:232	Confirmed
DQB1*03:113	Confirmed	DQB1*03:153	Unconfirmed	DQB1*03:193	Unconfirmed	DQB1*03:233	Confirmed
DQB1*03:114	Confirmed	DQB1*03:154	Unconfirmed	DQB1*03:194	Unconfirmed	DQB1*03:234	Unconfirmed
DQB1*03:115	Unconfirmed	DQB1*03:155	Unconfirmed	DQB1*03:195	Confirmed	DQB1*03:235	Unconfirmed
DQB1*03:116	Unconfirmed	DQB1*03:156	Unconfirmed	DQB1*03:196	Confirmed	DQB1*03:236	Unconfirmed
DQB1*03:117	Unconfirmed	DQB1*03:157	Unconfirmed	DQB1*03:197Q	Confirmed	DQB1*03:237N	Unconfirmed
DQB1*03:118N	Confirmed	DQB1*03:158	Confirmed	DQB1*03:198	Unconfirmed	DQB1*03:238	Unconfirmed
DQB1*03:119	Confirmed	DQB1*03:159	Unconfirmed	DQB1*03:199	Confirmed	DQB1*03:239	Confirmed
DQB1*03:120	Unconfirmed	DQB1*03:160	Unconfirmed	DQB1*03:200	Confirmed	DQB1*03:240	Confirmed
DQB1*03:121	Unconfirmed	DQB1*03:161	Unconfirmed	DQB1*03:201	Confirmed	DQB1*03:241	Unconfirmed
DQB1*03:122	Unconfirmed	DQB1*03:162	Confirmed	DQB1*03:202	Unconfirmed	DQB1*03:242	Unconfirmed
DQB1*03:123	Unconfirmed	DQB1*03:163	Unconfirmed	DQB1*03:203	Unconfirmed	DQB1*03:243	Confirmed
DQB1*03:124	Unconfirmed	DQB1*03:164	Unconfirmed	DQB1*03:204	Confirmed	DQB1*03:244	Unconfirmed
DQB1*03:125	Unconfirmed	DQB1*03:165	Unconfirmed	DQB1*03:205	Confirmed	DQB1*03:245	Confirmed
DQB1*03:126	Unconfirmed	DQB1*03:166	Unconfirmed	DQB1*03:206	Unconfirmed	DQB1*03:246	Confirmed
DQB1*03:127	Unconfirmed	DQB1*03:167	Unconfirmed	DQB1*03:207	Confirmed	DQB1*03:247	Confirmed
DQB1*03:128	Unconfirmed	DQB1*03:168	Unconfirmed	DQB1*03:208	Unconfirmed	DQB1*03:248	Unconfirmed
DQB1*03:129	Unconfirmed	DQB1*03:169	Unconfirmed	DQB1*03:209	Unconfirmed	DQB1*03:249	Unconfirmed
DQB1*03:130	Confirmed	DQB1*03:170	Confirmed	DQB1*03:210	Unconfirmed	DQB1*03:250	Unconfirmed
DQB1*03:131	Unconfirmed	DQB1*03:171	Unconfirmed	DQB1*03:211	Confirmed	DQB1*03:251	Unconfirmed
DQB1*03:132	Confirmed	DQB1*03:172	Unconfirmed	DQB1*03:212	Confirmed	DQB1*03:252	Unconfirmed
DQB1*03:133	Unconfirmed	DQB1*03:173	Unconfirmed	DQB1*03:213N	Unconfirmed	DQB1*03:253	Unconfirmed
DQB1*03:134	Unconfirmed	DQB1*03:174	Unconfirmed	DQB1*03:214	Unconfirmed	DQB1*03:254	Unconfirmed
DQB1*03:135	Unconfirmed	DQB1*03:175	Unconfirmed	DQB1*03:215	Unconfirmed	DQB1*03:255	Unconfirmed
DQB1*03:136	Unconfirmed	DQB1*03:176	Unconfirmed	DQB1*03:216	Unconfirmed	DQB1*03:256	Unconfirmed
DQB1*03:137	Unconfirmed	DQB1*03:177	Confirmed	DQB1*03:217	Unconfirmed		
DQB1*03:138	Confirmed	DQB1*03:178	Unconfirmed	DQB1*03:218	Unconfirmed		
DQB1*03:139	Confirmed	DQB1*03:179	Unconfirmed	DQB1*03:219	Unconfirmed		
DQB1*03:140	Unconfirmed	DQB1*03:180	Unconfirmed	DQB1*03:220	Unconfirmed		
DQB1*03:141	Unconfirmed	DQB1*03:181	Unconfirmed	DQB1*03:221	Confirmed		
DQB1*03:142	Unconfirmed	DQB1*03:182	Unconfirmed	DQB1*03:222	Unconfirmed		
DQB1*03:143	Unconfirmed	DQB1*03:183	Unconfirmed	DQB1*03:223	Unconfirmed		
DQB1*03:144	Unconfirmed	DQB1*03:184	Confirmed	DQB1*03:224	Unconfirmed		
DQB1*03:145	Unconfirmed	DQB1*03:185	Unconfirmed	DQB1*03:225	Confirmed		
DQB1*03:146	Unconfirmed	DQB1*03:186	Unconfirmed	DQB1*03:226	Unconfirmed		
DQB1*03:147	Unconfirmed	DQB1*03:187	Unconfirmed	DQB1*03:227	Unconfirmed		
DQB1*03:148	Unconfirmed	DQB1*03:188	Confirmed	DQB1*03:228	Confirmed		
DQB1*03:149	Unconfirmed	DQB1*03:189	Unconfirmed	DQB1*03:229	Unconfirmed		
DQB1*03:150	Confirmed	DQB1*03:190	Unconfirmed	DQB1*03:230	Unconfirmed		
DQB1*03:151	Unconfirmed	DQB1*03:191	Confirmed	DQB1*03:231	Unconfirmed		

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DQB1*04:01:01	Confirmed	DQB1*04:04	Unconfirmed	DQB1*04:24	Unconfirmed
DQB1*04:01:02	Unconfirmed	DQB1*04:05	Unconfirmed	DQB1*04:25N	Unconfirmed
DQB1*04:01:03	Unconfirmed	DQB1*04:06	Unconfirmed	DQB1*04:26	Unconfirmed
DQB1*04:01:04	Unconfirmed	DQB1*04:07	Unconfirmed	DQB1*04:27	Confirmed
DQB1*04:02:01:01	Confirmed	DQB1*04:08	Unconfirmed	DQB1*04:28	Unconfirmed
DQB1*04:02:01:02	Unconfirmed	DQB1*04:09	Unconfirmed	DQB1*04:29	Unconfirmed
DQB1*04:02:01:03	Unconfirmed	DQB1*04:10	Unconfirmed	DQB1*04:30	Unconfirmed
DQB1*04:02:02	Unconfirmed	DQB1*04:11	Confirmed	DQB1*04:31	Unconfirmed
DQB1*04:02:03	Confirmed	DQB1*04:12	Confirmed	DQB1*04:32	Confirmed
DQB1*04:02:04	Confirmed	DQB1*04:13	Unconfirmed	DQB1*04:33	Confirmed
DQB1*04:02:05	Unconfirmed	DQB1*04:14	Unconfirmed	DQB1*04:34	Unconfirmed
DQB1*04:02:06	Unconfirmed	DQB1*04:15	Unconfirmed	DQB1*04:35	Unconfirmed
DQB1*04:02:07	Confirmed	DQB1*04:16	Unconfirmed	DQB1*04:36N	Unconfirmed
DQB1*04:02:08	Unconfirmed	DQB1*04:17	Unconfirmed	DQB1*04:37	Unconfirmed
DQB1*04:02:09	Confirmed	DQB1*04:18	Confirmed	DQB1*04:38	Unconfirmed
DQB1*04:02:10	Unconfirmed	DQB1*04:19	Unconfirmed	DQB1*04:39	Unconfirmed
DQB1*04:02:11	Confirmed	DQB1*04:20	Unconfirmed		
DQB1*04:02:12	Unconfirmed	DQB1*04:21	Unconfirmed		
DQB1*04:03:01	Unconfirmed	DQB1*04:22	Unconfirmed		
DQB1*04:03:02	Unconfirmed	DQB1*04:23	Unconfirmed		

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DQB1*05:01:01:01	Confirmed	DQB1*05:03:12	Unconfirmed	DQB1*05:44	Confirmed	DQB1*05:92	Unconfirmed
DQB1*05:01:01:02	Unconfirmed	DQB1*05:03:13	Unconfirmed	DQB1*05:45	Unconfirmed	DQB1*05:93	Unconfirmed
DQB1*05:01:01:03	Confirmed	DQB1*05:03:14	Confirmed	DQB1*05:46	Unconfirmed	DQB1*05:94	Unconfirmed
DQB1*05:01:02	Unconfirmed	DQB1*05:03:15	Unconfirmed	DQB1*05:47	Unconfirmed	DQB1*05:95	Unconfirmed
DQB1*05:01:03	Unconfirmed	DQB1*05:03:16	Confirmed	DQB1*05:48	Unconfirmed	DQB1*05:96	Unconfirmed
DQB1*05:01:04	Confirmed	DQB1*05:03:17	Unconfirmed	DQB1*05:49	Confirmed	DQB1*05:97	Confirmed
DQB1*05:01:05	Confirmed	DQB1*05:04	Confirmed	DQB1*05:50	Unconfirmed	DQB1*05:98	Unconfirmed
DQB1*05:01:06	Unconfirmed	DQB1*05:05:01	Unconfirmed	DQB1*05:51	Confirmed	DQB1*05:99	Unconfirmed
DQB1*05:01:07	Unconfirmed	DQB1*05:05:02	Unconfirmed	DQB1*05:52	Confirmed	DQB1*05:100	Unconfirmed
DQB1*05:01:08	Confirmed	DQB1*05:06:01	Unconfirmed	DQB1*05:53	Unconfirmed	DQB1*05:101	Confirmed
DQB1*05:01:09	Unconfirmed	DQB1*05:06:02	Confirmed	DQB1*05:54	Unconfirmed	DQB1*05:102	Unconfirmed
DQB1*05:01:10	Unconfirmed	DQB1*05:07	Unconfirmed	DQB1*05:55	Unconfirmed	DQB1*05:103	Confirmed
DQB1*05:01:11	Unconfirmed	DQB1*05:08	Unconfirmed	DQB1*05:56	Unconfirmed	DQB1*05:104	Confirmed
DQB1*05:01:12	Unconfirmed	DQB1*05:09	Unconfirmed	DQB1*05:57	Confirmed	DQB1*05:105	Unconfirmed
DQB1*05:01:13	Unconfirmed	DQB1*05:10	Confirmed	DQB1*05:58	Unconfirmed	DQB1*05:106	Confirmed
DQB1*05:01:14	Unconfirmed	DQB1*05:11:01	Unconfirmed	DQB1*05:59	Unconfirmed	DQB1*05:107	Unconfirmed
DQB1*05:01:15	Unconfirmed	DQB1*05:11:02	Confirmed	DQB1*05:60	Unconfirmed	DQB1*05:108	Unconfirmed
DQB1*05:01:16	Unconfirmed	DQB1*05:12	Unconfirmed	DQB1*05:61	Unconfirmed	DQB1*05:109	Unconfirmed
DQB1*05:01:17	Unconfirmed	DQB1*05:13	Unconfirmed	DQB1*05:62	Unconfirmed	DQB1*05:110N	Unconfirmed
DQB1*05:01:18	Unconfirmed	DQB1*05:14	Unconfirmed	DQB1*05:63	Confirmed	DQB1*05:111	Confirmed
DQB1*05:01:19	Unconfirmed	DQB1*05:15	Unconfirmed	DQB1*05:64	Unconfirmed	DQB1*05:112	Confirmed
DQB1*05:01:20	Unconfirmed	DQB1*05:16	Unconfirmed	DQB1*05:65	Unconfirmed	DQB1*05:113	Unconfirmed
DQB1*05:01:21	Unconfirmed	DQB1*05:17	Unconfirmed	DQB1*05:66:01	Unconfirmed	DQB1*05:114	Unconfirmed
DQB1*05:01:22	Unconfirmed	DQB1*05:18	Confirmed	DQB1*05:66:02	Unconfirmed	DQB1*05:115	Unconfirmed
DQB1*05:01:23	Unconfirmed	DQB1*05:19	Confirmed	DQB1*05:67	Unconfirmed	DQB1*05:116	Unconfirmed
DQB1*05:02:01	Confirmed	DQB1*05:20	Confirmed	DQB1*05:68	Unconfirmed	DQB1*05:117	Unconfirmed
DQB1*05:02:02	Unconfirmed	DQB1*05:21	Unconfirmed	DQB1*05:69	Confirmed	DQB1*05:118	Unconfirmed
DQB1*05:02:03	Unconfirmed	DQB1*05:22	Unconfirmed	DQB1*05:70	Unconfirmed	DQB1*05:119	Unconfirmed
DQB1*05:02:04	Unconfirmed	DQB1*05:23	Confirmed	DQB1*05:71	Confirmed	DQB1*05:120	Unconfirmed
DQB1*05:02:05	Unconfirmed	DQB1*05:24	Unconfirmed	DQB1*05:72	Confirmed	DQB1*05:121	Unconfirmed
DQB1*05:02:06	Unconfirmed	DQB1*05:25	Unconfirmed	DQB1*05:73	Confirmed	DQB1*05:122	Unconfirmed
DQB1*05:02:07	Confirmed	DQB1*05:26	Confirmed	DQB1*05:74	Unconfirmed	DQB1*05:123	Unconfirmed
DQB1*05:02:08	Unconfirmed	DQB1*05:27	Confirmed	DQB1*05:75	Confirmed	DQB1*05:124	Unconfirmed
DQB1*05:02:09	Unconfirmed	DQB1*05:28	Confirmed	DQB1*05:76	Unconfirmed	DQB1*05:125	Unconfirmed
DQB1*05:02:10	Unconfirmed	DQB1*05:29	Unconfirmed	DQB1*05:77	Unconfirmed	DQB1*05:126	Unconfirmed
DQB1*05:02:11	Confirmed	DQB1*05:30	Unconfirmed	DQB1*05:78	Unconfirmed	DQB1*05:127	Unconfirmed
DQB1*05:02:12	Unconfirmed	DQB1*05:31	Unconfirmed	DQB1*05:79	Confirmed	DQB1*05:128N	Unconfirmed
DQB1*05:02:13	Unconfirmed	DQB1*05:32	Unconfirmed	DQB1*05:80	Confirmed	DQB1*05:129	Unconfirmed
DQB1*05:03:01:01	Confirmed	DQB1*05:33	Unconfirmed	DQB1*05:81	Unconfirmed	DQB1*05:130	Unconfirmed
DQB1*05:03:01:02	Unconfirmed	DQB1*05:34	Confirmed	DQB1*05:82	Unconfirmed	DQB1*05:131	Unconfirmed
DQB1*05:03:02	Confirmed	DQB1*05:35	Unconfirmed	DQB1*05:83	Unconfirmed	DQB1*05:132Q	Unconfirmed
DQB1*05:03:03	Unconfirmed	DQB1*05:36	Unconfirmed	DQB1*05:84	Unconfirmed	DQB1*05:133	Confirmed
DQB1*05:03:04	Unconfirmed	DQB1*05:37	Unconfirmed	DQB1*05:85	Unconfirmed	DQB1*05:134	Unconfirmed
DQB1*05:03:05	Unconfirmed	DQB1*05:38	Unconfirmed	DQB1*05:86	Unconfirmed	DQB1*05:135	Unconfirmed
DQB1*05:03:06	Unconfirmed	DQB1*05:39	Unconfirmed	DQB1*05:87Q	Unconfirmed	DQB1*05:136	Unconfirmed
DQB1*05:03:07	Unconfirmed	DQB1*05:40	Unconfirmed	DQB1*05:88	Confirmed	DQB1*05:137	Unconfirmed
DQB1*05:03:08	Confirmed	DQB1*05:41N	Unconfirmed	DQB1*05:89:01	Unconfirmed	DQB1*05:138	Unconfirmed
DQB1*05:03:09	Unconfirmed	DQB1*05:42	Unconfirmed	DQB1*05:89:02	Unconfirmed	DQB1*05:139	Unconfirmed
DQB1*05:03:10	Unconfirmed	DQB1*05:43:01	Confirmed	DQB1*05:90N	Unconfirmed	DQB1*05:140	Unconfirmed
DQB1*05:03:11	Unconfirmed	DQB1*05:43:02	Unconfirmed	DQB1*05:91	Unconfirmed	DQB1*05:141	Unconfirmed

Allele	Status ¹
DQB1*05:142	Unconfirmed
DQB1*05:143	Unconfirmed



Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DQB1*06:01:01	Confirmed	DQB1*06:03:08	Unconfirmed	DQB1*06:13:02	Confirmed	DQB1*06:53:01	Confirmed
DQB1*06:01:02	Confirmed	DQB1*06:03:09	Confirmed	DQB1*06:14:01	Confirmed	DQB1*06:53:02	Unconfirmed
DQB1*06:01:03	Unconfirmed	DQB1*06:03:10	Unconfirmed	DQB1*06:14:02	Unconfirmed	DQB1*06:54N	Unconfirmed
DQB1*06:01:04	Unconfirmed	DQB1*06:03:11	Confirmed	DQB1*06:14:03	Confirmed	DQB1*06:55	Confirmed
DQB1*06:01:05	Unconfirmed	DQB1*06:03:12	Confirmed	DQB1*06:15:01	Unconfirmed	DQB1*06:56	Unconfirmed
DQB1*06:01:06	Unconfirmed	DQB1*06:03:13	Unconfirmed	DQB1*06:15:02	Confirmed	DQB1*06:57	Confirmed
DQB1*06:01:07	Unconfirmed	DQB1*06:03:14	Confirmed	DQB1*06:16	Unconfirmed	DQB1*06:58	Confirmed
DQB1*06:01:08	Confirmed	DQB1*06:03:15	Unconfirmed	DQB1*06:17	Unconfirmed	DQB1*06:59	Confirmed
DQB1*06:01:09	Unconfirmed	DQB1*06:03:16	Unconfirmed	DQB1*06:18:01	Unconfirmed	DQB1*06:60	Unconfirmed
DQB1*06:01:10	Confirmed	DQB1*06:03:17	Unconfirmed	DQB1*06:18:02	Unconfirmed	DQB1*06:61	Unconfirmed
DQB1*06:01:11	Unconfirmed	DQB1*06:03:18	Unconfirmed	DQB1*06:19:01	Unconfirmed	DQB1*06:62	Confirmed
DQB1*06:01:12	Unconfirmed	DQB1*06:03:19	Unconfirmed	DQB1*06:19:02	Unconfirmed	DQB1*06:63	Unconfirmed
DQB1*06:01:13	Confirmed	DQB1*06:03:20	Unconfirmed	DQB1*06:20	Unconfirmed	DQB1*06:64	Confirmed
DQB1*06:01:14	Confirmed	DQB1*06:03:21	Confirmed	DQB1*06:21	Confirmed	DQB1*06:65	Unconfirmed
DQB1*06:01:15	Unconfirmed	DQB1*06:03:22	Unconfirmed	DQB1*06:22:01	Unconfirmed	DQB1*06:66	Unconfirmed
DQB1*06:02:01:01	Confirmed	DQB1*06:03:23	Confirmed	DQB1*06:22:02	Confirmed	DQB1*06:67	Unconfirmed
DQB1*06:02:01:02	Confirmed	DQB1*06:03:24	Unconfirmed	DQB1*06:22:03	Confirmed	DQB1*06:68	Confirmed
DQB1*06:02:01:03	Unconfirmed	DQB1*06:03:25	Unconfirmed	DQB1*06:23	Unconfirmed	DQB1*06:69	Confirmed
DQB1*06:02:02	Unconfirmed	DQB1*06:03:26	Unconfirmed	DQB1*06:24	Unconfirmed	DQB1*06:70	Unconfirmed
DQB1*06:02:03	Unconfirmed	DQB1*06:04:01	Confirmed	DQB1*06:25	Confirmed	DQB1*06:71	Confirmed
DQB1*06:02:04	Confirmed	DQB1*06:04:02	Unconfirmed	DQB1*06:26N	Confirmed	DQB1*06:72	Unconfirmed
DQB1*06:02:05	Confirmed	DQB1*06:04:03	Unconfirmed	DQB1*06:27:01	Unconfirmed	DQB1*06:73	Unconfirmed
DQB1*06:02:06	Unconfirmed	DQB1*06:04:04	Unconfirmed	DQB1*06:27:02	Confirmed	DQB1*06:74	Confirmed
DQB1*06:02:07	Confirmed	DQB1*06:04:05	Unconfirmed	DQB1*06:28	Unconfirmed	DQB1*06:75N	Confirmed
DQB1*06:02:08	Unconfirmed	DQB1*06:04:06	Confirmed	DQB1*06:29	Unconfirmed	DQB1*06:76	Unconfirmed
DQB1*06:02:09	Unconfirmed	DQB1*06:04:07	Unconfirmed	DQB1*06:30	Unconfirmed	DQB1*06:77N	Confirmed
DQB1*06:02:10	Unconfirmed	DQB1*06:04:08	Unconfirmed	DQB1*06:31	Unconfirmed	DQB1*06:78	Confirmed
DQB1*06:02:11	Confirmed	DQB1*06:04:09	Unconfirmed	DQB1*06:32:01	Unconfirmed	DQB1*06:79:01	Confirmed
DQB1*06:02:12	Unconfirmed	DQB1*06:04:10	Unconfirmed	DQB1*06:32:02	Unconfirmed	DQB1*06:79:02	Unconfirmed
DQB1*06:02:13	Unconfirmed	DQB1*06:05:01	Confirmed	DQB1*06:33	Confirmed	DQB1*06:80	Unconfirmed
DQB1*06:02:14	Confirmed	DQB1*06:05:02	Unconfirmed	DQB1*06:34	Confirmed	DQB1*06:81	Confirmed
DQB1*06:02:15	Unconfirmed	DQB1*06:06	Unconfirmed	DQB1*06:35	Unconfirmed	DQB1*06:82	Unconfirmed
DQB1*06:02:16	Unconfirmed	DQB1*06:07:01	Confirmed	DQB1*06:36	Confirmed	DQB1*06:83	Unconfirmed
DQB1*06:02:17	Confirmed	DQB1*06:07:02	Confirmed	DQB1*06:37	Confirmed	DQB1*06:84	Confirmed
DQB1*06:02:18	Unconfirmed	DQB1*06:08:01	Confirmed	DQB1*06:38	Unconfirmed	DQB1*06:85	Unconfirmed
DQB1*06:02:19	Unconfirmed	DQB1*06:08:02	Unconfirmed	DQB1*06:39	Unconfirmed	DQB1*06:86	Unconfirmed
DQB1*06:02:20	Confirmed	DQB1*06:08:03	Unconfirmed	DQB1*06:40	Unconfirmed	DQB1*06:87	Unconfirmed
DQB1*06:02:21	Unconfirmed	DQB1*06:09:01:01	Confirmed	DQB1*06:41	Confirmed	DQB1*06:88	Confirmed
DQB1*06:02:22	Confirmed	DQB1*06:09:01:02	Unconfirmed	DQB1*06:42	Unconfirmed	DQB1*06:89	Unconfirmed
DQB1*06:02:23	Confirmed	DQB1*06:09:02	Confirmed	DQB1*06:43	Unconfirmed	DQB1*06:90	Confirmed
DQB1*06:02:24	Unconfirmed	DQB1*06:09:03	Unconfirmed	DQB1*06:44	Confirmed	DQB1*06:91	Unconfirmed
DQB1*06:02:25	Unconfirmed	DQB1*06:09:04	Unconfirmed	DQB1*06:45	Unconfirmed	DQB1*06:92	Confirmed
DQB1*06:02:26	Unconfirmed	DQB1*06:09:05	Unconfirmed	DQB1*06:46	Confirmed	DQB1*06:93	Confirmed
DQB1*06:03:01	Confirmed	DQB1*06:09:06	Confirmed	DQB1*06:47	Unconfirmed	DQB1*06:94	Confirmed
DQB1*06:03:02	Confirmed	DQB1*06:10	Confirmed	DQB1*06:48	Confirmed	DQB1*06:95	Unconfirmed
DQB1*06:03:03	Unconfirmed	DQB1*06:11:01	Confirmed	DQB1*06:49	Confirmed	DQB1*06:96	Unconfirmed
DQB1*06:03:04	Unconfirmed	DQB1*06:11:02	Confirmed	DQB1*06:50	Confirmed	DQB1*06:97	Unconfirmed
DQB1*06:03:05	Confirmed	DQB1*06:11:03	Unconfirmed	DQB1*06:51:01	Unconfirmed	DQB1*06:98	Unconfirmed
DQB1*06:03:06	Unconfirmed	DQB1*06:12	Unconfirmed	DQB1*06:51:02	Unconfirmed	DQB1*06:99:01	Unconfirmed
DQB1*06:03:07	Confirmed	DQB1*06:13:01	Unconfirmed	DQB1*06:52	Unconfirmed	DQB1*06:99:02	Unconfirmed



Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
DQB1*06:100	Unconfirmed	DQB1*06:139	Unconfirmed	DQB1*06:178	Unconfirmed	DQB1*06:216N	Unconfirmed
DQB1*06:101	Unconfirmed	DQB1*06:140	Unconfirmed	DQB1*06:179N	Confirmed	DQB1*06:217	Confirmed
DQB1*06:102N	Unconfirmed	DQB1*06:141	Unconfirmed	DQB1*06:180	Unconfirmed	DQB1*06:218	Confirmed
DQB1*06:103	Confirmed	DQB1*06:142	Unconfirmed	DQB1*06:181	Unconfirmed	DQB1*06:219	Unconfirmed
DQB1*06:104	Unconfirmed	DQB1*06:143	Unconfirmed	DQB1*06:182	Unconfirmed	DQB1*06:221	Unconfirmed
DQB1*06:105	Unconfirmed	DQB1*06:144N	Unconfirmed	DQB1*06:183	Unconfirmed	DQB1*06:222	Unconfirmed
DQB1*06:106	Confirmed	DQB1*06:145	Unconfirmed	DQB1*06:184	Confirmed	DQB1*06:223	Unconfirmed
DQB1*06:107	Confirmed	DQB1*06:146:01	Unconfirmed	DQB1*06:185	Unconfirmed	DQB1*06:224	Confirmed
DQB1*06:108	Unconfirmed	DQB1*06:146:02	Unconfirmed	DQB1*06:186	Unconfirmed	DQB1*06:225	Unconfirmed
DQB1*06:109	Unconfirmed	DQB1*06:147	Unconfirmed	DQB1*06:187	Confirmed	DQB1*06:226	Unconfirmed
DQB1*06:110	Unconfirmed	DQB1*06:148	Unconfirmed	DQB1*06:188	Confirmed	DQB1*06:227	Unconfirmed
DQB1*06:111	Confirmed	DQB1*06:149	Confirmed	DQB1*06:189	Unconfirmed	DQB1*06:228	Unconfirmed
DQB1*06:112N	Unconfirmed	DQB1*06:150	Unconfirmed	DQB1*06:190:01	Confirmed	DQB1*06:229	Unconfirmed
DQB1*06:113	Unconfirmed	DQB1*06:151	Unconfirmed	DQB1*06:190:02	Unconfirmed	DQB1*06:230	Unconfirmed
DQB1*06:114	Unconfirmed	DQB1*06:152	Unconfirmed	DQB1*06:191	Confirmed		
DQB1*06:115	Unconfirmed	DQB1*06:153	Confirmed	DQB1*06:192	Unconfirmed		
DQB1*06:116	Unconfirmed	DQB1*06:154	Unconfirmed	DQB1*06:193N	Unconfirmed		
DQB1*06:117	Confirmed	DQB1*06:155	Confirmed	DQB1*06:194	Confirmed		
DQB1*06:118:01	Unconfirmed	DQB1*06:156	Unconfirmed	DQB1*06:195	Unconfirmed		
DQB1*06:118:02	Confirmed	DQB1*06:157	Unconfirmed	DQB1*06:196	Unconfirmed		
DQB1*06:119	Unconfirmed	DQB1*06:158N	Confirmed	DQB1*06:197	Unconfirmed		
DQB1*06:120	Confirmed	DQB1*06:159	Unconfirmed	DQB1*06:198	Unconfirmed		
DQB1*06:121	Unconfirmed	DQB1*06:160	Unconfirmed	DQB1*06:199	Unconfirmed		
DQB1*06:122	Confirmed	DQB1*06:161	Unconfirmed	DQB1*06:200	Unconfirmed		
DQB1*06:123	Confirmed	DQB1*06:162	Confirmed	DQB1*06:201	Unconfirmed		
DQB1*06:124	Unconfirmed	DQB1*06:163	Confirmed	DQB1*06:202	Confirmed		
DQB1*06:125	Confirmed	DQB1*06:164	Confirmed	DQB1*06:203	Unconfirmed		
DQB1*06:126	Confirmed	DQB1*06:165	Unconfirmed	DQB1*06:204	Unconfirmed		
DQB1*06:127	Unconfirmed	DQB1*06:166	Unconfirmed	DQB1*06:205	Unconfirmed		
DQB1*06:128	Unconfirmed	DQB1*06:167	Unconfirmed	DQB1*06:206:01	Unconfirmed		
DQB1*06:129	Unconfirmed	DQB1*06:168	Unconfirmed	DQB1*06:206:02	Unconfirmed		
DQB1*06:130	Unconfirmed	DQB1*06:169	Confirmed	DQB1*06:207	Unconfirmed		
DQB1*06:131	Unconfirmed	DQB1*06:170	Confirmed	DQB1*06:208	Unconfirmed		
DQB1*06:132	Unconfirmed	DQB1*06:171	Confirmed	DQB1*06:209	Unconfirmed		
DQB1*06:133	Unconfirmed	DQB1*06:172	Confirmed	DQB1*06:210	Confirmed		
DQB1*06:134	Unconfirmed	DQB1*06:173	Unconfirmed	DQB1*06:211	Unconfirmed		
DQB1*06:135	Unconfirmed	DQB1*06:174	Unconfirmed	DQB1*06:212	Confirmed		
DQB1*06:136	Unconfirmed	DQB1*06:175	Unconfirmed	DQB1*06:213	Unconfirmed		
DQB1*06:137	Unconfirmed	DQB1*06:176	Unconfirmed	DQB1*06:214	Unconfirmed		
DQB1*06:138	Unconfirmed	DQB1*06:177	Unconfirmed	DQB1*06:215	Unconfirmed		

¹Allele status “confirmed” or “unconfirmed” as listed on the IMGT/HLA web page 2017-January-20, release 3.27.0, www.ebi.ac.uk/imgt/hla.

RESOLUTION IN HOMO- AND HETEROZYGOTES

Results file with resolution in DQB1 high homo- and heterozygotes is available upon request.



Lot No.: **2F4**

Lot-specific information

SPECIFICITY TABLE

DQB1 high resolution for frequent alleles SSP typing

Specificities and sizes of the PCR products of the 95+1 primer mixes used for DQB1 high resolution for frequent alleles SSP typing

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified DQB1 alleles ³
1⁴	120 bp	515 bp	*02:01:01-02:02:04, 02:04-02:16, 02:18N-02:24, 02:26-02:39, 02:41-02:71, 02:73-02:76, 02:78-02:88
2^{4,7}	85 bp	430 bp	*02:03, 02:77
3⁶	145 bp	515 bp	*02:01:01-02:01:24, 02:04-02:05, 02:07:01-02:09, 02:13-02:25, 02:27-02:47, 02:49, 02:51-02:61, 02:63, 02:66-02:79, 02:81-02:83, 02:85-02:88, 03:01:01:01-03:23:02, 03:25:01-03:78, 03:80-03:96, 03:98-03:163, 03:166-03:256, 04:01:01-04:39, 05:01:01:01-05:13, 05:15-05:83, 05:85-05:143, 06:01:01-06:37, 06:39-06:85, 06:87-06:101, 06:105-06:217, 06:219, 06:221-06:226, 06:228-06:230
4	140 bp	430 bp	*02:02:01:01-02:03, 02:06, 02:10-02:12, 02:26, 02:50, 02:62, 02:64-02:65
	185 bp		*03:28, 03:50
5⁵	155 bp	430 bp	*02:04, 03:22, 03:96, 03:251, 04:23, 05:103, 06:04:01-06:09:06, 06:11:01-06:13:02, 06:15:01-06:32:02, 06:34-06:36, 06:38, 06:40-06:42, 06:45-06:46, 06:49-06:50, 06:51:02-06:54N, 06:56-06:61, 06:63-06:76, 06:79:02-06:85, 06:89, 06:91-06:97, 06:118:01-06:124, 06:126, 06:128-06:129, 06:131-06:147, 06:149-06:184, 06:186, 06:189-06:199, 06:201-06:204, 06:206:01-06:208, 06:210-06:215, 06:217, 06:229-06:230
	180 bp		*05:06:01-05:06:02, 05:07, 05:50
	210 bp		*05:106
	260 bp		*05:03:02, 05:51, 05:119, 05:125, 06:06, 06:129, 06:146:01
6	165 bp	515 bp	*03:01:01:01-03:01:07, 03:01:09-03:01:33, 03:01:35-03:01:37, 03:04:01-03:04:03, 03:09, 03:13, 03:16, 03:19:01-03:19:02, 03:21-03:22, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:60, 03:69, 03:71, 03:73-03:78, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108, 03:113-03:116, 03:118N-03:122, 03:127-03:131, 03:134-03:135, 03:139-03:140, 03:142-03:144, 03:147-03:148, 03:150-03:152, 03:157-03:160, 03:162-03:167, 03:169-03:173, 03:182, 03:184, 03:188, 03:191-03:194, 03:196-03:198, 03:201-03:202, 03:206-03:208, 03:216, 03:218-03:219, 03:231-03:232, 03:235-03:236, 03:241-03:243, 03:246, 03:252-03:256, 04:09
7	220 bp	515 bp	*03:01:01:01-03:01:33, 03:01:35-03:01:37, 03:04:01-03:04:03, 03:09-03:10:02, 03:13-03:14:02, 03:16, 03:19:01-03:19:02, 03:21-03:22, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:60, 03:69, 03:71, 03:73, 03:75-03:77, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108, 03:114-03:116, 03:118N-03:122, 03:127-03:131, 03:133-03:135, 03:138-03:140, 03:142-03:144, 03:147-03:148, 03:150, 03:152, 03:157-03:160, 03:162-03:167, 03:169-03:173, 03:180, 03:182-03:183, 03:186-03:188, 03:191-03:198, 03:201-03:202, 03:206-03:208, 03:216, 03:218-03:219, 03:231-03:232, 03:235-03:236, 03:241-03:243, 03:246, 03:252-03:256

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

8⁶	135 bp 225 bp	515 bp	*02:01:01-02:02:04, 02:04-02:16, 02:18N-02:36, 02:38-02:76, 02:78-02:88, 03:02:01:01-03:02:24, 03:07-03:08, 03:11, 03:18, 03:32, 03:37, 03:45, 03:62-03:64, 03:66N-03:68, 03:70, 03:81, 03:85, 03:106-03:107, 03:125, 03:146, 03:153, 03:161, 03:174-03:175, 03:178-03:179, 03:184-03:185, 03:189-03:190, 03:199, 03:203-03:205, 03:210-03:211, 03:213N-03:215, 03:220-03:221, 03:223-03:225, 03:228-03:229, 03:233, 03:237N, 03:240, 03:245, 03:247, 03:251, 06:03:13, 06:29, 06:63, 06:123, 06:139 *06:171
9	175 bp 220 bp	515 bp	*03:02:01:01-03:03:13, 03:05:01-03:05:04, 03:06 [?] -03:08 [?] , 03:11 [?] -03:15 [?] , 03:17:01 [?] -03:18 [?] , 03:19:01-03:19:02, 03:20 [?] , 03:23:01 [?] -03:23:02 [?] , 03:25:01, 03:25:02 [?] -03:26 [?] , 03:30-03:32, 03:33 [?] -03:34 [?] , 03:37 [?] , 03:38-03:39, 03:40 [?] , 03:41, 03:43, 03:45, 03:48 [?] , 03:52 [?] -03:71 [?] , 03:72, 03:74 [?] -03:78 [?] , 03:79, 03:81 [?] -03:82 [?] , 03:85-03:91Q, 03:95N-03:100, 03:101 [?] -03:104 [?] , 03:105, 03:106 [?] -03:112 [?] , 03:117, 03:118N [?] , 03:120 [?] -03:131 [?] , 03:132, 03:133 [?] -03:137 [?] , 03:138, 03:140 [?] -03:149 [?] , 03:150, 03:151 [?] -03:163 [?] , 03:165 [?] -03:167 [?] , 03:168, 03:170 [?] -03:179 [?] , 03:180-03:181, 03:183 [?] -03:185 [?] , 03:187 [?] -03:189 [?] , 03:190, 03:192 [?] -03:194 [?] , 03:195, 03:199-03:200, 03:201 [?] -03:205 [?] , 03:207 [?] -03:210 [?] , 03:211, 03:212 [?] -03:221 [?] , 03:222, 03:223 [?] -03:238 [?] , 03:239, 03:240 [?] , 03:244-03:245, 03:247-03:248, 03:250-03:251, 03:255-03:256, 04:01:01-04:39 *02:01:01-02:36, 02:38-02:88, 03:02:01:01-03:03:13, 03:06-03:08, 03:11-03:12, 03:15, 03:18, 03:20, 03:23:01-03:23:02, 03:25:01-03:26, 03:30-03:34, 03:37-03:41, 03:43, 03:45, 03:62-03:68, 03:70, 03:74, 03:79, 03:81, 03:85-03:91Q, 03:95N-03:99Q, 03:106-03:107, 03:110-03:113, 03:117, 03:123-03:126, 03:136-03:137, 03:141, 03:145-03:146, 03:149, 03:153, 03:155-03:156, 03:161, 03:168, 03:174-03:179, 03:184-03:185, 03:190, 03:199-03:200, 03:203-03:205, 03:209-03:215, 03:217, 03:220-03:225, 03:227-03:230, 03:233-03:234, 03:237N-03:240, 03:245, 03:247-03:249, 03:251, 04:03:01-04:03:02
10⁴	120 bp 145 bp 175 bp	430 bp	*02:11, 02:25 *02:26, 03:04:01-03:04:03, 03:14:01-03:14:02, 03:80, 03:138 *03:58
11⁴	95 bp 130 bp	430 bp	*03:20 *03:05:01-03:05:04, 03:17:01-03:17:02, 03:53, 03:61-03:62, 03:72, 03:100, 03:132, 03:181, 03:226, 03:244, 03:250, 04:09, 05:131
12^{4,7}	65 bp 105 bp 140 bp 195 bp	430 bp	*03:40 *03:03:06, 03:06, 03:25:01-03:25:02, 03:234, 03:239, 04:01:03 *03:08, 03:119, 03:137, 03:194, 06:02:02, 06:03:02, 06:210 *03:57
13⁴	110 bp 145 bp 190 bp	430 bp	*03:15 *02:62, 03:07, 03:16 *03:54, 03:162, 03:183, 03:204, 03:223, 03:233, 04:35
14	135 bp 180 bp 260 bp	430 bp	*03:09, 05:56 *03:55, 03:204, 03:250 *03:11, 03:26, 03:244
15	140 bp 195 bp	430 bp	*03:10:01-03:10:02, 03:12, 03:14:01-03:14:02, 03:70, 03:179, 03:183, 03:195, 06:01:01-06:01:15, 06:06, 06:43, 06:54N-06:58, 06:98, 06:99:02-06:105, 06:108, 06:120, 06:132, 06:140, 06:153, 06:157, 06:167-06:168, 06:177, 06:181, 06:194, 06:205, 06:209, 06:214, 06:229 *02:42, 02:78

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

	245 bp		*03:27, 03:82
16	260 bp	515 bp	*03:01:01:01-03:01:07, 03:01:09-03:01:33, 03:01:35-03:01:37, 03:04:01-03:04:03, 03:09-03:10:02, 03:12-03:14:02, 03:16, 03:19:01-03:19:02, 03:21-03:22, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:60, 03:70-03:71, 03:73-03:78, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108, 03:114-03:116, 03:118N-03:122, 03:127-03:131, 03:133-03:135, 03:139-03:140, 03:143-03:144, 03:147-03:148, 03:150-03:152, 03:157-03:160, 03:162-03:167, 03:169-03:173, 03:179, 03:182-03:183, 03:186-03:188, 03:191-03:198, 03:201-03:202, 03:206, 03:208, 03:216, 03:218-03:219, 03:231-03:232, 03:235-03:236, 03:241-03:243, 03:246, 03:252-03:256
17^{4,6}	110 bp 155 bp 185 bp 215 bp 300 bp	430 bp	*05:03:12, 05:43:02, 06:23, 06:82, 06:99:01, 06:156, 06:162, 06:169 *03:13, 04:32 *03:28, 03:232 *03:47 *06:34
18^{4,7}	80 bp 185 bp 220 bp 260 bp	430 bp	*03:188 *03:60, 03:73, 03:112, 03:118N, 03:226, 03:237N, 04:30, 05:73, 05:98, 05:116, 06:28, 06:56, 06:79:01-06:79:02, 06:144N *03:18 *02:36
19	135 bp	515 bp	*03:01:01:01-03:01:07, 03:01:09-03:01:33, 03:01:35-03:01:36, 03:02:01:01-03:15, 03:17:01-03:22, 03:24, 03:26-03:103, 03:106-03:108, 03:110-03:152, 03:155-03:171, 03:173-03:188, 03:190-03:204, 03:206-03:213N, 03:215-03:216, 03:218-03:233, 03:235-03:256, 04:09
20⁷	130 bp	430 bp	*02:03, 02:77, 03:01:01:01-03:01:02, 03:01:04-03:01:07, 03:01:09-03:01:12, 03:01:14, 03:01:16-03:01:24, 03:01:26-03:01:33, 03:01:35-03:01:37, 03:03:02:01-03:03:13, 03:06, 03:09-03:10:02, 03:12-03:13, 03:15-03:17:02, 03:19:01-03:24, 03:26-03:31, 03:33-03:36, 03:38-03:44, 03:46-03:60, 03:65, 03:69, 03:71-03:79, 03:82-03:84N, 03:86-03:103, 03:108, 03:111-03:124, 03:126-03:129, 03:131, 03:133-03:137, 03:139-03:145, 03:147-03:152, 03:155-03:160, 03:162-03:171, 03:173, 03:176-03:177, 03:180, 03:182-03:183, 03:186-03:188, 03:191-03:198, 03:200-03:201, 03:206-03:209, 03:212, 03:216-03:219, 03:222, 03:227, 03:230-03:232, 03:235-03:236, 03:238, 03:241-03:244, 03:246, 03:248-03:249, 03:252-03:256, 04:09-04:10
21^{4,7}	125 bp 205 bp	515 bp	*04:11 *04:01:01-04:01:04, 04:05-04:08, 04:14-04:17, 04:27, 04:38
22	140 bp 205 bp	430 bp	*04:14 *03:132, 04:02:01:01-04:02:01:03, 04:02:02-04:02:07, 04:02:09-04:03:01, 04:04, 04:09-04:13, 04:18-04:19, 04:21-04:37, 04:39
23	135 bp	430 bp	*05:01:01:01-05:01:23, 05:07, 05:11:01-05:12, 05:18, 05:20, 05:22, 05:25, 05:27, 05:29-05:32, 05:44-05:45, 05:48-05:49, 05:51, 05:54-05:55, 05:61-05:63, 05:68-05:76, 05:80-05:81, 05:84, 05:88-05:89:02, 05:92-05:93, 05:95, 05:99, 05:103-05:104, 05:107, 05:110N-05:112, 05:114-05:115, 05:120, 05:122, 05:124, 05:126-05:128N, 05:133, 05:137-05:139, 05:141
24⁴	115 bp	430 bp	*05:02:01-05:02:13, 05:05:01-05:05:02, 05:14, 05:17, 05:19, 05:26, 05:33-05:37, 05:46-05:47, 05:53, 05:57, 05:59, 05:65, 05:79, 05:83, 05:86-05:88, 05:90N, 05:94, 05:100, 05:102, 05:106, 05:113, 05:116-05:118, 05:123, 05:125, 05:135-05:136, 05:142

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Lot No.	Allele	Product	Allele	Product	Allele	Product
	215 bp					*05:72
25⁴	100 bp	515 bp	*03:226, 04:28, 05:01:09-05:01:10, 05:02:02, 05:03:01:01-05:03:11, 05:03:13, 05:03:15-05:03:17, 05:06:01-05:06:02, 05:08-05:10, 05:13, 05:15-05:16, 05:24, 05:28, 05:38-05:42, 05:50, 05:56, 05:58, 05:60, 05:64, 05:66:01-05:67, 05:78, 05:85, 05:91, 05:96, 05:98, 05:101, 05:108-05:109, 05:119, 05:121, 05:129-05:130, 05:140, 05:143			
	125 bp		*05:03:01:01-05:03:01:02, 05:03:03-05:03:17, 05:06:01-05:06:02, 05:08-05:10, 05:13, 05:15-05:16, 05:23-05:24, 05:28, 05:38-05:43:01, 05:50, 05:56, 05:58, 05:60, 05:64, 05:66:01-05:67, 05:78, 05:85, 05:91, 05:96, 05:98, 05:101, 05:108-05:109, 05:119, 05:121, 05:129-05:131, 05:134, 05:140			
26^{4,5}	220 bp	430 bp	*05:72, 05:128N			
	120 bp		*02:17, 02:32, 06:02:01:01-06:03:26, 06:06 [?] -06:07:01 [?] , 06:08:01 [?] -06:08:03 [?] , 06:09:02 [?] -06:09:03 [?] , 06:09:05 [?] -06:09:06 [?] , 06:10, 06:11:01 [?] -06:11:03 [?] , 06:13:01 [?] -06:13:02 [?] , 06:14:01-06:14:03, 06:15:01 [?] -06:18:01 [?] , 06:19:01 [?] -06:20 [?] , 06:22:01 [?] -06:32:02 [?] , 06:33, 06:35 [?] , 06:37, 06:39, 06:40 [?] , 06:44, 06:45 [?] -06:46 [?] , 06:47-06:48, 06:49 [?] -06:50 [?] , 06:51:02 [?] , 06:53:01 [?] -06:54N [?] , 06:56 [?] -06:57 [?] , 06:59 [?] -06:61 [?] , 06:62, 06:63 [?] -06:76 [?] , 06:77N-06:79:02, 06:80 [?] -06:83 [?] , 06:85 [?] , 06:87-06:88, 06:89 [?] , 06:90, 06:91 [?] -06:97 [?] , 06:106-06:117, 06:118:02 [?] -06:124 [?] , 06:125, 06:126 [?] , 06:127, 06:128 [?] , 06:130, 06:131 [?] -06:147 [?] , 06:148, 06:149 [?] -06:184 [?] , 06:185, 06:187-06:188, 06:190:01 [?] , 06:190:02, 06:191 [?] -06:199 [?] , 06:200, 06:201 [?] -06:204 [?] , 06:206:01 [?] , 06:207 [?] -06:208 [?] , 06:210 [?] -06:215 [?] , 06:216N, 06:218-06:219, 06:221-06:222, 06:224-06:230			
27	170 bp	430 bp	*03:67			
	210 bp		*06:01:15, 06:02:01:01-06:02:26, 06:05:02 [?] -06:06 [?] , 06:10-06:11:03, 06:13:01-06:13:02, 06:16, 06:18:01-06:20, 06:24, 06:29, 06:33, 06:35, 06:37, 06:43, 06:45, 06:47-06:51:02, 06:53:01-06:57, 06:68, 06:70-06:84, 06:95-06:109, 06:111-06:117, 06:120, 06:122-06:127, 06:130-06:132, 06:136-06:140, 06:147, 06:150-06:153, 06:157, 06:159, 06:161, 06:163, 06:166-06:167, 06:173-06:179N, 06:181-06:183, 06:188, 06:192, 06:197-06:198, 06:200-06:201, 06:205, 06:208-06:209, 06:211, 06:213-06:216N, 06:219, 06:224-06:229			
28	190 bp	430 bp	*04:10, 06:02:01:01-06:02:26, 06:14:01-06:16, 06:19:01-06:20, 06:23-06:24, 06:33, 06:37, 06:46-06:50, 06:51:02, 06:68, 06:70-06:84, 06:95, 06:97, 06:107, 06:109, 06:111-06:117, 06:122, 06:124-06:127, 06:136-06:138, 06:146:01-06:147, 06:150-06:152, 06:156, 06:159, 06:161-06:163, 06:166, 06:173-06:175, 06:178-06:179N, 06:182-06:183, 06:188, 06:192, 06:197-06:198, 06:200-06:201, 06:208, 06:211, 06:213, 06:215-06:216N, 06:219, 06:224-06:228			
29	130 bp	430 bp	*06:02:07, 06:03:01-06:03:09, 06:03:11-06:03:26, 06:07:01-06:07:02, 06:11:02-06:11:03, 06:14:01-06:14:03, 06:26N, 06:28, 06:30-06:32:02, 06:40-06:41, 06:44, 06:59-06:62, 06:64-06:65, 06:67, 06:87, 06:90-06:92, 06:110, 06:118:01-06:118:02, 06:128, 06:133-06:134, 06:141, 06:143-06:145, 06:148, 06:154, 06:165, 06:184-06:185, 06:187, 06:190:01-06:191, 06:195-06:196, 06:199, 06:203, 06:210, 06:218, 06:221-06:223, 06:230			
30^{4,5}	65 bp	430 bp	*02:67N			
	105 bp		*04:15-04:16			
	150 bp		*03:23:01, 05:03:12, 05:43:02, 06:03:01-06:03:03, 06:03:05-06:03:07, 06:03:11-06:03:25, 06:04:02, 06:07:01, 06:08:01, 06:09:02, 06:11:01-06:11:03, 06:26N, 06:28, 06:30-06:32:01, 06:40-06:41, 06:44, 06:59-06:62, 06:64-06:65, 06:67, 06:90-06:92, 06:94, 06:110, 06:118:01,			



Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

			06:128, 06:133-06:134, 06:141, 06:143-06:145, 06:148, 06:154, 06:161, 06:165, 06:169-06:170, 06:184-06:185, 06:187, 06:191, 06:195-06:196, 06:199, 06:203, 06:210, 06:218, 06:221-06:223, 06:230
	205 bp		*04:20, 04:27
	235 bp		*03:06, 03:25:01, 04:04-04:05, 04:18
31⁷	165 bp	515 bp	*06:02:07, 06:03:01-06:03:18, 06:03:20-06:03:26, 06:08:01-06:08:03, 06:11:02-06:12, 06:14:01-06:14:03, 06:21, 06:26N, 06:28, 06:31, 06:40-06:41, 06:44, 06:59, 06:61, 06:63-06:65, 06:67, 06:87, 06:90-06:92, 06:110, 06:128, 06:133-06:134, 06:141, 06:143-06:145, 06:148, 06:154, 06:165, 06:170, 06:184-06:185, 06:187, 06:190:01-06:191, 06:195-06:196, 06:199, 06:203, 06:206:01-06:206:02, 06:210, 06:218, 06:221-06:223, 06:230
32	210 bp	515 bp	*06:04:01-06:04:10, 06:07:01-06:07:02, 06:17, 06:21, 06:25, 06:34, 06:36, 06:38-06:39, 06:52, 06:58, 06:69, 06:85-06:86, 06:89, 06:92-06:93, 06:135, 06:149, 06:155, 06:158N, 06:160, 06:164, 06:172, 06:180, 06:193N, 06:202, 06:204, 06:217
33^{4,7}	110 bp 170 bp	430 bp	*03:133, 03:199, 03:200 *06:04:01-06:05:02, 06:06 ^w , 06:07:01-06:07:02, 06:09:01:01-06:09:03, 06:09:05-06:09:06, 06:18:01-06:18:02, 06:22:02, 06:25, 06:27:01-06:27:02, 06:32:01-06:32:02, 06:34, 06:36, 06:38-06:39, 06:52, 06:58, 06:66, 06:69, 06:85-06:86, 06:88-06:89, 06:93-06:94, 06:118:01-06:118:02, 06:121, 06:129, 06:135, 06:142, 06:155, 06:158N, 06:160, 06:164, 06:168, 06:171-06:172, 06:180, 06:186, 06:189, 06:193N, 06:202, 06:204, 06:207, 06:217
34⁴	95 bp 120 bp	430 bp	*05:07 *05:25, 05:75, 06:04:01-06:05:01, 06:06, 06:08:01-06:09:06, 06:12-06:13:02, 06:17-06:18:02, 06:21-06:22:03, 06:27:01-06:27:02, 06:34, 06:36, 06:38-06:39, 06:42, 06:52, 06:58, 06:69, 06:85-06:86, 06:88-06:89, 06:93, 06:119, 06:121, 06:129, 06:135, 06:142, 06:149, 06:155, 06:158N, 06:160, 06:164, 06:171, 06:180, 06:186, 06:189, 06:193N, 06:202, 06:204, 06:206:01-06:207, 06:217
35^{4,6,7}	180 bp 85 bp 155 bp 215 bp 255 bp	515 bp	*05:112 *05:103 *03:27, 03:48, 03:87, 06:35, 06:53:01-06:53:02, 06:145, 06:208 *05:101, 06:151, 06:212 *03:31, 05:73, 05:80, 05:98, 05:116, 06:05:01-06:06, 06:20, 06:31, 06:45, 06:156
36^{4,7}	95 bp 155 bp 180 bp 210 bp	430 bp	*03:90N *06:164 *03:68 *06:05:01, 06:05:02 [?] -06:06 [?] , 06:09:01:01-06:09:06, 06:12, 06:15:01-06:15:02, 06:22:01-06:22:03, 06:42, 06:46, 06:66, 06:88, 06:94, 06:118:01-06:119, 06:121, 06:142, 06:189, 06:207, 06:212 ^w
37⁵	130 bp 165 bp 185 bp	430 bp	*05:63 *05:39 *05:21, 06:06, 06:149
38^{4,6}	110 bp 140 bp 185 bp 235 bp	430 bp	*04:28, 05:38, 05:62, 05:79, 05:119-05:121, 06:14:01-06:14:03, 06:69, 06:156, 06:162, 06:168, 06:194, 06:206:01-06:206:02 *03:01:04, 05:04, 05:08, 05:18, 05:26, 05:52, 05:77, 05:132Q, 06:01:10, 06:221 *03:98, 04:33, 05:79, 06:10, 06:130, 06:188, 06:214 *04:33

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

39⁴	115 bp 140 bp 180 bp	430 bp	*03:88, 05:111 *03:165 *03:29-03:30, 03:52, 03:72, 03:100, 03:132, 03:179, 03:191, 03:215, 04:09, 04:13, 05:62, 05:120, 06:13:01-06:13:02, 06:16, 06:22:01-06:22:03, 06:29, 06:55, 06:69, 06:96, 06:119, 06:139, 06:202, 06:206:01-06:206:02, 06:209
40^{4,7}	110 bp 150 bp 205 bp	430 bp	*03:196, 03:245, 05:45, 06:40 *02:64, 02:79, 03:21-03:22, 03:96, 03:251, 04:23, 05:27, 05:35, 05:103, 06:04:01-06:09:06, 06:11:01-06:13:02, 06:15:01-06:32:02, 06:34-06:36, 06:38, 06:40-06:42, 06:45-06:46, 06:49-06:50, 06:51:02-06:54N, 06:56-06:61, 06:63-06:76, 06:79:02-06:85, 06:89, 06:91-06:97, 06:118:01-06:124, 06:126, 06:128-06:129, 06:131-06:147, 06:149-06:184, 06:186, 06:189-06:199, 06:201-06:204, 06:206:01-06:208, 06:210-06:215, 06:217, 06:229-06:230 *05:41N, 05:90N, 05:133
41⁶	165 bp 190 bp 230 bp 265 bp 300 bp	430 bp	*03:213N, 03:228, 06:17, 06:24, 06:30, 06:42 *03:45, 06:38 *02:47, 03:212 *06:127 *06:52
42^{4,6,7}	105 bp 135 bp 160 bp 190 bp 230 bp	515 bp	*03:41, 03:224, 04:06, 04:12, 04:16 *03:25:01-03:25:02, 03:66N, 03:234, 03:239, 04:03:01-04:03:02, 06:10, 06:25, 06:130, 06:193N *03:213N, 04:18 *02:08, 03:06, 03:25:01, 03:230, 04:03:01-04:03:02, 05:79 *04:18
43^{4,6,7}	85 bp 110 bp 155 bp 185 bp 215 bp 240 bp 270 bp	515 bp	*03:59, 06:37 *06:37 *02:05, 06:100 *06:26N, 06:100 *06:26N *03:71, 06:33, 06:143 *06:33, 06:143
44⁴	120 bp 165 bp	430 bp	*03:42, 03:89, 03:164 *03:94, 03:191, 04:02:07, 04:09, 05:62, 06:02:01:01-06:02:04, 06:02:06-06:02:26, 06:10, 06:13:01-06:16, 06:20, 06:22:01-06:24, 06:29, 06:33, 06:37, 06:46-06:51:02, 06:68-06:84, 06:96-06:97, 06:106-06:107, 06:109, 06:111-06:117, 06:119, 06:122, 06:124-06:127, 06:130, 06:136-06:138, 06:146:01-06:147, 06:150-06:151, 06:156, 06:159, 06:161-06:162, 06:166, 06:173-06:175, 06:178-06:179N, 06:182-06:183, 06:188, 06:192, 06:198, 06:200-06:201, 06:206:01-06:206:02, 06:208, 06:211, 06:213, 06:215-06:216N, 06:219, 06:224-06:228
45	130 bp 155 bp 225 bp 260 bp	515 bp	*06:07:01-06:07:02, 06:15:01-06:15:02, 06:36, 06:46, 06:66, 06:92, 06:118:01-06:118:02, 06:172, 06:216N, 06:224 *06:164 *06:212 *06:85
46	160 bp	515 bp	*06:03:01-06:03:14, 06:03:16-06:03:18, 06:03:20-06:03:26, 06:08:01-06:08:03, 06:14:01-06:14:03, 06:21, 06:28, 06:31, 06:40-06:41, 06:44, 06:59, 06:61, 06:63-06:65, 06:67, 06:87, 06:90-06:92, 06:110, 06:128, 06:133-06:134, 06:141, 06:143-06:145, 06:148, 06:154, 06:170, 06:184-06:185, 06:187, 06:190:01-06:191, 06:195-06:196,

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

	210 bp		06:199, 06:203, 06:206:01-06:206:02, 06:210, 06:218, 06:221-06:223, 06:230 *06:03:01-06:03:10, 06:03:12-06:03:26, 06:08:01-06:08:03, 06:14:01-06:14:03, 06:27:01-06:28, 06:30-06:32:02, 06:40-06:41, 06:44, 06:59-06:65, 06:67, 06:87, 06:90-06:91, 06:110, 06:128, 06:133-06:134, 06:141, 06:143-06:145, 06:148, 06:154, 06:168, 06:170, 06:184-06:185, 06:187, 06:190:01-06:191, 06:195-06:196, 06:206:01-06:206:02, 06:210, 06:218, 06:221-06:223, 06:230
47^{4,5}	110 bp 155 bp	430 bp	*03:42 *03:30, 03:72, 03:84N, 03:98, 03:100, 03:132, 03:215, 04:01:01-04:03:02, 04:06-04:10, 04:12-04:14, 04:16-04:36N, 04:38-04:39, 05:38, 05:119, 06:02:05, 06:02:16, 06:13:02, 06:14:03, 06:19:01-06:19:02, 06:22:01-06:22:02, 06:69, 06:95, 06:119, 06:123 ^w , 06:139, 06:152, 06:206:02 *06:02:04, 06:72
48	290 bp 170 bp 225 bp	430 bp	*03:32, 03:43, 04:07 *03:85
49^{4,6}	105 bp 140 bp 165 bp 240 bp	430 bp	*03:33, 03:196, 05:32, 05:42, 05:45, 06:228 *05:14, 05:87Q, 05:132Q *02:33, 02:41, 03:35 *03:246
50⁴	125 bp 155 bp 200 bp 230 bp	430 bp	*03:130 *03:80, 03:139 *03:37-03:38, 03:46, 03:162, 03:183, 03:223, 03:233, 03:252 *02:14:01-02:14:02, 03:34, 03:143, 03:201
51	125 bp 185 bp 270 bp	430 bp	*03:36, 03:39, 03:53, 03:62, 03:149 *03:51 *02:09
52⁴	95 bp	430 bp	*03:06 [?] -03:08 [?] , 03:11 [?] -03:15 [?] , 03:17:01 [?] -03:18 [?] , 03:19:01-03:19:02, 03:20 [?] , 03:23:01 [?] -03:23:02 [?] , 03:25:02 [?] -03:26 [?] , 03:37 [?] , 03:40 [?] , 03:48 [?] , 03:52 [?] -03:71 [?] , 03:74 [?] -03:78 [?] , 03:81 [?] -03:82 [?] , 03:101 [?] -03:104 [?] , 03:106 [?] -03:112 [?] , 03:118N [?] , 03:120 [?] -03:131 [?] , 03:133 [?] -03:137 [?] , 03:140 [?] -03:149 [?] , 03:150, 03:151 [?] -03:163 [?] , 03:165 [?] -03:167 [?] , 03:170 [?] -03:179 [?] , 03:183 [?] -03:185 [?] , 03:187 [?] -03:189 [?] , 03:192 [?] -03:194 [?] , 03:201 [?] -03:205 [?] , 03:207 [?] -03:210 [?] , 03:212 [?] -03:221 [?] , 03:223 [?] -03:238 [?] , 03:240 [?] , 03:244, 03:255-03:256, 04:08 *06:98, 06:126
53	220 bp 130 bp 160 bp 195 bp	430 bp	*03:162, 03:183, 03:223, 03:233, 06:153 *02:51 *02:07:02, 05:09, 05:29-05:30
54⁴	110 bp 150 bp 190 bp 245 bp	430 bp	*05:17, 05:23 *03:115, 06:109-06:110 *03:249, 05:10, 05:12, 05:47, 05:71, 06:105, 06:185, 06:225 *04:36N
55^{4,5}	115 bp	430 bp	*06:04:01-06:05:02, 06:06 [?] -06:07:01 [?] , 06:07:02, 06:08:01 [?] -06:08:03 [?] , 06:09:01:01-06:09:01:02, 06:09:02 [?] -06:09:03 [?] , 06:09:04, 06:09:05 [?] -06:09:06 [?] , 06:11:01 [?] -06:11:03 [?] , 06:12, 06:13:01 [?] -06:13:02 [?] , 06:14:02 [?] -06:18:01 [?] , 06:18:02, 06:19:01 [?] -06:20 [?] , 06:21, 06:22:01 [?] -06:32:02 [?] , 06:34, 06:35 [?] , 06:36, 06:38, 06:40 [?] , 06:41-06:42, 06:45 [?] -06:46 [?] , 06:49 [?] -06:50 [?] , 06:51:02 [?] , 06:52, 06:53:01 [?] -06:54N [?] , 06:56 [?] -06:57 [?] , 06:58, 06:59 [?] -06:61 [?] , 06:63 [?] -06:76 [?] , 06:79:02 [?] -06:83 [?] , 06:84, 06:85 [?] , 06:86, 06:89 [?] , 06:91 [?] -06:97 [?] , 06:118:01, 06:118:02 [?] -06:124 [?] , 06:126 [?] , 06:128 [?] , 06:129, 06:131 [?] -06:147 [?] , 06:149 [?] -06:184 [?] , 06:186, 06:189, 06:190:01 [?] , 06:190:02, 06:191 [?] -06:199 [?] , 06:201 [?] -06:204 [?] ,

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

			06:206:01 [?] , 06:206:02, 06:207 [?] -06:208 [?] , 06:210 [?] -06:215 [?] , 06:217, 06:229-06:230
56⁴	100 bp	430 bp	*06:44, 06:47
	140 bp		*03:177, 06:224
	185 bp		*03:114, 03:197Q
	220 bp		*06:43
	265 bp		*03:83
57	150 bp	430 bp	*02:54
	185 bp		*02:07:01, 03:01:01:01-03:01:35, 03:01:37, 03:04:01-03:04:03, 03:06 [?] -03:08 [?] , 03:09-03:10:02, 03:11 [?] -03:15 [?] , 03:16, 03:17:01 [?] -03:18 [?] , 03:20 [?] , 03:21-03:22, 03:23:01 [?] -03:23:02 [?] , 03:24, 03:25:02 [?] -03:26 [?] , 03:27-03:29, 03:33 [?] -03:34 [?] , 03:35-03:36, 03:37 [?] , 03:40 [?] , 03:42, 03:44, 03:46-03:47, 03:48 [?] , 03:49-03:51, 03:52 [?] -03:71 [?] , 03:73, 03:74 [?] -03:78 [?] , 03:80, 03:81 [?] -03:82 [?] , 03:83-03:84N, 03:92-03:94, 03:101 [?] -03:104 [?] , 03:106 [?] -03:112 [?] , 03:113-03:116, 03:118N [?] , 03:119, 03:120 [?] -03:131 [?] , 03:133 [?] -03:137 [?] , 03:139, 03:140 [?] -03:149 [?] , 03:151 [?] -03:163 [?] , 03:164, 03:165 [?] -03:167 [?] , 03:169, 03:170 [?] -03:179 [?] , 03:182, 03:183 [?] -03:185 [?] , 03:186, 03:187 [?] -03:189 [?] , 03:191, 03:192 [?] -03:194 [?] , 03:196-03:198, 03:201 [?] -03:205 [?] , 03:206, 03:207 [?] -03:210 [?] , 03:212 [?] -03:221 [?] , 03:223 [?] -03:238 [?] , 03:240 [?] , 03:241-03:244, 03:246, 03:252-03:256
58	145 bp	430 bp	*03:119, 03:228, 06:17, 06:24, 06:30, 06:42
	180 bp		*03:56, 05:05:01- 05:05:02, 05:11:01-05:11:02, 05:28, 05:66:01-05:66:02, 06:06, 06:129, 06:146:01
	215 bp 250 bp		*03:119, 03:228, 05:30, 06:17, 06:24, 06:30, 06:42 *03:56, 05:03:02, 05:51-05:52, 05:119, 05:125, 06:06, 06:129, 06:146:01
59^{4,5}	100 bp	430 bp	*03:44, 05:13, 05:117
	190 bp		*03:114, 06:54N
	220 bp		*05:24, 05:101
60⁴	110 bp	430 bp	*05:88
	150 bp		*05:40, 05:69
	180 bp		*05:15, 05:33, 05:49
61⁴	110 bp	430 bp	*05:88
	150 bp		*05:40, 05:69
	215 bp		*05:16, 05:34, 05:44, 05:73, 05:98, 05:116, 05:128N, 05:130, 06:156, 06:162, 06:169
62^{4,5}	90 bp	430 bp	*03:205, 03:214, 05:110N
	140 bp		*03:75, 06:97, 06:158N
	165 bp		*05:62, 06:02:19, 06:79:02, 06:95
	225 bp		*03:30, 03:72, 03:100, 03:215, 06:15:01-06:15:02, 06:22:01-06:22:03, 06:37, 06:48, 06:51:01-06:51:02, 06:69, 06:139
63⁴	95 bp	430 bp	*03:240
	180 bp		*06:50
	220 bp		*03:86
64^{6,7}	100 bp	430 bp	*02:16, 02:21
	160 bp		*02:51, 04:17
	230 bp		*02:20N
65⁴	115 bp	430 bp	*02:21, 05:36
	220 bp		*02:46
66⁴	110 bp	430 bp	*03:116, 05:57
	135 bp		*02:18N, 05:31, 05:46, 05:108, 06:111
	165 bp		*02:19

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

	205 bp		*05:41N, 05:90N, 05:133
	255 bp		*03:243
67⁴	100 bp	430 bp	*02:15, 02:27, 02:58N
	155 bp		*02:26, 02:30
68⁴	100 bp	430 bp	*02:23, 02:24, 02:27
	175 bp		*02:28
69⁴	105 bp	430 bp	*02:24, 02:31, 02:58N
	140 bp		*03:225
	185 bp		*02:29
70⁴	75 bp	430 bp	*03:61, 03:221, 06:209
71^{4,6}	50 bp	430 bp	*03:188
	135 bp		*03:63, 03:76, 03:139, 03:168, 03:211
	160 bp		*03:118N, 03:237N, 06:144N
	235 bp		*03:146, 03:207, 06:195
72	145 bp	430 bp	*03:64, 05:37, 05:104, 06:218
	180 bp		*02:50, 03:77, 03:118N, 03:237N
	210 bp		*03:99Q, 03:201
	270 bp		*03:150
73⁴	95 bp	430 bp	*02:53Q, 03:65
	125 bp		*03:23:01-03:23:02, 03:217, 04:10
	160 bp		*03:211
	225 bp		*03:78, 03:81, 03:159
	270 bp		*03:150
74⁴	115 bp	430 bp	*05:19, 05:22
	140 bp		*04:25N, 05:97, 05:143
	185 bp		*05:71
75^{4,7}	115 bp	430 bp	*05:75
	170 bp		*03:102, 03:197Q, 05:20, 05:122, 06:06, 06:55, 06:58, 06:202
	210 bp		*03:194, 06:122, 06:148, 06:161
76	150 bp	430 bp	*03:04:02, 03:10:01, 03:12, 03:14:02, 03:70, 03:179, 03:195, 06:01:01-06:01:06, 06:01:08-06:01:15, 06:43, 06:53:01-06:57, 06:98- 06:105, 06:108, 06:120, 06:132, 06:140, 06:153, 06:157, 06:167- 06:168, 06:177, 06:181, 06:194, 06:205, 06:209, 06:214, 06:229
	185 bp		*06:54N
77⁴	95 bp	430 bp	*06:155
	120 bp		*06:190:01-06:190:02
	190 bp		*06:60-06:61
78⁴	95 bp	430 bp	*03:217, 06:64
	130 bp		*03:232, 06:94, 06:170, 06:190:01-06:190:02
	170 bp		*06:62
79⁴	90 bp	430 bp	*06:65
	125 bp		*06:93
	200 bp		*06:116, 06:187
80^{4,7}	105 bp	430 bp	*06:70
	175 bp		*06:76, 06:163
81^{4,5}	120 bp	430 bp	*03:95N, 03:210
	170 bp		*03:30, 03:72, 03:100, 03:132, 03:158, 03:170, 03:215, 04:09, 06:29, 06:71, 06:96, 06:139, 06:163
82	170 bp	430 bp	*06:77N
	235 bp		*06:107
	280 bp		*06:73
83	165 bp	430 bp	*06:78

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

	235 bp		*06:107
	260 bp		*06:74
84⁴	115 bp	430 bp	*06:80, 06:217
	190 bp		*06:75N, 06:106
85	135 bp	430 bp	*06:184
	180 bp		*06:67, 06:174, 06:191
	215 bp		*06:68, 06:112N
	265 bp		*06:127
86⁴	95 bp	430 bp	*06:59, 06:90, 06:150, 06:191
	130 bp		*06:99:01-06:99:02, 06:112N, 06:210
	180 bp		*06:90, 06:191
	220 bp		*06:112N
87⁴	105 bp	430 bp	*03:08, 03:137, 03:194, 06:29, 06:59, 06:63, 06:87, 06:91, 06:96, 06:150
	140 bp		*06:179N, 06:184
88⁴	110 bp	430 bp	*06:125
	215 bp		*06:57, 06:81
89⁴	110 bp	430 bp	*06:125
	150 bp		*06:101
	185 bp		*06:144N
	260 bp		*06:83
90	130 bp	430 bp	*06:102N
	175 bp		*03:247
	240 bp		*03:117, 03:97, 06:86, 06:104, 06:117
91⁴	70 bp	430 bp	*03:10:02, 03:183, 04:02:07, 05:03:14, 05:43:01, 05:131, 06:01:12, 06:51:01, 06:66, 06:96, 06:99:01, 06:172
	100 bp		*06:01:12, 06:51:01-06:51:02, 06:99:01, 06:142
	185 bp		*03:158
92⁴	105 bp	430 bp	*02:22
	155 bp		*02:19, 02:54
	190 bp		*02:78
	230 bp		*02:12
93⁴	95 bp	430 bp	*03:91Q
	170 bp		*02:13, 02:41
	225 bp		*03:85
94⁷	170 bp	430 bp	*02:57, 02:77, 03:02:01:01-03:03:13, 03:05:01-03:08, 03:11, 03:15, 03:17:01-03:18, 03:20, 03:26, 03:30 ^w , 03:31-03:34, 03:37-03:41, 03:43, 03:45, 03:61-03:69, 03:72 ^w , 03:79, 03:81, 03:85-03:91Q, 03:95N-03:99Q, 03:100 ^w , 03:104-03:107, 03:110-03:113, 03:117, 03:123-03:126, 03:132 ^w , 03:136-03:138, 03:141, 03:145, 03:149, 03:153, 03:155-03:156, 03:161, 03:168, 03:174-03:178, 03:180-03:181, 03:184-03:185, 03:189-03:190, 03:199-03:200, 03:203-03:205, 03:209-03:213N, 03:215 ^w , 03:220-03:230, 03:233, 03:237N-03:238, 03:240, 03:244-03:245, 03:247-03:251, 04:09 ^w , 05:131, 06:29 ^w , 06:63, 06:66, 06:87, 06:96 ^w , 06:139 ^w , 06:172
95^{4,7}	95 bp	430 bp	*03:105-03:106, 03:109, 03:120
96⁸	-	-	Negative Control

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of DQB1 high resolution typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length.

Product Insert

101.221-12 – including *Taq* pol., IFU-01

Visit <https://labproducts.caredx.com> for

101.221-12u – without *Taq* pol., IFU-02

“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings. Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³For several DQB1 alleles 1st and/or 3rd exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁴HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁵Primer mixes 26, 30, 37, 47, 55, 59, 62 and 81 may give rise to a lower yield of HLA-specific PCR product than the other DQB1 primer mixes, most pronounced are primer mixes 26 and 55.

⁶Primer mixes 3, 8, 17, 35, 38, 41, 42, 43, 49, 64 and 71 have a tendency to giving rise to primer oligomer formation.

⁷Primer mixes 2, 12, 18, 20, 21, 31, 33, 35, 36, 40, 42, 43, 64, 75, 80 and 94 may have tendencies of unspecific amplifications.

⁸Primer mix 96 contains a negative control, which will amplify more than 95% of HLA amplicons as well as the amplicons generated by the control primer pairs matching the human growth hormone gene. HLA-specific PCR product sizes range from 75 to 200 base pairs and the PCR product generated by the HGH positive control primer pair is 430 base pairs.

‘?’ , nucleotide sequence information not available for the primer matching sequence.

‘w’ , might be weakly amplified.

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

PRIMER SPECIFICATION

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec. PCR product	120	85	145	140	155	165	220	135	175	120	95	65
				185	180			225	220	145	130	105
					210					175		140
					260							195
Length of int. pos. control ¹	515	430	515	430	430	515	515	515	515	430	430	430
5'-primer(s) ²	30(185) 5'-AAg 3'	57(266) 5'-TgA 3'	135(500) 5'-TgA 3'	26(173) 5'-TTA 3'	14(136) 5'-gCA 3'	45(230) 5'-ggA 3'	26(173) 5'-TTA 3'	26(173) 5'-TCT 3'	26(173) 5'-TCT 3'	23(163) 5'-AgT 3'	23(164) 5'-gCA 3'	38(208) 5'-ACT 3'
				102(400) 5'-TCT 3'	39(212) 5'-gCA 3'				140(514) 5'-CAA 3'	26(173) 5'-TTA 3'	26(173) 5'-ggg 3'	38(210) 5'-gCA 3'
					40(214) 5'-gCC 3'					30(184) 5'-gAg 3'	38(208) 5'-ACA 3'	
					40(216) 5'-TTg 3'					33(193) 5'-ATg 3'		
					113(435) 5'-AAA 3'							
					130(485) 5'-CCA 3'							
					136(502) 5'-ACg 3'							
3'-primer(s) ³	57(266) 5'-CAg 3'	71(309) 5'-CgT 3'	170(604) 5'-gAC 3'	77(325) 5'-CgC 3'	86(353) 5'-ACg 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	57(266) 5'-Cgg 3'	86(353) 5'-gCT 3'	57(266) 5'-Cgg 3'	55(260) 5'-gCg 3'	46(232) 5'-CAA 3'
	57(266) 5'-Cgg 3'		170(604) 5'-gAC 3'	135(500) 5'-ggC 3'	170(604) 5'-gAC 3'		86(353) 5'-gTT 3'	57(266) 5'-CAg 3'	185(650) 5'-CgA 3'	71(308) 5'-ggC 3'		56(263) 5'-CAA 3'
	58(270) 5'-TCA 3'			146(532) 5'-ggT 3'				57(266) 5'-Cgg 3'				62(282) 5'-CTA 3'
								58(270) 5'-TCA 3'				70(304) 5'-CCC 3'
								87(355) 5'-gTg 3'				75(320) 5'-CCC 3'
												89(362) 5'-TCA 3'
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Well No.	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec.	110	135	140	260	110	80	135	130	125	140	135	115
PCR product	145	180	195		155	185			205	205		215
	190	260	245		185	220						
					215	260						
					300							
Length of int. pos. control ¹	430	430	430	515	430	430	515	430	515	430	430	430
5'-primer(s) ²	34(197) 5'-CCA 3'	14(136) 5'-gCC 3'	12(132) 5'-AAA 3'	13(134) 5'-ggC 3'	26(173) 5'-ggg 3'	13(134) 5'-ggA 3'	55(260) 5'-gCC 3'	57(266) 5'-TgA 3'	23(164) 5'-gCT 3'	23(164) 5'-gCg 3'	26(173) 5'-ggg 3'	30(184) 5'-gAC 3'
	37(205) 5'-AgA 3'	37(205) 5'-AgA 3'	13(134) 5'-ggC 3'		26(173) 5'-TTA 3'	27(175) 5'-TTC 3'			24(166) 5'-gCT 3'	23(164) 5'-gCg 3'	26(173) 5'-ggg 3'	
	37(205) 5'-Agg 3'	41(217) 5'-TCC 3'	30(185) 5'-AAg 3'		102(400) 5'-TCT 3'	37(207) 5'-TAA 3'			49(241) 5'-ggA 3'	44(226) 5'-ACA 3'		
	48(239) 5'-CCA 3'	135(500) 5'-TgA 3'				38(209) 5'-CgT 3'						
	49(242) 5'-ggT 3'					39(211) 5'-CAT 3'						
	55(260) 5'-gCA 3'					74(316) 5'-CgA 3'						
	63(285) 5'-Agg 3'					140(514) 5'-CAA 3'						
3'-primer(s) ³	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	45(230) 5'-CCC 3'	86(353) 5'-gCT 3'	48(240) 5'-gCg 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	86(353) 5'-gCT 3'	77(327) 5'-ACg 3'	77(327) 5'-ACg 3'	57(266) 5'-CAA 3'	52(251) 5'-gCA 3'
		166(594) 5'-CgA 3'	79(332) 5'-TgT 3'		62(282) 5'-CTT 3'	87(356) 5'-ggA 3'						57(265) 5'-gCT 3'
		167(596) 5'-CAT 3'	83(343) 5'-gTg 3'		67(296) 5'-ggT 3'	170(604) 5'-gAT 3'						58(270) 5'-TCC 3'
					71(308) 5'-ggA 3'							88(358) 5'-gCC 3'
					77(325) 5'-CgC 3'							
					85(349) 5'-CAT 3'							
					142(520) 5'-AAT 3'							
					189(661) 5'-CCA 3'							
Well No.	13	14	15	16	17	18	19	20	21	22	23	24

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Well No.	25	26	27	28	29	30	31	32	33	34	35	36
Length of spec.	100	120	170	190	130	65	165	210	110	95	85	95
PCR product	125		210			105			170	120	155	155
	220					150				180	215	180
						205					255	210
						235						
Length of int. pos. control ¹	515	430	430	430	430	430	515	515	430	430	515	430
5'-primer(s) ²	30(184) 5'-gAC 3'	30(185) 5'-AAg 3'	30(184) 5'-gAT 3'	9(122) 5'-gTT 3'	27(177) 5'-gTA 3'	9(122) 5'-gTA 3'	27(177) 5'-gTA 3'	30(184) 5'-gAC 3'	27(177) 5'-gTA 3'	10(126) 5'-CAC 3'	14(136) 5'-gCC 3'	30(184) 5'-gAT 3'
		130(485) 5'-CCg 3'	30(184) 5'-gAT 3'		30(184) 5'-gAC 3'	15(140) 5'-gTA 3'	30(184) 5'-gAC 3'		44(226) 5'-ACA 3'	27(177) 5'-gTA 3'	25(170) 5'-gCA 3'	48(239) 5'-CCA 3'
						22(160) 5'-CCA 3'			48(239) 5'-TCA 3'	28(179) 5'-gAC 3'	45(230) 5'-ggA 3'	
						24(166) 5'-gCT 3'				30(185) 5'-ACC 3'	116(442) 5'-TgA 3'	
						52(251) 5'-gCT 3'				32(190) 5'-TCC 3'	146(533) 5'-CCA 3'	
						59(272) 5'-CgT 3'				39(212) 5'-gCA 3'		
3'-primer(s) ³	47(237) 5'-CgA 3'	57(265) 5'-ggT 3'	73(314) 5'-CCT 3'	57(266) 5'-CAT 3'	57(266) 5'-CAT 3'	16(144) 5'-AAC 3'	70(304) 5'-CCC 3'	86(353) 5'-ACC 3'	70(304) 5'-CCT 3'	57(266) 5'-CAA 3'	71(307) 5'-ggA 3'	48(238) 5'-Cg 3'
	47(237) 5'-CgA 3'	58(268) 5'-ggT 3'	86(353) 5'-ACg 3'	58(270) 5'-TCC 3'	58(270) 5'-TCA 3'	48(240) 5'-gCg 3'	70(304) 5'-CCC 3'	86(354) 5'-TAT 3'			81(338) 5'-TgC 3'	76(324) 5'-gTT 3'
	50(246) 5'-gTg 3'	154(558) 5'-AAA 3'	86(354) 5'-AAT 3'			77(327) 5'-ACg 3'					83(343) 5'-gTg 3'	86(353) 5'-ACC 3'
	57(266) 5'-CgT 3'		86(354) 5'-AAA 3'								86(353) 5'-ACC 3'	
	88(358) 5'-gCC 3'										87(356) 5'-ggA 3'	
	92(370) 5'-CTA 3'										130(485) 5'-ACT 3'	
											185(650) 5'-CgA 3'	
Well No.	25	26	27	28	29	30	31	32	33	34	35	36

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Well No.	37	38	39	40	41	42	43	44	45	46	47	48
Length of spec.	130	110	115	110	165	105	85	120	130	160	110	170
PCR product	165	140	140	150	190	135	110	165	155	210	155	225
	185	185	180	205	230	160	155		225		290	
		235			265	190	185		260			
					300	230	215					
							240					
							270					
Length of int.	430	430	430	430	430	515	515	430	515	515	430	430
pos. control ¹												
5'-primer(s) ²	25(170) 5'-gCA 3'	8(118) 5'-TCA 3'	9(121) 5'-TgC 3'	31(187) 5'-ACg 3'	26(173) 5'-TCT 3'	8(118) 5'-TCA 3'	10(126) 5'-CAT 3'	9(122) 5'-gTT 3'	13(134) 5'-ggA 3'	30(184) 5'-gAC 3'	9(122) 5'-gTT 3'	123(464) 5'-CTT 3'
	27(177) 5'-gTA 3'	9(121) 5'-TgC 3'	9(122) 5'-gTT 3'	112(432) 5'-CAA 3'	98(389) 5'-CAT 3'	15(140) 5'-gTA 3'	11(129) 5'-TTA 3'	140(514) 5'-CAA 3'	25(170) 5'-gCA 3'		140(514) 5'-CAA 3'	141(517) 5'-CTA 3'
	33(195) 5'-AAg 3'	9(122) 5'-gTT 3'	15(140) 5'-gTA 3'	118(449) 5'-CTA 3'	110(426) 5'-AAA 3'	26(173) 5'-TCT 3'	29(181) 5'-CCT 3'		48(239) 5'-CCA 3'			142(520) 5'-gCC 3'
	45(230) 5'-ggA 3'	24(166) 5'-gCT 3'	33(193) 5'-ATT 3'	130(485) 5'-CCA 3'	133(494) 5'-TCA 3'	55(259) 5'-ggT 3'	38(210) 5'-gTg 3'		57(266) 5'-TgA 3'			145(529) 5'-CCg 3'
		30(184) 5'-gAT 3'	140(514) 5'-CAA 3'	133(493) 5'-TTT 3'		59(272) 5'-CgT 3'	39(211) 5'-TgA 3'		154(558) 5'-ACT 3'			
		135(500) 5'-TgA 3'		139(512) 5'-gAT 3'			62(282) 5'-AAg 3'					
				145(530) 5'-CAT 3'								
3'-primer(s) ³				151(547) 5'-ACA 3'								
	74(317) 5'-CCg 3'	27(177) 5'-gTT 3'	55(260) 5'-gCg 3'	48(240) 5'-gCg 3'	64(286) 5'-CTA 3'	45(229) 5'-CCA 3'	77(326) 5'-CCg 3'	47(236) 5'-ggT 3'	86(353) 5'-ACC 3'	70(304) 5'-CCC 3'	47(237) 5'-CgA 3'	185(650) 5'-CgA 3'
		30(184) 5'-gTg 3'	55(260) 5'-gCg 3'	170(604) 5'-gAC 3'	71(307) 5'-ggC 3'	56(263) 5'-CAA 3'	86(353) 5'-ACg 3'	48(240) 5'-gCg 3'	185(650) 5'-CgT 3'	86(353) 5'-ACg 3'	47(237) 5'-CTA 3'	
		57(265) 5'-gCT 3'	57(266) 5'-CAA 3'		74(318) 5'-AAg 3'	57(265) 5'-gCT 3'		49(241) 5'-TgC 3'	186(653) 5'-CCg 3'		47(237) 5'-Cgg 3'	
		74(316) 5'-CgC 3'	60(274) 5'-gTT 3'		89(361) 5'-CgC 3'	61(278) 5'-TCT 3'		163(583) 5'-CAC 3'	186(654) 5'-TAA 3'		48(240) 5'-gCg 3'	
		167(596) 5'-CAT 3'	165(589) 5'-ggT 3'		185(650) 5'-Cgg 3'	64(286) 5'-CTA 3'		167(595) 5'-ACA 3'			50(244) 5'-CAT 3'	
		179(632) 5'-ggT 3'	172(611) 5'-AgC 3'			77(326) 5'-CCT 3'		170(604) 5'-gAA 3'			91(369) 5'-TTT 3'	
		179(633) 5'-AgT 3'	182(641) 5'-ggC 3'			77(327) 5'-ACg 3'		174(617) 5'-CgC 3'			163(583) 5'-CAC 3'	
			184(647) 5'-Tgg 3'					184(647) 5'-Tgg 3'			176(622) 5'-CTA 3'	
			187(656) 5'-ACA 3'								179(632) 5'-ggT 3'	
Well No.	37	38	39	40	41	42	43	44	45	46	47	48

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Well No.	49	50	51	52	53	54	55	56	57	58	59	60
Length of spec.	105	125	125	95	130	110	115	100	150	145	100	110
PCR product	140	155	185	220	160	150		140	185	180	190	150
	165	200	270		195	190		185		215	220	180
	240	230				245		220		250		
								265				
Length of int.	430	430	430	430	430	430	430	430	430	430	430	430
pos. control ¹												
5'-primer(s) ²	14(138)	20(154)	23(164)	13(134)	30(185)	26(173)	130(485)	13(134)	30(185)	14(136)	13(134)	30(184)
	5'-ATA 3'	5'-ACA 3'	5'-gCA 3'	5'-ggT 3'	5'-AAg 3'	5'-ggg 3'	5'-CCA 3'	5'-ggC 3'	5'-AAg 3'	5'-gCA 3'	5'-ggC 3'	5'-gAC 3'
	20(154)	20(154)	24(166)	23(163)	31(187)	135(500)		28(179)	140(514)	20(154)	14(136)	
	5'-ACT 3'	5'-ACC 3'	5'-gCA 3'	5'-AgT 3'	5'-ACC 3'	5'-TgA 3'		5'-gAC 3'	5'-CAA 3'	5'-ACA 3'	5'-gCC 3'	
	104(406)	27(175)	102(400)	30(185)	37(205)			154(558)		31(187)		
	5'-CAC 3'	5'-ATA 3'	5'-TCT 3'	5'-ATg 3'	5'-Agg 3'			5'-ACT 3'		5'-ACC 3'		
	129(481)	30(184)		167(596)						38(208)		
	5'-AAA 3'	5'-gAC 3'		5'-gCA 3'						5'-ACA 3'		
	133(494)	31(187)								38(210)		
	5'-TCA 3'	5'-ACg 3'								5'-gCg 3'		
	139(511)	36(204)								38(210)		
	5'-A 3'	5'-gAC 3'								5'-gCA 3'		
	146(533)	37(205)										
	5'-CCg 3'	5'-Agg 3'										
	147(535)	48(239)										
	5'-CCT 3'	5'-CCA 3'										
	151(547)	48(239)										
	5'-ACA 3'	5'-CCC 3'										
		58(268)										
		5'-ACT 3'										
3'-primer(s) ³	57(266)	86(353)	55(260)	86(353)	66(294)	51(248)	154(558)	60(274)	66(294)	71(307)	32(191)	52(251)
	5'-Cgg 3'	5'-gCT 3'	5'-gCg 3'	5'-ACg 3'	5'-ACT 3'	5'-gCT 3'	5'-AAA 3'	5'-gTA 3'	5'-ACT 3'	5'-ggC 3'	5'-TAC 3'	5'-gCA 3'
	170(604)		55(260)	185(650)	69(303)	74(316)		64(287)	81(338)	75(320)	61(279)	65(290)
	5'-gAC 3'		5'-gCg 3'	5'-CgA 3'	5'-CTg 3'	5'-Cgg 3'		5'-TCC 3'	5'-Tgg 3'	5'-CCC 3'	5'-TTT 3'	5'-CCC 3'
			127(477)		81(338)	78(328)		73(314)	185(650)	83(343)	64(287)	70(304)
			5'-TTC 3'		5'-TCg 3'	5'-CAA 3'		5'-CCA 3'	5'-Cgg 3'	5'-gTC 3'	5'-TCC 3'	5'-CCT 3'
			142(520)		86(353)	94(376)		88(359)		86(353)	71(307)	72(311)
			5'-AAT 3'		5'-ACg 3'	5'-TCA 3'		5'-TgA 3'		5'-ACg 3'	5'-ggA 3'	5'-CCT 3'
			150(544)			157(566)		174(618)			76(323)	74(317)
			5'-gTC 3'			5'-ggg 3'		5'-ACT 3'			5'-TgC 3'	5'-CCA 3'
			170(604)			168(599)		185(650)				80(335)
			5'-gAT 3'			5'-CTT 3'		5'-CgT 3'				5'-gTT 3'
			179(631)			175(619)						
			5'-gCg 3'			5'-CAT 3'						
						182(642)						
						5'-ggT 3'						
						186(652)						
						5'-CAT 3'						
Well No.	49	50	51	52	53	54	55	56	57	58	59	60



Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Well No.	61	62	63	64	65	66	67	68	69	70	71	72
Length of spec.	110	90	95	100	115	110	100	100	105	75	50	145
PCR product	150	140	180	160	220	135	155	175	140		135	180
	215	165	220	230		165			185		160	210
		225				205					235	270
						255						
Length of int. pos. control ¹	430	430	430	430	430	430	430	430	430	430	430	430
5'-primer(s) ²	30(184) 5'-gAC 3'	9(122) 5'-gTT 3'	30(184) 5'-gAT 3'	30(185) 5'-AAg 3'	30(185) 5'-AAg 3'	30(185) 5'-AAg 3'	17(147) 5'-TTA 3'	13(134) 5'-ggT 3'	9(122) 5'-gTg 3'	55(260) 5'-gCC 3'	13(133) 5'-AgA 3'	10(125) 5'-CCg 3'
		34(196) 5'-ACT 3'		146(533) 5'-CCT 3'	135(501) 5'-gAT 3'	99(391) 5'-CAA 3'	23(163) 5'-AgT 3'	36(204) 5'-gAT 3'	24(166) 5'-gCC 3'		37(207) 5'-TAA 3'	27(175) 5'-ATA 3'
		38(209) 5'-CgA 3'				112(432) 5'-CAA 3'	35(199) 5'-gAT 3'	39(212) 5'-gCT 3'	35(199) 5'-gAT 3'		44(227) 5'-CgC 3'	32(192) 5'-..T 3'
		52(251) 5'-gCT 3'				118(449) 5'-CTA 3'	41(217) 5'-TCT 3'	41(217) 5'-TCT 3'	36(204) 5'-gAT 3'		47(235) 5'-TgA 3'	37(207) 5'-TAA 3'
						141(517) 5'-CCA 3'	41(217) 5'-TCC 3'		38(209) 5'-CgC 3'		47(235) 5'-TgC 3'	38(208) 5'-TCA 3'
						145(530) 5'-CAA 3'					48(239) 5'-CCC 3'	41(219) 5'-gAg 3'
						147(535) 5'-CCg 3'					74(316) 5'-CgA 3'	48(239) 5'-TCC 3'
3'-primer(s) ³												135(501) 5'-gAT 3'
	52(251) 5'-gCA 3'	46(232) 5'-CAA 3'	48(238) 5'-CCA 3'	49(241) 5'-CgT 3'	51(248) 5'-gCC 3'	61(279) 5'-TTT 3'	57(266) 5'-Cgg 3'	57(266) 5'-Cgg 3'	57(266) 5'-Cgg 3'	66(294) 5'-ATg 3'	77(326) 5'-CCg 3'	86(353) 5'-gCT 3'
	65(290) 5'-CCC 3'	49(241) 5'-TgC 3'	76(322) 5'-gTg 3'	51(248) 5'-gCC 3'	89(362) 5'-TCT 3'	72(311) 5'-CCg 3'						167(595) 5'-ACA 3'
	70(304) 5'-CCT 3'	50(244) 5'-CAT 3'	90(365) 5'-AgT 3'	69(303) 5'-CTg 3'	163(585) 5'-gTT 3'	170(604) 5'-gAC 3'						172(610) 5'-ggC 3'
	82(341) 5'-AgC 3'	70(304) 5'-CCT 3'		94(376) 5'-TCA 3'								
	87(356) 5'-ggA 3'			185(650) 5'-CgA 3'								
	89(361) 5'-CCT 3'											
	92(370) 5'-CTA 3'											
Well No.	61	62	63	64	65	66	67	68	69	70	71	72

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Well No.	73	74	75	76	77	78	79	80	81	82	83	84
Length of spec.	95	115	115	150	95	95	90	105	120	170	165	115
PCR product	125	140	170	185	120	130	125	175	170	235	235	190
	160	185	210		190	170	200			280	260	
	225											
	270											
Length of int.	430	430	430	430	430	430	430	430	430	430	430	430
pos. control ¹												
5'-primer(s) ²	10(125) 5'-CCg ^{3'}	26(173) 5'-ggg ^{3'}	13(134) 5'-ggC ^{3'}	13(134) 5'-ggC ^{3'}	27(177) 5'-gTA ^{3'}	27(177) 5'-gTA ^{3'}	27(177) 5'-gTA ^{3'}	9(122) 5'-gTT ^{3'}	8(119) 5'-CgC ^{3'}	9(122) 5'-gTT ^{3'}	9(122) 5'-gTT ^{3'}	9(122) 5'-gTT ^{3'}
	23(163) 5'-AgA ^{3'}		15(140) 5'-gTA ^{3'}			38(210) 5'-gCA ^{3'}	100(394) 5'-TgT ^{3'}		9(122) 5'-gTT ^{3'}			154(558) 5'-ACT ^{3'}
	26(173) 5'-TTT ^{3'}		20(154) 5'-ACC ^{3'}				105(411) 5'-AgC ^{3'}		16(142) 5'-gCC ^{3'}			
	47(235) 5'-TgC ^{3'}		32(190) 5'-TCC ^{3'}						29(182) 5'-CAA ^{3'}			
	57(267) 5'-gAT ^{3'}											
	66(292) 5'-.g ^{3'}											
	70(306) 5'-AgC ^{3'}											
3'-primer(s) ³	86(353) 5'-gCT ^{3'}	49(242) 5'-CCT ^{3'}	52(250) 5'-CgC ^{3'}	47(237) 5'-CgA ^{3'}	46(233) 5'-ACC ^{3'}	48(238) 5'-gCA ^{3'}	43(224) 5'-Cgg ^{3'}	31(188) 5'-Agg ^{3'}	52(251) 5'-gCA ^{3'}	53(253) 5'-CTA ^{3'}	51(247) 5'-CgA ^{3'}	36(203) 5'-ACC ^{3'}
		54(257) 5'-gCT ^{3'}	57(266) 5'-CAA ^{3'}	53(254) 5'-CCT ^{3'}	53(254) 5'-CCA ^{3'}	53(254) 5'-CCA ^{3'}	55(260) 5'-gCA ^{3'}	52(251) 5'-gCA ^{3'}	53(253) 5'-CTT ^{3'}	74(317) 5'-ACC ^{3'}	74(317) 5'-ACC ^{3'}	59(271) 5'-CTT ^{3'}
		58(268) 5'-ggT ^{3'}	60(274) 5'-gTA ^{3'}	61(279) 5'-TTT ^{3'}	76(322) 5'-gTg ^{3'}	59(271) 5'-CTT ^{3'}	154(558) 5'-AAA ^{3'}	56(262) 5'-AgA ^{3'}	55(260) 5'-gCg ^{3'}	88(359) 5'-CgA ^{3'}	82(341) 5'-AgC ^{3'}	60(276) 5'-CAC ^{3'}
		59(271) 5'-CTA ^{3'}	70(304) 5'-CCC ^{3'}		77(326) 5'-CCA ^{3'}	71(308) 5'-ggA ^{3'}						178(629) 5'-TgC ^{3'}
		74(316) 5'-Cgg ^{3'}										
Well No.	73	74	75	76	77	78	79	80	81	82	83	84

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

Well No.	85	86	87	88	89	90	91	92	93	94	95
Length of spec.	135	95	105	110	110	130	70	105	95	170	95
PCR product	180	130	140	215	150	175	100	155	170		
	215	180			185	240	185	190	225		
	265	220			260			230			
Length of int. pos. control ¹	430	430	430	430	430	430	430	430	430	430	430
5'-primer(s) ²	10(124)	9(121)	34(196)	26(174)	14(136)	102(400)	8(119)	30(185)	14(138)	10(126)	71(309)
	5'-TCT 3'	5'-TgC 3'	5'-ACT 3'	5'-TAC 3'	5'-gCg 3'	5'-TCT 3'	5'-CgC 3'	5'-AAg 3'	5'-ATA 3'	5'-CAA 3'	5'-ACC 3'
	14(137)	10(124)	39(212)	31(187)	37(207)		47(237)		16(142)	13(134)	
	5'-CAC 3'	5'-TCT 3'	5'-gCA 3'	5'-ACT 3'	5'-TAA 3'		5'-TAC 3'		5'-gCA 3'	5'-ggg 3'	
	24(167)	24(166)	45(230)	63(285)	49(242)				123(464)		
	5'-CgA 3'	5'-gCA 3'	5'-ggC 3'	5'-Agg 3'	5'-ggA 3'				5'-CTT 3'		
	25(169)	24(167)	48(238)		63(285)				167(595)		
	5'-TgT 3'	5'-CgA 3'	5'-ACA 3'		5'-Agg 3'				5'-AgA 3'		
	39(212)	38(208)	55(260)								
	5'-gCA 3'	5'-ACT 3'	5'-gCC 3'								
110(426)	48(238)										
5'-AAA 3'	5'-ACA 3'										
3'-primer(s) ³	70(304)	38(209)	70(304)	87(356)	86(353)	131(488)	57(265)	51(248)	57(266)	55(260)	87(355)
	5'-CCC 3'	5'-gCg 3'	5'-CCC 3'	5'-ggA 3'	5'-ACg 3'	5'-ACT 3'	5'-gTC 3'	5'-gCA 3'	5'-Cgg 3'	5'-gCg 3'	5'-gAC 3'
	185(650)	70(304)				146(532)	66(294)	66(294)	185(650)	57(266)	88(358)
	5'-Cgg 3'	5'-CCC 3'				5'-ggC 3'	5'-ATg 3'	5'-ACT 3'	5'-CgA 3'	5'-Cgg 3'	5'-gCC 3'
						167(596)		72(311)			93(374)
						5'-CAA 3'		5'-CCg 3'			5'-gCA 3'
						167(596)		79(332)			
						5'-CAA 3'		5'-TgT 3'			
					172(611)		92(372)				
					5'-AgA 3'		5'-CgA 3'				
Well No.	85	86	87	88	89	90	91	92	93	94	95

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

CELL LINE VALIDATION SHEET																			
DQB1 high resolution primer set ²																			
				Well															
				1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
				201668201	201668202	201550003	201550004	201780305	201550006	201550007	201668208	201550009	201780310	201550011	201780312	201780313	201668214	201780315	201550016
	IHWC cell line ¹	DQB1		Production No.															
1	9001 SA	*05:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707	*06:01	*02:02	+	-	+	+	-	-	-	-	+	+	-	-	-	-	+	-
3	9011 E4181324	*06:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	+	-
4	9275 GU373	*02:01		+	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
5	9009 KAS011	*05:02		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM	*03:02	*06:01	-	-	+	-	-	-	-	-	+	+	-	-	-	-	+	-
7	9020 QBL	*02:01		+	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
8	9025 DEU	*03:01		-	-	+	-	-	+	+	-	-	-	-	-	-	-	-	+
9	9026 YAR	*03:02		-	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
10	9107 LKT3	*04:01		-	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
11	9051 PITOUT	*02:02		+	-	-	+	-	-	-	-	+	+	-	-	-	-	-	-
12	9052 DBB	*03:03		-	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
13	9004 JESTHOM	*05:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA	*04:02		-	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
15	9075 DKB	*03:03		-	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
16	9037 SWEIG007	*03:01		-	-	+	-	-	+	+	-	-	-	-	-	-	-	-	+
17	9282 CTM3953540	*02:01	*06:03	+	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
18	9257 32367	*06:02	*02:02	+	-	+	+	-	-	-	-	+	+	-	-	-	-	-	-
19	9038 BM16	*03:01		-	-	+	-	-	+	+	-	-	-	-	-	-	-	-	+
20	9059 SLE005	*06:04		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA	*03:01		-	-	+	-	-	+	+	-	-	-	-	-	-	-	-	+
22	9056 KOSE	*05:03	*06:04	-	-	+	-	W	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL	*05:03	*06:01	-	-	+	-	-	-	-	-	-	-	-	-	-	-	+	-
24	9035 JBUSH	*03:01		-	-	+	-	-	+	+	-	-	-	-	-	-	-	-	+
25	9049 IBW9	*02:02		+	-	-	+	-	-	-	-	+	+	-	-	-	-	-	-
26	9285 WT49	*02:01		+	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
27	9191 CH1007	*04:01	*05:01	-	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
28	9320 BEL5GB	*02:02	*03:01	+	-	+	+	-	+	+	+	+	+	-	-	-	-	-	+
29	9050 MOU	*02:02		+	-	-	+	-	-	-	-	+	+	-	-	-	-	-	-
30	9021 RSH	*04:02		-	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
31	9019 DUCAF	*02:01		+	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
32	9297 HAG	*03:01		-	-	+	-	-	+	+	-	-	-	-	-	-	-	-	+
33	9098 MT14B	*03:02		-	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
34	9104 DHIF	*03:01		-	-	+	-	-	+	+	-	-	-	-	-	-	-	-	+
35	9302 SSTO	*03:05		-	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
36	9024 KT17	*03:02		-	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
37	9065 HHKB	*06:03		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL	*03:01		-	-	+	-	-	+	+	-	-	-	-	-	-	-	-	+
39	9315 CML	*02:01	*03:01	+	-	+	-	-	+	+	+	+	+	-	-	-	-	-	+
40	9134 WHONP199	*02:02	*03:03	+	-	+	+	-	-	-	-	+	+	-	-	-	-	-	-
41	9055 H0301	*06:09		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089	*06:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	+	-
43	9076 T7526	*03:03		-	-	+	-	-	-	-	-	+	+	-	-	-	-	-	-
44	9057 TEM	*05:03		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO	*02:02		+	-	-	+	-	-	-	-	+	+	-	-	-	-	-	-
46	9013 SCHU	*06:02		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO	*03:01		-	-	+	-	-	+	+	-	-	-	-	-	-	-	-	+
48	9303 TER-ND	*05:01		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-



Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

CELL LINE VALIDATION SHEET																					
DQB1 high resolution primer set ²																					
				Production No.	Well																
					17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
					201780317	201780318	201550019	201550020	201780321	201550022	201550023	201668224	201780325	201550026	201550027	201550028	201550029	201780330	201550031	201550032	
	IHWC cell line ¹		DQB1																		
1	9001 SA		*05:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
2	9280 LK707		*06:01	*02:02	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
3	9011 E4181324		*06:01		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
4	9275 GU373		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011		*05:02		-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-
6	9353 SM		*03:02	*06:01	-	-	+	-	-	-	-	-	-	-	+	-	-	-	-	-	-
7	9020 QBL		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU		*03:01		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026 YAR		*03:02		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3		*04:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB		*03:03		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM		*05:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA		*04:02		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
15	9075 DKB		*03:03		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007		*03:01		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282 CTM3953540		*02:01	*06:03	-	-	-	-	-	-	-	-	-	+	-	-	+	+	+	-	-
18	9257 32367		*06:02	*02:02	-	-	-	-	-	-	-	-	-	+	+	+	-	-	-	-	-
19	9038 BM16		*03:01		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005		*06:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
21	9064 AMALA		*03:01		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE		*05:03	*06:04	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	+
23	9124 IHL		*05:03	*06:01	-	-	-	-	-	-	-	-	+	-	+	-	-	-	-	-	-
24	9035 JBUSH		*03:01		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007		*04:01	*05:01	-	-	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-
28	9320 BEL5GB		*02:02	*03:01	-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH		*04:02		-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-
31	9019 DUCAF		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG		*03:01		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B		*03:02		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF		*03:01		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO		*03:05		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17		*03:02		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB		*06:03		-	-	-	-	-	-	-	-	+	-	-	-	+	+	+	-	-
38	9099 LZL		*03:01		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML		*02:01	*03:01	-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134 WHONP199		*02:02	*03:03	-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055 H0301		*06:09		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089		*06:01		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
43	9076 T7526		*03:03		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM		*05:03		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
45	9239 SHJO		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU		*06:02		-	-	-	-	-	-	-	-	-	+	+	+	-	-	-	-	-
47	9045 TUBO		*03:01		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303 TER-ND		*05:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-



Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

CELL LINE VALIDATION SHEET																			
DQB1 high resolution primer set ²																			
				Well															
				33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48
				201780333	201780334	201780335	201668236	201550037	201780338	201780339	201780340	201780341	201780342	201550043	201780344	201780345	201550046	201550047	201550048
	IHWC cell line ¹	DQB1	Production No.																
1	9001 SA	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707	*06:01	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011 E4181324	*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275 GU373	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011	*05:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM	*03:02	*06:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
7	9020 QBL	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026 YAR	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
11	9051 PITOUT	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
15	9075 DKB	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282 CTM3953540	*02:01	*06:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
18	9257 32367	*06:02	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005	*06:04		+	+	-	-	-	-	-	-	W	-	-	-	-	-	-	-
21	9064 AMALA	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE	*05:03	*06:04	+	+	-	-	-	-	-	-	W	-	-	-	-	-	-	-
23	9124 IHL	*05:03	*06:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035 JBUSH	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007	*04:01	*05:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
28	9320 BEL5GB	*02:02	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
31	9019 DUCAF	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO	*03:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB	*06:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+
38	9099 LZL	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML	*02:01	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134 WHONP199	*02:02	*03:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055 H0301	*06:09		+	+	-	+	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089	*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM	*05:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303 TER-ND	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-



Product Insert

101.221-12 – including *Taq* pol., IFU-01
 101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
 “Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

CELL LINE VALIDATION SHEET																				
DQB1 high resolution primer set ²																				
				Well																
				49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
				Production No.	201780349	201780350	201780351	201550052	201668253	201780354	201668255	201780356	201668257	201780358	201668259	201668260	201780361	201780362	201780363	201668264
	IHWC cell line ¹		DQB1																	
1	9001 SA		*05:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707		*06:01	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011 E4181324		*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275 GU373		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011		*05:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM		*03:02	*06:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
7	9020 QBL		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU		*03:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
9	9026 YAR		*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB		*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA		*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075 DKB		*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007		*03:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
17	9282 CTM3953540		*02:01	*06:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257 32367		*06:02	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16		*03:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
20	9059 SLE005		*06:04		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-
21	9064 AMALA		*03:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
22	9056 KOSE		*05:03	*06:04	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-
23	9124 IHL		*05:03	*06:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035 JBUSH		*03:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
25	9049 IBW9		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007		*04:01	*05:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320 BEL5GB		*02:02	*03:01	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
29	9050 MOU		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH		*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019 DUCAF		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG		*03:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
33	9098 MT14B		*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF		*03:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
35	9302 SSTO		*03:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17		*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB		*06:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL		*03:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
39	9315 CML		*02:01	*03:01	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
40	9134 WHONP199		*02:02	*03:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055 H0301		*06:09		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-
42	9066 TAB089		*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526		*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM		*05:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU		*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO		*03:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
48	9303 TER-ND		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-



101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

CELL LINE VALIDATION SHEET																			
DQB1 high resolution primer set ²																			
				Well															
				65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80
				201668265	201780366	201780367	201550068	201780369	201550070	201780371	201780372	201780373	201780374	201780375	201780376	201780377	201780378	201780379	201780380
	IHWC cell line ¹		DQB1	Production No.															
1	9001 SA		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707		*06:01	*02:02	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
3	9011 E4181324		*06:01		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
4	9275 GU373		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011		*05:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM		*03:02	*06:01	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
7	9020 QBL		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026 YAR		*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3		*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB		*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA		*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075 DKB		*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282 CTM3953540		*02:01	*06:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257 32367		*06:02	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005		*06:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE		*05:03	*06:04	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL		*05:03	*06:01	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
24	9035 JBUSH		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007		*04:01	*05:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320 BEL5GB		*02:02	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH		*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019 DUCAF		*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B		*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO		*03:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17		*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB		*06:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML		*02:01	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134 WHONP199		*02:02	*03:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055 H0301		*06:09		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089		*06:01		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
43	9076 T7526		*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM		*05:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO		*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU		*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO		*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303 TER-ND		*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-



Product Insert

101.221-12 – including *Taq* pol., IFU-01
101.221-12u – without *Taq* pol., IFU-02

Visit <https://labproducts.caredx.com> for
“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

CELL LINE VALIDATION SHEET																					
DQB1 high resolution primer set ²																					
					Well																
					81	82	83	84	85	86	87	88	89	90	91	92	93	94	95		
					201668281	201550082	201550083	201780384	201780385	201780386	201780387	201780388	201668289	201780390	201780391	201780392	201668293	201780394	201550095		
	IHCW cell line ¹			DQB1	Production No.																
1	9001	SA	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
2	9280	LK707	*06:01	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
3	9011	E4181324	*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
4	9275	GU373	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
5	9009	KAS011	*05:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
6	9353	SM	*03:02	*06:01	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-		
7	9020	QBL	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
8	9025	DEU	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
9	9026	YAR	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-		
10	9107	LKT3	*04:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
11	9051	PITOUT	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
12	9052	DBB	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-		
13	9004	JESTHOM	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
14	9071	OLGA	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
15	9075	DKB	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-		
16	9037	SWEIG007	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
17	9282	CTM3953540	*02:01	*06:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
18	9257	32367	*06:02	*02:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
19	9038	BM16	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
20	9059	SLE005	*06:04		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
21	9064	AMALA	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
22	9056	KOSE	*05:03	*06:04	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
23	9124	IHL	*05:03	*06:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
24	9035	JBUSH	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
25	9049	IBW9	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
26	9285	WT49	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
27	9191	CH1007	*04:01	*05:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
28	9320	BEL5GB	*02:02	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
29	9050	MOU	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
30	9021	RSH	*04:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
31	9019	DUCAF	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
32	9297	HAG	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
33	9098	MT14B	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-		
34	9104	DHIF	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
35	9302	SSTO	*03:05		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-		
36	9024	KT17	*03:02		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-		
37	9065	HHKB	*06:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
38	9099	LZL	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
39	9315	CML	*02:01	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
40	9134	WHONP199	*02:02	*03:03	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-		
41	9055	H0301	*06:09		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
42	9066	TAB089	*06:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
43	9076	T7526	*03:03		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-		
44	9057	TEM	*05:03		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
45	9239	SHJO	*02:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
46	9013	SCHU	*06:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
47	9045	TUBO	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
48	9303	TER-ND	*05:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		



Product Insert

101.221-12 – including *Taq* pol., IFU-01

Visit <https://labproducts.caredx.com> for

101.221-12u – without *Taq* pol., IFU-02

“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 2, 5, 12 to 14, 17, 18, 37 to 43, 45, 48 to 51, 53, 54, 56, 58 to 75, 77 to 93 and 95 were available.

The specificities of the primers in primer solutions 2, 5, 12 to 14, 17, 18, 37 to 39, 41, 42, 45, 50, 51, 53, 58, 61, 62, 64, 66, 70, 72, 73, 75, 79, 81, 86, 87 and 91, were tested by separately adding additional 5'-primers, respectively additional 3'-primers.

In primer solution 40, 43, 48, 49, 67 to 69, 71, 85, 88, 89 and 93 it was only possible to test the 3'-primers, the 5'-primers were not possible to test.

In primer solution 54, 56, 59, 60, 63, 65, 74, 77, 78, 80, 82 to 84, 90, 92 and 95 it was only possible to test the 5'-primers, the 3'-primers were not possible to test.

In primer solutions 5, 10 to 15, 18, 21 to 23, 30, 33 to 39, 41, 42, 45, 50 to 53, 58, 62, 64, 66, 72, 73, 75, 79, 81, 86, 87, 91 and 94 one or more 5'-primers were not possible to test.

In primer solutions 1, 3, 4, 7, 8, 10, 12, 14, 15, 17, 18, 24 to 32, 35, 36, 38, 39, 41, 42, 44, 45, 47, 51, 53, 57, 58, 61, 62, 64, 66, 72, 75, 76, 79, 81 and 94 one or more 3'-primers were not possible to test.

Additional primers in primer solutions 4, 15, 25, 26, 29 to 31, 34, 35, 44, 46, 47, 52, 57 and 76 were tested by separately adding additional 5'-primers and/or 3'-primers.

Product Insert

101.221-12 – including *Taq* pol., IFU-01

Visit <https://labproducts.caredx.com> for

101.221-12u – without *Taq* pol., IFU-02

“Instructions for Use” (IFU)

Lot No.: **2F4**

Lot-specific information

ADDRESSES:

Manufacturer:

CareDx AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-508 939 00

Fax: +46-8-717 88 18

E-mail: orders-se@caredx.com

Web page: <https://labproducts.caredx.com/>

Distributed by:

CareDx GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

Tel: +43-1-710 15 00

Fax: +43-1-710 15 00 10

E-mail: orders-at@caredx.com

Web page: <https://labproducts.caredx.com/>

CareDx Lab Solutions Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-653-78171

Fax: 610-344-7989

E-mail: orders-us@caredx.com

Web page: <https://labproducts.caredx.com/>

For information on CareDx distributors worldwide, contact **CareDx GmbH**.