

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

Olerup SSP® HLA-C low resolution

Product number:	101.601-24/12 – including <i>Taq</i> pol. 101.601-24u/12u – without <i>Taq</i> pol.
Lot number:	2G8
Expiry date:	2020-12-01
Number of tests:	24 tests – Product No. 101.601-24/24u 12 tests – Product No. 101.601-12/12u
Number of wells per test:	31 + 1
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
- Adhesive PCR seals	RT
- Product Insert	RT

This Product Description is only valid for Lot No. 2G8.

Complete product documentation consists of generic Instructions for Use (IFU), lot specific Product Insert, Worksheet and Certificate.

CHANGES COMPARED TO THE PREVIOUS OLERUP SSP® HLA-C LOW RESOLUTION LOT (3F9)

The HLA-C low resolution specificity and interpretation tables have been updated for the HLA-C alleles described since the previous *Olerup SSP*® HLA-C low resolution lot was made (**Lot No. 3F9**). The kit design is based on IMGT/HLA database 3.30.0.

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
16	Added, modified	-	5'-primer added for the C*14:78 allele, 5'-primer modified for increased yield.
19	-	Exchanged, modified	3'-primer exchanged for the C*17:35 allele, 3'-primers modified for increased yield and improved HLA-specific amplification.

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

Well **32** contains Negative Control primer pairs, that will amplify more than 95% of the *Olerup* SSP[®] HLA Class I, DRB, DQB1, DPB1 and DQA1 amplicons as well as all the amplicons generated by the control primer pairs matching the human growth hormone gene.

HLA-specific PCR product sizes range from 75 to 200 base pairs.
The PCR product generated by the positive control primer pair is 430 base pairs.

Length of PCR product	105	200	105	80	75	80	85
5'-primer¹	164	340	440	45	45	43	36
	5'-CAC ^{3'}	5'-Agg ^{3'}	5'-TTA ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-TAC ^{3'}
							36
							5'-TAT ^{3'}
3'-primer²	231	2nd I	507	59	58	57	47
	5'-TgC ^{3'}	5'-AAA ^{3'}	5'-TTg ^{3'}	5'-CTC ^{3'}	5'-ggC ^{3'}	5'-CTC ^{3'}	5'-ACA ^{3'}
							48
							5'-gCA ^{3'}
							48
							5'-gCC ^{3'}
							52
							5'-TgT ^{3'}
A*	+	+	+				
B*	+	+	+				
C*	+	+	+				
DRB1				+	+		
DRB3				+	+		
DRB5				+			
DQB1					+		
DPB1						+	
DQA1							+

¹The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

²The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

PRODUCT DESCRIPTION

HLA-C-low resolution SSP subtyping

CONTENT

The primer set contains 5'- and 3'-primers for grouping the HLA-C*01:02 to C*18:10 alleles into the groups C*01:xx to C*18:xx.

PLATE LAYOUT

Each test consists of 32 PCR reactions in a 32 well cut PCR plate.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	NC

The 32 well cut PCR plate is marked with ‘HLA-C low’ in silver/gray ink.

Well No. 1 is marked with the Lot No. ‘2G8’.

Wells 1 to 31 – HLA-C low resolution primers.

Well 32 – Negative Control.

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded. The PCR plates are covered with a PCR-compatible foil.

Please note: When removing each 32 well PCR plate, make sure that the remaining plates stay sealed. Use a scalpel or a similar instrument to carefully cut the foil between the plates.

INTERPRETATION

Only HLA-C alleles will be amplified by the HLA-C low resolution typing kit, except that that a few HLA-A and HLA-B alleles will be amplified by some primer mixes. For further details see Specificity Table.

UNIQUELY IDENTIFIED ALLELES

All the HLA-C alleles, i.e. **C*01:02 to C*18:10**, recognized by the HLA Nomenclature Committee in October 2017^{1,2} will be amplified by the primers in the HLA-C low resolution SSP kit³. The HLA-C alleles will be grouped into the C*01xx to C*18xx groups.

¹HLA-C alleles listed on the IMGT/HLA web page 2017-October-27, release 3.30.0, www.ebi.ac.uk/imgt/hla.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

³The C*14:06, 14:09, 14:28:01-14:28:02 and 14:63 and the A*29:86, A*30:62, A*33:63, B*35:252 and B*39:114 alleles give rise to identical amplification patterns with the HLA-C low resolution primer set. These alleles are separated by the HLA-A low primer set and respective by the HLA-B low primer set.

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information
SPECIFICITY TABLE

HLA-C low resolution SSP typing

Specificities and sizes of the PCR products of the 31+1 primer mixes used for HLA-C low resolution SSP typing

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified HLA class I alleles ³
1¹²	95 bp, 150 bp	800 bp	*01:02:01:01-01:145N, 04:71, 07:316, 07:338, 07:579, 12:215, 14:58, 15:104, B*39:76, B*54:18
2^{4,14}	130 bp, 200 bp, 270 bp, 300 bp	800 bp	*01:10, 01:43, 02:02:01-02:02:03, 02:02:05-02:40:02, 02:42-02:57, 02:59-02:130, 04:32, 04:77, 04:265, 05:105, 06:08, 07:101, 07:148, 07:161, 07:583, 08:31, 12:119, 12:198, 14:25, 15:42, 15:111, 16:29, 16:50, 17:01:01-02-17:06, 17:08-17:16, 17:18-17:22, 17:24-17:32, 17:34-17:38, 18:03
3¹⁵	280 bp	800 bp	*02:02:01-02:02:03, 02:02:05-02:02:29, 02:02:31-02:03, 02:04 ^w , 02:05:01-02:13, 02:14:01 ^w , 02:15-02:25Q, 02:26:02-02:40:02, 02:42-02:86, 02:88-02:106, 02:108-02:130, 03:02:01-03:02:09, 03:02:11-03:02:17, 03:04:01:01-03:04:25, 03:04:27-03:10, 03:14-03:17:02, 03:19, 03:23-03:29, 03:32-03:38:02, 03:40:01-03:42, 03:44-03:48, 03:51, 03:54, 03:57:01-03:57:02, 03:60, 03:63-03:64:01, 03:65, 03:70-03:74, 03:77-03:78, 03:80:01-03:80:02, 03:82, 03:84, 03:87:01-03:87:02, 03:89-03:95, 03:98, 03:99:02-03:101, 03:104-03:111, 03:113:01-03:115, 03:117-03:118, 03:121N, 03:123, 03:125, 03:128-03:131, 03:134-03:136, 03:138-03:140, 03:143, 03:145-03:149, 03:153-03:155, 03:157, 03:159, 03:162-03:164, 03:169Q, 03:172-03:174, 03:178-03:181, 03:183-03:184:02, 03:186:01-03:186:02, 03:190-03:191, 03:193-03:194, 03:197-03:201N, 03:208N-03:213, 03:215-03:216, 03:218-03:219, 03:221-03:222, 03:224N-03:226, 03:232-03:236, 03:238-03:240, 03:244Q-03:250, 03:252, 03:255-03:261, 03:263:01-03:263:02, 03:265N-03:266, 03:269-03:270, 03:277N, 03:278 ^w , 03:279-03:283, 03:286-03:287, 03:292-03:294, 03:296-03:303, 03:305-03:306, 03:309-03:311, 03:313-03:315, 03:317-03:318N, 03:322-03:323N, 03:326, 03:328-03:335, 03:337-03:338, 03:340, 03:342-03:343, 03:347-03:350, 03:353-03:355, 03:358-03:359, 03:361-03:362, 03:365-03:366N, 04:03:01-04:03:03, 04:06, 04:16, 04:80, 04:107, 04:147, 04:160, 04:171, 04:256, 05:58:01, 05:58:03, 06:03:01, 06:132:01-06:132:02, 07:96:01-07:96:02, 07:314:02, 07:578, 12:03:23, 15:02:01:01-15:02:20, 15:02:22-15:02:24, 15:02:26-15:05:08, 15:05:10-15:09, 15:10:02-15:11, 15:13:01:01-15:13:01:02, 15:15-15:19, 15:21-15:22, 15:24-15:35, 15:37-15:60, 15:62, 15:64-15:106, 15:108-15:137, 15:139-15:143, 15:145N-15:146, 16:34, 16:70, B*44:02:44, B*56:01:09, B*58:74
4^{5,12}	85 bp, 170 bp, 275 bp	1070 bp	*01:67, 01:122, 01:131, 03:02:01-03:02:12, 03:02:14-03:04:06, 03:04:08-03:11:02, 03:13:01-03:15, 03:17:01-03:40:04, 03:42-03:57:02, 03:59-03:79, 03:80:02-03:85, 03:87:01-03:93, 03:95-03:98, 03:100-03:109, 03:111-03:112, 03:114-03:150, 03:152-03:164, 03:166-03:250, 03:252-03:259, 03:261-03:266, 03:268-03:285, 03:287-03:292, 03:294-03:331, 03:333-03:367, 04:212, 07:242
5	280 bp	800 bp	*03:03:01:01-03:03:20, 03:03:22-03:03:37, 03:11:01-03:11:02, 03:13:01-03:13:02, 03:18:02, 03:20N-03:22Q, 03:30-03:31,

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

			03:43-01-03:43:02, 03:49-03:50, 03:52-03:53, 03:55-03:56, 03:58-03:59, 03:61-03:62, 03:66, 03:67 ^w , 03:68-03:69, 03:75-03:76, 03:79, 03:81, 03:83, 03:85-03:86, 03:88, 03:96-03:97, 03:102-03:103, 03:112, 03:116:01-03:116:02, 03:119:01-03:120, 03:122, 03:124, 03:126-03:127, 03:132-03:133, 03:141-03:142, 03:150-03:152, 03:158, 03:160-03:161, 03:165, 03:167-03:168, 03:171, 03:175-03:177, 03:185, 03:187-03:189N, 03:192, 03:196, 03:202-03:207, 03:214, 03:217, 03:220, 03:223, 03:227-03:230, 03:237, 03:241-03:243, 03:251, 03:253-03:254, 03:262, 03:267-03:268, 03:271-03:276, 03:284-03:285, 03:288-03:291, 03:295 ^w , 03:304, 03:307-03:308, 03:312, 03:316N, 03:319-03:321, 03:324-03:325, 03:327, 03:336, 03:339, 03:341, 03:345-03:346, 03:351-03:352, 03:356-03:357, 03:360, 03:363N-03:364, 03:367, 15:12, 15:144
6¹⁵	135 bp, 330 bp	800 bp	*02:02:01-02:02:03, 02:02:05-02:02:12, 02:02:14-02:02:25, 02:02:27-02:02:28, 02:02:30-02:20, 02:22-02:25Q, 02:27:01-02:38N, 02:40:01-02:40:02, 02:42-02:44, 02:46-02:86, 02:88-02:100, 02:101 ^w , 02:102-02:106, 02:108-02:130, 03:231, 04:01:01:01-04:01:27, 04:01:29-04:01:60, 04:01:62, 04:01:64, 04:01:66-04:01:83, 04:03:01-04:15:03, 04:17-04:20, 04:23-04:41, 04:43-04:102, 04:104-04:139, 04:141-04:195, 04:197-04:212, 04:214-04:219, 04:221-04:235, 04:237-04:276, 04:278-04:279N, 05:26, 05:78:01-05:78:02, 07:02:09, 07:125, 07:356, 07:531, 07:583, 15:11, 15:36, 15:118, 16:34
7⁶	165 bp, 390 bp, 445 bp	1070 bp	*02:94, 04:129, 05:01:01:01-05:01:36, 05:03-05:154N, 06:05, 06:67, 08:10, 12:21, 12:33, 15:107, 16:40, 16:53, 16:110, 16:113, 17:05
8¹²	125 bp, 355 bp	800 bp	*01:90, 01:101-01:102, 01:113, 01:136, 01:145N, 02:06:01-02:06:02, 02:23, 02:36:01-02:36:02, 02:47, 04:108 ^w , 04:178, 05:132, 06:02:01:01-06:02:01:09, 06:02:03-06:02:38, 06:02:40-06:16N, 06:18-06:31, 06:33-06:76:01, 06:77-06:93, 06:96-06:98, 06:100-06:132:01, 06:133-06:179, 06:181-06:205, 08:113, 12:02:11, 12:03:09, 12:03:26, 12:15, 12:113, 12:208, 15:02:01:01-15:03, 15:05:01 ^w -15:05:12 ^w , 15:06:01-15:13:01:02, 15:15-15:19, 15:21, 15:22 ^w -15:23 ^w , 15:26-15:28, 15:29 ^w , 15:31-15:35, 15:36 ^w , 15:37-15:39, 15:41-15:45, 15:46 ^w , 15:47-15:53, 15:54 ^w , 15:55-15:58, 15:59 ^w , 15:60-15:64, 15:67-15:68, 15:69 ^w -15:70 ^w , 15:71-15:75, 15:78:01-15:89, 15:90 ^w , 15:91-15:101, 15:103, 15:104 ^w -15:105Q ^w , 15:106-15:107, 15:108 ^w , 15:109, 15:110 ^w -15:111 ^w , 15:112-15:114, 15:115N ^w -15:117 ^w , 15:118-15:124, 15:125 ^w , 15:126-15:139, 15:140 ^w , 15:141-15:146, 16:01:16, 16:20, 16:109, B*58:02:01
9⁷	245 bp, 425 bp	800 bp	*07:01:01:01-07:33N, 07:35-07:294, 07:296-07:600N, 12:181
10^{8,12,15}	110 bp, 165 bp, 390 bp	800 bp	*01:43, 02:87, 03:280, 07:101, 07:148, 07:161, 07:583, 08:01:01:01-08:63, 08:65-08:154, 12:127, 12:203
11	155 bp, 255 bp, 340 bp	1070 bp	*01:14, 01:59, 01:118, 02:02:01-02:02:03, 02:02:05-02:40:02, 02:42-02:86, 02:88-02:130, 03:07:01-03:07:02, 03:15, 03:45, 03:130, 03:140, 03:163, 03:231, 03:243, 03:248, 03:268, 03:297, 03:308, 04:01:01:01-04:01:83, 04:03:01-04:20, 04:23-04:276, 04:278-04:279N, 05:01:01:01-05:01:36, 05:03-05:128N, 05:130-05:154N, 06:02:01:01-06:02:01:09, 06:02:03-06:02:11, 06:02:13-06:10, 06:12-06:51, 06:53:01-06:121, 06:123, 06:125-06:146, 06:148-06:205, 07:07, 07:09, 07:49, 07:53, 07:76:01-07:76:02, 07:125, 07:210, 07:216, 07:238, 07:247, 07:315, 07:328, 07:356, 07:403, 07:406, 07:531, 07:598, 08:10, 12:04:01-12:05, 12:09, 12:21, 12:33, 12:41, 12:54, 12:60, 12:146, 12:188, 12:194, 14:04, 14:12, 14:49, 14:64,

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

			14:77, 15:02:01:01-15:05:10, 15:06:01-15:06:03, 15:08-15:13:01:02, 15:15-15:19, 15:22-15:24, 15:26-15:42, 15:44:01-15:70, 15:72-15:115N, 15:117-15:143, 15:145N-15:146, 16:01:01:01-16:02:13, 16:06-16:09, 16:11-16:28, 16:30N-16:32, 16:34, 16:36-16:39:02, 16:41, 16:43-16:44, 16:46-16:52, 16:54, 16:56-16:60, 16:62-16:65, 16:69-16:77N, 16:79-16:81, 16:83-16:108, 16:111-16:114, 17:01:01:02-17:21, 17:23-17:38, 18:01-18:10, A*24:174
12^{9,12,13}	100 bp, 155 bp, 220 bp	800 bp	*01:02:18, 01:04, 01:21, 06:02:38, 06:76:02, 12:02:01-12:02:13, 12:03:01:01-12:03:01:09, 12:03:03-12:03:07, 12:03:09-12:03:15, 12:03:17-12:03:19, 12:03:21-12:03:29, 12:03:31-12:03:33, 12:03:35-12:03:42, 12:04:02-12:08, 12:10:01-12:13, 12:14:02-12:18:01, 12:19-12:25, 12:27-12:32, 12:34-12:58, 12:60-12:143, 12:145-12:180, 12:182-12:184, 12:186-12:207, 12:209-12:216, 14:02:08, 16:01:17, 16:02:13, 16:15:02, B*57:01:24
13^{12,13}	120 bp, 250 bp	800 bp	*01:21, 02:12, 02:49, 02:55:01-02:55:02, 02:115, 04:01:01:01-04:01:83, 04:03:01-04:09N, 04:12-04:20, 04:23-04:35, 04:37-04:54, 04:56-04:152, 04:154-04:168, 04:170N-04:191N, 04:193, 04:195-04:209, 04:211-04:213, 04:216-04:226, 04:228-04:248, 04:250-04:276, 04:279N ^w , 05:42, 05:46, 05:112, 06:76:02, 07:02:09, 07:125, 07:356, 07:531, 07:583, 08:05, 08:21, 08:25, 08:137, 12:02:01-12:02:10, 12:02:12-12:03:03, 12:03:05-12:03:08, 12:03:10-12:03:12, 12:03:13 ^w , 12:03:14-12:03:23, 12:03:24 ^w , 12:03:25-12:03:33, 12:03:35-12:03:42, 12:04:02, 12:06-12:08, 12:10:01-12:20, 12:22-12:32, 12:34-12:48, 12:50-12:97, 12:99:01-12:122, 12:124-12:187, 12:189-12:194, 12:195:02-12:216, 15:03, 15:16, 15:25, 15:75, 16:01:01:01-16:02:13, 16:06-16:28, 16:30N-16:32, 16:34, 16:36-16:39:02, 16:41, 16:43-16:47, 16:49-16:52, 16:54, 16:56-16:60, 16:62-16:65, 16:67, 16:69-16:77N, 16:79-16:81, 16:83-16:87, 16:89N-16:90, 16:92-16:102, 16:104-16:108, 16:111-16:114, 17:01:04, B*67:02
14¹⁶	160 bp, 215 bp	800 bp	*01:04, 01:09, 02:05:01-02:05:03, 02:17, 04:42:01-04:42:02, 04:220, 06:02:01:01-06:02:01:09, 06:02:03-06:02:15, 06:02:17-06:02:42, 06:02:44-06:03:02, 06:07-06:13, 06:15-06:34:02, 06:36-06:39, 06:41-06:71, 06:73-06:78, 06:80, 06:82-06:100, 06:102:01-06:122, 06:124-06:126, 06:128N-06:135, 06:137-06:142, 06:145-06:152N, 06:154-06:196, 06:198-06:202, 06:205, 07:31:01-07:31:02, 07:125, 07:177, 07:356, 07:514, 07:531, 12:03:01:01-12:07, 12:11-12:13, 12:15, 12:19, 12:23, 12:25-12:26, 12:28-12:29, 12:31-12:35, 12:37-12:39N, 12:42Q-12:43, 12:45-12:48, 12:50-12:55, 12:57:01-12:63, 12:65-12:66, 12:70-12:71, 12:75-12:79, 12:81-12:82, 12:87-12:95, 12:97-12:102, 12:107-12:111, 12:113, 12:115-12:116, 12:119-12:122, 12:125, 12:129, 12:131, 12:133, 12:135, 12:138-12:141, 12:143-12:144, 12:147, 12:149-12:150, 12:152, 12:154, 12:156-12:160, 12:163, 12:165, 12:167, 12:170-12:174, 12:176, 12:178, 12:180, 12:182, 12:184-12:187, 12:189-12:192, 12:194-12:195:02, 12:197, 12:199, 12:201-12:203, 12:205-12:206, 12:209-12:211, 12:213, 12:215-12:216, 14:15-14:16, 14:87, 16:04:01:01-16:04:01:02, 16:04:03-16:04:05, 16:29, 16:33, 16:35, 16:42, 16:48, 16:55, 16:61, 16:66, 16:78, 16:82, 16:91, 16:109
15¹⁴	130 bp, 255 bp	800 bp	*02:02:01 ^w , 02:02:02:01-02:02:03, 02:02:05-02:13, 02:15-02:26:03, 02:28-02:40:02, 02:42-02:64, 02:66-02:86, 02:88-02:106, 02:108-02:114, 02:116-02:130, 03:07:01-03:07:02, 03:10, 03:15, 03:29, 03:45, 03:163, 03:268, 03:297, 03:302, 04:03:01-04:03:03, 04:06, 04:16, 04:37, 04:80, 04:103, 04:107, 04:147, 04:160, 04:171, 04:190, 04:256, 05:01:01:01-05:01:36, 05:03-05:19, 05:21-05:42, 05:44:01-05:77, 05:79-05:112, 05:114-05:154N, 06:03:01-06:03:02, 06:76:01-06:76:02, 06:132:01-06:132:02, 08:10, 12:04:01-12:05, 12:09, 12:21,

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

			12:33, 12:41, 12:60, 12:72, 12:135, 12:146, 12:154, 14:45, 15:02:01:01-15:06:03, 15:08-15:13:01:02, 15:15-15:19, 15:21-15:22, 15:24, 15:26-15:35, 15:37-15:42, 15:44:01-15:62, 15:64-15:84Q, 15:86-15:137, 15:139-15:143, 15:145N-15:146, 16:02:01-16:02:13, 16:09, 16:12, 16:19, 16:25, 16:37, 16:46-16:48, 16:60, 16:63, 16:69-16:70, 16:74, 16:77N, 16:84, 16:88-16:91, 16:99, 16:101-16:104, 16:107-16:108, 17:01:01:02-17:19, 17:21, 17:23-17:38, 18:09, B*15:277, B*35:222, B*42:22
16¹⁴	160 bp, 260 bp, 445 bp	1070 bp	*01:67, 01:122, 03:155, 03:231, 04:11, 04:29, 04:36, 04:55, 04:114, 04:172, 07:64, 07:402, 08:95, 12:55, 14:02:01:01-14:11, 14:13-14:48, 14:50-14:88, 16:40, 16:53, 16:110, 16:113, A*29:86, A*30:62, A*33:63, B*35:252, B*39:114
17^{10,12}	110 bp, 325 bp	800 bp	*01:90, 01:101-01:102, 01:113, 01:116, 02:06:01-02:06:02, 02:23, 02:36:01-02:36:02, 02:68, 03:81, 03:175, 03:199, 03:245, 03:317, 04:108, 04:112, 04:169, 04:178, 05:36, 05:132, 06:89, 07:123, 07:173, 07:294, 08:113, 12:08, 12:15, 12:81, 12:113, 12:188, 12:208, 15:02:01:01-15:13:01:02, 15:15-15:19, 15:21-15:24, 15:26-15:146, 16:20, 16:64, 16:70, 16:87, 16:109
18¹⁴	180 bp, 210 bp, 240 bp	1070 bp	*02:13, 02:18, 02:33, 02:49, 02:75, 02:115, 04:01:01:01-04:01:22, 04:01:24-04:01:73, 04:01:74 ^w , 04:01:75-04:01:83, 04:03:01-04:10, 04:12-04:20, 04:23-04:32, 04:34-04:106, 04:108-04:115N, 04:117-04:129, 04:131-04:168, 04:170N-04:171, 04:173N-04:230, 04:232-04:276, 04:278-04:279N, 05:17, 05:25, 05:42, 05:55, 05:68, 05:76, 05:79, 05:115, 05:134, 05:143, 06:05, 06:31, 06:76:02, 06:118, 07:02:09, 07:31:01-07:31:02, 07:154, 07:177, 07:514, 07:583, 08:01:01:01-08:01:17, 08:01:19-08:01:21, 08:03:01-08:03:04, 08:06, 08:08:01-08:11, 08:14, 08:16:01-08:16:02, 08:20-08:22, 08:24, 08:26N, 08:28, 08:36N, 08:38, 08:40-08:42, 08:44, 08:46-08:47, 08:50, 08:56-08:61, 08:65, 08:72:01-08:72:02, 08:78-08:89N, 08:91, 08:95-08:99, 08:101-08:102, 08:105-08:106, 08:109, 08:117-08:119, 08:121N-08:122, 08:124, 08:127N-08:131, 08:133, 08:135-08:137, 08:141Q, 08:143-08:145, 08:147-08:149, 08:153-08:154, 12:14:01-12:14:02, 12:28, 12:58, 12:132, 12:135, 12:146, 12:169, 12:176, 12:181, 14:10, 14:15, 14:17, 14:87, 15:12, 15:25, 15:31, 15:62, 15:65, 15:75, 15:144, 16:01:01:01-16:02:13, 16:04:01:01-16:04:01:02, 16:04:03-16:04:05, 16:06-16:39:02, 16:41-16:42, 16:44-16:52, 16:54-16:109, 16:111-16:112, 16:114
19¹³	225 bp, 250 bp, 290 bp	800 bp	*01:60, 02:40:01-02:40:02, 04:58, 04:160, 05:23, 05:62, 05:134, 05:143, 05:146, 05:151, 06:118, 08:07, 08:47, 08:104, 08:142, 12:14:01-12:14:02, 12:92, 12:176, 14:17, 15:25, 15:65, 17:01:01:02-17:38
20¹⁴	215 bp, 420 bp	800 bp	*01:02:01:01-01:20, 01:22-01:145N, 03:58, 03:86, 03:94, 03:99:01-03:99:02, 04:37, 04:230, 04:263, 05:16, 05:85, 05:107, 06:05-06:06, 06:17, 06:31, 06:204, 07:07, 07:09, 07:39-07:40, 07:49, 07:76:01-07:76:02, 07:177, 07:210, 07:238, 07:247, 07:315, 07:328, 07:335, 07:403, 07:406, 07:563, 07:598, 08:12, 12:09, 12:24, 12:176, 14:02:01:01-14:05, 14:07N, 14:10-14:14, 14:17-14:27, 14:29-14:52, 14:55-14:62, 14:64-14:76, 14:78-14:86, 14:88, 15:102, 16:04:01:01-16:04:01:02, 16:04:03-16:04:05, 16:29, 16:33, 16:42, 16:53, 16:55, 16:61, 16:66, 16:68, 16:78, 16:82, 16:88, 16:91, 16:109, 18:01-18:10, B*07:239, B*14:03
21	325 bp, 380 bp	1070 bp	*01:03, 01:24, 01:34, 01:49:01-01:50, 01:55, 01:78, 01:131, 02:22, 02:62, 02:118, 03:03:01:01-03:04:56, 03:06:01-03:11:02, 03:14, 03:18:01-03:24, 03:26, 03:28-03:32, 03:34, 03:37:01-03:59,

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

			03:61-03:70, 03:72-03:83, 03:85, 03:87:01-03:88, 03:90-03:93, 03:96, 03:98, 03:100-03:107, 03:109, 03:111-03:120, 03:122-03:131, 03:133-03:134, 03:136-03:138, 03:140-03:145, 03:147-03:153, 03:155-03:166, 03:168, 03:170-03:176, 03:179-03:189N, 03:191-03:193, 03:195-03:196, 03:202-03:215, 03:217-03:220, 03:223, 03:227-03:237, 03:239-03:244Q, 03:246-03:247, 03:250-03:257, 03:259-03:261, 03:262 [?] , 03:263:01-03:263:02, 03:265N-03:270, 03:272-03:291, 03:293-03:295, 03:297, 03:303-03:314, 03:316N-03:325, 03:327-03:328, 03:331, 03:333-03:334, 03:336-03:337, 03:339-03:348, 03:351-03:364, 03:366N-03:367, 04:01:01:01-04:01:15, 04:01:17-04:01:83, 04:03:01-04:20, 04:24-04:53, 04:55-04:71, 04:73-04:107, 04:109-04:276, 04:278-04:279N, 05:01:01:01-05:01:36, 05:03, 05:05:01-05:21, 05:23-05:102, 05:104-05:106:02, 05:108-05:122, 05:124-05:131, 05:133-05:134, 05:136-05:146, 05:148-05:154N, 06:09, 06:14, 06:35, 06:72, 06:143-06:144, 07:10, 07:28, 07:41, 07:43:01-07:43:02, 07:184, 07:196, 07:367, 07:536, 07:568, 08:01:01:01-08:08:02, 08:10, 08:12-08:47, 08:49-08:63, 08:65-08:82, 08:84-08:107, 08:109-08:112, 08:114-08:151, 08:153-08:154, 12:31, 12:44, 12:144, 14:54, 15:02:01:01-15:02:01:07, 15:02:02-15:13:01:02, 15:15-15:19, 15:21, 15:23-15:36, 15:38-15:54, 15:56-15:57, 15:59-15:64, 15:66-15:71, 15:73-15:101, 15:103-15:124, 15:126-15:146, 16:45, 17:01:01:02-17:38, 18:01-18:10, B*73:01^w
22¹⁷	135 bp	800 bp	*02:02:13, 02:02:29, 03:02:01-03:02:09, 03:02:11-03:03:20, 03:03:22-03:04:24, 03:04:27-03:11:02, 03:13:01-03:17:02, 03:18:02-03:38:02, 03:40:01-03:64:01, 03:65-03:66, 03:67 ^w , 03:68-03:98, 03:99:02-03:136, 03:138-03:143, 03:146-03:155, 03:157-03:165, 03:167-03:169Q, 03:171, 03:173-03:181, 03:183-03:194, 03:196-03:230, 03:232-03:247, 03:249-03:263:02, 03:265N-03:277N, 03:278 ^w , 03:279-03:294, 03:295 ^w , 03:296-03:323N, 03:325-03:341, 03:343, 03:345-03:367, 04:32, 04:77, 04:265, 06:03:01, 06:132:01-06:132:02, 07:96:01-07:96:02, 07:314:02, 07:578, 14:25, 15:02:10, 15:02:17, 15:43, 18:03, B*13:13:02, B*40:01:49
23	155 bp, 235 bp	1070 bp	*04:42:01-04:42:02, 04:220, 06:02:01:01-06:02:01:09, 06:02:03-06:02:50, 06:04:01-06:75, 06:78-06:131, 06:133-06:205, 07:01:01:01-07:01:40, 07:01:41 ^w , 07:01:42-07:02:07, 07:02:09-07:02:28, 07:02:30-07:02:40, 07:02:41 ^w , 07:02:42-07:25, 07:27:01-07:32N, 07:35-07:38:02, 07:41-07:63, 07:65-07:91, 07:93-07:95, 07:97-07:138, 07:140-07:151, 07:153-07:155, 07:157-07:176, 07:178-07:209, 07:211-07:237, 07:239-07:245, 07:247-07:266, 07:268-07:294, 07:297-07:313, 07:315-07:316, 07:318-07:321, 07:323-07:334, 07:336-07:350N, 07:352-07:401, 07:403-07:404, 07:406-07:414, 07:415 ^w -07:416 ^w , 07:417-07:448, 07:450-07:492, 07:494Q-07:516, 07:518-07:577, 07:579-07:582Q, 07:584-07:600N, 12:16, 12:147, 12:195:01-12:195:02, 17:20, 18:01-18:08, 18:10, B*08:145
24¹²	95 bp, 475 bp, 590 bp	800 bp	*02:02:01-02:02:03, 02:02:06-02:02:12, 02:02:15-02:02:33, 02:04, 02:05:02-02:05:03, 02:06:02-02:09, 02:11, 02:13, 02:15, 02:19-02:27:01, 02:28-02:31, 02:34-02:36:01, 02:37-02:40:01, 02:42-02:43:01, 02:44-02:50, 02:52N-02:55:01, 02:56-02:57, 02:59-02:71, 02:74-02:75, 02:77-02:88, 02:90-02:93, 02:96-02:108, 02:111-02:123, 02:126-02:130, 03:287, 05:16, 05:106:01-05:106:02, 06:05, 06:67, 08:01:19, 08:02:02, 08:12, 08:126, 12:03:17, 12:21, 12:33, 12:203
25^{11,14}	260 bp, 360 bp, 425 bp	1070 bp	*03:231, 03:248, 03:308, 04:01:01:01-04:01:83, 04:03:01-04:20, 04:23-04:212, 04:214-04:225N, 04:227-04:276, 04:278-04:279N, 05:78:01-05:78:02, 05:141, 06:120, 12:132, 12:146, 12:194, 16:34, 16:62

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

26	155 bp, 225 bp	1070 bp	*01:02:01:01-01:03, 01:06-01:07:01, 01:08, 01:10-01:20, 01:23-01:34, 01:37N-01:48, 01:51-01:54, 01:56N-01:78, 01:80-01:119, 01:121Q-01:130, 01:132-01:145N, 02:12, 02:49, 02:55:01-02:55:02, 02:115, 03:15, 03:27, 03:38:01-03:38:02, 03:58, 03:69, 03:86, 03:94, 03:99:01-03:99:02, 03:130, 03:136, 03:163, 03:246, 03:274, 03:297, 04:03:01, 04:03:03, 04:06, 04:16, 04:37, 04:80, 04:103, 04:107, 04:129, 04:147, 04:160, 04:171, 04:190, 04:230, 04:256, 04:263, 05:16, 05:42, 05:46, 05:85, 05:107, 06:03:01-06:03:02, 06:05-06:06, 06:67, 06:76:01-06:76:02, 06:132:01-06:132:02, 06:204, 07:26:01-07:26:02, 07:92, 07:96:01-07:96:02, 07:101, 07:148, 07:161, 07:314:01-07:314:03, 07:317, 07:351, 07:578, 07:583, 08:05, 08:12, 08:21, 08:25, 08:137, 12:02:01-12:03:12, 12:03:13 ^w , 12:03:14-12:04:02, 12:06-12:15, 12:17-12:20, 12:22-12:32, 12:34-12:48, 12:50-12:54, 12:56-12:97, 12:99:01-12:146, 12:148N-12:187, 12:189-12:194, 12:196-12:216, 14:02:01:01-14:05, 14:07N, 14:10-14:14, 14:17-14:27, 14:29-14:52, 14:55-14:62, 14:64-14:76, 14:78-14:86, 14:88, 15:03, 15:16, 15:25, 15:102, 16:15:01-16:15:02, 16:25, 16:53, 16:64, 16:68, 16:88, 17:01:01:02-17:19, 17:21-17:26, 17:27N ^w , 17:28-17:38, 18:09, B*07:13, B*07:15, B*07:160, B*07:239, B*14:03, B*67:02
27	140 bp, 245 bp	1070 bp	*01:02:34, 01:21, 02:12 ^w , 02:27:01-02:27:02, 02:87, 02:115, 02:126 ^w , 03:04:25, 04:11, 04:29, 04:36, 04:55, 04:172, 04:214 ^w , 07:02:09, 07:449, 07:583, 08:01:01:01-08:02:10, 08:02:12-08:09, 08:11-08:63, 08:65-08:94, 08:95 ^w , 08:96-08:154, 12:02:01-12:02:10, 12:02:12-12:03:03, 12:03:05-12:03:08, 12:03:10-12:03:23, 12:03:24 ^w , 12:03:25-12:03:33, 12:03:35-12:03:42, 12:06-12:08, 12:10:01-12:20, 12:22-12:26, 12:28-12:32, 12:34-12:40, 12:42Q-12:53, 12:55-12:59, 12:61-12:71, 12:72 ^w , 12:73-12:122, 12:124-12:134, 12:135 ^w , 12:136-12:145, 12:147-12:153, 12:154 ^w , 12:155Q-12:187, 12:189-12:194, 12:195:02-12:216, 14:02:03, 14:03, 14:08, 14:10, 14:22, 14:35N, 14:38, 14:41, 14:53-14:54, 14:61, 14:70, 14:78-14:79, 14:86, 15:07, 15:21 ^w , 15:25, 15:116 ^w , 15:144, 16:01:01:01-16:01:01:04, 16:01:03-16:01:15, 16:01:17-16:01:22, 16:04:01:01-16:04:01:02, 16:04:03-16:04:05, 16:06-16:08, 16:10-16:11, 16:13-16:18, 16:20-16:24, 16:26-16:36, 16:37 ^w , 16:38-16:45, 16:49-16:56, 16:58-16:59, 16:61-16:62, 16:64-16:68, 16:71 ^w , 16:72-16:73, 16:75-16:76, 16:78-16:83, 16:86-16:87, 16:92-16:93, 16:95-16:98, 16:100, 16:105, 16:109-16:114, B*08:01:36, B*35:08:02, B*35:08:05, B*67:02
28¹⁴	215 bp, 555 bp	1070 bp	*01:97, 01:114, 03:58, 03:86, 03:94, 03:99:01-03:99:02, 03:209, 03:302, 04:09N, 04:263, 05:85, 05:107, 12:178, 14:02:01:01-14:04, 14:07N, 14:10-14:14, 14:17-14:27, 14:29-14:44, 14:46:01-14:52, 14:55-14:62, 14:64-14:76, 14:78-14:86, 14:88, 15:37, 15:102, 16:18, B*54:38
29	175 bp	800 bp	*03:08, 03:29, 03:31, 03:246, 07:20, 07:96:01-07:96:02, 07:263, 15:02:01:01-15:09, 15:12-15:13:01:02, 15:15, 15:18-15:19, 15:21-15:24, 15:26, 15:28-15:42, 15:44:01-15:62, 15:64-15:106, 15:108-15:112, 15:114-15:123, 15:125-15:146, 16:70
30^{13,14}	160 bp, 370 bp, 540 bp	1070 bp	*01:21, 02:12, 02:49, 02:55:01-02:55:02, 02:115, 04:01:55, 04:14-04:15:01, 04:15:03, 04:68, 05:112, 06:02:37, 06:127:02, 07:01:07, 07:04:03, 07:14, 07:27:02, 07:31:02, 07:43:02, 07:53, 07:199:02, 07:216, 07:526, 12:02:01-12:02:02:03, 12:02:04-12:02:07, 12:02:09-12:02:11, 12:02:13-12:02:15, 12:03:29, 12:03:32, 12:08, 12:14:01-12:14:02, 12:16-12:18:02, 12:22, 12:27, 12:30, 12:36, 12:40-12:41, 12:56, 12:64, 12:67-12:69, 12:72-12:74, 12:80N, 12:83-12:86, 12:103-12:106, 12:112, 12:114, 12:117, 12:123-12:124, 12:126-12:128, 12:130, 12:132, 12:134, 12:136-12:137, 12:142, 12:145-12:146,

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

			12:148N, 12:151, 12:161, 12:164, 12:166, 12:168-12:169, 12:177, 12:179, 12:183, 12:193, 12:196, 12:198, 12:200, 12:204, 12:207, 12:212, 12:214, 15:75, 16:01:01:01-16:02:13, 16:06-16:28, 16:30N-16:32, 16:34, 16:36-16:39:02, 16:41, 16:43-16:47, 16:49-16:52, 16:54, 16:56-16:60, 16:62-16:65, 16:69-16:77N, 16:79-16:81, 16:83-16:87, 16:89N-16:90, 16:92-16:102, 16:104-16:108, 16:111-16:114, 17:01:10, A*24:96, A*24:106, A*24:174
31	520 bp	1070 bp	*01:14, 02:02:01-02:02:03, 02:02:05-02:02:33, 02:04-02:15, 02:17, 02:19-02:26:03, 02:28-02:40:02, 02:42-02:64, 02:66-02:71, 02:73-02:86, 02:88-02:114, 02:116-02:130, 04:01:01:01-04:01:76, 04:01:78-04:01:83, 04:03:01-04:03:03, 04:05, 04:07-04:10, 04:12, 04:15:01-04:20, 04:23-04:28, 04:30-04:33, 04:35, 04:37-04:54, 04:56-04:57, 04:59Q-04:67, 04:69-04:110, 04:112-04:113, 04:115N-04:121, 04:123N-04:159, 04:161-04:171, 04:173N-04:177, 04:179-04:187, 04:189-04:211, 04:213-04:248, 04:250-04:264, 04:266-04:276, 04:278-04:279N, 05:01:01:01-05:01:36, 05:03-05:10, 05:12-05:16, 05:18:01-05:19, 05:21-05:26, 05:28-05:31, 05:33-05:50, 05:52-05:61, 05:63-05:67, 05:69, 05:71-05:78:02, 05:80-05:111, 05:113N-05:114, 05:116-05:128N, 05:130-05:133, 05:135-05:145, 05:147-05:150, 05:152-05:154N, 06:02:01:01-06:02:01:09, 06:02:03-06:03:02, 06:05-06:10, 06:12-06:39, 06:41-06:51, 06:53:01-06:78, 06:80-06:81, 06:83-06:117, 06:119-06:146, 06:148-06:152N, 06:154-06:159, 06:161-06:182, 06:185-06:196, 06:198-06:205, 12:04:01-12:05, 12:09, 12:21, 12:33, 12:41, 12:60, 12:72, 12:135, 12:146, 12:154, 14:12, 14:49, 15:08, 15:74, 15:102, 15:143, 16:88, 16:91, 17:17, 18:01-18:10
32¹⁸	-	-	Negative Control

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-C low resolution SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherit feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³For several HLA Class I alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁴Primer mix 2 will for most C*02 alleles give rise to two specific PCR fragments.

⁵Primer mix 4 will for most C*03 alleles give rise to two specific PCR fragments.

⁶Primer mix 7 will for most C*05 alleles give rise to two specific PCR fragments.

⁷Primer mix 9 will for most C*07 alleles give rise to two specific PCR fragments.

⁸Primer mix 10 will for most C*08 alleles give rise to multiple specific PCR fragments.

⁹Primer mix 12 will for most C*12 alleles give rise to two specific PCR fragments.

¹⁰Primer mix 17 will for most C*15 alleles give rise to two specific PCR fragments.

¹¹Primer mix 25 will for most C*04 alleles give rise to two specific PCR fragments.

¹²HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

¹³Primer mixes 12, 13, 19 and 30 may have tendencies of unspecific amplifications.

¹⁴Primer mixes 2, 15, 16, 18, 20, 25, 28 and 30 have a tendency to giving rise to primer oligomer formation.

¹⁵Primer mixes 3, 6, 10 may give rise to a lower yield of HLA-specific PCR product than the other HLA-C low primer mixes.

¹⁶Primer mix 14 might faintly amplify most C*01 and the C*14 alleles.

¹⁷Primer mix 22 might generate a false band of about 500 base pairs. This band should be disregarded when interpreting HLA-C low resolution typings.

¹⁸Primer mix 32 contains a negative control, which will amplify more than 95% of HLA amplicons as well as the amplicons generated by the control primer pairs matching the human growth hormone gene. HLA-specific PCR product sizes range from 75 to 200 base pairs and the PCR product generated by the HGH positive control primer pair is 430 base pairs.

'w', might be weakly amplified.

'?', nucleotide sequence information not available for the primer matching sequence.

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

PRIMER SPECIFICATION

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec.	95	130	280	85	280	135	165	125	245	110	155	100
PCR product	150	200		170		330	390	355	425	165	255	155
		270		275			445			390	340	220
		300										
Length of int. pos. control ¹	800	800	800	1070	800	800	1070	800	800	800	1070	800
5'-primer(s) ²	89 5'-gAA 3'	47 5'-Agg 3'	105 5'-gCT 3'	98 5'-CTT 3'	105 5'-gCT 3'	47 5'-Agg 3'	1 st I 5'-CgA 3'	28 5'-TCA 3'	47 5'-Agg 3'	1 st I 5'-CgA 3'	1 st I 5'-CgA 3'	361 5'-AgT 3'
	89 5'-gAA 3'	89 5'-gAA 3'		355 5'-TCA 3'		112 5'-CCT 3'	176 5'-gCA 3'	2 nd I 5'-CCA 3'	648 5'-CAC 3'	176 5'-gCA 3'	419 5'-gTC 3'	419 5'-gTC 3'
	89 5'-gAA 3'	2 nd I 5'-CCA 3'		459 5'-gAT 3'		118 5'-CCA 3'	527 5'-TgA 3'			527 5'-TAC 3'		
	98 5'-CTT 3'	703 5'-CTA 3'										
3'-primer(s) ³	142 5'-TgA 3'	176 5'-ACT 3'	343 5'-C 3'	142 5'-TgA 3'	343 5'-T 3'	201 5'-CTT 3'	302 5'-ggT 3'	213 5'-Cgg 3'	302 5'-ggC 3'	302 5'-ggC 3'	118 5'-gCT 3'	474 5'-gCA 3'
	201 5'-CTT 3'	559 5'-CTC 3'	343 5'-g 3'	589 5'-CTT 3'		218 5'-gCT 3'	307 5'-CCA 3'	409 5'-ATg 3'	853 5'-CAT 3'	595 5'-CCC 3'	218 5'-gCT 3'	477 5'-gCA 3'
	201 5'-CTC 3'	861 5'-TCg 3'					3 rd I 5'-gCA 3'	420 5'-gCT 3'			302 5'-ggT 3'	538 5'-gCA 3'
											304 5'-CAA 3'	
											539 5'-TCT 3'	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Well No.	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec.	120	160	130	160	110	180	225	215	325	135	155	95
PCR product	250	210	255	260	325	210	250	420	380		235	475
				445		240	290					590
Length of int. pos. control ¹	800	800	800	1070	800	1070	800	800	1070	800	1070	800
5'-primer(s) ²	201 5'-CCA 3'	97 5'-TCg 3'	98 5'-CTA 3'	98 5'-CAC 3'	201 5'-CCA 3'	201 5'-CCA 3'	2 nd I 5'-CCA 3'	47 5'-Agg 3'	355 5'-TCA 3'	105 5'-gCT 3'	97 5'-TCg 3'	176 5'-gCA 3'
	218 5'-ggA 3'	361 5'-AgT 3'	368 5'-gTg 3'	98 5'-CTC 3'	409 5'-ggC 3'	2 nd I 5'-CCA 3'		117 5'-CCT 3'	412 5'-ATA 3'	459 5'-gAT 3'	97 5'-TTg 3'	486 5'-ACA 3'
	2 nd I 5'-CCA 3'			98 5'-CTT 3'		361 5'-AgT 3'		361 5'-AgT 3'				
				194 5'-CgT 3'								
				527 5'-TgA 3'								
3'-primer(s) ³	289 5'-AgC 3'	218 5'-gCT 3'	312 5'-AgT 3'	311 5'-ggT 3'	270 5'-TAg 3'	341 5'-CgT 3'	512 5'-CCA 3'	289 5'-AgT 3'	3 rd I 5'-CTC 3'	201 5'-CTC 3'	213 5'-Cgg 3'	361 5'-CCA 3'
	289 5'-AgC 3'	527 5'-CCg 3'	459 5'-AgA 3'	317 5'-CgT 3'	3 rd I 5'-CTC 3'	343 5'-T 3'	538 5'-gTC 3'	302 5'-ggT 3'		559 5'-CTC 3'	289 5'-AgC 3'	474 5'-gCA 3'
	291 5'-TCg 3'	538 5'-CCA 3'		3 rd I 5'-gCA 3'		343 5'-g 3'	580 5'-TCC 3'	527 5'-CCg 3'			289 5'-AgC 3'	538 5'-CCA 3'
	539 5'-TCT 3'	538 5'-gCA 3'				527 5'-CCg 3'		538 5'-CCg 3'				
						527 5'-CCg 3'						
						527 5'-CCg 3'						
						527 5'-CCg 3'						
						530 5'-CCA 3'						
Well No.	13	14	15	16	17	18	19	20	21	22	23	24

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

Well No.	25	26	27	28	29	30	31
Length of spec.	260	155	140	215	175	160	520
PCR product	360	225	245	555		370	
	425					540	
Length of int.	1070	1070	1070	1070	800	1070	1070
pos. control ¹							
5'-primer(s) ²	28	98	97	98	134	289	312
	5'-TCA 3'	5'-CTA 3'	5'-TCA 3'	5'-CTA 3'	5'-CCA 3'	5'-Agg 3'	5'-AAA 3'
	341	176	201	368		418	
	5'-ggA 3'	5'-gCA 3'	5'-CCA 3'	5'-gTT 3'		5'-Agg 3'	
		361		368		419	
		5'-AgT 3'		5'-gTT 3'		5'-gTC 3'	
				1018		419	
				5'-gTg 3'		5'-gTT 3'	
3'-primer(s) ³	112	289	302	361	270	369	539
	5'-CCA 3'	5'-AgC 3'	5'-ggC 3'	5'-CCA 3'	5'-TAg 3'	5'-CCg 3'	5'-TCC 3'
	118	289		538		539	
	5'-gCT 3'	5'-AgC 3'		5'-CCg 3'		5'-TCT 3'	
	218	538		1092			
	5'-gCT 3'	5'-CCg 3'		5'-TTA 3'			
	477						
	5'-gCA 3'						
Well No.	25	26	27	28	29	30	31

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

				Well																
				1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
				Prod. No.:	201786701	201786702	201786703	201786704	201786705	201786706	201786707	201786708	201786709	201786710	201786711	201786712	201550813	201786714	201888615	201890116
IHCW cell line ¹		C*																		
1	9001 SA	*07:02		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
2	9280 LK707	*07:01	*15:05	-	-	+	-	-	-	-	-	-	+	-	+	-	-	-	+	-
3	9011 E4181324	*12:02		-	-	-	-	-	-	-	-	-	-	-	-	+	+	-	-	-
4	9275 GU373	*03:04	*04:01	-	-	+	+	-	+	-	-	-	-	-	+	-	+	-	-	-
5	9009 KAS011	*06:02		-	-	-	-	-	-	-	-	+	-	-	+	-	-	+	-	-
6	9353 SM	*03:04	*07:02	-	-	+	+	-	-	-	-	-	+	-	-	-	-	-	-	-
7	9020 QBL	*05:01		-	-	-	-	-	-	-	+	-	-	-	+	-	-	-	+	-
8	9025 DEU	*04:01		-	-	-	-	-	-	+	-	-	-	-	+	-	+	-	-	-
9	9026 YAR	*12:03		-	-	-	-	-	-	-	-	-	-	-	+	-	+	+	-	-
10	9107 LKT3	*01:02		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT	*16:01		-	-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-
12	9052 DBB	*06:02		-	-	-	-	-	-	-	-	+	-	-	+	-	-	+	-	-
13	9004 JESTHOM	*01:02		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA	*01:02	*03:04	+	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075 DKB	*03:04		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007	*02:02		-	+	+	-	-	+	-	-	-	-	-	+	-	-	-	+	-
17	9282 CTM3953540	*03:03	*07:01	-	-	-	+	+	-	-	-	+	-	-	-	-	-	-	-	-
18	9257 32367	*01:02	*07:05	+	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
19	9038 BM16	*07:01		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
20	9059 SLE005	*03:04		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA	*03:03		-	-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	+	+	+	-	-
23	9124 IHL	*01:02	*15:02	+	-	+	-	-	-	-	+	-	-	+	-	-	-	-	+	-
24	9035 JBUSH	*12:03		-	-	-	-	-	-	-	-	-	-	-	+	+	+	-	-	-
25	9049 IBW9	*08:02		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
26	9285 WT49	*07:18		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
27	9191 CH1007	*07:04	*15:29	-	-	+	-	-	-	-	-	+	-	+	-	-	-	-	+	-
28	9320 BEL5GB	*05:01	*16:01	-	-	-	-	-	-	+	-	-	-	+	-	+	-	+	-	-
29	9050 MOU	*16:01		-	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
30	9021 RSH	*17:01		-	+	-	-	-	-	-	-	-	-	+	-	-	-	-	+	-
31	9019 DUCAF	*05:01		-	-	-	-	-	-	+	-	-	-	+	-	-	-	-	+	-
32	9297 HAG	*17:03		-	+	-	-	-	-	-	-	-	-	+	-	-	-	-	+	-
33	9098 MT14B	*03:04		-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF	*12:03		-	-	-	-	-	-	-	-	-	-	-	+	+	+	-	-	-
35	9302 SSTO	*05:01		-	-	-	-	-	-	+	-	-	-	+	-	-	-	-	+	-
36	9024 KT17	*03:03	*04:01	-	-	-	+	+	+	-	-	-	-	+	-	+	-	-	-	-
37	9065 HHKB	*07:02		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
38	9099 LZL	*03:03		-	-	-	+	+	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML	*02:02	*07:01	-	+	+	-	-	+	-	-	+	-	+	-	-	-	-	+	-
40	9134 WHONP199	*01:02	*06:02	+	-	-	-	-	-	-	+	-	-	+	-	-	-	+	-	-
41	9055 H0301	*08:02		-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-
42	9066 TAB089	*01:02		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526	*01:02	*08:01	+	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
44	9057 TEM	*12:03		-	-	-	-	-	-	-	-	-	-	-	-	+	+	+	-	-
45	9239 SHJO	*06:02	*17:01	-	+	-	-	-	-	-	+	-	-	+	-	-	-	+	+	-
46	9013 SCHU	*07:02		-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
47	9045 TUBO	*07:04	*15:02	-	-	+	-	-	-	-	+	-	+	-	-	-	-	-	+	-
48	9303 TER-ND	*04:01	*16:01	-	-	-	-	-	+	-	-	-	-	+	-	+	-	-	-	-

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

CELL LINE VALIDATION SHEET																			
HLA-C low resolution SSP primer set ²																			
				Well															
				17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
				Prod. No.:	201786717	201786718	201890119	201786720	201786721	201786722	201781723	201781724	201888625	201781726	201781727	201781728	201888629	201781730	201781731
	IHCW cell line ¹		C*																
1	9001	SA	*07:02		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
2	9280	LK707	*07:01	*15:05	+	-	-	-	+	-	+	-	-	-	-	+	-	-	-
3	9011	E4181324	*12:02		-	-	-	-	-	-	-	-	-	+	+	-	-	+	-
4	9275	GU373	*03:04	*04:01	-	+	-	-	+	+	-	-	+	-	-	-	-	-	+
5	9009	KAS011	*06:02		-	-	-	-	-	-	+	-	-	-	-	-	-	-	+
6	9353	SM	*03:04	*07:02	-	-	-	-	+	+	-	-	-	-	-	-	-	-	-
7	9020	QBL	*05:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	+
8	9025	DEU	*04:01		-	+	-	-	+	-	-	-	+	-	-	-	-	-	+
9	9026	YAR	*12:03		-	-	-	-	-	-	-	-	-	+	+	-	-	-	-
10	9107	LKT3	*01:02		-	-	-	+	-	-	-	-	-	+	-	-	-	-	-
11	9051	PITOUT	*16:01		-	+	-	-	-	-	-	-	-	-	+	-	-	+	-
12	9052	DBB	*06:02		-	-	-	-	-	-	+	-	-	-	-	-	-	-	+
13	9004	JESTHOM	*01:02		-	-	-	+	-	-	-	-	-	+	-	-	-	-	-
14	9071	OLGA	*01:02	*03:04	-	-	-	+	+	+	-	-	-	+	-	-	-	-	-
15	9075	DKB	*03:04		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-
16	9037	SWEIG007	*02:02		-	-	-	-	-	-	-	+	-	-	-	-	-	-	+
17	9282	CTM3953540	*03:03	*07:01	-	-	-	-	+	+	+	-	-	-	-	-	-	-	-
18	9257	32367	*01:02	*07:05	-	-	-	+	-	-	+	-	-	+	-	-	-	-	-
19	9038	BM16	*07:01		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
20	9059	SLE005	*03:04		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-
21	9064	AMALA	*03:03		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-
22	9056	KOSE	*12:03		-	-	-	-	-	-	-	-	-	+	+	-	-	-	-
23	9124	IHL	*01:02	*15:02	+	-	-	+	+	-	-	-	-	+	-	-	+	-	-
24	9035	JBUSH	*12:03		-	-	-	-	-	-	-	-	-	+	+	-	-	-	-
25	9049	IBW9	*08:02		-	-	-	-	+	-	-	-	-	-	+	-	-	-	-
26	9285	WT49	*07:18		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
27	9191	CH1007	*07:04	*15:29	+	-	-	-	+	-	+	-	-	-	-	-	+	-	-
28	9320	BEL5GB	*05:01	*16:01	-	+	-	-	+	-	-	-	-	-	+	-	-	+	+
29	9050	MOU	*16:01		-	+	-	-	-	-	-	-	-	-	+	-	-	+	-
30	9021	RSH	*17:01		-	-	+	-	+	-	-	-	-	+	-	-	-	-	-
31	9019	DUCAF	*05:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	+
32	9297	HAG	*17:03		-	-	+	-	+	-	-	-	-	+	-	-	-	-	-
33	9098	MT14B	*03:04		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-
34	9104	DHIF	*12:03		-	-	-	-	-	-	-	-	-	+	+	-	-	-	-
35	9302	SSTO	*05:01		-	-	-	-	+	-	-	-	-	-	-	-	-	-	+
36	9024	KT17	*03:03	*04:01	-	+	-	-	+	+	-	-	+	-	-	-	-	-	+
37	9065	HHKB	*07:02		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
38	9099	LZL	*03:03		-	-	-	-	+	+	-	-	-	-	-	-	-	-	-
39	9315	CML	*02:02	*07:01	-	-	-	-	-	-	+	+	-	-	-	-	-	-	+
40	9134	WHONP199	*01:02	*06:02	-	-	-	+	-	-	+	-	-	+	-	-	-	-	+
41	9055	H0301	*08:02		-	-	-	-	+	-	-	-	-	-	+	-	-	-	-
42	9066	TAB089	*01:02		-	-	-	+	-	-	-	-	-	+	-	-	-	-	-
43	9076	T7526	*01:02	*08:01	-	+	-	+	+	-	-	-	-	+	+	-	-	-	-
44	9057	TEM	*12:03		-	-	-	-	-	-	-	-	-	+	+	-	-	-	-
45	9239	SHJO	*06:02	*17:01	-	-	+	-	+	-	+	-	-	+	-	-	-	-	+
46	9013	SCHU	*07:02		-	-	-	-	-	-	+	-	-	-	-	-	-	-	-
47	9045	TUBO	*07:04	*15:02	+	-	-	-	+	-	+	-	-	-	-	-	+	-	-
48	9303	TER-ND	*04:01	*16:01	-	+	-	-	+	-	-	-	+	-	+	-	-	+	+



101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

Additional 5'-primers and 3'-primers in primer solutions 1, 2, 4, 7, 14 to 16, 20, 22, 24 to 26 and 28 were tested by separately adding one or two additional 3'-primers, respectively one or two additional 5'-primers. Additional 3'-primers in primer solution 8, 11, 18 and 19 were tested by separately adding one or two additional 5'-primers. Additional 5'-primers in primer solutions 6, 13, 27 and 30 were tested by separately adding one or two additional 3'-primers.

In primer solutions 1, 16, 23 and 28 one or two 5'-primers were not possible to test, and in primer solutions 3, 7, 11 to 16 and 18 one or two 3'-primers were not possible to test.

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

101.601-24/12 – including *Taq* pol., IFU-01
101.601-24u/12u – without *Taq* pol., IFU-02

Visit www.olerup.com for
“Instructions for Use” (IFU)

Lot No.: **2G8**

Lot-specific information

ADDRESSES:

Manufacturer:

Olerup SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-717 88 27

Fax: +46-8-717 88 18

E-mail: olerup-se@caredx.com

Web page: <http://www.olerup.com>

Distributed by:

Olerup GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

Tel: +43-1-710 15 00

Fax: +43-1-710 15 00 10

E-mail: olerup-at@caredx.com

Web page: <http://www.olerup.com>

Olerup Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-OLERUP1

Fax: 610-344-7989

E-mail: olerup-us@caredx.com

Web page: <http://www.olerup.com>

For information on *Olerup* distributors worldwide, contact **Olerup GmbH**.