

Olerup SSP[®] HLA-B*27 - bulk

Product number:	101.532-48 – including <i>Taq</i> polymerase 101.532-48u – without <i>Taq</i> polymerase
Lot number:	38G
Expiry date:	2011-October-01
Number of tests:	48

CHANGES COMPARED TO THE PREVIOUS HLA-B*27 – BULK LOT (54F):
HLA-B*27 bulk primer set is unchanged.

THE NUMBER OF PRIMER MIXES is unchanged.

ALLELE COVERAGE:

B*2701 to B*2759N, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2009-October-19, release 2.27.0.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:

The B*3702, B*4704 and B*4705 alleles are amplified by primer mix 2 of the HLA-B*27 bulk kit.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.