

## **Olerup SSP® HLA-B\*27 - bulk**

|                  |   |
|------------------|---|
| Product number:  | 101.532-48 – including <i>Taq</i> polymerase<br>101.532-48u – without <i>Taq</i> polymerase |
| Lot number:      | 83K   |
| Expiry date:     | 2013-July-01  |
| Number of tests: | 48  |

### **CHANGES COMPARED TO THE PREVIOUS HLA-B\*27 – BULK LOT (38G):**

HLA-B\*27 bulk primer set is unchanged.

THE NUMBER OF PRIMER MIXES is unchanged.

### **ALLELE COVERAGE:**

HLA-B\*27:01 to B\*27:73, i.e. all the currently recognized HLA-B\*27 alleles, will be amplified by the B\*27 primer mix; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2010-October-15, release 3.2.0.

### **INFLUENCE ON THE INTERPRETATION OF HLA-B\*27 TYPINGS BY NON-HLA-B\*27 ALLELES:**

The B\*44:97 is amplified by primer mix 1 and the B\*37:02, B\*47:04 and B\*47:05 alleles are amplified by primer mix 2 of the HLA-B\*27 kit.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.