

Olerup SSP® HLA-B*27 - bulk

Product number:	101.532-48 – including <i>Taq</i> polymerase 101.532-48u – without <i>Taq</i> polymerase
Lot number:	30S
Expiry date:	2015-April-01
Number of tests:	48

CHANGES COMPARED TO THE PREVIOUS HLA-B*27 – BULK LOT (20N):

HLA-B*27 bulk primer set is unchanged.

THE NUMBER OF PRIMER MIXES is unchanged.

ALLELE COVERAGE:

HLA-B*27:01 to B*27:100, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2013-January-11, release 3.11.0.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:

The B*44:97 is amplified by primer mix 1 and the B*37:02, 47:04-47:05 and 48:26 alleles are amplified by primer mix 2 of the HLA-B*27 kit.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.