

101.532-48 – including *Taq* polymerase101.532-48u – without *Taq* polymerase

Lot No.: 16V

**Olerup SSP® HLA-B\*27 - bulk**

**Product number:** 101.532-48 – including *Taq* polymerase  
101.532-48u – without *Taq* polymerase

**Lot number:** 16V

**Expiry date:** 2016-January-01

**Number of tests:** 48

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*27 – BULK LOT (30S):**

HLA-B\*27 bulk primer set is unchanged.

THE NUMBER OF PRIMER MIXES is unchanged.

**ALLELE COVERAGE:**HLA-B\*27:01 to B\*27:107, i.e. all the currently recognized HLA-B\*27 alleles, will be amplified by the B\*27 primer mix; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2013-October-11, release 3.14.0.**INFLUENCE ON THE INTERPRETATION OF HLA-B\*27 TYPINGS BY NON-HLA-B\*27 ALLELES:**

Due to the sharing of sequence motifs between HLA-B alleles a few non-HLA-B\*27 alleles will be amplified by both primer mixes.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.