

## **Olerup SSP® HLA-B\*27 - bulk**

<b>Product number:</b>	101.532-48 – including <i>Taq</i> polymerase 101.532-48u – without <i>Taq</i> polymerase
<b>Lot number:</b>	46Y
<b>Expiry date:</b>	2017-May-01
<b>Number of tests:</b>	48

### **CHANGES COMPARED TO THE PREVIOUS HLA-B\*27 – BULK LOT (15X):**

HLA-B\*27 bulk primer set is unchanged.

THE NUMBER OF PRIMER MIXES is unchanged.

### **ALLELE COVERAGE:**

HLA-B\*27:01 to B\*27:135, i.e. all the currently recognized HLA-B\*27 alleles, will be amplified by the B\*27 primer mix; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2015-January-19, release 3.19.0.

### **INFLUENCE ON THE INTERPRETATION OF HLA-B\*27 TYPINGS BY NON-HLA-B\*27 ALLELES:**

Due to the sharing of sequence motifs between HLA-B alleles a few non-HLA-B\*27 alleles will be amplified by both primer mixes.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.