

Olerup SSP® HLA-B*27 - bulk

| | |
|-------------------------|---|
| Product number: | 101.532-48 – including <i>Taq</i> polymerase 101.532-48u – without <i>Taq</i> polymerase |
| Lot number: | 4D4 |
| Expiry date: | 2017-December-01 |
| Number of tests: | 48 |

CHANGES COMPARED TO THE PREVIOUS HLA-B*27 – BULK LOT (46Y):

HLA-B*27 bulk primer set is unchanged.

THE NUMBER OF PRIMER MIXES is unchanged.

ALLELE COVERAGE:

HLA-B*27:01 to B*27:142, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2015-August-11, release 3.21.1.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:

Due to the sharing of sequence motifs between HLA-B alleles a few non-HLA-B*27 alleles will be amplified by both primer mixes.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.