

101.532-48 – including *Taq* polymerase

101.532-48u – without *Taq* polymerase

Lot No.: 2F8

Olerup SSP[®] HLA-B*27 - bulk

Product number: 101.532-48 – including *Taq* polymerase
101.532-48u – without *Taq* polymerase
Lot number: 2F8
Expiry date: 2019-05-01
Number of tests: 48

CHANGES COMPARED TO THE PREVIOUS HLA-B*27 – BULK LOT (4D4):

Well	5'-primer	3'-primer	rationale
1	Modified	Added	Strength of control band has been optimized, 5'-primer modified for improved HLA-specific amplification, 3'-primer added for the B*27:153 allele.

THE NUMBER OF PRIMER MIXES is unchanged.

ALLELE COVERAGE:

HLA-B*27:01 to B*27:156, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2016-October-14, release 3.26.0.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:

Due to the sharing of sequence motifs between HLA-B alleles a few non-HLA-B*27 alleles will be amplified by both primer mixes.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.