

## Olerup SSP<sup>®</sup> HLA-B\*27 - bulk

Product number:	101.532-48 – including <i>Taq</i> polymerase 101.532-48u – without <i>Taq</i> polymerase
Lot number:	2G7
Expiry date:	2021-10-01
Number of tests:	48

### CHANGES COMPARED TO THE PREVIOUS HLA-B\*27 – BULK LOT (2F8):

HLA-B\*27 bulk primer set is unchanged.

THE NUMBER OF PRIMER MIXES is unchanged.

### ALLELE COVERAGE:

HLA-B\*27:01 to B\*27:164, i.e. all the currently recognized HLA-B\*27 alleles, will be amplified by the B\*27 primer mix; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2017-October-27, release 3.30.0.

### INFLUENCE ON THE INTERPRETATION OF HLA-B\*27 TYPINGS BY NON-HLA-B\*27 ALLELES:

Due to the sharing of sequence motifs between HLA-B alleles a few non-HLA-B\*27 alleles will be amplified by both primer mixes.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.

Change in revision R01 compared to R00:

1. The expiration date has been altered due to extension of shelf-life.