

Olerup SSP[®] HLA-B*27 - bulk

Product number:	101.532-48 – including <i>Taq</i> polymerase 101.532-48u – without <i>Taq</i> polymerase
Lot number:	4L0
Expiry date:	2022-08-01
Number of tests:	48

CHANGES COMPARED TO THE PREVIOUS HLA-B*27 – BULK LOT (5H6):

The HLA-B*27 bulk primer set is unchanged compared to the previous lot.

THE NUMBER OF PRIMER MIXES is unchanged.

ALLELE COVERAGE:

HLA-B*27:01 to B*27:226, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2020-April-20, release 3.40.0.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:

Due to the sharing of sequence motifs between HLA-B alleles a few non-HLA-B*27 alleles will be amplified by both primer mixes.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.