

## Olerup SSP<sup>®</sup> HLA-B\*27 - bulk

Product number:	101.532-48 – including <i>Taq</i> polymerase 101.532-48u – without <i>Taq</i> polymerase
Lot number:	2N4
Expiry date:	2023-11-01
Number of tests:	48

### Changes compared to the previous HLA-B\*27 – BULK Lot (4L0):

The HLA-B\*27 bulk primer set is unchanged compared to the previous lot.

THE NUMBER OF PRIMER MIXES is unchanged.

### ALLELE COVERAGE:

HLA-B\*27:01 to B\*27:248, i.e. all the currently recognized HLA-B\*27 alleles, will be amplified by the B\*27 primer mix; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2021-July-12, release 3.45.0.

### INFLUENCE ON THE INTERPRETATION OF HLA-B\*27 TYPINGS BY NON-HLA-B\*27 ALLELES:

Due to the sharing of sequence motifs between HLA-B alleles a few non-HLA-B\*27 alleles will be amplified by both primer mixes.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.



0197

For *In Vitro* Diagnostic Use  
MA125 v01 SSP Release Note Template  
Date: January 2022, Rev. No: 00