

101.532-48 – including *Taq* polymerase

101.532-48u – without *Taq* polymerase

Lot No.: 5S9

Olerup SSP® HLA-B*27 - bulk

Product number:	101.532-48 – including <i>Taq</i> polymerase 101.532-48u – without <i>Taq</i> polymerase
Lot number:	5S9
Expiry date:	2025-10-01
Number of tests:	48

Changes compared to the previous HLA-B*27 – BULK Lot (7R4):

The HLA-B*27 bulk primer set is unchanged compared to the previous lot.

THE NUMBER OF PRIMER MIXES is unchanged.

ALLELE COVERAGE:

HLA-B*27:01 to B*27:267, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2023-July-12, release 3.53.0.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:

Due to the sharing of sequence motifs between HLA-B alleles a few non-HLA-B*27 alleles will be amplified by both primer mixes.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.



0197