

101.532-48 – including *Taq* polymerase101.532-48u – without *Taq* polymerase

Lot No.: 7R4

Olerup SSP[®] HLA-B*27 - bulk

Product number: 101.532-48 – including *Taq* polymerase
101.532-48u – without *Taq* polymerase

Lot number: 7R4

Expiry date: 2025-02-01

Number of tests: 48

Changes compared to the previous HLA-B*27 – BULK Lot (2N4):

The HLA-B*27 bulk primer set is unchanged compared to the previous lot.

THE NUMBER OF PRIMER MIXES is unchanged.

ALLELE COVERAGE:

HLA-B*27:01 to B*27:264, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2022-October-12, release 3.50.0.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:

Due to the sharing of sequence motifs between HLA-B alleles a few non-HLA-B*27 alleles will be amplified by both primer mixes.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.

**0197**For *In Vitro* Diagnostic Use

MA125 v01 SSP Release Note Template

Date: March 2023, Rev. No: 00