

HLA-B*27 – bulk
101.531-48 – licensed for PCR
101.531-48u – not licensed for PCR
Lot No.: **X19**

1

Release note *Olerup SSP*TM HLA-B*27 – bulk – Lot Nr. **X19**

Product number: 101.532-48 – licensed for PCR
101.532-48u – not licensed for PCR
Lot number: X19
Expiry date: 2009-June-01
Number of tests: 48

CHANGES COMPARED TO THE PREVIOUS HLA-B*27 – BULK LOT (V37):
HLA-B*27 bulk primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*2701 to B*2737, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2007-April-12, release 2.17.0.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:

The B*3702, B*4704 and B*4705 alleles are amplified by primer mix 2 of the HLA-B*27 bulk kit.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.