

## **Olerup SSP<sup>®</sup> HLA-B\*27 – unit dose**

Product number:	101.531-48 – including <i>Taq</i> polymerase 101.531-48u – without <i>Taq</i> polymerase
Lot number:	77G
Expiry date:	2011-December-01
Number of tests:	48

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*27 – UNIT DOSE LOT (53F):**  
HLA-B\*27 unit dose primer set is unchanged.

**THE NUMBER OF PRIMER MIXES** is unchanged.

**ALLELE COVERAGE:**

B\*2701 to B\*2759N, i.e. all the currently recognized HLA-B\*27 alleles, will be amplified by the B\*27 primer mix; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2009-October-19, release 2.27.0.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*27 TYPINGS BY NON-HLA-B\*27 ALLELES:**

The B\*3702, B\*4704 and B\*4705 alleles are amplified by primer mix 2 of the HLA-B\*27 unit dose kit.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.