

Olerup SSP[®] HLA-B*27 – unit dose

Product number: 101.531-48 – including *Taq* polymerase
101.531-48u – without *Taq* polymerase
Lot number: 84K
Expiry date: 2013-May-01
Number of tests: 48

CHANGES COMPARED TO THE PREVIOUS HLA-B*27 – UNIT DOSE LOT (77G):
HLA-B*27 unit dose primer set is unchanged.

THE NUMBER OF PRIMER MIXES is unchanged.

ALLELE COVERAGE:

B*27:01 to B*27:73, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2010-October-15, release 3.2.0.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:

The B*44:97 is amplified by primer mix 1 and the B*37:02, B*47:04 and B*47:05 alleles are amplified by primer mix 2 of the HLA-B*27 kit.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.