

Olerup SSP® HLA-B*27 – unit dose

Product number:	101.531-48 – including <i>Taq</i> polymerase 101.531-48u – without <i>Taq</i> polymerase
Lot number:	16X
Expiry date:	2017-March-01
Number of tests:	48

CHANGES COMPARED TO THE PREVIOUS HLA-B*27 – UNIT DOSE LOT (65S):

HLA-B*27 unit dose primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*27:01 to B*27:125, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2014-April-14, release 3.16.0.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:

Due to the sharing of sequence motifs between HLA-B alleles, a few non-HLA-B*27 alleles will be amplified by primer mixes 1 and 2.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.