

101.531-48 – including *Taq* polymerase101.531-48u – without *Taq* polymerase

Lot No.: 57Y

Olerup SSP® HLA-B*27 – unit dose

Product number: 101.531-48 – including *Taq* polymerase
101.531-48u – without *Taq* polymerase

Lot number: 57Y

Expiry date: 2017-December-01

Number of tests: 48

CHANGES COMPARED TO THE PREVIOUS HLA-B*27 – UNIT DOSE LOT (16X):

HLA-B*27 unit dose primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:B*27:01 to B*27:138, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2015-April-17, release 3.20.0.**INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:**

Due to the sharing of sequence motifs between HLA-B alleles, a few non-HLA-B*27 alleles will be amplified by primer mixes 1 and 2.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.