

## **Olerup SSP® HLA-B\*27 – unit dose**

<b>Product number:</b>	101.531-48 – including <i>Taq</i> polymerase 101.531-48u – without <i>Taq</i> polymerase
<b>Lot number:</b>	7D3
<b>Expiry date:</b>	2018-09-01
<b>Number of tests:</b>	48

### **CHANGES COMPARED TO THE PREVIOUS HLA-B\*27 – UNIT DOSE LOT (57Y):**

HLA-B\*27 unit dose primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

### **ALLELE COVERAGE:**

B\*27:01 to B\*27:145, i.e. all the currently recognized HLA-B\*27 alleles, will be amplified by the B\*27 primer mix; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2015-October-10, release 3.22.0.

### **INFLUENCE ON THE INTERPRETATION OF HLA-B\*27 TYPINGS BY NON-HLA-B\*27 ALLELES:**

Due to the sharing of sequence motifs between HLA-B alleles, a few non-HLA-B\*27 alleles will be amplified by primer mixes 1 and 2.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.