

## **Olerup SSP® HLA-B\*27 – unit dose**

**Product number:** 101.531-48 – including *Taq* polymerase  
101.531-48u – without *Taq* polymerase  
**Lot number:** 5F1  
**Expiry date:** 2020-02-01  
**Number of tests:** 48

### **CHANGES COMPARED TO THE PREVIOUS HLA-B\*27 – UNIT DOSE LOT (7D3):**

<b>Well</b>	<b>5'-primer</b>	<b>3'-primer</b>	<b>rationale</b>
1	Modified	Added	Strength of control band has been optimized, 5'-primer modified for improved HLA-specific amplification, 3'-primer added for the B*27:153 allele.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

B\*27:01 to B\*27:161, i.e. all the currently recognized HLA-B\*27 alleles, will be amplified by the B\*27 primer mix; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2017-January-20, release 3.27.0.

### **INFLUENCE ON THE INTERPRETATION OF HLA-B\*27 TYPINGS BY NON-HLA-B\*27 ALLELES:**

Due to the sharing of sequence motifs between HLA-B alleles, a few non-HLA-B\*27 alleles will be amplified by primer mixes 1 and 2.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.