

**101.531-48 – including *Taq* polymerase****101.531-48u – without *Taq* polymerase****Lot No.: 3G2****Olerup SSP® HLA-B\*27 – unit dose**

|                         |   |
|-------------------------|---|
| <b>Product number:</b>  | <b>101.531-48 – including <i>Taq</i> polymerase</b><br><b>101.531-48u – without <i>Taq</i> polymerase</b> |
| <b>Lot number:</b>      | <b>3G2</b>  |
| <b>Expiry date:</b>     | <b>2020-10-01</b>   |
| <b>Number of tests:</b> | <b>48</b>   |

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*27 – UNIT DOSE LOT (5F1):**

HLA-B\*27 unit dose primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

**ALLELE COVERAGE:**B\*27:01 to B\*27:164, i.e. all the currently recognized HLA-B\*27 alleles, will be amplified by the B\*27 primer mix; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2017-October-27, release 3.30.0.**INFLUENCE ON THE INTERPRETATION OF HLA-B\*27 TYPINGS BY NON-HLA-B\*27 ALLELES:**

Due to the sharing of sequence motifs between HLA-B alleles, a few non-HLA-B\*27 alleles will be amplified by primer mixes 1 and 2.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.