

**101.531-48 – including *Taq* polymerase****101.531-48u – without *Taq* polymerase****Lot No.: 2K6**

## **Olerup SSP<sup>®</sup> HLA-B\*27 – unit dose**

**Product number:** 101.531-48 – including *Taq* polymerase  
101.531-48u – without *Taq* polymerase

**Lot number:** 2K6

**Expiry date:** 2023-09-01

**Number of tests:** 48

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*27 – UNIT DOSE LOT (3H1):**  
HLA-B\*27 unit dose primer set is unchanged.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

B\*27:01 to B\*27:197, i.e. all the currently recognized HLA-B\*27 alleles, will be amplified by the B\*27 primer mix; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2019-April-17, release 3.36.0.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*27 TYPINGS BY NON-HLA-B\*27 ALLELES:**

Due to the sharing of sequence motifs between HLA-B alleles, a few non-HLA-B\*27 alleles will be amplified by primer mixes 1 and 2.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.