

101.531-48 – including *Taq* polymerase101.531-48u – without *Taq* polymerase

Lot No.: 2L2

Olerup SSP[®] HLA-B*27 – unit dose

Product number: 101.531-48 – including *Taq* polymerase
101.531-48u – without *Taq* polymerase
Lot number: 2L2
Expiry date: 2024-05-01
Number of tests: 48

CHANGES COMPARED TO THE PREVIOUS HLA-B*27 – UNIT DOSE LOT (2K6):

The HLA-B*27 unit dose primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*27:01 to B*27:223N, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2020-January-20, release 3.39.0.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:

Due to the sharing of sequence motifs between HLA-B alleles, a few non-HLA-B*27 alleles will be amplified by primer mixes 1 and 2.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.