

101.531-48 – including *Taq* polymerase101.531-48u – without *Taq* polymeraseLot No.: **1R9****Olerup SSP<sup>®</sup> HLA-B\*27 – unit dose**

**Product number:** 101.531-48 – including *Taq* polymerase  
101.531-48u – without *Taq* polymerase

**Lot number:** 1R9

**Expiry date:** 2026-06-01

**Number of tests:** 48

**Changes compared to the previous HLA-B\*27 – UNIT DOSE Lot (5N7):**

The HLA-B\*27 unit dose primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

B\*27:01 to B\*27:248, i.e. all the currently recognized HLA-B\*27 alleles, will be amplified by the B\*27 primer mix; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2022-April-19, release 3.48.0.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*27 TYPINGS BY NON-HLA-B\*27 ALLELES:**

Due to the sharing of sequence motifs between HLA-B alleles, a few non-HLA-B\*27 alleles will be amplified by primer mixes 1 and 2.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.

**0197**

For *In Vitro* Diagnostic Use  
MA125 v01 SSP Release Note Template  
Date: June 2022, Rev. No: 00