

HLA-B*27 – unit dose
101.531-48 – licensed for PCR
101.531-48u – not licensed for PCR
Lot No.: Y02

1

Release Note - Olerup SSP™ HLA-B*27 – Lot No. Y02

Product number: 101.531-48 – licensed for PCR
101.531-48u – not licensed for PCR
Lot number: Y02
Expiry date: 2009-June-01
Number of tests: 48
Number of tubes per test: 2

CHANGES COMPARED TO THE PREVIOUS HLA-B*27 – UNIT DOSE LOT (X18):

Tube	5'-primer	3'-primer	rationale
3	Modified	-	Increased yield of specific PCR product.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*2701 to B*2738, i.e. all the currently recognized HLA-B*27 alleles, will be amplified by the B*27 primer mix; www.ebi.ac.uk/imgt/hla, 2007-July-09, release 2.18.0.

INFLUENCE ON THE INTERPRETATION OF HLA-B*27 TYPINGS BY NON-HLA-B*27 ALLELES:

The B*3702, B*4704 and B*4705 alleles are amplified by primer mix 2 of the HLA-B*27 kit.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.