

Olerup SSP[®] DPA1

Product number: 101.331-24/06 – including *Taq* pol.
101.331-24u/06u – without *Taq* pol.
Lot number: 11E
Expiry date: 2010-February-01
Number of tests: 24 test – Product No. 101.331-24
6 tests – Product No. 101.331-06
Number of Wells per test: 16

CHANGES COMPARED TO THE PREVIOUS DPA1 LOT (X42):

Tube	5'-primer	3'-primer	rationale
6	-	Added	Primer added for the DPA1*0204 allele.
16	Modified	-	Increased specificity of specific primer pair.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DPA1*0103 to 0109, DPA1*0201 to 0204, DPA1*0301 to 0303 and DPA1*0401, i.e. all the currently recognized DPA1 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-January-11, release 2.20.0.

RESOLUTION IN DPA1 HOMO- AND HETEROZYGOTES:

Excellent. All frequently occurring genotypes give rise to unique amplification patterns.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.