

## Olerup SSP<sup>®</sup> DPA1

**Product number:** 101.331-24/06 – including *Taq* pol.  
101.331-24u/06u – without *Taq* pol.  
**Lot number:** 20F  
**Expiry date:** 2010-October-01  
**Number of tests:** 24 test – Product No. 101.331-24  
6 tests – Product No. 101.331-06  
**Number of Wells per test:** 16

### CHANGES COMPARED TO THE PREVIOUS DPA1 LOT (11E):

Well	5'-primer	3'-primer	rationale
6	Added	-	Primer added for the DPA1*0110 allele.

**THE NUMBER OF WELLS** is unchanged.

### ALLELE COVERAGE:

DPA1\*0103 to 0110, DPA1\*0201 to 0204, DPA1\*0301 to 0303 and DPA1\*0401, i.e. all the currently recognized DPA1 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2008-October-10, release 2.23.0.

### RESOLUTION IN DPA1 HOMO- AND HETEROZYGOTES:

Excellent. All frequently occurring genotypes give rise to unique amplification patterns.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.